

**LongAmp®  
Hot Start Taq  
DNA Polymerase**



1-800-632-7799  
info@neb.com  
www.neb.com



**M0534S**



**500 units      2,500 U/ml      Lot: 0031212**  
**RECOMBINANT    Store at -20°C    Exp: 12/14**

**Description:** LongAmp Hot Start *Taq* DNA Polymerase is a unique blend of aptamer-based Hot Start *Taq* and Deep Vent<sub>r</sub>™ DNA Polymerases. The aptamer-based inhibitor binds reversibly to the enzyme, inhibiting polymerase activity at temperatures below 45°C, but releases the enzyme during normal PCR cycling conditions. This permits assembly of PCR reactions at room temperature. An advantage of the aptamer-based hot start mechanism is that it does not require a separate high temperature incubation step to activate the enzyme. The 3'→5' exonuclease activity of Deep Vent<sub>r</sub> DNA Polymerase increases the fidelity and robust amplification of Hot Start *Taq* DNA Polymerase (1). LongAmp Hot Start *Taq* DNA Polymerase offers two-fold higher fidelity than Hot Start *Taq* DNA Polymerase alone. A wide range of PCR products can be generated; up to 30 kb from lambda or human genomic DNA.

**Source:** An *E. coli* strain that carries the *Taq* DNA Polymerase gene from *Thermus aquaticus* YT-1 and an *E. coli* strain that carries the Deep Vent<sub>r</sub> DNA Polymerase gene from *Pyrococcus* species GB-D.

**Applications:**

- High-specificity Long Range PCR

Supplied in: 100 mM KCl, 10 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 1 mM dithiothreitol, 0.5% Tween® 20, 0.5% IGEPAL® CA-630 and 50% glycerol.

**Reagents Supplied with Enzyme:**

5X LongAmp *Taq* Reaction Buffer

**Reaction Conditions:** 1X LongAmp *Taq* Reaction Buffer, DNA template, primers, 300 μM dNTPs (not included) and 1–5 units of LongAmp Hot Start *Taq* DNA Polymerase in a total reaction volume of 50 μl.

**1X LongAmp Taq Reaction Buffer:**

60 mM Tris-SO<sub>4</sub> (pH 9.0 @ 25°C)  
20 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>  
2 mM MgSO<sub>4</sub>  
3% glycerol  
0.06% IGEPAL CA-630  
0.05% Tween 20

**Unit Definition:** One unit is defined as the amount of enzyme that will incorporate 10 nmol of dNTP into acid insoluble material in 30 minutes at 75°C.

**Unit Assay Conditions:** 1X ThermoPol® Reaction Buffer, 200 μM dNTPs including [<sup>3</sup>H]-dTTP and 200 μg/ml activated Calf Thymus DNA.

**Heat Inactivation:** No

**Quality Control Assays**


**Long Amplicon PCR (Lambda DNA):** 28 cycles of PCR amplification in a 25 μl reaction containing 1 ng of Lambda DNA with 2.5 units of LongAmp Hot Start *Taq* DNA Polymerase in the presence of 300 μM dNTPs and 0.4 μM primers in LongAmp *Taq* Reaction Buffer results in the expected 20 and 30 kb products.

**Long Amplicon PCR (Genomic DNA):** 28 cycles of PCR amplification in a 25 μl reaction containing 500 ng of human genomic DNA with 2.5 units of LongAmp Hot Start *Taq* DNA Polymerase in the presence of 300 μM dNTPs and 0.4 μM primers in LongAmp *Taq* Reaction Buffer results in the expected 20 and 30 kb products.

**High Sensitivity PCR:** 35 cycles of PCR amplification of 2 ng of human genomic DNA with 5 units of LongAmp Hot Start *Taq* DNA Polymerase in the presence of 200 μM dNTPs and 0.2 μM primers in 1X LongAmp *Taq* Reaction Buffer results in the hot start-specific expected 306 bp product after pre-incubation at room temperature for 1 hour.

**Inhibition Assay:** > 95% inhibition is observed after a 16 hour incubation at 25°C in a 50 μl primer extension assay containing 10 units of LongAmp Hot Start *Taq* DNA Polymerase in 1X ThermoPol Reaction Buffer with 200 μM dNTPs including [<sup>3</sup>H]-dTTP and 15 nM primed M13 DNA template.

**Note:** Product specifications for individual components in the LongAmp Hot Start *Taq* DNA Polymerase are available separately.



Amplification of a 6 kb amplicon from varying amounts of human genomic DNA using LongAmp Hot Start *Taq* DNA Polymerase. Starting template amounts are indicated below the gel. Marker M is the NEB 2-Log DNA Ladder (NEB #N3200).

**PCR**

The Polymerase Chain Reaction (PCR) is a powerful and sensitive technique of DNA amplification (2). *Taq* DNA Polymerase is an enzyme widely used in PCR (3). LongAmp Hot Start *Taq* DNA Polymerase allows for greater PCR sensitivity and permits room temperature PCR set-up. The following guidelines are provided to ensure successful PCR using New England Biolabs' LongAmp Hot Start *Taq* DNA Polymerase. These guidelines cover routine PCR reactions. Amplification of templates with high GC content, high secondary structure or low template concentrations may require further optimization.

**Reaction Setup:**

Due to the hot-start nature of the enzyme, reactions can be assembled on the bench at room temperature and transferred to a thermocycler. No separate activation step is required.

COMPONENT	25 μl REACTION	50 μl REACTION	FINAL CONCENTRATION
5X LongAmp <i>Taq</i> Reaction Buffer	5 μl	10 μl	1X
10 mM dNTPs	0.75 μl	1.5 μl	300 μM
10 μM Forward Primer	1 μl	2 μl	0.4 μM (0.05–1 μM)
10 μM Reverse Primer	1 μl	2 μl	0.4 μM (0.05–1 μM)
Template DNA	variable	variable	<1,000 ng
LongAmp Hot Start <i>Taq</i> DNA Polymerase	1 μl	2 μl	5 units/50 μl PCR
Nuclease-free water	to 25 μl	to 50 μl	

Notes: Gently mix the reaction. Avoid pipetting samples containing target DNA when amplicons above 20 kb are desired. Collect all liquid to the bottom of the tube by a quick spin if necessary. Overlay the sample with mineral oil if using a PCR machine without a heated lid.

Transfer PCR tubes to a PCR machine and begin thermocycling:

**Thermocycling Conditions for a Routine PCR:**

STEP	TEMP	TIME
Initial Denaturation	94°C	30 seconds
30 Cycles	94°C	10–30 seconds
	45–65°C	15–60 seconds
	65°C	50 seconds/kb
Final Extension	65°C	10 minutes
Hold	4–10°C	

**General Guidelines:**

1. **Template:**  
The quality of the DNA template is essential for long-range PCR amplification. Recommended amounts of DNA template for a 50 μl reaction are as follows:

DNA	Up to 15 kb	Above 15 kb
Genomic	1 ng–500 ng	10 ng–1 μg
Plasmid or Viral	1 pg–1 ng	10 pg–10 ng

Successful amplification above 20 kb largely depends on the quality of DNA templates and the primer sequences.

2. **Primers:**  
Oligonucleotide primers are generally 20–40 nucleotides in length and ideally have a GC content of 40–60%. Computer programs such as Primer3 (<http://frodo.wi.mit.edu/primer3>) can be used to design or analyze primers. For amplicons larger than 20 kb, it is desirable to have primers with GC content above 50%, matched T<sub>m</sub>s above 60°C and at least 24 nucleotides in length. The final concentration of each primer in a PCR reaction may be 0.05–1 μM, typically 0.1–0.5 μM.
3. **Mg<sup>++</sup> and additives:**  
Mg<sup>++</sup> concentration of 1.0–2.0 mM is optimal for most PCR products generated with LongAmp Hot Start *Taq* DNA Polymerase. The final Mg<sup>++</sup> concentration in 1X LongAmp *Taq* Reaction Buffer is 2 mM. This supports satisfactory amplification of most amplicons. However, Mg<sup>++</sup> can be further optimized in 0.5 or 1.0 mM increments using MgSO<sub>4</sub> (not provided).

Amplification of some difficult targets, like GC-rich sequences, may be improved with additives, such as DMSO (4) or formamide (5).

(see other side)

4. Deoxynucleotides:  
The recommended final concentration of dNTPs for long-range PCR is 300  $\mu$ M of each deoxy-nucleotide.

5. LongAmp Hot Start *Taq* DNA Polymerase concentration:  
We generally recommend using LongAmp Hot Start *Taq* DNA Polymerase at a concentration of 200 units/ml (5 units/50  $\mu$ l reaction). However, the optimal concentration of LongAmp Hot Start *Taq* DNA Polymerase may vary in specialized applications.

6. Denaturation:  
No separate activation step is required to release the hot start inhibitor from the enzyme.  
An initial denaturation of 30 seconds at 94°C is sufficient for most amplicons from pure DNA templates. For difficult templates such as GC-rich sequences, a longer denaturation of 2–4 minutes at 94°C is recommended prior to PCR cycling to fully denature the template. With colony PCR, an initial 5 minute denaturation at 94°C is recommended.

During thermocycling, a 10-30 second denaturation at 94 °C is recommended.

7. Annealing:  
The annealing step is typically 15–60 seconds. Annealing temperature is based on the  $T_m$  of the primer pair and is typically 45–65°C. Annealing temperatures can be optimized by doing a temperature gradient PCR starting 5°C below the calculated  $T_m$ . We recommend using NEB's  $T_m$  Calculator, available at [www.neb.com/TmCalculator](http://www.neb.com/TmCalculator) to determine appropriate annealing temperatures for PCR.

When primers with annealing temperatures above 60°C are used, a 2-step PCR protocol is possible (see #10).

8. Extension:  
The recommended extension temperature is 65°C. Extension times are generally 50 seconds per kb. A final extension of 10 minutes at 65°C is recommended.

9. Cycle Number:  
Generally, 25–35 cycles yields sufficient product. Up to 45 cycles may be required to detect low copy number targets.

10. 2-step PCR:  
When primers with annealing temperatures above 60°C are used, a 2-step thermocycling protocol is possible.

#### Thermocycling Conditions for a Routine 2-Step PCR:

STEP	TEMP	TIME
Initial Denaturation	94°C	30 seconds
30 Cycles	94°C	10–30 seconds
	60–65°C	50 seconds/kb
Final Extension	60–65°C	10 minutes
Hold	4–10°C	

11. PCR product:  
The majority of the PCR products generated using LongAmp Hot Start *Taq* DNA Polymerase contain dA overhangs at the 3'-end; therefore the PCR products can be ligated to the dT/dU-overhang vectors.

#### FAQs:

1. *What is the recommended enzyme amount when using LongAmp Hot Start?*  
In general, we recommend 5 units of LongAmp Hot Start *Taq* DNA Polymerase in a 50  $\mu$ l reaction. For amplicons < 8 kb, we recommend 1-2.5 units per 50  $\mu$ l reaction for higher fidelity.

2. *What is the fidelity of the LongAmp Hot Start Taq DNA Polymerase compared to Taq DNA Polymerase?*  
The LongAmp Hot Start *Taq* DNA Polymerase offers two-fold higher fidelity than *Taq*.

3. *Can the extension step be carried out at 72°C when using LongAmp Hot Start?*  
Yes, LongAmp Hot Start *Taq* DNA Polymerase can be used at 72°C. However, extension at 65–68°C is a better choice for most amplicons.

4. *What is the extension rate when using LongAmp Hot Start?*  
We recommend 50 seconds per kb for maximum yields. Extension rates such as 30 seconds per kb can be used for targets up to 4 kb using a 3-step PCR protocol. Shorter extension rates, such as 15 seconds per kb, can be used for targets up to 2 kb using a 3-step PCR protocol on a fast PCR machine.

5. *What type of DNA ends result from a primer extension reaction or a PCR using LongAmp Hot Start Taq DNA Polymerase?*  
The majority of the PCR products generated using LongAmp Hot Start *Taq* DNA Polymerase contain dA overhangs at the 3'-end; therefore, the PCR products can be ligated to dT/dU-overhang vectors.

6. *Why is the product a smear when visualized on an agarose gel?*  
When PCR conditions are not optimal, a smear or high level of background is often observed. Try one or more of the following suggestions:

- use lower amount of enzymes
- use 65°C for extension
- raise annealing temperature
- try 2-step cycling protocols

7. *Can LongAmp Hot Start Taq DNA Polymerase be used to amplify GC-rich amplicons?*  
Yes. The addition of DMSO up to 10% helps amplify GC-rich amplicons.

#### References:

1. Barnes, W.M. (1994) *Proc. Natl. Acad. Sci. USA*, 91, 2216–2220.
2. Saiki R.K. et al. (1985) *Science*, 230, 1350–1354.
3. Powell, L.M. et al. (1987) *Cell*, 50, 831–840.
4. Sun, Y., Hegamyer, G. and Colburn, N. (1993) *Biotechniques*, 15, 372–374.
5. Sarkar, G., Kapelner, S. and Sommer, S.S. (1990) *Nucleic Acids Res.*, 18, 7465.

#### Companion Products Sold Separately:

- LongAmp *Taq* (Mg-free) Reaction Buffer Pack  
#B0322S 6.0 ml
- LongAmp *Taq* Reaction Buffer Pack  
#B0323S 6.0 ml
- Crimson LongAmp *Taq* Reaction Buffer Pack  
#B0326S 6.0 ml
- Magnesium Sulfate (MgSO<sub>4</sub>) Solution  
#B1003S 6.0 ml
- Diluent F  
#B8006S 4.0 ml
- LongAmp Hot Start *Taq* 2X Master Mix  
#M0533S 100 Reactions  
#M0533L 500 Reactions

LongAmp *Taq* PCR Kit  
#E5200S 100 Reactions

LongAmp *Taq* 2X Master Mix  
#M0287S 100 Reactions  
#M0287L 500 Reactions

Crimson LongAmp *Taq* DNA Polymerase  
#M0326S 250 units  
#M0326L 1,250 units

Deoxynucleotide Solution Set  
#N0446S 25  $\mu$ mol of each

Deoxynucleotide Solution Mix  
#N0447S 8  $\mu$ mol of each  
#N0447L 40  $\mu$ mol of each



Notice to Purchaser: Nucleic acid-based aptamers for use with thermophilic DNA polymerases are licensed exclusively by New England Biolabs, Inc. from SomaLogic, Inc. (See Patent Nos. 5,475,096; 5,670,637; 5,696,249; 5,874,557; and 5,693,502). New England Biolabs, Inc. gives the Buyer/User a non-exclusive license to use the aptamer-based LongAmp Hot Start *Taq* DNA Polymerase for RESEARCH PURPOSES ONLY. Commercial use of the aptamer-based LongAmp Hot Start *Taq* DNA Polymerase requires a license from New England Biolabs, Inc. Please contact [busdev@neb.com](mailto:busdev@neb.com) for more information.

Limited Label License: This product is covered by certain patents outside the United States (the "Non-US Patents") corresponding to the expired U.S. Patent No. 5,436,149 (claims 6-16). The purchase of this product conveys to the purchaser only the limited, non-transferable right under the Non-US Patents to use the purchased quantity of the product for the purchaser's own research by the purchaser only. No rights are granted to the purchaser with respect to the Non-US Patents to sell, modify for resale or otherwise transfer this product, either alone or as a component of another product, to any third party. Takara Bio reserves all other rights to the Non-US Patents, and this product may not be used in any manner other than as provided herein. For information on obtaining a license with respect to the Non-US Patents to use this product for purposes other than research, please contact Takara Bio Inc., Seta 3-4-1, Otsu, Shiga 520-2193, Japan (Fax +81-77-543-9254).

LONGAMP®, NEW ENGLAND BIOLABS® and THERMOPOL® are registered trademarks of New England Biolabs, Inc. DEEP VENT™ is a trademark of New England Biolabs, Inc.

IGEPAL® is a registered trademark of Rhodia Operations. TWEEN® is a registered trademark of Uniqema Americas LLC.