



NEW ENGLAND
BioLabs[®] Inc.

be INSPIRED
drive DISCOVERY
stay GENUINE

2021•22



In loving memory

of our Founder, Donald G. Comb
(1927 – 2020)



Watch a video tribute
to New England Biolabs
Founder, Don Comb

Since our establishment in 1974, New England Biolabs has been different. From our founding principles — placing the advancement of science and stewardship of the environment as our highest priorities — to our unique corporate culture, NEB's philosophy can be distilled down to three core values: passion, humility and being genuine. This was a priority for our Founder, Don Comb, and these values continue to guide us today.

Don's provocative curiosity, boundless enthusiasm and commitment to the environment will continue to serve as inspiration to all of us. His vision for how a company could be run was unconventional but highly effective — people and passion over process and profit, the supply of high quality reagents at a fair price, the co-mingling of basic and applied research, and the prioritization of protecting the planet and advancing science. His generosity and the faith he had in each of us resulted in a family-like culture where everyone has a voice and all employees feel valued. We are all honored and humbled to have been a part of his life and look forward to continuing his legacy.

When thinking about everything that has happened in 2020 —from the pandemic, deep political polarization and graphic images of racial inequity — I can say that this has been the most trying, challenging and emotionally draining year of my thirty-seven at NEB. I'm guessing that many of you feel the same way. At the same time, it has also been one of the most exhilarating, selfless and hopeful —with science taking center stage and the promise of SARS-CoV-2 vaccines, still believing in the inherent goodness in people and the glimmers of hope for peace and unity across the globe. We have always known that science is more than a vocation. It embodies an ethos that inspires acts of compassion, brilliance, originality and a sense of community. Let's continue to care for, support and appreciate one another. I know we will get through all of this and be stronger than ever.

Thank you for your ongoing trust and support. If there is anything you believe we could be doing differently, please share your thoughts with us. We wish you continued success in your research.

Jim Ellard
CEO, New England Biolabs, Inc.

Want to experience the NEB catalog in augmented reality?



Download

the NEB Augmented Reality (AR) app for iPhone® or iPad® at the Apple® App Store or for Android™ on Google Play™.



Explore

the NEB Catalog & Technical Reference and keep an eye out for the augmented reality icon.



Discover

videos, tutorials and immersive experiences that can be accessed by scanning the icon using the AR app.



be
INSPIRED

drive
DISCOVERY

stay
GENUINE

At New England Biolabs, we are motivated by a set of core values that are still as true today as they were when the company was founded 45 years ago. These principles continue to guide us both as a company and as individuals.



Take a tour of
New England Biolabs



Advancement of Science

We believe that basic research and the cultivation of scientific knowledge is critical for us to stay connected with our customers and to drive scientific breakthroughs. At NEB, over 30 labs participate in research projects, which are aided by post-doctoral fellows and students in Masters and Ph.D. programs. NEB researchers have authored or co-authored over 1,300 publications to date, most of which are in peer-reviewed journals.

Environmental Stewardship

We continuously strive to promote ecologically sound practices and environmental sustainability in order to protect our natural resources, both locally and globally. It is our goal to continuously improve our business processes to minimize our impact on the environment. Additionally, NEB supports other organizations that are advancing environmental research and stewardship. We are excited to announce that NEB has become a Certified B Corporation™, a recognition awarded to organizations with the highest standards for social, and environmental performance, transparency and accountability. To learn more, visit www.neb.com/environmentalphilosophy.



Learn more about
our commitment
to the environment

Social Responsibility

We see opportunities where science can be used to improve lives, and we continue to be guided by our responsibility to each other and our community to work towards a kinder and more just world. This philosophy lies behind NEB's longstanding commitment to its parasitology research program, which contributes to the understanding and treatment of poorly-funded and understudied tropical diseases. NEB also supports several organizations devoted to humanitarian efforts. Further, we recognize that we must work together to build a more equitable society and improve diversity, equity and inclusion in our workplace. To learn more, visit www.neb.com/corporateresponsibility.

Delivering the Highest Quality Product

It is our goal to deliver best-in-class product quality and technical support. With a reliance on recombinant technologies, the majority of our products are designed and manufactured in our ISO 13485:2016 and ISO 9001:2015 certified facility in Ipswich, MA, USA. We are constantly improving the stringency and range of our quality controls to ensure that our products will perform to your expectations, every time. In addition, our manufacturing facility in Rowley, MA produces GMP-grade* materials for customers requiring an enhanced level of quality documentation and support.



*See page 7 for more details.

A Unique Approach to Wastewater Treatment

Our state of the art Solar Aquatics System™ utilizes and accelerates processes found in streams and wetlands, and employs proven wastewater treatment practices. This innovative method can treat up to 27,500 gallons of effluent per day through a series of filtration beds and aerated translucent tanks that host plant communities and aerobic microorganisms. The entire wastewater stream from our facility is treated, and the clean effluent is used for groundwater recharge and irrigation.

Solar Aquatics System™ is a trademark of Ecological Engineering Associates.

The NEB Facility

Our research and production facility, located in Ipswich, MA, USA, is LEED® certified, which is awarded based on environmentally-focused standards that include:

- Water efficiency
- Energy conservation
- Atmospheric protection
- Sustainable building materials and resources
- Indoor environmental quality
- Innovation and building design

LEED® is a registered trademark of the U.S. Green Building Council Corporation.

Diversity, Equity & Inclusion

We recognize that there are areas that we need to grow and are taking steps to raise our level of consciousness to injustices that have been overlooked and underestimated. Our Diversity, Equity & Inclusion Team works to address these issues, and is divided into four subgroups focusing on:

- Social Justice Philanthropy
- STEM Education & Mentorship
- Social Justice Outreach
- Diversity, Equity & Inclusion at NEB



Partnering with NEB

NEB has over 45 years of experience in the discovery, development and manufacture of molecular biology reagents. These are essential components in a vast array of genomic and proteomic technologies that continue to transform our understanding of the world we live in, and ultimately the diagnosis and treatment of disease. With experience in fields as diverse as next generation sequencing, RNA biology, qPCR and protein engineering, NEB is ready to work with you to develop custom reagent solutions, and to help bring your technologies to market. Further, our global distribution network can help to ensure that your products have worldwide reach.



Customized Solutions

From development to commercialization, NEB provides the technical expertise, consistent scalable manufacturing, quality systems and a global distribution network to enable a successful long-term partnership. Our dedicated team is ready to work with you to develop novel, high performance enzymes tailored to your application, optimize these enzymes in your workflow, enable small to large scale production, generate quality controls and customize packaging. With our ISO 13485 and ISO 9001 certified manufacturing processes, as well as the ability to manufacture GMP-grade* products, you can be confident in our robust process, documentation and risk mitigation for the product you need. For more information, contact custom@neb.com.

Global Development

The ability of NEB to successfully operate as both a research institute and a commercial enterprise in service of our customers is amplified by the extent of our global reach. The Global Development (GD) team at NEB operates worldwide to generate sustainable growth through an exceptional network of commercial operations that includes eight wholly owned subsidiaries located in Australia, Canada, China, France, Germany, Japan, Singapore and the United Kingdom, and over 60 distribution partners covering our involvement in more than 90 countries. By leveraging the talents and assets of NEB, including scientific and commercial resources, we ensure that our customers are serviced by a stable, ethical and engaged global network. More information can be found on the inside back cover or contact globaldev@neb.com.

NEBnow Freezer Program Network

With NEBnow® on-site freezers, enjoy convenient and affordable access to NEB's high quality reagents, anytime. Our NEBnow Freezer Program Team works closely with your institution to customize inventory best suited to your research program. Save time and avoid shipping fees with consolidated shipments. For more information, contact freezers@neb.com.

Enzymes for Innovation

The NEB catalog highlights a wide variety of enzyme functionalities found in nature or engineered for specific purposes. However, in molecular biology, new tools can often lead to new discoveries. Taking advantage of the enzymology expertise at NEB, we now offer a growing selection of novel enzymes with interesting and unique activities for manipulating DNA, RNA, proteins and glycans, even if specific applications for them have yet to be discovered. If you are looking for an enzyme functionality that it is not currently available, visit www.enzymesforinnovation.com or contact enzymesforinnovation@neb.com.



What are Enzymes
for Innovation?

**See page 7 for more details.*

“When it comes time to choose a source of enzymes, my first choice is NEB. In addition to having a wide selection of enzymes, I’ve been impressed with their rigorous test procedures and overall quality of their products. The staff is responsive, knowledgeable, customer-focused and a pleasure to work with.”

– Senior Fellow, Analytics & Knowledge Transfer,
Molecular Biology Reagents Provider

“The NEBnow freezer program is an amazing addition; having 24/7 access to the NGS reagents has been extremely helpful. NEB allows you to customize the freezer to hold the items you use most, saving us valuable time by cutting out the ordering process!”

– Assistant Director, Genomic Sequencing and Analysis Facility,
University of Texas, Austin





Practicing Ethical Science

New England Biolabs (NEB) is committed to practicing ethical science – we believe it is our job as researchers to ask the important questions that, when answered, help preserve our quality of life and the world that we live in. However, this research should always be done in safe and ethical manner. Learn more at www.neb.com/neb-ethics.

Supporting COVID-19 Research

NEB is supplying and supporting customers who are working diligently to develop better diagnostic tools and vaccines for the SARS-CoV-2 virus, including diagnostics for at home or point-of-care settings. Visit www.neb.com/COVID19 to learn more, and find products and publications related to the applications below.



RNA Extraction

Extraction of RNA is often the first step in viral detection assays and has become an increasingly important laboratory technique in the midst of this global pandemic. Though NEB does not offer a kit specifically designed for viral RNA extraction, both the Monarch® Total RNA Miniprep Kit and the Monarch RNA Cleanup Kits have been used successfully to extract total RNA, including viral RNA, from clinically-relevant samples. See pages 137-138 for more information.

Virus Detection

Virus detection often utilizes nucleic acid amplification to identify the presence of specific sequences in SARS-CoV-2 viral RNA. These methods include RT-qPCR for real time quantitation and loop-mediated isothermal amplification (LAMP) for robust amplification performed at a single reaction temperature. NEB has several research use only (RUO) products that can be used for detection of SARS-CoV-2 viral RNA, including our Luna® qPCR and RT-qPCR kits, which have been mentioned in several emergency use authorization (EUA) protocols. See pages 72-78 for more information.

Sequencing & Epidemiology

During the COVID-19 pandemic, sequencing has enabled scientists, often in unprecedented collaborations, to make rapid advances in epidemiology and surveillance, basic and clinical research, and diagnostics. A fast-growing number of methods are being developed to address these applications, including targeted SARS-CoV-2 sequencing protocols. Our NEBNext ARTIC kits—developed for long and short read sequencing—were based on the original work of the ARTIC Network (1) and are available for both Oxford Nanopore Technologies® and Illumina® Platforms. See pages 161-162 for more information.

Viral Biology

Proteolytic cleavage between the S1 and S2 unit of the spike protein of SARS-CoV-2 is an important part of the process to allow viral entry into the host cell after binding to the surface receptor. A subsequent cleavage on the S2 unit further triggers membrane fusion. Furin is one of the key players in a family of proteases responsible for these steps. The susceptibility of the spike protein to furin is essential to understand not only the viral biology, but also COVID-19's transmissibility.

The spike protein is also extensively glycosylated. These glycans have profound implications in modulating cell entry and membrane fusion, as well as host immune response. These glycosylations can be analyzed and profiled by PNGase F, Trypsin and Endoproteinase LysC. See pages 256-257 and 265-269 for more information.

Vaccine Development

mRNA-based vaccines are an emerging alternative to conventional vaccine approaches. Unlike protein-based vaccines, mRNA vaccines involve introduction of an mRNA sequence encoding disease-specific antigens. Once delivered into the cells, they get translated by the cellular machinery, resulting in the synthesis of protein antigens. These antigens are recognized by the immune system, and immune responses are mounted. Several mRNA vaccine candidates are either already being administered or are in development for SARS-CoV-2. mRNA vaccines can be rapidly and scalably produced *in vitro* using enzymes to transcribe and modify mRNA. See pages 199-205 for more information.

(1) Josh Quick 2020. nCoV-2019 sequencing protocol v2 (Gunlt). protocols.io <https://dx.doi.org/10.17504/protocols.io.bdp7i5m>.

* "GMP-grade" is a branding term NEB uses to describe reagents manufactured at NEB's Rowley facility. The Rowley facility was designed to manufacture reagents under more rigorous infrastructure and process controls to achieve more stringent product specifications and customer requirements. Reagents manufactured at NEB's Rowley facility are manufactured in compliance with ISO 9001 and ISO 13485 quality management system standards. However, at this time, NEB does not manufacture or sell products known as Active Pharmaceutical Ingredients (APIs), nor does NEB manufacture its products in compliance with all of the Current Good Manufacturing Practice regulations.

GMP Capabilities

With over 45 years in enzymology expertise, NEB is a world leader in the discovery and production of reagents for the life science industry. This expertise effectively positions us to supply reagents for the synthesis of high-quality RNA – from template generation and transcription, to capping, tailing and cleanup after synthesis. These products are designed and manufactured based on decades of molecular biology experience, so that you can be confident they will work for your application.

To better serve the needs of customers in regulated markets, NEB has opened a state-of-the-art, 43,000 sq. ft. production facility in Rowley, MA for the manufacture of GMP-grade* materials. This facility includes Quality Control and Production functions ranging from a shipping/receiving area and dedicated warehouse, to separate inoculation preparation, fermentation, purification and filling suites. Visit www.neb.com/GMP to learn more.




Supporting Molecular Diagnostics

The COVID-19 pandemic has highlighted many of the challenges faced by clinical labs, requiring new and innovative solutions to address them. Technology development is happening faster than ever before, exposing weaknesses in diagnostic supply chains, increasing the need for innovation and thinking differently about how diagnostics should be developed, manufactured and deployed.

Many scientists know NEB as a trusted reagent provider to the life science community. What many do not know is that we also offer a portfolio of products that serve as critical components for a wide array of diagnostic products and services. Extensive molecular biology and enzymology experience provide NEB with the unique ability to help customers solve the challenges inherent in technology development and ultimately in scale-up and commercialization. Visit www.neb.com/MDx to learn more.

Supporting Non-Profits and Foundations

New England Biolabs has played a role in the establishment of several organizations that are advancing social responsibility and environmental stewardship worldwide.



The New England Biolabs Foundation is a private independent foundation whose mission is to foster community-based conservation of landscapes and seascapes, and the bio-cultural diversity found in these places. The foundation supports projects in selected countries of Central America, Andean South America, West Africa, and in coastal communities on the North Shore of Massachusetts. Learn more at NEBF.org.



Creative Action Institute catalyzes community-driven solutions that advance gender equality and build a sustainable planet. Through our experiential trainings, convenings and coaching, we develop leadership, build networks and support grassroots advocacy. Visit CreativeActionInstitute.org to learn more.



The Ocean Genome Legacy Center of New England Biolabs is a non-profit research center dedicated to the conservation of marine genome diversity and maintaining a repository of genomic DNA from marine organisms around the world. Learn more at northeastern.edu/oglc.

Ordering Information

Placing an Order

Orders can be placed:
9:00 AM to 8:00 PM EST (Monday – Friday)
Phone: 1-800-632-5227 (1-800-NEB-LABS)
Fax: 1-800-632-7440
Online: www.neb.com

To inquire about an existing order, email customerservice@neb.com or chat with Live Agent on the US Customer and Order Support page at www.neb.com.

NEB's main office and its subsidiaries in Canada, Germany and the UK all support e-mail ordering. For US customers only, email neborders@neb.com. For international orders, please contact the appropriate office.

Electronic Ordering (EDI, cXML)

Please inquire to our EDI administrator: edi@neb.com

Credit Card Purchases

NEB accepts Mastercard, Visa and American Express at the point of sale.

OEM & Customized Products

To inquire about OEM and customized products:
Phone: 1-800-NEB-LABS, ext. 7275
Email: custom@neb.com
Online: www.neb.com/neb-custom

For more information, see page 4 of this Catalog.

Shipping Rates

NEB offers free shipping to all U.S. available locations for all orders totaling \$350 (net) and over.

U.S. orders under \$350 (net) are currently charged \$31 per order for shipping and handling. Charges are prepaid and added to your invoice.

Please inquire for international shipping charge quotation.

Shipping terms and policies are subject to change at NEB's discretion.

Delivery

Orders received Monday through Thursday by 7:30 PM EST are shipped the same day for next business day delivery (within carrier's overnight delivery area)*. After 7:30 PM but before 8:00 PM EST, please place orders with Customer Service at 1-800-NEB-LABS, who can confirm if next business day delivery is still possible.

Orders placed by 7:30 PM EST on Friday will be shipped (with extra ice packs, as necessary, to NEB's courier, who stores the products at the appropriate storage temperature, prior to Sunday shipment, for delivery the following Monday (or next business day) nationwide. This service is available for the same cost as our standard shipping service.

PLEASE NOTE: Monday delivery service is not available for most products that have storage temperatures other than -20°C or require shipment on dry ice.

Please contact Customer Service at 1-800-NEB-LABS or customerservice@neb.com with any inquiries regarding NEB holiday deliveries or weather-related delays.

**"Local" orders received Monday through Friday before 10:00 AM EST can be delivered the same day to most of the Boston and Cambridge area. Please contact Customer Service, 1-800-NEB-LABS or customerservice@neb.com to arrange for this special service at no extra charge.

Terms

Net 30 days in U.S. dollars.

Discounts

In the U.S., NEB passes along discounts to its customers based on account history. Other discount opportunities include bulk (OEM) and new lab discounts. Please note: discounts are not cumulative.

Research Use Only

The products in this Catalog are intended for research purposes only. The products are not intended to be used for therapeutic or diagnostic purposes in humans or animals.

Terms of Sale

The purchase, acceptance, and payment of and for NEB's products is pursuant to NEB's Terms of Sale at www.neb.com/support/terms-of-sale. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

Legal Information

The products in this Catalog are covered by one or more patents, trademarks and/or copyrights owned or controlled by NEB.


While NEB develops and validates its products for various applications, the use of any product may require the buyer to obtain additional third party intellectual property rights for certain applications.

For more information about commercial rights, please contact NEB's Global Business Development team at gbd@neb.com.

Notice to Buyer

The information in this Catalog is provided in good faith based on NEB's knowledge as of the publication date. The Catalog is only a guide and may not be up-to-date. For current information, including a .pdf of current price list, please see NEB's website at www.neb.com. NEB's terms, conditions, policies and product information are subject to change at NEB's discretion.

Copyright ©2021
New England Biolabs, Inc. All rights reserved.



The employee photo has been a tradition in the catalog for over 30 years. Given the current situation, we went with a more creative approach to bring the NEB family together.

Table of Contents

Corporate Philosophy	2	GMP Capabilities	7
Partnering with NEB	4	Supporting Molecular Diagnostics	7
Supporting COVID-19 Research	6	Supporting Non-Profits & Foundations	8
		Ordering Information	8

Restriction Endonucleases

AatII	BceAI	BspCNI	DraIII-HF	HpyCH4III	NsiI	ScrFI
AbaSI	BcgI	BspDI	DrdI	HpyCH4IV	NsiI-HF	SexAI
AccI	BciVI	BspEI	EaeI	HpyCH4V	NspI	StaNi
Acc65I	BclI	BspHI	EagI-HF	KasI	PacI	StcI
Acil	BclII-HF	BspMI	EarI	KpnI	PaeR7I	SfiI
AcII	BcoDI	BspQI	Ecil	KpnI-HF	PaqCI	SfoI
AcuI	Bfal	BsrI	Eco53kI	LpnPI	PciI	SgrAI
AfeI	BfuAI	BsrBI	EcoNI	MboI	PfiFI	SmaI
AfIII	BglI	BsrDI	EcoO109I	MboII	PfiMI	SmlI
AfIII	BglII	BsrFI-v2	EcoP15I	MfeI	PleI	SnaBI
AgeI	BlpI	BsrGI	EcoRI	MfeI-HF	PluTI	SpeI
AgeI-HF	BmgBI	BsrGI-HF	EcoRI-HF	MluI	PmeI	SpeI-HF
AhdI	BmrI	BssHII	EcoRV	MluI-HF	PmlI	SphI
AleI-v2	BmtI	BssSI-v2	EcoRV-HF	MluCI	PpuMI	SphI-HF
AluI	BmtI-HF	BstAPI	Esp3I	MlyI	PshAI	SrfI
AlwI	BpmI	BstBI	FatI	MmeI	Psil-v2	SspI
AlwNI	Bpu10I	BstEI	FauI	MniI	PspGI	SspI-HF
Apal	BpuEI	BstEII-HF	Fnu4HI	MscI	PspOMI	StuI
ApalI	Bsal-HFv2	BstNI	FokI	MseI	PspXI	Styl-HF
ApeKI	BsaAI	BstUI	FseI	MsiI	PstI	StyDI
ApoI	BsaBI	BstXI	FspI	MspI	PstI-HF	Swal
ApoI-HF	BsaHI	BstYI	FspEI	MspA1I	PvuI	TaqI-v2
Ascl	BsaJI	BstZ17I-HF	HaeII	MspJI	PvuI-HF	TfiI
Asel	BsaWI	Bsu36I	HaeIII	MwoI	PvuII	TseI
AsiSI	BsaXI	BtgI	HgaI	NaeI	PvuII-HF	Tsp45I
AvaI	BseRI	BtgZI	HhaI	NarI	RsaI	TspMI
Avall	BseYI	BtsI-v2	HincII	NciI	RsrI	TspRI
AvrII	BsgI	BtsIMutI	HindIII	NcoI	SacI	Tth111I
BaeI	BsiEI	BtsCI	HindIII-HF	NcoI-HF	SacI-HF	XbaI
BaeGI	BsiHKAII	Cac8I	HinfI	NdeI	SacII	XcmI
BamHI	BsiWI	Clal	HinPI	NgoMIV	Sall	XhoI
BamHI-HF	BsiWI-HF	CspCI	HpaI	NheI-HF	Sall-HF	XmaI
BanI	BslI	CviAI	HpaII	NlaI	SapI	XmnI
BanII	BsmI	CviKI-1	HphI	NlaI	Sau3AI	ZrI
BbsI	BsmAI	CviQ1	Hpy99I	NmeAIII	Sau96I	
BbsI-HF	BsmBI-v2	DdeI	Hpy166II	NotI	SbfI	
BbvI	BsmFI	DpnI	Hpy188I	NotI-HF	SbfI-HF	
BbvCI	BsoBI	DpnII	Hpy188III	NruI	Scal-HF	
BccI	Bsp1286I	DraI	HpyAV	NruI-HF		

High-Fidelity Restriction Enzymes	21, 316
--	---------

Methylation Sensitive Restriction Enzymes for Epigenetics Studies	
AbaSI, DpnI, DpnII, FspEI, HpaII, LpnPI, McrBC, MspI, MspJI	280–281

Nicking Endonucleases	
Nb.BbvCI, Nb.BsmI, Nb.BsrDI, Nb.BssSI, Nb.BtsI, Nt.AlaI, Nt.BbvCI, Nt.BsmAI, Nt.BspQI, Nt.BstNBI, Nt.CviPII	55–57

Homing Endonucleases	
I-CeuI, I-SceI, PI-PspI, PI-SceI	57–58

NEBuffers, Diluents, Gel Loading Dyes, BSA & Recombinant Albumin, Molecular Biology Grade, NEB Tube Opener	58–59
---	-------

Polymerases & Amplification Technologies

Amplification-based Molecular Diagnostics	64	Specialty PCR	
PCR Polymerase Selection Chart	65	LongAmp Taq DNA Polymerase & Master Mix	70
High Fidelity PCR		LongAmp Hot Start Taq DNA Polymerase & 2X Master Mix	70
Q5 High-Fidelity DNA Polymerase & Master Mix	66	LongAmp Taq PCR Kit, Hemo KlenTaq	70
Q5 Hot Start High-Fidelity DNA Polymerase & Master Mix	66	EpiMark Hot Start Taq DNA Polymerase	71
Q5U Hot Start High-Fidelity DNA Polymerase	66	Other PCR Polymerases	
Q5 High-Fidelity PCR Kit	66	Vent & Vent (exo-) DNA Polymerases	71
Q5 Site-Directed Mutagenesis Kit (with or without competent cells)	98	Deep Vent & Deep Vent (exo-) DNA Polymerases	71
Phusion High-Fidelity DNA Polymerase & Master Mixes	67	qPCR & RT-qPCR	
Phusion Hot Start Flex DNA Polymerase & Master Mix	67	Luna qPCR & RT-qPCR	72
Phusion High-Fidelity PCR Kit	67	Find the Right Luna Product for Your Application	72
Routine PCR		Luna Universal qPCR Master Mix	73
One Taq DNA Polymerase & Master Mixes	68	Luna Universal Probe qPCR Master Mix	73
One Taq RT-PCR & One-step RT-PCR Kits	68	LunaScript RT SuperMix Kit	73
One Taq Hot Start DNA Polymerase & Master Mixes	68	LunaScript RT Master Mix (Primer-free)	73
Taq DNA Polymerase w/ThermoPol or Std. Buffer	69	Luna Probe One-Step RT-qPCR Kit (No ROX)	74
Taq DNA Polymerase w/Standard Taq (Mg-free) Buffer	69	Luna Probe One-Step RT-qPCR 4X Mix with UDG	74
Taq Master Mixes (2X, 5X, Quick-Load)	69	Luna Universal One-Step RT-qPCR Kit	74
Taq PCR Kit	69	Luna Universal Probe One-Step RT-qPCR Kit	74
Multiplex PCR 5X Master Mix	69	Luna Cell Ready One-Step RT-qPCR Kit	75
Hot Start Taq DNA Polymerase & 2X Master Mix	69	Luna Cell Ready Probe One-Step RT-qPCR Kit	75
		Luna Cell Ready Lysis Module	75

SARS-CoV-2 Detection	
SARS-CoV-2 Rapid Colorimetric LAMP Assay Kit	76
SARS-CoV-2 Positive Control (N gene)	76
WarmStart Fluorescent LAMP/RT-LAMP Kit w/UDG	77
WarmStart Multi-Purpose LAMP/RT-LAMP 2X Master Mix (with UDG)	77
LAMP Fluorescent Dye	77
Luna SARS-CoV-2 RT-qPCR Multiplex Assay Kit	77

Isothermal Amplification and Strand Displacement	
WarmStart Colorimetric LAMP/RT-LAMP 2X Master Mix (with or without UDG)	78
WarmStart Fluorescent LAMP/RT-LAMP Kit	78
Bst DNA Polymerase, (Large Frag. and Full-Length)	79
Bst 2.0 & Bst 2.0 WarmStart DNA Polymerases	79
Bst 3.0 DNA Polymerase	79
IsoAmp II Universal tHDA Kit	80
phi29 DNA Polymerase	80

RT-PCR	
One Taq One-Step RT-PCR Kit, One Taq RT-PCR Kit	81

Polymerases for DNA Manipulation (labeling, blunting, etc.)	
PreCR Repair Mix	82
Sulfolobus DNA Polymerase IV	82
Therminator DNA Polymerase	82
DNA Polymerase I (<i>E. coli</i>)	83
DNA Polymerase I, Large (Klenow) Fragment	83
Klenow Fragment (3'→5' exo-)	83
T4 DNA Polymerase	83
T7 DNA Polymerase (unmodified)	84
Bsu DNA Polymerase, Large Fragment	84
Terminal Transferase	84

Nucleotide Solutions & Buffers	
Reaction Buffers	85
Acyclonucleotide Set	85
Deoxynucleotide (dNTP) Solution Mix/Set	85
Ribonucleotide Solution Mix/Set	85
7-deaza-dGTP	85
Adenosine-5' Triphosphate (ATP)	85
5-methyl-dCTP	85
dATP, dUTP, dGTP Solutions	85

cDNA Synthesis	
cDNA Synthesis Selection Chart	86
LunaScript RT SuperMix Kit	73
LunaScript RT Master Mix Kit (Primer-free)	73
ProtoScript II First Strand cDNA Synthesis Kit	209
ProtoScript First Strand cDNA Synthesis Kit	209
Template Switching RT Enzyme Mix	208
ProtoScript II Reverse Transcriptase	206
M-MuLV, AMV Reverse Transcriptases	207
WarmStart RTx Reverse Transcriptase	207

PCR Cleanup	
Monarch PCR & DNA Cleanup Kit (5 µg)	87
Exo-CIP Rapid PCR Kit	87

DNA Modifying Enzymes & Cloning Technologies

Cloning & Synthetic Biology	
Cloning Workflow Comparison	92
NEBuilder HiFi DNA Assembly Master Mix & Cloning Kit	94
NEBuilder HiFi DNA Assembly Bundle for Large Fragments	94
Gibson Assembly Master Mix/Cloning Kit	95
Synthetic Biology/DNA Assembly Selection Chart	95
NEB Golden Gate Assembly Kits (BsaI-HFv2, BsmBI-v2)	96
BioBrick Assembly Kit	97
NEB PCR Cloning Kit (with and without competent cells)	97
USER & Thermolabile USER II Enzymes	97
Q5 Site Directed Mutagenesis Kit (with and without competent cells)	98
Quick Blunting Kit	98

DNA Ligases	
DNA Ligase Selection Chart	99
Substrate-based Ligase Selection Chart	100
T4 DNA Ligase	101
Instant Sticky-end Ligase Master Mix	101

Blunt T/A Ligase Master Mix	101	Monarch gDNA Buffers (Tissue Lysis, Cell Lysis, Blood Lysis, Binding, Wash, Elution)	135	Ultra II RNA	
Quick Ligation Kit	101	Monarch RNase A	135	NEBNext Ultra II Directional RNA Library Prep Kits +/- Sample Purification Beads	153
ElectroLigase	101	Proteinase K, Molecular Biology	135	NEBNext Ultra II RNA Library Prep Kits +/- Sample Purification Beads	153
Immobilized T4, Hi-T4, Salt-T4 DNA Ligases	101	HMW DNA Extraction		Reagents for RNA Depletion	
T3, T7, <i>E. coli</i> DNA Ligases	102	Monarch HMW DNA Extraction Kit for Tissue	136	NEBNext rRNA Depletion Kits v2 (Human/Mouse/Rat) +/- Sample Purification Beads	155
<i>Taq</i> , HiFi <i>Taq</i> DNA Ligases	103	Monarch HMW DNA Extraction Kit for Cells & Blood	136	NEBNext rRNA Depletion Kits (Bacteria) +/- Sample Purification Beads	155
9°N DNA Ligase, SplintR Ligase	214	Companion Products		NEBNext Globin RNA & rRNA Depletion Kits (H/M/R) +/- Sample Purification Beads	156
DNA Kinases		Monarch Pestle Set, Monarch 2 ml Tubes	136	NEBNext Globin RNA & rRNA Depletion Kit (Human/Mouse/Rat) with RNA Sample Purification Beads	156
T4 Polynucleotide Kinase	104	Monarch DNA Capture Beads	136	NEBNext Poly(A) mRNA Magnetic Isolation Module	156
T4 Polynucleotide Kinase (3' phosphatase minus)	104	Monarch Bead Retainers	136	Customized RNA Depletion	
5-hydroxymethyluridine Kinase	105	Monarch gDNA Nuclei Prep & Lysis Buffer Pack	136	NEBNext RNA Depletion Core Reagent Set (with & without RNA Sample Purification Beads)	157
DNA Phosphatases & Sulfurylases		Monarch RBC Lysis Buffer	136	Single Cell RNA	
Phosphatase Selection Chart	105	Monarch gDNA Elution Buffer II	136	NEBNext Single/Cell Low Input RNA Library Prep Kit for Illumina	158
Quick CIP	105	Monarch HMW gDNA Tissue Lysis Buffer	136	NEBNext Single/Cell Low Input cDNA Synthesis & Amplification Module	158
Shrimp Alkaline Phosphatase (rSAP)	106	Monarch Protein Separation Solution	136	NEBNext Single Cell Lysis Module	158
Antarctic Phosphatase	106	Monarch Precipitation Enhancer	136	Small RNA	
Pyrophosphatase, Inorganic (<i>E. coli</i> , yeast)	203	RNA Purification		NEBNext Multiplex Small RNA Library Prep Sets for Illumina (Sets 1 & 2)	159
Thermostable Inorganic Pyrophosphatase	203	Monarch Total RNA Miniprep Kit	137	NEBNext Multiplex Small RNA Library Prep Kit for Illumina (Index Primers 1-48)	159
NudC Pyrophosphatase	203	Monarch RNA Cleanup Kit (10 µg, 50 µg, 500 µg)	138	NEBNext Small RNA Library Prep Set for Illumina (Multiplex Compatible)	159
Apyrase, <i>Tte</i> UvrD Helicase	107	Companion Products		Modules & Enzymes	
Nucleases		Monarch RNA Purification Columns	137	NEBNext Ultra II FS DNA Module	150
Properties of Exonucleases and Non-specific Endonucleases	108	Monarch gDNA Removal Columns	137	Ultra II RNA First Strand Synthesis Module	154
Common Apps for Exonucleases & Endonucleases	109	Monarch DNA/RNA Protection Reagent	137	Ultra II Directional RNA Second Strand Synthesis Module	154
DNase I (RNase-Free)	217	Monarch RNA Buffers (Lysis, Priming, Wash, Cleanup Binding, Cleanup Wash)	137, 138	Ultra II Non-directional RNA Second Strand Synthesis Module	154
Lambda Exonuclease	110	Monarch Total RNA Miniprep Enzyme Pack	137	NEBNext Ultra II End Repair/dA-Tailing Module	150, 154
Exonuclease I, Thermolabile Exonuclease I	110	Monarch RNA Cleanup Columns (10 µg, 50 µg, 500 µg)	138	NEBNext Ultra II Ligation Module	150, 154
Exonuclease III, VII, Exonuclease V (RecBCD)	111	Monarch Collection Tubes II	137, 138	NEBNext Ultra II Q5 Master Mix	150
Exonuclease VIII, truncated, Exonuclease T	112	Nuclease-free Water	137, 138	Adaptors & Primers	
Thermostable FEN1	112	DNA Cleanup		NEBNext Multiplex Oligos (Unique Dual Index UMI Adaptors)	160
Micrococcal Nuclease, Mung Bean Nuclease	113	Monarch DNA Gel Extraction Kit	140	NEBNext Multiplex Oligos (96 Unique Dual Index Primer Pairs Sets 1-4)	160
Nuclease BAL-31	113	Monarch PCR & DNA Cleanup Kit (5 µg)	139	NEBNext Multiplex Oligos for Enzymatic Methyl-seq (Unique Dual Index Primer Pairs)	160
Nuclease P1, RecJ _h , T5 Exonuclease	114	Companion Products		NEBNext Multiplex Oligos (Dual Index Primers Sets 1 & 2)	160
T7 Exonuclease, Nucleoside Digestion Mix	115	Monarch DNA Cleanup Columns (5 µg)	139, 140	NEBNext Multiplex Oligos (Index Primers Sets 1-4)	160
DNA Repair Enzymes		Monarch Gel Dissolving Buffer	140	NEBNext Multiplex Oligos (96 Index Primers)	160
Properties of DNA Repair Enzymes and Structure-specific Endonucleases	116	Monarch DNA Cleanup Binding Buffer	139	NEBNext Multiplex Oligos (Methylated Adaptor, Index Primers Set 1)	160
DNA Repair Enzymes on Damaged and Non-Standard Bases	117	Monarch DNA Wash Buffer	139, 140	NEBNext Adaptor Dilution Buffer	160
APE 1, Endonuclease IV, <i>Tth</i> Endonuclease IV	118	Monarch DNA Elution Buffer	139, 140	NEBNext ARTIC Products for SARS-CoV-2 Sequencing	
Endonuclease III (Nth), <i>Tma</i> Endonuclease III	119	Plasmid Purification		NEBNext ARTIC SARS-CoV-2 FS Library Prep Kit (Illumina)	161
Endonuclease V	119	Monarch Plasmid Miniprep Kit	140	NEBNext ARTIC SARS-CoV-2 Library Prep Kit (Illumina)	161
Endonuclease VIII	120	Companion Products		NEBNext ARTIC SARS-CoV-2 Companion Kit (Oxford Nanopore Technologies)	161
Fpg, hAAG	120	Monarch Plasmid Miniprep Columns	140	Reagents for Oxford Nanopore Technologies Sequencing	
T4 PDG (T4 Endonuclease V)	121	Monarch Plasmid Buffers (Resuspension (B1), Lysis (B2), Neutralization (B3), Wash 1, Wash 2)	140	NEBNext Companion Module for Oxford Nanopore Technologies Ligation Sequencing	162
Uracil-DNA Glycosylase (UDG) & <i>Afu</i> (UDG)	121	Monarch DNA Elution Buffer	140	NEBNext ARTIC SARS-CoV-2 Companion Kit (Oxford Nanopore Technologies)	162
Antarctic Thermolabile UDG	122	Accessories		Target Enrichment	
hSMUG1, T7 Endonuclease I	122	Monarch Microfuge Tube EcoRack	141	NEBNext Direct Custom Ready Panels	164
PreCR Repair Mix	82	Monarch Pestle Set, Monarch 2 ml Tubes	136	NEBNext Direct Genotyping Solution Sample Prep Kits	165
Other Enzymes		Monarch Collection Tubes II	137	Immune Sequencing Kits (Human & Mouse)	
Cre Recombinase	123	NEBNext Reagents for Next Generation Sequencing		NEBNext Immune Sequencing Kits (Human)	166
TelIN Protelomerase, Topoisomerase I (<i>E. coli</i>), β-Agarase I	123	Reagents for Illumina Sequencing		NEBNext Immune Sequencing Kits (Mouse)	166
DNA Methyltransferases		DNA Product Selection Chart	146	Library Quantitation	
CpG & GpC Methyltransferases	124	RNA Product Selection Chart	147	NEBNext Library Quant Kit for Illumina	170
AluI, BamHI Methyltransferase	124	Ultra II DNA		NEBNext Library Dilution Buffer	170
<i>dam</i> , EcoGII, EcoRI, HaeIII, HhaI Methyltransferase	125	NEBNext Ultra II DNA Library Prep Kits +/- Sample Purification Beads	149		
HpaII, MspI, <i>TaqI</i> Methyltransferase	126	NEBNext Ultra II FS DNA Library Prep Kits +/- Sample Purification Beads	149		
Human Dnmt1 Methyltransferase	286	NEBNext Ultra II DNA & FS DNA PCR-Free Library Prep Kits +/- Sample Purification Beads	149		
ssDNA Binding Proteins		Methylation Analysis			
RecA, RecA ₄ , T4 Gene 32 Protein	127	NEBNext Enzymatic Methyl-seq Kit	151		
ET SSB (Extreme Thermostable ssBinding Protein)	127	NEBNext Enzymatic Methyl-seq Conversion Module	151		
Cloning Plasmids & DNAs		NEBNext Q5U Master Mix	151		
Cloning Plasmids (pBR322, pUC19)	128	NEBNext Multiplex Oligos for Enzymatic Methyl-seq (Unique Dual Index Primer Pairs)	151		
DNAs (M13mp18, Lambda, φX174)	128				
Other					
M13K07 Helper Phage, <i>Tth</i> Argonaute (TiAgo)	128				
Competent Cell Selection Chart for Cloning	129				
Monarch Nucleic Acid Purification Kits	129				
Nucleic Acid Purification					
Genomic DNA Purification					
Monarch Genomic DNA Purification Kit	135				
Companion Products					
Monarch gDNA Purification Columns	135				
Monarch Collection Tubes II	135				

Table of Contents (Cont.)

Reagents for Ion Torrent

NEBNext Fast DNA Library Prep Set for Ion Torrent	171
NEBNext Fast DNA Fragmentation & Library Prep Set for Ion Torrent	171

Reagents Suitable for Any Platform

DNA Repair

NEBNext FFPE DNA Repair Mix	152
-----------------------------	-----

DNA Separation

NEBNext Microbiome DNA Enrichment Kit	167
---------------------------------------	-----

Modules & Enzymes

NEBNext Magnesium RNA Fragmentation Module	154
NEBNext dsDNA Fragmentase	168
NEBNext dsDNA Fragmentase Reaction Buffer v2	168
NEBNext Ultra II Q5 Master Mix	169

Magnetic Separation

NEBNext Magnetic Separation Rack	159
----------------------------------	-----

Markers and Ladders (DNA, RNA & Protein)

DNA Markers & Ladders

Size Ranges of DNA Ladders	180
Gel Loading Dye, Purple (6X) (with or without SDS)	180

DNA Ladders

1 kb, 100 bp, 1 kb Plus, 50 bp DNA Ladders	181
Low Molecular Weight DNA Ladder	181
PCR Marker	181
Supercoiled DNA Ladder	184

Quick-Load Purple Formats

Quick-Load Purple 1 kb Plus DNA Ladder	182
Quick-Load Purple 1 kb, 100 bp, 50 bp DNA Ladders	182
Quick-Load Purple Low Molecular Weight DNA Ladder	182

Other Ready-to-Load Formats

1 kb Plus DNA Ladder for Safe Stains	182
Fast DNA Ladder	182

TriDye Formats

TriDye DNA Ladders	182
TriDye Ultra Low Range DNA Ladder	182

Quick-Load Formats

Quick-Load DNA Ladders	182
Quick-Load 1 kb Extend DNA Ladder	182

PFG Ladders

Lambda PFG, MidRange PFG Ladders	182
Lambda DNA-Mono Cut Mix	182

Conventional DNA Markers

Lambda DNA (HindIII Digest, BstEII Digest)	183
φX174 DNA-HaeIII Digest	183
pBR322 DNA (BstNI Digest, MspI Digest)	183

RNA Markers & Ladders

dsRNA, ssRNA, Low Range ssRNA Ladder	184
microRNA Marker	184

Protein Standards

Unstained Protein Standard, Broad Range (10–200 kDa)	185
Color Prestained Protein Standard, Broad Range (10–250 kDa)	185
Blue Prestained Protein Standard, Broad Range (11–250 kDa)	185

Genome Editing

Featured Products Supporting CRISPR Workflows

EnGen Spy Cas9 NLS	190, 191
EnGen Mutation Detection Kit	190, 192
EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i>	190, 193
EnGen Spy Cas9 Nickase	190, 191
EnGen Spy dCas9 (SNAP-tag)	190, 191
EnGen Lba Cas12a (Cpf1)	190, 191
EnGen Sau Cas9	190, 191
<i>Th</i> Argonaute (TiAgo)	191
Cas9 Nuclease, <i>S. pyogenes</i>	190, 191
Monarch Total RNA Miniprep Kit	137
Monarch RNA Cleanup Kit (50 µg)	138
Q5 Site-Directed Mutagenesis Kit	

(with or without competent cells)	98
Q5 High-Fidelity DNA Polymerases	66
NEBuilder HiFi DNA Assembly Master Mix	94
NEBuilder HiFi DNA Assembly Cloning Kit	94
HiScribe T7 ARCA mRNA Kit (with or without tailing)	200
HiScribe T7 High Yield RNA Synthesis Kit	199
HiScribe T7 Quick High Yield RNA Synthesis Kit	199
T7 Endonuclease I	122

RNA Reagents

Avoiding Ribonuclease Contamination	198
-------------------------------------	-----

RNA Synthesis

HiScribe T7 Quick High Yield RNA Synthesis Kit	199
HiScribe T7 High Yield RNA Synthesis Kit	199
HiScribe SP6 RNA Synthesis Kit	199
HiScribe T7 ARCA mRNA Kit (with or without tailing)	200
EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i>	193
Recommended HiScribe Kits by Application	201
T3, T7, SP6 RNA Polymerases	202
Hi-T7 RNA Polymerase (Std. & High Conc.)	202
<i>E. coli</i> Poly(A) Polymerase	202
Poly(U) Polymerase	202
<i>E. coli</i> RNA Polymerase, (Core and Holoenzyme)	203
Pyrophosphatase, Inorganic (<i>E. coli</i> , yeast)	203
Thermostable Inorganic Pyrophosphatase	203
NudC Pyrophosphatase	203
Ribonucleotide Solution Set, Solution Mix	203
Sce Pst1	203
Vaccinia Capping System	204
mRNA Cap 2'- <i>O</i> -Methyltransferase	204
RNA Cap Analog Selection Chart	204
3'-Desthiobiotin-GTP, 3'-Biotin-GTP	205
yDcpS, mRNA Decapping Enzyme	205

cDNA Synthesis, RT-PCR & RT-qPCR

cDNA Synthesis Selection Chart	206
ProtoScript II Reverse Transcriptase	206
LunaScript RT SuperMix Kit	73
LunaScript RT SuperMix Kit (Primer-free)	73
M-MuLV Reverse Transcriptase	207
AMV Reverse Transcriptase	207
WarmStart RTx Reverse Transcriptase	207
Template Switching RT Enzyme Mix	208
Primers for cDNA Synthesis	208
ProtoScript II First Strand cDNA Synthesis Kit	209
ProtoScript First Strand cDNA Synthesis Kit	209
Luna Probe One-Step RT-qPCR Kit (No ROX)	74
Luna Probe One-Step RT-qPCR 4X Mix with UDG	74
Luna Cell Ready One-Step RT-qPCR Kit	75
Luna Cell Ready Probe One-Step RT-qPCR Kit	75
Luna Universal One-Step RT-qPCR Kit	74
Luna Universal Probe One-Step RT-qPCR Kit	74
One Taq One-Step RT-PCR Kit, One Taq RT-PCR Kit	81
Luna SARS-CoV-2 RT-qPCR Multiplex Assay Kit	76

RNA Ligases & Modifying Enzymes

RNA Ligase Activity Chart	210
RNA Ligase Selection Chart	211
T4 RNA Ligase 1, 2	212
T4 RNA Ligase 2, truncated	212
T4 RNA Ligase 2, truncated K227Q and KQ	213
RtcB Ligase, Thermostable 5' AppDNA/RNA Ligase	213
5' DNA Adenylation Kit	214
SplintR Ligase	214
RNA 5' Pyrophosphohydrolase (RppH)	214
5' Deadenylase	214
RNase I _h , RNase H, Thermostable RNase H	215
RNase HII	215
Quick CIP	105
Antarctic Phosphatase	106
Shrimp Alkaline Phosphatase (rSAP)	106
T4 Polynucleotide Kinase	104
ShortCut RNase III, XRN-1	216

RNase Control

Exonuclease T, Nucleoside Digestion Mix	217
DNase I (RNase-free)	217
RNase Inhibitor, Murine	217

RNase Inhibitor, Human Placenta	218
Ribonucleoside Vanadyl Complex	218

NEBNext Reagents for RNA Library Preparation

NEBNext Ultra II Directional RNA Library Prep Kit for Illumina	153
NEBNext Ultra II Directional RNA Library Prep with Sample Purification Beads	153
NEBNext Ultra II RNA Library Prep Kit for Illumina	153
NEBNext Ultra II RNA Library Prep with Sample Purification Beads	153
NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina	158
NEBNext Single Cell/Low Input cDNA Synthesis & Amplification Module	158

RNA Cleanup & Isolation

Magnetic mRNA Isolation Kit	219
Oligo d(T) ₂₅ Magnetic Beads	219
Streptavidin Magnetic Beads	219
EpiMark N6-Methyladenosine Enrichment Kit	220
Monarch Total RNA Miniprep Kit	137
Monarch RNA Cleanup Kits (10, 50, 500 µg)	138

RNA Markers & Ladders

dsRNA Ladder	221
microRNA Marker	221
ssRNA, Low Range ssRNA Ladder	221
RNA Loading Dye (2X)	221
Universal miRNA Cloning Linker	221

Protein Expression & Purification Technologies

Protein Expression Overview	226
-----------------------------	-----

Expression Systems

Cell-free Expression

NEBExpress Cell-free <i>E. coli</i> Protein Synthesis System	227
NEBExpress GamS Nuclease Inhibitor	227
PURExpress <i>In Vitro</i> Protein Synthesis Kit	228
PURExpress Δ Ribosome, Δ (aa, tRNA), Δ RF123 Kits	228

PURExpress Companion Products:

PURExpress Disulfide Bond Enhancer	228
<i>E. coli</i> Ribosome	228

E. coli

NEBExpress MBP Fusion and Purification System	229
---	-----

NEBExpress Companion Products:

Amylose Resin, Anti-MBP Monoclonal Antibody	229
pMAL-c6T Vector, TEV Protease	229

IMPACT Kit, IMPACT Vector Comparison	230
--------------------------------------	-----

IMPACT Companion Products:

Chitin Resin, Anti-CBD Monoclonal Antibody	230
IMPACT Vectors	230

Yeast

<i>K. lactis</i> Protein Expression Kit	231
---	-----

K. lactis Companion Products:

SacII	231
Yeast Carbon Base Medium Powder	231
<i>K. lactis</i> GG79 Competent Cells	231

Competent Cells

Competent Cells for Protein Expression	232, 245–249
--	--------------

Purification

Purification Beads, Columns & Resins Selection Chart	233
--	-----

Polyhistidine-tagged Protein Purification

NEBExpress Ni-NTA Magnetic Beads	234
NEBExpress Ni Resin, NEBExpress Ni Spin Column	234
TEV Protease	234

Maltose Binding Protein (MBP) Purification

Amylose Resin, Amylose Resin High Flow	234
Amylose, Anti-MBP Magnetic Beads	234
Anti-MBP Monoclonal Antibody	234

Chitin Binding Domain (CBD) Purification

Chitin Resin, Chitin Magnetic Beads	234
Anti-CBD Monoclonal Antibody	234

Magnetic Bead Purification Products	
Oligo d(T) ₂₅ and Streptavidin Magnetic Beads	235
Magnetic mRNA Isolation Kit	235
Hydrophilic Streptavidin Magnetic Beads	235
Protein A, Protein G Magnetic Beads	235
Goat Anti-Mouse, Anti-Rabbit, Anti-Rat IgG Magnetic Beads	235
6-Tube, 50 ml, 12-Tube Magnetic Separation Racks	235
96-Well Microtiter Plate Magnetic Separation Rack	235
NEBNext Magnetic Separation Rack	235

Competent Cells

Strain Properties Chart	240
-------------------------	-----

Cloning Strains

NEB Cloning Competent <i>E. coli</i> Sampler	242
NEB Turbo Competent <i>E. coli</i> (High Efficiency)	242
NEB 10-beta Competent <i>E. coli</i> (High Efficiency)	243
Electrocompetent	243
NEB 5-alpha Competent <i>E. coli</i> (High Efficiency)	243
(Subcloning Efficiency)	243
Electrocompetent	243
NEB 5-alpha F ['] /r Competent <i>E. coli</i> (High Efficiency)	243
NEB Stable Competent <i>E. coli</i> (High Efficiency)	244
dam-/dcm- Competent <i>E. coli</i>	244
NEB Tube Opener	244

Companion Products:

SOC Outgrowth Medium	242
NEB 10-beta/Stable Outgrowth Medium	243

Protein Expression Strains

BL21, BL21(DE3), Lemo21(DE3) Competent <i>E. coli</i>	245
NiCo21(DE3) Competent <i>E. coli</i>	246
NEBExpress Competent <i>E. coli</i> (High Efficiency)	246
NEBExpress I ^r Competent <i>E. coli</i> (High Efficiency)	246
T7 Express Competent <i>E. coli</i> (High Efficiency)	247
T7 Express <i>lysY</i> Competent <i>E. coli</i> (High Efficiency)	247
T7 Express <i>lysY</i> /r Competent <i>E. coli</i> (High Efficiency)	247
SHuffle Express Competent <i>E. coli</i>	248
SHuffle T7 Express Competent <i>E. coli</i>	248
SHuffle T7 Competent <i>E. coli</i>	249
SHuffle T7 Express <i>lysY</i> Competent <i>E. coli</i>	249
How do I express my protein in SHuffle cells?	249

Glycobiology & Protein Tools

Endoglycosidases

Endo F2, Endo F3	254
Endo D, Endo H, Endo H ₁ , Endo S	255
PNGase A, PNGase F and PNGase F, recombinant, RNase B	256
Rapid PNGase F and Rapid PNGase F (non-reducing format)	257
Rapid PNGase F Antibody Standard	257
Remove-iT PNGase F	257
O-Glycosidase	258
Protein Deglycosylation Mix II	258
Fetuin	258

Exoglycosidases

α-N-Acetylglucosaminidase, β-N-Acetylglucosaminidase S	259
β-N-Acetylhexosaminidase, α1-2 Fucosidase	259
α1-3,4 Fucosidase, α1-2,3,4,6 Fucosidase	260
α1-2,4,6 Fucosidase O, α1-3,6 Galactosidase	260
α1-3,4,6 Galactosidase, β1-3 Galactosidase	261
β1-3,4 Galactosidase, β1-4 Galactosidase S	261
α1-2,3 Mannosidase, α1-2,3,6 Mannosidase	262
α1-6 Mannosidase, α2-3,6,8 Neuraminidase	262
α2-3,6,8,9 Neuraminidase A, α2-3 Neuraminidase S	263

Heparin Lyases

<i>Bacteroides</i> Heparinase I, II, III	264
--	-----

Proteases

Endoglycoceramidase I (EGCase I)	265
IdeZ Protease, O-Glycoprotease	265
Trypsin-ultra, Mass Spectrometry Grade	266
α-Lytic Protease	266
Endoproteinase LysC, GluC, AspN	267
Trypsin-digested BSA MS Standard (CAM Modified)	267
Proteinase K, Molecular Biology Grade	268
Thermolabile Proteinase K, TEV Protease	268

Factor Xa Protease	269
Enterokinase, light chain	269
Furin	269

Protein Phosphatases & Kinases

Lambda Protein Phosphatase (Lambda PP)	270
p-Nitrophenylphosphate (PNPP)	270
Sodium Orthovanadate	270
cAMP-dependent Protein Kinase (PKA)	271
Casein Kinase II (CK2)	271
Adenosine 5' Triphosphate (ATP)	271

Phage Display

Ph.D.-7, -12, -C7C Phage Display Peptide Library Kits	272
Ph.D.-12 Phage Display Peptide Library	272
Ph.D. Peptide Display Cloning System	273
Anti-M13 pIII Monoclonal Antibody	273

Epigenetics

Epigenetics

EpiMark 5-hmC and 5-mC Analysis Kit	278
T4 Phage β-glucosyltransferase	278
EpiMark Methylated DNA Enrichment Kit	279
EpiMark N6-Methyladenosine Enrichment Kit	220, 279
EpiMark Bisulfite Conversion Kit	279
EpiMark Hot Start Taq DNA Polymerase	71, 279

Methylation-Dependent Restriction Enzymes

AbaSI, FspEI, LpnPI, MspJI	280
McrBC	281

Additional Restriction Enzymes for Epigenetic Analysis

DpnI, DpnII	36
HpaII	40
MspI	43

Histones

EpiMark Nucleosome Assembly Kit	283
Nucleosome Control DNA	283
Histone H3.1/H4 Tetramer Human, Recombinant	283
Histone H2A/H2B Dimer Human, Recombinant	283
Histone H1 ⁹ Human, Recombinant	283
Histone H2A, H2B, H3.1 Human, Recombinant	284
Histone H3.2, H3.3, H4 Human, Recombinant	285

DNA Methyltransferases

CpG Methyltransferase (M.SssI)	124
GpC Methyltransferase (M.CviPI)	124
AluI, BamHI Methyltransferases	124
dam, EcoGII, EcoRI, HaeIII, HhaI Methyltransferases	125
HpaII, MspI, TaqI Methyltransferases	126
Human DNA (cytosine-5) Methyltransferase (Dnmt1)	286

Control DNA

HeLa Genomic DNA	287
5-methyl-dCTP	287

NEBNext Reagents for ChIP-Seq Library Prep

NEBNext Ultra II DNA Library Prep Kit for Illumina	149
NEBNext Ultra II DNA Library Prep with Sample Purification Beads	149
NEBNext Multiplex Oligos	160
NEBNext Ultra II End Repair/da-Tailing Module	150
NEBNext Ultra II Ligation Module	150
NEBNext Ultra II Q5 Master Mix	150
NEBNext Enzymatic Methyl-seq Kit	151, 289
NEBNext Enzymatic Methyl-seq Conversion Module	151, 289
NEBNext Q5U Master Mix	289
NEBNext Multiplex Oligos for Enzymatic Methyl-seq (Unique Dual Index Primer Pairs)	289

Cellular Analysis

Cellular Imaging & Analysis

Overview of SNAP- and CLIP-tag Technologies	294
Comparison of SNAP-tag/CLIP-tag Technology to GFP	295
Fluorescent Substrates	296
Blocking Agents, SNAP-tag Antibodies, Purified Protein	297
Cloning Vectors, Biotin Labels	298
SNAP-Capture, Building Blocks	299

Reference Appendix

Online Interactive Tools	304
Databases, Mobile Apps & Freezer Program Locator	305

Restriction Endonucleases

Optimizing Restriction Enzyme Reactions	306
Double Digestion	307
Types of Restriction Enzymes	307
Troubleshooting Guide	308

Performance

Performance Chart for Restriction Enzymes	309–314
Activity of Enzymes at 37°C	315
Activity of DNA Modifying Enzymes in CutSmart Buffer	315
Tips for Avoiding Star Activity	316
High-Fidelity (HF) Restriction Enzymes	316

Properties of Restriction Enzymes

Isoschizomers	317
Survival in a Reaction	317–318
Cleavage of Supercoiled DNA	318

Generating New Cleavage Sites

Compatible Cohesive Ends	319–320
--------------------------	---------

Methylation Sensitivity

Dam, Dcm and CpG Methylation	321
------------------------------	-----

PCR, qPCR & RT-qPCR

General Guidelines for PCR Optimization	322
PCR Troubleshooting Guide	323
Optimization Tips for Luna qPCR	324
Optimization Tips for Luna One-Step RT-qPCR	325
Luna qPCR Troubleshooting Guide	326
Luna One-Step RT-qPCR Troubleshooting Guide	327
Cleavage Close to the End of DNA Fragments	328
Activity of Restriction Enzymes in PCR Buffers	329–330

Cloning & DNA Assembly

Getting Started with Molecular Cloning	331
Optimization Tips for Golden Gate Assembly	332
Optimization Tips for NEBuilder HiFi DNA Assembly	333
Traditional Cloning Quick Guide	334–335
Troubleshooting Guide for Cloning	336–339
Optimization Tips for Your Cloning Reactions	339–340

Nucleic Acid Purification

Troubleshooting Guide for DNA Cleanup & Plasmid Purification Using Monarch Kits	341
Choosing Sample Input Amounts When Using the Monarch Genomic DNA Purification Kit	342
Troubleshooting Guide for Genomic DNA Purification Using Monarch Kit	343
Choosing Sample Input Amounts When Using the Monarch Total RNA Miniprep Kit	344
Troubleshooting Guide for Total RNA Extraction & Purification using Monarch Kits	345
Guidelines for Choosing Sample Input Amounts When Using the Monarch HMW DNA Kit for Cells & Blood	346
Guidelines for Choosing Sample Input Amounts When Using the Monarch HMW DNA Extraction Kit for Tissue, Bacteria and Other Samples	347

Transformation

Genetic Markers	348–349
Enhancing Transformation Efficiency	350
Protein Expression with T7 Express Strains	351

NGS Library Preparation

Troubleshooting for Guide for NEBNext DNA Lib. Prep Kits	352
Guidelines for NGS Library Prep	353

Protein Expression

Labeling with SNAP-tag Technology Troubleshooting Guide	354
DNA Maps	355–365

Odds and Ends

General Molecular Biology Information	366–369
Index	370

Application Reference Guide

DNA Cloning

A-tailing	84
Alkaline Elution	118–120
Alkaline Unwinding	118–120
Chromosomal DNA Removal	109–113
Cleavage Close to the Ends	328
Cloning Optimization	339–340
Cloning Plasmids	128
Cloning Quick Guide	334–345
Cloning Troubleshooting Guide	336–339
Cloning Workflow	92
DNA Assembly	94–97
DNA Cleanup	87, 129, 139–140
DNA Digestion	18–59
DNA Extraction from Agarose	123, 129, 140
DNA Ligase Selection Chart	99–100
DNA Maps	355–365
DNA Markers and Ladders	178–185
DNA Repair	
Enzymes	82, 116–122
Polymerases	82
Removal	110, 217
Repair Mix	82, 152
Electrophoresis	
Loading Dye	59, 180
Markers & Ladders	178–185
End Treatments	
Blunting	83–84, 98, 112–113
dA-Tailing	84
Dephosphorylation	105–106
End Labeling	104
Phosphorylation	104
Genomic DNA Purification	135–136
Gibson Assembly	94–95
Golden Gate Assembly	96
HMW DNA Extraction	136
Library Construction	144–175, 218, 288
Ligase Chain Reaction	103–104
Ligase Detection Reaction	103–104
Ligation	99–104, 214
Selection Chart	99–100
Ligation Independent Cloning	92
Methyltransferases	124–126, 286
Nick Repair	101–104
Nicking Endonucleases	19, 55–57
Nucleic Acid Purification	87, 129, 132–141
Padlock Probes	103
PCR Cloning	92, 97
Plasmid Purification	129, 140
Polymerases	66–71
Recombination	123, 127
Recombinational Cloning	92
Restriction Digests	
Activity in PCR Buffers	329–330
Double Digestion	307
High Fidelity Restriction Enzymes	19, 316
Methylation Dependent	280–281
Methylation Sensitivity	309–314, 321
Optimization	306
Restriction Enzymes	18–59
Time Saver Qualified Restriction Enzymes	309–314, 306
Troubleshooting	308
Seamless Cloning	92, 94–97
Single Cell Gel Electrophoresis	119–121
Site Directed Mutagenesis	94, 98, 111

Site Mapping	113
SNP Analysis	82, 106, 109
Synthetic Biology	94–97
T/A Cloning	92, 101–104
Transcriptional Mapping	113
Transformation Tips	350–351
Transformation/Competent cells	129, 238–249
USER Cloning	97
Vectors	128
Vector Maps	355–365

DNA Amplification & PCR

cDNA Synthesis	86, 206–209
Chain Termination	82, 85
Colorimetric LAMP	76–78
DNA Polymerases	63–87
Damaged Templates	82, 152
DNA Removal	105–106, 109–115
DNA Repair	
Polymerases	82
Repair Mix	82, 152
DNA Sequencing	144–175
Genome Amplification	80
Helicase Dependent Amplification	80
Incorporation of Modified Nucleotides	82
Isothermal Amplification	76–80
LAMP	76–78
Lesion Bypass	82
Library Construction	144–175
Molecular Diagnostics-based	64
Next Generation Sequencing	144–175
Nick Translation	82–83
Nucleic Acid Purification	87, 132–141
Nucleotides	85
One-Step RT-PCR	81
PCR	
AT-Rich Targets	65–68, 71
Cleanup	87, 139
Colony PCR	65, 68–69
Extraction-free PCR	65, 70
Fast	65, 66–67
GC-Rich Targets	65, 66, 68, 70
High-Fidelity	65, 66–67
High Throughput	65–69, 71
Long Amplicon	65, 66–67, 70
Multiplex PCR	65, 66–69
Optimization	322
Routine	65–70
Selection Chart	65
Troubleshooting	323
Probe Labeling	83
qPCR	72–73
Optimization	324
Troubleshooting	326
RT-LAMP	78
RT-PCR	81
RT-qPCR	72, 73–75
Optimization	325
Troubleshooting	327
SARS-CoV-2 Detection	76–77
Second Strand cDNA Synthesis	83
SNP Analysis	82, 287
ssDNA Binding Proteins	127
Strand Displacement Amplification	79, 80, 84
Nicking Enzymes	19, 55–57

Site Directed Mutagenesis	94, 98, 111
Two Step RT-PCR	73, 81, 208
Whole Genome Amplification	80

Genome Editing

Cas9-based Genome Editing	190–191
Mutation Detection	190, 192
sgRNA Synthesis	190, 193
Programmable Nucleases	191

Epigenetics

5-mC & 5-hmC Detection, Enrichment & Quantitation	278–279, 289
5-hmC Labeling	278
Bisulfite Sequencing	279, 289
ChIP Assay	283
ChIP-Seq Sample Preparation	149–150, 288
Combined Bisulfide Restriction Analysis (COBRA)	287
Control DNA Substrates	287
Glucosylation of 5-hmC	278
HAT Assay	283
Histone Modification Studies	282–285
Immunodetection of 5-hmC	278
m6A Enrichment	220, 279
Methyl-Seq	151, 289
Methylation Dependent Cleavage	280–281
Methylation Sensitive Restriction Enzyme Cleavage	280, 309–314, 321
Methylation Sensitive Single Nucleotide Primer Extension (Ms-SNuPE)	287
Methylation Specific PCR (MSP)	287
Methyltransferases	124–126, 286
Nucleosome Assembly	282–285
SNP Analysis	287
Southern Blotting	287

RNA Analysis

Adaptor Ligation	210–213
Blunting	217
Capping	24–205
cDNA Synthesis	86, 206–209
Dephosphorylation	105–106, 216
m6A Enrichment	219–220
Gene Silencing	216
In vitro Translation	199–203
Isolation and Enrichment	219–220
Labeling	202, 204, 212
Library Construction	147, 153, 158, 218
Ligation	210–214
Activity Chart	210
Selection Chart	211
Loading Dye	221
mRNA Isolation	219–220
Phosphorylation	104, 216
Poly(A) Tailing	202
Poly(U) Tailing	202
Probe Generation	202
RACE	86, 206, 208
Reverse Transcription	86, 206–209
RNAi	216
RNA Cleanup	137–138, 220
RNA Labeling	202
RNA Library Preparation	147, 153, 158, 218

RNA Markers and Ladders	184, 221
RNA Miniprep	137, 220
RNA Removal	217–218
RNA Synthesis	199–203
RNase Inhibitors	217–218
RT-PCR	81
RT-qPCR	72, 73–75
Sequencing	147, 153, 158, 218
siRNA	216
Tailing	202

Sample Prep for Next Gen Sequencing

Adaptors & Primers	160
Adenylation of ssDNA Linkers	214
ARTIC Workflow	161
Barcodes	160
Enzymatic Fragmentation	149, 150, 151, 161, 289
Enzymatic Methyl-seq	151, 289
FFPE-treated DNA	152
Library Construction	
Amplification	169
ChIP-Seq	149–150, 288
DNA	146, 148–152, 162, 167–171

dsDNA Fragmentation	149, 150, 161, 168
mRNA	143, 153–157, 161
rRNA Depletion	155–157
Single Cell RNA	158
Small RNA	159
Library Quantitation	170
Magnetic Separation	159
Microbiome Enrichment	167
Oxford Nanopore Sequencing	162
Sample Preparation	144–175
SARS-CoV-2 Sequencing	161–162
Target Enrichment	163–165

Protein Expression & Analysis

Expression	224–235
Bacterial	229–230
Cell-free	227–228
Difficult to Express	
Proteins	227–229, 232, 245–249
Yeast	231
Glycosylation	252–263
Endoglycosidases	254–258
Exoglycosidases	259–263
MS Standards	267

Peptide Ligation	230
Phage Display	272–273
Proteases	265–269
Protein Labeling	236, 292–299
Protein Standards	182
Protein Purification	
Chitin Binding Domain (CBD) Purification	234
Expression & Purification Kits	226–231
Magnetic Matrices and Racks	235
Maltose Binding Protein (MBP) Purification	234
Polyhistidine-tagged Purification	234
Pull-Down Assays	299
Transformation	
Competent Cells	232, 245–249
Troubleshooting & Tips	350–351
Yeast	231
Strain Properties	240

Cellular Analysis

Protein Labeling (e.g., SNAP-tag)	292–299
Protein Localization	292–299
Pulse Chase	292–299
Receptor Internalization	292–299
Protein Pull-Down Assays	299



View of the moody sky over Thailand at twilight.
Credit: Boy Anupong, Getty images

Protecting a Legacy

Each edition of the New England Biolabs Catalog contains a collection of mini-reviews that addresses various scientific, environmental and/or humanitarian topics. This year, we are dedicating the Catalog to our Founder, Donald G. Comb, and sharing some of the values that he was passionate about.


16 By Scientists for Scientists	176 Be Adventurous	250 Mostly Mozart
60 Conserving Biodiversity	186 Advancement of Science	274 Promoting Environmental Sustainability
88 Ask Big Questions	194 Recognizing Those Making a Difference	290 Giving Back
130 Putting Recycling to Practice	222 The Intersection of Art and Science	300 And the NEB Catalog
142 Neglected Tropical Diseases	236 Sustainable Growth	

Mini-reviews were authored by Joanne Gibson, Lydia Morrison and Deana Martin. Joanne is a Marketing Communications Writer at NEB. She received her Ph.D. in Molecular Biology from the University of Sydney, Australia. Lydia is also a Marketing Communications Writer at NEB, and manages our social media program. She received her M.S. in Biochemistry and Biophysics from the University of North Carolina at Chapel Hill School of Medicine. Deana is the Associate Director of Marketing Communications. She received her Ph.D. in Biochemistry from Boston College, Chestnut Hill, MA.

NEB's first home on Cabot St.
in Beverly, MA



NEB's Founder,
Donald Comb
(1927-2020)



NEB group photo, 2019



NEB laboratory facility
in Ipswich, MA



NEW ENGLAND
Biolabs
281 CABOT ST. BEVERLY, MA 01915 (978) 927-5054

RESTRICTION ENDONUCLEASES
NUCLEIC ACIDS
RNA's
DNA LIGASE
BACTERIAL CELLS

1975 - 1976
CATALOG

By Scientists for Scientists

In 1974, a scientist named Don Comb left his associate professorship in the Biochemistry Department at Harvard Medical School — disenthralled with the traditional avenues of securing research funding and the rigidity of academia. He left with the goal of bringing together like-minded scientists to form a cooperative company — a company that would provide tools for molecular biologists. Further, the company would use its revenue to fund independent research projects, an idea that was unprecedented at the time. And thus, New England Biolabs® Inc. (NEB®) was founded, with the goal of enabling life science research, both within and outside the company.

In the early days, NEB's focus was to provide high quality restriction enzymes at a fair price. Over the years, NEB researchers have discovered a wide range of enzyme specificities — as a result NEB now offers >285 restriction enzymes, the largest commercially-available selection. We have also expanded our product offerings into areas related to PCR, gene expression, sample preparation for next generation sequencing, synthetic biology, glycobiology and RNA analysis.

Don hired scientists who were passionate and curious. He encouraged them to explore their own research interests while carrying out their job responsibilities. In doing so, NEB scientists became experts in the functionality and biophysical properties of a wide range of enzyme activities, resulting in an in-depth understanding of product performance. This enabled the collegial exchange of ideas, protocols and techniques with customers — the beginnings of a more personal approach to technical support that still continues today.


This expertise also enabled NEB to produce recombinant enzymes at commercial scale, dramatically improving quality and yield — greatly advancing the field of molecular biology. As yields increased, NEB was able to lower prices and pass savings on to customers.

NEB remains committed to research and has maintained an active program for over 45 years. Currently, more than 30 labs participate in basic research in areas ranging from molecular biology to parasitology, and this work has resulted in over 1,300 publications to date. Basic research and publications help to ensure that NEB scientists are staying curious, driving discovery and thinking like their customers. Additionally, NEB scientists are invited to speak at renowned academic institutions and conferences around the world, yet equally value the importance of advocating for science at local high schools. This emphasis on scientific collaboration and advancement has led to the sharing of innovative products and ideas worldwide.

Don's aspiration was for NEB to remain a small, private business that focused on steady, sustainable growth. His devotion to the advancement of science, stewardship of the environment and altruistic philanthropy has been a priority since day one. NEB employees are encouraged to pursue their passions, whether it be getting involved in groundbreaking research, helping out at a local science fair, sharing ideas to improve the sustainability of our business practices, or volunteering at a local food bank — everyone feels a responsibility to each other and to the community. NEB's CEO Jim Ellard said "Don's generosity and the faith that he had in each employee resulted in a family-like culture, where everyone has a voice and all employees feel valued."

A culture built by scientists, for scientists.

Learn about
NEB and the
photos shown.



NEB group photo, taken
at the Tozer Road Facility
in Beverly, MA, 1997



First NEB group photo, 1980/81



Entrance to
NEB laboratory facility
in Ipswich, MA

Restriction Endonucleases

The leader in the discovery & production of restriction enzymes.

Having supplied restriction enzymes to the research community for over 45 years, NEB has earned the reputation of being a leader in enzyme technologies. Working continuously to be worthy of that distinction, NEB strives to develop enzymes of the highest purity and unparalleled quality.


NEB scientists continue to improve our existing portfolio, as well as explore the utility of NEB reagents in new technologies. As a result, NEB scientists continue to publish scientific papers and be awarded grants in this area. With the industry's largest research and development group dedicated to restriction enzymes, we are proud to have been there first: the first to commercialize a recombinant enzyme, the first to introduce a nicking enzyme, and the first to supply a true restriction enzyme master mix. In addition, NEB has a continuing history of innovation by engineering restriction enzymes with altered specificities and improved performance. Through ongoing research in these areas, we are committed to driving the innovations that allow us to offer maximum performance and convenience.

Featured Tools and Resources

309 Performance/Activity Chart for Restriction Enzymes

306 Tips for Restriction Enzyme Optimization


308 Restriction Enzyme Troubleshooting Guide


 Visit NEBrestrictionenzymes.com to find additional online tools, video tech tips and tutorials to help you in setting up restriction enzyme reactions.





What are restriction enzymes?


Icon Descriptions


 The gene encoding this enzyme was cloned at NEB.

 This enzyme is purified from a recombinant source.

 This enzyme has been engineered for maximum performance.


 Time-Saver qualified enzymes will digest 1 µg of substrate DNA in 5–15 minutes using 1 µl of enzyme under recommended reaction conditions. These enzymes can also be used overnight with no loss of sample.


 Indicates that the restriction enzyme requires two or more sites for cleavage.







Indicates which reaction buffer is supplied with the enzyme for optimal activity. Enzymes with buffer requirements not met by one of the four standard NEBuffers are supplied with their own unique NEBuffer (NEB U). NEBuffers are color-coded (NEB r1.1 & 1.1 - yellow, NEB r2.1 & 2.1 - blue, NEB r3.1 & 3.1 - red, rCutSmart & CutSmart - green) and supplied as 10X stocks with each enzyme. For more information, consult the Performance Chart on pages 309–314.

 This enzyme is EpiMark validated for epigenetics studies.


 This enzyme is supplied with a separate tube of S-adenosylmethionine (SAM). To obtain 100% activity, SAM should be added to the 1X reaction mix as indicated. When required, a concentrated stock of SAM is supplied with the enzyme.


 This restriction enzyme is sensitive to dam, dcm, or CpG methylation. (Note that CpG methylation is applicable to eukaryotic genomic DNA only.) For more information, see page 321.

 Indicates whether or not the enzyme can be heat inactivated. Enzymes are first tested by incubation at 65°C for 20 minutes; any enzyme not inactivated at 65°C is then tested by incubation at 80°C for 20 minutes. If an enzyme can be heat inactivated, the temperature is indicated in the icon.



Indicates the enzyme's optimal incubation temperature.

 Indicates which diluent buffer (A, B or C) is recommended for making dilutions of restriction enzymes. For more information see pages 309–314.

 Enzymes for Innovation.

Note: Starting April 2021, NEB will begin the transition to supply all restriction enzymes with Recombinant Albumin-containing reaction buffer. For more details, see page 20 or visit www.neb.com/BSA-free

AatII	BfuAI	BstAPI	FspI	NciI	SapI
AbaSI	BglI	BstBI	FspEI	NcoI	Sau3AI
AccI	BglII	BstEII	HaeII	NcoI-HF	Sau96I
Acc65I	BipI	BstEII-HF	HaeIII	NdeI	SbfI
Acil	BmgBI	BstNI	HgaI	NgoMIV	SbfI-HF
AcII	Bmrl	BstUI	HhaI	NheI-HF	ScaI-HF
AcuI	BmtI	BstXI	HincII	NlaIII	ScrFI
AfeI	BmtI-HF	BstYI	HindIII	NlaIV	SexAI
AflII	BpmI	BstZ17I-HF	HindIII-HF	NmeAIII	SfaNI
AflIII	Bpu10I	Bsu36I	Hinfl	NotI	SfiI
AgeI	BpuEI	BtgI	HinPII	NotI-HF	SfiI
AgeI-HF	BsaI-HFv2	BtgZI	HpaI	NruI	SfiI
AhdI	BsaAI	BtsI-v2	HpaII	NruI-HF	SgrAI
AleI-v2	BsaBI	BtsIMutI	HphI	NsiI	SmaI
AluI	BsaHI	BtsCI	Hpy99I	Nsil-HF	SmlI
AlwI	BsaJI	Cac8I	Hpy166II	NspI	SnaBI
AlwNI	BsaWI	ClaI	Hpy188I	PacI	SpeI
Apal	BsaXI	CspCI	Hpy188III	PaeR7I	SpeI-HF
ApalI	BseRI	CviAII	HpyAV	PaqCI	SphI
ApeKI	BseYI	CviKI-1	HpyCH4III	PciI	SphI-HF
ApoI	BsgI	CviQ1	HpyCH4IV	PfIF	SrfI
ApoI-HF	BsiEI	DdeI	HpyCH4V	PHIM	SspI
AscI	BsiHKA	DpnI	KasI	PleI	SspI-HF
Asel	BsiWI	DpnII	KpnI	PluTI	StuI
AsiSI	BsiWI-HF	DraI	KpnI-HF	PmeI	StyI-HF
AvaI	BsII	DraIII-HF	LpnPI	PmlI	StyD4I
AvaII	BsmI	DrdI	MboI	PpuMI	Swal
AvrII	BsmAI	EaeI	MbolI	PshAI	TaqI-v2
BaeI	BsmBI-v2	EagI-HF	MfeI	Psil-v2	TfiI
BaeGI	BsmFI	EarI	MfeI-HF	PspGI	TseI
BamHI	BsoBI	EciI	MluI	PspOMI	Tsp45I
BamHI-HF	Bsp1286I	Eco53kI	MluI-HF	PspXI	TspMI
BanI	BspCNI	EcoNI	MluCI	PstI	TspRI
BanII	BspDI	EcoO109I	MlyI	PstI-HF	Tth111I
BbsI	BspEI	EcoP15I	MmeI	PvuI	XbaI
BbsI-HF	BspHI	EcoRI	MnlI	PvuI-HF	XcmI
BbvI	BspMI	EcoRI-HF	MscI	PvuII	XhoI
BbvCI	BspQI	EcoRV	MseI	PvuII-HF	XmaI
BccI	BsrI	EcoRV-HF	MslI	RsaI	XmnI
BceAI	BsrBI	Esp3I	MspI	RsrII	ZraI
BcgI	BsrDI	FatI	MspA1I	SacI	
BciVI	BsrFI-v2	FauI	MspJI	SacI-HF	
BclI	BsrGI	Fnu4HI	MwoI	SacII	
BclI-HF	BsrGI-HF	FokI	NaeI	Sall	
BcoDI	BssHII	FseI	NarI	Sall-HF	
Bfal	BssSI-v2				

High-Fidelity Restriction Enzymes

AgeI-HF, ApoI-HF, BamHI-HF, BbsI-HF, BclI-HF, BmtI-HF, BsaI-HFv2, BsiWI-HF, BsrGI-HF, BstEII-HF, BstZ17I-HF, DraIII-HF, EagI-HF, EcoRI-HF, EcoRV-HF, HindIII-HF, KpnI-HF, MfeI-HF, NcoI-HF, NheI-HF, NotI-HF, NruI-HF, Nsil-HF, PstI-HF, PvuI-HF, PvuII-HF, SacI-HF, Sall-HF, SbfI-HF, ScaI-HF, SpeI-HF, SphI-HF, SspI-HF, Styl-HF
21, 316

Methylation Sensitive Restriction Enzymes for Epigenetics Studies

AbaSI, DpnI, DpnII, FspEI, HpaII, LpnPI, McrBC, MspI, MspJI
280–281

Nicking Endonucleases

Nb.BbvCI, Nb.BsmI, Nb.BsrDI, Nb.BssSI, Nb.BtsI, Nt.AlwI, Nt.BbvCI, Nt.BsmAI, Nt.BspQI, Nt.BstNBI, Nt.CviPII
55–57

Homing Endonucleases

I-CeuI, I-SceI, PI-PspI, PI-SceI
57–58

NEBuffers, Diluents, Gel Loading Dyes, BSA & Recombinant Albumin, Molecular Biology Grade, NEB Tube Opener
58–59

Enzymes in green are all 100% active in CutSmart Buffer (see page 309–314).

One or more of these products are covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at gbd@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications. Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at <https://www.neb.com/support/terms-of-sale>. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

Looking to bring convenience to your workflow?

Speed up digestions with Time-Saver™ Qualified Restriction Enzymes


There are > 195 NEB restriction enzymes that can digest DNA in 5–15 minutes; many of which are supplied with rCutSmart or CutSmart Buffer or are High-Fidelity (HF®) Restriction Enzymes. If you prefer, you can also digest overnight with no unwanted star activity. All of our enzymes are rigorously tested for nuclease contamination. Only NEB can offer enzymes with the power to digest in 5–15 minutes, and the flexibility to withstand overnight digestions with no loss of sample.

For more information, visit www.neb.com/timesaver



Improve your analysis with our Purple Gel Loading Dye

Our Gel Loading Dye, Purple (6X), which is supplied with most restriction enzymes and all HF enzymes, sharpens bands and eliminates the UV shadow seen with other dyes. This solution contains SDS, which often results in sharper bands, as some restriction enzymes are known to remain bound to DNA following cleavage.



The Gel Loading Dye, Purple (6X) (Lane 1) included in the Quick-Load Purple 1 kb DNA Ladder does not cast a UV shadow over the underlying bands, unlike the Gel Loading Dye, Blue (6X) (Lane 2).

Simplify reaction setup and double digestion with rCutSmart™ Buffer

Over 215 enzymes are 100% active in a single buffer, rCutSmart Buffer, making it significantly easier to set up double digest reactions. Since rCutSmart Buffer includes Recombinant Albumin, there are fewer tubes and pipetting steps to worry about. Additionally, many DNA modifying enzymes are 100% functionally active in rCutSmart Buffer, eliminating the need for subsequent purification (see page 315). For more information, visit www.NEBCutSmart.com

Same high performance, now with BSA-free reaction buffer

BSA is commonly used to stabilize some proteins during reactions, and can also prevent adhesion to reaction tubes and pipette surfaces. We understand that there is an increased need for comparable performance using BSA-free reagents — beginning April 2021, we will be switching our current BSA-containing reaction buffers (NEBuffer 1.1, 2.1, 3.1 and CutSmart® Buffer) to Recombinant Albumin-containing buffers (NEBuffer r1.1, r2.1, r3.1 and rCutSmart™ Buffer). These buffers have been rigorously tested, and there is no difference in performance when using either buffer system. We anticipate this switch may take up to 6 months to complete. Over that time, you may receive product with BSA- or rAlbumin-containing buffer with your product — either will work for your reactions.

For more details visit www.neb.com/BSA-free


Looking to optimize performance in your reaction?

Choose a High-Fidelity (HF®) Restriction Enzyme

NEB High-Fidelity (HF) restriction enzymes have the same specificity as native enzymes, with the added benefits of reduced star activity, rapid digestion (5–15 minutes), and 100% activity in rCutSmart or CutSmart Buffer. Enjoy the improved performance of our engineered enzymes at the same price as the native enzymes!

For more information, visit www.neb.com/HF

EcoRI-HF (NEB #R3101) shows no star activity in overnight digests, even when used at higher concentrations. 50 µl rxns were set up using 1 µg of Lambda DNA, the indicated amount of enzyme and the recommended reaction buffer. Rxns were incubated overnight at 37°C. Marker M is the 1 kb DNA Ladder (NEB #N3232).




Benefit from industry-leading quality controls

NEB's reputation as a leader in enzyme technologies stems from the quality and reliability of our restriction enzymes. All of our restriction enzymes undergo stringent quality control testing, ensuring the highest levels of purity and lot-to-lot consistency.

For more information on quality at NEB, visit www.neb.com/quality

Restriction Enzyme Competitor Study: Nuclease Contamination



EcoRI, NotI, and BamHI from multiple suppliers were tested in reactions containing a fluorescent labeled single stranded, double stranded blunt, 3' overhang or 5' overhang containing oligonucleotides. The percent degradation is determined by capillary electrophoresis and peak analysis. The resolution is at the single nucleotide level.

Learn about the benefits of HF enzymes.



AatII

rCutSmart  RR  dII B 37°  CpG

#R0117S 500 units\$64
 #R0117L 2,500 units\$262

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity <10 50* 50 100

5'... G A C G T C ... 3'
 3'... C T G C A G ... 5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 20,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Note: *May exhibit star activity in this buffer.

AbaSI

rCutSmart  RR  Epi dII C 25° 

#R0665S 1,000 units\$113

See page 280 for more information.

AccI

rCutSmart  RR  dII A 37°  CpG

#R0161S 1,000 units\$75
 #R0161L 5,000 units\$308

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 50 50 10 100

5'... G T M K A C ... 3'
 3'... C A K M T G ... 5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

Acc65I

 RR  NEB r3.1 dII A 37°  dcm CpG

#R0599S 2,000 units\$70
 #R0599L 10,000 units\$287

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 10 75* 100 25

5'... G G T A C C ... 3'
 3'... C C A T G G ... 5'

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

Methylation Sensitivity: Blocked by some combinations of overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

Note: *May exhibit star activity in this buffer.

AciI

rCutSmart  RR  dII A 37°  CpG

#R0551S 200 units\$70
 #R0551L 1,000 units\$287

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity <10 25 100 100

5'... C C G G C ... 3'
 3'... G G C C G ... 5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

AclI

rCutSmart  RR  dII B 37°  CpG

#R0598S 300 units\$70
 #R0598L 1,500 units\$287

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity <10 <10 <10 100

5'... A A C G T T ... 3'
 3'... T T G C A A ... 5'

Reaction Conditions: rCutSmart Buffer, 37°C

Concentration: 5,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

AclI

rCutSmart  RR  dII B 37° 

#R0641S 300 units\$67
 #R0641L 1,500 units\$274

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 50 100 50 100

5'... C T G A A G (N)₁₆ ... 3'
 3'... G A C T T C (N)₁₄ ... 5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Star Activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

AfeI

rCutSmart  RR  dII B 37°  CpG

#R0652S 200 units\$72
 #R0652L 1,000 units\$296

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 25 100 25 100

5'... A G C G C T ... 3'
 3'... T C G C G A ... 5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

AflII

rCutSmart  RR  dII A 37° 

#R0520S 2,000 units\$67
 #R0520L 10,000 units\$274

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 50 100 10 100

5'... C T T A A G ... 3'
 3'... G A A T T C ... 5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 20,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

AflIII



#R0541S	250 units\$70
#R0541L	1,250 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	50	100	50

5'...A[▼]CRYGT...3'
3'...TGYRCA...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

AgeI



#R0552S	300 units\$74
#R0552L	1,500 units\$304

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	75	25	75

5'...A[▼]CCGGT...3'
3'...TGGCA...5'

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: NEBuffer r1.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Note: Star activity may result from extended digestion.

Concentration: 5,000 units/ml

AgeI-HF[®]

#R3552S	300 units\$74
#R3552L	1,500 units\$304

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	50	10	100

5'...A[▼]CCGGT...3'
3'...TGGCA...5'

Concentration: 20,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

AhdI



#R0584S	1,000 units\$75
#R0584L	5,000 units\$308

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	25	10	100

5'...GACNN[▼]NGTC...3'
3'...CTGNNNN[▼]CAG...5'

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by some combinations of overlapping CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

AleI-v2



#R0685S	500 units\$70
#R0685L	2,500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	<10	<10	100

5'...CACNN[▼]NGTG...3'
3'...GTGNN[▼]NCAC...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Impaired by overlapping CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

AluI



#R0137S	1,000 units\$70
#R0137L	5,000 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	100	50	100

5'...AG[▼]CT...3'
3'...TCGA...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

AlwI



#R0513S	500 units\$70
#R0513L	2,500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	50	10	100

5'...GGATC(N)₄[▼]...3'
3'...CCTAG(N)₅[▼]...5'

Methylation Sensitivity: Blocked by *dam* methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C

Concentration: 10,000 units/ml

Note: Star activity may result from extended digestion, high enzyme concentration or glycerol concentration of > 5%.

AlwNI



#R0514S	500 units\$70
#R0514L	2,500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	100	50	100

5'...CAGNN[▼]CTG...3'
3'...GTCNN[▼]NGAC...5'

Methylation Sensitivity: Blocked by overlapping *dcm* methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 10,000 units/ml

Note: Conditions of low ionic strength, high enzyme concentration, glycerol concentration > 5%, or pH > 8.0 may result in star activity.

ApaI

rCutSmart    25°  dcm CpG

#R0114S 5,000 units\$72
 #R0114L 25,000 units\$296

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 25 25 <10 100

5'... G G G C C C ...3'
 3'... C C C G G G ...5'

Activity at 37°C: 100% However, the half-life of ApaI at 37°C is only 30 minutes.

Methylation Sensitivity: Blocked by overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 25°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 50,000 units/ml

ApaLI

rCutSmart    37°  CpG

#R0507S 2,500 units\$67
 #R0507L 12,500 units\$274
 for high (5X) concentration
 #R0507M 12,500 units\$274

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 100 100 10 100



5'... G T G C A C ...3'
 3'... C A C G T G ...5'

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C

ApeKI

   NEB r3.1  75°  CpG

#R0643S 250 units\$72
 #R0643L 1,250 units\$296

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 25 50 100 10

5'... G C W G C ...3'
 3'... C G W C A G ...5'

Activity at 37°C: 10%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

Reaction Conditions: NEBuffer r3.1, 75°C

Concentration: 5,000 units/ml

ApoI

   NEB r3.1  50° 

#R0566S 1,000 units\$72

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 10 75 100 75

5'... R A A T T Y ...3'
 3'... Y T T A A R ...5'

Activity at 37°C: 50%

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: NEBuffer r3.1, 50°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 10,000 units/ml

ApoI-HF®

rCutSmart      37° 

#R3566S 1,000 units\$72
 #R3566L 5,000 units\$296

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity <10 100 10 100

5'... R A A T T Y ...3'
 3'... Y T T A A R ...5'

Concentration: 20,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

AseI

rCutSmart    37°  CpG

#R0558S 500 units\$70
 #R0558L 2,500 units\$287

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity <10 10 10 100

5'... G G C G C G C C ...3'
 3'... C C G C G C A G G ...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

AseI

   NEB r3.1  37° 

#R0526S 2,000 units\$67
 #R0526L 10,000 units\$274
 for high (5X) concentration
 #R0526M 10,000 units\$274

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity <10 50* 100 10

5'... A T T A A T ...3'
 3'... T A A T A ...5'

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Note: Star activity may result from a glycerol concentration of > 5%. *May exhibit star activity in this buffer.

AsiSI

rCutSmart     37°  CpG

#R0630S 500 units\$70
 #R0630L 2,500 units\$287

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 100 100 25 100

5'... G C G A T C G C ...3'
 3'... C G C A T A G C G ...5'

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Note: Star activity may result from extended digestion.

Concentration: 10,000 units/ml

AvaI

rCutSmart  RRI  dIIA 37°  CpG

#R0152S	2,000 units	\$67
#R0152L	10,000 units	\$274

for high (5X) concentration

#R0152T	2,000 units	\$67
#R0152M	10,000 units	\$274

5'...C[▼]Y C G R G...3'
3'...G R G C Y[▲]C...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	100	25	100

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

AvaII

rCutSmart  RRI  dIIA 37°  dcm CpG

#R0153S	2,000 units	\$67
#R0153L	10,000 units	\$274

for high (5X) concentration

#R0153M	10,000 units	\$274
---------	--------------	-------	-------

5'...G[▼]G W C C...3'
3'...C C W G[▲]G...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	75	10	100

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Blocked by overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

AvrII

rCutSmart  RRI  dIB 37° 

#R0174S	100 units	\$72
#R0174L	500 units	\$296

5'...C[▼]C T A G G...3'
3'...G G A T C[▲]C...5'

Reaction Conditions: rCutSmart Buffer, 37°C

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	50	50	100

Concentration: 5,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

BaeI

rCutSmart  RRI  SAM dIIA 25°  CpG

#R0613S	250 units	\$67
---------	-----------	-------	------

5'...¹⁰(N)₁AC(N)₁GTAYC(N)₁₂³'
3'...¹⁵(N)₁TG(N)₁CATRG(N)₂[▲]5'

Reaction Conditions: rCutSmart Buffer + SAM, 25°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	50	100

Activity at 37°C: 100%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

BaeGI

 RRI  r3.1 dIIA 37° 

#R0708S	500 units	\$67
---------	-----------	-------	------

5'...G K G C M[▼]C...3'
3'...C[▲]M C G K G...5'

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	75	100	25

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

BamHI

 RRI  r3.1 dIIA 37° 

#R0136S	10,000 units	\$62
#R0136L	50,000 units	\$254

for high (5X) concentration

#R0136T	10,000 units	\$62
#R0136M	50,000 units	\$254

5'...G[▼]G A T C C...3'
3'...C C T A G[▲]G...5'

Reaction Conditions: NEBuffer r3.1, 37°C

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75*	100*	100	100*

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Star activity may result from a glycerol concentration of > 5%. *May exhibit star activity in this buffer.

BamHI-HF[®]rCutSmart  RRI  e  dIIA 37° 

#R3136S	10,000 units	\$62
#R3136L	50,000 units	\$254

for high (5X) concentration

#R3136T	10,000 units	\$62
#R3136M	50,000 units	\$254

5'...G[▼]G A T C C...3'
3'...C C T A G[▲]G...5'

Reaction Conditions: rCutSmart Buffer, 37°C

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	50	10	100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

BanI

rCutSmart  RRI  dIIA 37°  dcm CpG

#R0118S	5,000 units	\$70
---------	-------------	-------	------

5'...G[▼]G Y R C C...3'
3'...C C R Y G[▲]G...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 20,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	25	<10	100

Methylation Sensitivity: Blocked by some combinations of overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

BanII

rCutSmart   37° 

#R0119S 2,000 units\$70

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	50	100

5'... G R G C Y^Y C ... 3'
3'... C Y C G R G ... 5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Star activity may result from extended digestion.

BbsI

  37° 

#R0539S 300 units\$70

#R0539L 1,500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	25	75

5'... G A A G A C (N)₂ ... 3'
3'... C T T C T G (N)₆ ... 5'

Reaction Conditions: NEBuffer r2.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Store at -80°C.

BbsI-HF[®]

High-Fidelity
rCutSmart   37° 

#R3539S 300 units\$70

#R3539L 1,500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	10	<10	100

5'... G A A G A C (N)₂ ... 3'
3'... C T T C T G (N)₆ ... 5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 20,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

BbvCI

rCutSmart   37°  CpG

#R0601S 100 units\$67

#R0601L 500 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	100	50	100

5'... C C^T T C A G C ... 3'
3'... G G A G T C G ... 5'

Reaction Conditions: rCutSmart Buffer, 37°C.

Concentration: 2,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by overlapping CpG methylation (see p. 321).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

BccI

rCutSmart   37° 

#R0704S 1,000 units\$70

#R0704L 5,000 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	50	10	100

5'... C C A T C (N)₄ ... 3'
3'... G G T A G (N)₆ ... 5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Star activity may result from a glycerol concentration of > 5%.

BceAI

  37°  CpG

#R0623S 50 units\$67

#R0623L 250 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100*	100*	100	100*

5'... A C G G C (N)₁₂ ... 3'
3'... T G C C G (N)₁₄ ... 5'

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 2,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%. *May exhibit star activity in this buffer.

BbvI

rCutSmart   37° 

#R0173S 300 units\$70

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	25	100

5'... G C A G C (N)₈ ... 3'
3'... C G T C G (N)₁₂ ... 5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 2,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Star activity may result from a glycerol concentration of > 5%.

BcgI

  37°  *dam* CpG

#R0545S 250 units\$70

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	75*	100	50*

5'...₁₀(N) C G A (N)₈ T G C (N)₁₂ ... 3'
3'...₁₂(N) G C T (N)₈ A C G (N)₁₀ ... 5'

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 2,000 units/ml

Methylation Sensitivity: Impaired by overlapping *dam* methylation. Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

Note: *May exhibit star activity in this buffer.

BciVI

#R0596S	200 units	\$67
#R0596L	1,000 units	\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	25	<10	100

5'...GTATCC(N)₆...3'
3'...CATAGG(N)₅...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

BclI

#R0160S	3,000 units	\$62
---------	-------------	-------	------

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	100	75

5'...TGATCA...3'
3'...ACTAGT...5'

Activity at 37°C: 50%

Methylation Sensitivity: Blocked by *dam* methylation (see p. 321).

Reaction Conditions: NEBuffer r3.1, 50°C

Concentration: 10,000 units/ml

BclI-HF®

#R3160S	3,000 units	\$62
#R3160L	15,000 units	\$254

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	10	100

5'...TGATCA...3'
3'...ACTAGT...5'

Concentration: 20,000 units/ml

Methylation Sensitivity: Blocked by *dam* methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

BcoDI

#R0542S	1,000 units	\$70
---------	-------------	-------	------

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	75	75	100

5'...GTCTC(N)₁...3'
3'...CAGAG(N)₅...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by some combinations of overlapping CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C

BfaI

#R0568S	500 units	\$75
#R0568L	2,500 units	\$308

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	10	<10	100

5'...CTAG...3'
3'...GATC...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Note: Star activity may result from extended digestion. Store at -80°C.

BfuAI

#R0701S	250 units	\$67
#R0701L	1,250 units	\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	25	100	10

5'...ACCTGC(N)₄...3'
3'...TGAGC(N)₈...5'

Activity at 37°C: 25%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by overlapping CpG methylation (see p. 321).

Reaction Conditions: NEBuffer r3.1, 50°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

Note: Star activity may result from a glycerol concentration of > 5%.

BglI

#R0143S	2,000 units	\$62
#R0143L	10,000 units	\$254

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	25	100	10

5'...GCCNNNNNGGC...3'
3'...CGGNNNNCCG...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

BglII

#R0144S	2,000 units	\$62
#R0144L	10,000 units	\$254

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	10	100	<10

for high (5X) concentration

#R0144M	10,000 units	\$254
---------	--------------	-------	-------

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

5'...AGATCT...3'
3'...TCTAGA...5'

Reaction Conditions: NEBuffer r3.1, 37°C.

BlnI



#R0585S 500 units\$75
#R0585L 2,500 units\$308

5'... G C T N A G C ...3'
3'... C G A N T C G ...5'


Reaction Conditions: rCutSmart Buffer, 37°C

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	10	100

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

BmgBI



#R0628S 500 units\$72
#R0628L 2,500 units\$296

5'... C A C G T C ...3'
3'... G T G C A G ...5'

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 65°C for 20 minutes.


Concentration: 10,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	10	100	10

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Note: Star activity may result from a glycerol concentration of > 5%.

BmrI



#R0600S 100 units\$70

5'... A C T G G G (N)₅ ...3'
3'... T G A C C C (N)₄ ...5'

Reaction Conditions: NEBuffer r2.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	100	75	100*

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: *May exhibit star activity in this buffer.

BmtI



#R0658S 300 units\$70

5'... G C T A G C ...3'
3'... C G A T C G ...5'

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	100	100*

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Star activity may result from extended digestion. *For added flexibility, NEB offers an isoschizomer or HF enzyme, supplied with rCutSmart Buffer.

BmtI-HF®



#R3658S 300 units\$70
#R3658L 1,500 units\$287

5'... G C T A G C ...3'
3'... C G A T C G ...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	10	100

Concentration: 20,000 units/ml.

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Store at -80°C.

BpmI



#R0565S 100 units\$70
#R0565L 500 units\$287

5'... C T G G A G (N)₁₆ ...3'
3'... G A C C T C (N)₁₄ ...5'

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 2,000 units/ml.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	100	100	100*

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Star activity may result from a extended digestion. *May exhibit star activity in this buffer.

Bpu10I



#R0649S 200 units\$70

5'... C C T N A G C ...3'
3'... G G A N T C G ...5'

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 5,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	25	100	25

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Star activity may result from a glycerol concentration of > 5%.

BpuEI



#R0633S 500 units\$70

5'... C T T G A G (N)₁₆ ...3'
3'... G A A C T C (N)₁₄ ...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50*	100	50*	100

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: *May exhibit star activity in this buffer.

BsaI

BsaI has been replaced by BsaI-HFv2 (engineered for performance).

BsaI-HF[®]v2



#R3733S 1,000 units\$70
#R3733L 5,000 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	100	100

5'... G G T C T C (N)₁... 3'
3'... C C A G A G (N)₅... 5'

Methylation Sensitivity: Impaired by some combinations of overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 20,000 units/ml

BsaAI



#R0531S 500 units\$70

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	100	100

5'... Y A C G T R ... 3'
3'... R T G C A Y ... 5'

Concentration: 5,000 units/ml

Reaction Conditions: rCutSmart Buffer, 37°C

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

BsaBI



#R0537S 2,000 units\$70

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	75	100

5'... G A T N N N N A T C ... 3'
3'... C T A N N N N T A G ... 5'

Reaction Conditions: rCutSmart Buffer, 60°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 10,000 units/ml

Activity at 37°C: 25%

Methylation Sensitivity: Blocked by overlapping *dam* methylation. Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

Note: Star activity may result from extended digestion.

BsaHI



#R0556S 2,000 units\$67

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	100	100

5'... G R C G Y C ... 3'
3'... C Y G C R G ... 5'

Concentration: 10,000 units/ml

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Methylation Sensitivity: Blocked by some combinations of overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

BsaJI



#R0536S 1,000 units\$70

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	100	100

5'... C C N N G G ... 3'
3'... G G N N C C ... 5'

Reaction Conditions: rCutSmart Buffer, 60°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 10,000 units/ml

Activity at 37°C: 25%

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

BsaWI



#R0567S 250 units\$75

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	100	50	100

5'... W C C G G W ... 3'
3'... W G G C C W ... 5'

Reaction Conditions: rCutSmart Buffer, 60°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 10,000 units/ml

Activity at 37°C: 50%

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

BsaXI



#R0609S 100 units\$70

#R0609L 500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50*	100*	10	100

5'...₉(N) A C (N)₁ C T C C (N)₁₀... 3'
3'...₁₂(N) T G (N)₈ G A G G (N)₇... 5'

Reaction Conditions: rCutSmart Buffer, 37°C

Concentration: 2,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: *May exhibit star activity in this buffer.

BseRI



#R0581S 200 units\$75

#R0581L 1,000 units\$308

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	75	100

5'... G A G G A G (N)₁₀... 3'
3'... C T C C T C (N)₈... 5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 5,000 units/ml.

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

BseYI



#R0635S	100 units\$70
#R0635L	500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	50	100	50

5'...C⁺CCAGC...3'
3'...GGGT⁺C_A...5'

Concentration: 5,000 units/ml.

Methylation Sensitivity:

Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

BsgI



#R0559S	50 units\$70
#R0559L	250 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	50	25	100

5'...GTGCAG(N)₁₆...3'
3'...CACGTC(N)₁₄...5'

Concentration: 5,000 units/ml

Methylation Sensitivity: Not

sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

BsiEI



#R0554S	1,000 units\$70
---------	-------------	-----------

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	50	<10	100

5'...CGRY⁺CG...3'
3'...GCY⁺RGC...5'

Activity at 37°C: 10%

Methylation Sensitivity: Cleavage

of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 60°C.

Concentration: 10,000 units/ml

BsiHKA1



#R0570S	1,000 units\$67
---------	-------------	-----------

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	100	100	100

5'...GWGCW⁺C...3'
3'...C_AWCGWG...5'

Activity at 37°C: 10%


Methylation Sensitivity: Not

sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 65°C.

Concentration: 10,000 units/ml.

BsiWI



#R0553S	300 units\$70
#R0553L	1,500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	50*	100	25

5'...C⁺G T A C G...3'
3'...G C A T G⁺C...5'

Activity at 37°C: 25%

Methylation Sensitivity: Cleavage

of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: NEBuffer r3.1, 55°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

Note: *May exhibit star activity in this buffer.

BsiWI-HF[®]

#R3553S	300 units\$70
#R3553L	1,500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	10	100

5'...C⁺G T A C G...3'
3'...G C A T G⁺C...5'

Concentration: 20,000 units/ml

Methylation Sensitivity: Cleavage

of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C.

BslI



#R0555S	1,000 units\$67
#R0555L	5,000 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	75	100	100

5'...CCNNNNNN⁺NNGG...3'
3'...GGNNNNNN⁺CC...5'

Activity at 37°C: 50%

Methylation Sensitivity: Blocked

by some combinations of overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 55°C

Concentration: 10,000 units/ml

BsmI



#R0134S	500 units\$70
#R0134L	2,500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	100	<10	100

5'...GAATGCN⁺...3'
3'...CTTAC_AGN...5'

Concentration: 10,000 units/ml

Activity at 37°C: 10%

Methylation Sensitivity: Not

sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 65°C. Heat inactivation: 80°C for 20 minutes.



Recommended Buffer



Cloned at NEB



Recombinant Enzyme



EpiMark Validated



Engineered for Performance



Multiple Recognition Site



Time-Saver Qualified



Diluent Buffer



Incubation Temperature



Heat Inactivation



Methylation Sensitivity

BsmAI

rCutSmart  RR             

BspHI

rCutSmart  RR             

BsrGI-HF[®]

#R3575S	1,000 units\$70
#R3575L	5,000 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	100	100	100

5'...T[↓]GTACA...3'
3'...ACATG[↓]T...5'

Concentration: 20,000 units/ml.

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

BssHII

#R0199S	500 units\$70
#R0199L	2,500 units\$287

for high (5X) concentration

#R0199M	2,500 units\$287
---------	-------------	------------

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	100	100

5'...G[↓]CGCGC...3'
3'...CGCGC[↓]G...5'

Concentration: 5,000 and 25,000 units/ml**Activity at 37°C:** 100%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 50°C. Heat inactivation: 65°C for 20 minutes.

BssSI-v2

#R0680S	200 units\$75
#R0680L	1,000 units\$308

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	25	<10	100

5'...C[↓]ACGAG...3'
3'...GTGCT[↓]C...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C

BstAPI

#R0654S	200 units\$95
#R0654L	1,000 units\$389

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	25	100

5'...GCANNN[↓]NTGC...3'
3'...CGTNNNNACG...5'

Activity at 37°C: 10%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 60°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 5,000 units/ml**BstBI**

#R0519S	2,500 units\$64
#R0519L	12,500 units\$262

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	100	10	100

5'...T[↓]TCGAA...3'
3'...AAGCT[↓]T...5'

Activity at 37°C: 10%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 65°C

Concentration: 20,000 units/ml**BstEII**

#R0162S	2,000 units\$64
#R0162L	10,000 units\$262

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	75*	100	75*

5'...G[↓]GTNACC...3'
3'...CCANTG[↓]G...5'

Activity at 37°C: 10%

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: NEBuffer r3.1, 60°C

Concentration: 10,000 units/ml

Note: Star activity may result from a glycerol concentration of > 5%. *May exhibit star activity in this buffer.

BstEII-HF[®]

#R3162S	2,000 units\$64
#R3162L	10,000 units\$262

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	10	<10	100

for high (5X) concentration

#R3162M	10,000 units\$262
---------	--------------	------------

Concentration: 20,000 and 100,000 units/ml.

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C

BstNI

#R0168S	3,000 units\$62
#R0168L	15,000 units\$254

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	100	100	75

5'...CC[↓]WGG...3'
3'...GGW[↓]CC...5'

Activity at 37°C: 25%

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: NEBuffer r3.1, 60°C

Concentration: 10,000 units/ml

BstUI

rCutSmart   60°  

#R0518S 1,000 units\$64
#R0518L 5,000 units\$262

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	25	100

5'...CGTCG...3'
3'...GCGC...5'

Activity at 37°C: 25%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 60°C

Concentration: 10,000 units/ml

BstXI

    37°  

#R0113S 1,000 units\$70
#R0113L 5,000 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	50	100	25

5'...CCANNNNNTGG...3'
3'...GGTNNNNNACC...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Blocked by some combinations of overlapping *dcm* methylation (see p. 321).

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

Note: Star activity may result from a glycerol concentration of > 5%.

BstYI

    60°  

#R0523S 2,000 units\$67

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	100	75	100

5'...RGATCY...3'
3'...YCTAGR...5'

Activity at 37°C: 10%

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: NEBuffer r2.1, 60°C

Concentration: 10,000 units/ml

BstZ17I-HF®

rCutSmart      37°  

#R3594S 1,000 units\$70
#R3594L 5,000 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	10	100

5'...GTATAC...3'
3'...CATATG...5'

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C.

Concentration: 20,000 units/ml

Bsu36I

rCutSmart     37°  

#R0524S 1,000 units\$67
#R0524L 5,000 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	100	100	100

5'...CCTNAGG...3'
3'...GGANTCC...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

BtgI

rCutSmart     37°  

#R0608S 1,000 units\$67

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	100	100

5'...CRYGG...3'
3'...GGYRC...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

BtgZI

rCutSmart    60°  

#R0703S 100 units\$70
#R0703L 500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	25	<10	100

5'...GCGATG(N)₁₀...3'
3'...CGCTAC(N)₁₄...5'

Activity at 37°C: 50%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 60°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 5,000 units/ml

Note: Star activity may result from a glycerol concentration of > 5%.

BtsI-v2

rCutSmart      37°  

#R0667S 500 units\$70
#R0667L 2,500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	25	100

5'...GCAGTGNN...3'
3'...CGTCAcNN...5'

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C.

Concentration: 10,000 units/ml

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

BtsIMutIrCutSmart      

#R0664S 100 units\$107

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	50	10	100

5'...CAGTGN[▼]...3'
3'...GTCAC[▲]NN...5'

Concentration: 1,000 units/ml**Activity at 37°C:** 50%**Methylation Sensitivity:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.**Reaction Conditions:** rCutSmart Buffer, 55°C. Heat inactivation: 80°C for 20 minutes.**BtsCI**rCutSmart      

#R0647S 2,000 units\$72

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	100	25	100

5'...GGATGN[▼]...3'
3'...CCTAC[▲]NN...5'

Concentration: 20,000 units/ml**Activity at 37°C:** 25%**Methylation Sensitivity:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.**Reaction Conditions:** rCutSmart Buffer, 50°C. Heat inactivation: 80°C for 20 minutes.**Cac8I**rCutSmart      

#R0579L 500 units\$308

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	75	100	100

5'...GCN[▼]NGC...3'
3'...CGN[▲]CGC...5'

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).**Note:** Store at -80°C.**Reaction Conditions:** rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.**Concentration:** 5,000 units/ml**ClaI**rCutSmart        

#R0197S 1,000 units\$67

#R0197L 5,000 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	50	50	100

5'...AT[▼]CGAT...3'
3'...TAG[▲]CTA...5'

Concentration: 10,000 units/ml**Methylation Sensitivity:** Blocked by overlapping *dam* methylation. Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).**Reaction Conditions:** rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.**CspCI**rCutSmart       

#R0645S 500 units\$69

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	100	10	100

5'...¹⁰⁻¹¹(N)CAA(N)₂GTTG(N)₁₀₋₁₃...3'
3'...¹⁰⁻¹³(N)GTT(N)₂CACC(N)₁₀₋₁₁...5'

Note: The cleavage point may shift one base pair depending on the DNA sequence context before and after the recognition site. For a given sequence, one site will predominate. For details, see www.neb.com.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.**Concentration:** 5,000 units/ml**Methylation Sensitivity:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.**CviAI**rCutSmart      

#R0640S 200 units\$67

#R0640L 1,000 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	50	10	100

5'...C[▼]ATG...3'
3'...GTA[▲]C...5'

Concentration: 10,000 units/ml**Activity at 37°C:** 10%**Methylation Sensitivity:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.**Reaction Conditions:** rCutSmart Buffer, 25°C. Heat inactivation: 65°C for 20 minutes.**CviKI-1**rCutSmart      

#R0710S 250 units\$72

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	100	100	100

5'...RG[▼]CY...3'
3'...YCG[▲]R...5'

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.**Reaction Conditions:** rCutSmart Buffer, 37°C**Concentration:** 5,000 units/ml**Note:** Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.**CviQI**      

#R0639S 2,000 units\$72

#R0639L 10,000 units\$296

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	100*	100	75*

5'...G[▼]TAC...3'
3'...CAT[▲]G...5'

Activity at 37°C: 25%**Methylation Sensitivity:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.**Reaction Conditions:** NEBuffer r3.1, 25°C**Concentration:** 10,000 units/ml**Note:** *May exhibit star activity in this buffer.

DdeI

rCutSmart     37°  

#R0175S	1,000 units\$70
#R0175L	5,000 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	100	100	100

Concentration: 10,000 units/ml**Methylation Sensitivity:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

5'...CTNAG...3'
3'...GANTC...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

DpnI

rCutSmart     37°  

#R0176S	1,000 units\$67
#R0176L	5,000 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	75	100

Concentration: 20,000 units/ml**Methylation Sensitivity:** DpnI cleaves only when its recognition site is methylated. DNA purified from a *dam*⁺ strain will be a substrate for DpnI. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

CH₃
5'...GATC...3'
3'...CTAG...5'
CH₃

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

DpnII

    37°  

#R0543S	1,000 units\$72
#R0543L	5,000 units\$296

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	25	100*	25

for high (5X) concentration

#R0543T	1,000 units\$72
#R0543M	5,000 units\$296

Concentration: 10,000 and 50,000 units/ml**Methylation Sensitivity:** Blocked by *dam* methylation (see p. 321).**Note:** *May exhibit star activity in this buffer.

5'...GATC...3'
3'...CTAG...5'

Reaction Conditions: NEBuffer DpnII, 37°C. Heat inactivation: 65°C for 20 minutes.

DraI

rCutSmart     37°  

#R0129S	2,000 units\$67
#R0129L	10,000 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	75	50	100

Concentration: 20,000 units/ml.**Methylation Sensitivity:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

5'...TTTAAA...3'
3'...AAATT...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.DraIII-HF[®]rCutSmart     37°   

#R3510S	1,000 units\$70
#R3510L	5,000 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	50	10	100

Concentration: 20,000 units/ml**Methylation Sensitivity:** Cleavage of mammalian genomic DNA is impaired by some combinations of overlapping CpG methylation (see p. 321).

5'...CACNNNGTG...3'
3'...GTGNNNCAC...5'

Reaction Conditions: rCutSmart Buffer, 37°C

DrdI

rCutSmart     37°  

#R0530S	300 units\$70
#R0530L	1,500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	50	10	100

5'...GACNNNNNGTC...3'
3'...CTGNNNNNCAG...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.**Concentration:** 5,000 units/ml**Methylation Sensitivity:** Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).**Note:** Star activity may result from a glycerol concentration of > 5%.

EaeI

rCutSmart     37°   

#R0508S	200 units\$67
#R0508L	1,000 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	50	<10	100

5'...YGGCCR...3'
3'...RCGGY...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.**Concentration:** 5,000 units/ml**Methylation Sensitivity:** Blocked by overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

EagI

EagI has been replaced by EagI-HF (engineered for performance).

High-Fidelity

EagI-HF®

rCutSmart  RR                       

EcoRI-HF®rCutSmart     High-Fidelity

#R3101S	10,000 units\$62
#R3101L	50,000 units\$254

for high (5X) concentration

#R3101T	10,000 units\$62
#R3101M	50,000 units\$254

5'...GAATTC...3'
3'...CTTAAG...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	100	<10	100

Concentration: 20,000 and 100,00 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

EcoRV    NEB r3.1 dIIA

#R0195S	4,000 units\$62
#R0195L	20,000 units\$254

for high (5X) concentration

#R0195T	4,000 units\$62
#R0195M	20,000 units\$254

5'...GATATC...3'
3'...CTAATG...5'

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	50	100	10

Concentration: 20,000 and 100,000 units/ml.

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by some combinations of overlapping CpG methylation (see p. 321).

EcoRV-HF®rCutSmart     High-Fidelity

#R3195S	4,000 units\$62
#R3195L	20,000 units\$254

for high (5X) concentration

#R3195T	4,000 units\$62
#R3195M	20,000 units\$254

5'...GATATC...3'
3'...CTAATG...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	100	100	100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by some combinations of overlapping CpG methylation (see p. 321).

Esp3IrCutSmart     

#R0734S	300 units\$77
#R0734L	1,500 units\$350

5'...CGTCTC(N)₁...3'
3'...GCAGAG(N)₅...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	<10	100

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

FatI     NEB r2.1 dIIA

#R0650S	50 units\$95
#R0650L	250 units\$389

5'...CATG...3'
3'...GTAC...5'

Reaction Conditions: NEBuffer r2.1, 55°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 2,000 units/ml

Activity at 37°C: 100%

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	100	50	50

FauIrCutSmart     

#R0651S	200 units\$72
---------	-----------	-----------

5'...CCCGC(N)₄...3'
3'...GGGCG(N)₆...5'

Reaction Conditions: rCutSmart Buffer, 55°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

Activity at 37°C: 20%

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Note: Star activity may result from a glycerol concentration of > 5%.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	50	10	100

Fnu4HIrCutSmart    

#R0178S	200 units\$70
#R0178L	1,000 units\$287

5'...G^NCGC...3'
3'...CG^NCG...5'

Reaction Conditions: rCutSmart Buffer, 37°C

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	<10	<10	100

FokIrCutSmart      2'-site

#R0109S	1,000 units\$70
#R0109L	5,000 units\$287

5'...GGATG(N)₉...3'
3'...CTAC(N)₁₃...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

Methylation Sensitivity: Impaired by overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is impaired by overlapping CpG methylation (see p. 321).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	75	100

FseI

rCutSmart  RR    B 37°   dcm CpG#R0588S 100 units\$75
#R0588L 500 units\$308NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 100 75 <10 1005'...GGCCGGC...3'
3'...CCGGCCGG...5'**Methylation Sensitivity:** Impaired by some combinations of overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).**Reaction Conditions:** rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.**Note:** Store at -80°C.**Concentration:** 2,000 units/ml

FspI

rCutSmart  RR    C 37°   CpG#R0135S 500 units\$70
#R0135L 2,500 units\$287NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 10 100 10 1005'...TGCGCA...3'
3'...ACGCGT...5'**Concentration:** 10,000 units/ml.**Methylation Sensitivity:** Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).**Reaction Conditions:** rCutSmart Buffer, 37°C

FspEI

rCutSmart  RR    B 37°  

#R0662S 200 units\$112

See page 280 for more information.

HaeII

rCutSmart  RR    A 37°   CpG#R0107S 2,000 units\$64
#R0107L 10,000 units\$262NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 25 100 10 1005'...RGCGCY...3'
3'...YCGCGR...5'**Concentration:** 20,000 units/ml**Methylation Sensitivity:** Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).**Reaction Conditions:** rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

HaeIII

rCutSmart  RR    A 37°  #R0108S 3,000 units\$65
#R0108L 15,000 units\$266NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 50 100 25 100

for high (5X) concentration

#R0108T 3,000 units\$65
#R0108M 15,000 units\$266**Concentration:** 10,000 and 50,000 units/ml**Methylation Sensitivity:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.5'...GGCC...3'
3'...CCGG...5'**Reaction Conditions:** rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

HgaI

 RR    A 37°  CpG




#R0154S 100 units\$70

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 100 100 25 100*5'...GACGC(N)₅...3'
3'...CTGCG(N)₁₀...5'**Reaction Conditions:** NEBuffer r1.1, 37°C. Heat inactivation: 65°C for 20 minutes.**Concentration:** 2,000 units/ml**Methylation Sensitivity:** Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).**Note:** Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%. *May exhibit star activity in this buffer.






HhaI

rCutSmart  RR    A 37°  CpG#R0139S 2,000 units\$62
#R0139L 10,000 units\$254NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 25 100 100 1005'...GCGC...3'
3'...CGCG...5'**Concentration:** 20,000 units/ml**Reaction Conditions:** rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.**Methylation Sensitivity:** Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

HincII

 RR     B 37°  CpG#R0103S 1,000 units\$64
#R0103L 5,000 units\$262NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 25 100 100 1005'...GTY^YRAC...3'
3'...CARYT^AGTG...5'**Concentration:** 10,000 units/ml**Reaction Conditions:** NEBuffer r3.1, 37°C. Heat inactivation: 65°C for 20 minutes.**Methylation Sensitivity:** Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

HindIII

 RR     B 37° #R0104S 10,000 units\$62
#R0104L 50,000 units\$254NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 25 100 50 50

for high (5X) concentration

#R0104T 10,000 units\$62
#R0104M 50,000 units\$254**Concentration:** 20,000 and 100,000 units/ml5'...A^AAGCTT...3'
3'...TTCGA^A...5'**Methylation Sensitivity:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.**Reaction Conditions:** NEBuffer r2.1, 37°C. Heat inactivation: 80°C for 20 minutes.**Note:** Star activity may result from extended digestion.

HindIII-HF®rCutSmart       High-Fidelity

#R3104S	10,000 units	\$62
#R3104L	50,000 units	\$254

for high (5X) concentration

#R3104T	10,000 units	\$62
#R3104M	50,000 units	\$254

5'...A[▼]AGCTT...3'
3'...TTCGA[▲]A...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	100	10	100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

HinfIrCutSmart      

#R0155S	5,000 units	\$64
#R0155L	25,000 units	\$262

for high (5X) concentration

#R0155T	5,000 units	\$64
#R0155M	25,000 units	\$262

5'...G[▼]ANTC...3'
3'...CTNAG[▲]...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	100	100

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

HinP1IrCutSmart      

#R0124S	2,000 units	\$70
---------	-------------	-------	------

5'...G[▼]CGC...3'
3'...CGC[▲]G...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	100	100

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

HpaIrCutSmart      

#R0105S	500 units	\$65
#R0105L	2,500 units	\$266

5'...GTT[▼]AAC...3'
3'...CAA[▲]TG...5'

Reaction Conditions: rCutSmart Buffer, 37°C

Concentration: 5,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	75*	25	100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%. *May exhibit star activity in this buffer.

HpaIIrCutSmart      

#R0171S	2,000 units	\$67
#R0171L	10,000 units	\$274

for high (5X) concentration

#R0171M	10,000 units	\$274
---------	--------------	-------	-------

5'...C[▼]CGG...3'
3'...GGC[▲]...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	50	<10	100

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

HphIrCutSmart       

#R0158S	1,000 units	\$67
#R0158L	5,000 units	\$274

5'...GGTGA(N)[▼]...3'
3'...CCACT(N)[▲]...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	50	<10	100

Methylation Sensitivity: This enzyme is blocked by *dam* and *dcm* methylation (see p. 321).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

Hpy99IrCutSmart      

#R0615S	100 units	\$70
#R0615L	500 units	\$287

5'...CGWCG[▼]...3'
3'...G[▲]CGWGC...5'

Reaction Conditions: rCutSmart Buffer Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	10	<10	100

Concentration: 2,000 units/ml.

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Hpy166IIrCutSmart      

#R0616S	1,000 units	\$72
---------	-------------	-------	------

5'...GTN[▼]NAC...3'
3'...CAN[▲]TG...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	50	100

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

Hpy188I

rCutSmart  RRII dIIA 37° 65° dam

#R0617S 1,000 units\$70
 #R0617L 5,000 units\$287

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 25 100 50 100

5'...TCN^NGA...3'
 3'...AGN^NCT...5'

Methylation Sensitivity: Blocked by overlapping *dam* methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

Hpy188III

rCutSmart  RRII dII B 37° 65° dam CpG

#R0622S 500 units\$70
 #R0622L 2,500 units\$287

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 100 100 10 100

5'...T^NCNNGA...3'
 3'...AGNN^NCT...5'

Methylation Sensitivity: Blocked by overlapping *dam* methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

Note: Star activity may result from a glycerol concentration of > 5%.

HpyAV

rCutSmart  RRII dII B 37° 65° CpG

#R0621S 100 units\$72
 #R0621L 500 units\$296

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 100 100 25 100

5'...CCTTC(N)₆...3'
 3'...GGAAG(N)₆...5'

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by overlapping CpG methylation. (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 2,000 units/ml

Note: Star activity may result from a glycerol concentration of > 5%.

HpyCH4III

rCutSmart  RRII dIIA 37° 65°

#R0618S 250 units\$70
 #R0618L 1,250 units\$287

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 100 25 <10 100

5'...ACN^NGT...3'
 3'...T^NGNCA...5'

Concentration: 5,000 units/ml

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

HpyCH4IV

rCutSmart  RRII dIIA 37° 65° CpG

#R0619S 500 units\$70
 #R0619L 2,500 units\$287

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 100 50 25 100

5'...A^NCGT...3'
 3'...TGCA...5'

Concentration: 10,000 units/ml

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

HpyCH4V

rCutSmart  RRII dIIA 37° 65°

#R0620S 100 units\$70
 #R0620L 500 units\$287

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 50 50 25 100

5'...T^NGCA...3'
 3'...ACGT...5'

Concentration: 5,000 units/ml

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

KasI

rCutSmart  RRII dII B 37° 65° CpG

#R0544S 250 units\$70
 #R0544L 1,250 units\$287

NEBuffer r1.1 r2.1 r3.1 rCutSmart
 % Activity 50 100 50 100

5'...G^NGCGCC...3'
 3'...CCGCG^NG...5'

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

Note: Star activity may result from a glycerol concentration of > 5%. Store at -80°C.

KpnI

RRII  NEB r1.1 dIIA 37° 65°

#R0142S 4,000 units\$67
 #R0142L 20,000 units\$274
 for high (5X) concentration
 #R0142M 20,000 units\$274

NEBuffer **r1.1** r2.1 r3.1 rCutSmart
 % Activity **100** 75 <10 50

5'...GGTACC...3'
 3'...CATG...5'

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: NEBuffer r1.1, 37°C

Concentration: 10,000 and 50,000 units/ml

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

KpnI-HF®rCutSmart       High-Fidelity

#R3142S	4,000 units\$67
#R3142L	20,000 units\$274
for high (5X) concentration		
#R3142M	20,000 units\$274

5'...GGTAC[▼]C...3'
3'...[▲]CATGG...5'

Reaction Conditions: rCutSmart Buffer, 37°C

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	25	<10	100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

LpnPIrCutSmart     

#R0663S	200 units\$112
---------	-----------	------------

See page 280 for more information.

MboIrCutSmart      

#R0147S	500 units\$75
#R0147L	2,500 units\$308
for high (5X) concentration		
#R0147M	2,500 units\$308

5'...[▼]GATC...3'
3'...CTAG[▲]...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	100	100	100

Concentration: 5,000 and 25,000 units/ml

Methylation Sensitivity: Blocked by *dam* methylation. Its isoschizomer Sau3AI is not. Cleavage of mammalian genomic DNA is impaired by overlapping CpG methylation (see p. 321).

MboIIrCutSmart      

#R0148S	300 units\$70
#R0148L	1,500 units\$287

5'...GAAGA(N)₈[▼]...3'
3'...CTTCT(N)₇[▲]...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100*	100	50	100

Concentration: 5,000 units/ml

Methylation Sensitivity: Blocked by overlapping *dam* methylation (see p. 321).

Note: *May exhibit star activity in this buffer.

MfeIrCutSmart    

#R0589S	500 units\$75
#R0589L	2,500 units\$308

5'...CAATTG...3'
3'...GTTAA[▲]C...5'

Reaction Conditions: rCutSmart Buffer, 37°C

Concentration: 10,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	50	10	100

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Star activity may result from extended digestion.

MfeI-HF®rCutSmart       High-Fidelity

#R3589S	500 units\$75
#R3589L	2,500 units\$308








5'...C[▼]AATTG...3'
3'...GTTAA[▲]C...5'

Reaction Conditions: rCutSmart Buffer, 37°C

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	25	<10	100

Concentration: 20,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

MluI      

#R0198S	1,000 units\$67
#R0198L	5,000 units\$274

5'...A[▼]CGCGT...3'
3'...TGC[▲]GCA...5'

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	50	100	25

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

MluI-HF®rCutSmart        High-Fidelity

#R3198S	1,000 units\$67
#R3198L	5,000 units\$274

5'...A[▼]CGCGT...3'
3'...TGC[▲]GCA...5'

Reaction Conditions: rCutSmart Buffer, 37°C

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	100	100	100

Concentration: 20,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

MluCIrCutSmart     

#R0538S	1,000 units\$67
#R0538L	5,000 units\$274

5'...AATT...3'
3'...TTAA[▲]...5'

Reaction Conditions: rCutSmart Buffer, 37°C

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	10	10	100

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

MlyI



#R0610S	1,000 units\$67
#R0610L	5,000 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	50	10	100

5'...GAGTC(N)₅...3'
3'...CTCAG(N)₅...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

MmeI



#R0637S	100 units\$70
#R0637L	500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	50	100

5'...TCCRAC(N)₂₀...3'
3'...AGGYTG(N)₁₈...5'

Concentration: 2,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

MnlI



#R0163S	500 units\$70
#R0163L	2,500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	100	50	100

5'...CCTC(N)₇...3'
3'...GGAG(N)₆...5'

Concentration: 5,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

MscI



#R0534S	250 units\$72
#R0534L	1,250 units\$296

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	100	100	100

for high (5X) concentration

#R0534M	1,250 units\$296
---------	-------------	------------

Concentration: 5,000 and 25,000 units/ml

Methylation Sensitivity: Blocked by overlapping *dcm* methylation (see p. 321). The single MscI site in pBR322 overlaps a *dcm* methylation site; consequently, pBR322 which has been grown in a *dcm*⁻ host should be used for cloning.

5'...TGG[▼]CCA...3'
3'...ACCGGT...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

MseI



#R0525S	500 units\$70
#R0525L	2,500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	100	75	100

for high (5X) concentration

#R0525M	2,500 units\$287
---------	-------------	------------

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

5'...T[▼]TAA...3'
3'...AAT[▲]T...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

MslI



#R0571S	500 units\$72
#R0571L	2,500 units\$296

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	50	<10	100

5'...CAYNN[▼]NNRTG...3'
3'...GTRNN[▲]NYAC...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

MspI



#R0106S	5,000 units\$67
#R0106L	25,000 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	100	50	100

for high (5X) concentration

#R0106T	5,000 units\$67
---------	-------------	-----------

#R0106M	25,000 units\$274
---------	--------------	------------

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

5'...C[▼]CGG...3'
3'...GGC[▲]C...5'

Reaction Conditions: rCutSmart Buffer, 37°C

MspAII



#R0577S	500 units\$72
#R0577L	2,500 units\$296

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	50	10	100

5'...CMG[▼]CKG...3'
3'...GKCGMC...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

MspJI

rCutSmart    37° #R0661S 200 units\$112
#R0661L 1,000 units\$459

See page 280 for more information.

MwoI

rCutSmart    60°  CpG#R0573S 500 units\$72
#R0573L 2,500 units\$296NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity <10 100 100 1005'...GCNNNNNNNGC...3'
3'...CGNNNNNNNCG...5'**Activity at 37°C:** 25%**Methylation Sensitivity:** Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).**Reaction Conditions:** rCutSmart Buffer, 60°C**Concentration:** 5,000 units/ml

NaeI

rCutSmart    37°  CpG#R0190S 500 units\$64
#R0190L 2,500 units\$262NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 25 25 <10 1005'...GCCGGC...3'
3'...CGGCCG...5'**Concentration:** 10,000 units/ml**Methylation Sensitivity:** Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).**Reaction Conditions:** rCutSmart Buffer, 37°C

NarI

rCutSmart    37°  CpG#R0191S 500 units\$70
#R0191L 2,500 units\$287NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 100 100 10 1005'...GGCGCC...3'
3'...CCGCGG...5'**Concentration:** 5,000 units/ml**Methylation Sensitivity:** Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).**Reaction Conditions:** rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NciI

rCutSmart    37°  CpG#R0196S 2,000 units\$67
#R0196L 10,000 units\$274NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 100 25 10 1005'...CCSGG...3'
3'...GGSCC...5'**Concentration:** 20,000 units/ml**Methylation Sensitivity:** Cleavage of mammalian genomic DNA is impaired by overlapping CpG methylation (see p. 321).**Reaction Conditions:** rCutSmart Buffer, 37°C

NcoI

rCutSmart    37° #R0193S 1,000 units\$65
#R0193L 5,000 units\$266NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 100 100 100 100*

for high (5X) concentration

#R0193T 1,000 units\$65
#R0193M 5,000 units\$266**Concentration:** 10,000 and 50,000 units/ml5'...CCATGG...3'
3'...GGTACC...5'**Methylation Sensitivity:**Not sensitive to *dam*, *dcm* or mammalian CpG methylation. *For added flexibility, NEB offers an isoschizomer or HF enzyme, supplied with rCutSmart Buffer.**Reaction Conditions:**

NEBuffer r3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

NcoI-HF®

rCutSmart    37°  High-Fidelity#R3193S 1,000 units\$65
#R3193L 5,000 units\$266NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 50 100 10 100

for high (5X) concentration

#R3193M 5,000 units\$266

Concentration: 20,000 and 100,000 units/ml5'...CCATGG...3'
3'...GGTACC...5'**Methylation Sensitivity:** Notsensitive to *dam*, *dcm* or mammalian CpG methylation.**Reaction Conditions:** rCutSmart

Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NdeI

rCutSmart    37° #R0111S 4,000 units\$67
#R0111L 20,000 units\$274NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 75 100 100 1005'...CATATG...3'
3'...GTATAC...5'**Concentration:** 20,000 units/ml**Methylation Sensitivity:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.**Reaction Conditions:** rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NgoMIV

rCutSmart     37°  CpG#R0564S 1,000 units\$67
#R0564L 5,000 units\$274NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 100 50 10 1005'...GCCGGC...3'
3'...CGGCCG...5'**Methylation Sensitivity:** Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 334).**Reaction Conditions:** rCutSmart Buffer, 37°C**Concentration:** 10,000 units/ml**Note:** Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

NheI

NheI has been replaced by NheI-HF (engineered for performance).

NmeAIII

rCutSmart  RR 2' site  **dnB** 37° 

#R0711S 250 units\$72

NEBuffer r1.1 r2.1 r3.1 **rCutSmart**
% Activity 10 10 <10 **100**

5'...GCCGAG(N)₂₀₋₂₁...3'
3'...CGGCTC(N)₁₈₋₁₉...5'

Note: The cleavage point may shift one base pair depending on the DNA sequence context between the recognition site and the position of cleavage. For a given sequence, generally one site will predominate. For details, see www.neb.com.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 2,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

NheI-HF®

High-Fidelity
rCutSmart  RR  **e**  **dnC** 37°  **180°**  CpG

#R3131S 1,000 units\$70

#R3131L 5,000 units\$287

for high (5X) concentration

#R3131M 5,000 units\$287

NEBuffer r1.1 r2.1 r3.1 **rCutSmart**
% Activity 100 25 <10 **100**

Concentration: 20,000 and 100,000 units/ml

5'...GCTAGC...3'
3'...CGATCG...5'

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NotI

 RR  **dnB** 37°  CpG

#R0189S 500 units\$75

#R0189L 2,500 units\$308

for high (5X) concentration

#R0189M 2,500 units\$308

NEBuffer r1.1 r2.1 **r3.1** **rCutSmart**
% Activity <10 50 **100** 25

Concentration: 10,000 and 50,000 units/ml

5'...GCGGCCGC...3'
3'...CGCCGGCG...5'

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

NlaIII

rCutSmart  RR  **dnB** 37° 

#R0125S 500 units\$70

#R0125L 2,500 units\$287

NEBuffer r1.1 r2.1 r3.1 **rCutSmart**
% Activity <10 <10 <10 **100**

Concentration: 10,000 units/ml

5'...CATG...3'
3'...GTAC...5'

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Note: Store at -80°C.

NotI-HF®

High-Fidelity
rCutSmart  RR  **e**  **dnA** 37°  **180°**  CpG

#R3189S 500 units\$75

#R3189L 2,500 units\$308

for high (5X) concentration

#R3189M 2,500 units\$308

NEBuffer r1.1 r2.1 r3.1 **rCutSmart**
% Activity 25 100 25 **100**

Concentration: 20,000 and 100,000 units/ml

5'...GCGGCCGC...3'
3'...CGCCGGCG...5'

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NlaIV

rCutSmart  RR  **dnB** 37°  **180°**  **dcm**  CpG

#R0126S 200 units\$70

NEBuffer r1.1 r2.1 r3.1 **rCutSmart**
% Activity 10 10 10 **100**

Concentration: 2,000 units/ml

5'...GGNNCC...3'
3'...CCNNGG...5'

Methylation Sensitivity: Blocked by overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NruI

 RR  **dnB** 37°  **180°**  **dam**  CpG

#R0192S 1,000 units\$62

#R0192L 5,000 units\$254

NEBuffer r1.1 r2.1 **r3.1** **rCutSmart**
% Activity <10 10 **100** 10

Concentration: 10,000 units/ml

5'...TCGCGA...3'
3'...AGCGCT...5'

Methylation Sensitivity: Blocked by overlapping *dam* methylation. Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: NEBuffer r3.1, 37°C

NruI-HF[®]
High-Fidelity

#R3192S	1,000 units\$62
#R3192L	5,000 units\$254

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	0	25	50	100

5'...T C G C G A...3'
3'...A G C G C T...5'

Concentration: 20,000 units/ml

Methylation Sensitivity: Blocked by overlapping *dam* methylation. Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C.

NsiI

#R0127S	1,000 units\$67
#R0127L	5,000 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	75	100	25

5'...A T G C A T...3'
3'...T A C G T A...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

NsiI-HF[®]
High-Fidelity

#R3127S	1,000 units\$67
#R3127L	5,000 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	20	<10	100

5'...A T G C A T...3'
3'...T A C G T A...5'

Concentration: 20,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NspI

#R0602S	250 units\$70
#R0602L	1,250 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	<10	100

5'...R C A T G Y...3'
3'...Y G T A C R...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

PacI

#R0547S	250 units\$71
#R0547L	1,250 units\$291

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	75	10	100

5'...T T A A T T A A...3'
3'...A A T T A A T T...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

PaeR7I

#R0177S	2,000 units\$62
---------	-------------	-----------

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	100	10	100

5'...C T C G A G...3'
3'...G A G C T C...5'

Concentration: 20,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C

NEW
PaqCI[™]

#R0745S	200 units\$100
#R0745L	1,000 units\$410

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	100	10	100

5'...C A C C T G C (N)₄...3'
3'...G T G G A C G (N)₄...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%. Use with PaqCI Activator.

PciI

#R0655S	200 units\$72
#R0655L	1,000 units\$296

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	75	100	50*

5'...A C A T G T...3'
3'...T G T A C A...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

Note: *May exhibit star activity in this buffer.

PflFIrCutSmart    37° 

#R0595S 2,000 units\$67

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	100	25	100

5'...GACNNNGTC...3'
3'...CTGNNNCAG...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

PflMI   37°  

#R0509S 1,000 units\$67

#R0509L 5,000 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	0	100	100	50

5'...CCANNNNNTGG...3'
3'...GGTNNNNNACC...5'

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

Methylation Sensitivity: Blocked by overlapping *dcm* methylation (see p. 321).

Note: Star activity may result from a glycerol concentration of > 5%.

PleIrCutSmart    37°  

#R0515S 1,000 units\$67

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	50	25	100

5'...GAGTC(N)_A...3'
3'...CTCAG(N)_A...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

PluTIrCutSmart    37°  

#R0713S 500 units\$69

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	25	<10	100

5'...GGCGC_T...3'
3'...CCGCGG...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

PmeIrCutSmart    37°  

#R0560S 500 units\$75

#R0560L 2,500 units\$308

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	50	10	100

5'...GTTTAAAC...3'
3'...CAAAATTTG...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

PmlIrCutSmart    37°  

#R0532S 2,000 units\$70

#R0532L 10,000 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	50	<10	100

5'...CACGTG...3'
3'...GTGCAC...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 20,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

PpuMIrCutSmart    37°  

#R0506S 500 units\$70

#R0506L 2,500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	<10	<10	100

5'...RG_TGWCCY...3'
3'...YCCW_TGR...5'

Reaction Conditions: rCutSmart Buffer, 37°C

Concentration: 10,000 units/ml

Methylation Sensitivity: Blocked by overlapping *dcm* methylation (see p. 321).

PshAIrCutSmart    37°  

#R0593S 1,000 units\$70

#R0593L 5,000 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	50	10	100

5'...GACNNNGTC...3'
3'...CTGNNNNCAG...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).

PsiI

PsiI has been replaced by PsiI-v2 (optimized for performance).

NEW

PsiI-v2

rCutSmart RR e dclB 37° 165°

#R0744S 400 units\$119
#R0744L 2,000 units\$482

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	50	10	100

5'...TTATAA...3'
3'...AATATT...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Note: Star activity may result from a glycerol concentration of > 5%.

PspGI

rCutSmart RR dclA 75° 165° dcm

#R0611S 1,000 units\$67

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	100	50	100

5'...CCWGG...3'
3'...GGWCC...5'

Activity at 37°C: 25%

Methylation Sensitivity: Blocked by *dcm* methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 75°C

Concentration: 10,000 units/ml

Note: Star activity may result from a glycerol concentration of > 5%.

PspOMI

rCutSmart RR dclB 37° 165° dcm CpG

#R0653S 1,500 units\$64
#R0653L 7,500 units\$262

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	10	<10	100

5'...GGGCC...3'
3'...CCCGG...5'

Concentration: 20,000 units/ml

Methylation Sensitivity: Impaired by some combinations of overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

PspXI

rCutSmart RR dclB 37° 165° CpG

#R0656S 200 units\$72
#R0656L 1,000 units\$296

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	100	25	100

5'...VCTCGAGB...3'
3'...BGAGCTCV...5'

Concentration: 5,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by CpG methylation (see page 321).

Reaction Conditions: rCutSmart Buffer, 37°C

PstI

RR e dclC 37° 165°

#R0140S 10,000 units\$67
#R0140L 50,000 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	75	100	50*

for high (5X) concentration

#R0140T 10,000 units\$67
#R0140M 50,000 units\$274

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: *May exhibit star activity in this buffer.

5'...CTGCA⁺G...3'
3'...GACGTC...5'

Reaction Conditions: NEBuffer r3.1, 37°C. Heat inactivation: 80°C for 20 minutes.

PstI-HF[®]

High-Fidelity
rCutSmart RR e dclC 37° 165°

#R3140S 10,000 units\$67
#R3140L 50,000 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	75	50	100

for high (5X) concentration

#R3140T 10,000 units\$67
#R3140M 50,000 units\$274

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

5'...CTGCA⁺G...3'
3'...GACGTC...5'

Reaction Conditions: rCutSmart Buffer, 37°C

PvuI

RR e dclB 37° 165° CpG

#R0150S 500 units\$74
#R0150L 2,500 units\$304

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	25	100	<10

5'...CGAT⁺CG...3'
3'...GCTAGC...5'

Concentration: 10,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: NEBuffer r3.1, 37°C.

PvuI-HF[®]

High-Fidelity
rCutSmart RR e dclB 37° 165° CpG

#R3150S 500 units\$74
#R3150L 2,500 units\$304

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	100	100	100

5'...CGAT⁺CG...3'
3'...GCTAGC...5'

Concentration: 20,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C

SalI-HF®

rCutSmart    **37°**   **High-Fidelity**

#R3138S	2,000 units	\$65
#R3138L	10,000 units	\$266
for high (5X) concentration			
#R3138T	2,000 units	\$65
#R3138M	10,000 units	\$266

5'...GTCGAC...3'
3'...CAGCTG...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	100	100	100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

SapI

rCutSmart    **37°** 

#R0569S	250 units	\$70
#R0569L	1,250 units	\$287

5'...GCTCTTC(N)₁...3'
3'...CGAGAAAG(N)₄...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	50	>10	100

Concentration: 10,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Sau3AI

   **37°**  

#R0169S	200 units	\$64
#R0169L	1,000 units	\$262

5'...GATC...3'
3'...CTAG...5'

Reaction Conditions: NEBuffer r1.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	50	10	100*

Methylation Sensitivity: Unlike DpnII and MboI, Sau3AI is not blocked by *dam* methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

Note: *For added flexibility, NEB offers an isoschizomer or HF enzyme, supplied with rCutSmart Buffer.

Sau96I

rCutSmart   **37°**   

#R0165S	1,000 units	\$65
---------	-------------	-------	------

5'...G⁺GNCC...3'
3'...CCNGG...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	100	100

Methylation Sensitivity: Blocked by overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

SbfI

rCutSmart    **37°** 

#R0642S	500 units	\$75
#R0642L	2,500 units	\$308

5'...CCTGCA⁺GG...3'
3'...GGACGTCC...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 10,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	25	<10	100

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Star activity may result from a glycerol concentration of > 5%.

SbfI-HF®

rCutSmart    **37°**  **High-Fidelity**

#R3642S	500 units	\$75
#R3642L	2,500 units	\$308

5'...CCTGCA⁺GG...3'
3'...GGACGTCC...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	25	<10	100

Concentration: 20,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

ScaI-HF®

rCutSmart    **37°**  **High-Fidelity**

#R3122S	1,000 units	\$67
#R3122L	5,000 units	\$274
for high (5X) concentration			
#R3122M	5,000 units	\$274

5'...AGT⁺ACT...3'
3'...TCATGA...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	10	100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

ScrFI

rCutSmart   **37°**   

#R0110S	1,000 units	\$70
---------	-------------	-------	------

5'...CC⁺NGG...3'
3'...GGNCC...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	100	100

Methylation Sensitivity: Blocked by overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by overlapping CpG methylation (see p. 321).

Note: Star activity may result from extended digestion.



Recommended Buffer

Cloned at NEB

Recombinant Enzyme

EpiMark Validated

Engineered for Performance

Multiple Recognition Site

Time-Saver Qualified

Diluent Buffer

Incubation Temperature

Heat Inactivation

Methylation Sensitivity

SexAI

rCutSmart  RR  dil A 37°  dcm

#R0605S 200 units\$70
 #R0605L 1,000 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	75	50	100

5'...A[↓]CCWGGT...3'
 3'...TGGWCC[↓]A...5'

Concentration: 5,000 units/ml**Methylation Sensitivity:** Blocked by *dcm* methylation (see p. 321).**Reaction Conditions:** rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.**Note:** Star activity may result from a glycerol concentration of > 5%.

SfaNI

 RR  NEB r3.1  dil B 37°  CpG

#R0172S 300 units\$124

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	75	100	25

5'...GCATC(N)₆...3'
 3'...CGTAG(N)₉...5'

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by some combinations of overlapping CpG methylation (see p. 321).**Reaction Conditions:** NEBuffer r3.1, 37°C. Heat inactivation: 65°C for 20 minutes.**Concentration:** 2,000 units/ml**Note:** Star activity may result from a glycerol concentration of > 5%.

SfcI

rCutSmart  RR  dil B 37°  CpG

#R0561S 200 units\$67
 #R0561L 1,000 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	50	25	100

5'...C[↓]TRYAG...3'
 3'...GAYRT[↓]C...5'

Concentration: 10,000 units/ml**Methylation Sensitivity:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.**Reaction Conditions:** rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.**Note:** Star activity may result from a glycerol concentration of > 5%.

SfiI

rCutSmart  RR  2nd site  dil C 50°  dcm  CpG

#R0123S 3,000 units\$72
 #R0123L 15,000 units\$296

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	100	50	100

5'...GGCCNNNN[↓]NGGCC...3'
 3'...CCGGN[↓]NNNNCCGG...5'

Activity at 37°C: 10%**Methylation Sensitivity:** Impaired by overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by some combinations of overlapping CpG methylation (see p. 321).**Reaction Conditions:** rCutSmart Buffer, 50°C**Concentration:** 20,000 units/ml

SfoI

rCutSmart  RR  dil B 37°  dcm  CpG

#R0606S 500 units\$67
 #R0606L 2,500 units\$274

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	100	100

5'...GGC[↓]GCC...3'
 3'...CCG[↓]CGG...5'

Methylation Sensitivity: Blocked by some combinations of overlapping *dcm* methylation. Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).**Reaction Conditions:** rCutSmart Buffer, 37°C**Concentration:** 10,000 units/ml

SgrAI

rCutSmart  RR  2nd site  dil A 37°  CpG

#R0603S 1,000 units\$70
 #R0603L 5,000 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	10	100

5'...C[↓]RCGGYG...3'
 3'...GYGGC[↓]RC...5'

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).**Reaction Conditions:** rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.**Concentration:** 10,000 units/ml**Note:** Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

SmaI

rCutSmart  RR  dil B 25°  CpG

#R0141S 2,000 units\$65
 #R0141L 10,000 units\$266

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	<10	<10	100

5'...CCC[↓]GCG...3'
 3'...GGC[↓]CCC...5'

Concentration: 20,000 units/ml**Activity at 37°C:** 25%**Reaction Conditions:** rCutSmart Buffer, 25°C. Heat inactivation: 65°C for 20 minutes.**Methylation Sensitivity:** Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

SmlI

rCutSmart  RR  dil A 55°  CpG

#R0597S 500 units\$70
 #R0597L 2,500 units\$287

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	75	25	100

5'...C[↓]TRYAG...3'
 3'...GARYT[↓]C...5'

Activity at 37°C: 10%**Reaction Conditions:** rCutSmart Buffer, 55°C**Concentration:** 10,000 units/ml**Methylation Sensitivity:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

SnaBI

rCutSmart  RRII dIIA 37°  CpG

#R0130S	500 units\$62
#R0130L	2,500 units\$254
for high (5X) concentration		
#R0130M	2,500 units\$254

5'...TACGTA...3'
3'...ATGCA...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 5,000 and 25,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50*	50	10	100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Note: Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%. *May exhibit star activity in this buffer.

SpeI

rCutSmart  RRII dIIIC 37°  CpG

#R0133S	500 units\$70
#R0133L	2,500 units\$287
for high (5X) concentration		
#R0133M	2,500 units\$287

5'...A[▼]CTAGT...3'
3'...TGATCA...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	75	100	25	100

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

SpeI-HF[®]High-Fidelity
rCutSmart  RRII  dIIIC 37°  CpG

#R3133S	500 units\$70
#R3133L	2,500 units\$287
for high (5X) concentration		
#R3133M	2,500 units\$287

5'...A[▼]CTAGT...3'
3'...TGATCA...5'



Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	50	<10	100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

SphI

 RRII NEB r2.1 dIIB 37° 

#R0182S	500 units\$70
#R0182L	2,500 units\$287
for high (8X) concentration		
#R0182M	2,500 units\$287

5'...GCATG[▼]C...3'
3'...CGTACG...5'

Reaction Conditions: NEBuffer r2.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 and 80,000 units/ml

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	100	100	50	100*

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Star activity may result from extended digestion. *For added flexibility, NEB offers an isoschizomer or HF enzyme, supplied with rCutSmart Buffer.

SphI-HF[®]High-Fidelity
rCutSmart  RRII  dIIB 37°  CpG

#R3182S	500 units\$70
#R3182L	2,500 units\$287
for high (5X) concentration		
#R3182M	2,500 units\$287

5'...GCATG[▼]C...3'
3'...CGTACG...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	25	10	100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

SrfI

rCutSmart  RRII  dIIB 37°  CpG

#R0629S	500 units\$80
#R0629L	2,500 units\$363

5'...GCCCCG[▼]GGC...3'
3'...CGGGCCCG...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	50	0	100

Concentration: 20,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

SspI

 RRII  NEB dIIC 37° 

#R0132S	1,000 units\$72
#R0132L	5,000 units\$296
for high (5X) concentration		
#R0132M	5,000 units\$296

5'...AATA[▼]TTT...3'
3'...TTATAA...5'

Reaction Conditions: NEBuffer EcoRI/SspI, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	50	100	50	50

Concentration: 5,000 and 25,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

SspI-HF[®]High-Fidelity
rCutSmart  RRII  dIIB 37° 

#R3132S	1,000 units\$72
#R3132L	5,000 units\$296
for high (5X) concentration		
#R3132M	5,000 units\$296

5'...AATA[▼]TTT...3'
3'...TTATAA...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	25	100	<10	100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.



Recommended Buffer



StuI

rCutSmart     37°  

#R0187S 1,000 units\$65
#R0187L 5,000 units\$266

for high (10X) concentration

#R0187M 5,000 units\$266

5'...AGGCTT...3'
3'...TCCGGA...5'

Reaction Conditions: rCutSmart
Buffer, 37°C

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 50 100 50 100

Concentration: 10,000 and
100,000 units/ml

Methylation Sensitivity: Blocked
by overlapping *dcm* methylation
(see p. 321).

TaqI-v2

rCutSmart      65°  

#R0149S 4,000 units\$64
#R0149L 20,000 units\$262

for high (5X) concentration

#R0149T 4,000 units\$64
#R0149M 20,000 units\$262

5'...TCGA...3'
3'...AGCT...5'

Reaction Conditions: rCutSmart
Buffer, 65°C.

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 50 100 50 100

Concentration: 20,000 and
100,000 units/ml

Activity at 37°C: 10%

Methylation Sensitivity: Blocked
by overlapping *dam* methylation (see
p. 321).

StyI

StyI has been replaced by StyI-HF (engineered for performance).

StyI-HF®

 rCutSmart      37°  

#R3500S 3,000 units\$72
#R3500L 15,000 units\$296

5'...C'CW'WG...3'
3'...G'GW'WC...5'

Reaction Conditions: rCutSmart
Buffer, 37°C. Heat inactivation: 65°C
for 20 minutes.

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 25 100 25 100

Concentration: 20,000 units/ml

Methylation Sensitivity: Not
sensitive to *dam*, *dcm* or mammalian
CpG methylation.

TfiI

rCutSmart     65°  

#R0546S 500 units\$70

5'...G'AWTC...3'
3'...CTWAG...5'

Reaction Conditions: rCutSmart
Buffer, 65°C

Concentration: 10,000 units/ml

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 50 100 100 100

Activity at 37°C: 10%

Methylation Sensitivity: Cleavage
of mammalian genomic DNA is
blocked by some combinations of
overlapping CpG methylation (see
p. 321).

StyD4I

rCutSmart     37°   

#R0638S 200 units\$67

5'...CCNGG...3'
3'...GGNCC...5'

Reaction Conditions: rCutSmart
Buffer, 37°C. Heat inactivation: 65°C
for 20 minutes.

Concentration: 5,000 units/ml

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 10 100 100 100

Methylation Sensitivity: Blocked
by overlapping *dcm* methylation.
Cleavage of mammalian genomic
DNA is impaired by overlapping CpG
methylation. (see p. 321).

TseI

rCutSmart     65°  

#R0591L 500 units\$274

5'...G'CWGC...3'
3'...CGWCG...5'

Reaction Conditions: rCutSmart
Buffer, 65°C

Concentration: 10,000 units/ml.

Activity at 37°C: 10%

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 75 100 100 100

Methylation Sensitivity: Cleavage
of mammalian genomic DNA is
blocked by some combinations of
overlapping CpG methylation (see
p. 321).

Note: Star activity may result from a
glycerol concentration of > 5%.

SwaI

     25° 

#R0604S 2,000 units\$70
#R0604L 10,000 units\$287

5'...ATTTAAAT...3'
3'...TAAATTTA...5'

Reaction Conditions: NEBuffer r3.1,
25°C. Heat inactivation: 65°C for
20 minutes.

Concentration: 10,000 units/ml

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 10 10 100 10

Activity at 37°C: 25%

Methylation Sensitivity: Not
sensitive to *dam*, *dcm* or mammalian
CpG methylation.

Tsp45I

rCutSmart    65° 

#R0583S 400 units\$72
#R0583L 2,000 units\$296

5'...GTSAC...3'
3'...CASTG...5'

Reaction Conditions: rCutSmart
Buffer, 65°C

Concentration: 10,000 units/ml

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 100 50 <10 100

Activity at 37°C: 10%

Methylation Sensitivity: Not
sensitive to *dam*, *dcm* or mammalian
CpG methylation.

TspMI

rCutSmart   75°  CpG

#R0709S 200 units\$72

5'...C[▼]CCGGG...3'
3'...GGGCC[▲]C...5'

Reaction Conditions: rCutSmart Buffer, 75°C

Concentration: 5,000 units/ml

Activity at 37°C: 10%

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 50* 75* 50* 100

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Note: *May exhibit star activity in this buffer.

TspRI

rCutSmart    65° 

#R0582S 1,000 units\$76

5'...NNCASTGNN[▼]...3'
3'...[▲]NGTSACNN...5'

Reaction Conditions: rCutSmart Buffer, 65°C

Concentration: 10,000 units/ml

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 25 50 25 100

Activity at 37°C: 10%

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Tth111I

rCutSmart    65° 

#R0185S 400 units\$62

5'...GACN[▼]NGTC...3'
3'...CTGN[▲]NCAG...5'

Reaction Conditions: rCutSmart Buffer, 65°C

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 25 100 25 100

Concentration: 5,000 units/ml

Activity at 37°C: 10%

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

XbaI

rCutSmart    37°  

#R0145S 3,000 units\$72
#R0145L 15,000 units\$296

for high (5X) concentration

#R0145T 3,000 units\$72
#R0145M 15,000 units\$296

5'...T[▼]CTAGA...3'
3'...AGATC[▲]...5'





Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity <10 100 75 100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Blocked by overlapping *dam* methylation (see p. 321).

XcmI

   37° 

#R0533S 1,000 units\$70
#R0533L 5,000 units\$287

5'...C CANNNN[▼]NNNTGG...3'
3'...GGTNNNN[▲]NNNNNAACC...5'

Reaction Conditions: NEBuffer r2.1, 37°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 5,000 units/ml

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 10 100 25 100*

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Note: Star activity may result from extended digestion. *May exhibit star activity in this buffer.

XhoI

rCutSmart    37°  CpG

#R0146S 5,000 units\$72
#R0146L 25,000 units\$296
for high (5X) concentration
#R0146M 25,000 units\$296

5'...C[▼]TGAG...3'
3'...GAGCT[▲]C...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 75 100 100 100

Concentration: 20,000 and 100,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by CpG methylation (see p. 321).

XmaI

rCutSmart    37°  CpG

#R0180S 500 units\$70
#R0180L 2,500 units\$287
for high (5X) concentration
#R0180M 2,500 units\$287

5'...C[▼]CGGG...3'
3'...GGGCC[▲]C...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 25 50 <10 100

Concentration: 10,000 and 50,000 units/ml

Methylation Sensitivity: Cleavage of mammalian genomic DNA is impaired by CpG methylation (see p. 321).

Note: Star activity may result from a glycerol concentration of > 5%.

XmnI

rCutSmart    37° 

#R0194S 1,000 units\$65
#R0194L 5,000 units\$266

5'...GAANN[▼]NTTC...3'
3'...CTTN[▲]NNAAG...5'

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 50 75 <10 100

Concentration: 20,000 units/ml

Methylation Sensitivity: Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

ZraI

rCutSmart  RR  dmB 37°  CpG#R0659S 200 units\$67
#R0659L 1,000 units\$274NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 100 25 10 1005'...GAC^TGTC...3'
3'...CTG^ACAG...5'**Concentration:** 10,000 units/ml**Reaction Conditions:** rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.**Methylation Sensitivity:** Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Nicking Endonucleases

As a rule, when restriction endonucleases bind to their recognition sequences in DNA, they cleave both strands of the duplex at the same time. Two independent hydrolytic reactions proceed in parallel, most often driven by the presence of two catalytic sites within each enzyme, one for hydrolyzing each strand. One of the focuses of our research program is to engineer restriction enzymes so that they hydrolyze only one strand of the duplex, generating DNA molecules that are “nicked”, rather than cleaved. These conventional nicks (3'-hydroxyl, 5'-phosphate) can serve as initiation points for a variety of further enzymatic reactions such as replacement DNA synthesis, strand-displacement amplification (1), exonucleolytic degradation or the creation of small gaps (2).

Nicking endonucleases (NEases) are as simple to use as restriction endonucleases. Since the nicks generated by 6- or 7-base nicking endonucleases do not fragment DNA, their activities are monitored by conversion of supercoiled plasmids to open circles. Alternatively, substrates with nicking sites close enough on opposite strands to create a double-stranded cut can be used instead.

The uses of nicking endonucleases are still being explored. NEases can generate nicked or gapped duplex DNA for DNA mismatch repair studies and for diagnostic applications. The long overhangs that nicking enzymes make can be used in DNA fragment assembly. Nt.BbvCI has been used to generate long and non-complementary overhangs when used with XbaI in the USER® cloning protocol from NEB. Nicking endonucleases are also useful for isothermal DNA amplification, which relies on the production of site-specific nicks. For example, isothermal DNA amplification using Nt.BstNBI in concert with Vent® (exo-) DNA Polymerase (NEB #M0257) (EXPAR) has been reported for detection of a specific DNA sequence in a sample (3). Another isothermal DNA amplification technique [Nicking Endonuclease Mediated- DNA Amplification (NEMDA)] (4) has been described using the 3-base cutter Nt.CviPII and Bst DNA Polymerase I. Frequent cutting NEases can generate short partial duplex DNA fragments from genomic DNA. These fragments can be used for cloning or used as probes for hybridization-based applications. Nicking enzymes have also been used for genome mapping.

NEB continues to engineer more nicking enzymes, particularly in response to specific customer needs and applications.

References:

- (1) Walker, G.T. et al. (1992) *Proc. Natl. Acad. Sci. USA*, 89, 392–396.
- (2) Wang, H. and Hays, J.B. (2000) *Mol. Biotechnol.*, 15, 97–104.
- (3) Van Ness, J. et al. (2003) *Proc. Natl. Acad. Sci. USA*, 89, 4504–4509.
- (4) Chan, S.H. et al. (2004) *Nucl. Acids Res.*, 32, 6187–6199.

Nb.BbvCI

rCutSmart  RR  dmA 37° #R0631S 1,000 units\$72
#R0631L 5,000 units\$290NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity 25 100 100 1005'...CCTCAGC...3'
3'...GGAGT^ACG...5'**Concentration:** 10,000 units/ml**Reaction Conditions:** rCutSmart Buffer, 37°C. Heat inactivation: 80°C for 20 minutes.**Methylation Sensitivity:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Nb.BsmI

 RR  NEB r3.1  dmA 65° 

#R0706S 1,000 units\$72

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity <10 50 100 105'...GAATGCN...3'
3'...CTTAC^AGN...5'**Concentration:** 10,000 units/ml**Reaction Conditions:** NEBuffer r3.1, 65°C. Heat inactivation: 80°C for 20 minutes.**Activity at 37°C:** 25%**Methylation Sensitivity:** Not sensitive to *dam*, *dcm* or mammalian CpG methylation.

Nb.BsrDIrCutSmart                           

Nt.CviPII

rCutSmart   37°  CpG

#R0626S 40 units\$132

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	100	25	100

5'...CCD...3'
3'...GGH...5'

Concentration: 2,000 units/ml

D = A or G or T (not C)
H = A or C or T (not G)

Methylation Sensitivity: Cleavage of mammalian genomic DNA is blocked by CpG methylation (see p. 321).

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Homing Endonucleases

Homing endonucleases are double stranded DNases that have large, asymmetric recognition sites (12–40 base pairs) and coding sequences that are usually embedded in either introns or inteins (1). Introns are spliced out of precursor RNAs, while inteins are spliced out of precursor proteins (2,3). Homing endonucleases are named using conventions similar to those of restriction endonucleases with intron-encoded endonucleases containing the prefix, “I-” and intein endonucleases containing the prefix, “PI-”(1,7).

Homing endonuclease recognition sites are extremely rare. For example, an 18 base pair recognition sequence will occur only once in every 7×10^{10} base pairs of random sequence. This is equivalent to only one site in 20 mammalian-sized genomes (4). However, unlike standard restriction endonucleases, homing endonucleases tolerate some sequence degeneracy within their recognition sequence (5,6). As a result, their observed sequence specificity is typically in the range of 10–12 base pairs.

Homing endonucleases do not have stringently-defined recognition sequences in the way that restriction enzymes do. That is, single base changes do not abolish cleavage but reduce its efficiency to variable extents. The precise boundary of required bases is generally not known. The recognition sequence listed is one site that is known to be recognized and cleaved.

References:

- (1) Belfort, M. and Roberts, R.J. (1997) *Nucleic Acids Res.*, 25, 3379–3388.
- (2) Dujon, B. et al. (1989) *Gene*, 82, 115–118.
- (3) Perler, F.B. et al. (1994) *Nucleic Acids Res.*, 22, 1125–1127.
- (4) Jasin, M. (1996) *Trends in Genetics*, 12, 224–228.
- (5) Gimble, F.S. and Wang, J. (1996) *J. Mol. Biol.*, 263, 163–180.
- (6) Argast, M.G. et al. (1998) *J. Mol. Biol.*, 280, 345–353.
- (7) Roberts, R.J. et al. (2003) *Nucleic Acids Res.*, 31, 1805–1812.

I-CeuI

rCutSmart   37° 

#R0699S 500 units\$75

#R0699L 2,500 units\$308

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	10	10	<10	100

(Supplied with 5 µg of plasmid DNA)

Concentration: 5,000 units/ml

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Note: Homing endonucleases do not have stringently-defined recognition sequences.

Specificity: The homing or recognition site for this endonuclease is shown below:

5'...TAACATAACGGCTCTAAGGTAGCGAA...3'
3'...ATTGATATTGCCAGATTCCATCGCTT...5'

I-SceI

rCutSmart   37° 

#R0694S 500 units\$72

#R0694L 2,500 units\$296

NEBuffer	r1.1	r2.1	r3.1	rCutSmart
% Activity	<10	50	25	100

(Supplied with 5 µg of plasmid DNA)

Concentration: 5,000 units/ml

Reaction Conditions: rCutSmart Buffer, 37°C. Heat inactivation: 65°C for 20 minutes.

Note: Homing endonucleases do not have stringently-defined recognition sequences. Store at -80°C.

Specificity: The homing or recognition site for this endonuclease is shown below:

5'...TAGGGATAACAGGGTAAT...3'
3'...ATCCCTATTGTCCCATTA...5'

PI-PspI



#R0695S 500 units\$72

(Supplied with 5 µg of plasmid DNA)

Reaction Conditions: NEBuffer
PI-PspI + Recombinant Albumin, 65°C

Specificity: The homing or recognition site for this endonuclease is shown below:

5'...TGGCAACAGCTATTATGGGTATTATGGGT...3'
3'...ACCGTTTGTGCTAATACCCATAATACCA...5'

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity <10 10 10 10

Concentration: 5,000 units/ml

Activity at 37°C: 5%

Note: Homing endonucleases do not have stringently-defined recognition sequences.

Recombinant Albumin, Molecular Biology Grade

#B9200S 12 mg\$44

Companion Product:

BSA, Molecular Biology Grade

#B9000S 12 mg\$32

Reaction Conditions: Recombinant Albumin, Molecular Biology Grade, is a non-bovine derived albumin that can serve as an alternative to Bovine Serum Albumin (BSA). Like BSA, it has been shown to prevent adhesion of enzymes to reaction tubes and pipette surfaces. It also stabilizes some proteins during incubation. Choose Recombinant Albumin, when there is a need to avoid BSA.

PI-SceI



#R0696S 250 units\$72

(Supplied with 5 µg of plasmid DNA)

Reaction Conditions: NEBuffer
PI-SceI + Recombinant Albumin, 37°C.
Heat inactivation: 65°C for 20 minutes.

Specificity: The homing or recognition site for this endonuclease is shown below:

5'...ATCTATGTGGGTGCGGAGAAAGAGGTAATGAAATGG...3'
3'...TAGATACAGCCACGCCCTCTTCTCCATTACTTACC...5'

NEBuffer r1.1 r2.1 r3.1 rCutSmart
% Activity <10 <10 <10 <10


Concentration: 5,000 units/ml

Note: Homing endonucleases do not have stringently-defined recognition sequences.

NEB Tube Opener

#C1008S 2 units\$21

Description: Use to open a variety of microcentrifuge tubes. Can be used for snap-on caps or screw-on caps.



Reaction Buffers

NEBuffer™ 1.1 (10X)
#B7201S 5.0 ml\$25

NEBuffer 2.1 (10X)
#B7202S 5.0 ml\$25

NEBuffer 3.1 (10X)
#B7203S 5.0 ml\$25

CutSmart Buffer (10X)
#B7204S 5.0 ml\$25

NEBuffer Set (EcoRI/SspI, DpnII)
#B7006S 2.5 ml of each\$25

NEBuffer Set (1.1, 2.1, 3.1 & CutSmart)
#B7200S 1.25 ml of each\$25

NEW
NEBuffer Set (r1.1, r2.1, r3.1 & rCutSmart™)
#B7030S 1.25 ml of each\$25

S-adenosylmethionine (SAM)
#B9003S (32 mM) 0.5 ml\$40

Nuclease-free Water
#B1500S 25 ml\$27
#B1500L 100 ml\$63

Description: New England Biolabs provides a color-coded 10X NEBuffer with each restriction endonuclease to ensure optimal (100%) activity. Most of our enzymes are supplied with one of four standard NEBuffers. Occasionally, an enzyme has specific buffer requirements not met by one of the four standard NEBuffers, in which case the enzyme is supplied with its own unique NEBuffer.

The NEBuffer Sets contain multiple vials of buffers, which are indicated in the product name. NEBuffer Set (1.1, 2.1, 3.1 & CutSmart) is formulated with BSA, while NEBuffer Set (r1.1, r2.1, r3.1 and rCutSmart) is formulated with Recombinant Albumin.

Nuclease-free Water is ideal for the preparation of reagents and for use in enzymatic reactions. No toxic agents, such as DEPC, are used in the production of this water, so as to avoid inhibition in enzymatic reactions.

For more information about NEBuffer compatibility, please turn to page 309.

Reaction Buffer Compositions: Visit www.neb.com for details.

Diluent Buffers

Diluent A #B8001S	5.0 ml \$34
Diluent B #B8002S	5.0 ml \$34
Diluent C #B8003S	5.0 ml \$34

Description: Diluent Buffers are recommended for making dilutions of restriction endonucleases. When necessary, we recommend diluting enzymes just prior to use and suggest that the final concentration of diluted enzymes be at least 1,000 units/ml. Diluent preference for each restriction endonuclease is listed with its catalog entry; a complete listing of all restriction endonucleases and their diluents can be found on page 309.

Storage Conditions:

Store at -20°C.

Diluent Buffer Compositions: Visit www.neb.com for details.

Gel Loading Dyes

Gel Loading Dye, Blue (6X) #B7021S	4 ml \$50
Gel Loading Dye, Orange (6X) #B7022S	4 ml \$50
Gel Loading Dye, Purple (6X) #B7024S	4 ml \$45
Gel Loading Dye, Purple (6X), no SDS #B7025S	4 ml \$45

Description: Pre-mixed loading dye solutions are available with a choice of tracking dyes, for agarose and non-denaturing polyacrylamide gel electrophoresis. Three of the solutions contain SDS, which often results in sharper bands, as some restriction enzymes remain bound to their DNA substrates following cleavage. Some interference may be observed when using SYBR® or GelRed® as precast dyes in the presence of increased concentrations of SDS. When using these dyes as precast dyes, NEB recommends using our Gel Loading Dye, Purple (6X), no SDS (NEB #B7025). Gel Loading Dyes contain EDTA to chelate up to 10 mM magnesium, thereby stopping the reaction. Bromophenol blue is the standard tracking dye for electrophoresis. It migrates at approximately 300 bp on a standard 1% TBE agarose gel. Orange G will not appear in gel photographs, and runs ahead of all but the smallest restriction fragments, migrating at approximately 50 bp on a standard 1% TBE agarose gel. NEB also offers a unique purple dye which migrates similarly to Bromophenol blue. Specifically chosen, this dye does not leave a shadow under UV light.

Gel Loading Dye Compositions: Visit www.neb.com for details.

Note: Use 5 µl of Gel Loading Dye per 25 µl reaction, or 10 µl per 50 µl reaction. Mix well before loading gel. Store at room temperature.

SYBR® is a registered trademark of Molecular Probes, Inc.
GELRED® is a registered trademark of Biotium.

Scott has been with NEB for 2 years as a Purification Scientist in our Rowley Production Facility. He primarily works on downstream processing of products, but also helps with enzyme purification. In his spare time, Scott enjoys cooking and spending time with friends and family outdoors.





Aim AR camera at this page
to see image in motion.



Conserving Biodiversity

Biodiversity is a broad concept that refers to the “web of life”. It describes all the organisms on Earth (species diversity), the variety of organisms in an ecosystem (ecosystem diversity), and the genetic make-up of those organisms (genetic diversity). Biodiversity is responsible for sustaining life as we know it — from producing oxygen, to pollinating our food sources and cleaning our water supplies.

It is estimated that 80% of the world’s biodiversity exists underwater, within the just over 70% of Earth’s surface that is covered by ocean. And, if you have had the opportunity to explore a coastline or to gaze underwater, this likely comes as no surprise. Teeming with creatures from microscopic to the largest known on Earth (the blue whale), ocean life is amazingly diverse.

It is easy to understand why Don Comb’s curiosity was piqued by the marine ecosystems he visited as he pursued his Ph.D. at the University of Michigan. His graduate research first took him to the Bermuda Biological Research Station, where he studied sea urchin development. Firsthand experience observing rich and diverse marine organisms and their interactions within the ecosystem inspired Don to make protecting biodiversity a priority.

In the 1970s, Don supported the Conservation Law Foundation’s efforts to stop offshore drilling on George’s Bank, an area recognized for its diverse marine biology. Not long after, he helped establish a marine sanctuary in the waters off the island of Saint Barthélemy through the New England Biolabs Foundation (NEBF), which was established shortly after NEB became profitable. This sanctuary protects breeding grounds for local fish, as well as a diversely inhabited coral reef ecosystem.

In 2000 Don read an opinion piece in *Science*, which asked the question: Why aren’t we saving DNA from the world’s threatened and endangered species? Don realized that it wouldn’t require new technology to do so, just someone with the initiative. In 2001, he founded The Ocean Genome Legacy (OGL), a non-profit environmental research organization and DNA bank dedicated to promoting new methods for studying and conserving marine species through preservation and genome analysis. As part of its mission, OGL created the Ocean Genome Resource Collection, a publicly accessible genome biorepository containing samples of DNA from species around the world, abundant and rare, who call the ocean home. In 2013, this priceless collection was moved to Northeastern University’s Marine Science Center, located in Nahant, MA. Here it continues to help uncover some of our oceans’ deepest mysteries and reveal genomic information that can help cure diseases, improve the sustainability of global food and energy supplies, and protect the environment.

Don challenged everyone to consider their impact on the planet’s biodiversity. In the 2002/2003 NEB Catalog, he included a personal note, urging readers to “learn to live well on less and remember that small is beautiful. We all need to recognize that the earth’s biodiversity is our future.” Through the NEBF, Don has encouraged individuals to do just that, all over the world, by providing support and funding to communities in biodiversity-rich areas and the organizations working to conserve and protect them.

Sea Urchin moving its parts.
Credit: Uw.art, Adobe Stock

View the calming
motions of the
sea urchin.



DNA Polymerases & Amplification Technologies



NEB has pursued the discovery & development of DNA polymerases for over 25 years.

As the first company to sell *Taq* DNA polymerase and the first to discover a PCR-stable high-fidelity DNA polymerase, NEB established a foundation in developing innovative, high quality tools for the research and diagnostic community.

PCR & qPCR

NEB's product portfolio features a large selection of polymerases for PCR. Q5 High-Fidelity DNA Polymerase offers fidelity 280 times higher than *Taq*, along with superior performance with minimal optimization. One *Taq* is ideal for robust amplification in routine PCR applications. Both are available in hot start formulations, utilizing a novel aptamer-based approach that does not require an activation step.

Fluorescence-based quantitative real-time PCR (qPCR) is the gold standard for detection and quantitation of nucleic acids due to its specificity and sensitivity. Luna qPCR & RT-qPCR products feature an inert blue tracking dye for easy reaction setup and are available for intercalating dye or fluorescent probe-based detection methods.


Isothermal Amplification

Sequence specific isothermal amplification approaches eliminate the need for temperature cycling, providing advantages for certain point-of-care diagnostic needs. NEB's broad suite of reagents continue to enable advancement in isothermal amplification. Nicking enzymes, WarmStart enzymes, strand-displacing DNA polymerases, and RNA polymerases offer flexibility to assemble and design an isothermal amplification platform.

Featured Products

- 66** Q5® High-Fidelity DNA Polymerase
- 68** One *Taq*® DNA Polymerase
- 72** Luna® qPCR & RT-qPCR Products
- 78** WarmStart® Colorimetric LAMP/RT-LAMP 2X Master Mix

Featured Tools & Resources

- 65** PCR Polymerase Selection Chart
 - 322** Guidelines for PCR Optimization
 - 323** PCR Troubleshooting Guide
-  Visit www.neb.com/PCR to find additional online tools, video tech tips and tutorials to help you in setting up your PCR experiments.



Find an overview of PCR.



Amplification-based Molecular Diagnostics	64
PCR Polymerase Selection Chart	65

High Fidelity PCR

Q5 High-Fidelity DNA Polymerase	66
Q5 High-Fidelity 2X Master Mix	66
Q5 Hot Start High-Fidelity DNA Polymerase	66
Q5 Hot Start High-Fidelity 2X Master Mix	66
Q5U Hot Start High-Fidelity DNA Polymerase	66
Q5 High-Fidelity PCR Kit	66
Q5 Site-Directed Mutagenesis Kit	98
Q5 Site-Directed Mutagenesis Kit (without competent cells)	98
Phusion High-Fidelity DNA Polymerase	67
Phusion High-Fidelity PCR Master Mix (HF/GC) Buffer	67
Phusion Hot Start Flex DNA Polymerase	67
Phusion Hot Start Flex 2X Master Mix	67
Phusion High-Fidelity PCR Kit	67

Routine PCR

OneTaq DNA Polymerase (Std., Quick-Load)	68
OneTaq Master Mixes (Std., Quick-Load)	68
OneTaq RT-PCR & One-step RT-PCR Kits	68
OneTaq Hot Start DNA Polymerase	68
OneTaq Hot Start Master Mixes (Std., Quick-Load)	68
Taq DNA Polymerase w/ThermoPol Buffer	69
Taq DNA Polymerase w/Standard Taq Buffer	69
Taq DNA Polymerase w/Standard Taq (Mg-free) Buffer	69
Taq Master Mixes (2X, 5X, Quick-Load)	69
Taq PCR Kit	69
Multiplex PCR 5X Master Mix	69
Hot Start Taq DNA Polymerase	69
Hot Start Taq 2X Master Mix	69

Specialty PCR

LongAmp Taq DNA Polymerase	70
LongAmp Hot Start Taq DNA Polymerase	70
LongAmp Taq 2X Master Mix	70
LongAmp Hot Start Taq 2X Master Mix	70
LongAmp Taq PCR Kit	70
Hemo KlenTaq	70
EpiMark Hot Start Taq DNA Polymerase	71

Other PCR Polymerases

Vent/Vent (exo-) DNA Polymerase	71
Deep Vent/Deep Vent (exo-) DNA Polymerase	71

qPCR & RT-qPCR

Luna qPCR & RT-qPCR	72
Find the Right Luna Product for Your Application	72
Luna Universal qPCR Master Mix	73
Luna Universal Probe qPCR Master Mix	73
LunaScript RT SuperMix Kit	73
LunaScript RT Master Mix (Primer-free)	73
Luna Probe One-Step RT-qPCR Kit (No ROX)	74
Luna Probe One-Step RT-qPCR 4X Mix with UDG	74
Luna Universal One-Step RT-qPCR Kit	74
Luna Universal Probe One-Step RT-qPCR Kit	74
Luna Cell Ready One-Step RT-qPCR Kit	75
Luna Cell Ready Probe One-Step RT-qPCR Kit	75
Luna Cell Ready Lysis Module	75

SARS-CoV-2 Detection

SARS-CoV-2 Rapid Colorimetric LAMP Assay Kit	76
SARS-CoV-2 Positive Control (N gene)	76
WarmStart Fluorescent LAMP/RT-LAMP Kit w/UDG	77
WarmStart Multi-Purpose LAMP/RT-LAMP 2X Master Mix (with UDG)	77
LAMP Fluorescent Dye	77
Luna SARS-CoV-2 RT-qPCR Multiplex Assay Kit	77

Isothermal Amplification and Strand Displacement

WarmStart Colorimetric LAMP/RT-LAMP 2X Master Mix	78
WarmStart Colorimetric LAMP/RT-LAMP 2X Master Mix with UDG	78
WarmStart Fluorescent LAMP/RT-LAMP Kit	78
Bst DNA Polymerase, (Large Frag. and Full-Length)	79
Bst 2.0 DNA Polymerase	79
Bst 2.0 WarmStart DNA Polymerase	79
Bst 3.0 DNA Polymerase	79
IsoAmp II Universal tHDA Kit	80
phi29 DNA Polymerase	80

RT-PCR

OneTaq One-Step RT-PCR Kit	81
OneTaq RT-PCR Kit	81

Polymerases for DNA Manipulation (labeling, blunting, etc.)

PreCR Repair Mix	82
Sulfolobus DNA Polymerase IV	82
Therminator DNA Polymerase	82
DNA Polymerase I (<i>E. coli</i>)	83
DNA Polymerase I, Large (Klenow) Fragment	83
Klenow Fragment (3'→5' exo-)	83
T4 DNA Polymerase	83
T7 DNA Polymerase (unmodified)	84
Bsu DNA Polymerase, Large Fragment	84
Terminal Transferase	84

Nucleotide Solutions & Buffers

Reaction Buffers	85
Acyclonucleotide Set	85
Deoxynucleotide (dNTP) Solution Mix/Set	85
Ribonucleotide Solution Mix/Set	85
7-deaza-dGTP	85
Adenosine-5' Triphosphate (ATP)	85
5-methyl-dCTP	85
dATP/dUTP/dGTP Solutions	85

cDNA Synthesis

cDNA Synthesis Selection Chart	86
LunaScript RT SuperMix Kit	73
LunaScript RT Master Mix Kit (Primer-free)	73
ProtoScript II First Strand cDNA Synthesis Kit	209
ProtoScript First Strand cDNA Synthesis Kit	209
Template Switching RT Enzyme Mix	208
ProtoScript II Reverse Transcriptase	206
M-MuLV Reverse Transcriptase	207
AMV Reverse Transcriptase	207
WarmStart RTx Reverse Transcriptase	207

PCR Cleanup

Monarch PCR & DNA Cleanup Kit (5 µg)	87
Exo-CIP Rapid PCR Kit	87

Recombinant Enzyme

One or more of these products are covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at qbd@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications. Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at <https://www.neb.com/support/terms-of-sale>. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

Amplification-based Molecular Diagnostics

NEB has a long history in the development of reliable and convenient tools for amplification, and offers a large selection of products for PCR, qPCR, RT-qPCR and isothermal amplification. Our extensive expertise in this area has allowed us to develop optimized enzymes for a variety of applications, including incorporation into diagnostics. The table below summarizes products available from NEB for diagnostic applications. Bulk and/or custom formats are available for all products, with more specific customization details included below. Learn more at www.neb.com/MDx.

APPLICATION	PRODUCTS	PRODUCT NOTES	CUSTOM FORMULATIONS AVAILABLE
PCR APPLICATIONS	DNA, Dye Luna Universal qPCR Master Mix (NEB #M3003) DNA, Probe Luna[®] Universal Probe qPCR Master Mix (NEB #M3004)	<ul style="list-style-type: none"> Sensitive, reproducible and reliable performance Compatible with automated liquid handling and reaction miniaturization Room temperature stable for ≥ 24 hours 	<ul style="list-style-type: none"> ROX-free Blue-dye-free Lyo-compatible
	qPCR/RT-qPCR RNA (1-step), Dye Luna Universal One-Step RT-qPCR Kit (NEB #E3005) RNA (1-step), Probe Luna Probe One-Step RT-qPCR 4X Mix with UDG (NEB #M3019) Luna Universal Probe One-Step RT-qPCR Kit (NEB #E3006) Luna Probe One-Step RT-qPCR Kit (No ROX) (NEB #E3007) Luna SARS-CoV-2 RT-qPCR Multiplex Assay Kit (NEB #E3019)	<ul style="list-style-type: none"> Luna WarmStart RT paired with Hot Start <i>Taq</i> increases reaction specificity and robustness Compatible with automation and reaction miniaturization Room temperature stable for ≥ 24 hours High conc. ideal for viral targets (NEB #M3019) Includes carryover prevention (NEB #M3019) 	<ul style="list-style-type: none"> Rox-free Blue-dye-free Lyo-compatible
	RNA (2-step) LunaScript[®] RT SuperMix Kit (NEB #E3010)	<ul style="list-style-type: none"> Novel thermostable RT Single-tube format 13-minute protocol 	<ul style="list-style-type: none"> Blue-dye-free
	PCR/RT-PCR Master Mixes Q5[®] Hot Start High-Fidelity 2X Master Mix (NEB #M0494) Q5 High-Fidelity 2X Master Mix (NEB #M0492) Standalone Enzyme & Buffer Q5 Hot Start High-Fidelity DNA Polymerase (NEB #M0493) Q5 High-Fidelity DNA Polymerase (NEB #M0491)	<ul style="list-style-type: none"> ~280X fidelity of <i>Taq</i> Consistent, fast, reliable performance Compatible with automation and reaction miniaturization Room temperature stable for ≥ 24 hours 	<ul style="list-style-type: none"> High conc. Glycerol-free
	Hemo KlenTaq (NEB #M0332)	<ul style="list-style-type: none"> Amplification direct from blood 	<ul style="list-style-type: none"> +/- Hot Start High conc.
	Hot Start <i>Taq</i> DNA Polymerase (NEB #M0495) Hot Start <i>Taq</i> 2X Master Mix (NEB #M0496)	<ul style="list-style-type: none"> Unique aptamer-based enzyme control supports fast protocols Compatible with automation and reaction miniaturization 	<ul style="list-style-type: none"> High conc. Glycerol-free
	WarmStart[®] Colorimetric LAMP 2X Master Mix (DNA & RNA) (NEB #M1800) WarmStart Colorimetric LAMP 2X Master Mix with UDG (NEB #M1804)	<ul style="list-style-type: none"> Fast, clear pink-to-yellow visible detection of amplification Results in approximately 30 minutes 	
	SARS-CoV-2 Rapid Colorimetric LAMP Assay Kit (NEB #E2019)	<ul style="list-style-type: none"> Simple, colorimetric detection of amplification of SARS-CoV-2 nucleic acid Automation-compatible when coupled with absorbance plate reader 	
ISOTHERMAL APPLICATIONS	LAMP WarmStart Fluorescent LAMP/RT-LAMP Kit (NEB #E1700) WarmStart Fluorescent LAMP/RT-LAMP Kit (with UDG) (NEB #E1708) WarmStart Multi-Purpose LAMP/RT-LAMP 2X Master Mix (with UDG) (NEB #M1708)	<ul style="list-style-type: none"> Master mix for LAMP and RT-LAMP workflows Supports multiple detection methods, including fluorescence and turbidity Automation compatible 	<ul style="list-style-type: none"> Lyo-compatible High conc.
	Bsf 2.0 WarmStart DNA Polymerase (NEB #M0538) Bsf 2.0 DNA Polymerase (NEB #M0537)	<ul style="list-style-type: none"> Improved reaction properties compared to wild-type <i>Bsf</i> DNA Polymerase Increased dUTP tolerance enables carryover prevention 	<ul style="list-style-type: none"> Glycerol-free High conc.
	Bsf 3.0 DNA Polymerase (NEB #M0374)	<ul style="list-style-type: none"> DNA binding domain fusion supports robust performance Significantly increased RT activity up to 72°C enables single enzyme RT-LAMP 	<ul style="list-style-type: none"> Glycerol-free High conc.
	WarmStart RTx Reverse Transcriptase (NEB #M0380)	<ul style="list-style-type: none"> <i>In-silico</i> designed RT for RT-LAMP with reversibly-bound aptamer that inhibits activity below 40°C 	<ul style="list-style-type: none"> Glycerol-free High conc.
	Strand Displacement Nt.BstNBI (NEB #R0607)	<ul style="list-style-type: none"> High purity, high quality nicking endonuclease 	<ul style="list-style-type: none"> Glycerol-free High conc.
	Helicase-dependent Amplification Tte UvrD Helicase (NEB #M1202) IsoAmp II Universal tHDA Kit (NEB #H0110)	<ul style="list-style-type: none"> Thermostable Improves specifically of problematic fluorescent LAMP reactions Requires only two primer Produces short, discrete DNA products 	<ul style="list-style-type: none"> High conc.
	Other Bsu DNA Polymerase Large Fragment (NEB #M0330) T4 Gene 32 Protein (NEB #M0300) Deoxynucleotide (dNTP) Solution Mix (NEB #N0447)	<ul style="list-style-type: none"> Enables low temperature isothermal applications Can increase yield and efficiency of amplification reactions Highly pure Individual mixes available 	<ul style="list-style-type: none"> High conc. Glycerol-free Custom conc.
	Antarctic Thermolabile UDG (NEB #M0372)	<ul style="list-style-type: none"> Unique thermolabile version is completely inactivated in typical isothermal and RT-qPCR workflows 	<ul style="list-style-type: none"> High conc.
	Proteinase K, Molecular Biology Grade (NEB #P8107)		<ul style="list-style-type: none"> Custom conc.
	Thermolabile Proteinase K (NEB #P8111)	<ul style="list-style-type: none"> Unique thermolabile version is completely inactivated in typical isothermal and RT-qPCR workflows 	<ul style="list-style-type: none"> Custom conc.

PCR Polymerase Selection Chart

For over 45 years, New England Biolabs, Inc. has been a world leader in the discovery and production of reagents for the life science industry. NEB offers a wide range of DNA polymerases, and through our commitment to research, ensures the development of innovative and high quality tools for PCR and related applications. The following table simplifies the selection of a polymerase that best suits your experiment.

High-fidelity polymerases benefit from a Tm+3 annealing temp. Use the **NEB Tm Calculator** to ensure successful PCR at **TmCalculator.neb.com**

★ indicates recommended choice for application

	STANDARD PCR		HIGH-FIDELITY PCR			SPECIALTY PCR			
	One Taq®/ One Taq Hot Start	Taq / Hot Start Taq	HIGHEST FIDELITY Q5®/Q5 Hot Start	Phusion®(4)/ Phusion Flex	MODERATE FIDELITY Vent®/ Deep Vent™	LONG AMPLICONS LongAmp®/ LongAmp Hot Start Taq	dU TOLERANCE Epimark® Hot Start Taq	Q5U™	BLOOD DIRECT PCR Hemo KlenTaq®
PROPERTIES									
Fidelity vs. Taq	2X	1X	~280X ⁽²⁾	> 39X	5–6X	2X	1X	ND	ND
Amplicon Size	< 6 kb	≤ 5 kb	≤ 20 kb	≤ 20 kb	≤ 6 kb	≤ 30 kb	≤ 1 kb	app-specific	≤ 2 kb
Extension Time	1 kb/min	1 kb/min	6 kb/min	4 kb/min	1 kb/min	1.2 kb/min	1 kb/min	2 kb/min	0.5 kb/min
Resulting Ends	3' A/Blunt	3' A	Blunt	Blunt	Blunt	3' A/Blunt	3' A	Blunt	3' A
3'→5' exo	Yes	No	Yes	Yes	Yes	Yes	No	Yes	No
5'→3' exo	Yes	Yes	No	No	No	Yes	Yes	No	No
Units/50 µl Reaction	1.25	1.25	1.0	1.0	0.5–1.0	5.0	1.25	1.0	N/A
Annealing Temperature	Tm-5	Tm-5	Tm+3	Tm+3	Tm-5	Tm-5	Tm-5	Tm+3	Tm-5
APPLICATIONS									
Routine PCR	★	•	•	•	•	•			
Colony PCR	★	•							
Enhanced Fidelity	•		★	•	•	•			
High Fidelity			★	•					
High Yield	★	•	★	•					
Fast			★	•					
Long Amplicon			★	•		★			
GC-rich Targets	★		★		•	•			
AT-rich Targets	★	•	★	•		•	•	★	
High Throughput	•	•	•	•			•	★	
Multiplex PCR	•	★ ⁽¹⁾	•	•					
Extraction-free PCR									•
DNA Labeling		★							
Site-directed Mutagenesis			★	•					
Carryover Prevention							•	★	
USER® Cloning							•	★	
NGS APPLICATIONS									
NGS Library Amplification			★ ⁽³⁾	•			★ ⁽⁵⁾		
FORMATS									
Hot Start Available	•	•	•	•		•	•	•	
Kit		•	•	•		•		•	
Master Mix Available	•	•	•	•		•		•	
Direct Gel Loading	•	•							

(1) Use Multiplex PCR 5X Master Mix.

(2) Due to the very low frequency of misincorporation events being measured, the error rate of high-fidelity enzymes like Q5 is challenging to measure in a statistically significant manner. We continue to investigate improved assays to characterize Q5's very low error rate to ensure that we present the most robust accurate fidelity data possible (Popatov, V. and Ong, J.L. (2017) PLoS One, 12(1):e0169774. doi 10.1371/journal.pone.0169774).

(3) Use NEBNext High-Fidelity 2X PCR Master Mix.

(4) Phusion DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific.

(5) Use NEBNext Enzymatic Methyl-seq Kit (EM-seq™) for Illumina.

Why choose
Q5 for your
PCR?



Q5[®] High-Fidelity DNA Polymerase

RR NEB U HiFi PCR Tm³

Q5 Hot Start High-Fidelity DNA Polymerase

RR NEB U HiFi PCR Tm³

Q5 High-Fidelity DNA Polymerase	
#M0491S 100 units\$113
#M0491L 500 units\$450

Q5 High-Fidelity 2X Master Mix	
#M0492S 100 rxns (50 µl vol)\$185
#M0492L 500 rxns (50 µl vol)\$739

Q5 Hot Start High-Fidelity DNA Polymerase	
#M0493S 100 units\$140
#M0493L 500 units\$554

Q5 Hot Start High-Fidelity 2X Master Mix	
#M0494S 100 rxns (2 x 1.25 ml)\$217
#M0494L 500 rxns (10 x 1.25 ml)\$862
#M0494X 500 rxns (1 x 12.5 ml)\$862

NOW AVAILABLE: NEBNext Ultra II Q5 Master Mix – see page 169.

Q5 POLYMERASE DETAILS	
Extension Rate	6 kb/min
Amplicon Size	≤ 20 kb
Fidelity	~ 280X <i>Taq</i>
Units / 50 µl rxn	1 unit
Resulting Ends	Blunt
3'→5' Exonuclease Activity	Yes
5'→3' Exonuclease Activity	No
Supplied Buffer	Q5 Reaction Buffer
Supplied Enhancer	Q5 High GC Enhancer
Extraction-free PCR	Yes
PRODUCT FORMATS	
Hot Start Available	Yes
Activation Required	No
Master Mix Available	Yes
PCR Kit Available	Yes
NGS Version Available	Yes
APPLICATIONS	
High-Fidelity PCR	Yes
Difficult PCR	Yes
High GC PCR	Yes
T/A, U/A Cloning	No
Colony PCR	No
Blunt Cloning	Yes
Multiplex PCR	Yes
USER Cloning	Yes (Q5U)
Carryover Prevention	Yes (Q5U)

APPLICATIONS	
High-Fidelity PCR	Yes
Difficult PCR	Yes
High GC PCR	Yes
T/A, U/A Cloning	No
Colony PCR	No
Blunt Cloning	Yes
Multiplex PCR	Yes
USER Cloning	Yes (Q5U)
Carryover Prevention	Yes (Q5U)

Tm Calculator
High-Fidelity polymerases benefit from a Tm³ annealing temp. Use the **NEB Tm Calculator** to ensure successful PCR at TmCalculator.neb.com

Q5U [™] Hot Start High-Fidelity DNA Polymerase	
#M0515S 100 units\$171
#M0515L 500 units\$660

Q5 High-Fidelity PCR Kit	
#E0555S 50 rxns (50 µl vol)\$109
#E0555L 200 rxns (50 µl vol)\$334

Q5 Site-Directed Mutagenesis Kit	
#E0554S 10 rxns\$201

Q5 Site-Directed Mutagenesis Kit (without competent cells)	
#E0552S 10 rxns\$139

For a complete listing of Deoxynucleotide Solutions, see page 85.

For a complete listing of Reaction Buffers, see pages 85.

Description: Q5 High-Fidelity DNA Polymerase sets a new standard for both fidelity and performance. With the highest fidelity amplification available (~280 times higher than *Taq*), Q5 DNA Polymerase results in ultra-low error rates. Q5 DNA Polymerase is composed of a novel polymerase that is fused to the processivity-enhancing Sso7d DNA binding domain, improving speed, fidelity and reliability of performance.

The Q5 buffer system is designed to provide superior performance with minimal optimization across a broad range of amplicons, regardless of GC content. For routine or complex amplicons up to ~65% GC content, Q5 Reaction Buffer provides reliable and robust amplification. For amplicons with high GC content (>65% GC), addition of the Q5 High GC Enhancer ensures continued maximum performance.


Q5 Hot Start DNA Polymerase: In contrast to chemically-modified or antibody-based hot start polymerases, NEB's Q5 Hot Start utilizes a unique synthetic aptamer. This structure binds to the polymerase through non-covalent interactions, blocking activity during the reaction setup. The polymerase is activated during normal cycling conditions, allowing reactions to be set up at room temperature. Q5 Hot Start does not require a separate high temperature activation step, reducing the potential for sample degradation, shortening reaction times and increasing ease-of-use. Q5 Hot Start is an ideal choice for high specificity amplification and provides robust amplification of a wide variety of amplicons, regardless of GC content.

Q5U Hot Start High-Fidelity DNA Polymerase: A modified version of Q5 High-Fidelity DNA Polymerase capable of incorporating dUTP for carryover prevention. Q5U is also compatible with USER cloning methods and enables the amplification of bisulfite treated/deaminated DNA.

Additional Formats: For added convenience, Q5 DNA Polymerase is available in master mix format or as a kit. Master mix formulations include dNTPs, Mg⁺⁺ and all necessary buffer components. The Q5 High-Fidelity PCR Kit contains the Q5 High-Fidelity 2X Master Mix, nuclease-free water and the Quick-Load Purple 1 kb Plus DNA Ladder. For information on the Q5 Site-Directed Mutagenesis Kit, with or without competent cells, see page 98.


Concentration: 2,000 units/ml

Visit Q5PCR.com for more information.





Robust amplification with Q5 (A) and Q5 Hot Start (B) High-Fidelity DNA Polymerases, regardless of GC content.

Amplification of a variety of human genomic amplicons from low-to-high GC content using either Q5 or Q5 Hot Start High-Fidelity DNA Polymerase. Reactions using Q5 Hot Start were set up at room temperature. All reactions were conducted using 30 cycles of amplification and visualized by microfluidic LabChip[®] analysis.



The stand-alone enzyme comes with a reaction buffer that supports robust amplification of high AT to routine targets. Addition of the High GC Enhancer allows amplification of GC rich and difficult targets. For added convenience, the master mix formats allow robust amplification of a broad range of targets with a single formulation.

LABCHIP[®] is a registered trademark of Caliper Life Sciences, Inc.



Phusion® High-Fidelity DNA Polymerase

RR HiFi PCR  Tm³

Phusion Hot Start Flex DNA Polymerase

RR HiFi PCR  Tm³

Phusion High-Fidelity DNA Polymerase	
#M0530S 100 units\$111
#M0530L 500 units\$446

Phusion High-Fidelity PCR Master Mix with HF Buffer	
#M0531S 100 rxns (50 µl vol)\$190
#M0531L 500 rxns (50 µl vol)\$736

Phusion High-Fidelity PCR Master Mix with GC Buffer	
#M0532S 100 rxns (50 µl vol)\$186
#M0532L 500 rxns (50 µl vol)\$751

Phusion Hot Start Flex DNA Polymerase	
#M0535S 100 units\$142
#M0535L 500 units\$566

Phusion Hot Start Flex 2X Master Mix	
#M0536S 100 rxns (50 µl vol)\$221
#M0536L 500 rxns (50 µl vol)\$880

Phusion High-Fidelity PCR Kit	
#E0535S 50 rxns (50 µl vol)\$74
#E0535L 200 rxns (50 µl vol)\$253

For a complete listing of Deoxynucleotide Solutions, see page 85.

For a complete listing of Reaction Buffers, see pages 85.

PHUSION POLYMERASE DETAILS	
Extension Rate	4 kb/min
Amplicon Size	≤ 20 kb
Fidelity	> 50X Taq
Units / 50 µl rxn	1 unit
Resulting Ends	Blunt
3'→5' Exonuclease Activity	Yes
5'→3' Exonuclease Activity	No
Supplied Buffers	- 5X Phusion HF Buffer - 5X Phusion GC Buffer
Supplied Enhancer	100% DMSO
PRODUCT FORMATS	
Hot Start Available	Yes
Activation Required	No
Master Mix Available	Yes
PCR Kit Available	Yes
APPLICATIONS	
High-Fidelity PCR	Yes
T/A, U/A Cloning	No
Colony PCR	No
Blunt Cloning	Yes

Description: DNA polymerases with high fidelity are important for applications in which the DNA sequence needs to be correct after amplification. Manufactured and quality controlled at New England Biolabs, Thermo Scientific® Phusion High-Fidelity DNA Polymerase offers both high fidelity and robust performance, and thus can be used for all PCR applications. Its unique structure, a novel *Pyrococcus*-like enzyme fused with a processivity-enhancing domain, increases fidelity and speed. Phusion Hot Start Flex DNA Polymerase is available as a standalone enzyme or in a master mix format, and enables high specificity amplification. Phusion DNA Polymerase is an ideal choice for cloning and can be used for long amplicons.

Additional Formats: Phusion and Phusion Hot Start Flex DNA Polymerases are also available in master mix format. Phusion master mixes are available with HF or GC Buffer.


The Phusion PCR Kit contains Phusion Polymerase, Phusion HF and GC buffers, deoxynucleotides, MgCl₂, DMSO and DNA size standards.

Concentration: 2,000 units/ml

* Phusion DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific. Phusion® is a registered trademark and property of Thermo Fisher Scientific.


Tm Calculator

High-Fidelity polymerases benefit from a Tm³ annealing temp. Use the **NEB Tm Calculator** to ensure successful PCR at TmCalculator.neb.com



Phusion DNA Polymerase generates robust amplification even with short extension times.

A 1.2 kb *C. elegans* genomic amplicon was analyzed using different polymerases. Reactions were set up in duplicate and visualized by microfluidic LabChip analysis. All reactions were conducted according to manufacturer's instructions using 30 cycles and extension times of either 30 seconds or 2 minutes as indicated.



Phusion Hot Start Flex DNA Polymerase delivers robust amplification. All amplicons are from human Jurkat template except for the 1.3 kb *C. elegans* amplicon. Amplicon sizes are indicated above gel. All reactions were conducted in duplicate according to manufacturer's instructions using 30 cycles and visualized after microfluidic LabChip analysis.

OneTaq® DNA Polymerase



OneTaq Hot Start DNA Polymerase



OneTaq Products

OneTaq DNA Polymerase		
#M0480S	200 units\$44
#M0480L	1,000 units\$174
#M0480X	5,000 units\$684

OneTaq 2X Master Mix with Standard Buffer		
#M0482S	100 rxns (50 µl vol)\$44
#M0482L	500 rxns (50 µl vol)\$192

OneTaq Quick-Load® DNA Polymerase		
#M0509S	100 units\$24
#M0509L	500 units\$90
#M0509X	2,500 units\$352

OneTaq Quick-Load 2X Master Mix with Standard Buffer		
#M0486S	100 rxns (50 µl vol)\$44
#M0486L	500 rxns (50 µl vol)\$192

OneTaq RT-PCR Kit		
#E5310S	30 rxns\$149

OneTaq One-step RT-PCR Kit		
#E5315S	30 rxns\$168

OneTaq Hot Start Products

OneTaq Hot Start DNA Polymerase		
#M0481S	200 units\$91
#M0481L	1,000 units\$362
#M0481X	5,000 units\$1,444

OneTaq Hot Start 2X Master Mix with Standard Buffer		
#M0484S	100 rxns (50 µl vol)\$73
#M0484L	500 rxns (50 µl vol)\$331

OneTaq Hot Start 2X Master Mix with GC Buffer		
#M0485S	100 rxns (50 µl vol)\$73
#M0485L	500 rxns (50 µl vol)\$331

OneTaq Hot Start Quick-Load 2X Master Mix with Standard Buffer		
#M0488S	100 rxns (50 µl vol)\$73
#M0488L	500 rxns (50 µl vol)\$331

OneTaq Hot Start Quick-Load 2X Master Mix with GC Buffer		
#M0489S	100 rxns (50 µl vol)\$73
#M0489L	500 rxns (50 µl vol)\$331

For a complete listing of Deoxynucleotide Solutions, see page 85.

For a complete listing of Reaction Buffers, see pages 85.

Description: OneTaq DNA Polymerase is an optimized blend of *Taq* and Deep Vent® DNA polymerases for use with routine and difficult PCR experiments. The 3'→5' exonuclease activity of Deep Vent DNA Polymerase increases the fidelity and robust amplification of *Taq* DNA Polymerase. The OneTaq reaction buffers and High GC Enhancer have been formulated for robust yields with minimal optimization, regardless of a template's GC content.

OneTaq Hot Start DNA Polymerase: OneTaq Hot Start DNA Polymerase utilizes an aptamer-based inhibitor. The hot start formulation combines convenience with decreased interference from primer-dimers and secondary products.

The aptamer-based inhibitor binds reversibly, blocking polymerase activity at temperatures below 45°C, allowing reactions to be set up at room temperature. OneTaq Hot Start DNA Polymerase is activated during normal cycling conditions, eliminating the need for a separate high temperature incubation step to activate the enzyme. OneTaq Hot Start DNA Polymerase can therefore be substituted into typical or existing *Taq*-based protocols.

OneTaq and OneTaq Hot Start are provided with Standard Reaction Buffer, GC Reaction Buffer and High GC Enhancer. Recommendations for buffer selection, based on % GC content, are shown below.

Quick-Load formats offer direct loading of PCR products onto gels, eliminating the need for a PCR clean-up step.

Additional Formats: For added convenience, OneTaq and OneTaq Hot Start DNA Polymerases are available in master mix format. Master mixes are available with Standard or GC Buffer. High GC Enhancer is also provided with master mixes containing GC Buffer. For direct gel loading, Quick-Load versions of master mixes are also available. For information on the OneTaq RT-PCR Kit or OneTaq One-step RT-PCR Kit, see page 81.

Concentration: 5,000 units/ml

ONE TAg POLYMERASE DETAILS

Extension Rate	1 kb/min
Amplicon Size	≤ 6 kb
Fidelity	> 2X <i>Taq</i>
Units / 50 µl rxn	1.25 units
Resulting Ends	3' A/Blunt
3'→5' Exonuclease Activity	Yes
5'→3' Exonuclease Activity	Yes
Supplied Buffers	- OneTaq Std Rxn Buffer - OneTaq GC Rxn Buffer
Supplied Enhancer	OneTaq High GC Enhancer

PRODUCT FORMATS


Hot Start Available	Yes
Activation Required	No
Master Mix Available	Yes
Direct Gel-loading Available	Yes
PCR Kit Available	Yes

APPLICATIONS


Routine PCR	Yes
T/A, U/A Cloning	Yes
Colony PCR	Yes

OneTaq Buffer Recommendations

AMPLICON % GC CONTENT	RECOMMENDED DEFAULT BUFFER	OPTIMIZATION NOTES
< 50% GC	OneTaq Standard Reaction Buffer	Adjust annealing temperature, primer/template concentration, etc., if needed.
50–65% GC	OneTaq Standard Reaction Buffer	OneTaq GC Reaction Buffer can be used to enhance performance of difficult amplicons.
> 65% GC	OneTaq GC Reaction Buffer	OneTaq GC Reaction Buffer with 10–20% OneTaq High GC Enhancer can be used to enhance performance of difficult amplicons.



Amplification of a selection of sequences with varying GC content from human and *C. elegans* genomic DNA using OneTaq DNA Polymerase. GC content is indicated above gel. Marker M is the 1 kb DNA Ladder (NEB #N3232).



Taq DNA Polymerase



Hot Start Taq DNA Polymerase



Taq Polymerases

Taq DNA Polymerase with ThermoPol® Buffer

#M0267S	400 units\$64
#M0267L	2,000 units\$254
#M0267X	4,000 units\$462
#M0267E	20,000 units\$1,249

Taq DNA Polymerase with Std. Taq Buffer

#M0273S	400 units\$65
#M0273L	2,000 units\$259
#M0273X	4,000 units\$469
#M0273E	20,000 units\$1,249

Taq DNA Polymerase with Std. Taq (Mg-free) Buffer

#M0320S	400 units\$65
#M0320L	2,000 units\$263

Taq PCR Kit

#E5000S	200 rxns\$108
---------	----------	------------

Taq 2X Master Mix

#M0270L	500 rxns (50 µl vol)\$146
---------	----------------------	------------

Quick-Load Taq 2X Master Mix

#M0271L	500 rxns (50 µl vol)\$238
---------	----------------------	------------

Taq 5X Master Mix

#M0285L	500 rxns (50 µl vol)\$146
---------	----------------------	------------

Multiplex PCR 5X Master Mix

#M0284S	100 rxns (50 µl vol)\$246
---------	----------------------	------------

Hot Start Taq Products

Hot Start Taq DNA Polymerase

#M0495S	200 units\$68
#M0495L	1,000 units\$271

Hot Start Taq 2X Master Mix

#M0496S	100 rxns (50 µl vol)\$60
#M0496L	500 rxns (50 µl vol)\$238

For a complete listing of Deoxynucleotide Solutions, see page 85.

For a complete listing of Reaction Buffers, see pages 85.

Description: *Taq* DNA Polymerase is a thermostable DNA polymerase that possesses a 5'→3' polymerase activity and a 5' flap endonuclease activity. It is the most widely used enzyme for PCR. To accommodate a variety of PCR applications, *Taq* is available with different reaction buffers. Standard *Taq* Reaction Buffer is designed to support existing PCR platforms, and is an ideal choice for DHPLC and high-throughput applications. ThermoPol Reaction Buffer was designed at NEB, and is formulated to promote high product yields, even under demanding conditions. For additional convenience, *Taq* DNA Polymerase is also available in kit and master mix formats. For direct gel loading, a Quick-Load version of the *Taq* 2X Master Mix is also available.

Hot Start Taq DNA Polymerase: With value pricing and attractive commercial terms, Hot Start *Taq* from NEB is an ideal choice for molecular diagnostics and other applications. In contrast to chemically modified or antibody-based hot start polymerases, NEB's Hot Start *Taq* utilizes an aptamer-based technology. The unique aptamer binds to the polymerase through non-covalent interactions, inhibiting polymerization at non-permissive temperatures. This novel method eliminates the need for an activation step, reducing the potential for sample degradation and decreasing overall reaction time.

Additional Formats: For added convenience, *Taq* and Hot Start *Taq* DNA Polymerase are available in master mix format. For direct gel loading, a Quick-Load version of the *Taq* 2X Master Mix is also available. The *Taq* PCR Kit contains *Taq* DNA Polymerase, dNTP Mix, Buffer, MgCl₂ and the Quick-Load Purple 1 kb Plus DNA Ladder.

The Multiplex PCR 5X Master Mix formulation has been specifically optimized for enhanced performance in multiplex PCR reactions.

Concentration: 5,000 units/ml

Taq Buffer Selection Chart

CHOICE OF BUFFER	AVAILABLE PRODUCTS	NEB #
ThermoPol Reaction Buffer: Designed for optimal yield and specificity	<i>Taq</i> DNA Polymerase with ThermoPol Buffer	M0267
Standard <i>Taq</i> Reaction Buffer: Detergent-free and designed to be compatible with existing assay systems	<i>Taq</i> DNA Polymerase with Standard <i>Taq</i> Buffer	M0273
	<i>Taq</i> DNA Polymerase with Standard <i>Taq</i> (Mg-free) Buffer	M0320

Taq DNA POLYMERASE DETAILS	
Extension Rate	1 kb/min
Amplicon Size	≤ 5 kb
Units / 50 µl rxn	1.25 units
Resulting Ends	3' A
3'→5' Exonuclease Activity	No
5'→3' Exonuclease Activity	Yes
Supplied Buffer (product dependent)	- Std <i>Taq</i> Rxn Buffer or - ThermoPol Rxn Buffer
PRODUCT FORMATS	
Hot Start Available	Yes
Activation Required	No
Master Mix Available	Yes
Direct Gel-loading Available	Yes
PCR Kit Available	Yes
APPLICATIONS	
Routine PCR	Yes
SNP Detection	Yes
T/A, U/A Cloning	Yes
Colony PCR	Yes

LongAmp® *Taq* DNA Polymerase



LongAmp Hot Start *Taq* DNA Polymerase



LongAmp *Taq* DNA Polymerase

#M0323S	500 units\$87
#M0323L	2,500 units\$349

LongAmp Hot Start *Taq* DNA Polymerase

#M0534S	500 units\$113
#M0534L	2,500 units\$450

LongAmp *Taq* 2X Master Mix

#M0287S	100 rxns (50 µl vol)\$138
#M0287L	500 rxns (50 µl vol)\$548

LongAmp Hot Start *Taq* 2X Master Mix

#M0533S	100 rxns (50 µl vol)\$179
#M0533L	500 rxns (50 µl vol)\$719


LongAmp *Taq* PCR Kit

#E5200S	100 rxns (50 µl vol)\$120
---------	----------------------	------------

For a complete listing of Deoxynucleotide Solutions, see page 85.

For a complete listing of Reaction Buffers, see pages 85.

LONGAMP <i>Taq</i> POLYMERASE DETAILS	
Extension Rate	1.2 kb/min
Amplicon Size	≤ 30 kb
Fidelity	2X <i>Taq</i> DNA Polymerase
Units / 50 µl rxn	5 units
Resulting Ends	3' A/Blunt
3'→5' Exonuclease Activity	Yes
5'→3' Exonuclease Activity	Yes
Supplied Buffer (product dependent)	LongAmp <i>Taq</i> Rxn Buffer
PRODUCT FORMATS	
Hot Start Available	Yes
Activation Required	No
Master Mix Available	Yes
Direct Gel-loading Available	No
PCR Kit Available	Yes
APPLICATIONS	
Long Amplicons	Yes
Routine PCR	Yes
T/A, U/A Cloning	Yes
Colony PCR	Yes



Amplification of longer templates with LongAmp *Taq*.
Amplification of specific sequences from human genomic DNA using LongAmp *Taq* DNA Polymerase. Amplicon sizes are indicated below gel. Ladder (M) is NEB 1 kb DNA Ladder (NEB #N3232).

Hemo KlenTaq®




#M0332S	200 rxns (25 µl vol)\$111
#M0332L	1,000 rxns (25 µl vol)\$441

For a complete listing of Deoxynucleotide Solutions, see page 85.

HEMO KLENTAQ DETAILS	
Extension Rate	0.5 kb/min
Amplicon Size	≤ 2 kb
Units / 50 µl rxn	4 units
Resulting Ends	3' A
3'→5' Exonuclease Activity	No
5'→3' Exonuclease Activity	No
Supplied Buffer	Hemo KlenTaq Rxn Buffer
APPLICATIONS	
Extraction-free PCR	Yes
T/A, U/A Cloning	Yes
Colony PCR	Yes

Description: Hemo KlenTaq is a truncated version of *Taq* DNA Polymerase, lacking the first 280 amino acids. Hemo KlenTaq also contains mutations that make it resistant to inhibitors present in whole blood. It can amplify DNA from whole blood samples from humans and mice, without the need for DNA extraction. Hemo KlenTaq tolerates up to 20% whole blood in a 25 µl reaction (30% in a 50 µl reaction). Hemo KlenTaq works well with most common anticoagulants, including heparin, citrate and EDTA.




Source: An *E. coli* strain that carries a mutant *Taq* DNA polymerase gene. The protein lacks the N-terminal 5' flap endonuclease domain and the gene has three internal point mutations.

Reaction Buffer: 1X Hemo KlenTaq Reaction Buffer

KLENTAQ® is a registered trademark of Wayne M. Barnes.

Amplification of human whole blood with Hemo KlenTaq. Lane 1: 5% blood + Na-EDTA; Lane 2: 5% blood + K-EDTA; Lane 3: 5% blood + Na-Heparin; Lane 4: 5% blood + Na-Citrate; Lane 5: 1.2 mm² FTA Guthrie Card containing dried human blood + Na-Heparin; Lane 6: 1.2 mm² PTA Guthrie Card containing dried human blood + Na-Heparin (washed with 50 µl H₂O at 50°C for 5 min.). Ladder (M) is the 1 kb Plus DNA Ladder (NEB #N3200).



EpiMark® Hot Start *Taq* DNA Polymerase

#M0490S 100 rxns (50 µl vol)\$66
#M0490L 500 rxns (50 µl vol)\$265

For a complete listing of Deoxynucleotide Solutions, see page 85.

EPIMARK POLYMERASE DETAILS	
Extension Rate	1 kb/min
Amplicon Size	≤ 1 kb
Units / 50 µl rxn	1.25 units
Resulting Ends	3' A
3' → 5' Exonuclease Activity	No
5' → 3' Exonuclease Activity	Yes
Supplied Buffer	EpiMark Hot Start <i>Taq</i> Rxn Buffer
APPLICATIONS	
AT-rich Targets	Yes
Bisulfite-converted DNA	Yes
Routine PCR	Yes
T/A, U/A Cloning	Yes

Description: EpiMark Hot Start *Taq* DNA Polymerase is an excellent choice for use on bisulfite-converted DNA. It is a mixture of *Taq* DNA Polymerase and a temperature-sensitive, aptamer-based inhibitor. The inhibitor binds reversibly to the enzyme, inhibiting polymerase activity at temperatures below 45°C, but releases the enzyme during normal PCR cycling conditions. This permits room temperature reaction assembly with no separate high-temperature incubation step to activate the enzyme.

Companion Product:

EpiMark Bisulfite Conversion Kit
#E3318S 48 reactions\$175



Source: An *E. coli* strain that carries the *Taq* DNA Polymerase gene from *Thermus aquaticus* YT-1.

Concentration: 5,000 units/ml

Vent® & Deep Vent® DNA Polymerases

Vent DNA Polymerase
#M0254S 200 units\$73
#M0254L 1,000 units\$294

Vent (exo⁻) DNA Polymerase
#M0257S 200 units\$73
#M0257L 1,000 units\$294

Deep Vent DNA Polymerase
#M0258S 200 units\$106
#M0258L 1,000 units\$423

Deep Vent (exo⁻) DNA Polymerase
#M0259S 200 units\$106
#M0259L 1,000 units\$423

For a complete listing of Deoxynucleotide Solutions, see page 85.

VENT/DEEP VENT POLYMERASES DETAILS	
Extension Rate	1 kb/min
Amplicon Size	≤ 6 kb
Fidelity	5-6X <i>Taq</i>
Resulting Ends	Blunt
3' → 5' Exonuclease Activity	Yes (M0254, M0258)
5' → 3' Exonuclease Activity	No
Supplied Buffer	ThermoPol Rxn Buffer

Description: Vent DNA Polymerase was the first high-fidelity thermophilic DNA polymerase available for PCR. The fidelity is 5-fold higher than that observed for *Taq* DNA Polymerase, and is derived in part from an integral 3' → 5' proofreading exonuclease activity. Greater than 90% activity remains following a 1 hour incubation at 95°C.

Deep Vent DNA Polymerase is a modified version of Vent DNA Polymerase. Deep Vent has similar fidelity with even more stability than Vent at temperatures of 95 to 100°C, with a half-life of 8 hours at 100°C.


Vent (exo⁻) DNA Polymerase has been genetically engineered to eliminate the 3' → 5' proofreading exonuclease activity associated with Vent DNA Polymerase. The fidelity of polymerization by this form is reduced to a level about 2-fold higher than that of *Taq* DNA Polymerase. Likewise, Deep Vent (exo⁻) DNA Polymerase has been genetically engineered to eliminate the 3' → 5' proofreading exonuclease activity associated with Deep Vent DNA Polymerase.



Source: Vent DNA Polymerase is purified from a strain of *E. coli* that carries the Vent DNA Polymerase gene from the archaea *Thermococcus litoralis*. Vent (exo⁻) is purified from an *E. coli* strain that carries the Vent (D141A/E143A) DNA Polymerase gene, a genetically engineered form of the native DNA polymerase.

Deep Vent DNA Polymerase is purified from a strain of *E. coli* that carries the Deep Vent DNA Polymerase gene from *Pyrococcus* species GB-D. Deep Vent (exo⁻) is purified from an *E. coli* strain that carries the Deep Vent (D141A/E143A) DNA Polymerase gene, a genetically engineered form of the native polymerase.

Concentration: 2,000 units/ml



Amplification of Jurkat genomic DNA with Vent DNA Polymerase.
Amplicon Sizes are indicated below gel. Marker (M) is the 1 kb DNA Ladder (NEB #N3232).

Luna® qPCR and RT-qPCR

- Convenient master mix and supermix formats with user-friendly protocols simplify reaction setup
- Non-interfering, visible tracking dye helps to eliminate pipetting errors
- Novel, thermostable reverse transcriptase (RT) improves performance
- One-Step RT-qPCR kits feature Luna WarmStart RT paired with Hot Start Taq for increased reaction specificity and robustness

Fluorescence-based quantitative real-time PCR (qPCR) is the gold standard for the detection and quantification of nucleic acids due to its sensitivity and specificity. Luna products from NEB are optimized for qPCR or RT-qPCR, and are available for either intercalating dye or probe-based detection methods. Each Hot Start Taq-based Luna qPCR master mix has been formulated with a unique passive reference dye that is compatible across a wide variety of thermal cyclers, regardless of ROX requirements. No additional components are required to ensure compatibility. For two-step RT-qPCR, the LunaScript® RT SuperMix Kit offers a fast (less than 15 minute), robust, and sensitive option for cDNA synthesis upstream of your Luna qPCR experiment.

The Luna Probe One-Step RT-qPCR Mix with UDG is supplied at a 4X concentration and enables higher amounts of sample input, which is relevant for applications where increased sensitivity is needed, such as molecular diagnostics. Performance in multiplexing applications has been optimized, with linear detection achieved for up to 5 targets across a range of inputs.

Find the right Luna product for your application


1 Select your target



2 Select your detection method

	Dye-based	Probe-based
Genomic DNA or cDNA	Luna® Universal qPCR Master Mix (NEB #M3003)	Luna Universal Probe qPCR Master Mix (NEB #M3004)
Purified RNA	Luna Universal One-Step RT-qPCR Kit (NEB #E3005)	Luna Universal Probe One-Step RT-qPCR: • Kit (NEB #E3006) • Kit (no ROX) (NEB #E3007) • 4X Mix with UDG (NEB #M3019)
Two-Step RT-qPCR	LunaScript® RT SuperMix Kit (NEB #E3010) + Luna Universal qPCR Master Mix (NEB #M3003)	LunaScript RT SuperMix Kit (NEB #E3010) + Luna Universal Probe qPCR Master Mix (NEB #M3004)
RNA from cell lysate	Luna Cell Ready One-Step RT-qPCR Kit (NEB #E3030)	Luna Cell Ready Probe One-Step RT-qPCR Kit (NEB #E3031)


Avoiding pipetting errors with Luna



Find an overview of qPCR.



Cloned at NEB



Recombinant Enzyme



High-Fidelity



PCR Enzyme



Hot Start/WarmStart



Requires BSA



Heat Inactivation



Annealing Temperature

Luna Universal qPCR & Probe qPCR Master Mixes



Luna Universal qPCR Master Mix

#M3003S	200 reactions\$120
#M3003L	500 reactions\$270
#M3003X	1,000 reactions\$469
#M3003E	2,500 reactions\$1,040

Luna Universal Probe qPCR Master Mix

#M3004S	200 reactions\$104
#M3004L	500 reactions\$232
#M3004X	1,000 reactions\$403
#M3004E	2,500 reactions\$895

Description: The Luna Universal qPCR Master Mix is an optimized 2X reaction mix for real-time qPCR detection and quantitation of target DNA sequences using the SYBR®/FAM channel of most real-time qPCR instruments.

The Luna Universal Probe qPCR Master Mix is a 2X reaction mix optimized for real-time qPCR detection and quantitation of target DNA sequences using hydrolysis probes.

Each Hot Start Taq-based Luna qPCR master mix has been formulated with a unique passive reference dye

that is compatible across a wide variety of instrument platforms, including those that require a ROX reference signal. This means that no additional components are required to ensure machine compatibility. The mixes also contain dUTP, enabling carryover prevention when reactions are treated with Antarctic Thermolabile UDG. A blue visible dye assists in tracking the reagents when pipetting into clear, multi-welled PCR plates. These features, combined with rapid, sensitive and precise real-time qPCR performance, make Luna the universal choice for all your qPCR experiments.

SYBR® is a registered trademark of Thermo Fisher Scientific.


Companion Product:

Antarctic Thermolabile UDG

#M0372S	100 units\$78
#M0372L	500 units\$316

- Convenient master mix formats and user-friendly protocols simplify reaction setup
- Non-interfering, visible tracking dye helps to eliminate pipetting errors
- Rigorously tested to optimize specificity, sensitivity, accuracy and reproducibility
- Unique passive reference dye for compatibility across wide range of instruments

Learn more about our comprehensive qPCR/RT-qPCR testing and "dots in boxes" data visualization at LUNAqPCR.com.



Extensive performance evaluation of commercially available dye-based qPCR reagents demonstrates the robustness and specificity of Luna. qPCR reagents from NEB and other manufacturers were tested across 16–18 qPCR targets varying in abundance, length and %GC, using either Jurkat genomic DNA or Jurkat-derived cDNA as input (10 genomic DNA targets and 8 cDNA targets on Bio-Rad real-time instrument, 9 genomic and 7 cDNA targets on ABI instrument). For each testing condition, data was collected by 2 users and according to manufacturer's specifications. Results were evaluated for efficiency, low input detection and lack of non-template amplification (where ΔCq = average Cq of non-template control – average Cq of lowest input). In addition, consistency, reproducibility and overall curve quality were assessed (Quality Score). Bar graph indicates % of targets that met acceptable performance criteria (indicated by green box on dot plot and Quality Score > 3). NEB's Luna Universal qPCR Master Mix outperformed all other reagents tested.

LunaScript® RT SuperMix Kit & LunaScript RT Master Mix Kit (Primer-free)



LunaScript RT SuperMix Kit

#E3010S	25 reactions\$141
#E3010L	100 reactions\$437

NEW

LunaScript RT Master Mix Kit (Primer-free)

#E3025S	25 reactions\$117
#E3025L	100 reactions\$364

- Less than 15 minute first-strand cDNA synthesis protocol
- Combine with Luna qPCR master mixes for robust RT-qPCR results

Description: The LunaScript RT SuperMix Kit is an optimized master mix containing all the necessary components for first strand cDNA synthesis in the context of a two-step RT-qPCR workflow. It features the thermostable Luna Reverse Transcriptase, which supports cDNA synthesis at elevated temperatures. Murine RNase inhibitor is also included to protect template RNA from degradation. The LunaScript RT SuperMix Kit is available with or without random hexamer and poly-dT primers, which allow for even coverage across the length of the RNA targets.

LunaScript RT Master Mix Kit (Primer-free) is the same formulation as the LunaScript RT SuperMix Kit but does not contain primers. It includes all components needed for carrying out first strand cDNA synthesis except for RNA template and user-supplied primers. It provides flexible options for using different primers (dT primers, random

primers, gene-specific primers, modified primers, etc.) for first strand cDNA synthesis.

The LunaScript RT SuperMix Kit Includes:

- LunaScript RT SuperMix
- No-RT Control Mix
- Nuclease-free Water

The LunaScript RT Master Mix Kit (Primer-free) Includes:

- LunaScript RT Master Mix (Primer-free)
- No-RT Control Mix (Primer-free)
- Nuclease-free Water

How can we ensure best in class performance with Luna?



Luna One-Step RT-qPCR Products



NEW

Luna Probe One-Step RT-qPCR Kit (No ROX)
#E3007E 2,500 reactions\$1,908

NEW

Luna Probe One-Step RT-qPCR 4X Mix
with UDG
#M3019S 200 reactions\$310
#M3019L 500 reactions\$698
#M3019X 1,000 reactions\$1,241
#M3019E 2,000 reactions\$2,190

Luna Universal One-Step RT-qPCR Kit
#E3005S 200 reactions\$244
#E3005L 500 reactions\$549
#E3005X 1,000 reactions\$952
#E3005E 2,500 reactions\$2,116

Luna Universal Probe One-Step RT-qPCR Kit
#E3006S 200 reactions\$220
#E3006L 500 reactions\$495
#E3006X 1,000 reactions\$867
#E3006E 2,500 reactions\$1,908

Companion Product:

Antarctic Thermolabile UDG
#M0372S 100 units\$78
#M0372L 500 units\$316

- Novel, thermostable reverse transcriptase (RT) improves performance
- WarmStart RT paired with Hot Start Taq increases reaction specificity and robustness
- Non-interfering, visible tracking dye helps to eliminate pipetting errors
- Products perform consistently across a wide variety of sample sources

Description: The Luna RT-qPCR kits contain a novel, *in silico*-designed reverse transcriptase (RT) engineered for improved performance. Both the Luna WarmStart Reverse Transcriptase and Hot Start Taq DNA Polymerase, included in these kits, utilize a temperature-sensitive, reversible aptamer, which inhibits activity below 45°C. This enables room temperature reaction setup and prevents undesired non-specific activity. Furthermore, the WarmStart RT has increased thermostability, improving performance at higher reaction temperatures.

The Luna Probe One-Step RT-qPCR 4X Mix with UDG (NEB #M3019) supports robust, sensitive detection and quantitation of up to 5 targets in a multiplexed reaction. It is supplied at a 4X concentration and enables higher amounts of sample input, which is relevant for applications where RNA present in low abundance is of interest, such as pathogen detection. For simplified reaction setup, the single tube master mix format consolidates components for the one-step RT-qPCR reaction. It also includes dUTP and UDG in the mix for reduced risk of carryover contamination.

The NEB Luna Universal One-Step RT-qPCR Kit is optimized for dye-based real-time quantitation of target RNA sequences via the SYBR/FAM fluorescence channel of most real-time instruments.

The NEB Luna Universal Probe One-Step RT-qPCR Kit is optimized for real-time quantitation of target RNA sequences using hydrolysis probes.

For instruments that do not utilize ROX normalization, the Luna Probe One-Step RT-qPCR Kit (No ROX) (NEB #E3007) contains no reference dye. If ROX normalization is desired, ROX can be added; this is only necessary with the E3007 product.

The other Luna products contain dUTP and enable carryover prevention when reactions are treated with Antarctic Thermolabile UDG (NEB #M0372). A blue visible dye assists in tracking the reagents when pipetting into clear, multi-welled PCR plates. The reverse transcriptase, featured in the Luna One-Step RT-qPCR products is a novel, engineered WarmStart enzyme developed for robust performance and increased thermostability. These features, combined with rapid, sensitive and precise real-time qPCR performance, make Luna the universal choice for all your qPCR and RT-qPCR experiments.

The Luna Probe One-Step RT-qPCR Kit (No ROX) Includes:

- Luna Universal Probe One-Step Reaction Mix (No ROX)
- Luna WarmStart RT Enzyme Mix
- Nuclease-free Water

The Luna Universal One-Step RT-qPCR Kit Includes:

- Luna Universal One-Step Reaction Mix
- Luna WarmStart RT Enzyme Mix
- Nuclease-free Water

The Luna Universal Probe One-Step RT-qPCR Kit Includes:

- Luna Universal Probe One-Step Reaction Mix
- Luna WarmStart RT Enzyme Mix
- Nuclease-free Water

	Luna Universal Probe One-Step RT-qPCR Kit (NEB #E3006)		Luna Probe RT-qPCR 4X Mix with UDG (NEB #M3019)
	20X RT Enzyme Mix	2X Reaction Mix	4X Master Mix
Number of tubes	2		1
Sample-limiting concentration	2X		4X
dUTP included	✓		✓
UDG included	✗		✓
Universal ROX included	✓		✓

A new product choice for one-step RT-qPCR assays

Luna Cell Ready One-Step RT-qPCR Kit



Luna Cell Ready Probe One-Step RT-qPCR Kit

NEW

Luna Cell Ready One-Step RT-qPCR Kit
#E3030S 100 reactions\$990

NEW

Luna Cell Ready Probe One-Step RT-qPCR Kit
#E3031S 100 reactions\$939

NEW

Luna Cell Ready Lysis Module
#E3032S 100 rxns (50 µl vol)\$505

Companion Product:

Antarctic Thermolabile UDG
#M0372S 100 units\$78
#M0372L 500 units\$316

- Go direct from cells to RNA quantitation without purification
- Coordinated cell lysis, RNA release, and genomic DNA removal in a fast 15-minute protocol
- Effective cell lysis preparation from 10 to 100,000 cells across numerous cell lines
- Features Luna WarmStart RT paired with Hot Start Taq for increased thermostability and room temperature setup

Description: The Luna Cell Ready One-Step RT-qPCR Kit provides all the necessary components for direct, dye-based RNA detection and quantitation, bypassing the need for RNA extraction and purification.

The Luna Cell Ready Probe One-Step RT-qPCR Kit provides all the necessary components for direct, probe-based RNA detection and quantitation, bypassing the need for RNA extraction and purification.


Cell cultures are often analyzed for gene expression or treatment responses as a proxy for a living organism. Traditionally, RNA is extracted and purified from treated cells via column-based or chemical methods. Coordinating the actions of DNase I and the Luna Cell Ready Protease, the Luna Cell Ready Lysis Module offers a simple alternative workflow resulting in effective cell lysis, RNA release, and genomic DNA removal simultaneously in a 15-minute protocol. The Lysis Module includes a unique Luna Cell Ready RNA Protection Reagent that maintains RNA integrity during cell lysis. The lysis capacity spans 10–100,000 cells in a 50 µl lysis reaction. Up to 2 µl of lysate (equivalent to RNA from 0.2–4,000 cells) can be transferred into 20 µl downstream RT-qPCR reactions. Similar to other Luna products, the lysis buffer includes an inert blue tracking dye for visual assistance throughout the workflow.

The Luna Cell Ready One-Step RT-qPCR Kit Includes:

- Luna Cell Ready Lysis Module
- Luna Universal One-Step RT-qPCR Kit (NEB #E3005)


The Luna Cell Ready Probe One-Step RT-qPCR Kit Includes:

- Luna Cell Ready Lysis Module
- Luna Universal Probe One-Step RT-qPCR Kit (NEB #E3006)



The Luna Cell Ready One-Step RT-qPCR Kit provides all the necessary components for direct RNA detection and quantitation from cultured cells (up to 100,000 cells per 50 µl lysis reaction). By coordinating the actions of DNase I and the Luna Cell Ready Protease, the Luna Cell Ready Lysis Module offers a simple workflow resulting in effective cell lysis, RNA release and genomic DNA removal simultaneously in a 15-minute protocol. Up to 2 µl lysate (equivalent to RNA from 0.2 - 4,000 cells) can be transferred into 20 µl downstream RT-qPCR reactions.

SARS-CoV-2 Detection



Virus detection often utilizes nucleic acid amplification to identify the presence of specific sequences in SARS-CoV-2 viral RNA. These methods include RT-qPCR for real time quantitation and loop-mediated isothermal amplification (LAMP) for robust amplification performed at a single reaction temperature.

NEW

SARS-CoV-2 Rapid Colorimetric LAMP Assay Kit



#E2019S 96 reactions\$750

Companion Product:


SARS-CoV-2 Positive Control (N gene)
#N2117S 0.05 ml\$70

- *Colorimetric LAMP enables simple, visual detection (pink-to-yellow) of amplification of SARS-CoV-2 nucleic acid*
- *Set up reactions quickly and easily, using a simple heat source and unique WarmStart technology*
- *Reduce risk of carryover contamination, with UDG and dUTP included in the master mix*
- *Assay targets N and E regions of the SARS-CoV-2 genome, for optimized sensitivity and specificity*
- *Bring confidence to your results using the provided controls*

Description: The SARS-CoV-2 Rapid Colorimetric LAMP Assay Kit utilizes Loop-Mediated Isothermal Amplification (LAMP) to detect SARS-CoV-2 nucleic acid. The kit is available for research use only and includes WarmStart Colorimetric LAMP 2X Master Mix with UDG and a primer mix targeting the N and E regions of the viral genome. Controls are provided to verify assay performance, and include an internal control primer set and a positive control template. Guanidine hydrochloride has been found to increase the speed and sensitivity of the RT-LAMP reaction and is also included.

The SARS-CoV-2 Rapid Colorimetric LAMP Assay Kit Includes:

- WarmStart Colorimetric LAMP 2X Master Mix with UDG
- SARS-CoV-2 Positive Control (N gene)
- Internal Control LAMP Primer Mix
- SARS-CoV-2 LAMP Primer Mix (N/E)
- Nuclease-free Water
- Guanidine Hydrochloride




In the SARS-CoV-2 Colorimetric LAMP Assay Kit, three control reactions are run with each test sample. All reactions should be pink prior to incubation. Using the provided protocol, a Non-Template Control (NTC) reaction will contain all materials of the test sample (master mix, primers, etc) except for the test input nucleic acid itself and serves as a measure of reaction contamination and primer-based mis-amplification. The NTC sample should not amplify and should stay pink throughout the experiment. The Positive Control (PC) will contain master mix, a plasmid that contains the SARS-CoV-2 N-gene (GenBank: MN908947.3) and primers that will amplify this sequence. Amplification should be observed in the PC and the sample should become yellow after incubation. The Internal Control (IC) will contain master mix, test input nucleic acid, and LAMP primers for *rActin*, an endogenous housekeeping gene. If reagents are active and samples have been handled appropriately, the IC should become yellow after incubation.



View our loop mediated isothermal amplification tutorial.



Cloned at NEB



Recombinant Enzyme



High-Fidelity



PCR Enzyme



Hot Start/WarmStart



BSA Requires BSA



Heat Inactivation



Annealing Temperature

NEW

WarmStart Fluorescent LAMP/RT-LAMP Kit (with UDG)



#E1708S	100 reactions\$273
#E1708L	500 reactions\$1,092

Companion Products:

LAMP Fluorescent Dye	
#B1700S	0.25 ml\$45

WarmStart Multi-Purpose LAMP/RT-LAMP 2X Master Mix (with UDG)	
#M1708S	100 reactions\$270
#M1708L	500 reactions\$1,080

Description: The WarmStart Fluorescent LAMP/RT-LAMP Kit (with UDG) is designed to provide a simple, one-step solution for Loop-Mediated Isothermal Amplification (LAMP) of DNA or RNA (RT-LAMP) targets. This kit is supplied with the WarmStart Multi-Purpose LAMP/RT-LAMP 2X Master Mix (with UDG), which contains a blend of *Bst* 2.0 WarmStart DNA Polymerase, WarmStart RTx Reverse Transcriptase, Antarctic Thermolabile UDG and dUTP in an optimized LAMP buffer solution. Both *Bst* 2.0 WarmStart DNA Polymerase and WarmStart RTx Reverse Transcriptase have been engineered for improved performance in LAMP and RT-LAMP reactions. The inclusion of dUTP and UDG in the master mix reduces

the possibility of carryover contamination between reactions. A fluorescent dye is also supplied to enable real-time fluorescence measurement of LAMP. The WarmStart Fluorescent LAMP/RT-LAMP Kit (with UDG) is compatible with multiple detection methods, including turbidity detection, real-time fluorescence detection (when used with LAMP fluorescent dye) and end-point visualization.

The WarmStart Fluorescent LAMP/RT-LAMP Kit (with UDG) Includes:

- WarmStart Multi-Purpose LAMP/RT-LAMP 2X Master Mix (with UDG)
- LAMP Fluorescent Dye

- LAMP amplification of DNA or RNA targets
- Reduce risk of carryover contamination with UDG and dUTP included in the master mix
- Improve LAMP specificity and sensitivity with optimized master mixes
- Set up reactions at room temperature with our unique dual WarmStart formulation
- WarmStart Multi-Purpose LAMP/RT-LAMP 2X Master Mix (with UDG) (NEB #M1708) supports a variety of detection methods including turbidity, visual detection, and electrophoresis

NEW

Luna SARS-CoV-2 RT-qPCR Multiplex Assay Kit



#E3019S	96 reactions\$295
#E3019L	480 reactions\$1,288

Description: The Luna SARS-CoV-2 RT-qPCR Multiplex Assay Kit is optimized for real-time qualitative detection of SARS-CoV-2 nucleic acid using hydrolysis probes. It features the Luna Probe One-Step RT-qPCR 4X Mix with UDG, an optimized SARS-CoV-2 Primer/Probe mix containing primers and probes specific to two regions of the SARS-CoV-2 virus N-gene, and a positive control template. The probes have been modified to contain different fluorophores (N1, HEX; N2, FAM) to enable simultaneous observation on two different channels of a real-time instrument. To ensure the integrity of the input material and absence of inhibition, an internal control (IC) primer and probe set, designed to amplify the human RNase P gene, is also provided in the primer mix.

The reverse primer of this target has been modified from the CDC design to target an exon/exon boundary to reduce background amplification from possible contaminating genomic DNA. Amplification of the IC is observed in the Cy5 channel.

The Luna SARS-CoV-2 RT-qPCR Multiplex Assay Kit Includes:

- Luna Probe One-Step RT-qPCR 4X Mix with UDG (NEB #M3019)
- SARS-CoV-2 Positive Control (N gene)
- SARS-CoV-2 Primer/Probe Mix (N1/N2/RP)
- Nuclease-free Water

- Multiplex detection of N1 and N2 targets
- Enables sample pooling of purified RNA
- Internal controls include a redesigned RNase P reverse primer for reduced background amplification

For the isolation of viral RNA, try our Monarch RNA purification kits. See pages 137-138 for details.



WarmStart® Colorimetric LAMP/RT-LAMP 2X Master Mix



NEW

WarmStart Colorimetric LAMP/RT-LAMP 2X Master Mix (with UDG)



WarmStart Colorimetric LAMP/RT-LAMP 2X Master Mix

#M1800S	100 reactions\$213
#M1800L	500 reactions\$849


WarmStart Colorimetric LAMP/RT-LAMP 2X Master Mix (with UDG)

#M1804S	100 reactions\$270
#M1804L	500 reactions\$1,080

- Simple, one-step solution for Loop-Mediated Isothermal Amplification (LAMP) of DNA or RNA (RT-LAMP) targets
- Fast, clear, pink-to-yellow visible detection of amplification
- WarmStart feature inhibits enzyme activity at room temperature allowing for flexible reaction setup
- Reduce risk of carryover contamination with UDG and dUTP included in the master mix (NEB #M1804)

Description: The WarmStart Colorimetric LAMP/RT-LAMP 2X Master Mix is an optimized formulation of *Bst* 2.0 WarmStart DNA Polymerase and WarmStart RTx in a special low-buffer reaction solution containing a visible pH indicator for rapid and easy detection of Loop Mediated Isothermal Amplification (LAMP) and RT-LAMP reactions. WarmStart Colorimetric LAMP/RT-LAMP 2X Master Mix (with UDG) contains dUTP and UDG in the master mix, which reduces the possibility of carryover contamination between reactions.

This system is designed to provide a fast, clear visual detection of amplification based on the production of protons and subsequent drop in pH that occurs from the extensive DNA polymerase activity in a LAMP reaction, producing a change in solution color from pink to yellow (an overview of LAMP and primer design can be found in the Featured Videos section). This mix can be used for any LAMP or RT-LAMP reaction and requires only a heated chamber and samples, with readout of positive amplification judged by eye in 15–40 minutes.



WarmStart Fluorescent LAMP/RT-LAMP Kit



#E1700S	100 reactions\$216
#E1700L	500 reactions\$861

Companion Product:

LAMP Fluorescent Dye		
#B1700S	0.25 ml\$45
<i>Tte</i> UvrD Helicase		
#M1202S	50 reactions\$73

- LAMP amplification of RNA and DNA targets
- Improve LAMP specificity and sensitivity with optimized master mixes
- Set up reactions at room temperature with our unique dual WarmStart formulation
- Use with a variety of detection methods including fluorescence, turbidity, visual detection and electrophoresis


Description: The WarmStart Fluorescent LAMP/RT-LAMP Kit is designed to provide a simple, one-step solution for Loop-Mediated Isothermal Amplification (LAMP) of DNA or RNA (RT-LAMP) targets. LAMP and RT-LAMP are commonly used isothermal amplification techniques that provides rapid detection of a target nucleic acid using LAMP-specific primers (supplied by the user) and a strand-displacing DNA polymerase. This kit is supplied with the WarmStart LAMP 2X Master Mix, which contains a blend of *Bst* 2.0 WarmStart DNA Polymerase and WarmStart RTx Reverse Transcriptase in an optimized LAMP buffer solution. Both *Bst* 2.0 WarmStart DNA Polymerase and WarmStart RTx Reverse Transcriptase have been engineered for improved performance in LAMP and RT-LAMP reactions. A fluorescent dye is also supplied to enable real-time fluorescence measurement of LAMP. The WarmStart Fluorescent LAMP/RT-LAMP Kit is compatible with multiple detection methods, including turbidity detection, real-time fluorescence detection (when used with LAMP fluorescent dye) and end-point visualization.

The WarmStart Fluorescent LAMP/RT-LAMP Kit Includes:

- WarmStart LAMP 2X Master Mix
- LAMP Fluorescent Dye (50X)



How is colorimetric LAMP used in point of care?



Bst DNA Polymerase-based Products for Isothermal DNA Amplification

	5'→3' EXO ACTIVITY	AMPLIFICATION SPEED	ROOM TEMPERATURE SETUP	REVERSE TRANSCRIPTASE ACTIVITY	INHIBITOR TOLERANCE	APPLICATIONS
Bst DNA Polymerase, Full Length	★★	N/A	N/A	N/A	★	<ul style="list-style-type: none"> Nick translation reactions at elevated temperatures Primer extension
Bst DNA Polymerase, Large Fragment	N/A	★	N/A	★	★	<ul style="list-style-type: none"> General strand-displacement reactions
Bst 2.0 DNA Polymerase	N/A	★★	N/A	★★	★	<ul style="list-style-type: none"> Improved LAMP, SDA, and other amplification reactions Minimal effect of substitution of dTTP with dUTP
Bst 2.0 WarmStart DNA Polymerase	N/A	★★	★★★	★★	★★	<ul style="list-style-type: none"> Consistent, room-temperature, and high-throughput amplification assays Minimal effect of substitution of dTTP with dUTP
Bst 3.0 DNA Polymerase	N/A	★★★	★★	★★★	★★★	<ul style="list-style-type: none"> Fused to novel nucleic acid binding domain for enhanced performance Fastest, most robust LAMP and RT-LAMP reactions High reverse transcriptase activity up to 72°C Strand displacement DNA synthesis

★★★ Optimal, recommended product for selected application
 ★★ Works well for selected application

★ Will perform selected application, but is not recommended
 N/A Not applicable to this application

Bst DNA Polymerases



Bst DNA Polymerase, Large Fragment

#M0275S 1,600 units \$72
 #M0275L 8,000 units \$289
 for high (15X) concentration
 #M0275M 8,000 units \$289

Bst DNA Polymerase, Full Length

#M0328S 500 units \$72

Bst 2.0 DNA Polymerase

#M0537S 1,600 units \$72
 #M0537L 8,000 units \$289
 for high (15X) concentration
 #M0537M 8,000 units \$289

Bst 2.0 WarmStart DNA Polymerase

#M0538S 1,600 units \$81
 #M0538L 8,000 units \$319
 for high (15X) concentration
 #M0538M 8,000 units \$319

Bst 3.0 DNA Polymerase

#M0374S 1,600 units \$71
 #M0374L 8,000 units \$284
 for high (15X) concentration
 #M0374M 8,000 units \$284

Companion Products:

WarmStart RTx Reverse Transcriptase
 #M0380S 50 reactions \$69
 #M0380L 250 reactions \$277

Tte UvrD Helicase

#M1202S 0.5 µg \$73

For a complete listing of Deoxynucleotide Solutions, see page 85.

Description: *Bst* DNA Polymerase, Large Fragment, is the portion of the *Bacillus stearothermophilus* DNA Polymerase protein that contains the 5'→3' polymerase activity, but lacks 5'→3' exonuclease activity.

Bst DNA Polymerase, Full Length is the full length polymerase from *Bacillus stearothermophilus*. It has 5'→3' polymerase and double-strand specific 5'→3' exonuclease activities, but lacks 3'→5' exonuclease activity.

Bst 2.0 DNA Polymerase is an *in silico* designed homologue of *Bst* DNA Polymerase I, Large Fragment. It contains 5'→3' DNA polymerase activity and strong strand displacement activity but lacks 5'→3' exonuclease activity. It has improved amplification speed, yield, salt tolerance and thermostability compared to wild-type *Bst* DNA Polymerase, Large Fragment.

Bst 2.0 WarmStart DNA Polymerase utilizes aptamer technology to inhibit activity at non-permissive temperatures (< 50°C). Like "Hot Start" PCR polymerases, the WarmStart feature enables room temperature set up and prevents non-templated addition of nucleotides, increasing reaction efficiencies. Additionally, no separate activation step is required to release the aptamer from the enzyme. *Bst* 2.0 WarmStart DNA Polymerase permits reaction temperatures from 60–72°C.


Bst 3.0 DNA Polymerase is a similarly designed *in silico* homologue engineered and fused to a novel nucleic acid binding domain for improved isothermal amplification performance and increased reverse transcription activity. *Bst* 3.0 DNA Polymerase contains 5'→3' DNA polymerase activity with either DNA or RNA templates but lacks 5'→3' and 3'→5' exonuclease activity. It demonstrates robust performance in the presence of inhibitors and significantly increase reverse transcriptase activity compared to *Bst* DNA Polymerase.

Concentration: *Bst* DNA Polymerase, Full Length: 5,000 units/ml. All others: 8,000 and 120,000 units/ml

Heat Inactivation: 80°C for 20 minutes

Usage Notes: No *Bst* DNA Polymerase-based products can be used for thermal cycle sequencing or PCR.

Bst 2.0 WarmStart DNA Polymerase permits reaction temperatures from 60–72°C. Generally, reaction temperatures above 72°C are not recommended for any *Bst* DNA Polymerase-based product.



Benefits of Bst 2.0 WarmStart: Identical LAMP reactions were run either immediately after setup (solid line) or after a 2 hour incubation at 25°C. Without the protection from Bst 2.0 WarmStart, this room temperature incubation results in variable LAMP performance. Bst 2.0 WarmStart provides more consistent amplification reaction and enables room-temperature and high-throughput setup.

IsoAmp® II Universal tHDA Kit

#H0110S 50 reactions\$477


- Easy-to-use for assay development
- Helicase eliminates need for thermocycler
- Reactions performed at constant temp
- Amplify & detect short DNA sequences (70–120 bp)
- Use with a variety of templates (microbial genomic DNA, viral DNA, plasmid DNA and cDNA)
- Amplify a single copy of target DNA by tHDA when optimized primers and buffer are used

Description: Thermophilic Helicase-Dependent Amplification (tHDA) is a novel method for isothermal amplification of nucleic acids. Like PCR, the tHDA reaction selectively amplifies a target sequence defined by two primers. However, unlike PCR, tHDA uses an enzyme called a helicase to separate DNA, rather than heat. This allows DNA amplification without the need for thermocycling. The tHDA reaction can also be coupled with reverse transcription for RNA analysis.

IsoAmp II Universal tHDA Kit is based on a second-generation thermophilic Helicase-Dependent Amplification platform. The reactions supported by IsoAmp II Universal tHDA Kit include tHDA, reverse transcription HDA (RT-HDA), real-time quantitative HDA (qHDA) and real-time quantitative RT-HDA (qRT-HDA), from a single reaction buffer.

The IsoAmp II Universal tHDA Kit Includes:

- IsoAmp dNTP solution and IsoAmp Enzyme Mix
- 10X Annealing Buffer II, 100 mM MgSO₄ and 500 mM NaCl
- Control template and amplification primers



HDA technology. Helicase Dependent Amplification:
Step 1: Helicase unwinding and primer binding. Step 2: DNA polymerization. Step 3: DNA amplification.

Developed by BioHelix Corporation a NEB-affiliated company, now part of Quidel Corporation. ISOAMP® is a registered trademark of BioHelix Corporation

phi29 DNA Polymerase

#M0269S 250 units\$59
#M0269L 1,250 units\$233

For a complete listing of Deoxynucleotide Solutions, see page 85.

- Extreme processivity
- Extreme strand displacement
- Replication requiring a high degree of strand displacement and/or processive synthesis
- High-fidelity replication at moderate temperatures


Description: phi29 DNA Polymerase is the replicative polymerase from the *Bacillus subtilis* phage phi29 (φ29). This polymerase has exceptional strand displacement and processive synthesis properties. The polymerase has an inherent 3' → 5' proofreading exonuclease activity.



Concentration: 10,000 units/ml

Heat Inactivation: 65°C for 10 minutes

As a member of our Marketing Team, Michelle provides support to our various ordering channels. This can include our Freezer Programs, as well as other third-party platforms used at various academic institutions. Prior to taking this role, Michelle was a member of our Technical Support Group.



OneTaq One-Step RT-PCR Kit



#E5315S 30 reactions\$168

Companion Products:

ProtoScript II First Strand cDNA Synthesis Kit	
#E6560S	30 reactions\$171
#E6560L	150 reactions\$685
ProtoScript II Reverse Transcriptase	
#M0368S	4,000 units\$92
#M0368L	10,000 units\$184
#M0368X	40,000 units\$656
RNase Inhibitor, Murine	
#M0314S	3,000 units\$73
#M0314L	15,000 units\$292

For a complete listing of OneTaq products, see page 68.

- Combine cDNA synthesis and PCR in a single reaction
- Detect at little as 0.1 pg of a GAPH target
- Robust amplification from 100 bp to 9 kb
- Faster protocols with less hands-on time
- Quick-Load Reaction Mix allows instant gel loading

Description: The OneTaq One-Step RT-PCR Kit offers sensitive and robust end-point detection of RNA templates. cDNA synthesis and PCR amplification steps are performed in a single reaction using gene-specific primers, resulting in a streamlined RT-PCR protocol.


The kit combines three optimized mixes: OneTaq One-Step Enzyme Mix, OneTaq One-Step Reaction Mix and OneTaq One-Step Quick-Load Reaction Mix. The Enzyme Mix combines ProtoScript II Reverse Transcriptase, Murine RNase Inhibitor and OneTaq Hot Start DNA Polymerase. ProtoScript II Reverse Transcriptase is a mutant M-MuLV reverse transcriptase with reduced RNase H activity and increased thermostability. OneTaq Hot Start DNA Polymerase is mixture of a Hot Start Taq DNA Polymerase combined with a proof-reading DNA polymerase, resulting in high-yield amplification with minimal optimization. The OneTaq One-Step RT-PCR Kit is capable of amplifying long transcripts up to 9 kb in length.

Two optimized reaction mixes are included, OneTaq One-Step Reaction Mix and Quick-Load OneTaq One-Step Reaction Mix. The reaction mixes offer robust conditions for both cDNA synthesis and PCR amplification. The unique Quick-Load OneTaq One-Step Reaction Mix contains additional dyes, offering color indication for reaction setup as well as direct gel loading.

Both total RNA and mRNA can be used as template. The kit can detect a GAPDH target as low as 0.1 pg per reaction. It can routinely detect RNA targets up to 9 kb. The OneTaq One-Step RT-PCR Kit is capable of multiplex detection of two or three targets.

The OneTaq One-Step RT-PCR Kit Includes:

- OneTaq One-Step Enzyme Mix
- OneTaq One-Step Reaction Mix
- Quick-Load OneTaq One-Step Reaction Mix
- Nuclease-free Water



Detection of RNA templates of different length. About 100 ng of Jurkat total RNA was used in 50 µl reactions following the standard protocol. The target sizes were Lane 1: 0.7 kb, Lane 2: 1.1 kb, Lane 3: 1.9 kb, Lane 4: 2.3 kb, Lane 5: 2.5 kb, Lane 6: 5.5 kb, Lane 7: 7.6 kb and Lane 8: 9.3 kb. The marker lane (M) contains Quick-Load 1 kb Plus DNA Ladder (NEB #N0469).

OneTaq RT-PCR Kit



#E5310S 30 reactions\$149

Companion Products:

RNase Inhibitor, Murine	
#M0314S	3,000 units\$73
#M0314L	15,000 units\$292
M-MuLV Reverse Transcriptase	
#M0253S	10,000 units\$73
#M0253L	50,000 units\$294
OneTaq Hot Start 2X Master Mix with Standard Buffer	
#M0484S	100 rxns (50 µl vol)\$73
#M0484L	500 rxns (50 µl vol)\$331

For a complete listing of OneTaq products, see page 68.


- Robust RT-PCR reactions
- Convenient master mix formats
- OneTaq DNA Polymerase offers robust amplification for a wide range of templates

Description: OneTaq RT-PCR Kit combines two powerful mixes, M-MuLV Enzyme Mix and OneTaq Hot Start 2X Master Mix with Standard Buffer for 2-step RT-PCR applications. The two mixes require minimal handling during reaction setup and yet offer consistent and robust RT-PCR reactions.

The first strand cDNA synthesis is achieved by using two optimized mixes, M-MuLV Enzyme Mix and M-MuLV Reaction Mix. M-MuLV Enzyme Mix combines M-MuLV Reverse Transcriptase and RNase Inhibitor, Murine while M-MuLV Reaction Mix contains dNTPs and an optimized buffer. The kit also contains two optimized primers for reverse transcription and nuclease-free water. An anchored Oligo-dT primer [d(T)₂₃VN] forces the primer to anneal to the beginning of the polyA tail. The optimized Random Primer Mix provides random and consistent priming sites covering the entire RNA templates including both mRNAs and non-polyadenylated RNAs.

The OneTaq RT-PCR Kit Includes:

- 10X M-MuLV Enzyme Mix
- 2X M-MuLV Reaction Mix
- OneTaq Hot Start 2X Master Mix with Standard Buffer
- Random Primer Mix (60 µM)
- Oligo d(T)₂₃VN Primer (50 µM)**
- Nuclease-free Water



First strand cDNA synthesis. Synthesis was carried out in the presence of 1X M-MuLV Enzyme Mix at 42°C using 0.5 µg of human spleen total RNA in the presence of d(T)₂₃VN (lanes 1, 4 and 7) or Random Hexamer Mix (lanes 2, 5 and 8). No-RT controls were lanes 3, 6 and 9. OneTaq Hot Start Master Mix was used to amplify a 1.5 kb fragment of beta-actin gene, a 0.6 kb fragment of GAPDH gene, and a 5.5 kb fragment from p532 gene in 35 cycles. The marker lane (M) contains 1 kb Plus DNA Ladder (NEB #N3200).

** Oligo d(T)₂₃VN and Random Primer Mix contain 1 mM dNTP

PreCR® Repair Mix

#M0309S	30 reactions	\$164
#M0309L	150 reactions	\$654

Companion Product:

beta-Nicotinamide adenine dinucleotide (NAD+)	
#B9007S	0.2 ml\$34

- Repair DNA prior to its use in DNA-related technologies
- Easy-to-use protocols
- Does not harm template

Need to repair FFPE-treated DNA prior to next gen sequencing?

See page 152 for more information.

Description: The PreCR Repair Mix is an enzyme cocktail formulated to repair damaged template DNA prior to its use in the polymerase chain reaction (PCR), microarrays, or other DNA technologies. The PreCR Repair Mix is active on a broad range of DNA damage, including those that block PCR (e.g., apurinic/apyrimidinic sites, thymidine dimers, nicks and gaps) and those that are mutagenic (e.g., deaminated cytosine and 8-oxo-guanine). In addition, it will remove a variety of moieties from the 3' end of DNA leaving a hydroxyl group. The PreCR Repair Mix will not repair all damage that inhibits/interferes with PCR. It can be used in conjunction with any thermophilic polymerase.

Applications:

- Repair DNA prior to its use as a template in PCR or other DNA technologies

Reagents Supplied:


1X PreCR Repair Mix
10X ThermoPol Reaction Buffer
100X NAD⁺ solution
Control Template (UV damaged λ DNA)
PCR primers for control template
Purified BSA

Repair of different types of DNA damage with the

PreCR Repair Mix. The gel shows amplification of damaged DNA that was either not treated (–) or treated (+) with the PreCR Repair Mix. Type of DNA damage is shown. Note: heat treated DNA is incubated at 99°C for 3 minutes. Marker (M) is the 1 kb Plus DNA Ladder (NEB #N3200).

Types of DNA Damage

DNA DAMAGE	CAUSE	REPAIRED BY PRECR REPAIR MIX?
abasic sites	hydrolysis	yes
nicks	hydrolysis nucleases shearing	yes
thymidine dimers	UV radiation	yes
blocked 3'-ends	multiple	yes
oxidized guanine	oxidation	yes
oxidized pyrimidines	oxidation	yes
deaminated cytosine	hydrolysis	yes
fragmentation	hydrolysis nucleases shearing	no
protein-DNA crosslinks	formaldehyde	no



Sulfolobus DNA Polymerase IV

#M0327S	100 units	\$124
---------	-----------------	-------

For a complete listing of Deoxynucleotide Solutions, see page 85.

- Synthesis of DNA through DNA lesions (lesion bypass)
- DNA Repair

Description: *Sulfolobus* DNA Polymerase IV is a thermostable Y-family lesion-bypass DNA Polymerase that efficiently synthesizes DNA across a variety of DNA template lesions.



Concentration: 2,000 units/ml

Therminator™ DNA Polymerase

#M0261S	200 units	\$106
#M0261L	1,000 units	\$423

For a complete listing of Deoxynucleotide Solutions, see page 85.

- Incorporation of modified nucleotides
- DNA sequencing by partial ribosubstitution
- DNA sequencing or SNP analysis using dideoxy or acyclo chain terminators

Description: Therminator DNA Polymerase is a 9°N™ DNA Polymerase variant with an enhanced ability to incorporate modified substrates such as dideoxynucleotides, ribonucleotides and acyclo nucleotides.

Source: An *E. coli* strain that carries the 9°N (D141A / E143A / A485L) DNA Polymerase gene, a genetically engineered form of the native DNA polymerase from *Thermococcus* species 9°N-7.

Concentration: 2,000 units/ml

Usage Notes: Amplification of extended regions may require optimization of reaction conditions.

DNA Polymerase I (*E. coli*)



#M0209S	500 units\$70
#M0209L	2,500 units\$281

Description: DNA Polymerase I (*E. coli*) is a DNA-dependent DNA polymerase with inherent 3'→5' and 5'→3' exonuclease activities. The 5'→3' exonuclease activity removes nucleotides ahead of the growing DNA chain, allowing nick translation.

Concentration: 10,000 units/ml

Heat Inactivation: 75°C for 20 minutes

Usage Notes: DNase I is not included with this enzyme and must be added for nick translation reactions.

- Nick translation of DNA
- Second strand cDNA synthesis

DNA Polymerase I, Large (Klenow) Fragment



#M0210S	200 units\$63
#M0210L	1,000 units\$252
for high (10X) concentration		
#M0210M	1,000 units\$252

For a complete listing of Deoxynucleotide Solutions, see page 85.

Description: DNA Polymerase I, Large (Klenow) Fragment was originally derived as a proteolytic product of *E. coli* DNA Polymerase I that retains polymerase and 3'→5' exonuclease activity, but lacks 5'→3' exonuclease activity. Klenow retains the polymerization fidelity of the holoenzyme without degrading 5' termini.

Concentration: 5,000 and 50,000 units/ml

Heat Inactivation: 75°C for 20 minutes

Usage Notes: Elevated temperatures, excessive amounts of enzyme, failure to supplement with dNTPs or long reaction times will result in recessed ends due to the 3'→5' exonuclease activity of the enzyme.

- Generates probes using random primers
- Removal of 3' overhangs or fill-in of 5' overhangs to form blunt ends
- Second strand cDNA synthesis

Source: An *E. coli* strain that contains the *E. coli* *polA* gene that has had its 5'→3' exonuclease domain removed.

Klenow Fragment (3'→5' exo⁻)



#M0212S	200 units\$63
#M0212L	1,000 units\$252
for high (10X) concentration		
#M0212M	1,000 units\$252

For a complete listing of Deoxynucleotide Solutions, see page 85.

Description: Klenow Fragment (3'→5' exo⁻) is an N-terminal truncation of DNA Polymerase I that retains polymerase activity, but has lost the 5'→3' exonuclease activity, and has mutations (D355A, E357A) that abolish the 3'→5' exonuclease activity.

Concentration: 5,000 and 50,000 units/ml

Heat Inactivation: 75°C for 20 minutes

Usage Notes: Klenow Fragment (3'→5' exo⁻) is not suitable for generating blunt ends because it lacks the 3'→5' exonuclease activity necessary to remove non-templated 3' additions.

- Generates probes using random primers
- Random priming labeling
- Second strand cDNA synthesis

Source: An *E. coli* strain containing a plasmid with a fragment of the *E. coli* *polA* (D355A, E357A) gene starting at codon 324.

T4 DNA Polymerase



#M0203S	150 units\$69
#M0203L	750 units\$276

Companion Product:

Quick Blunting™ Kit		
#E1201S	20 reactions\$86
#E1201L	100 reactions\$347

For a complete listing of Deoxynucleotide Solutions, see page 85.

Description: T4 DNA Polymerase catalyzes the synthesis of DNA in the 5'→3' direction and requires the presence of template and primer. This enzyme has a 3'→5' exonuclease activity which is much more active than that found in *E. coli* DNA Polymerase I. Unlike DNA Polymerase I, T4 DNA Polymerase does not have a 5'→3' exonuclease function.

Heat Inactivation: 75°C for 20 minutes

Usage Notes: Elevated temperatures, excessive amounts of enzyme, failure to supplement with dNTPs or long reaction times will result in recessed ends due to the 3'→5' exonuclease activity of the enzyme.

Note: This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

- Gap filling (no strand displacement activity)
- Removal of 3' overhangs or fill-in of 5' overhangs to form blunt ends
- Probe labeling using replacement synthesis
- Single-strand deletion subcloning

Concentration: 3,000 units/ml

T7 DNA Polymerase (unmodified)

#M0274S	300 units	\$71
#M0274L	1,500 units	\$286

For a complete listing of Deoxynucleotide Solutions, see page 85.

- Gap-filling reaction (no strand displacement)

Description: T7 DNA Polymerase catalyzes the replication of T7 phage DNA during infection. The protein dimer has two catalytic activities: DNA polymerase activity and strong 3' → 5' exonuclease. The high fidelity and rapid extension rate of the enzyme make it particularly useful in copying long stretches of DNA template.

Source: T7 DNA Polymerase consists of two subunits: T7 gene 5 protein (84 kDa) and *E. coli* thioredoxin (12 kDa). Each protein is cloned and overexpressed in a T7 expression system in *E. coli*.

Reaction Conditions: 1X T7 DNA Polymerase Reaction Buffer. Supplement with BSA and dNTPs (not included). Incubate at 37°C. Heat inactivation: 75°C for 20 minutes.



Unit Definition: One unit is defined as the amount of enzyme that will incorporate 10 nmol of dNTP into acid insoluble material in 30 minutes at 37°C.

Concentration: 10,000 units/ml

Usage Notes: The high polymerization rate of the enzyme makes long incubations unnecessary. T7 DNA Polymerase is not suitable for DNA sequencing.

Bsu DNA Polymerase, Large Fragment

#M0330S	200 units	\$68
#M0330L	1,000 units	\$272

For a complete listing of Deoxynucleotide Solutions, see page 85.

- Random primer labeling
- Second strand cDNA synthesis
- Single dA tailing
- Strand displacement DNA synthesis

Description: *Bsu* DNA Polymerase I, Large Fragment retains the 5' → 3' polymerase activity of the *Bacillus subtilis* DNA polymerase I, but lacks the 5' → 3' exonuclease domain. This large fragment naturally lacks 3' → 5' exonuclease activity.

Source: An *E. coli* strain that contains the *Bacillus subtilis* DNA polymerase I gene (starting from codon 297 thus lacking the 5' → 3' exonuclease domain).

Concentration: 5,000 units/ml

Heat Inactivation: 75°C for 20 minutes



Usage Notes: *Bsu* DNA Polymerase, Large Fragment is not suitable for generating blunt ends because it lacks the 3' → 5' exonuclease activity necessary to remove non-templated 3' additions.

Bsu DNA Polymerase, Large Fragment retains 50% activity at 25°C and is twice as active as Klenow Fragment (3' → 5' exo-) at this temperature.

Terminal Transferase

#M0315S	500 units	\$73
#M0315L	2,500 units	\$294

For a complete listing of Deoxynucleotide Solutions, see page 85.

- Addition of homopolymer tails to the 3' ends of DNA
- Labeling the 3' ends of DNA with modified nucleotides (e.g., ddNTP, DIG-dUTP)
- TUNEL assay (in situ localization of apoptosis)
- TdT dependent PCR

Description: Terminal Transferase (TdT) is a template independent polymerase that catalyzes the addition of deoxynucleotides to the 3' hydroxyl terminus of DNA molecules. Protruding, recessed or blunt-ended double or single-stranded DNA molecules serve as a substrate for TdT. The 58.3 kDa enzyme does not have 5' or 3' exonuclease activity. The addition of Co²⁺ in the reaction makes tailing more efficient.



Concentration: 20,000 units/ml

Heat Inactivation: 75°C for 20 minutes

Polymerase Reaction Buffers

Q5

Q5 Reaction Buffer Pack [Q5 Reaction Buffer (5X), Q5 High GC Enhancer (5X)]
#B9027S 6.0 ml \$29

Phusion*

Phusion HF Buffer Pack [Phusion HF Reaction Buffer (5X), MgCl₂ (50 mM), DMSO]
#B0518S 6.0 ml \$25

Phusion GC Buffer Pack [Phusion GC Reaction Buffer (5X), MgCl₂ (50 mM), DMSO]
#B0519S 6.0 ml \$25

Taq

Standard Taq Reaction Buffer
[Standard Taq Reaction Buffer (10X), MgCl₂ (25 mM)]
#B9014S 6.0 ml \$23

Standard Taq (Mg-free) Reaction Buffer Pack
[Standard Taq (Mg-free) Reaction Buffer (10X), MgCl₂ (25 mM)]
#B9015S 6.0 ml \$23

Other

ThermoPol Reaction Buffer Pack
[ThermoPol Reaction Buffer (10X), MgSO₄ (100 mM)]
#B9004S 6.0 ml \$23

Isothermal Amplification Buffer Pack
[Isothermal Amplification Buffer (10X)]
#B0537S 6.0 ml \$29

Isothermal Amplification Buffer II Pack
[Isothermal Amplification Buffer II (10X)]
#B0374S 6.0 ml \$28

* Phusion DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific. Phusion® is a registered trademark and property of Thermo Fisher Scientific.

Description: Q5 Reaction Buffer and High GC Enhancer are provided with both Q5 and Q5 Hot Start High-Fidelity DNA Polymerases.

Phusion High-Fidelity DNA Polymerase is supplied with 5X Phusion HF Buffer, 5X Phusion GC Buffer, DMSO, and 50 mM MgCl₂.

Standard Taq Reaction Buffer is provided with Taq DNA Polymerase as an alternative to the ThermoPol Reaction Buffer.

ThermoPol Reaction Buffer is provided with Taq, Vent, Deep Vent, Bst Full Length and Bst Large Fragment, Sulfolobus IV and Terminator DNA Polymerases; this buffer contains 2 mM MgSO₄ when the buffer is diluted to its final 1X concentration.

Isothermal Amplification Buffer is supplied with Bst 2.0 and Bst 2.0 WarmStart DNA Polymerases.

Isothermal Amplification Buffer II is supplied with Bst 3.0 DNA Polymerase.

Nucleotides

Acyclonucleotide Set
#N0460S 0.5 µmol of each \$72

Deoxynucleotide (dNTP) Solution Set
#N0446S 25 µmol of each \$176

Deoxynucleotide (dNTP) Solution Mix
#N0447S 8 µmol of each \$64
#N0447L 40 µmol of each \$255

Ribonucleotide Solution Set
#N0450S 10 µmol of each \$72
#N0450L 50 µmol of each \$289

Ribonucleotide Solution Mix
#N0466S 10 µmol of each \$75
#N0466L 50 µmol of each \$302

7-deaza-dGTP
#N0445S 0.3 µmol of each \$70
#N0445L 1.5 µmol of each \$280

Adenosine 5'-Triphosphate (ATP)
#P0756S 1.0 ml \$37
#P0756S 5.0 ml \$143

5-methyl-dCTP
#N0356S 1 µmol \$73

dATP Solution
#N0440S 25 µmol \$56

NEW
dUTP Solution
#N0459S 25 µmol \$62

NEW
dGTP Solution
#N0442S 25 µmol \$61

Description:

Deoxynucleotide Solution Set:
Four separate solutions of ultrapure deoxynucleotide (dATP, dCTP, dGTP and dTTP). Each deoxynucleotide is supplied at a concentration of 100 mM.

Deoxynucleotide Solution Mix:
An equimolar solution of ultrapure dATP, dCTP, dGTP and dTTP. Each deoxynucleotide is supplied at a concentration of 10 mM.

Ribonucleotide Solution Set:
Four separate solutions of ATP, CTP, GTP and UTP, pH 7.5, as sodium salts.

Ribonucleotide Solution Mix:
A buffered equimolar solution of ribonucleotide triphosphates: rATP, rCTP, rGTP and rUTP. Each nucleotide is supplied at a concentration of 25 mM (total rNTP concentration equals 100 mM).

7-deaza-dGTP:
7-deaza contains a 5 mM solution of 7-deaza-dGTP as a dilithium salt.

5-methyl-dCTP:
dm5CTP supplied as 10 mM solution at pH 7.0.

dATP Solution:
dATP Solution contains a 100 mM solution of dATP as a sodium salt at pH 7.4.

dUTP Solution:
dUTP Solution contains a 100 mM solution of dUTP as a sodium salt at pH 7.5.

dGTP Solution:
dGTP Solution contains a 100 mM solution of dGTP as a sodium salt at pH 7.5.

Acyclonucleotide Set:
Four separate tubes of acyNTPs (acyATP, acyCTP, acyGTP and acyTTP). Acyclonucleotides are supplied as a dry powder. Addition of 50 µl of distilled or de-ionized (Milli-Q®) water will result in a final concentration of 10 mM acyNTP.

Acyclonucleotides (acyNTPs) act as chain terminators and are thus useful in applications that normally employ dideoxynucleotides such as DNA sequencing and SNP detection. AcyNTPs are especially useful in applications with archaeal DNA Polymerases, more specifically with Terminator DNA Polymerase. Terminator DNA Polymerase is an engineered enzyme with an increased capacity to incorporate analogs with altered sugars, such as ribonucleotides, dideoxynucleotides, 2' deoxynucleotides and especially acyclo-base analogs.

For RNA Cap Analogs, see page 204.

MILLI-Q® is a registered trademark of Millipore, Inc.

Products for cDNA Synthesis

cDNA SYNTHESIS	FEATURES	SIZE	PRICE
KITS			
LunaScript RT SuperMix Kit (NEB #E3010)	<ul style="list-style-type: none"> Ideal for cDNA synthesis in a two-step RT-qPCR workflow Single tube supermix contains random hexamer and oligo-dT primers, dNTPs, Murine RNase Inhibitor, and Luna Reverse Transcriptase Visible blue tracking dye for easy reaction setup Fast 13-minute protocol 	25/100 rxns	\$141/\$437
NEW LunaScript RT Master Mix Kit (Primer-free) (NEB #E3025)	<ul style="list-style-type: none"> Ideal for first strand cDNA synthesis Master mix contains dNTPs, Murine RNase Inhibitor and Luna Reverse Transcriptase Visible blue tracking dye for easy reaction setup Flexible choice of primers 	25/100 rxns	\$117/\$364
ProtoScript® II First Strand cDNA Synthesis Kit (NEB #E6560)	<ul style="list-style-type: none"> Generates cDNA at least 10 kb in length Contains ProtoScript II Reverse Transcriptase, an enzyme with increased thermostability and reduced RNase H activity Convenient 2-tube kit includes dNTPs, Oligo-dT primer and Random Primer Mix 	30/150 rxns	\$171/\$685
ProtoScript First Strand cDNA Synthesis Kit (NEB #E6300)	<ul style="list-style-type: none"> Generates cDNA at least 5 kb in length Contains M-MuLV Reverse Transcriptase Convenient 2-tube kit includes dNTPs, Oligo-dT primer and Random Primer Mix 	30/150 rxns	\$162/\$645
Template Switching RT Enzyme Mix (NEB #M0466)	<ul style="list-style-type: none"> Incorporates a universal adaptor sequence at the 3' end of cDNA during the RT reaction Enzyme mix and buffer are optimized for efficient template switching RT enzyme mix includes RNase Inhibitor High sensitivity for cDNA amplification – enables transcriptome analysis by RNA-seq from single cells or as low as 2 pg of human total RNA Robust and simple workflow for 5' Rapid Amplification of cDNA Ends (RACE) Retains the complete 5' end of transcripts for 2nd Strand cDNA Synthesis 	20/100 rxns	\$90/\$360
STANDALONE REAGENTS			
ProtoScript II Reverse Transcriptase (NEB #M0368) An alternative to SuperScript® II	<ul style="list-style-type: none"> RNase H⁻ mutant of M-MuLV Reverse Transcriptase with increased thermostability and reduced RNase H activity Increased reaction temperatures (37–50°C) 	4,000/10,000/40,000 units	\$92/\$184/\$656
M-MuLV Reverse Transcriptase (NEB #M0253)	<ul style="list-style-type: none"> Robust reverse transcriptase for a variety of templates Standard reaction temperatures (37–45°C) 	10,000/50,000 units	\$73/\$294
AMV Reverse Transcriptase (NEB #M0277)	<ul style="list-style-type: none"> Robust reverse transcriptase for a broad temperature range (37–52°C) Can be used for templates requiring higher reaction temperatures 	200/1,000 units	\$74/\$298
WarmStart RTx Reverse Transcription (NEB #M0380)	<ul style="list-style-type: none"> Permits room temperature reaction setup Increased reaction temperatures (50–65°C) Optimized for RT-LAMP isothermal detection 	50/250 rxns	\$69/\$277

More information about our reverse transcriptases and cDNA synthesis kits can be found in the RNA analysis chapter, pages 206–209.

Monarch PCR & DNA Cleanup Kit (5 µg)

#T1030S	50 preps\$104
#T1030L	250 preps\$475

Companion Products:

Monarch DNA Cleanup Columns (5 µg)		
#T1034L	100 columns + tubes\$131
Monarch DNA Cleanup Binding Buffer		
#T1031L	175 ml\$100
Monarch DNA Wash Buffer		
#T1032L	25 ml\$35
Monarch DNA Elution Buffer		
#T1016L	25 ml\$35
Monarch Plasmid Miniprep Kit		
#T1010S	50 preps\$85
#T1010L	250 preps\$361
Monarch DNA Gel Extraction Kit		
#T1020S	50 preps\$104
#T1020L	250 preps\$475

Description: The Monarch PCR & DNA Cleanup Kit (5 µg) is a rapid and reliable method for the purification and concentration of up to 5 µg of high-quality, double-stranded DNA from enzymatic reactions such as PCR, restriction digestion, ligation and reverse transcription. This method employs a bind/wash/elute workflow with minimal incubation and spin times, resulting in purification in less than 5 minutes. DNA Cleanup Binding Buffer is used to dilute the samples and ensure they are compatible for loading onto the proprietary silica matrix under high salt conditions. The DNA Wash Buffer ensures enzymes, short primers (≤ 40 nt), detergents and other low-molecular weight reaction components (e.g., nucleotides, DMSO, betaine) are removed, thereby allowing low-volume elution of concentrated, high-purity DNA. Eluted DNA is ready for use in restriction digests, DNA sequencing, ligation and other enzymatic manipulations. The unique column design ensures no buffer retention and no carryover of contaminants, allowing elution of sample in volumes as low as 6 µl. A slight protocol modification enables purification of small DNA and oligonucleotides.

Applications:

- PCR cleanup
- Enzymatic reaction cleanup
- cDNA cleanup
- Labeling cleanup
- Plasmid cleanup
- Oligo cleanup

The Monarch PCR & DNA Cleanup Kit Includes:

- Monarch DNA Cleanup Columns (5 µg)
- Monarch DNA Cleanup Binding Buffer
- Monarch DNA Wash Buffer
- Monarch DNA Elution Buffer
- Monarch Collection Tubes

With Monarch PCR & DNA Cleanup Kit, you can purify your DNA in as little as 5 minutes.

- Elute in as little as 6 µl
- Prevent buffer retention and salt carry-over with optimized column design
- Purify small DNA and oligos with a slight protocol modification
- Save time with fast, user-friendly protocol
- Purchase optimized kit formats or buffers & columns separately for your convenience

Exo-CIP™ Rapid PCR Cleanup Kit


#E1050S	100 reactions\$94
#E1050L	400 reactions\$312

- 5 minute protocol for enzymatic cleanup of primers and dNTPs
- Improves sequencing results, allowing longer reads


Description: The Exo-CIP Rapid PCR Cleanup Kit contains optimized formulations of thermolabile Exo-nuclease I and thermolabile Calf Intestinal Phosphatase and is used to rapidly degrade residual PCR primers and dephosphorylate excess dNTPs after amplification. Degradation occurs in only 4 minutes at 37°C, and is immediately followed by rapid inactivation of the enzymes by heating for 1 minute at 80°C. In just 5 minutes, the PCR product is ready for downstream analysis such as Sanger sequencing, SNP detection, or library preparation for NGS. The Exo-CIP Rapid PCR Cleanup Kit is compatible with all commonly-used reaction buffers.

The Exo-CIP Rapid PCR Cleanup Kit Includes:

- Exo-CIP Tube A (thermolabile Exo I)
- Exo-CIP Tube B (thermolabile CIP)



Exo-CIP Rapid PCR Cleanup Kit workflow. 1 µl of Exo-CIP Tube A (thermolabile Exo I) and 1 µl of Exo-CIP Tube B (thermolabile CIP) are added to the PCR product to degrade excess primers and dNTPs. The mixture is incubated at 37°C for 4 minutes, followed by a 1 minute incubation at 80°C to irreversibly inactivate both enzymes. The cleaned PCR product is ready for downstream applications or analysis.



Ask Big Questions

Voltaire said, “Judge a man by his questions rather than by his answers.”

Asking questions is easy. What’s for dinner? Where’s the library? But asking big questions is much harder. It takes courage to challenge the status quo, and to ask questions that take time, energy and exploration to answer. But aren’t these questions, the tough ones with no definitive answer, the ones worth asking? These are the questions that drive discovery and innovation.

Questions can be powerful tools. The art of asking impactful questions can be innate in curious and emotionally intelligent individuals, but it can also be learned. The first step toward asking better questions is to avoid “yes” or “no” questions — these often lead to short answers with little detail. Anyone with young children knows this well! Instead, ask questions that start with “what,” “why,” “how,” or “if” — these types of questions are open-ended, and can lead to a more thoughtful discussion and new ideas to explore. This brings us to the second step, asking follow-up questions that dig deeper. Scientists know that asking one great question can get the thought process rolling, and that one great question often leads to two more, which can easily lead to four more, and so on.

Don Comb’s passion was reflected in the intensity of his questions. New England Biolabs’ CEO Jim Ellard recalls a dinner some 30 years ago when Don leaned across the table and asked, “What if you could come up with a way to keep algae, seaweed and barnacles from sticking to the bottom of boats? You would reduce drag, conserve fuel and revolutionize the shipping industry!” — always the inquisitor!

Over the years, New England Biolabs’ culture has been shaped by great questions that Don asked early on — for example, “How can we run a company that does great things without impacting the environment?” A company can take responsibility for the inevitable footprint it creates by offsetting negative impacts through the use of solar energy, managing treatment of its own wastewater, supporting the reforestation of land to improve global carbon capture, practicing extensive recycling and composting, and investing in conservation efforts around the world. Another big question Don asked was, “How can we help underdeveloped countries that are struggling with filariasis?” For 40 years, NEB scientists have been dedicated to research on neglected tropical diseases caused by parasitic worms, such as lymphedema and elephantiasis. These diseases are the second leading cause of disability worldwide, and our scientists continue to improve understanding and detection, and freely share their findings with the research community.

Asking provocative questions has stimulated many exciting conversations at NEB, and has helped nurture the curiosity of its employees. In 2012, NEB introduced the Innovation Awards to emphasize the importance of “blue-sky” thinking, the courage to take risks, and the drive to implement new ideas. These awards are given out annually, and recognize employees who contribute and bring to fruition novel ideas that help make NEB a better place and help to carry on Don’s legacy of provocative curiosity.

View of the moody sky
over Thailand at twilight.
Credit: Boy Anupong, Getty Images

Learn more about
NEB’s Innovation
Awards.



DNA Modifying Enzymes & Cloning Technologies



The trusted source for DNA-modifying enzymes & cloning technologies.

Molecular biology, which also includes synthetic biology, is a fundamental area of research and development. Central to these advances has been the use of DNA modifying enzymes and novel cloning technologies. Some common examples of DNA-modifying enzymes include kinases, ligases, methylases, exonucleases and endonucleases, while newer cloning technologies include NEBuilder® HiFi DNA Assembly and Golden Gate Assembly.

Over 45 years as a leader in enzyme technologies gives you confidence in the products and support you'll receive. NEB continues to serve the scientific community by providing the tools to carry out the most innovative research, from start to finish. All NEB products pass stringent quality control assays to ensure the highest level of functionality and purity.

NEB offers several online tools to aid in your cloning experiments, including:

- **NEBcloner®** – find the right products and protocols for each step of your traditional cloning experiment, including double digests and mutagenesis
- **NEBioCalculator®** – use this tool for your scientific calculations and conversions
- **NEBuilder Assembly Tool** – use this tool for help with your DNA assembly primer design
- **ExoSelector** – find the right exonuclease for your workflows
- **NEB® Golden Gate Assembly Tool** – use this tool for help with construct design for Golden Gate Assembly
- **Ligase Fidelity Tools** – utilize ligation preferences for the design of high-fidelity Golden Gate Assembly

To view the full list of online tools available, see page 304.

Featured Products

- 94** NEBuilder HiFi DNA Assembly Master Mix & Cloning Kit
- 111** Exonuclease V
- 96** NEB Golden Gate Assembly Kit (BsaI-HFv2, BsmBI-v2)
- 98** Q5® Site-Directed Mutagenesis Kit

Featured Tools & Resources

- 92** Cloning Overview
- 331** Getting Started with Cloning
- 336** Troubleshooting Guide for Cloning



Visit www.neb.com/FastCloning to accelerate your cloning experiments with reagents from NEB.



Visit ClonewithNEB.com to view our online tutorials explaining each of the steps in the cloning workflow.



Cloning & Synthetic Biology

Cloning Workflow Comparison	92
NEBuilder HiFi DNA Assembly Master Mix	94
NEBuilder HiFi DNA Assembly Cloning Kit	94
NEBuilder HiFi DNA Assembly Bundle for Large Fragments	94
Gibson Assembly Master Mix/Cloning Kit	95
Synthetic Biology/DNA Assembly Selection Chart	95
NEB Golden Gate Assembly Kits (BsaI-HFv2, BsmBI-v2)	96
BioBrick Assembly Kit	97
NEB PCR Cloning Kit (with and without competent cells)	97
USER & Thermolabile USER II Enzymes	97
Q5 Site Directed Mutagenesis Kit (with and without competent cells)	98
Quick Blunting Kit	98

DNA Ligases

DNA Ligase Selection Chart	99
Substrate-based Ligase Selection Chart	100
T4 DNA Ligase	101
Instant Sticky-end Ligase Master Mix	101
Blunt T/A Ligase Master Mix	101
Quick Ligation Kit	101
ElectroLigase	101
Immobilized T4 DNA Ligase	101
Hi-T4 DNA Ligase	101
Salt-T4 DNA Ligase	101
T3 DNA Ligase	102
T7 DNA Ligase	102
<i>E. coli</i> DNA Ligase	102
HiFi <i>Taq</i> DNA Ligase	103
<i>Taq</i> DNA Ligase	103
9°N DNA Ligase	104
SplintR Ligase	214

DNA Kinases

T4 Polynucleotide Kinase	104
T4 Polynucleotide Kinase (3' phosphatase minus)	104
5-hydroxymethyluridine Kinase	105

DNA Phosphatases & Sulfurylases

Phosphatase Selection Chart	105
Quick CIP	105
Shrimp Alkaline Phosphatase (rSAP)	106
Antarctic Phosphatase	106
Pyrophosphatase, Inorganic (<i>E. coli</i> , yeast)	203
Thermostable Inorganic Pyrophosphatase	203
NudC Pyrophosphatase	203
Apyrase	107
<i>Tte</i> UvrD Helicase	107

Nucleases

Properties of Exonucleases and Non-specific Endonucleases	108
Common Apps for Exonucleases & Endonucleases	109
DNase I (RNase-Free)	217
Lambda Exonuclease	110
Exonuclease I	110
Thermolabile Exonuclease I	110
Exonuclease III	111
Exonuclease V (RecBCD)	111
Exonuclease VII	111
Exonuclease VIII, truncated	112
Exonuclease T	112
Thermostable FEN1	112

Micrococcal Nuclease	113
Mung Bean Nuclease	113
Nuclease BAL-31	113
Nuclease P1	114
RecJ ₁	114
T5 Exonuclease	114
T7 Exonuclease	115
Nucleoside Digestion Mix	115

DNA Repair Enzymes

Properties of DNA Repair Enzymes and Structure-specific Endonucleases	116
DNA Repair Enzymes on Damaged and Non-Standard Bases	117
APE 1	118
Endonuclease IV	118
<i>Tth</i> Endonuclease IV	118
Endonuclease III (Nth)	119
<i>Tma</i> Endonuclease III	119
Endonuclease V	119
Endonuclease VIII	120
Fpg	120
hAAG	120
T4 PDG (T4 Endonuclease V)	121
Uracil-DNA Glycosylase (UDG)	121
<i>Afu</i> Uracil-DNA Glycosylase (UDG)	121
Antarctic Thermolabile UDG	122
hSMUG1	122
T7 Endonuclease I	122
PreCR Repair Mix	82

Other Enzymes

Cre Recombinase	123
TelN Protelomerase	123
Topoisomerase I (<i>E. coli</i>)	123
β -Agarase I	123

DNA Methyltransferases

CpG Methyltransferase (M.SssI)	124
GpC Methyltransferase (M.CviPI)	124
AluI Methyltransferase	124
BamHI Methyltransferase	124
<i>dam</i> Methyltransferase	125
EcoGII Methyltransferase	125
EcoRI Methyltransferase	125
HaeIII Methyltransferase	125
HhaI Methyltransferase	125
HpaII Methyltransferase	126
MspI Methyltransferase	126
<i>TaqI</i> Methyltransferase	126
Human Dnmt1 Methyltransferase	286

ssDNA Binding Proteins

RecA, RecA ₁	127
T4 Gene 32 Protein	127
ET SSB (Extreme Thermostable ssBinding Protein)	127

Cloning Plasmids & DNAs

Cloning Plasmids (pBR322, pUC19)	128
DNAs (M13mp18, Lambda, ϕ X174)	128

Other


M13K07 Helper Phage	128
<i>Tth</i> Argonaute (TiAgo)	128
Competent Cell Selection Chart for Cloning	129
Monarch Nucleic Acid Purification Kits	129

Recombinant Enzyme

One or more of these products are covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at gbd@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications. Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at <https://www.neb.com/support/terms-of-sale>. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

Cloning Workflow Comparison

The figure below compares the steps for the various cloning methodologies, along with the time needed for each step in the workflows.






* Note that times are based on estimates for moving a gene from one plasmid to another. If the source for gene transfer is gDNA, add 2 hours to calculation for the traditional cloning method. Total time does not include transformation, isolation or analysis.

** 70 minutes for recombination occurs on second day




For help with choosing the right product for each step in the cloning workflow, try NEBcloner at NEBcloner.neb.com

VECTOR PREPARATION



DNA Isolation (Plasmid Purification)



DNA Analysis



Protein Expression

Functional Analysis

Site-Directed Mutagenesis

NEBuilder® HiFi DNA Assembly Master Mix & Cloning Kit



NEBuilder HiFi DNA Assembly Master Mix	
#E2621S	10 reactions \$159
#E2621L	50 reactions \$630
#E2621X	250 reactions \$2,520

NEBuilder HiFi DNA Assembly Cloning Kit	
#E5520S	10 reactions \$188

NEBuilder HiFi DNA Assembly Bundle for Large Fragments	
#E2623S	20 reactions \$495

- Simple and fast seamless cloning in as little as 15 minutes
- Use one system for both "standard-size" cloning and larger gene assembly products (up to 11 fragments and 20 kb)
- DNA can be used immediately for transformation or as template for PCR or RCA
- Adapts easily for multiple DNA manipulations, including site-directed mutagenesis
- Enjoy less screening/re-sequencing of constructs, with virtually error-free, high-fidelity assembly
- Use NEBuilder HiFi in successive rounds of assembly, as it removes 5' and 3' end mismatches
- Bridge two ds-fragments with a synthetic ssDNA oligo for simple and fast construction (e.g., linker insertion or gRNA library)
- No licensing fee requirements from NEB for NEBuilder products
- NEBuilder HiFi DNA Assembly Cloning Kit includes the NEBuilder HiFi DNA Assembly Master Mix and NEB 5-alpha Competent *E. coli*
- NEBuilder HiFi DNA Assembly Bundle for Large Fragments includes the NEBuilder HiFi DNA Assembly Master Mix and NEB 10-beta Competent *E. coli* for assemblies larger than 15 kb

To learn how simple NEBuilder HiFi is, visit NEBuilderHiFi.com

Description: NEBuilder HiFi DNA Assembly Master Mix was developed to improve the efficiency and accuracy of DNA assembly. This method allows for seamless assembly of multiple DNA fragments, regardless of fragment length or end compatibility. This method has been used to assemble either single-stranded oligonucleotides or different sizes of DNA fragments with varied overlaps (15–80 bp). It has utility for the synthetic biology community, as well as those interested in one-step cloning of multiple fragments due to its ease of use, flexibility and simple master-mix format. The reaction features different enzymes that perform in the same buffer:

- Exonuclease creates single-stranded 3' overhangs that facilitate the annealing of fragments that share complementarity at one end (the overlap region)
- The polymerase fills in gaps within each annealed fragment
- The DNA ligase seals nicks in the assembled DNA

The end result is a double-stranded, fully sealed DNA molecule that can serve as template for PCR, RCA or a variety of other molecular biology applications, including direct transformation of *E. coli*.

The NEBuilder HiFi DNA Assembly Cloning Kit combines the power of the NEBuilder HiFi DNA Assembly Master Mix with NEB 5-alpha Competent *E. coli*, enabling fragment assembly and transformation in just under 2 hours.

NEBuilder HiFi Kits can be purchased with NEB 5-alpha Competent *E. coli* (Cloning Kit, NEB #E5520) or as a bundle with NEB 10-beta Competent *E. coli* (Bundle for Large Fragments, NEB #E2623). NEB 5-alpha competent cells are excellent for routine assemblies of 15 kb or less. NEB recommends NEB 10-beta competent cells for assemblies larger than 15 kb.

The NEBuilder HiFi DNA Assembly Master Mix Includes:

- NEBuilder HiFi DNA Assembly Master Mix
- NEBuilder Positive Control

The NEBuilder HiFi DNA Assembly Cloning Kit Includes:


- NEBuilder HiFi DNA Assembly Master Mix
- NEBuilder Positive Control
- NEB 5-alpha Competent *E. coli* (High Efficiency)
- SOC Outgrowth Medium
- pUC19 Control DNA

The NEBuilder HiFi DNA Assembly Bundle for Large Fragments Includes:

- NEBuilder HiFi DNA Assembly Master Mix
- NEBuilder Positive Control
- NEB 10-beta Competent *E. coli* (High Efficiency)
- NEB 10-beta/Stable Outgrowth Medium
- pUC19 Control DNA



Speed up your experimental design with our primer design tool at NEBuilder.neb.com



Overview of the NEBuilder HiFi DNA Assembly Cloning Method.



How does NEBuilder HiFi DNA Assembly work?

Gibson Assembly® Master Mix & Cloning Kit

Gibson Assembly Master Mix

#E2611S	10 reactions\$164
#E2611L	50 reactions\$649

Gibson Assembly Cloning Kit

#E5510S	10 reactions\$199
---------	--------------	------------

- Increased number of successful assembly products, particularly for longer or greater numbers of fragments
- Flexible sequence design with no need to engineer cloning sites
- Assemble multiple DNA fragments and transform in just under 2 hours
- Clone into any vector with no additional sequence added (scarless)
- No PCR clean-up step required

Description: Gibson Assembly Master Mix was developed by Dr. Daniel Gibson and his colleagues at the J. Craig Venter Institute and licensed to NEB by Synthetic Genomics, Inc. It allows for successful assembly of multiple DNA fragments, regardless of fragment length or end compatibility.

Gibson Assembly efficiently joins multiple overlapping DNA fragments in a single-tube isothermal reaction. The Gibson Assembly Master Mix includes three different enzymatic activities that perform in a single buffer:

- The exonuclease creates a single-stranded 3' overhang that facilitates the annealing of fragments that share complementarity at one end
- The polymerase fills in gaps within each annealed fragment
- The DNA ligase seals nicks in the assembled DNA

The Gibson Assembly Cloning Kit combines the power of the Gibson Assembly Master Mix with NEB 5-alpha Competent *E. coli*, enabling fragment assembly and transformation in just under 2 hours.

The Gibson Assembly Cloning Kit has been optimized for the assembly and cloning of up to 6 fragments.



The Gibson Assembly Master Mix Includes:

- Gibson Assembly Master Mix
- NEBuilder Positive Control

The Gibson Assembly Cloning Kit Includes:

- Gibson Assembly Master Mix
- NEBuilder Positive Control
- NEB 5-alpha Competent *E. coli* (High Efficiency)
- SOC Outgrowth Medium
- pUC19 Control DNA

GIBSON ASSEMBLY® is a registered trademark of Synthetic Genomics Inc.



Speed up your experimental design with our primer design tool at NEBuilder.neb.com

To learn how simple Gibson Assembly is, view our online tutorials at NEBGibson.com

Synthetic Biology/DNA Assembly Selection Chart

	NEBuilder® HiFi DNA Assembly (NEB #E2621) (NEB #E5520) (NEB #E2623)	NEB Gibson® Assembly (NEB #E5510) (NEB #E2611)	NEB® Golden Gate Assembly Kit (BsaI-HFv2, BsmBI-v2) (NEB #E1601) (NEB #E1602)	USER® Enzyme (NEB #M5505) Thermolabile USER II Enzyme (NEB #M5508)
PROPERTIES				
Removes 5' or 3' End Mismatches	★★★	★	N/A	N/A
Assembles with High Fidelity at Junctions	★★★	★★	★★★	★★★
Tolerates Repetitive Sequences at Ends	★	★	★★★	★★★
Generates Fully Ligated Product	★★★	★★★	★★★	NR
Joins dsDNA with Single-stranded Oligo	★★★	★★	NR	NR
Assembles with High Efficiency with Low Amounts of DNA	★★★	★★	★★	★★
Accommodates Flexible Overlap Lengths	★★★	★★★	★	★★
APPLICATIONS				
Simple Cloning (1-2 Fragments)	★★★	★★★	★★★	★★★
4-6 Fragment Assembly (one pot)	★★★	★★★	★★★	★★★
7-11 Fragment Assembly (one pot)	★★★	★★	★★★	★★★
12-35 Fragment Assembly (one pot) ⁽¹⁾	★	★	★★★	NR
Template Construction for <i>In vitro</i> Transcription	★★★	★★★	★★★	★★★
Synthetic Whole Genome Assembly	★★★	★	★	★
Multiple Site-directed Mutagenesis	★★★	★★	★★	★★
Library Generation	★★★	★★★	★★★	★★
Metabolic Pathway Engineering	★★★	★★	★★★	★★★
TALENs	★★	★★	★★★	★★
Short Hairpin RNA Cloning (shRNA)	★★★	★★	★	★
gRNA Library Generation	★★★	★★	★	★
Large Fragment (> 10 kb) Assembly	★★★	★★★	★★★	★★
Small Fragment (< 100 bp) Assembly	★★★	★	★★★	★★★
Use in Successive Rounds of Restriction Enzyme Assembly	★★★	★	NR	★

- ★★★ Optimal, recommended product for selected application
 ★★ Works well for selected application
 ★ Will perform for selected application, but is not recommended
 N/A Not applicable to this application
 NR Not recommended

(1) Please visit www.neb.com/GoldenGate for more information.

NEB Golden Gate Assembly Kits

NEB Golden Gate Assembly Kit (Bsal-HFv2)		
#E1601S	20 reactions\$153
#E1601L	100 reactions\$408

NEW

NEB Golden Gate Assembly Kit (BsmBI-v2)		
#E1602S	20 reactions\$153
#E1602L	100 reactions\$408

Companion Products:

NEB 5-alpha Competent <i>E. coli</i> (High Efficiency)		
#C2987H	20 x 0.05 ml\$199
#C2987I	6 x 0.2 ml\$153
#C2987P	1 x 96 well plate\$503
#C2987R	1 x 384 well plate\$1,100
#C2987U	96 x 50 µl/tube\$748

NEB 10-beta Competent <i>E. coli</i> (High Efficiency)		
#C3019H	20 x 0.05 ml\$233
#C3019I	6 x 0.2 ml\$179

NEB Cloning Competent <i>E. coli</i> Sampler		
#C1010S	8 tubes\$122

Q5 Hot Start High-Fidelity 2X Master Mix		
#M0494S	100 rxns (50 µl vol)\$217
#M0494L	500 rxns (50 µl vol)\$862
#M0494X	500 rxns (50 µl vol)\$862

- Choose from Bsal-HFv2 or BsmBI-v2
- Seamless cloning – no scar remains following assembly
- Includes destination plasmid with T7/SP6 promoters
- Ordered assembly of multiple fragments (2-35+) in a single reaction*
- Can also be used for cloning of single inserts and library preparations
- Efficient with regions of high GC content and areas of repeats
- Compatible with a broad range of fragment sizes (< 100 bps to > 15 kb)

Type IIS Enzymes used in Golden Gate:

- BbsI (NEB #R0539)
- BbsI-HF (NEB #R3539)
- Bsal-HFv2 (NEB #R3733)
- BsmBI-v2 (NEB #R0739)
- BspQI (NEB #R0712)
- BtgZI (NEB #R0703)
- Esp3I (NEB #R0734)
- PaeCI (NEB #R0745)
- SapI (NEB #R0569)

Description: The NEB Golden Gate Assembly Kits (Bsal-HFv2 and BsmBI-v2) contain an optimized mix of restriction enzyme and T4 DNA Ligase. Together these enzymes can direct the assembly of multiple inserts/modules using the Golden Gate approach. Also included is the pGGAselect destination plasmid, which provides a backbone for your assembly, features convenient restriction enzyme sites for subcloning, and has T7/SP6 promoter sequences to enable *in vitro* transcription.

The efficient and seamless assembly of DNA fragments, is commonly referred to as Golden Gate Assembly. In this workflow, Type IIS restriction enzymes bind to their recognition sites but cut the DNA downstream from that site at a positional, not sequence-specific, cut site. Thus, a single Type IIS restriction enzyme can be used to generate DNA fragments with unique overhangs. As an example, Bsal has a recognition site of GGTCTC(N1/N5), where the GGTCTC represents the recognition/binding site, and the N1/N5 indicates the cut site is one base downstream on the top strand, and five bases downstream on the bottom strand. Assembly of digested fragments proceeds through annealing of complementary four base overhangs on adjacent fragments. The digested fragments and the final assembly no longer contain Type IIS restriction enzyme recognition sites, so no further cutting is possible. The assembly product accumulates with time.

***Note:** Complexities up to 24 fragments have been routinely achieved with both precloned and amplicon insert test systems. Complexities of 35+ fragments have only used amplicon inserts to date, as we have not yet constructed a 35 fragment precloned insert test system.

NEB Golden Gate Assembly Tool

Speed up your experimental design with our assembly tool at GoldenGate.neb.com

While particularly useful for multi-fragment assemblies such as Transcription Activator Like Effectors (TALEs) and TALEs fused to a FokI nuclease catalytic domain (TALENs), the Golden Gate method can also be used for cloning of single inserts and inserts from diverse populations that enable library creation.

Advances in Ligase Fidelity: Research at NEB has led to increased understanding of ligase fidelity, including the development of a comprehensive method for profiling end-joining ligation fidelity in order to predict which overhangs have improved fidelity. This research allows careful choice of overhang sets, which is especially important for achieving complex Golden Gate Assemblies. To learn more, visit www.neb.com/goldengate.

The NEB Golden Gate Assembly Kit (Bsal-HFv2) Includes:

- NEB Golden Gate Assembly Mix
- T4 DNA Ligase Reaction Buffer (10X)
- pGGAselect Destination Plasmid

The NEB Golden Gate Assembly Kit (BsmBI-v2) Includes:


- NEB Golden Gate Assembly Mix
- T4 DNA Ligase Reaction Buffer (10X)
- pGGAselect Destination Plasmid

Ligase Fidelity Tools

For the design of high-fidelity Golden Gate Assemblies

- Ligase Fidelity Viewer™ (v2) – visualize overhang ligation junctions
- GetSet™ – predict high-fidelity junction sets
- SplitSet™ – split DNA sequence for scarless high-fidelity assembly

Access these tools at www.neb.com/research/nebata-tools.



In its simplest form, Golden Gate Assembly requires a Type IIS recognition site, in this case, Bsal-HFv2 (GGTCTC), or BsmBI-v2 (CGTCTC) added to both ends of a dsDNA fragment. After digestion, these sites are left behind, with each fragment bearing the designed 4-base overhangs that direct the assembly.*



How does
Golden Gate
Assembly work?

BioBrick® Assembly Kit

#E0546S 50 reactions\$284

The BioBrick Assembly Kit was developed in partnership with Ginkgo BioWorks. For more details and for technical questions, please see: ginkgobioworks.com/support

BIOBRICK® is a registered trademark of The BioBricks Foundation. (<http://www.biobricks.org>)

Description: The BioBrick Assembly Kit provides a streamlined method for assembly of BioBrick parts into multi-component genetic systems. BioBrick parts are DNA sequences that encode a defined biological function and can be readily assembled with any other BioBrick part. The process for assembling any two BioBrick parts is identical and results in a new composite BioBrick part.

Please refer to the individual datacards for each reagent's recommended use and storage conditions.

Note: This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.



The BioBrick Assembly Kit Includes:

- EcoRI-HF
- XbaI
- SpeI
- PstI
- 10X NEBuffer r2.1
- T4 DNA Ligase
- 10X T4 DNA Ligase Buffer

NEB® PCR Cloning Kit (with or without competent cells)

NEB PCR Cloning Kit
#E1202S 20 reactions\$289

NEB PCR Cloning Kit
(without competent cells)
#E1203S 20 reactions\$168

- *In vitro* transcription with both SP6 & T7 promoters
- Easy cloning of all PCR products, including blunt and TA ends
- Fast cloning with 5-minute ligation step
- Simplified screening with low/no colony background and no blue/white selection
- Save time by eliminating purification steps
- More flanking restriction sites available for easy subcloning, including choice of two single digest options
- BsaI site removed to allow cloning of Golden Gate modules

Description: The NEB PCR Cloning Kit contains optimized Cloning Master Mixes with a proprietary ligation enhancer and a linearized vector that uses a novel mechanism for background colony suppression to give a low background. It allows simple and quick cloning of any PCR amplicon, whether the amplification reactions are performed with proofreading DNA polymerases, such as Q5® or Phusion® which produce blunt ends; or nonproofreading DNA polymerases, such as *Taq* or *Taq* mixes (One *Taq*, LongAmp *Taq*) which produce single-base overhangs. This is possible due to "invisible" end polishing components in the master mix that are active during the ligation step. The kit also allows direct cloning from amplification reactions without purification, and works well whether or not the primers used in the PCR possess 5'-phosphate groups.

- Provided analysis primers allow for downstream colony PCR screening or sequencing
- Ready-to-use kit components include 1 kb control amplicon, linearized cloning vector and single-use competent *E. coli* (NEB #E1202 only)
- Longer shelf life (12 months), as compared to some commercially available products



The PCR Cloning Kit Includes:

- Linearized pMiniT™ 2.0 Vector
- Cloning Mix 1
- Cloning Mix 2
- Amplicon Cloning Control (1 kb)
- Cloning Analysis Forward Primer
- Cloning Analysis Reverse Primer
- NEB 10-beta Competent *E. coli* (Cloning Efficiency) (NEB #E1202 only)
- NEB 10-beta/Stable Outgrowth Medium (NEB #E1202 only)
- pUC19 Control DNA

PHUSION® is a registered trademark of Thermo Fisher Scientific.

USER® Enzyme



Thermolabile USER II Enzyme



USER Enzyme
#M5505S 50 units\$75
#M5505L 250 units\$303

Thermolabile USER II Enzyme
#M5508S 50 units\$98
#M5508L 250 units\$392

- USER Cloning
- Directional RNA-Seq
- NEBNext adaptor cleavage

Description: USER (Uracil-Specific Excision Reagent) Enzyme generates a single nucleotide gap at the location of a uracil. USER Enzyme is a mixture of Uracil DNA glycosylase (UDG) and the DNA glycosylase-lyase Endonuclease VIII. UDG catalyzes the excision of a uracil base, forming an abasic (apyrimidinic) site while leaving the phosphodiester backbone intact. The lyase activity of Endonuclease VIII breaks the phosphodiester backbone at the 3' and 5' sides of the abasic site so that base-free deoxyribose is released.

Thermolabile Uracil DNA Glycosylase (UDG) catalyzes the excision of a uracil base, forming an abasic (apyrimidinic) site while leaving the phosphodiester backbone intact. The lyase activity of Endonuclease III breaks the phosphodiester backbone at the 3' and 5' sides of the abasic site. In addition to generating a different 3'-terminus than USER Enzyme, Thermolabile USER II Enzyme can also be completely heat inactivated after 10 minutes at 65°C.

Reaction Conditions: rCutSmart Reaction Buffer. Incubate at 37°C. Heat Inactivation of Thermolabile USER II Enzyme: 65°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to nick 10 pmol of a 34-mer oligonucleotide duplex containing a single uracil base, in 15 minutes at 37°C in a total reaction volume of 10 µl. Unit assay conditions can be found at www.neb.com.

Concentration: 1,000 units/ml

Note: This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

How does
the NEB PCR
Cloning Kit work?



Q5® Site Directed Mutagenesis Kit (with or without competent cells)



Q5 Site Directed Mutagenesis Kit
#E0554S 10 reactions\$201

Q5 Site Directed Mutagenesis Kit
(without competent cells)
#E0552S 10 reactions\$139

- Robust exponential amplification generates high yields of desired mutations from a wide range of templates.
- Low error rate of Q5 High-Fidelity DNA Polymerase reduces screening time.
- Room temperature reaction setup
- Use of standard primers eliminates need for phosphorylated or purified oligos
- Easy-to-use master mix format

For help with primer design, try
NEBaseChanger at NEBaseChanger.neb.com


Description: The Q5 Site-Directed Mutagenesis Kit allows rapid site-specific mutagenesis of double-stranded plasmid DNA in less than 2 hours. The kit utilizes the robust Q5 Hot Start High-Fidelity DNA Polymerase along with custom mutagenic primers to create substitutions, deletions and insertions in a wide variety of plasmids. Transformation into high-efficiency NEB 5-alpha Competent *E. coli*, provided with (NEB #E0554), ensures robust results with plasmids up to 14 kb in length.

Applications:

- Generation of mutations, insertions or deletions in plasmid DNA

The Q5 Site Directed Mutagenesis Kit Includes:

- Q5 Hot Start High-Fidelity Master Mix (2X)
- KLD Enzyme Mix (10X) and Reaction Buffer (2X)
- Control Primer Mix (10 µM each) and Template DNA (5 ng/µl)
- NEB 5-alpha Competent *E. coli* (High Efficiency), pUC19 Transformation Control Plasmid (5 pg/µl), SOC Outgrowth Medium (NEB #E0554 only)



Q5 Site-Directed Mutagenesis Kit Overview. The first step is an exponential amplification using Q5 Hot Start High-Fidelity DNA Polymerase. The second step is a unique enzyme mix containing a kinase, ligase and DpnI. Together, these enzymes allow for rapid circularization of the PCR product and removal of the template DNA. The last step is a high-efficiency transformation into chemically competent cells.

Quick Blunting™ Kit

#E1201S 20 reactions\$86
#E1201L 100 reactions\$347

Special Offer:

Quick Blunting and Quick Ligation Kits
#E0542S 20 reactions\$164

See page 101 for details on the
Quick Ligation Kit.

- Restriction enzyme-digested DNA is blunted in less than 30 minutes
- Reactions are performed at room temperature in a ready-to-use mix
- Suitable for restriction enzyme-digested DNA, sheared or nebulized DNA or PCR product

Description: The Quick Blunting Kit is used to convert DNA with incompatible 5' or 3' overhangs to 5' phosphorylated, blunt-ended DNA for efficient blunt-end ligation into DNA cloning vectors. DNA is blunted using T4 DNA Polymerase (NEB #M0203) which has both 3' → 5' exonuclease activity and 5' → 3' polymerase activity. T4 Polynucleotide Kinase (NEB #M0201) is included in the enzyme mix for phosphorylation of the 5' ends of blunt-ended DNA for subsequent ligation into a cloning vector. This kit is optimized for blunting up to 5 µg of DNA in a single reaction.

Applications:

- Prepare sheared, nebulized or restriction enzyme-digested DNA for blunt-end ligation into a plasmid, cosmid, fosmid or BAC vector
- Prepare PCR products for efficient blunt-end cloning


The Quick Blunting Kit Includes:

- Blunting Enzyme Mix
- 10X Blunting Buffer
- 1 mM Deoxynucleotide (dNTP) Solution Mix



Heat Inactivation: 70°C for 10 minutes.

Notes: PCR generated DNA must be purified before blunting by using a commercial purification kit, such as Monarch PCR & DNA Cleanup Kit (NEB #T1030) phenol extraction/ethanol precipitation or gel electrophoresis. Restriction enzyme digested DNA can be blunted directly without purification.



Learn more about
the benefits of
the Q5 SDM Kit.

DNA Ligase Selection Chart

NEB offers a variety of ligases for DNA research. Many of these enzymes are recombinant, and all offer the quality and value you have come to expect from our products. While more than one ligase may work for your application, the following selection chart presents our recommendations for optimal performance.

Visit NEBStickTogether.com for more information on DNA Ligases.

	Blunt/TA Ligase Master Mix	Instant Sticky-end Ligase Master Mix	Electro Ligase®	T4 DNA Ligase	Quick Ligation™ Kit	Salt-T4® DNA Ligase	Hi-T4™ DNA Ligase	Immo- bilized T4 DNA Ligase	T3 DNA Ligase	T7 DNA Ligase	HiFi Taq DNA Ligase	E. coli DNA Ligase	Taq DNA Ligase	9°N™ DNA Ligase	SplintR® Ligase
DNA APPLICATIONS															
Ligation of sticky ends	★★	★★★	★★	★★	★★★	★★	★★	★★★	★★	★★	★	★	★	★	
Ligation of blunt ends	★★★	★	★★	★★	★★★	★★	★★	★	★★						
T/A cloning	★★★	★	★★	★★	★★	★★	★★		★	★					
Electroporation			★★★	★★			★★								
Ligation of sticky ends only										★★★					
Repair of nicks in dsDNA	★★	★★	★★	★★★	★★	★★★	★★★		★★	★★	★★	★★	★★	★★	★★
High complexity library cloning	★★	★★	★★	★★★	★★										
Adaptor Ligation	★★★	★★	★★	★	★★			▲	★						
Ligation-Dependent DNA Sequence & SNP Detection (LCR, LDR & related methods)											★★★		★★	★★	
Ligation-Dependent RNA Sequence & SNP Detection				★											★★★
Ligation of adjacent ssDNAs on an RNA Splint															★★★
NGS APPLICATIONS															
NGS Library Prep dsDNA-dsDNA (ligation)	▲			▲					▲						
FEATURES															
Salt tolerance (> 2X that of T4 DNA Ligase)						✓			✓						
Ligation in 15 min. or less	✓	✓		✓	✓	✓	✓		✓	✓	✓		✓	✓	✓
Master Mix Formulation	✓	✓													
Thermostable											✓		✓	✓	
Thermotolerant							✓				✓		✓	✓	
Recombinant	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓

★★★ Optimal, recommended ligase for selected application

★★ Works well for selected application

★ Will perform selected application, but is not recommended

▲ Please consult the specific NGS protocol to determine the optimal enzyme for your needs

Helpful Online Tools:

For help with choosing the best ligase for your experiment, try [NEBcloner](http://NEBcloner.neb.com) at NEBcloner.neb.com

For help with estimating incubation temperature when using thermostable ligases, try the [Thermostable Ligase Reaction Temperature Calculator](http://LigaseCalc.neb.com) at LigaseCalc.neb.com

For help with scientific calculations and conversions, try [NEBioCalculator](http://NEBioCalculator.neb.com) at NEBioCalculator.neb.com

Find an
overview of
ligation.



Substrate-based Ligase Selection Chart

— DNA
 ~ RNA

This chart provides our recommendation for a choice of ligase to use in a reaction, based upon the type of substrate present. DNA, RNA and hybrid substrates are represented and require specific enzymes to achieve the highest-efficiency ligation.

	RECOMMENDED LIGASE	COMMENTS
NICKED DNA/RNA		
5' ~ OH p ~ 3' 3' ~ 5'	T4 RNA Ligase 2	
5' ~ OH p ~ 3' 3' — 5'	T4 RNA Ligase 2	
5' ~ OH p — 3' 3' ~ 5'	T4 RNA Ligase 2	
5' ~ OH p — 3' 3' — 5'	T4 DNA Ligase	
5' — OH p ~ 3' 3' ~ 5'	N/A	No ligase optimized for this activity
5' — OH p ~ 3' 3' — 5'	N/A	No ligase optimized for this activity
5' — OH p — 3' 3' ~ 5'	SplintR Ligase	100 – 1,000-fold higher efficiency than T4 DNA Ligase
5' — OH p — 3' 3' — 5'	T4 DNA Ligase	For high temperatures, we recommend <i>Taq</i> DNA Ligase. For highest fidelity, we recommend HiFi <i>Taq</i> DNA Ligase.
ssDNA/RNA		
5' — OH p — 3'	N/A	See CircLigase™
5' — OH p ~ 3'	N/A	No ligase optimized for this activity
5' ~ OH p — 3'	T4 RNA Ligase 1	Supplement with ATP
5' ~ OH p ~ 3'	T4 RNA Ligase 1	
5' ~ OH App — 3'	T4 RNA Ligase 2 Truncated KQ	
5' ~ OH App ~ 3'	T4 RNA Ligase 2 Truncated KQ	
5' — OH App — 3'	Thermostable 5' App DNA/RNA Ligase	We recommend a Proteinase K cleanup
5' — OH App ~ 3'	Thermostable 5' App DNA/RNA Ligase	We recommend a Proteinase K cleanup
5' — p OH ~ 3'	RtcB Ligase	Supplement with GTP and Mn ²⁺
5' ~ p OH — 3'	RtcB Ligase	Supplement with GTP and Mn ²⁺
5' ~ OH pNp 3'	T4 RNA Ligase 1	
5' ~ OH pdNp 3'	T4 RNA Ligase 1	Reported to work, but ligates inefficiently. Consider pdCp.
5' — OH pNp 3'	T4 RNA Ligase 1	
5' — OH pdNp 3'	T4 RNA Ligase 1	
dsDNA/RNA		
5' — OH p — 3' 3' — p OH — 5'	Blunt T/A Ligase Master Mix	
5' — A-OH p — 3' 3' — p OH-T — 5'	Blunt T/A Ligase Master Mix	
5' — OH p — 3' 3' — p OH — 5'	Quick Ligation Kit or Instant Sticky-end Ligase Master Mix	For ligating ends under high salt conditions, we recommend T3 DNA Ligase. For ligation of cohesive ends ONLY, we recommend T7 DNA Ligase.
5' — OH p — 3' 3' — p OH — 5'	Quick Ligation Kit or Instant Sticky-end Ligase Master Mix	For ligating ends under high salt conditions, we recommend T3 DNA Ligase. For ligation of cohesive ends ONLY, we recommend T7 DNA Ligase.

CIRCLIGASE™ is a trademark of EpiCentre Technologies Corp.

T4 DNA Ligase Products

For help with molar ratio calculations, try
NEBioCalculator at [NEBioCalculator.neb.com](https://www.neb.com/nebiocalculator)

T4 DNA Ligase catalyzes the formation of a phosphodiester bond between juxtaposed 5' phosphate and 3' hydroxyl termini in duplex DNA or RNA. This enzyme will join blunt-end and cohesive-end termini, as well as repair single stranded nicks in duplex DNA, RNA or DNA/RNA hybrids. T4 DNA Ligase is available in a variety of formulations and variants. The table below lists products available from NEB.

PRODUCT	NEB #	FEATURES	REACTION CONDITIONS	SIZE	PRICE
T4 DNA Ligase	M0202S/L/T/M	<ul style="list-style-type: none"> Standalone enzyme ideal for a variety of ligation reactions Ligation can be performed in supplied buffer, or in any of the four restriction endonuclease NEBuffers, or in T4 Polynucleotide Kinase Buffer if supplemented with 1 mM ATP 	1X T4 DNA Ligase Reaction Buffer. Incubate at 16°C. Heat inactivate at 65°C for 10 minutes.	S/L Regular Concentration (400,000 cohesive end units/ml) 20,000/100,000 units T/M High Concentration (2,000,000 cohesive end units/ml) 20,000/100,000 units	\$65/\$260/ \$65/\$260
MASTER MIXES					
Instant Sticky-end Ligase Master Mix	M0370S/L	<ul style="list-style-type: none"> Ready-to-use 2X solution of T4 DNA Ligase and a proprietary ligation enhancer in an optimized reaction buffer Specifically formulated to rapidly ligate cohesive-end (2-4 bp) substrates and improve transformation 	1X Instant Sticky-end Ligase Master Mix with DNA substrates in a 10 µl reaction volume	50/250 rxns	\$98/\$392
Blunt/TA Ligase Master Mix	M0367S/L	<ul style="list-style-type: none"> Ready-to-use 2X solution of T4 DNA Ligase and a proprietary ligation enhancer in an optimized reaction buffer Specifically formulated to improve ligation and transformation of both-blunt-end and single-base overhang substrates 	1X Blunt/TA Ligase Master Mix with DNA substrates in a 10 µl reaction volume	50/250 rxns	\$100/\$400
FORMULATIONS					
Quick Ligation Kit	M2200S/L	<ul style="list-style-type: none"> Ligation of cohesive- or blunt-end DNA fragments in 5 minutes at room temperature (25°C) 	1X Quick Ligation Reaction Buffer. Incubate at room temperature (25°C).	30/150 rxns	\$99/\$388
ElectroLigase	M0369S	<ul style="list-style-type: none"> Combines T4 DNA Ligase and an optimized, ready-to-use 2X reaction buffer containing a proprietary ligation enhancer and no PEG Promotes robust ligation of all types of DNA ends (blunt, sticky, TA) and is suitable for electroporation, without desalting or purification 	1X ElectroLigase Reaction Buffer with DNA substrates and 1 µl ElectroLigase in an 11 µl reaction volume incubated at 25°C	50 rxns	\$123
NEW Immobilized T4 DNA Ligase	M0569S	<ul style="list-style-type: none"> Enzyme is covalently linked to a magnetic bead, and can be removed from a reaction and reused Enables ligated product to be used directly with no heat inactivation step 	1X T4 DNA Ligase Reaction Buffer. Incubate at 25°C. Place on magnet for 3 minutes to remove.	1.1 mg	\$275
VARIANTS					
NEW Hi-T4™ DNA Ligase	M2622S/L	<ul style="list-style-type: none"> Enables ligation with improved thermostability Active in temperatures up to 50°C 	1X T4 DNA Ligase Reaction Buffer. Incubate at 25°C. Heat inactivate at 65°C for 10 minutes.	20,000/100,000 units	\$65/\$261
NEW Salt-T4® DNA Ligase	M0467S/L	<ul style="list-style-type: none"> Enables ligation with improved salt tolerance Active in reactions as high as 300 mM salt with no loss in activity 	1X T4 DNA Ligase Reaction Buffer. Incubate at 25°C. Heat inactivate at 65°C for 10 minutes.	20,000/100,000 units	\$65/\$261

T3 DNA Ligase

#M0317S	100,000 units	\$70
#M0317L	750,000 units	\$281

- Ligation of sticky or blunt ends
- Increased salt tolerance
- Repair of nicks in dsDNA

Description: T3 DNA Ligase is an ATP-dependent ds-DNA ligase from bacteriophage T3. Cohesive ends, blunt ends, and nick sealing can all be efficiently catalyzed by T3 DNA Ligase. Blunt end ligation is enhanced by the addition of PEG 6000 to the reaction. T3 DNA Ligase exhibits a higher tolerance (2-fold) for NaCl in the reaction compared to T4 DNA Ligase, making the enzyme a versatile choice for *in vitro* molecular biology protocols requiring DNA ligase activity.

Reaction Conditions: 1X StickTogether™ DNA Ligase Buffer. Incubate at 25°C.



Unit Definition: One unit is defined as the amount of enzyme required to give 50% ligation of 100 ng HindIII fragments of λ DNA in a total reaction volume of 20 µl in 1 minute at 25°C in 1X StickTogether DNA Ligase Buffer.

Concentration: 3,000,000 units/ml

Notes: ATP is an essential cofactor for the reaction.

T3 DNA Ligase is also active in buffers without PEG 6000, including T4 DNA Ligase Buffer for applications in which PEG 6000 is detrimental. Supplement with 1 mM ATP (final concentration). In these buffers T3 DNA Ligase exhibits ~10-fold reduction in activity. In applications where a high concentration of NaCl needs to be maintained, we suggest using a reaction buffer without PEG 6000.

T7 DNA Ligase

#M0318S	100,000 units	\$70
#M0318L	750,000 units	\$281

- Ligation of sticky ends only
- Repair of nicks in dsDNA

Description: T7 DNA Ligase is an ATP-dependent ds-DNA ligase from bacteriophage T7. It will catalyze the formation of a phosphodiester bond between adjacent 5' phosphate and 3' hydroxyl groups of duplex DNA. Cohesive end ligation and nick sealing can be efficiently catalyzed by T7 DNA Ligase. However, unlike T4 and T3 DNA Ligases, blunt end ligation is not efficiently catalyzed by T7 DNA Ligase, making it a good choice for applications in which blunt and sticky ends of DNA are present but only the sticky ends are to be joined.

Reaction Conditions: 1X StickTogether DNA Ligase Buffer. Incubate at 25°C.



Unit Definition: One unit is defined as the amount of enzyme required to give 50% ligation of 100 ng HindIII fragments of λ DNA in a total reaction volume of 20 µl in 30 minutes at 25°C in 1X StickTogether DNA Ligase Buffer.

Concentration: 3,000,000 units/ml

Notes: ATP is an essential cofactor for the reaction.

T7 DNA Ligase is also active in buffers without PEG 6000, including T4 DNA Ligase Buffer for applications in which PEG 6000 is detrimental. Supplement the reaction with 1 mM ATP (final concentration). Using these buffers, the activity of T7 DNA Ligase is reduced ~10-fold.

E. coli DNA Ligase

#M0205S	200 units	\$63
#M0205L	1,000 units	\$252

- Selective ligation of nicks in dsDNA without significant joining of dsDNA fragments regardless of end type
- cDNA synthesis

Description: *E. coli* DNA Ligase catalyzes the formation of a phosphodiester bond between the 5'-phosphate and 3'-hydroxyl of two adjacent DNA strands in duplex DNA with cohesive ends. It is not appreciably active on blunt-ended substrates. *E. coli* DNA Ligase uses NAD as a cofactor and can be heat-inactivated. *E. coli* DNA Ligase is active at a range of temperatures (4–37°C).

Reaction Conditions: 1X *E. coli* DNA Ligase Reaction Buffer. Optimal ligation occurs at 16°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to give 50% ligation of HindIII



fragments of λ DNA (5' DNA termini concentration of 0.12 µM, 300 µg/ml) in a total reaction volume of 20 µl in 30 minutes at 16°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

Usage Notes:

Requires NAD⁺ (nicotinamide adenine dinucleotide) as a cofactor.

Ligation of blunt-ended fragments is extremely inefficient. For ligation of blunt-ended fragments we recommend the Blunt/TA Ligase Master Mix (NEB #M0367) or the Quick Ligation Kit (NEB #M2200).

HiFi Taq DNA Ligase

#M0647S 50 reactions\$155


- High fidelity, thermostable
- Repair of nicks in dsDNA
- Allele-specific gene detection using ligase-dependent methods, including the Ligase Chain Reaction (LCR) and Ligase Detection Reaction (LDR)
- Ligation of padlock probes

Description: An optimized blend of a thermostable DNA Ligase and a proprietary additive, HiFi Taq DNA Ligase efficiently seals nicks in DNA with unmatched high fidelity. The formation of a phosphodiester bond between juxtaposed 5' phosphate and 3' hydroxyl termini of two adjacent oligonucleotides that are hybridized to a complementary target DNA is enhanced in the improved reaction buffer and mismatch ligation is dramatically reduced. The improved formulation allows higher

resolution discrimination between ligation donors and acceptors, enabling precise detection of SNPs and other allele variants. HiFi Taq DNA Ligase is active at elevated temperatures (37–75°C).

Reaction Conditions: 1X HiFi Taq DNA Ligase Reaction Buffer. Incubate at 25°C.

For help with calculating ligation temp, try our **Thermostable Ligase Reaction Temperature Calculator at LigaseCalc.neb.com**



HiFi Taq DNA Ligase displays increased fidelity. (A) Schematic of multiplexed substrate pools. Each substrate pool contained a single splint with a defined NN at the ligation junction (e.g., AA, AC, AG...) along with all four upstream probes and all four FAM-labeled downstream probes. Each probe that encodes the base at the ligation junction is of unique length allowing for separation and analysis by capillary electrophoresis. A total of 16 substrate pools were prepared, one for each unique splint. (B) Comparison of the ligation fidelity of Ampligase (Epicentre), Taq DNA Ligase and HiFi Taq DNA Ligase. Fidelity measurements were performed using 1 µl of ligase in a 50 µl reaction mixture in the supplied buffers at 1X concentration. Reactions were incubated 30 min at 55°C, using multiplexed substrate pools as outlined in (A). Rows represent a single template sequence, while columns indicate a particular ligation product resulting from a specific pair of probes ligating with the indicated bases at the ligation junction. A dot indicates detection of a product (see legend above). The diagonal from the top left to the bottom right represents Watson-Crick ligation products; all other spaces indicate mismatch ligation products. While Taq DNA Ligase and Ampligase perform similarly under these conditions, with a range of mismatch products detectable, HiFi Taq DNA Ligase shows dramatically fewer mismatch products while maintaining high yields.

Thermus aquaticus (Taq) DNA Ligase

#M0208S 2,000 units\$81
#M0208L 10,000 units\$324

- Thermostable
- Repair of nicks in dsDNA
- Used in Gibson Assembly method
- Allele-specific gene detection using Ligase Detection Reaction and Ligase Chain Reaction

Description: Taq DNA Ligase is a thermostable ligase that catalyzes the formation of a phosphodiester bond between the 5'-phosphate and 3'-hydroxyl termini of two adjacent DNA strands. The strands to be ligated need to be hybridized and accurately paired, with no gap, to a complementary DNA strand; allowing resolution of single nucleotide variants. Taq DNA Ligase uses NAD as a cofactor and is active at elevated temperatures (37–75°C).

Reaction Conditions: 1X Taq DNA Ligase Reaction Buffer. Incubate at 45°C.

For help with calculating ligation temp, try our **Thermostable Ligase Reaction Temperature Calculator at LigaseCalc.neb.com**

Requires NAD⁺ as a cofactor. NAD⁺ is supplied in the 10X Taq DNA Ligase Reaction Buffer; the buffer should be stored at –70°C to extend the half life of the NAD⁺ cofactor.

Unit Definition: (Cohesive End Unit)

One unit is defined as the amount of enzyme required to give 50% ligation of the 12-base pair cohesive ends of 1 µg of BstEII-digested λ DNA in a total reaction volume of 50 µl in 15 minutes at 45°C. Unit assay conditions can be found at www.neb.com.

Concentration: 40,000 units/ml

Notes: Will not ligate short 4-base overlaps (typical of restriction enzyme digests), while it efficiently ligates 12 base pair overlaps.

9°N™ DNA Ligase

#M0238S 2,500 units \$81

- Repair of nicks in DNA while incubating at high temperatures
- Thermostable
- Allele-specific gene detection using Ligase Detection Reaction and Ligase Chain Reaction

Description: 9°N DNA Ligase is a thermostable ligase that catalyzes the formation of a phosphodiester bond between the 5'-phosphate and 3'-hydroxyl of two adjacent DNA strands that are hybridized and accurately paired, with no gap, to a complementary DNA strand. 9°N DNA Ligase uses ATP as a cofactor and it is active at elevated temperatures (45–70°C).

Reaction Conditions: 1X 9°N DNA Ligase Reaction Buffer. Incubate at 45°C.

  45° 

Unit Definition: (Cohesive End Unit)

One unit is defined as the amount of enzyme required to give 50% ligation of the 12-base pair cohesive ends of 1 µg of BstEII-digested λ DNA in a total reaction volume of 50 µl in 15 minutes at 45°C. Unit assay conditions can be found at www.neb.com.

Concentration: 40,000 units/ml

Notes: Will not ligate short 4-base overlaps (typical of restriction enzyme digests), while it efficiently ligates 12 base pair overlaps.

SplintR® Ligase

#M0375S 1,250 units \$98
#M0375L 6,250 units \$439

SplintR Ligase, also known as PBCV-1 DNA Ligase or *Chlorella* virus DNA Ligase, efficiently

  25° 

See page 214 for more information.

catalyzes the ligation of adjacent, single-stranded DNA oligonucleotides splinted by a complementary RNA strand.

T4 Polynucleotide Kinase & T4 Polynucleotide Kinase (3' phosphatase minus)

  37° 

T4 Polynucleotide Kinase
#M0201S 500 units \$57
#M0201L 2,500 units \$228

T4 Polynucleotide Kinase
(3' phosphatase minus)
#M0236S 200 units \$111
#M0236L 1,000 units \$445

- 5' phosphorylation of DNA/RNA for subsequent ligation
- End labeling DNA or RNA for probes and DNA sequencing
- Removal of 3' phosphoryl groups with T4 Polynucleotide Kinase (NEB #M0201)
- T4 PNK (3' phosphatase minus) (NEB #M0236) can be used for the 5' phosphorylation of 3' phosphorylated mononucleotide to generate a substrate (pNp) that can be added to the 3' end of DNA or RNA
- 5' end labeling of 3' phosphorylated oligos with T4 PNK (3' phosphatase minus) (NEB #M0236)

Description: T4 Polynucleotide Kinase catalyzes the transfer and exchange of P_i from the γ position of ATP to the 5' hydroxyl terminus of polynucleotides (double- and single-stranded DNA and RNA), as well as nucleoside 3' monophosphates. The T4 Polynucleotide Kinase (NEB #M0201) also catalyzes the removal of 3' phosphoryl groups from 3' phosphoryl polynucleotides, deoxynucleoside 3' monophosphates and deoxynucleoside 3' diphosphates. The modified version (NEB #M0236) exhibits full kinase activity with no 3' phosphatase activity.

Reaction Conditions: 1X T4 Polynucleotide Kinase Reaction Buffer. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Usage Notes:

Fresh buffer is required for optimal activity (in older buffers, loss of DTT due to oxidation lowers activity).

CTP, GTP, TTP, UTP, dATP or dTTP can be substituted for ATP as a phosphate donor.

Protocols for phosphorylation (radioactive and non-radioactive) of DNA & RNA can be found at www.neb.com.

The efficiencies of blunt and recessed 5' end phosphorylation can be improved by heating to 70°C for 5 minutes, then chilling on ice prior to kinase addition and by adding PEG-8,000 to 5% (w/v).

T4 Polynucleotide Kinase requires ATP for activity, but the supplied reaction buffer does not contain ATP to allow for high specific activity radiolabeling reactions.

Often, a kinase reaction is followed by a ligation reaction. In such cases, the T4 PNK reaction is performed in ligation buffer at 37°C for 30 minutes. The product of this reaction can be used directly in the ligation reaction without a buffer change or heat inactivation UNLESS there is a need to keep other DNA fragments dephosphorylated during ligation. When this is desirable, PNK should be heat inactivated prior to ligation.

Unit Definition: One Richardson unit is defined as the amount of enzyme catalyzing the incorporation of 1 nmol of acid-insoluble [³²P] in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

5-hydroxymethyluridine DNA Kinase



#M0659S 1,000 units\$73


This is an **Enzyme for Innovation (EFI)**. EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit www.neb.com/EnzymesforInnovation to view the full list.

Description: 5-hydroxymethyluridine DNA Kinase (5-HMUDK) transfers the gamma phosphate from ATP to the hydroxymethyl moiety of 5-hydroxymethyluridine in polymeric DNA.

Reaction Conditions: 1X T4 DNA Ligase Reaction Buffer. Incubate at 37°C. Heat inactivation: 80°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 µg of *Bacillus subtilis* bacteriophage SP8 genomic DNA in 30 minutes at 37°C in a total reaction volume of 20 µl against subsequent cleavage by NcoI-HF restriction endonuclease.

Concentration: 20,000 units/ml



5-HMUDK activity.

Phosphatase Selection Chart

	Quick CIP	Recombinant Shrimp Alkaline Phosphatase (rSAP)	Antarctic Phosphatase
FEATURES			
100% heat inactivation	2 minutes/80°C	5 minutes/65°C	2 minutes/80°C
High specific activity	•	•	
Improved stability	•	•	
Works directly in NEBuffers	•	•	•
Requires additive			• (Zn ²⁺)
Quick Protocol	•		

NEW Quick CIP

#M0525S 1,000 reactions\$90
#M0525L 5,000 reactions\$430

- Dephosphorylation of 5' and 3' ends of DNA and RNA
- Dephosphorylation of cloning vector DNA to prevent recircularization during ligation
- Removal of dNTPs in PCR reactions prior to sequencing or SNP analysis
- Dephosphorylation of DNA prior to end-labeling using T4 Polynucleotide Kinase

Description: Quick CIP is a heat-labile recombinant version of calf intestinal alkaline-phosphatase (CIP). Quick CIP nonspecifically catalyzes the dephosphorylation of 5' and 3' ends of DNA and RNA phosphomonoesters. It also hydrolyses ribo- and deoxyribonucleoside triphosphates (NTPs and dNTPs). Quick CIP is useful in many applications that require the dephosphorylation of DNA or RNA ends. In cloning, dephosphorylation prevents re-ligation of linearized plasmid DNA. The enzyme can quickly dephosphorylate 5' protruding, 5' recessed, and blunt ends in just 10 minutes. Quick CIP may also be used to degrade unincorporated dNTPs in PCR reactions to prepare templates for DNA sequencing or SNP analysis.

Quick CIP is completely and irreversibly inactivated by heating it at 80°C for 2 minutes, unlike wild type CIP, which is not heat-inactivatable. This makes removal of Quick CIP prior to ligation or end-labeling unnecessary.



Reaction Conditions: 1X rCutSmart Reaction Buffer. Incubate at 37°C. Heat inactivation: 80°C for 2 minutes.

Unit Definition: One unit is defined as the amount of enzyme that hydrolyzes 1 µmol of *p*-Nitrophenyl Phosphate, PNPP (NEB #P0757) in a total reaction volume of 1 ml in 1 minute at 37°C.

Concentration: 5,000 units/ml

Note: This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

The Quick CIP Includes:

- rCutSmart Buffer
- Quick CIP

Find an overview of dephosphorylation.



Shrimp Alkaline Phosphatase (rSAP)

rCutSmart RR 37° 65°

#M0371S	500 units \$63
#M0371L	2,500 units \$253

- Dephosphorylation of DNA and RNA
- Dephosphorylation of cloning vector DNA to prevent recircularization during ligation
- Removal of dNTPs and pyrophosphate from PCR reactions prior to sequencing or SNP analysis
- Dephosphorylation of DNA prior to 5' end-labeling using T4 Polynucleotide Kinase

Description: Shrimp Alkaline Phosphatase (rSAP) is a heat labile alkaline phosphatase purified from a recombinant source. rSAP is identical to the native enzyme, and contains no affinity tags or other modifications. rSAP nonspecifically catalyzes the dephosphorylation of 5' and 3' ends of DNA and RNA phosphomonoesters. Also, rSAP hydrolyzes ribo-, as well as deoxyribonucleoside triphosphates (NTPs and dNTPs). rSAP is useful in many applications that require the dephosphorylation of DNA or RNA ends. In cloning, dephosphorylation prevents religation of linearized plasmid DNA. rSAP may also be used to degrade unincorporated dNTPs in PCR reactions to prepare templates for DNA sequencing

or SNP analysis. rSAP is completely and irreversibly inactivated by heating at 65°C for 5 minutes, thereby making removal of rSAP prior to ligation or end-labeling unnecessary.

Reaction Conditions: 1X rCutSmart Reaction Buffer. Incubate at 37°C. Heat inactivation: 65°C for 5 minutes.

Unit Definition: One unit is defined as the amount of enzyme that hydrolyzes 1 µmol of *p*-Nitrophenyl Phosphate, PNPP (NEB #P0757) in a total reaction volume of 1 ml in 1 minute at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 1,000 units/ml

Note: This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

Antarctic Phosphatase

RR NEBU 37° 60°

#M0289S	1,000 units \$73
#M0289L	5,000 units \$295

- Dephosphorylation of DNA and RNA
- Dephosphorylation of cloning vector DNA to prevent recircularization during ligation
- Removal of dNTPs and pyrophosphate from PCR reactions prior to sequencing or SNP analysis
- Dephosphorylation of DNA prior to 5' end-labeling using T4 Polynucleotide Kinase

Description: Antarctic Phosphatase catalyzes the dephosphorylation of 5' and 3' ends of DNA and RNA phosphomonoesters. Antarctic Phosphatase also hydrolyzes ribo-, as well as deoxyribonucleoside triphosphates (NTPs and dNTPs). Antarctic Phosphatase is useful in many applications that require the dephosphorylation of DNA or RNA ends. In cloning, dephosphorylation prevents religation of linearized plasmid DNA. The enzyme acts on 5' protruding, 5' recessed and blunt ends. Antarctic Phosphatase may also be used to degrade unincorporated dNTPs in PCR reactions to prepare templates for DNA sequencing. The enzyme is completely and irreversibly inactivated by heating at 70°C for 5 minutes, thereby making removal of Antarctic Phosphatase prior to ligation or end-labeling unnecessary.

Reaction Conditions: 1X Antarctic Phosphatase Reaction Buffer. Incubate at 37°C. Heat inactivation: 80°C for 2 minutes.

Unit Definition: One unit is defined as the amount of enzyme that will dephosphorylate 1 µg of pUC19 vector DNA cut with a restriction enzyme generating 5' recessed ends in 30 minutes at 37°C. Dephosphorylation is defined as > 95% inhibition of recirculation in a self-ligation reaction and is measured by transformation into *E. coli*. Unit assay conditions can be found at www.neb.com.

Concentration: 5,000 units/ml

Jessica serves as a Technician in our Production Group, helping to prepare buffers for purification processes, assisting with enzyme formulation and fermentation as needed. When not at work, Jessica spends time rollerblading and listening to vinyl records.



Inorganic & Thermostable Inorganic Pyrophosphatases



See page 203 for more information.

Pyrophosphatase, Inorganic (*E. coli*)

#M0361S	10 units	\$64
#M0361L	50 units	\$257

Pyrophosphatase Inorganic (Yeast)

#M2403S	10 units	\$69
#M2403L	50 units	\$278

Thermostable Inorganic Pyrophosphatase

#M0296S	250 units	\$74
#M0296L	1,250 units	\$299

NEW

NudC Pyrophosphatase

#M0607S	250 pmol	\$72
---------	----------	-------	------

Apyrase

#M0398S	10 units	\$74
#M0398L	50 units	\$299

- Highly efficient degradation of ATP to AMP
- Removal of deoxynucleotides in DNA pyrosequencing between cycles
- Conversion of 5' triphosphorylated RNA to ligatable monophosphorylated
- Conversion of 5' triphosphorylated RNA to 5' exonuclease *XRN-1* (NEB #M0338) sensitive monophosphorylated RNA
- Supplied at 10-fold higher concentration

Description: Apyrase (recombinant, *E. coli*) is a highly active ATP-diphosphohydrolase that catalyses the sequential hydrolysis of ATP to ADP and ADP to AMP releasing inorganic phosphate. It is a recombinant version of one of several isoforms of apyrase. It can also hydrolyse 5' tri- and diphosphate ribonucleosides and deoxyribonucleosides to their respective 5' monophosphates. Apyrase can catalyse the conversion of 5' triphosphorylated RNA to 5' monophosphorylated RNA by sequential removal of γ and β phosphates.

Reaction Conditions: 1X Apyrase Reaction Buffer. Heat inactivation: 65°C for 20 minutes.



Unit Definition: One unit is defined as the amount of enzyme that catalyses the release of 1 μ mol of inorganic phosphate from ATP (1 mM, NEB #P0756) in 1X Apyrase Reaction Buffer in 1 minute at 30°C in a total reaction volume of 50 μ l. Unit assay conditions can be found at www.neb.com.

Concentration: 500 units/ml

Notes: Apyrase has a higher ratio of activity for ATP:ADP (14:1).

Apyrase is a calcium-activated enzyme. It is approximately 50% active when Mg^{2+} substitutes Ca^{2+} in Apyrase Reaction Buffer.

As a metal-dependent enzyme, Apyrase can be inhibited by EGTA and EDTA.

The activity of Apyrase is approximately 30% in NEBuffers 1.1, 2.1, 3.1 and CutSmart Buffer.

Apyrase does not remove 5' caps from eukaryotic mRNA.

Tte UvrD Helicase

#M1202S	0.5 μ g	\$73
---------	-------------	-------	------

This is an **Enzyme for Innovation (EFI)**. EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit www.neb.com/EnzymesforInnovation to view the full list.


- Unwinds double-stranded DNA
- Thermostable to 65°C
- Reduces non-specific product formation in isothermal amplification (e.g. LAMP)

Description: *Tte* UvrD Helicase is a repair helicase from the thermophilic organism *Thermoanaerobacter tengcongensis*. It is capable of unwinding double-stranded DNA without a requirement for a specific flap or overhang structure. *Tte* UvrD Helicase is active on a wide range of DNA substrates and, along with its thermostability (active to 70°C), *Tte* UvrD Helicase has been demonstrated to be a useful additive for improving specificity of isothermal amplification reactions.



Reaction Conditions: 1X Isothermal Amplification Buffer. Supplement with 1 mM ATP. Heat inactivation: 80°C for 20 minutes.

Concentration: 20 μ g/ml



Tte UvrD Helicase activity.

Properties of Exonucleases and Non-specific Endonucleases

ENZYME	POLARITY	ACTIVITY ON ssDNA		ACTIVITY ON dsDNA ¹ :					PARTIAL DIGESTION TO GENERATE ss EXTENSION ²	PRODUCTS PRODUCED ³	INHIBITION BY PHOSPHOROTHIOATE ⁴	NOTES
		LINEAR	CIRCULAR	LINEAR 5' ext	LINEAR 3' ext	LINEAR BLUNT	NICKED CIRCULAR/ LINEAR	CIRCULAR (SUPER-COILED)				
Exonuclease I (<i>E. coli</i>)	3'→5'	+	—	—	— ¹⁵	— ⁵	—	—	No	dNMP, dinucleotide ⁵	Yes	15, 5, 6
Thermolabile Exonuclease I	3'→5'	+	—	—	— ¹⁵	— ⁵	—	—	No	dNMP, dinucleotide ⁵	Yes	15, 5, 6
Exonuclease T	3'→5'	+	—	—	— ⁷	— ⁵	—	—	No	dNMP, dinucleotide, short oligo	Yes	5, 7
Exonuclease VII	both	+	+ ⁸	— ¹⁷	— ¹⁷	—	—	—	No	short oligos	No	8
RecJ₁	5'→3'	+	—	— ¹⁶	—	— ⁵	—	—	No	dNMP, ssDNA	Yes	5, 15
Mung Bean Nuclease	endonuclease	+	+	—	—	—	—	—	No	dNMP, dsDNA	No	
Nuclease P1	endonuclease	+	+	—	—	—	—	—	No	5' mononucleotides	No	
Exonuclease III (<i>E. coli</i>)	3'→5'	+/- ¹⁶	—	+	+/- ¹⁴	+	+	—	5'	dNMP, ssDNA	Yes	14
T7 Exonuclease	5'→3'	—	—	+/-	+	+	+	—	3'	dNMP, dinucleotide, ssDNA ⁹	Yes	9
Exonuclease V (RecBCD)	both	+	+	+	+	+	—	—	Yes	short oligos	No	
Exonuclease VIII, truncated	5'→3'	+/- ¹⁰	—	+	+	+	—	—	3'	dNMP, ssDNA	No	10
Lambda Exonuclease	5'→3'	+/- ¹⁰	—	+/- ¹¹	+	+	+/- ¹¹	—	3'	dNMP, dinucleotide, ssDNA	Yes	10, 11
Nuclease BAL-31	both ¹²	+	+	+	+	+	+	—	Yes	dsDNA, dNMP	No	12
T5 Exonuclease	5'→3'	+	+	+	+	+	+	—	3'	dNMP to 6-mer	No	
DNase I (RNase-free)	endonuclease	+	+	+	+	+	+	+	NA	dinucleotides, trinucleotides, oligonucleotides, ssDNA, dsDNA	No	
Micrococcal Nuclease	endonuclease	+	+	+	+	+	+	+	NA	diphosphonucleotides, ssDNA, dsDNA 3'-monophosphonucleotides ¹³	No	13

This table is intended to be used as a guide. Not all reported activities and properties for each exonuclease or endonuclease are listed. The amount of enzyme, substrate and time of incubation can have a dramatic effect upon the desired outcome of the experiment.

Table Legend:

- + activity; preferred substrate
- no significant activity
- +/- activity greatly reduced relative to preferred substrate
- NA not applicable
- ss single-stranded
- ds double-stranded
- ext extension
- dNMP deoxyribonucleoside monophosphate

Footnotes:

- The ability to act on short extensions, blunt ends and nicks distinguishes these enzymes; some of these ends are conveniently generated by restriction digestion. The 5' and 3' extensions tested were 4 nt in length.
- Partial digestion of dsDNA by Lambda Exonuclease, T7 Exonuclease and Exonuclease III will produce dsDNA products with ss extensions. Complete digestion produces ssDNA as products.
- Complete hydrolysis of the preferred substrate will generate the listed products.

- To inhibit exonucleases, use of at least 5 phosphorothioate (pt) bonds in a row is recommended. These bonds must be placed at the end of the DNA corresponding to the Polarity of the enzyme; 5' end for 5'→3' nucleases, the 3' end for 3'→5' nucleases, and at both ends if the nucleases cannot initiate at both ends. Endonucleases cannot be inhibited by pt bonds unless the entire sequence has pt bonds between all nucleotides.
- Depending upon the DNA sequence and amount of exonuclease, RecJ₁, Thermolabile Exonuclease I, Exonuclease I and Exonuclease T may remove a few nucleotides from blunt termini.
- Thermolabile Exonuclease I, Exonuclease I release dNMP from ssDNA, except from the last hydrolytic step where a dinucleotide is produced.
- Exonuclease T can be used to make 3' extensions blunt, however, the yield is low.
- Exonuclease VII will not be able to digest circular ssDNA when EDTA is present in the reaction. In the absence of Mg⁺⁺ the enzyme will act as a pure exonuclease.
- It has been reported that the initial first product hydrolyzed from dsDNA by T7 Exonuclease is a dinucleotide. Subsequent hydrolytic cleavage releases dNMP.
- Lambda Exonuclease and Exonuclease VIII, truncated only cut ssDNA if the 5' contains a phosphate.

- Lambda Exonuclease has a strong preference for initiating on dsDNA containing a 5' phosphate. Thus if linear dsDNA has a 5' phosphate at one end and lacks a 5' phosphate on the other end, then Lambda Exonuclease will preferentially degrade the DNA that contains the phosphorylated end.
- BAL-31 Nuclease has been reported as having both ss endonuclease activity as well as 3' to 5' exonuclease activity. Thus any linear DNA is substrate for this enzyme.
- Products of Micrococcal Nuclease degradation have 3' phosphates. Also cuts RNA whereas DNaseI does not.
- Exonuclease III will be inhibited by overhangs > 4 nucleotides.
- RecJ₁ is not suitable for making 5' extensions blunt. Thermolabile Exonuclease I and Exonuclease I are not suitable for making 3' extensions blunt. These enzymes require longer length ssDNA extensions to initiate than those generated by restriction enzymes.
- Exonuclease III exhibits 5-10X less activity on linear ssDNA versus linear dsDNA.
- For information on removing ssDNA extensions from dsDNA see the Blunting Selection chart.

Common Applications for Exonucleases and Endonucleases

APPLICATION	RECOMMENDED ENZYME(S)	NOTES
Removal of 3' overhangs	Quick Blunting Kit	
5' overhang fill in treatment:	Quick Blunting Kit	
Removal of single-stranded primers for nested PCR reactions	Thermolabile Exonuclease I	
Removal of primers post PCR prior to DNA sequencing or SNP detection	Exonuclease I Thermolabile Exonuclease I Exonuclease VII	Quick heat inactivation versus Exonuclease I For 3' chemically modified primers
Mapping positions of introns in genomic DNA	Exonuclease VII	
Removal of primers with or without 3' or 5' terminal phosphorothioate bonds	Exonuclease VII	
Generating ssDNA from linear dsDNA: If 5' → 3' polarity required	Lambda Exonuclease	Strand targeted for removal requires one 5' end with phosphate
If 3' → 5' polarity required	Exonuclease III	Strand targeted for removal requires a 5' overhang, a blunt end, or a 3' overhang (with less than 4 bases)
Preparation of nested deletions in double-stranded DNA	Exonuclease III (<i>E. coli</i>) plus Exonuclease VII	
Site-directed mutagenesis	Exonuclease III (<i>E. coli</i>) T7 Exonuclease	Removes nicked-strand DNA from 3' to 5' Removes nicked-strand DNA from 5' to 3'
Nick-site extension	T7 Exonuclease	
Degradation of denatured DNA from alkaline-based plasmid purification methods for improving DNA cloning	T5 Exonuclease	
Removal of chromosomal/linear DNA in plasmid preparations	T5 Exonuclease Exonuclease V (RecBCD)	Degrades linear ss + dsDNA, nicked DNA Degrades linear ss + dsDNA: PREFERRED as Exo V will save nicked plasmids resulting in higher yields especially for low-copy number plasmid prep
Removal of unligated products (linear dsDNA) from ligated circular double-stranded DNA	T5 Exonuclease Exonuclease V (RecBCD)	Only the un-nicked form of ligated circular double-stranded remains Both nicked and unnicked-form of ligated circular double-stranded DNA remains
Removal of residual gDNA after purification of low copy plasmid	Exonuclease V (RecBCD)	
Removal of contaminated DNA from RNA samples	DNase I	
Conversion of single-stranded DNA or RNA to 5'-mononucleotides	Nuclease P1	
Analysis of base composition, potential damage and modification of nucleic acids	Nuclease P1	
Progressive shortening of both ends of double-stranded DNA	Nuclease BAL-31	
Preparation of double-stranded DNA fragments with 5'-OH and 3'-phosphate	Micrococcal Nuclease	
Degradation of nucleic acids (both DNA and RNA) in crude cell-free extracts	Micrococcal Nuclease	
Preparation of rabbit reticulocyte	Micrococcal Nuclease	
Chromatin Immunoprecipitation (ChIP) analysis	Micrococcal Nuclease	



DNase I (RNase-Free)




See page 217 for more information.

#M0303S	1,000 units	\$70
#M0303L	5,000 units	\$282

DNase I (RNase-free) is a DNA-specific endonuclease that degrades ds- and ss-DNA to release short oligos with 5'-phosphorylated and 3'-hydroxylated ends.

Note: GMP-grade reagent also available. See page 7.

- Removal of contaminating genomic DNA from RNA samples
- Degradation of DNA templates in transcription reactions



Lambda Exonuclease



#M0262S	1,000 units	\$71
#M0262L	5,000 units	\$286

- Conversion of linear dsDNA to single-stranded DNA via preferred activity on 5'-phosphorylated ends


Description: Lambda exonuclease is a highly processive DNA-specific exonuclease that catalyzes the removal of nucleotides from linear or nicked dsDNA in a 5' → 3' direction. The preferred substrate is 5' phosphorylated dsDNA, although it will also degrade ssDNA and non-phosphorylated substrates at a reduced rate.

Unit Definition: One unit is defined as the amount of enzyme required to produce 10 nmol of acid-soluble deoxyribonucleotide from double-stranded substrate in a total reaction volume of 50 µl in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Reaction Conditions: 1X Lambda Exonuclease Reaction Buffer. Incubate at 37°C. Heat inactivation: 75°C for 10 minutes.

Concentration: 5,000 units/ml

Note: 5'-OH ends are digested 10X slower than 5'-PO₄ ends. ssDNA is digested 100X slower than dsDNA.



Exonuclease I (*E. coli*)



#M0293S	3,000 units	\$71
#M0293L	15,000 units	\$286

- Removal of single-stranded primers in PCR reactions prior to Sanger DNA sequencing or SNP analysis
- Removal of single-stranded primers for nested PCR reactions
- Removal of linear ssDNA, leaving behind dsDNA in the sample


Description: Exonuclease I is a DNA-specific exonuclease that catalyzes the removal of nucleotides from linear single-stranded DNA in the 3' to 5' direction.

Unit Definition: One unit is defined as the amount of enzyme that will catalyze the release of 10 nmol of acid-soluble nucleotide in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.


Reaction Conditions: 1X Exonuclease I Reaction Buffer. Incubate at 37°C. Heat inactivation: 80°C for 20 minutes.

Concentration: 20,000 units/ml

Note: Thermolabile version also available (NEB #M0568).



Thermolabile Exonuclease I



#M0568S	3,000 units	\$77
#M0568L	15,000 units	\$312

- Removal of single-stranded primers in PCR reactions prior to DNA sequencing or SNP analysis
- Removal of single-stranded primers for nested PCR reactions
- Removal of linear ssDNA, leaving behind dsDNA in the sample

Description: Thermolabile Exonuclease I is a DNA-specific exonuclease that catalyzes the removal of nucleotides from single-stranded DNA in the 3' to 5' direction, and can be heat inactivated at 80°C in 1 minute.

soluble nucleotide in a total reaction volume of 100 µl in 6 minutes at 37°C. Unit assay conditions can be found at www.neb.com.


Reaction Conditions: 1X NEBuffer r3.1. Incubate at 37°C. Heat inactivation: 80°C for 1 minute.

Concentration: 20,000 units/ml

Unit Definition: One unit is defined as the amount of enzyme that will catalyze the release of 2 nmol of acid-

Note: Try Exo-CIP Rapid PCR Cleanup Kit (NEB #E1050).

This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.



Exonuclease III (*E. coli*)

#M0206S	5,000 units\$63
#M0206L	25,000 units\$252

- Site-directed mutagenesis
- Preparation of ssDNA for dideoxy sequencing
- Preparation of nested deletions in dsDNA

Description: Exonuclease III is a dsDNA-specific exonuclease that catalyzes the removal of nucleotides from linear or nicked dsDNA in the 3' → 5' direction. A limited number of nucleotides are removed during each binding event, resulting in coordinated progressive deletions within the population of DNA molecules.

Initiation occurs at the 3' termini of linear double-stranded DNA with 5' overhangs or blunt ends and 3' overhangs containing less than four bases.

Exonuclease III has also been reported to have RNase H, 3' phosphatase and AP-endonuclease activities.


RR NEB 37° 100

Reaction Conditions: 1X NEBuffer 1. Incubate at 37°C. Heat inactivation: 70°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to produce 1 nmol of acid-soluble total nucleotide in a total reaction volume of 50 µl in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 100,000 units/ml

Usage Notes: Phosphorothioate linkages are not cleaved by Exonuclease III. Unidirectional deletions can also be created by protecting one terminus by incorporation of an α-phosphorothioate-containing nucleotide.



Exonuclease V (RecBCD)

#M0345S	1,000 units\$76
#M0345L	5,000 units\$307

- Degradation of contaminating linear DNA in plasmid samples
- Removal of residual gDNA after purification of low copy plasmid


Description: Exonuclease V, (RecBCD) is a DNA-Specific exonuclease that also acts as an endonuclease on ssDNA. Activity initiates at both the 5' and 3' ends and is processive, generating oligonucleotides. Activity requires ATP and divalent cations. Mg²⁺ is required for the exonuclease activity, while Ca²⁺ inhibits the exonuclease activity and allows dsDNA unwinding (helicase activity).

RR NEB 37° 100

Reaction Conditions: 1X NEBuffer 4. Supplement with 1 mM ATP. Incubate at 37°C. Heat inactivation: 70°C for 30 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to produce 1 nmol of acid-soluble deoxyribonucleotide from double-stranded DNA in 30 minutes at 37°C in a total reaction volume of 50 µl. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml



Exonuclease VII

#M0379S	200 units\$163
#M0379L	1,000 units\$653

- Removal of primers with or without 3' or 5' terminal phosphorothioate bonds
- Mapping positions of introns in genomic DNA
- Removal of ssDNA, leaving behind the dsDNA in a sample

Description: Exonuclease VII, (Exo VII) is a DNA-specific exonuclease that cleaves linear ssDNA in both 5' → 3' and 3' → 5' direction. The preferred substrate is linear ssDNA.

Reaction Conditions: 1X Exonuclease VII Reaction Buffer. Incubate at 37°C. Heat inactivation: 95°C for 10 minutes.

RR NEB 37° 100

Unit Definition: One unit is defined as the amount of enzyme that will catalyze the release of 1 nmol of acid-soluble nucleotide in a total reaction volume of 50 µl in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml



Exonuclease VIII, truncated

#M0545S 1,000 units \$74

- Removal of linear dsDNA, leaving behind circular DNA in the sample


Description: Exonuclease VIII, truncated, is a dsDNA-specific exonuclease. Exonuclease VIII, truncated initiates nucleotide removal from the 5' termini of linear double-stranded DNA in the 5' to 3' direction. The enzyme does not degrade supercoiled dsDNA and circular ssDNA.

Reaction Conditions: 1X NEBuffer 4. Incubate at 37°C. Heat inactivation: 70°C for 15 minutes.



Unit Definition: One unit is defined as the amount of enzyme required to produce 1 nmol of acid-soluble deoxyribonucleotide from double-stranded DNA in a total reaction volume of 50 µl in 30 minutes at 37°C in 1X NEBuffer 4 with 0.15 mM sonicated duplex [³H] DNA. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml



Exonuclease T (RNase T)

#M0265S 250 units \$74
#M0265L 1,250 units \$299

- Removal of 3' overhangs of dsDNA to generate blunt-ends (sequence-dependent)
- Removal of single-stranded primers in PCR reactions prior to DNA sequencing or SNP analysis
- Removal of single-stranded primers for nested PCR reactions
- Removal of ssDNA, leaving behind dsDNA in the sample

Description: Exonuclease T (Exo T), also known as RNase T, is a single-stranded RNA or DNA specific nuclease that requires a free 3' terminus and removes nucleotides in the 3'→5' direction. Exonuclease T can be used to generate blunt ends from RNA or DNA molecules that have 3' extensions.

Reaction Conditions: 1X NEBuffer 4. Incubate at 25°C. Heat inactivation: 65°C for 20 minutes.


Unit Definition: One unit is defined as the amount of enzyme required to produce 0.1 nmol of TCA



soluble nucleotides from 1 nmol of [³H]-labeled polythymidine in a total reaction volume of 100 µl in 30 minutes at 25°C in 1X NEBuffer 4 with 1 nmol [³H]-labeled polythymidine DNA.

Concentration: 5,000 units/ml

Usage Note: Exo T has different activity on RNA vs. DNA. For RNA, 1 unit of Exo T is required to completely digest 1.0 pmol of rA20 under standard reaction conditions.



Thermostable FEN1

#M0645S 1,600 units \$75

This is an **Enzyme for Innovation (EFI)**. EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit www.neb.com/EnzymesforInnovation to view the full list.

Description: Thermostable Flap Endonuclease 1, FEN1, is a thermostable DNA and RNA endonuclease that catalyzes the cleavage of 5' DNA flaps from branched dsDNA substrates, creating a 5' phosphate terminus. FEN1 products can be ligated by DNA ligase to create dsDNA.

Reaction Conditions: 1X ThermoPol Reaction Buffer. Incubate at 65°C.



Unit Definition: One unit is defined as the amount of enzyme required to cleave 10 pmol of 5' flap containing oligonucleotide substrate in a total reaction volume of 10 µl in 10 minutes at 65°C. Unit assay conditions can be found at www.neb.com.

Concentration: 32,000 units/ml



Micrococcal Nuclease

#M0247S 320,000 gel units \$76

- Preparation of double-stranded DNA fragments with 5'-OH and 3'-phosphate
- Studies of chromatin structure
- Degradation of nucleic acids in crude cell-free extracts
- Preparation of rabbit reticulocyte

Description: Micrococcal Nuclease is a DNA and RNA endonuclease that degrades ds- and ss-DNA and RNA. Both DNA and RNA are degraded to 3' phosphomononucleotides and dinucleotides.

Reaction Conditions: 1X Micrococcal Nuclease Reaction Buffer. Supplement with 100 µg/ml BSA. Incubate at 37°C.

Unit Definition: (Agarose Gel Unit) One gel unit is defined as the amount of enzyme required to digest 1 μ g of


RR NEB U BSA 37° No

lambda DNA in 15 minutes at 37°C into molecular DNA fragments (100–400 base pairs) on a 1.2% agarose gel.

Note: 10,000 Gel Units are approximately equal to 1,000 Kunitz Units.

Concentration: 2×10^6 gel units/ml

Notes: 1–5 mM Ca^{2+} is required for activity. The enzyme is active in the pH range 7–10, with optimal activity at pH 9.2, as long as salt concentration is less than 100 mM. Enzyme can be inactivated by addition of excess EGTA.



Mung Bean Nuclease

#M0250S	1,500 units\$67
#M0250L	7,500 units\$270

- Removal of both 3' and 5' single-stranded overhangs from dsDNA to create blunt ends
- Cleavage of ssDNA and RNA
- Cleavage of the single-stranded region in a DNA hairpin
- Mapping of RNA transcripts

Description: Mung Bean Nuclease is a single-strand specific DNA and RNA endonuclease which will degrade single-stranded extensions from the ends of DNA and RNA molecules, leaving blunt, ligatable ends.

Reaction Conditions: 1X Mung Bean Nuclease Reaction Buffer. Incubate at 30°C.


Unit Definition: One unit is defined as the amount of enzyme required to produce 1 μg of acid-soluble total nucleotide in a total reaction volume of 50 μl in 1 minute

NEB U 30° 

at 37°C in 1X Mung Bean Nuclease Reaction Buffer with 0.5 mg/ml denatured Calf Thymus DNA. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

Usage Note: Do not attempt to heat inactivate, DNA will "breathe" before enzyme inactivates, causing undesirable degradation.



Nuclease BAL-31

#M0213S 50 units \$69

- Progressive shortening of dsDNA

Description: Nuclease BAL-31 is a DNA and RNA exonuclease that degrades linear dsDNA from both 5' and 3' termini. The enzyme is also a highly-specific single-stranded endonuclease which cleaves at nicks, gaps and single-stranded regions of duplex DNA and RNA.

Reaction Conditions: 1X Nuclease BAL-31 Reaction Buffer. Incubate at 30°C. Heat inactivation: 65°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to remove 200 base pairs from each end of linearized double-stranded ϕ X174 DNA (40 μ g/ml)


NEB U 30° 65°

in 50 μ l of 1X Nuclease BAL-31 Reaction Buffer in 10 minutes at 30°C.

Concentration: 1,000 units/ml

Usage Notes: Duplex products of the exonuclease are a mixture of blunt and staggered ends. This mixture can be cloned directly, although maximal ligation efficiency requires repairing the staggered ends with a suitable DNA polymerase.

If necessary, the enzyme may be diluted in reaction buffer just prior to use.



Nuclease P1

#M0660S 10,000 units \$52

- Conversion of ssDNA or RNA to 5'-mono-nucleotides
- Analysis of the base composition of nucleic acids
- Studies of the potential damage and modification in DNA

Description: Nuclease P1 (from *P. citrinum*) is a zinc-dependent ssDNA or RNA specific endonuclease which hydrolyzes 3' → 5' phosphodiester bonds in RNA and ssDNA with no base specificity. Nuclease P1 also exhibits 3'-phosphomonoesterase activity.

Although a single-strand specific nuclease, Nuclease P1 does display some activity toward double-stranded DNA (dsDNA) in Nuclease P1 Reaction Buffer. If preferentially degrading single-stranded nucleic acids (ssDNA or RNA) in the presence of dsDNA, we recommend using Nuclease P1 in 1X NEBuffer 1.1, to limit activity on dsDNA while maintaining single-strand nuclease activity.

NEB 37° 15'


Reaction Conditions: 1X Nuclease P1 Reaction Buffer. Incubate at 37°C. Heat inactivation: 75°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to liberate 1.0 µg of acid soluble nucleotides from *Torula* Yeast total RNA per minute at 37°C in 1X Nuclease P1 Reaction Buffer.

Concentration: 100,000 units/ml

Usage Notes: Substrate specificity for Nuclease P1 is as follows: 3' AMP > RNA > ssDNA >> dsDNA.

The rate of hydrolysis of 2' AMP is 3,000-fold less than that of 3' AMP.



RecJ_f

#M0264S 1,000 units \$72
#M0264L 5,000 units \$290

- Degradation of single-stranded DNA from the 5'-end
- Removal of 5' protruding single-stranded termini at the ends of linear dsDNA (blunt ends are not exclusively created)

Description: RecJ_f is a ssDNA-specific exonuclease that catalyzes the removal of nucleotides from linear ssDNA in the 5' → 3' direction. The preferred substrate is dsDNA with 5' single-stranded overhangs > 6 nucleotides long.

DNA substrate containing a 22 base 5' extension results in products that are a mixture of DNA fragments that have blunt-ends, 5' extensions (1–5 nucleotides) and recessed 5' ends (1–8 nucleotides). RecJ_f does not require a 5' phosphate.

RR NEB 2 37° 65'

Reaction Conditions: 1X NEBuffer 2. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to produce 0.05 nmol TCA soluble deoxyribonucleotide in a total reaction volume of 50 µl in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 30,000 units/ml



T5 Exonuclease

#M0663S 1,000 units \$67
#M0663L 5,000 units \$270

- Removal of incomplete ligation products from ligated circular dsDNA
- Degradation of denatured DNA from alkaline-based plasmid purification methods for improved DNA cloning
- Degradation of contaminating linear and nicked DNA in plasmid samples


Description: T5 Exonuclease is a dsDNA-specific exonuclease and ssDNA endonuclease. It initiates at the 5' termini of linear or nicked dsDNA, and cleaves in the 5' → 3' direction. T5 Exonuclease is able to initiate nucleotide removal from the 5' termini or at gaps and nicks of linear or circular dsDNA. However, the enzyme does not degrade supercoiled dsDNA.

RR NEB 4 37° 15'

Reaction Conditions: 1X NEBuffer 4. Incubate at 37°C.

Unit Definition: One unit of T5 Exonuclease is defined as the amount of enzyme required to cause the change of 0.00032 A₂₆₀ nm/min at 37°C in CutSmart Buffer.

Concentration: 10,000 units/ml



T7 Exonuclease



#M0263S	1,000 units\$66
#M0263L	5,000 units\$266


- *Site-directed mutagenesis*
- *Nick-site extension*

Description: T7 Exonuclease is a dsDNA specific exonuclease that catalyzes removal of nucleotides from linear or nicked dsDNA in the 5' to 3' direction. It initiates at the 5' termini or at gaps and nicks of double-stranded DNA. It will degrade both 5' phosphorylated or 5' dephosphorylated DNA. It has also been reported to degrade RNA and DNA from RNA/DNA hybrids in the 5' to 3' direction, but it is unable to degrade either ds- or ssRNA.

Reaction Conditions: 1X NEBuffer 4. Incubate at 25°C.

Unit Definition: One unit is defined as the amount of enzyme required to produce 1 nmol of acid-soluble deoxyribonucleotide from double-stranded DNA in a total reaction volume of 50 µl in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml



Nucleoside Digestion Mix

37°

See page 217 for more information.

#M0649S	50 reactions\$121
---------	--------------	------------

- *Convenient one-step protocol*
- *Digests both DNA and RNA to single nucleosides*
- *Low-glycerol formulation significantly reduces glycerol-induced ion suppression during MS analysis*

The Nucleoside Digestion Mix is a mixture of enzymes that provides a convenient one-step method to generate single nucleosides from DNA or RNA. Optimized for

quantitative analysis by liquid chromatography-mass spectrometry (LC-MS), this reagent eliminates the need for sequential multi-step, time-consuming digestion protocols.

Annie has been with NEB for over 12 years as a member of our Marketing Communications Team. She helps coordinate and ship materials for vendor tabletops and tradeshow and supports our Sales Team. You may have spoken with her directly, as she sometimes works at the front desk!



Properties of DNA Repair Enzymes and Structure-specific Endonucleases

NEB carries an array of reliable DNA repair enzymes, for use in multiple applications.

ENZYME	MAJOR SUBSTRATE ^{1,2}	CLEAVAGE SITE	PRODUCT(S) PRODUCED	TERMINI CREATED FROM CLEAVAGE		MAJOR ACTIVITY	THERMOSTABLE
				5'-terminus	3'-terminus		
APE 1	AP site	1st phosphodiester bond 5' to AP site	1nt gap	dR5P	OH	Endonuclease	
Endo III	AP site, damaged pyrimidines, Tg	N-glycosidic bond, 1st phosphodiester bond 3' to AP site	1nt gap	P	PA	Glycosylase & AP lyase	
Tma Endo III	AP site, damaged pyrimidines, Tg	N-glycosidic bond, 1st phosphodiester bond 3' to AP site	1nt gap	P	PA	Glycosylase & AP lyase	Yes
Endo IV	AP site	1st phosphodiester bond 5' to AP site	1nt gap	dR5P	OH	Endonuclease	
T7 Endo IV	AP site	1st phosphodiester bond 5' to AP site	1nt gap	dR5P	OH	Endonuclease	Yes
Endo V	dl ⁴ , dU, AP site	2nd phosphodiester bond 3' to dl	nick	P	OH	Endonuclease	
Endo VIII	AP site ⁴	phosphodiester bond 3' and 5' to AP site	1nt gap	P	P	AP lyase	
Fpg	8oxoG, oxidized purines	N-glycosidic bond, phosphodiester bond 3' and 5' to AP site	AP site, 1nt gap	P	P	Glycosylase & AP lyase	
hAAG	3mA, 7mG, dl, dX	N-glycosidic bond	AP site	N/A	N/A	Glycosylase	
hSMUG1	dU ⁴ , 5-hmU, 5-hoU, 5fU	N-glycosidic bond	AP site	N/A	N/A	Glycosylase	
RNaseHIII	rN in dsDNA	phosphodiester bond 5' to ribo	nick	P	OH	Endonuclease	
T4 PDG	CPD, AP site	N-glycosidic bond, phosphodiester bond 3' to AP site	AP site, 1nt gap	P ⁶		Glycosylase & AP lyase	
UDG	dU ⁴	N-glycosidic bond	AP site	N/A	N/A	Glycosylase	
Afu UDG	dU ⁴	N-glycosidic bond	AP site	N/A	N/A	Glycosylase	Yes
Antarctic Thermolabile UDG	dU ⁴	N-glycosidic bond	AP site	N/A	N/A	Glycosylase	
T7 Endo I	Cruciforms, mismatches, Holliday junctions, across DNA nicks	phosphodiester bond 5' to structure	nick	P	OH	Endonuclease	
Thermostable FEN1	5' DNA flap ³	phosphodiester bond at base of flap	nick	P	OH (on flap)	Endonuclease	Yes

Table Legend:

3mA	3-methyladenine	dR5P	deoxyribose-5'-phosphate
5fU	5-formyluridine	dU	deoxyuridine
5-hmU	5-hydroxymethyluridine	dX	deoxyxanthosine
5-hoU	5-hydroxyuridine	NA	not applicable
7mG	7-methylguanine	OH	Hydroxyl
8oxoG	8-oxo-7, 8-dihydroguanine	P	Phosphate
AP	apurinic/aprimidinic sites	PA	3' phospho- α , β -unsaturated aldehyde
CPDs	Cyclobutane pyrimidine dimers	rN	ribonucleotides
dl	deoxyinosine	Tg	Thymine Glycol

Footnotes:

1. Activity is on dsDNA unless noted otherwise.
2. Minor activities, substrates, and references can be found at www.neb.com.
3. 5' flaps of 1-40 nt in length have been confirmed substrates.
4. Enzyme has robust activity on ssDNA in addition to dsDNA.
5. Antarctic Thermolabile UDG can be heat inactivated.
6. CPD still covalently attached.



What are endonucleases and their applications?

DNA Repair Enzymes on Damaged and Non-standard Bases

NEB carries an array of endonucleases and glycosylases for Base-excision repair (BER) for use in multiple applications. The following table indicates the level of repair on either double-stranded or single-stranded DNA oligos for various damaged and non-standard bases.

		DOUBLE-STRANDED DNA OLIGOS (34-MERS)										
ENZYME		AP:A	DHT:A	5-hmU:A	5-hmU:G	I:T	6-MeA:T	8-OG:C	8-OG:G	U:A	U:G	THYMINE GLYCOL:A
REPAIR ENZYMES	APE 1	++++	+	-	-	-	-	-	-	-	-	-
	Endo III	++++	+	-	-	-	-	-	-	-	-	+
	Tma Endo III	++++	++	-	-	-	-	+	+	-	-	++
	Endo IV	++++	+	-	-	-	-	-	-	-	-	-
	Tth Endo IV	++++	+	-	-	-	-	-	+	-	-	-
	Endo V*	+++	+	+	+	++++	+	++	+	+	++++	++
	Endo VIII	++++	++	-	-	-	-	-	-	-	-	+++
	Fpg	+	+	-	-	-	-	++++	++++	-	-	+
	hAAG	-	-	-	-	++++	-	-	-	-	-	-
	hNEIL1	++++	++	-	-	-	-	+	+	-	-	++
	hOGG1	++	-	-	-	-	-	++++	+	-	-	-
	T4 PDG	++++	-	-	-	-	-	-	-	-	-	-
	UDG	N/A	-	-	-	-	-	-	-	-	++++	+
Afu UDG	N/A	-	-	-	-	-	-	-	-	++++	+	-
hSMUG1	N/A	-	+++	+++	-	-	-	-	-	++++	++++	-

Standard reaction conditions were used to titer the enzymes with the alternate base.

*Nicks only, does not remove damage

		SINGLE-STRANDED DNA OLIGOS (34-MERS)							
	ENZYME	AP	DHT	5-hmU	I	6-MeA	8-OG	U	THYMINE GLYCOL:A
REPAIR ENZYMES	APE 1	++	–	–	–	–	–	–	–
	Endo III	++	–	–	–	–	–	–	–
	Tma Endo III	++	+	–	–	–	–	–	–
	Endo IV	–	–	–	–	–	–	–	–
	Tth Endo IV	–	–	–	–	–	–	–	–
	Endo V	+	–	–	++++	–	+	–	–
	Endo VIII	+++	–	–	–	–	–	–	–
	Fpg	+	+	–	–	–	+	–	+
	hAAG	–	–	–	+	–	–	–	–
	hNEIL1	+	+	–	–	–	–	–	+
	hOGG1	++	–	–	–	–	+	–	–
	T4 PDG	–	–	–	–	–	–	–	–
	UDG	N/A	–	–	–	–	–	+++	–
Afu UDG	N/A	–	–	–	–	–	+++	–	
hSMUG1	N/A	–	++	–	–	–	+++	–	

Table Legend:

AP	apurinic/aprymidinic site. The AP site is created by treating a uracil containing oligo with UDG.
DHT	5,6-dihydrothymine
5-hmU	5-hydroxymethyluracil
I	inosine
6-MeA	6-methyladenine
8-OG	8-oxoguanine
U	uridine
AP:A	apurinic/aprymidinic site base paired with adenine
DHT:A	5,6 dihydrothymine base paired with an adenine
5-hmU:A	5-hydroxymethyluracil base paired with an adenine
5-hmU:G	5-hydroxymethyluracil base paired with a guanine
I:T	inosine base paired with a thymine
6-MeA:T	6-methyladenine base paired with a thymine
8-OG:C	8-oxoguanine base paired with a cytosine
8-OG:G	8-oxoguanine base paired with a guanine
U:A	uridine paired with an adenine
U:G	uridine paired with a guanine

Activity Level:

++++	100%
+++	50%
++	10% – 25%
+	< 10%
–	no detectable enzyme activity (< 0.7%)
N/A	not applicable

APE 1

#M0282S	1,000 units	\$74
#M0282L	5,000 units	\$299


- Single-cell gel electrophoresis (Comet assay)
- Alkaline elution
- Alkaline unwinding
- Modified nick translation

Description: Human apurinic/aprimidinic (AP) endonuclease, APE 1, also known as HAP 1 or Ref-1, shares homology with *E. coli* Exonuclease III. APE 1 catalyzes the cleavage of DNA phosphodiester backbone at AP sites via hydrolysis leaving a 1 nucleotide gap with 3'-hydroxyl and 5'-deoxyribose phosphate termini. APE 1 has also been reported to have weak DNA 3'-diesterase, 3' to 5' exonuclease and RNase H activities.

Reaction Conditions: 1X NEBuffer 4. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to cleave 20 pmol of a 34-mer oligonucleotide duplex containing a single AP site in a total reaction volume of 10 µl in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml



Endonuclease IV

#M0304S	1,000 units	\$77
#M0304L	5,000 units	\$312


- Single-cell gel electrophoresis (Comet assay)
- Alkaline elution
- Alkaline unwinding

Description: Endonuclease IV can act on several types of oxidative damage in DNA. The enzyme is an apurinic/aprimidinic (AP) endonuclease that catalyzes the cleavage of DNA phosphodiester backbone at AP sites via hydrolysis leaving a 1 nucleotide gap with 3'-hydroxyl and 5'-deoxyribose phosphate (dRP) termini. The enzyme has 3'-diesterase activity which can remove 3' phosphate, 3'-α, β-unsaturated aldehyde, phosphoglycoaldehyde, and other 3' blocking groups.

Reaction Conditions: 1X NEBuffer 3. Incubate at 37°C. Heat inactivation: 85°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to cleave 1 pmol of a 34-mer oligonucleotide duplex containing a single AP site in a total reaction volume of 10 µl in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml



Tth Endonuclease IV

#M0294S	500 units	\$77
---------	-----------------	------


- Thermostable
- Alkaline elution
- Alkaline unwinding

Description: *Tth* Endonuclease IV is a thermostable apurinic/aprimidinic (AP) endonuclease that catalyzes the cleavage of DNA phosphodiester backbone at AP sites via hydrolysis, leaving a 1 nucleotide gap with 3'-hydroxyl and 5'-deoxyribose phosphate (dRP) termini. The enzyme also has a 3'-diesterase activity that can remove 3' phosphate, 3'-α, β-unsaturated aldehyde, phosphoglycoaldehyde, and other 3' blocking groups.

Reaction Conditions: 1X ThermoPol Reaction Buffer. Incubate at 65°C.

Unit Definition: One unit is defined as the amount of enzyme required to cleave 1 pmol of a 60-mer oligonucleotide duplex containing a single AP site in a total reaction volume of 10 µl in 1 hour at 65°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml



Endonuclease III (Nth)

#M0268S 1,000 units\$74

- Single-cell gel electrophoresis (Comet assay)
- Alkaline elution
- Alkaline unwinding

Description: Endonuclease III (Nth) protein from *E. coli* acts both as a *N*-glycosylase and an AP-lyase. The *N*-glycosylase activity releases damaged pyrimidines, including thymine glycol and 5, 6-dihydroxythymine, generating an AP site. The AP lyase activity cleaves an AP site via β -elimination, creating a 1-nucleotide gap with 3'- α , β -unsaturated aldehyde and 5'-phosphate termini.


Some of the damaged bases recognized and removed by Endonuclease III (Nth) include urea, 5, 6 dihydroxythymine, thymine glycol, 5-hydroxy-5 methylhydantoin, uracil glycol, 6-hydroxy-5, 6-dihydrothymine and methyltartronylurea.



Reaction Conditions: 1X Endonuclease III (Nth) Reaction Buffer. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to cleave 1 pmol of a 34-mer oligonucleotide duplex containing a single AP site in a total reaction volume of 10 μ l in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml



Tma Endonuclease III

#M0291S 500 units\$77

- Alkaline elution
- Alkaline unwinding

Description: A thermostable homolog of the *E. coli* Endonuclease III (Nth) that acts as an *N*-glycosylase and an AP-lyase. The *N*-glycosylase activity releases damaged pyrimidines, including thymine glycol and 5,6-dihydroxythymine, generating an AP site. The AP lyase activity cleaves an AP site via β -elimination, creating a 1-nucleotide gap with 3'- α , β -unsaturated aldehyde and 5'-phosphate termini.


Tma Endonuclease III recognizes abasic sites, 5,6 dihydroxythymine and thymine glycol in DNA.



Reaction Conditions: 1X ThermoPol Reaction Buffer. Incubate at 65°C.

Unit Definition: One unit is defined as the amount of enzyme required to cleave 1 pmol of a 60-mer oligonucleotide duplex containing a single AP site in a total reaction volume of 10 μ l in 1 hour at 65°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml



Endonuclease V

#M0305S 250 units\$77

- Cleavage of oligonucleotides containing deoxyinosines
- Mismatch cleavage

Description: Endonuclease V is a repair enzyme found in *E. coli* that recognizes deoxyinosine, a deamination product of deoxyadenosine in DNA. Endonuclease V, often called Deoxyinosine 3' Endonuclease, recognizes DNA containing deoxyinosines (paired or not) on dsDNA, ssDNA with deoxyinosines and, to a lesser degree, DNA containing abasic sites (AP) or urea, base mismatches, insertion/deletion mismatches, hairpin or unpaired loops, flaps and pseudo-Y structures.


Endonuclease V catalyzes cleavage of the second phosphodiester bond 3' to the mismatch of deoxyinosine, leaving a nick with 3'-hydroxyl and 5'-phosphate.



Reaction Conditions: 1X NEBuffer 4. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to cleave 1 pmol of a 34-mer oligonucleotide duplex containing a single deoxyinosine site in a total reaction volume of 10 μ l in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml



Endonuclease VIII

#M0299S	1,000 units	\$77
#M0299L	5,000 units	\$312

- Single-cell gel electrophoresis (Comet assay)
- Alkaline elution
- Alkaline unwinding

Description: Endonuclease VIII acts as both an *N*-glycosylase and an AP-lyase. The *N*-glycosylase activity releases damaged pyrimidines, including thymine glycol and uracil glycol. The AP lyase activity cleaves DNA phosphodiester backbone at AP sites via β and δ -elimination, creating a 1 nucleotide DNA gap with 5' and 3' phosphate termini.

Damaged bases recognized and removed by Endonuclease VIII include urea, 5, 6- dihydroxythymine, thymine glycol, 5-hydroxy-5- methylhydantoin, uracil glycol, 6-hydroxy-5, 6-dihydrothymine and methyltartronylurea. While Endonuclease VIII is similar to Endonuclease III, Endonuclease




VIII has β and δ lyase activity while Endonuclease III has β lyase activity.

Reaction Conditions: 1X Endonuclease VIII Reaction Buffer. Incubate at 37°C. Heat inactivation: 75°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to cleave 1 pmol of a 34 mer oligonucleotide duplex containing a single AP site in a total reaction volume of 10 μ l in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml



Fpg

#M0240S	500 units	\$76
#M0240L	2,500 units	\$307

- Single-cell gel electrophoresis (Comet assay)
- Alkaline elution
- Alkaline unwinding

Description: Fpg (formamidopyrimidine [fapy]-DNA glycosylase), also known as 8-oxoguanine DNA glycosylase, acts both as an *N*-glycosylase and an AP-lyase. *N*-glycosylase activity releases damaged purines, including 2,6-diamino-4-hydroxy-5-formamidopyrimidine (FapyG) and 8-oxo-7,8-dihydroguanine (8oxoG), generating an AP site. The AP lyase activity cleaves an AP site, via β and δ -elimination, creating a 1 nucleotide DNA gap with 5' and 3' phosphate termini.


Some of the damaged bases recognized and removed by Fpg include 7, 8-dihydro-8-oxoguanine (8-oxoguanine), 8-oxoadenine, fapy-guanine, methy-fapy-guanine, fapy-adenine, aflatoxin B1-fapy-guanine, 5-hydroxy-cytosine and 5-hydroxy-uracil.



Reaction Conditions: 1X NEBuffer 1. Supplement with 100 μ g BSA. Incubate at 37°C. Heat inactivation: 60°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to cleave 10 pmol of a 34-mer oligonucleotide duplex containing a single 8-oxoguanine base paired with a cytosine in a total reaction volume of 10 μ l in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 8,000 units/ml



hAAG

#M0313S	500 units	\$77
---------	-----------------	------

- Single-cell gel electrophoresis (Comet Assay)
- Alkaline elution
- Alkaline unwinding


Description: Human Alkyladenine DNA Glycosylase (hAAG) excises alkylated and oxidative DNA damaged sites, including 3-methyladenine, 7-methylguanine, 1,N⁶-ethenoadenine and hypoxanthine. hAAG catalyzes the hydrolysis of the *N*-glycosidic bond to release the damaged base. hAAG is also known as methylpurine DNA glycosylase (MPG) or 3-methyladenine-DNA glycosylase (ANPG).

Reaction Conditions: 1X ThermoPol Reaction Buffer. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.



Unit Definition: One unit is defined as the amount of enzyme required to create an AP site from 1 pmol of a 34-mer oligonucleotide duplex containing a single deoxyinosine site in a total reaction volume of 10 μ l in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml



T4 PDG (T4 Endonuclease V)

#M0308S 2,000 units\$77

- DNA damage studies
- Single-cell gel electrophoresis (Comet assay)

Description: T4 PDG (pyrimidine dimer glycosylase) has both DNA glycosylase and AP-lyase activity. The *N*-glycosylase activity releases cis-syn cyclobutane pyrimidine dimers, including T^{AT}, T^{AC} and C^{AC}, generating an AP site. The AP lyase activity cleaves an AP site via β -elimination, creating a 1 nucleotide DNA gap with 3'- α , β -unsaturated aldehyde and 5'-phosphate termini.

Reaction Conditions: 1X T4 PDG Reaction Buffer. Supplement with BSA. Incubate at 37°C.

RR NEB BSA 37° 100%

Unit Definition: One unit is defined as the amount of enzyme that catalyzes the conversion of 0.5 μ g of UV irradiated supercoiled pUC19 DNA to > 95% nicked plasmid in a total reaction volume of 20 μ l in 30 minutes at 37°C. Nicking is assessed by agarose gel electrophoresis. Irradiated plasmid contains an average of 3–5 pyrimidine dimers. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

Note: Incubation time should be \leq 30 minutes for best results.



Uracil-DNA Glycosylase (UDG)

#M0280S 1,000 units\$76
#M0280L 5,000 units\$307

Companion Product:

Uracil Glycosylase Inhibitor (UGI)

#M0281S 200 units\$74
#M0281L 1,000 units\$299

- Eliminates PCR carry-over contamination
- Release of uracil from ss- or ds- DNA

Description: *E. coli* Uracil-DNA Glycosylase (UDG) is a monofunctional DNA glycosylase that catalyzes the release of uracil from uracil-containing DNA. UDG efficiently hydrolyzes uracil from single-stranded or double-stranded DNA, but not from oligomers (6 or fewer bases).

Reaction Conditions: 1X UDG Reaction Buffer. Incubate at 37°C.


Unit Definition: One unit is defined as the amount of enzyme that catalyzes the release of 60 pmol of uracil

RR NEB 37° 100%

per minute from double-stranded, uracil-containing DNA. Activity is measured by release of [³H]-uracil in a 50 μ l reaction containing 0.2 μ g DNA (10^4 – 10^5 cpm/ μ g) in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 5,000 units/ml

Usage Notes: UDG is active over a broad pH range with an optimum at pH 8.0, does not require divalent cation, and is inhibited by high ionic strength (> 200 mM).



Afu Uracil-DNA Glycosylase (UDG)

#M0279S 200 units\$77

- Thermostable
- Release of uracil from ss- or ds- DNA

Description: A thermostable homolog of the *E. coli* Uracil-DNA Glycosylase (UDG) from *Archaeoglobus fulgidus*. *Afu* UDG catalyzes the release of uracil from uracil-containing DNA. *Afu* UDG efficiently hydrolyzes uracil from ss- or ds-DNA.

Reaction Conditions: 1X ThermoPol II (Mg-free) Reaction Buffer. Incubate at 65°C.


Unit Definition: One unit is defined as the amount of enzyme that catalyzes the release of 60 pmol of uracil per minute from double-stranded uracil-containing

RR NEB 65° 100%

DNA. Activity is measured by release of [³H]-uracil in a 50 μ l reaction containing 0.2 μ g DNA (10^4 – 10^5 cpm/ μ g) in 30 minutes at 65°C. Unit assay conditions can be found at www.neb.com.

Concentration: 2,000 units/ml

Usage Notes: *Afu* UDG retains 50% activity in the presence of 150 mM NaCl. *Afu* UDG retains less than 1% activity after boiling for 30 minutes in standard reaction conditions. Under standard reaction conditions, uracil glycosylase inhibitor (UGI) does not inhibit *Afu* UDG.



Antarctic Thermolabile UDG

#M0372S	100 units	\$78
#M0372L	500 units	\$316


- Eliminates PCR carry-over contamination
- Release of uracil from ss- or ds- DNA
- Component of Thermolabile User II Enzyme (NEB #M5508)

Description: Antarctic Thermolabile UDG (Uracil-DNA Glycosylase) is a monofunctional DNA glycosylase that catalyzes the release of free uracil from uracil-containing ss- or ds-DNA. The resulting abasic sites are susceptible to the hydrolytic cleavage at the elevated temperature and high pH. This enzyme is sensitive to heat and can be rapidly and completely inactivated at temperatures above 50°C.

Reaction Conditions: 1X Standard *Taq* Reaction Buffer. Incubate at 37°C. Heat inactivation: 50°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme that catalyzes the release of 60 pmol of uracil per minute from double-stranded, uracil-containing DNA. Activity is measured by release of [³H]-uracil in a 50 µl Standard *Taq* reaction buffer containing 0.2 µg DNA (10⁴-10⁵ cpm/µg) in 30 minutes at 37°C.

Concentration: 1,000 units/ml



hSMUG1

#M0336S	500 units	\$77
---------	-----------------	------

- Oxidative DNA damage studies
- Single-cell gel electrophoresis (Comet assay)

Description: Human single-strand-selective monofunctional uracil-DNA Glycosylase (hSMUG1) excises deoxyuracil and deoxyuracil-derivatives bearing an oxidized group at C5, such as 5-hydroxyuracil, 5-hydroxymethyluracil and 5-formyluracil in ssDNA and dsDNA. Major substrates include uracil, 5-hydroxyuracil, 5-hydroxymethyluracil, and 5-formyluracil.

Reaction Conditions: 1X NEBuffer 1. Supplement with 100 µg BSA. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to excise 1 pmol of deoxyuracil from a 34 mer oligonucleotide duplex containing a single dU site in a total reaction volume of 10 µl in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 5,000 units/ml

Usage Notes: hSMUG1 has 50% activity on 5-hydroxymethyluracil when compared to uracil. hSMUG1 has 50% activity on ssDNA compared to dsDNA.



T7 Endonuclease I

#M0302S	250 units	\$71
#M0302L	1,250 units	\$284

- Recognition of mismatched DNA
- Resolve four-way junction or branched DNA
- Detection or cleavage of heteroduplex and nicked DNA
- Random cleavage of linear DNA for shotgun cloning
- Key enzyme for genome editing mutation detection
- Also available: EnGen® Mutation Detection Kit (NEB #E3321)

Description: T7 Endonuclease I is a DNA endonuclease that catalyzes the cleavage of DNA mismatches and non-β DNA structures, including Holliday junctions and cruciform, leaving 3'-OH and 5' phosphate. It is best at C mismatches and does not recognize all DNA mismatches, and to a lesser extent cleaves across a nick in dsDNA.

Source: An *E. coli* strain that carries a fusion of maltose binding protein (MBP) and T7 Endo I.

Reaction Conditions: 1X NEBuffer 2. Incubate at 37°C.


Unit Definition: One unit is defined as the amount of enzyme required to convert > 90% of 1 µg of super-

coiled cruciform pUC(AT)* to > 90% linear form in a total reaction volume of 50 µl in 1 hour at 37°C. Unit assay conditions can be found at www.neb.com.

*pUC(AT) is derived from pUC19 with a modification of the polylinker between the EcoRI and PstI sites.

Concentration: 10,000 units/ml

Usage Note: It is important to control the amount of enzyme and the reaction time used for cleavage of a particular substrate. Temps. > 42°C increase nonspecific nuclease activity.



PreCR® Repair Mix

#M0309S	30 reactions	\$164
#M0309L	150 reactions	\$654

The PreCR Repair Mix is an enzyme cocktail formulated to repair damaged template DNA prior to its use in the

polymerase chain reaction (PCR), microarrays or other DNA technologies.



See page 82 for more information.

Cre Recombinase



#M0298S	50 units\$71
#M0298L	250 units\$286
for high (15X) concentration		
#M0298M	250 units\$286

- *Excision of DNA between two loxP sites*
- *Fusion of DNA molecules containing loxP sites*
- *Inversion of DNA between loxP sites*

Description: Cre Recombinase is a Type I topoisomerase from bacteriophage P1 that catalyzes the site-specific recombination of DNA between *loxP* sites. The enzyme requires no energy cofactors, and Cre-mediated recombination quickly reaches equilibrium between substrate and reaction products. The *loxP* recognition element is a 34 base pair (bp) sequence comprised of two 13 bp inverted repeats flanking an 8 bp spacer region which confers directionality. Recombination products depend on the location and relative orientation of the *loxP* sites. Two DNA species containing single *loxP* sites will be fused. DNA between repeated *loxP* sites will be excised in circular form while DNA between opposing *loxP* sites will be inverted with respect to external sequences.

Reaction Conditions: 1X Cre Recombinase Reaction Buffer. Incubate at 37°C. Heat inactivation: 70°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme necessary to produce maximal site-specific recombination of 0.25 µg pLox2+ control DNA in a total reaction volume of 50 µl in 30 minutes at 37°C. Maximal recombination is determined by agarose gel analysis and by transformation of reactions followed by selection on ampicillin plates.

Concentration: 1,000 units/ml and 15,000 units/ml

TelN Protelomerase



#M0651S	250 units\$74
---------	-----------	-----------

This is an **Enzyme for Innovation (EFI)**. EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit www.neb.com/EnzymesforInnovation to view the full list.

Description: TelN Protelomerase, from phage N15, cuts dsDNA at a TelN recognition sequence (56 bp) and leaves covalently closed ends at the site of cleavage.

Reaction Conditions: 1X ThermoPol Reaction Buffer. Incubate at 30°C. Heat inactivation: 75°C for 5 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to cleave 0.5 µg of pMiniT-TelN BsaI-linearized control plasmid (313 fmol of TelN recognition sites) in a total reaction volume of 50 µl in 30 minutes at 30°C in 1X ThermoPol Reaction Buffer.

Concentration: 5,000 units/ml

5'...TATCAGCACACAATTGCCATTATACGCGGTATAATGGACTATTGTGTGCTGATA...3'
3'...ATAGTCGTGTGTTAACGGGTAATATGCGCATATTACCTGATAACACACGACTAT...5'

↓ TelN Protelomerase

5'...TATCAGCACACAATTGCCATTATACGCGGTATAATGGACTATTGTGTGCTGATA...3'
3'...ATAGTCGTGTGTTAACGGGTAATATGCGCATATTACCTGATAACACACGACTAT...5'

Topoisomerase I (*E. coli*)



#M0301S	100 units\$73
#M0301L	500 units\$295

- *Recognition of mismatched DNA*
- *Catalyzes relaxation of negatively-supercoiled DNA*

Description: Topoisomerase I (*E. coli*) catalyzes the relaxation of negatively-supercoiled DNA. Topoisomerase I has also been implicated in knotting and unknotting DNA, and in linking complementary rings of single-stranded DNA into double-stranded rings. The intact holoenzyme is a 97 kDa protein.

Reaction Conditions: 1X rCutSmart Buffer. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme that catalyzes the relaxation of > 95% of 0.5 µg of negatively supercoiled pUC19 RF I DNA in a total reaction volume of 25 µl in 15 minutes at 37°C. DNA supercoiling is assessed by agarose gel electrophoresis in the absence of ethidium bromide.

Concentration: 5,000 units/ml

Note: This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

β-Agarase I



#M0392S	100 units\$82
#M0392L	500 units\$328

- *Extraction of DNA from agarose gels*
- *DNase and RNase free*

Description: β-Agarase cleaves the agarose subunit, unsubstituted neoagarobiose [3,6-anhydro-α-L-galactopyranosyl-1-3-D-galactose] to neoagaro-oligosaccharides. β-Agarase I digests agarose, releasing trapped DNA and producing carbohydrate molecules which can no longer gel. The remaining carbohydrate molecules and β-Agarase I will not, in general, interfere with subsequent DNA manipulations such as restriction endonuclease digestion, ligation and transformation.

Reaction Conditions: 1X β-Agarase I Reaction Buffer. Incubate at 42°C.


Unit Definition: One unit is defined as the amount of enzyme required to digest 200 µl of molten low melting point or NuSieve agarose to nonprecipitable neoagaro-oligosaccharides in 1 hour at 42°C.

Concentration: 1,000 units/ml

Heat Inactivation: Incubation at 95°C for 2 minutes or incubation at 65°C for 15 minutes inactivates 50 units of β-Agarase I. β-Agarase I retains activity for several hours at 40–45°C and is stabilized by the presence of agarose in the reaction.

CpG Methyltransferase (M.SssI)

#M0226S	100 units \$75
#M0226L	500 units\$302
for high (5X) concentration		
#M0226M	500 units\$302



- Blocking restriction enzyme cleavage
- Studying of CpG methylation-dependent gene expression
- Probing sequence-specific contacts within the major groove of DNA
- Altering the physical properties of DNA
- Uniform [³H]-labeling of DNA
- Decreasing the number of RE cut sites, yielding an apparent increase in specificity

Description: The CpG Methyltransferase (M.SssI) methylates all cytosine residues (C⁵) within the double-stranded dinucleotide recognition sequence 5'...CG...3'.

Reaction Conditions: 1X NEBuffer 2 + SAM. Supplement with 160 μM S-adenosylmethionine (supplied). Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.



Note: MgCl₂ is not required as a cofactor. In the presence of Mg²⁺, methylation becomes distributive rather than processive and exhibits topoisomerase activity.


Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μg of λ DNA in a total reaction volume of 20 μl in 1 hour at 37°C against cleavage by BstUI restriction endonuclease.

Concentration: 4,000 units/ml and 20,000 units/ml. Assayed on λ DNA.

See page 321 for more information on CpG methylation.

GpC Methyltransferase (M.CviPI)

#M0227S	200 units \$78
#M0227L	1,000 units\$314



- Blocking restriction enzyme cleavage
- Altering the physical properties of DNA
- Uniform [³H]-labeling of DNA

Description: The GpC Methyltransferase (M.CviPI) methylates all cytosine residues (C⁵) within the double-stranded dinucleotide recognition sequence 5'...GC...3'.

Reaction Conditions: 1X GC Reaction Buffer. Supplement with 160 μM SAM (supplied). Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Note: MgCl₂ is not required as a cofactor.




Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μg of λ DNA in a total reaction volume of 20 μl in 1 hour at 37°C against cleavage by HaeIII restriction endonuclease.

Concentration: 4,000 units/ml. Assayed on λ DNA.

For more information on products to study DNA methylation, see pages 276–289.

AluI Methyltransferase

#M0220S	100 units \$77
---------	-----------	------------



Description: AluI Methyltransferase modifies the cytosine residue (C⁵) in the sequence to the left.

Reaction Conditions: 1X AluI Methyltransferase Reaction Buffer + SAM. Supplement with 80 μM S-adenosylmethionine (supplied). Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.




Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μg λ DNA in 1 hour at 37°C in a total reaction volume of 10 μl against cleavage by AluI restriction endonuclease.

Concentration: 5,000 units/ml

BamHI Methyltransferase

#M0223S	100 units \$77
---------	-----------	------------



Description: BamHI Methyltransferase modifies the internal cytosine residue (N⁴) in the sequence to the left.

Reaction Conditions: 1X BamHI Methyltransferase Reaction Buffer + SAM. Supplement with 80 μM S-adenosylmethionine (supplied). Incubate at 37°C.




Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μg λ DNA in 1 hour at 37°C in a total reaction volume of 10 μl against cleavage by BamHI restriction endonuclease.

Concentration: 4,000 units/ml

dam Methyltransferase

RR NEB U SAM 37° 65°

#M0222S 500 units\$79
 #M0222L 2,500 units\$318



Description: *dam* Methyltransferase modifies the adenine residue (N⁶) in the sequence to the left.

Reaction Conditions: 1X *dam* Methyltransferase Reaction Buffer + SAM. Supplement with 80 μM S-adenosylmethionine (supplied). Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.


Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μg λ DNA in 1 hour at 37°C in a total reaction volume of 10 μl against cleavage by MboI restriction endonuclease.

Concentration: 8,000 units/ml

EcoGII Methyltransferase

rCutSmart RR NEB U SAM 37° 65°

#M0603S 200 units\$76



This is an **Enzyme for Innovation (EFI)**. EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit www.neb.com/EnzymesforInnovation to view the full list.

Description: EcoGII Methyltransferase is a non-specific methyltransferase that modifies adenine residues (N⁶) in any sequence context.

Reaction Conditions: 1X rCutSmart Buffer + SAM. Supplement with 160 μM S-adenosylmethionine (supplied). Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to protect 100 ng FAM-labeled dsDNA in 30 minutes at 37°C in a total reaction volume of 20 μl against cleavage by MboI restriction endonuclease.

Concentration: 5,000 units/ml


Note: SAM is unstable at 37°C (pH 7.5) and should be replenished in reactions incubated > 4 hours. For use of methylation reaction SAM should be diluted 1:200 to a final concentration of 160 μM. EcoGII Methyltransferase is sensitive to salt. Make sure DNA solution is low in salt concentration or that it makes up only a small % of the final reaction volume.

This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

EcoRI Methyltransferase

RR NEB U SAM 37° 65°

#M0211S 10,000 units\$69



Description: EcoRI Methyltransferase modifies the internal adenine residue (N⁶) in the sequence to the left.

Reaction Conditions: 1X EcoRI Methyltransferase Reaction Buffer + SAM. Supplement with 80 μM S-adenosylmethionine (supplied). Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μg λ DNA in 1 hour at 37°C

in a total reaction volume of 10 μl against cleavage by EcoRI restriction endonuclease.

Concentration: 40,000 units/ml


Note: EcoRI Methyltransferase is inhibited by MgCl₂. Only 50% activity is retained at a concentration of 4 mM MgCl₂.



HaeIII Methyltransferase

RR NEB U SAM 37° 65°

#M0224S 500 units\$77



Description: HaeIII Methyltransferase modifies the internal cytosine residue (C⁵) in the sequence to the left.

Reaction Conditions: 1X HaeIII Methyltransferase Reaction Buffer + SAM. Supplement with 80 μM S-adenosylmethionine (supplied). Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μg λ DNA in 1 hour at 37°C in a total reaction volume of 10 μl against cleavage by HaeIII restriction endonuclease.


Concentration: 10,000 units/ml

Note: HaeIII Methyltransferase protects DNA against cleavage by NotI.

HhaI Methyltransferase

rCutSmart RR NEB U SAM 37° 65°

#M0217S 1,000 units\$77



Description: HhaI Methyltransferase modifies the internal cytosine residue (C⁵) in the sequence to the left.

Reaction Conditions: 1X rCutSmart Buffer + SAM. Supplement with 80 μM S-adenosylmethionine (supplied). Incubate at 37°C.

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μg λ DNA in 1 hour at 37°C in a total reaction volume of 10 μl against cleavage by HhaI restriction endonuclease.


Concentration: 25,000 units/ml

Note: This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

HpaII Methyltransferase

rCutSmart  RR SAM 37° 

#M0214S 100 units\$72



Description: HpaII Methyltransferase recognizes the same sequence as the MspI Methyltransferase, but modifies the internal cytosine residue (C⁵) in the sequence to the left.

Reaction Conditions: 1X rCutSmart Buffer + SAM. Supplement with 80 μM S-adenosylmethionine (supplied). Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μg λ DNA in 1 hour at 37°C in a total reaction volume of 10 μl against cleavage by HpaII restriction endonuclease.


Concentration: 4,000 units/ml

Note: This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

MspI Methyltransferase

 RR NEB U SAM 37° 

#M0215S 100 units\$77



Description: MspI Methyltransferase recognizes the same sequence as the HpaII Methyltransferase, but modifies the external cytosine residue (C⁵) in the sequence to the left.

Reaction Conditions: 1X MspI Methyltransferase Reaction Buffer + SAM. Supplement with 80 μM S-adenosylmethionine (supplied). Incubate at 37°C.


Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μg λ DNA in 1 hour at 37°C in a total reaction volume of 10 μl against cleavage by MspI restriction endonuclease.

Concentration: 5,000 units/ml

TaqI Methyltransferase

rCutSmart  RR SAM 65° 

#M0219S 1,000 units\$77



Description: TaqI Methyltransferase modifies the adenine residue (N⁶) in the sequence to the left.

Reaction Conditions: 1X rCutSmart Buffer + SAM. Supplement with 80 μM S-adenosylmethionine (supplied). Incubate at 65°C.

Note: Activity at 37°C is 25%.

Unit Definition: One unit is defined as the amount of enzyme required to protect 1 μg λ DNA in 1 hour at 65°C in a total reaction volume of 20 μl against cleavage by TaqI restriction endonuclease.

Concentration: 10,000 units/ml

Note: This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

Human DNA (cytosine-5) Methyltransferase (Dnmt1) RR Epi NEB U SAM BSA 37°

See page 286 for more information.

#M0230S 50 units\$130
#M0230L 250 units\$520

Dnmt1 methylates cytosine residues in hemimethylated DNA at 5'...CG...3'. Mammalian Dnmt1 is believed to be involved in carcinogenesis, embryonic development

and several other biological functions. The bulk of the methylation takes place during DNA replication in the S-phase of the cell cycle.

For the past 4 years, Greg has worked as a Fermentation Scientist, helping to produce various enzymes and proteins. When not at work, he enjoys reading, traveling and attending sporting events.



RecA & RecA_f



RecA		
#M0249S	200 µg\$73
#M0249L	1,000 µg\$294
RecA _f		
#M0355S	200 µg\$73

- Visualization of DNA structures with electron microscopy
- D-loop mutagenesis
- Screening libraries using RecA-coated probes
- Cleavage of DNA at a single predetermined site
- RecA mediated affinity capture for full length cDNA cloning

Description: *E. coli* RecA is necessary for genetic recombination, reactions involving DNA repair and UV-induced mutagenesis. RecA promotes the autodigestion of the LexA repressor, umuD protein and lambda repressor. Cleavage of LexA derepresses more than 20 genes. *In vitro* studies indicate that in the presence of ATP, RecA promotes the strand exchange of single-strand DNA fragments with homologous duplex DNA. The reaction has three distinct steps: (i) RecA polymerizes on the single-strand DNA, (ii) the nucleoprotein filament binds the duplex DNA and searches for a homologous region, (iii) the strands are exchanged.

RecA_f is a N-terminal 6X His tagged recombinant protein.

Reaction Conditions: 1X RecA Reaction Buffer. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Molecular Weight: RecA: 37,973 daltons.
RecA_f: 38,907 daltons.

Concentration: 2 mg/ml

T4 Gene 32 Protein



#M0300S	100 µg\$80
#M0300L	500 µg\$324

- Increase yield and processivity of reverse transcriptase during RT-PCR
- Increase yield and specificity of PCR products from soil samples
- Stabilization and marking of ssDNA structures

Description: T4 Gene 32 Protein is a single-stranded DNA (ssDNA) binding protein required for bacteriophage T4 replication and repair. It cooperatively binds to and stabilizes transiently formed regions of ssDNA and plays an important structural role during T4 phage replication. It also has been used extensively to stabilize and mark regions of ssDNA for electron microscopic examination of intracellular DNA structures. Recently, it has been shown to improve restriction enzyme digestion, improve the yield and efficiency of reverse transcription (RT) reactions during RT-PCR, enhance T4 DNA polymerase activity, as well as increase the yield of PCR products.

Reaction Conditions: 1X NEBuffer 4. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Unit Definition: Sold by mass of pure protein as determined by OD₂₈₀ (A₂₈₀=1.184 at 1 mg/ml, 1 cm)

Molecular Weight: 33,485 daltons.

Concentration: 10 mg/ml

ET SSB (Extreme Thermostable Single-Stranded Binding Protein)



#M2401S	50 µg\$165
---------	-------	------------

- Improve the processivity of DNA polymerase
- Stabilization and marking of ssDNA structure
- Increase the yield and specificity of PCR
- Increase the yield and processivity of RT during RT-PCR
- Improve DNA sequencing through regions with strong secondary structure

Description: ET SSB (Extreme Thermostable Single-Stranded DNA Binding Protein) is a single-stranded DNA binding protein isolated from a hyperthermophilic microorganism, which remains fully active after incubation at 95°C for 60 minutes. Due to the extreme thermostability, ET SSB can be used in applications that require extremely high temperature conditions, such as nucleic acid amplification and sequencing.

Unit Definition: Sold by mass of pure protein as determined by OD₂₈₀ (A₂₈₀ = 0.774 at 1 mg/ml, 1 cm).

Concentration: 500 µg/ml

Molecular Weight: 16 kDa

Usage Note: ET SSB is active in any polymerase buffer. Add 200 ng of ET SSB per 50 µl reaction.

Cloning Plasmids and DNAs

CLONING PLASMID/DNA	NEB #	FEATURES	CONCENTRATION	MW/SIZE	SIZE	PRICE
pBR322 Vector	N3033S/L	• Commonly used cloning vectors	1,000 µg/ml	2.83 x 10 ⁶ Da/4,361 bp	50/250 µg	\$78/\$310
pUC19 Vector	N3041S/L	• Amp resistance	1,000 µg/ml	1.75 x 10 ⁶ Da/2,686 bp	50/250 µg	\$76/\$308
M13mp18 RF I DNA	N4018S	• Phage vectors derived from bacteriophage M13 • DNA, covalently closed circular • 13 Unique RE sites with β-gal gene • Blue/white selection	100 µg/ml	7,249 bp	10 µg	\$89
M13mp18 Single-stranded DNA	N4040S		250 µg/ml	7,249 bp	10 µg	\$43
Lambda DNA	N3011S/L	• Commonly used DNA substrate	500 µg/ml	31.5 x 10 ⁶ Da/48,502 bp	250/1,250 µg	\$71/\$280
Lambda DNA (N⁺-methyladenine-free)	N3013S/L		500 µg/ml	31.5 x 10 ⁶ Da/48,502 bp	250/1,250 µg	\$72/\$282
φX174 RF I DNA	N3021S/L	• Commonly used DNA substrate • Covalently closed circular form of φX174	1,000 µg/ml	3.5 x 10 ⁶ Da/5,386 bp	30/150 µg	\$77/\$305
φX174 RF II DNA	N3022L	• Commonly used DNA substrate • Double-stranded nicked circular form of φX174	1,000 µg/ml	3.5 x 10 ⁶ Da/5,386 bp	150 µg	\$310
φX174 Virion DNA	N3023S/L	• Single-stranded viral DNA	1,000 µg/ml	1.7 x 10 ⁶ Da/5,386 bp	50/250 µg	\$86/\$330

NEB offers a selection of common cloning plasmids and DNAs for use as substrates.

Additional information for many of these DNAs can be found in the technical reference section of this catalog or at www.neb.com.

M13KO7 Helper Phage

#N0315S 1.8 ml\$114

- Production of single-stranded phagemid DNA for sequencing and mutagenesis

Description: M13KO7 is a derivative of M13 phage with the origin of replication from P15A and the kanamycin resistance gene from Tn903 both inserted within the M13 origin of replication, which is able to replicate in the absence of phagemid DNA. In the presence of phagemid bearing a wild-type M13 or f1 origin, single-stranded phagemid is packaged preferentially and secreted into the culture medium. This allows easy production of single-stranded phagemid DNA for mutagenesis or sequencing. M13KO7 carries the kanamycin resistance marker.

Source: M13KO7 phage supernatant is isolated from infected *E. coli* ER2738 by a standard procedure.

Concentration: 1.0 x 10¹¹ pfu/ml

Note: NEB does not recommend the use of M13KO7 as a cloning vector. For cloning peptides displayed on M13 phage, we recommend the Ph.D.[™] Peptide Display Cloning System (see page 273).

Programmable Nucleases

Site-specific gene modification and highly-specific *in vitro* cutting is enabled by nucleases that can be easily programmed with nucleic acids. In addition to programmed with nucleic acids. In addition to RNA-guided Cas enzymes, *Tth* Argonaute can be programmed with DNA, further expanded the range of available tools.

NEW

Tth Argonaute (TtAgo)

#M0665S 50 pmol\$75

- Short 16-18 nucleotide 5'-phosphorylated ssDNA guides are cost effective and can be phosphorylated with T4 Polynucleotide Kinase
- Guide/target sequence selection is not limited by the requirement of an adjacent sequence motif
- Highly active on ssDNA and most dsDNA substrates (generates a nick in dsDNA substrates), with mild activity on ssRNA substrates

Description: *Thermus thermophilus* argonaute (TtAgo) is a programmable DNA-endonuclease which requires a short 5'-phosphorylated single-stranded DNA guide to target its activity to a specific corresponding sequence on a substrate. TtAgo introduces one break in the phosphodiester backbone of the complementary substrate between positions 10 and 11 of the DNA guide.

Source: *Thermus thermophilus* argonaute (TtAgo) is purified from an *E. coli* strain that carries a cloned gene from the Gram-negative, thermophilic bacterium *Thermus thermophilus* which is expressed as a recombinant N-terminal 6X His-tagged fusion.

Concentration: 1 µM

Note: Visit www.neb.com/M0665 for usage guidelines.

Competent Cell Selection Chart for Cloning

	NEB 5-alpha Competent <i>E. coli</i> (#C2987)	NEB Turbo Competent <i>E. coli</i> (#C2984)	NEB 5-alpha F' / ⁺ Competent <i>E. coli</i> (#C2992)	NEB 10-beta Competent <i>E. coli</i> (#C3019)	<i>dam</i> ⁻ / <i>dcm</i> ⁻ Competent <i>E. coli</i> (#C2925)	NEB Stable Competent <i>E. coli</i> (#C3040)
FEATURES						
Versatile	•					•
Fast growth (< 8 hours)		•				
Toxic gene cloning			•			•
Large plasmid/BAC cloning				•		•
Dam/Dcm-free plasmid growth					•	
Retroviral/lentiviral vector cloning						•
<i>recA</i> ⁻	•		•	•		•
<i>endA</i> ⁻	•	•	•	•	•	•
FORMATS						
Chemically competent	•	•	•	•	•	•
Electrocompetent	•			•		
Subcloning	•					
96-well format*	•			•		
384-well format*	•					
12 x 8-tube strips*	•					

* Other strains are available upon request. For more information, contact custom@neb.com.

Monarch[®] Nucleic Acid Purification Kits

Need ultra HMW DNA? Check out our new kits for HMW DNA extraction, see page 136.

Monarch kits provide fast and reliable purification of high quality DNA and RNA from a variety of sources using best-in-class silica-column technology. DNA and RNA purified with Monarch kits is highly pure and suitable for use in a wide range of applications, including sequencing, cloning, PCR and other enzymatic

manipulations. Monarch kits are all designed with sustainability in mind; wherever possible, they use less plastic and are packaged in responsibly-sourced, recyclable material. For convenience, all Monarch kit components are available separately. Learn more at NEBmonarch.com.

PRODUCT	NEB #	SIZE	PRICE
Monarch Plasmid Miniprep Kit	T1010S/L	50/250 preps	\$85/\$361
Monarch DNA Gel Extraction Kit	T1020S/L	50/250 preps	\$104/\$475
Monarch PCR & DNA Cleanup Kit (5 µg)	T1030S/L	50/250 preps	\$104/\$475
Monarch Genomic DNA Purification Kit	T3010S/L	50/150 preps	\$156/\$402
Monarch Total RNA Miniprep Kit	T2010S	50 preps	\$260
Monarch RNA Cleanup Kit (10 µg)	T2030S/L	10/100 preps	\$54/\$290
Monarch RNA Cleanup Kit (50 µg)	T2040S/L	10/100 preps	\$54/\$290
Monarch RNA Cleanup Kit (500 µg)	T2050S/L	10/100 preps	\$60/\$480
COLUMNS AVAILABLE SEPARATELY			
Monarch Plasmid Miniprep Columns	T1017L	100 columns + tubes	\$90
Monarch DNA Cleanup Columns (5 µg)	T1034L	100 columns + tubes	\$131
Monarch RNA Purification Columns	T2007L	100 columns + tubes	\$155
Monarch gDNA Purification Columns	T3017L	100 columns + tubes	\$132
Monarch RNA Cleanup Columns (10 µg)	T2037L	100 columns + tubes	\$204
Monarch RNA Cleanup Columns (50 µg)	T2047L	100 columns + tubes	\$204
Monarch RNA Cleanup Columns (500 µg)	T2057L	100 columns + tubes	\$346



NEW ENGLAND

BioLabs[®] Inc.





Putting Recycling to Practice

New England Biolabs has always placed stewardship of the environment as one of our highest priorities. We continuously work to promote ecologically sound practices and environmental sustainability, in order to protect our natural resources. This was instilled in us by our Founder, Don Comb, who believed that all companies should share the responsibility of protecting the environment. Don had a passion for ecologically sustainable practices and encouraged employees to find ways to make a difference, both on the NEB campus and in their personal lives.

An undeniable burden that results from life science research is the extensive use of plastics. Its impact on the environment is increasingly at the forefront of the conscientious consumer's mind, particularly as the toxicity to both terrestrial and marine life becomes more evident. Since 1950, nine million tons of plastic have been produced, and it is predicted that only 9% of that colossal amount is recycled, while 12% is incinerated and 79% ends up in landfills. And while researchers are not solely responsible for the overuse of plastic, one cannot deny that we are a major contributor — the increase in precision, speed and convenience makes single-use microfuge and conical tubes, cell culture dishes, pipettes and gloves much more desirable over the alternatives.

Sustainable solutions focusing on reduction can be more effective than relying on the end-user to recycle or reuse. As one example, NEB carefully designed its Monarch® nucleic acid purification kits to reduce the amounts of plastic used. This was achieved by making the walls of the collection tubes thinner and by using appropriately-sized, lighter-weight bottles.

Cold chain shipping has long been a challenge from an environmental standpoint. Many life science reagents require shipment on ice or dry ice, and maintaining shipping temperature is critical. Expanded polystyrene (EPS), commonly referred to as Styrofoam®, has always been the gold standard for cold shipping — it is light, durable, and well-known for its insulative properties. That said, it is difficult to recycle and often ends up in landfills. For over 40 years, NEB has encouraged customers to use the free return label and send their shipping box back to NEB. More recently, we introduced a 100% recyclable ClimaCell® cooler as an alternative to EPS. The packaging liner is made of starch and paper and is designed with thousands of tiny air pockets that maintain the required temperature for the duration of product delivery.

There are other, smaller changes to business practices that can also help. In 1975, NEB printed its first catalog on recycled paper, and uses recycled, Forest Stewardship Council® (FSC®) certified paper for all its printing needs. In 2015, the NEB Catalog became the first carbon neutral-certified catalog produced in the U.S. At NEB headquarters, we manage an extensive recycling program for paper, plastics, EPS, glass, aluminum, batteries and electronics. These are just some examples of how NEB is working to minimize its impact.

In order to continue to address these types of challenges, we realize it is important to communicate and learn from each other. To that end, NEB supports Labconscious (Labconscious.com) — a community of biologists interested in sustainable laboratory practices that will protect nature and enhance lab productivity. Learn more at Labconscious.com.

NEB shipping boxes being prepared for use.
Credit: Mileidy Rodriguez, NEB

STYROFOAM® is a registered trademark of the Dow Chemical Company.
CLIMACELL® is a registered trademark of Temperpack Technologies, Inc.
FOREST STEWARDSHIP COUNCIL® and FSC® are registered trademarks of The Forest Stewardship Council.

Explore the
ClimaCell Cooler.



Nucleic Acid Purification



Time for change.

Nucleic acid purification is an important step in molecular biology workflows and there are many commercially-available solutions from which to choose. Our Research and Development team spent time with customers to better understand what could be done to improve upon current nucleic acid purification kits. This feedback helped us develop our line of Monarch® Nucleic Acid Purification kits, which have been optimized for maximum performance and minimal environmental impact.

Monarch kits are available for DNA and RNA extraction and cleanup, plasmid purification, and gel extraction. This year, we introduced a novel glass bead-based solution for extraction of high molecular weight DNA to support long read sequencing. Our other Monarch kits utilize unique column designs, which enable the isolation of highly-pure nucleic acids, free from contaminants and often in low volumes. Monarch kits are supported by a variety of validated, user-friendly protocols to support multiple workflows and applications.

We know that it can be difficult to be environmentally friendly in the lab, where sterility and convenience are of utmost importance. At times, it may seem that sustainability and benchwork are at odds with each other. Although we can't completely solve this problem, we can make changes to our product design to help move toward the goal of a greener lab, and that's exactly what we did with the design of our Monarch kits. Wherever possible, these kits use less plastic, as well as responsibly-sourced and recyclable packaging. The columns have thinner walls, reducing total plastic usage without affecting performance. All bottles were carefully chosen to minimize plastic usage, and the kit boxes are made from 100% post consumer content and are designed to be reused.

Let's work together to clean up the world of nucleic acid purification, one prep at a time.

Featured Products

- 136** Monarch HMW DNA Extraction Kits
- 135** Monarch Genomic DNA Purification Kit
- 137** Monarch Total RNA Miniprep Kit
- 140** Monarch Plasmid Miniprep Kit

Featured Tools & Resources



Visit [NEBMonarch.com](https://www.nebmonarch.com) to access our full line of products for nucleic acid purification as well as associated protocols, videos, FAQs and troubleshooting guides



GENOMIC DNA PURIFICATION

Products

Monarch Genomic DNA Purification Kit	135
--------------------------------------	-----

Companion Products

Monarch gDNA Purification Columns	135
Monarch Collection Tubes II	135
Monarch gDNA Tissue Lysis Buffer	135
Monarch gDNA Cell Lysis Buffer	135
Monarch gDNA Blood Lysis Buffer	135
Monarch gDNA Binding Buffer	135
Monarch gDNA Wash Buffer	135
Monarch gDNA Elution Buffer	135
Monarch RNase A	135
Proteinase K, Molecular Biology	135

HMW DNA EXTRACTION

Products

Monarch HMW DNA Extraction Kit for Tissue	136
Monarch HMW DNA Extraction Kit for Cells & Blood	136

Companion Products

Monarch Pestle Set	136
Monarch 2 ml Tubes	136
Monarch DNA Capture Beads	136
Monarch Bead Retainers	136
Monarch gDNA Nuclei Prep & Lysis Buffer Pack	136
Monarch RBC Lysis Buffer	136
Monarch gDNA Elution Buffer II	136
Monarch HMW gDNA Tissue Lysis Buffer	136
Monarch Protein Separation Solution	136
Monarch Precipitation Enhancer	136

RNA PURIFICATION

Products

Monarch Total RNA Miniprep Kit	137
Monarch RNA Cleanup Kit (10 µg)	138
Monarch RNA Cleanup Kit (50 µg)	138
Monarch RNA Cleanup Kit (500 µg)	138

Companion Products

Monarch RNA Purification Columns	137
Monarch gDNA Removal Columns	137
Monarch DNA/RNA Protection Reagent	137
Monarch RNA Lysis Buffer	137
Monarch Total RNA Miniprep Enzyme Pack	137
Monarch RNA Priming Buffer	137
Monarch RNA Wash Buffer	137

Monarch RNA Cleanup Columns (10 µg)	138
Monarch RNA Cleanup Columns (50 µg)	138
Monarch RNA Cleanup Columns (500 µg)	138
Monarch Collection Tubes II	137, 138
Monarch RNA Cleanup Binding Buffer	138
Monarch RNA Cleanup Wash Buffer	138
Nuclease-free Water	137, 138

DNA CLEANUP

Products

Monarch DNA Gel Extraction Kit	140
Monarch PCR & DNA Cleanup Kit (5 µg)	139

Companion Products

Monarch DNA Cleanup Columns (5 µg)	139, 140
Monarch Gel Dissolving Buffer	140
Monarch DNA Cleanup Binding Buffer	139
Monarch DNA Wash Buffer	139, 140
Monarch DNA Elution Buffer	139, 140

PLASMID PURIFICATION

Products

Monarch Plasmid Miniprep Kit	140
------------------------------	-----

Companion Products

Monarch Plasmid Miniprep Columns	140
Monarch Plasmid Resuspension Buffer (B1)	140
Monarch Plasmid Lysis Buffer (B2)	140
Monarch Plasmid Neutralization Buffer (B3)	140
Monarch Plasmid Wash Buffer 1	140
Monarch Plasmid Wash Buffer 2	140
Monarch DNA Elution Buffer	140

ACCESSORIES

Monarch Microfuge Tube EcoRack	141
Monarch Pestle Set	136
Monarch Collection Tubes II	137
Monarch 2 ml Tubes	136

Make the right choice and migrate to Monarch®

Workflows for detecting, analyzing, amplifying or manipulating DNA and RNA often require extraction and purification from a biological sample and/or enzymatic reactions. Monarch nucleic acid purification kits provide fast and reliable purification of high-quality DNA and RNA from a variety of sources using best-in-class silica column technology and a novel glass-bead based workflow for HMW DNA extraction. DNA and RNA purified with Monarch Kits is highly-pure and suitable for use in a wide variety of downstream applications including sequencing, cloning, PCR and other enzymatic manipulations. Monarch kits are all designed with sustainability in mind; they use less plastic whenever possible and are packaged in responsibly-sourced, recyclable material. For convenience, all Monarch kit components are also available separately.

Experience exceptional performance and streamlined workflows

- Efficient extraction of high quality DNA and RNA from a variety of samples
- Novel glass bead-based workflow for extraction of HMW DNA
- Simplified DNA and RNA cleanup in low elution volumes
- Enhanced column designs for improved performance
- Fast, user-friendly protocols
- Optimized buffer systems

Choose Monarch Kits for pure value

- Buffers and columns available separately
- No additional shipping or handling charges**
- No hazardous materials fees**
- Competitive pricing

Reduce your impact on the environment

- Less plastic used in product design
- Responsibly sourced and recyclable packaging
- Packaging and protocol cards are printed with water and soy-based inks
- Reusable kit boxes made from post-consumer content

** In the US and select subsidiary locations. Contact your local distributor for shipping policies.

Designed for sustainability – Monarch kits'...



have up to
44%
less plastic

could eliminate
>140
tons of plastic
each year







use
recyclable
packaging
materials

use
boxes made from
100%
post-consumer
paper



* Visit www.NEBMonarchPackaging.com for details.

To learn more, visit NEBMonarch.com



Monarch Genomic DNA Purification Kit

#T3010S	50 preps\$156
#T3010L	150 preps\$402

Companion Products:

Monarch HMW DNA Extraction Kit for Tissue		
#T3060S	5 preps\$75
#T3060L	50 preps\$450
Monarch HMW DNA Extraction Kit for Cells & Blood		
#T3050S	5 preps\$70
#T3050L	50 preps\$400
Monarch gDNA Purification Columns		
#T3017L	100 columns + tubes\$132
Monarch Collection Tubes II		
#T2018L	100 tubes\$29
Monarch gDNA Tissue Lysis Buffer		
#T3011L	34 ml\$41
Monarch gDNA Cell Lysis Buffer		
#T3012L	20 ml\$41
Monarch gDNA Blood Lysis Buffer		
#T3013L	20 ml\$41
Monarch gDNA Binding Buffer		
#T3014L	65 ml\$112
Monarch gDNA Wash Buffer		
#T3015L	60 ml\$35
Monarch gDNA Elution Buffer		
#T3016L	34 ml\$35
Monarch RNase A		
#T3018L	1 ml\$64
Proteinase K, Molecular Biology		
#P8107S	2 ml\$85

- Purify high quality gDNA from a wide variety of sample types (cells, blood, tissues, and more)
- Experience extremely low residual RNA contamination (typically < 1%)
- Isolate high molecular weight gDNA (peak size typically ≥ 50 kb)
- Take advantage of user-friendly protocols with fast and efficient lysis steps
- Additionally, use the kit to clean up genomic DNA
- Enjoy the flexibility to purchase kit components separately

The Monarch Genomic DNA Purification Kit is an excellent complement to the NEBNext Library Preparation products for NGS.


Description: The Monarch Genomic DNA Purification Kit is a comprehensive solution for cell lysis, RNA removal, and purification of intact genomic DNA (gDNA) from a wide variety of biological samples, including cultured cells, blood, and mammalian tissues. Additionally, bacteria and yeast can be processed with extra steps to enhance lysis in these tough-to-lyse samples. Protocols are also included to enable purification from clinically-relevant samples, such as saliva and cheek swabs, as well as rapid cleanup of previously extracted gDNA. Purified gDNA has high quality metrics, including $A_{260/280} > 1.8$ and $A_{260/230} > 2.0$, high DIN scores and minimal residual RNA. The purified gDNA is suitable for downstream applications, such as endpoint PCR, qPCR and library prep for next generation sequencing (NGS). Typical peak size is 50–70 kb, making this kit an excellent choice upstream of long-read sequencing platforms.

The Monarch gDNA Purification Kit Includes:

- gDNA Purification Columns
- Collection Tubes II
- gDNA Tissue Lysis Buffer
- gDNA Cell Lysis Buffer
- gDNA Blood Lysis Buffer
- gDNA Binding Buffer
- gDNA Wash Buffer
- gDNA Elution Buffer
- RNase A
- Proteinase K, Molecular Biology

SPECIFICATIONS	
Input	<ul style="list-style-type: none"> • Cultured mammalian cells: up to 5×10^6 cells • Mammalian whole blood: 100 μl • Tissue: up to 25 mg, depending on tissue type • Bacteria: up to 2×10^9 • Yeast: up to 5×10^7 • Saliva: up to 500 μl • Buccal swabs • Genomic DNA requiring cleanup
Binding Capacity	30 μ g genomic DNA
Yield	Varies depending on sample type*
Genomic DNA Size	Peak size > 50 kb for most sample types; may be lower for saliva and buccal swabs
RNA Content	< 1% (with included RNase A treatment)
Purity	$A_{260/280} \geq 1.8$, $A_{260/230} \geq 2.0$

*See "Guidelines for Choosing Sample Inputs" (page 342) or at www.neb.com/MonarchgDNAInputs.



The Monarch Genomic DNA Purification Kit efficiently purifies high-quality, high molecular weight gDNA from a variety of sample types. 100 ng of genomic DNA from each sample was loaded on a 0.75% agarose gel. gDNA was isolated following the standard protocols for blood, cultured cells and tissue, and the supplemental protocols for buccal swabs, saliva, Gram- and Gram+ bacteria. Starting material used: 1×10^6 HeLa cells, 100 μ l human blood, 10 μ l bird blood, 10 mg frozen tissue powder, 1 buccal swab, 500 μ l saliva and $\sim 1 \times 10^9$ bacterial cells. Lambda DNA-Hind III digest (NEB #N3012) was used as a marker in the last lane (M). Purified gDNA samples were analyzed using a Genomic DNA ScreenTape® on an Agilent Technologies® 4200 TapeStation®. Samples typically yield peak sizes 50–70 kb and DINs of ~9. The cell fractions processed in the buccal swab and saliva preps contain dead cells, as expected, causing a smear like pattern with typical low molecular weight apoptotic bands.

NEW

Monarch HMW DNA Extraction Kits

Monarch HMW DNA Extraction Kit for Tissue
#T3060S 5 preps\$75
#T3060L 50 preps\$450

Monarch HMW DNA Extraction Kit
for Cells & Blood
#T3050S 5 preps\$70
#T3050L 50 preps\$400

Companion Products:

Monarch Pestle Set
#T3000L 100 sets\$115

Monarch 2 ml Tubes
#T3003L 100 tubes\$60

Monarch DNA Capture Beads
#T3005L 200 beads\$60

Monarch Bead Retainers
#T3004L 100 retainers\$29

Monarch gDNA Nuclei Prep & Lysis Buffer Pack
#T3054L 1 pack\$55

Monarch RBC Lysis Buffer
#T3051L 160 ml\$70

Monarch gDNA Elution Buffer II
#T3056L 24 ml\$35

Monarch HMW gDNA Tissue Lysis Buffer
#T3061L 62 ml\$60

Monarch Protein Separation Solution
#T3062L 36 ml\$60


Monarch Precipitation Enhancer
#T3055L 10 ml\$35

Description: Monarch HMW DNA Extraction Kits provide a rapid and reliable process for extracting high molecular weight DNA (HMW DNA) from biological samples including cells, blood, tissue, bacteria and other sample types. Utilizing an optimized process that combines gentle cell lysis with a tunable fragment length generation, followed by precipitation of the extracted DNA onto the surface of large glass beads, the prep proceeds rapidly. DNA size ranges from 50-250 kb for the standard protocol and into the Mb range on several sample types when the lowest agitation speeds are used. Purified DNA is recovered in high yield with excellent purity, including nearly complete removal of RNA. Purified HMW gDNA is easy to dissolve and is suitable for a variety of downstream applications including long-read sequencing (e.g., PacBio® and Oxford Nanopore Technologies®).

The Monarch HMW DNA Extraction Kits Include:


- DNA Capture Beads & Bead Retainers
- 2 ml Tubes & Collection Tubes II
- RNase A
- Proteinase K, Molecular Biology Grade
- RBC Lysis Buffer (NEB #T3050 only)
- gDNA Nuclei Prep & Nuclei Lysis Buffers (NEB #T3050 only)
- Precipitation Enhancer (NEB #T3050 only)
- Protein Separation Solution (NEB #T3060 only)
- Pestles & Pestle Tubes (NEB #T3060 only)
- HMW gDNA Tissue Lysis Buffer (NEB #T3060 only)
- gDNA Wash Buffer & Elution Buffer II

Workflow for cells.



- Fast workflow (cells: 30 min, blood: 60 min, tissue/bacteria: 90 min)
- Extract DNA into the megabase (Mb) size range with cells, blood, soft organ tissues, and bacteria
- Tune DNA size based on agitation speed during lysis
- Achieve best-in-class yields and purity
- Consistently achieve reproducible results
- Effectively remove RNA
- Elute DNA easily and completely

Try the NEBNext® Companion Module for Oxford Nanopore Technologies Ligation Sequencing (NEB #E7180) for library prep after extraction. See page 162 for details.



DNA fragment size is tunable based on agitation speed during lysis. Preps were performed on duplicate aliquots of 1×10^6 HEK 293 cells and 500 μl fresh human blood. Samples were agitated at the indicated speed during the lysis step to control the fragmentation of the DNA. Equal amounts of DNA from the replicates (cells: 500 ng; blood: 650 ng) were resolved by PFGE (1% agarose gel, 6 V/cm, 13°C for 20 hours, switch times ramped from 0.5–94 seconds on a BioRad CHEF-DR III System). Yield and purity ratios of the individual preps are shown in the accompanying tables. Lambda PFG Ladder and Lambda DNA-Hind III Digest (NEB #N0341 and #N3012) were used as molecular weight standards. Yield, purity ratios and DINs of the individual preps are shown in the accompanying tables.

Monarch Total RNA Miniprep Kit

#T2010S 50 preps\$260

Companion Products:

Monarch RNA Purification Columns		
#T2007L	100 columns + tubes\$155
Monarch gDNA Removal Columns		
#T2017L	100 columns + tubes\$128
Monarch Collection Tubes II		
#T2018L	100 tubes\$29
Monarch DNA/RNA Protection Reagent		
#T2011L	56 ml\$95
Monarch RNA Lysis Buffer		
#T2012L	100 ml\$63
Monarch Total RNA Miniprep Enzyme Pack		
#T2019L	1 pack\$65
Monarch RNA Priming Buffer		
#T2013L	56 ml\$63
Monarch RNA Wash Buffer		
#T2014L	50 ml\$63
Nuclease-free Water		
#B1500S	25 ml\$27
#B1500L	100 ml\$63

- Use with a wide variety of sample types
- Purify RNA of all sizes, including miRNA & small RNAs > 20 nucleotides
- Includes DNase I, gDNA removal columns, Proteinase K, and a stabilization reagent
- Efficiently remove contaminating genomic DNA
- Compatible with Qiacube® and KingFisher™ Flex automation platforms
- Protocols available for RNA fractionation and RNA cleanup
- Save money with value pricing for an all-in-one kit

Description: The Monarch Total RNA Miniprep Kit is a comprehensive solution for sample preservation, cell lysis, gDNA removal, and purification of total RNA from a wide variety of biological samples, including cultured cells, blood, and mammalian tissues. Additionally, tough-to-lyse samples, such as bacteria, yeast, and plant, can be processed with additional steps that enhance lysis. Total RNA, including viral RNA, can also be extracted from clinically-relevant sample like saliva, and swabs (buccal, nasopharyngeal, etc.). Cleanup of enzymatic reactions or purification of RNA from TRIzol®-extracted samples is also possible using this kit. Purified RNA has high quality metrics, including A_{260/280} and A_{260/230} ratios > 1.8, high RIN scores, and minimal residual gDNA. Captured RNA ranges in size from full-length rRNAs down to intact miRNAs. Additionally, differential binding conditions allow selective capture or exclusion of the sub-200 nucleotide RNA pool that includes miRNA, 5S rRNA, and tRNA. Purified RNA is suitable for downstream applications, such as RT-qPCR, cDNA synthesis, RNA-seq, Northern blot analysis, etc.

The Monarch Total RNA Miniprep Kit Includes:


- gDNA Removal Columns
- RNA Purification Columns
- Collection Tubes II
- DNA/RNA Protection Reagent (2X)
- RNA Lysis Buffer
- Proteinase K
- Proteinase K Resuspension Buffer
- Proteinase K Reaction Buffer
- DNase I
- DNase I Reaction Buffer
- RNA Priming Buffer
- RNA Wash Buffer (5X)
- Nuclease-free Water

SPECIFICATIONS

Binding Capacity	100 µg RNA
RNA Size	≥ 20 nt
Purity	A _{260/280} and A _{260/230} usually ≥ 1.8
Input Amount	up to 10 ⁷ cells or 50 mg tissue*
Elution Volume	30–100 µl
Yield	varies depending on sample type
Compatible downstream applications	RNA Library prep for NGS, RT-PCR, RT-qPCR, Northern blots

*See "Guidelines for Choosing Sample Input Amounts" (page 344) or at www.neb.com/MonarchRNAInputs.

BIOANALYZER® is a registered trademark of Agilent Technologies, Inc.
KINGFISHER™ is a trademark of ThermoFisher Scientific.
QIAGEN® and RNEASY® are registered trademarks of Qiagen, Inc.
QIACUBE® and TRIzol® are registered trademarks of Molecular Research Center, Inc.
TRINEAN® and DROPSENSE® are registered trademarks of Trinean NV/SA.



Monarch-purified RNA is high-quality and compatible with a wide variety of downstream applications. Total RNA from a broad array of sample types was purified using the Monarch Total RNA Miniprep Kit. Aliquots were run on an Agilent Bioanalyzer® 2100 using the Nano 6000 RNA chip (S. cerevisiae RNA was run using a plant Nano assay). RIN values and O.D. ratios confirm the overall integrity and purity of the RNA. To demonstrate compatibility with downstream applications, samples were subsequently used for RT-PCR (+/- RT) (A) for detection of 4 different RNA species using Protoscript® II Reverse Transcriptase/LongAmp Taq DNA Polymerase, NGS library prep (B) using NEBNext® Ultra™ II RNA Library Prep Kit and RT-qPCR (C) using Luna® One-Step RT-qPCR Reagents.



View tips for successful RNA purification using Monarch Total RNA Miniprep Kit.

Monarch RNA Cleanup Kits

Monarch RNA Cleanup Kit (10 µg)
 #T2030S 10 preps\$54
 #T2030L 100 preps\$290

Monarch RNA Cleanup Kit (50 µg)
 #T2040S 10 preps\$54
 #T2040L 100 preps\$290

Monarch RNA Cleanup Kit (500 µg)
 #T2050S 10 preps\$60
 #T2050L 100 preps\$480

Companion Products:

Monarch RNA Cleanup Columns (10 µg)
 #T2037L 100 columns + tubes\$204

Monarch RNA Cleanup Columns (50 µg)
 #T2047L 100 columns + tubes\$204

Monarch RNA Cleanup Columns (500 µg)
 #T2057L 100 columns + tubes\$346

Monarch Collection Tubes II
 #T2018L 100 tubes\$29

Monarch RNA Cleanup Binding Buffer
 #T2041S 80 ml\$97

Monarch RNA Cleanup Wash Buffer
 #T2042S 40 ml\$35

Nuclease-free Water
 #B1500S 25 ml\$27
 #B1500L 100 ml\$63

- Choose from 3 different binding capacities and flexible elution volumes
- Quickly and easily purify large quantities of high-quality RNA from *in vitro* transcription (IVT) reactions
- Efficiently remove unincorporated nucleotides from your RNA sample
- Compatible with Qiacube® and KingFisher™ Flex automation platforms

Great for RNA cleanup following *in vitro* transcription with HiScribe™ Kits, see pages 185–187.

Description: The Monarch RNA Cleanup Kits provide a fast and simple silica spin column-based solution for RNA cleanup and concentration after any enzymatic reaction (including *in vitro* transcription, DNase I treatment, capping and labeling) and after other purification methods such as phenol/chloroform extraction. These kits can also be used to extract total RNA from cells, saliva and swabs (buccal, nasopharyngeal, etc.). The Monarch RNA Cleanup Kits are available in 3 different binding capacities: 10 µg, 50 µg and 500 µg. Each kit contains unique columns, all designed to prevent buffer retention and ensure no carryover of contaminants, enabling low-volume elution of highly-pure RNA. Following the standard protocol, RNA ≥ 25 nt is purified with this kit; however, a modified protocol is available to enable the binding of RNA as small as 15 nt (including miRNAs).


Applications:

- RNA cleanup and concentration (including from the TRIzol aqueous phase)
- Enzymatic reaction cleanup
- *In vitro* transcription cleanup
- Total RNA extraction from some samples
- RNA gel extraction
- RNA fractionation

The Monarch RNA Cleanup Kits Include:

- RNA Cleanup Columns (10, 50 or 500 µg)
- RNA Cleanup Binding Buffer
- RNA Cleanup Wash Buffer
- Collection Tubes II
- Nuclease-free Water

Monarch RNA Cleanup Kit	NEB #T2030 (10 µg)	NEB #T2040 (50 µg)	NEB #T2050 (500 µg)
Binding Capacity	10 µg	50 µg	500 µg
RNA Size Range	≥ 25 nt (≥ 15 nt with modified protocol)		
Typical Recovery	70–100%		
Elution Volume	6–20 µl	20–50 µl	50–100 µl
Purity	$A_{260/280} > 1.8$ and $A_{260/230} > 1.8$		
Protocol Time	5 minutes of spin and incubation time		10–15 minutes of spin and incubation time



Recovery of RNA from Monarch RNA Cleanup Kits with Varying Elution Volumes. rRNA (10, 50 or 500 µg, respectively of 16S and 23S Ribosomal Standard from *E. coli*, Sigma) was purified using a Monarch RNA Cleanup Kit (10 µg, NEB #T2030) (50 µg, NEB #T2040) (500 µg, NEB #T2050). Nuclease-free water was used to elute the RNA. The percent recovery of the RNA was calculated from the resulting A_{260} as measured using a Trinean® DropSense® 16. ~80% of RNA can be efficiently recovered in 6 µl from the Monarch RNA Cleanup Kit (10 µg, NEB #T2030), 20 µl from the Monarch RNA Cleanup Kit (50 µg, NEB #T2040), and 50 µl from the Monarch RNA Cleanup Kit (500 µg, NEB #T2050).



View tips for successful purification using the Monarch RNA Cleanup Kits.

Monarch PCR & DNA Cleanup Kit (5 µg)

#T1030S	50 preps\$104
#T1030L	250 preps\$475

Companion Products:

Exo-CIP Rapid PCR Cleanup Kit		
#E1050S	100 reactions\$94
#E1050L	400 reactions\$312

Q5 High-Fidelity DNA Polymerase		
#M0491S	100 units\$113
#M0491L	500 units\$450

Monarch DNA Cleanup Columns (5 µg)		
#T1034L	100 columns + tubes\$131

Monarch DNA Cleanup Binding Buffer		
#T1031L	175 ml\$100

Monarch DNA Wash Buffer		
#T1032L	25 ml\$35

Monarch DNA Elution Buffer		
#T1016L	25 ml\$35

Monarch Plasmid Miniprep Kit		
#T1010S	50 preps\$85
#T1010L	250 preps\$361

Monarch DNA Gel Extraction Kit		
#T1020S	50 preps\$104
#T1020L	250 preps\$475

Monarch Genomic DNA Purification Kit		
#T3010S	50 preps\$156
#T3010L	150 preps\$402


- Elute in as little as 6 µl
- Prevent buffer retention and salt carry-over with optimized column design
- Purify small DNA and oligos with a slight protocol modification
- Save time with fast, user-friendly protocol
- Purchase optimized kit formats or buffers & columns separately for your convenience

With the Monarch PCR & DNA Cleanup Kit, you can purify your DNA in as little as 5 minutes.

Prefer an enzymatic cleanup approach? Check out the Exo-CIP Rapid PCR Cleanup Kit (page 87).

SYBR® is a registered trademark of Molecular Probes, Inc.

Description: The Monarch PCR & DNA Cleanup Kit (5 µg) is a rapid and reliable method for the purification and concentration of up to 5 µg of high-quality, double-stranded DNA from enzymatic reactions such as PCR, restriction digestion, ligation and reverse transcription. This method employs a bind/wash/elute workflow with minimal incubation and spin times, resulting in purification in less than 5 minutes. DNA Cleanup Binding Buffer is used to dilute the samples and ensure they are compatible for loading onto the proprietary silica matrix under high salt conditions. The DNA Wash Buffer ensures enzymes, short primers (≤ 40 nt), detergents and other low-molecular weight reaction components (e.g., nucleotides, DMSO, betaine) are removed, thereby allowing low-volume elution of concentrated, high-purity DNA. Eluted DNA is ready for use in restriction digests, DNA sequencing, ligation and other enzymatic manipulations. The unique column design ensures no buffer retention and no carryover of contaminants, allowing elution of sample in volumes as low as 6 µl. A slight protocol modification enables purification of small DNA and oligonucleotides.



Monarch PCR & DNA Cleanup Kit (5 µg) removes low molecular weight primers from dsDNA samples. Three independent amplicons (267 bp, 520 bp, 1003 bp) were spiked with two oligonucleotides (16-mer, 24-mer) to a final concentration of 1 µM. Half of each mix was purified with the Monarch PCR & DNA Cleanup Kit (5 µg) following the included protocol. Equivalent fractions of the original mixture and the eluted material were resolved on a 20% TBE acrylamide gel at 100V for one hour and stained with SYBR® Green II.


Applications:

- PCR cleanup
- Enzymatic reaction cleanup
- cDNA cleanup
- Labeling cleanup
- Plasmid cleanup
- Oligonucleotide cleanup

The Monarch PCR & DNA Cleanup Kit Includes:

- DNA Cleanup Columns (5 µg)
- DNA Cleanup Binding Buffer
- DNA Wash Buffer
- DNA Elution Buffer
- Collection Tubes

SPECIFICATIONS	
Binding Capacity	5 µg
DNA Size Range	~50 bp–25 kb DNA ≥ 15 bp to 25 kb (dsDNA) and DNA ≥ 18 nt to 10 kb (ssDNA) can also be purified using the Oligonucleotide Cleanup Protocol
Elution Volume	≥ 6 µl
Typical Recovery	DNA 50 bp–10 kb 70–90% DNA 11–23 kb 50–70% ssDNA ≥ 18 nt and dsDNA ≥ 15 bp 70–85%
Protocol Time:	5 minutes



Recovery of ssDNA and dsDNA oligonucleotides (1 µg) using the Oligonucleotide Cleanup Protocol and the Monarch PCR & DNA Cleanup Kit. Synthesized ssDNA and dsDNA oligonucleotides (1 µg in 50 µl H₂O) of varying lengths (6–20 nt) were purified using the Oligonucleotide Cleanup Protocol and the Monarch PCR & DNA Cleanup Kit (NEB #T1030) and were eluted in 50 µl water. The average percent recovery (n=3) of the oligonucleotides was calculated from the resulting A260 as measured using a Trinean DropSense™ 16. Use of the Oligonucleotide Cleanup Protocol and the Monarch PCR & DNA Cleanup Kit (NEB #T1030) results in the efficient removal of small oligonucleotides (6–12 nt) and >70% recovery and cleanup of oligonucleotides ≥ 15 bp (dsDNA) or ≥ 18 nt (ssDNA).



View tips for using the Monarch PCR & DNA Cleanup Kit.

Monarch DNA Gel Extraction Kit

#T1020S	50 preps\$104
#T1020L	250 preps\$475

Companion Products:

β-Agarase I		
#M0392S	100 units\$82
#M0392L	500 units\$328

Monarch DNA Cleanup Columns (5 µg)		
#T1034L	100 columns + tubes\$131

Monarch Gel Dissolving Buffer		
#T1021L	235 ml\$102

Monarch DNA Wash Buffer		
#T1032L	25 ml\$35

Monarch DNA Elution Buffer		
#T1016L	25 ml\$35

Monarch Plasmid Miniprep Kit		
#T1010S	50 preps\$85
#T1010L	250 preps\$361

Monarch PCR & DNA Cleanup Kit (5 µg)		
#T1030S	50 preps\$104
#T1030L	250 preps\$475

Monarch Genomic DNA Purification Kit		
#T3010S	50 preps\$156
#T3010L	150 preps\$402

- Elute in as little as 6 µl
- Prevent buffer retention and salt carry-over with optimized column design
- Save time with fast, user-friendly protocol
- Purchase optimized kit formats or buffers & columns separately for your convenience

With Monarch DNA Cleanup Columns, DNA can be eluted in as little as 6 µl.

Description: The Monarch DNA Gel Extraction Kit rapidly and reliably purifies up to 5 µg of concentrated high-quality, double-stranded DNA from agarose gels. This method employs a bind/wash/elute workflow with minimal incubation and spin times, resulting in purification in less than 15 minutes. The Monarch Gel Dissolving Buffer is used to digest the agarose gel slice and ensure the sample is compatible for loading the DNA onto the proprietary silica matrix under high salt conditions. The wash buffer ensures trace amounts of DNA binding dyes, electrophoresis buffer salts and gel loading buffer components are removed. Low-volume elution produces concentrated, highly pure DNA ready for use in restriction digests, DNA sequencing, ligation, and other enzymatic manipulations. The unique column design ensures no buffer retention and no carryover of contaminants, allowing elution of sample in volumes as low as 6 µl.

The Monarch DNA Gel Extraction Kit Includes:


- DNA Cleanup Columns (5 µg)
- Gel Dissolving Buffer
- DNA Wash Buffer
- DNA Elution Buffer
- Collection Tubes

SPECIFICATIONS

Binding Capacity	5 µg
DNA Size Range	50 bp–25 kb
Elution Volume	≥ 6 µl
Typical Recovery	DNA 50 bp–10 kb 70–90% DNA 11–25 kb 50–70%
Protocol Time:	~10 minutes



Monarch DNA Gel Extraction Kit reproducibly recovers DNA over a broad range of molecular weights. A mixture of 7 DNA fragments ranging from 10 kb down to 0.5 kb was prepared and one-half of the mixture was resolved on a 1% gel. Each fragment was manually excised from the agarose gel and processed using the Monarch DNA Gel Extraction Kit. The entire elution of each fragment was resolved on a new gel with the remainder of the original mixture for comparison.



Our optimized column design supplied with the Monarch DNA Gel Extraction and PCR & DNA Cleanup Kits enables elution in as little as 6 µl, and eliminates buffer retention



View tips for using the Monarch Gel Extraction Kit.

Monarch Plasmid Miniprep Kit

#T1010S	50 preps\$85
#T1010L	250 preps\$361

Companion Products:


Monarch Genomic DNA Purification Kit		
#T3010S	50 preps\$156
#T3010L	150 preps\$402
Exonuclease V (RecBCD)		
#M0345S	1,000 units\$76
#M0345L	5,000 units\$307
Monarch Plasmid Miniprep Columns		
#T1017L	100 columns + tubes\$90
Monarch Plasmid Resuspension Buffer (B1)		
#T1011L	55 ml\$36
Monarch Plasmid Lysis Buffer (B2)		
#T1012L	2 x 27 ml\$36
Monarch Plasmid Neutralization Buffer (B3)		
#T1013L	110 ml\$57
Monarch Plasmid Wash Buffer 1		
#T1014L	2 x 27 ml\$35
Monarch Plasmid Wash Buffer 2		
#T1015L	30 ml\$35
Monarch DNA Elution Buffer		
#T1016L	25 ml\$35
Monarch DNA Gel Extraction Kit		
#T1020S	50 preps\$104
#T1020L	250 preps\$475
Monarch PCR & DNA Cleanup Kit (5 µg)		
#T1030S	50 preps\$104
#T1030L	250 preps\$475

- Prevent buffer retention and salt carry-over with optimized column design
- Monitor completeness of certain steps using colored buffer system
- No need to add RNase before starting
- Elute in low volumes
- Purchase optimized kit formats or buffers & columns separately for your convenience
- Easily label columns using tab and frosted surfaces

Description: The Monarch Plasmid Miniprep Kit is a rapid and reliable method for the purification of high-quality plasmid DNA. This method employs standard cell resuspension, alkaline lysis and neutralization steps, with the additional benefit of color indicators at certain steps to easily monitor completion. After clarification of the lysate by centrifugation, the DNA is bound to the proprietary silica matrix under high salt conditions. Unique wash buffers ensure salts, proteins, RNA and other cellular components (endotoxins) are removed, allowing low-volume elution of concentrated, highly pure DNA, ready for use in restriction digests, DNA sequencing, PCR and other enzymatic manipulations.

The Monarch Plasmid Miniprep Kit Includes:

- Plasmid Miniprep Columns
- Plasmid Resuspension Buffer (B1)
- Plasmid Lysis Buffer (B2)
- Plasmid Neutralization Buffer (B3)
- Plasmid Wash Buffer 1
- Plasmid Wash Buffer 2
- DNA Elution Buffer
- Collection Tubes



Monarch Plasmid Miniprep Kits consistently yield more concentrated plasmid DNA with equivalent purity and functionality as the leading supplier. Preps were performed according to recommended protocols using 1.5 ml aliquots of the same overnight culture. One microliter of each prep was digested with HindIII-HF (NEB #R3104) to linearize the vector and the digests were resolved on a 1% w/v agarose gel.



Monarch Microfuge Tube EcoRack


#T5020S	2 racks\$34
---------	---------	-----------

Description: The Monarch Microfuge Tube EcoRack is a bench-top tube rack made from plastic recovered during the manufacture of Monarch Nucleic Acid Purification Columns. Plastic that would otherwise be discarded during the injection molding process is recovered and re-molded into this useful lab accessory that can hold up to 48 tubes each side. One side can accommodate tubes 1.5-2 mls and the other can accommodate 0.5 ml tubes.



View tips for using the Monarch Plasmid Miniprep Kit.





Aim AR camera at this page
to see image in motion.



Neglected Tropical Diseases

Neglected Tropical Diseases (NTDs) is the term given to 17 understudied diseases that are collectively the second leading cause of disability worldwide after mental health issues.

Among the NTDs, the diseases caused by filarial nematodes garner the least attention — they are the most understudied of the neglected diseases, yet they threaten the well-being of about 20% of the world's population. Insect vectors harbor eight species of these parasitic worms and transmit them to humans. They cause a range of devastating diseases such as lymphatic filariasis and onchocerciasis, commonly known as elephantiasis and river blindness, respectively. Because they are found in tropical regions and do not affect more developed countries, very little funding is dedicated to basic research or development of diagnostics, vaccines and drugs to ease suffering.

Over the years, Don Comb traveled extensively and saw firsthand the suffering caused by filariasis. Compelled to find a tangible way to make a difference in these burdened communities, he embarked on an altruistic research venture that has had a tremendous impact over the course of 40 years, and continues to forge new ground in an otherwise overlooked health crisis. Don used what he knew best to make a difference — science. He established the NEB parasitology research program in 1980, funding the initiative through profits from NEB's restriction enzyme sales.

NEB's parasitology group is uniquely positioned to identify the gaps needed to move the field forward. They have published hundreds of research papers shedding light on the biology of filarial parasites and their *Wolbachia* endosymbionts, describing advances in genomics, glycomics, vaccine and drug development, and developing diagnostics that can be used in low-resource settings.

There is a great need for simple tools that empower local agencies to perform diagnosis and surveillance, in an effort to control transmission of these devastating diseases. One example of NEB's contribution is the development of loop-mediated isothermal amplification (LAMP) tests for several filarial diseases. LAMP uses a simplified, single-temperature, robust nucleic acid amplification method that does not require sophisticated laboratory equipment. Successful amplification of filarial nucleic acid is determined by an easily visible, pH-based color change. It is suited for use in the field — no specialized lab equipment is required. Researchers in Africa are currently using LAMP-based tests for various human filarial infections to determine communities at risk, identify people who need treatment and provide information on the impact of the drugs used to control transmission. The NEB parasitology group also played a significant role in developing new antibiotic regimens, in use or in clinical trials, which target the *Wolbachia* endosymbiont of filarial parasites.

Don was a humanitarian at heart and wanted to see change happen where it is needed most. NEB's parasitology research program has bridged the gap between the bench and the field and has made significant contributions that continue to help make a difference to those suffering in neglected communities.



NEBNext® Reagents for Next Generation Sequencing

Leading the way in sample preparation for next generation sequencing.

Library preparation is a critical part of the next generation sequencing workflow; successful sequencing requires the generation of high quality libraries of sufficient yield and quality.

As sequencing technologies continue to improve and capacities expand, the need for high performance sample prep is greater than ever, from decreasing input quantities and samples of lower quality, for an expanding range of applications.

To meet these growing challenges, the NEBNext suite of products continues to evolve to support next generation sequencing with sample preparation tools that expand applications, streamline workflows, minimize inputs, and improve library quality and yields.

NEBNext reagents are available for sample preparation for DNA, RNA, ChIP, FFPE, small RNA, single cell and microbiome samples, for use with Illumina®, Oxford Nanopore Technologies®, Ion Torrent™ and other sequencing platforms. Kits for specialized applications including methylome analysis, virus sequencing and immune repertoire sequencing are also included in the NEBNext line.

Products are in user-friendly formats including kits and modules. A fast-growing range of adaptors and primers is available separately, for maximized flexibility. Use of NEBNext products has been cited in thousands of peer-reviewed publications.

For additional convenience and cost-effectiveness in high-throughput workflows, NEBNext reagents are also available in bulk and customized formats. For more information, contact Custom@neb.com.



Find an
overview of NGS
library preparation.

Featured Products

- 151** NEBNext Enzymatic Methyl-seq
- 157** Customizable RNA Depletion – NEBNext RNA Depletion Core Reagents Set & the online NEBNext Custom RNA Depletion Design Tool
- 155** NEBNext rRNA Depletion for Human/Mouse/Rat v2
- 155** NEBNext rRNA Depletion for Bacteria
- 160** NEBNext Adaptors & Primers
- 149** NEBNext Ultra II FS DNA Library Prep Kits
- 153** NEBNext Ultra II RNA Library Prep Kits

Featured Tools & Resources



Visit NEBNextSelector.neb.com for help with selecting products.



Visit the [NEBNext Custom RNA Depletion Design Tool](#) to obtain custom probe sequences.



Visit www.NEBNext.com to keep up to date on everything NEBNext.

REAGENTS FOR ILLUMINA SEQUENCING

Product Selection Charts

DNA Product Selection Chart	146
RNA Product Selection Chart	147

Ultra II DNA

NEBNext Ultra II DNA Library Prep Kits	
+/- Sample Purification Beads	149
NEBNext Ultra II FS DNA Library Prep Kits	
+/- Sample Purification Beads	149
NEBNext Ultra II DNA & FS DNA PCR-Free Library Prep Kits +/- Sample Purification Beads	149

Methylation Analysis

NEBNext Enzymatic Methyl-seq Kit	151
NEBNext Enzymatic Methyl-seq Conversion Module	151
NEBNext Q5U Master Mix	151
NEBNext Multiplex Oligos for Enzymatic Methyl-seq (Unique Dual Index Primer Pairs)	151

Ultra II RNA

NEBNext Ultra II Directional RNA Library Prep Kits	
+/- Sample Purification Beads	153
NEBNext Ultra II RNA Library Prep Kits	
+/- Sample Purification Beads	153

Reagents for RNA Depletion

NEBNext rRNA Depletion Kits v2 (Human/Mouse/Rat)	
+/- Sample Purification Beads	155
NEBNext rRNA Depletion Kits (Bacteria)	
+/- Sample Purification Beads	155
NEBNext Globin RNA & rRNA Depletion Kits (Human/Mouse/Rat)	
+/- Sample Purification Beads	156
NEBNext Poly(A) mRNA	
Magnetic Isolation Module	156

Customized RNA Depletion

NEBNext RNA Depletion Core Reagent Set	157
NEBNext RNA Depletion Core Reagent Set with RNA Sample Purification Beads	157

Single Cell RNA

NEBNext Single/Cell Low Input RNA Library Prep Kit for Illumina	158
NEBNext Single/Cell Low Input cDNA Synthesis & Amplification Module	158
NEBNext Single Cell Lysis Module	158

Small RNA

NEBNext Multiplex Small RNA Library Prep Sets for Illumina (Sets 1 & 2)	159
NEBNext Multiplex Small RNA Library Prep Kit for Illumina (Index Primers 1-48)	159
NEBNext Small RNA Library Prep Set for Illumina (Multiplex Compatible)	159

Modules & Enzymes

NEBNext Ultra II FS DNA Module	150
Ultra II RNA First Strand Synthesis Module	154
Ultra II Directional RNA Second Strand Synthesis Module	154
Ultra II Non-directional RNA Second Strand Synthesis Module	154
NEBNext Ultra II End Repair/	
dA-Tailing Module	150, 154
NEBNext Ultra II Ligation Module	150, 154
NEBNext Ultra II Q5 Master Mix	150

Adaptors & Primers

NEBNext Multiplex Oligos	
(Unique Dual Index UMI Adaptors)	160
NEBNext Multiplex Oligos	
(96 Unique Dual Index Primer Pairs Sets 1-4)	160
NEBNext Multiplex Oligos for Enzymatic Methyl-seq (Unique Dual Index Primer Pairs)	160
NEBNext Multiplex Oligos	
(Dual Index Primers Sets 1 & 2)	160
NEBNext Multiplex Oligos (Index Primers Sets 1-4)	160
NEBNext Multiplex Oligos (96 Index Primers)	160
NEBNext Multiplex Oligos	
(Methylated Adaptor, Index Primers Set 1)	160
NEBNext Adaptor Dilution Buffer	160

NEBNext ARTIC PRODUCTS FOR SARS-COV-2 SEQUENCING

NEBNext ARTIC SARS-CoV-2 FS Library Prep Kit (Illumina)	161
NEBNext ARTIC SARS-CoV-2 Library Prep Kit (Illumina)	161
NEBNext ARTIC SARS-CoV-2 Companion Kit (Oxford Nanopore Technologies)	161

REAGENTS FOR OXFORD NANOPORE TECHNOLOGIES SEQUENCING

NEBNext Companion Module for Oxford Nanopore Technologies Ligation Sequencing	162
NEBNext ARTIC SARS-CoV-2 Companion Kit (Oxford Nanopore Technologies)	162

Target Enrichment

NEBNext Direct Custom Ready Panels	164
NEBNext Direct Genotyping Solution Sample Prep Kits	165

Immune Sequencing Kits (Human & Mouse)

NEBNext Immune Sequencing Kits (Human)	166
NEBNext Immune Sequencing Kits (Mouse)	166

Library Quantitation

NEBNext Library Quant Kit for Illumina	170
NEBNext Library Dilution Buffer	170

REAGENTS FOR ION TORRENT

NEBNext Fast DNA Library Prep Set for Ion Torrent	171
NEBNext Fast DNA Fragmentation & Library Prep Set for Ion Torrent	171

REAGENTS SUITABLE FOR ANY PLATFORM

DNA Repair

NEBNext FFPE DNA Repair Mix	152
-----------------------------	-----

DNA Separation

NEBNext Microbiome DNA Enrichment Kit	167
---------------------------------------	-----

Modules & Enzymes

NEBNext Magnesium RNA Fragmentation Module	154
NEBNext dsDNA Fragmentase	168
NEBNext dsDNA Fragmentase Reaction Buffer v2	168
NEBNext Ultra II Q5 Master Mix	169

Magnetic Separation

NEBNext Magnetic Separation Rack	159
----------------------------------	-----

One or more of these products are covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. For more information, please email us at gbd@neb.com. The use of these products may require you to obtain additional third party intellectual property rights for certain applications.


Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at <https://www.neb.com/support/terms-of-sale>. NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

ILLUMINA® is a registered trademark of Illumina, Inc.
ION TORRENT™ is a trademark of Life Technologies, Inc.

Illumina® DNA Product Selection Chart

Use the following chart to determine the best NEBNext® products for your Illumina DNA library prep needs. For the most up-to-date product and pricing information, visit [NEBNext.com](https://www.nebnext.com).

For help selecting products, try our online product selection tool at NEBNextSelector.neb.com




* Module and EM-seq Oligos also available.

Reagents for original Ultra workflow and standard workflow are also available. See ordering information.

Illumina RNA Product Selection Chart

Use the following chart to determine the best NEBNext products for your Illumina RNA sequencing needs. For the most up-to-date product and pricing information, visit NEBNext.com.

For help selecting products, try our online product selection tool at NEBNextSelector.neb.com



Reagents for original Ultra workflow and standard workflow are also available. See ordering information.


The heart of the matter – NEBNext® Ultra™ II Workflow

As sequencing technologies continue to improve and applications expand, the need for compatibility with ever-decreasing input amounts and sub-optimal sample quality grows. Reliability and high performance are critical, along with faster turnaround, higher throughput and automation compatibility.


The NEBNext Ultra II workflow lies at the heart of NEB's portfolio for next generation sequencing library preparation. NEBNext Ultra II kits and modules for Illumina are the perfect combination of reagents, optimized formulations and simplified workflows, enabling you to generate DNA or RNA libraries of the highest quality and yield, even when starting from extremely low input amounts.

- Learn one central workflow and apply it to a suite of different applications
- Save time with streamlined modular workflows, reduced hands-on time and automation compatibility
- Benefit from low input amount requirements, fewer PCR cycles and uniform GC coverage in all applications

The Ultra II workflow is central to many NEBNext solutions, including:




- Hands on time: < 15 minutes
- Total time: ~ 2:30 – 3:00 hrs



Ultra II DNA FS input amounts:
100 pg – 0.5 µg unsheared DNA




- Hands on time: < 15 minutes
- Total time: ~ 2:30 hrs



Ultra II RNA input amounts:
5 ng – 1 µg total RNA (rRNA depletion)
10 ng – 1 µg total RNA (poly(A) mRNA)




- Hands on time: < 30 minutes
- Total time: ~ 5:30 – 6:30 hrs



EM-seq input amounts:
10 ng – 200 ng



- Total time: 6 hrs – 7 hrs



Single Cell/Low Input RNA input amounts:
2 pg – 200 ng total RNA



- Hands on time: < 30 minutes
- Total time: 6 hrs – 7 hrs

NEBNext Ultra II DNA, FS and PCR-free DNA Library Prep Kits for Illumina

NEBNext Ultra II DNA Library Prep Kit
for Illumina

#E7645S	24 reactions\$590
#E7645L	96 reactions\$2,245

NEBNext Ultra II DNA Library Prep
with Sample Purification Beads

#E7103S	24 reactions\$645
#E7103L	96 reactions\$2,465

NEBNext Ultra II FS DNA Library Kit
for Illumina

#E7805S	24 reactions\$635
#E7805L	96 reactions\$2,400

NEBNext Ultra II FS DNA Library Prep
with Sample Purification Beads

#E6177S	24 reactions\$695
#E6177L	96 reactions\$2,625

NEW

NEBNext Ultra II DNA PCR-Free Library Prep
Kit for Illumina

#E7410S	24 reactions\$560
#E7410L	96 reactions\$2,132

NEW

NEBNext Ultra II DNA PCR-Free Library Prep
with Sample Purification Beads

#E7415S	24 reactions\$612
#E7415L	96 reactions\$2,341

NEW

NEBNext Ultra II FS DNA PCR-Free Library
Prep Kit for Illumina

#E7430S	24 reactions\$2,280
#E7430L	96 reactions\$603

NEW

NEBNext Ultra II FS DNA PCR-Free Library
Prep with Sample Purification Beads

#E7435S	24 reactions\$660
#E7435L	96 reactions\$2,493

See ordering information for
NEBNext Ultra II DNA modules.

- Get more of what you need, with the highest library yields
- Generate high quality libraries even with limited amounts of DNA, as low as 500 pg
- Prepare libraries from ALL of your samples, including GC-rich targets and FFPE DNA samples
- Save time with streamlined workflows, reduced hands-on time, and automation compatibility
- Access reliable and easy-to-use, scalable enzymatic DNA fragmentation, integrated into the workflow with the FS kit
- Enjoy the flexibility and reliability of the gold standard SPRIselect size selection and clean-up beads, supplied in just the amounts you need


Visit NEBNextUltraII.com for more information, including our technical notes and protocol videos

Description: NEBNext Ultra II DNA Library Prep Kits for Illumina meet the challenge of constructing high quality libraries from ever-decreasing input quantities, enabling high yield preparation of high quality libraries from 500 picograms to 1 microgram of input DNA. Ultra II kits use a fast, streamlined, automatable workflow and enable use of fewer PCR cycles while also improving GC coverage. The kit is also effective with challenging samples such as FFPE DNA.


The Ultra II FS DNA Library Prep Kit combines robust enzymatic DNA fragmentation with end repair and dA-tailing, integrated into a streamlined library prep workflow.

PCR-free kits are now available for both the Ultra II DNA and Ultra II FS DNA workflows.

All Ultra II kits are available with or without SPRIselect® beads.



NEBNext Ultra II provides uniform GC coverage for microbial genomic DNA over a broad range of GC composition and input amounts. Libraries were made using 500 pg, 1 ng and 100 ng of the genomic DNAs shown and the Ultra II DNA Library Prep Kit and sequenced on an Illumina MiSeq®. Reads were mapped using Bowtie 2.2.4 and GC coverage information was calculated using Picard's CollectGCBiasMetrics (v1.117). Expected normalized coverage of 1.0 is indicated by the horizontal grey line, the number of 100 bp regions at each GC% is indicated by the vertical grey bars, and the colored lines represent the normalized coverage for each library.





NEBNext Ultra II provides the highest and most uniform coverage of difficult sequence regions.

A: Indexed libraries were prepared from 100 ng of Human NA19240 genomic DNA using a PCR-free workflow or the library prep kits shown, following manufacturers' recommendations. The PCR-free library was prepared using NEBNext Ultra II. Libraries were sequenced on the Illumina NextSeq® 500. 420 million reads were randomly extracted from each dataset, to produce an average coverage of 10X. Reads were mapped to the GRCh37 reference genome using Bowtie 2.2.4. Reads on each region were counted using bedtools v2.19.1. The number of reads overlapping distinct difficult, low-coverage regions of the human genome (1) are shown for each library. Ultra II provides the highest and most uniform coverage of these difficult regions, and provides the coverage closest to that obtained with a PCR-free protocol.

B: From the 420 million 75 bp reads randomly extracted from each dataset, 10X coverage was expected. The percentage of difficult regions covered at ≥ 10X is shown for each library prep kit and for the PCR-free workflow. Ultra II provides the highest percentage of reads at ≥ 10X coverage and also provides the coverage closest to that obtained with a PCR-free protocol. (1) Aird, D. et al. (2011). Analyzing and minimizing PCR amplification bias in Illumina sequencing libraries. Genome Biology 12(2), R18.

NEBNext Ultra II DNA Reagents for Illumina Sequencing

NEBNext Ultra II Kits for DNA are available with or without integrated enzymatic DNA fragmentation. Note that adaptors and primers are supplied separately. In addition to stringent QC's on individual components, the NEBNext DNA kits are also functionally validated by preparation of a library, followed by Illumina sequencing.

Input							
Ultra II DNA Workflow: 500 pg – 1 µg (PCR); 250 ng – 1,000 ng (PCR-free) Ultra II FS DNA Workflow: 500 pg – 1 µg (PCR); 50 ng – 500 ng (PCR-free)							
Ultra II DNA Library Prep Kits	Fragmentation	End Repair/dA-Tailing	Adaptor Ligation	Clean Up/ Size Selection	PCR Enrichment	Clean Up	Total Workflow
	NEBNext Ultra II DNA Library Prep (NEB #E7645) – with Sample Purification Beads (NEB #E7103)						 Hands-On (not including fragmentation) 12 – 13 min
	<ul style="list-style-type: none">Ultra II End Prep Enzyme MixUltra II End Prep Reaction Buffer (10X)	<ul style="list-style-type: none">Ultra II Ligation Master MixLigation Enhancer	<ul style="list-style-type: none">Sample Purification Beads (SPRIselect) (NEB #E7103 only)	<ul style="list-style-type: none">NEBNext Ultra II Q5 Master Mix	<ul style="list-style-type: none">Sample Purification Beads (SPRIselect) (NEB #E7103 only)		
	NEBNext Ultra II DNA PCR-free Library Prep (NEB #E7410) – with Sample Purification Beads (NEB #E7415)						Total 1.7 – 3.2 hrs
	<ul style="list-style-type: none">Ultra II End Prep Enzyme MixUltra II End Prep Reaction Buffer (10X)	<ul style="list-style-type: none">Ultra II Ligation Master MixLigation Enhancer	<ul style="list-style-type: none">Sample Purification Beads (SPRIselect) (NEB #E7415 only)				
Ultra II DNA Modules	NEBNext Ultra II FS DNA Library Prep (NEB #E7805) – with Sample Purification Beads (NEB #E6177)						 Hands-On (including fragmentation) 12 – 13 min
	<ul style="list-style-type: none">Ultra II FS Enzyme MixUltra II FS Reaction Buffer	<ul style="list-style-type: none">Ultra II Ligation Master MixLigation Enhancer	<ul style="list-style-type: none">Sample Purification Beads (SPRIselect) (NEB #E6177 only)	<ul style="list-style-type: none">Ultra II Q5 Master Mix	<ul style="list-style-type: none">Sample Purification Beads (SPRIselect) (NEB #E6177 only)		
	NEBNext Ultra II FS DNA PCR-free Library Prep (NEB #E7430) – with Sample Purification Beads (NEB #E7435)						Total 1.4 – 3.2 hrs
	<ul style="list-style-type: none">Ultra II FS Enzyme MixUltra II FS Reaction Buffer	<ul style="list-style-type: none">Ultra II Ligation Master MixLigation Enhancer	<ul style="list-style-type: none">Sample Purification Beads (SPRIselect) (NEB #E7435 only)				
	NEBNext Ultra II FS DNA Module (NEB #E7810)						
<ul style="list-style-type: none">Ultra II FS Enzyme MixUltra II FS Reaction Buffer							
NEBNext dsDNA Fragmentase® (NEB #M0348)	NEBNext Ultra II End Repair/dA-Tailing Module (NEB #E7546)	NEBNext Ultra II Ligation Module (NEB #E7595)		NEBNext Ultra II Q5 Master Mix (NEB #M0544)			
<ul style="list-style-type: none">dsDNA FragmentaseReaction Buffer v2Magnesium Chloride	<ul style="list-style-type: none">Ultra II End Prep Enzyme MixUltra II End Prep Reaction Buffer (10X)	<ul style="list-style-type: none">Ultra II Ligation Master MixLigation Enhancer		<ul style="list-style-type: none">Ultra II Q5 Master Mix			

NEBNext Enzymatic Methyl-seq

NEBNext Enzymatic Methyl-seq Kit

#E7120S	24 reactions\$895
#E7120L	96 reactions\$3,360

NEBNext Enzymatic Methyl-seq Conversion Module

#E7125S	24 reactions\$185
#E7125L	96 reactions\$675

NEBNext Q5U™ Master Mix

#M0597S	50 reactions\$125
#M0597L	250 reactions\$495

NEBNext Multiplex Oligos for Enzymatic Methyl-seq (Unique Dual Index Primer Pairs)


#E7140S	24 reactions\$135
#E7140L	96 reactions\$535

- Superior sensitivity of detection of 5mC and 5hmC
- Larger library insert sizes
- More uniform GC coverage
- Greater mapping efficiency
- High-efficiency library preparation


NEBNext Enzymatic Methyl-seq is an enzymatic alternative to bisulfite conversion with superior performance. For more information, including extensive performance data, visit NEBNext.com.

Description: While bisulfite sequencing has been the gold standard for the study of DNA methylation, this conversion treatment is damaging to DNA, resulting in DNA fragmentation, loss and GC bias. The NEBNext Enzymatic Methyl-seq Kit (EM-seq™) provides an enzymatic alternative to whole genome bisulfite sequencing (WGBS), combined with high-efficiency streamlined library preparation suitable for Illumina sequencing.

The highly effective EM-seq enzymatic conversion minimizes damage to DNA and, in combination with the supplied NEBNext Ultra II library preparation workflow reagents, results in high quality libraries that enable superior detection of 5mC and 5hmC from fewer sequencing reads.



NEBNext Enzymatic Methyl-seq libraries have larger insert sizes 50 ng Human NA12878 genomic DNA was sheared to 300 bp using the Covaris® S2 instrument and used as input into EM-seq and WGBS protocols. For WGBS, NEBNext Ultra II DNA was used for library construction, followed by the Zymo Research EZ DNA Methylation-Gold™ kit for bisulfite conversion. Libraries were sequenced on an Illumina MiSeq (2 x 76 bases) and insert sizes were determined using Picard 2.18.14. The normalized frequency of each insert size was plotted, illustrating that library insert sizes are larger for EM-seq than for WGBS, and indicating that EM-seq does not damage DNA as bisulfite treatment does in WGBS.



EM-seq identifies more CpGs than WGBS, at lower sequencing coverage depth with superior uniformity of GC coverage. 10, 50 and 200 ng Human NA12878 genomic DNA was sheared to 300 bp using the Covaris S2 instrument and used as input into EM-seq and WGBS protocols. For WGBS, NEBNext Ultra II DNA was used for library construction, followed by the Zymo Research EZ DNA Methylation-Gold Kit for bisulfite conversion. Libraries were sequenced on an Illumina NovaSeq® 6000 (2 x 100 bases). Reads were aligned to hg38 using bwa-meth 0.2.2.

A: Coverage of CpGs with EM-seq and WGBS libraries was analyzed using 324 million paired end reads, and each top and bottom strand CpGs were counted independently, yielding a maximum of 56 million possible CpG sites. EM-seq identifies more CpGs at lower depth of sequencing.

B: GC coverage was analyzed using Picard 2.17.2 and the distribution of normalized coverage across different GC contents of the genome (0-100%) was plotted. EM-seq libraries have significantly more uniform GC coverage, and lack the AT over-representation and GC under-representation typical of WGBS libraries.



NEBNext FFPE DNA Repair Mix

#M6630S	24 reactions	\$165
#M6630L	96 reactions	\$575

- Construct high-quality NGS libraries from FFPE DNA samples
- Use upstream of library prep for any NGS platform
- No alteration of DNA sequence
- Rely on NEB's NGS validation process for FFPE DNA library prep

Description: Archiving of clinical materials as Formalin-Fixed, Paraffin-Embedded (FFPE) samples is a common practice. However, the methods used for fixation and storage significantly damage and compromise the quality of nucleic acids from these samples. As a result, it can be challenging to obtain useful information, including high quality sequence data, especially when sample amounts are limited. The NEBNext FFPE DNA Repair Mix is a cocktail of enzymes formulated to repair DNA, and specifically optimized and validated for repair of FFPE DNA samples. The FFPE DNA Repair Mix increases library yield and overall library success rates, without introduction of bias.

Functional Validation: Each lot is functionally validated by repair of FFPE DNA followed by library construction and Illumina sequencing.

For the most up-to-date product and pricing information, visit NEBNext.com

FFPE DAMAGE TYPE	REPAIRED BY THE FFPE DNA ENZYME REPAIR MIX?
Deamination of cytosine to uracil	Yes
Nicks and gaps	Yes
Oxidized bases	Yes
Blocked 3' ends	Yes
DNA fragmentation	No
DNA-protein crosslinks	No

Table 1: Types of FFPE DNA damage and ability to be repaired by the NEBNext FFPE DNA Repair Mix.




Figure 1: Effect of FFPE DNA Repair Mix on library yields. An example of Agilent® Bioanalyzer® traces of libraries prepared from stomach tumor FFPE DNA that was treated with the FFPE DNA Repair Mix, or was untreated, before library construction. Yield improvements of 101% to 458% have been observed.

NEBNext Ultra II Library Prep Kits for RNA

NEBNext Ultra II Directional

RNA Library Prep Kit for Illumina

#E7760S	24 reactions\$935
#E7760L	96 reactions\$3,175

NEBNext Ultra II Directional RNA

Library Prep with Sample Purification Beads

#E7765S	24 reactions\$1,035
#E7765L	96 reactions\$3,525

NEBNext Ultra II RNA Library Prep Kit

for Illumina

#E7770S	24 reactions\$890
#E7770L	96 reactions\$3,025


NEBNext Ultra II RNA Library Prep

with Sample Purification Beads

#E7775S	24 reactions\$985
#E7775L	96 reactions\$3,350

See ordering information for
NEBNext Ultra II RNA modules.

Do you need increased sensitivity and specificity from your RNA-seq experiments? Do you have ever-decreasing amounts of input RNA? To address these challenges, our latest generation of RNA library prep kits generate several fold higher yields of high quality libraries and enable use of lower input amounts and fewer PCR cycles. The kits have streamlined, automatable workflows and are available for directional (strand-specific, using the "dUTP method") and non-directional library prep, with the option of SPRISelect beads for size-selection and clean-up steps.




NEBNext Ultra II Directional RNA produces the highest yields, from a range of input amounts.

Poly(A)-containing mRNA was isolated from 10 ng, 100 ng and 1 µg of Universal Human Reference RNA (Agilent #740000) and libraries were made using the NEBNext Ultra II Directional RNA kit, Kapa Stranded mRNA-Seq kit, Kapa mRNA HyperPrep kit and Illumina TruSeq Stranded mRNA Kit. The input RNA amount and number of PCR cycles are indicated. Library yields from an average of three replicates are shown. Error bars indicate standard deviation. Library yields were assessed using the Agilent Bioanalyzer.


- Get more of what you need, with the highest library yields
- Generate high quality libraries with limited amounts of RNA:
 - 10 ng–1 µg Total RNA (polyA mRNA workflow)
 - 5 ng–1 µg (rRNA depletion workflow)
- Minimize bias, with fewer PCR cycles required
- Maximize the flexibility to order reagents for your specific workflow needs
 - Directional (strand-specific, using the "dUTP method") and non-directional workflow options available
 - rRNA Depletion and poly(A) mRNA Isolation reagents are available separately
 - Adaptors and primers for multiplexing, are available separately
- Save time with streamlined workflows, reduced hands-on time, and automation compatibility
- Enjoy the reliability of the gold standard SPRIselect size selection and clean-up beads, supplied in just the amounts you need
- Rely on robust performance, even with low quality RNA

Visit UltraII.RNA.com to learn more
and to view performance data





Uniformity of Coverage across the DAM1 transcript. Poly(A)-containing mRNA was isolated from 10 ng, 100 ng and 1 µg of Universal Human Reference RNA (Agilent #740000) and libraries were made using the NEBNext Ultra II Directional RNA kit, Kapa Stranded mRNA-Seq kit, Kapa mRNA HyperPrep kit and Illumina TruSeq Stranded mRNA Kit. Coverage across transcript ENST00000369541.3 (DAM1) was assessed by mapping reads directly to the transcriptome (Hisat 2.0.3) and assessing coverage using bedtools cov in 100 bins along the transcript length. Libraries prepared using the NEBNext Ultra II Directional RNA Kit provided superior coverage across the transcript at 100 ng and 10 ng input amounts.

View the NEBNext
Ultra II Directional
RNA Workflow.



NEBNext Ultra II RNA Reagents for Illumina Sequencing

NEBNext Ultra II RNA Kits are available for directional (strand-specific) and non-directional library preparation, and for bulk RNA and single cell samples. These kits utilize streamlined workflows and have been designed for performance with input amounts as low as 5 ng. Note that reagents for rRNA depletion and poly(A) mRNA enrichment are supplied separately, as are adaptors and primers. In addition to stringent QC's on individual components, the NEBNext RNA kits are functionally validated by preparation of a library, followed by Illumina sequencing.

Input		Poly(A) mRNA Workflow: 10 ng – 1 µg		rRNA Depletion Workflow: 5 ng – 1 µg																			
mRNA Isolation/ rRNA Depletion		mRNA Fragmentation		First Strand cDNA Synthesis		Second Strand cDNA Synthesis		End Repair/ dA Tailing		Adaptor Ligation		Size Selection		PCR Enrichment		Clean Up		Total Workflow					
Ultra II Directional Kits		NEBNext Ultra II Directional RNA Library Prep Kit for Illumina (NEB #E7760) – with Sample Purification Beads (NEB #7765)																 Hands-On 27 min Total 5.5 – 5.7 hrs* 6.6 – 6.8 hrs**					
		• First Strand Synthesis Reaction Buffer • First Strand Synthesis Enzyme Mix • Random Primers • Strand Specificity Reagent		• Second Strand Synthesis Reaction Buffer with dUTP Mix • Second Strand Synthesis Enzyme Mix • Nuclease-free Water		• End Prep Enzyme Mix • End Repair Reaction Buffer		• Ultra II Ligation Master Mix • Ligation Enhancer • Adaptor Dilution Buffer		• Sample Purification Beads (SPRIselect) – NEB #E7765 only		• Ultra II Q5 Master Mix		• Sample Purification Beads (SPRIselect) – NEB #E7765 only									
Ultra II Non-directional Kits		NEBNext Ultra II RNA Library Prep Kit for Illumina (NEB #E7770) – with Sample Purification Beads (NEB #E7775)																 Hands-On 27 min Total 5.5 – 5.7 hrs* 6.6 – 6.8 hrs**					
		• First Strand Synthesis Reaction Buffer • First Strand Synthesis Enzyme Mix • Random Primers • Strand Specificity Reagent		• Second Strand Synthesis Reaction Buffer • Second Strand Synthesis Enzyme Mix • Nuclease-free Water		• End Prep Enzyme Mix • End Repair Reaction Buffer		• Ultra II Ligation Master Mix • Ligation Enhancer • Adaptor Dilution Buffer		• Sample Purification Beads (SPRIselect) – NEB #E7775 only		• Ultra II Q5 Master Mix		• Sample Purification Beads (SPRIselect) – NEB #E7775 only									
rRNA Depletion Kit v2 (Human/Mouse/Rat) (NEB #E7400, #E7405)		Magnesium RNA Fragmentation Module (NEB #E6150)	Ultra II RNA First Strand Synthesis Module (NEB #E7771)	Ultra II Directional RNA Second Strand Synthesis Module (NEB #E7550)	NEBNext Ultra II End Repair/dA-Tailing Module (NEB #E7546)	NEBNext Ultra II Ligation Module (NEB #E7595)											NEBNext Ultra II Q5 Master Mix (NEB #M0544)	* Including poly(A) mRNA isolation ** Including rRNA depletion					
• RNase H/RNase H Reaction Buffer • rRNA Depletion Solution • Probe Hybridization Buffer • DNase I/DNase I Reaction Buffer • Nuclease-free Water • RNA Sample Purification Beads (NEB #E7405 only)		• RNA Fragmentation Buffer • RNA Fragmentation Stop Solution	• First Strand Synthesis Reaction Buffer • First Strand Synthesis Enzyme Mix • Random Primers • Strand Specificity Reagent	• Second Strand Synthesis Enzyme Mix • Second Strand Synthesis Reaction Buffer with dUTP	• Ultra II End Prep Enzyme Mix • Ultra II End Prep Reaction Buffer	• Ultra II Ligation Master Mix • Ligation Enhancer											• Ultra II Q5 Master Mix						
				Ultra II Non-directional RNA Second Strand Synthesis Module (NEB #E6111)																			
				• Second Strand Synthesis Enzyme Mix • Second Strand Synthesis Reaction Buffer																			
mRNA Isolation/rRNA Depletion (continued)																							
Globin & rRNA Depletion Kit (Human/Mouse/Rat) (NEB #E7750, #E7755)		NEBNext RNA Depletion Core Reagent Set (NEB #E7865, #E7870)		Poly(A) mRNA Magnetic Isolation Module (NEB #E7490)																			
• RNase H/RNase H Reaction Buffer • Globin & rRNA Depletion Solution • Probe Hybridization Buffer • DNase I/DNase I Reaction Buffer • Nuclease-free Water • RNA Sample Purification Beads (NEB #E7775 only)		• Probe Hybridization Buffer • Thermostable RNase H • RNase H & DNase I Reaction Buffers • DNase I • Nuclease-free Water • RNA Sample Purification Beads (NEB #E7870 only)		• Oligo d(T)25 beads • RNA Binding Buffer • Wash Buffer • Nuclease-free Water																			

NEBNext rRNA Depletion Kits (Human/Mouse/Rat and Bacteria)

NEW

NEBNext rRNA Depletion Kit v2
(Human/Mouse/ Rat)

#E7400S	6 reactions\$297
#E7400L	24 reactions\$1,080
#E7400X	96 reactions\$3,885

NEW

NEBNext rRNA Depletion Kit v2
(Human/Mouse/Rat) with RNA Sample
Purification Beads

#E7405S	6 reactions\$305
#E7405L	24 reactions\$1,125
#E7405X	96 reactions\$4,050

NEW

NEBNext rRNA Depletion Kit (Bacteria)

#E7850S	6 reactions\$297
#E7850L	24 reactions\$1,080
#E7850X	96 reactions\$3,885

NEW

NEBNext rRNA Depletion Kit (Bacteria) with
RNA Sample Purification Beads

#E7860S	6 reactions\$305
#E7860L	24 reactions\$1,125
#E7860X	96 reactions\$4,050

- Suitable for low-quality (e.g., FFPE) and high-quality RNA
- Compatible with a broad range of input amounts: 10 ng–1 µg
- Superior depletion of abundant RNAs, with retention of RNAs of interest
- Fast workflow: 2 hours, with less than 10 minutes hands-on time
- Depleted RNA is suitable for RNA-seq, random-primed cDNA synthesis, or other downstream RNA analysis applications
- Available with optional Agencourt® RNAClean® XP Beads


Description: The NEBNext rRNA Depletion Kit (Human/Mouse/Rat) employs an RNase H-based method (1,2) to deplete cytoplasmic (5S, 5.8S, 18S, 28S, human ITS, ETS) and mitochondrial (12S and 16S) rRNA from human total RNA preparations.

Specific enrichment of bacterial mRNAs is challenging due to their lack of poly(A) tails, precluding the use of oligo d(T)-based enrichment methods. For these


samples, specific removal of bacterial rRNAs is an efficient way to enrich for RNAs of interest.

The NEBNext rRNA Depletion Kit (Bacteria) targets removal of rRNA (5S, 16S and 23S) from gram-positive and gram-negative organisms, from monocultures or samples with mixed bacterial species.

- (1) Adiconis, X. et al (2013) *Nature Methods*, 10, 623–629.
- (2) Morlon, J.D. et al (2012) *PLoS One*, 77 e42882.



NEBNext rRNA Depletion Kit Workflow



Depletion of ribosomal RNA enriches for RNAs of interest, and maintains expression correlation, across a mock community of bacterial species and a range of input amounts. Total RNA was extracted from a lyophilized pool of 20 different bacterial organisms (ATCC® #MSA-2002). Ribosomal RNA was depleted using the NEBNext rRNA Depletion Kit (Bacteria). RNA-seq libraries were prepared from untreated and depleted RNA using the NEBNext Ultra II Directional RNA Library Prep Kit for Illumina followed by paired-end sequencing (2 x 75 bp). 4 Million read pairs were sampled (seqtk) from each library, mapped to a composite genome (Bowtie 2.3.2) before counting reads on genes (htseq-count) and correlating their levels. Effective depletion of sequences overlapping with annotated rRNA regions was observed at 100 ng and 10 ng of input RNA for most of the organisms. Correlation analysis of the transcripts indicates consistent transcript expression regardless of treatment or input amount.

NEBNext Globin & rRNA Depletion Kits

NEBNext Globin & rRNA Depletion Kit
(Human/Mouse/Rat)

#E7750S	6 reactions\$340
#E7750L	24 reactions\$1,235
#E7750X	96 reactions\$4,450

NEBNext Globin & rRNA Depletion Kit
(Human/Mouse/Rat) with RNA Sample
Purification Beads

#E7755S	6 reactions\$350
#E7755L	24 reactions\$1,290
#E7755X	96 reactions\$4,635

- Efficient, specific depletion of globin mRNA and rRNA
- Suitable for low and high-quality RNA
- Broad range of input amounts: 10 ng–1 µg
- Superior deletion of abundant RNAs
- Fast workflow: 2 hours, with less than 10 minutes hands-on time

Description: The great majority of RNA in blood samples is comprised of globin mRNA as well as cytoplasmic and mitochondrial ribosomal RNAs (rRNA). These highly abundant RNA species can conceal the biological significance of less abundant transcripts, and so their efficient and specific removal is desirable.


The NEBNext Globin & rRNA Depletion Kit (Human/Mouse/Rat) employs the NEBNext RNase H-based RNA depletion workflow to deplete the following:

- Globin mRNA (HBA1/2, HBB, HBD, HBM, HBG1/2, HBE1, HBQ1 and HBZ)
- Cytoplasmic rRNA (5S, 5.8S, 18S, 28S, ITS and ETS)
- Mitochondrial rRNA (12S and 16S)

The kit is effective with human, mouse and rat total RNA preparations, both intact and degraded. The resulting depleted RNA is suitable for RNA-seq, random-primed cDNA synthesis, or other downstream RNA analysis.

This kit can also be used following poly(A) mRNA enrichment (e.g., using the NEBNext poly(A) mRNA Magnetic Isolation Module, NEB #E7490), so that the final depleted RNA contains only mRNA of interest and no non-coding RNA.

This kit is available with or without RNAClean beads.



NEBNext Globin & rRNA Depletion Kit efficiently removes Globin mRNA and rRNA. Ribosomal RNA (rRNA) and globin mRNA were depleted from Human, Mouse, and Rat Whole Blood Total RNA (100 ng) using the NEBNext Globin & rRNA Depletion Kit (Human/Mouse/Rat) (NEB #E7750) or the Globin-Zero® Gold rRNA Removal Kit (Illumina #GZG1224). Libraries were prepared from the depleted RNA using the NEBNext Ultra II Directional RNA Library Prep Kit for Illumina and sequenced on an Illumina NextSeq® instrument (2 x 75 bp). Reads were down sampled to 20 million reads per sample for analysis, and were identified as globin mRNA (A) or rRNA (B) using Mirabait (6 or more, 25-mers). The data represents an average of 3-4 replicates. Error bars indicate standard error.

NEBNext Poly(A) mRNA Magnetic Isolation Module

#E7490S	24 reactions\$69
#E7490L	96 reactions\$247

Description: The NEBNext Poly(A) mRNA Magnetic Isolation Module is designed to isolate intact poly(A)+ RNA from previously-isolated total RNA. The technology is based on the coupling of Oligo d(T)₂₅ to 1 µm paramagnetic beads which are then used as the solid support for the direct binding of poly(A)+ RNA. Thus, the procedure permits the manual processing of multiple samples and can be adapted for automated high-throughput applications. Additionally, magnetic separation technology permits elution of intact mRNA in small volumes eliminating the need for precipitating the

poly(A)+ transcripts in the eluent. Intact poly(A)+ RNA which is fully representative of the mRNA population of the original sample can be obtained in less than one hour.

The mRNA Magnetic Isolation Module Includes:

- NEBNext Oligo d(T)₂₅ Beads
- NEBNext RNA Binding Buffer (2X)
- NEBNext Wash Buffer
- Nuclease-free Water
- NEBNext Tris Buffer

Customized Depletion of Unwanted RNA

NEW

NEBNext RNA Depletion Core Reagent Set		
#E7865S	6 reactions	\$254
#E7865L	24 reactions	\$924
#E7865X	96 reactions	\$3,326

NEW

NEBNext RNA Depletion Core Reagent Set with RNA Sample Purification Beads		
#E7870S	6 reactions	\$262
#E7870L	24 reactions	\$963
#E7870X	96 reactions	\$3,465

- Compatible with a broad range of input amounts: 10 ng - 1 µg
- Suitable for low-quality or high-quality RNA
- Fast workflow: 2 hours, with less than 10 minutes hands-on time



Design oligos for depletion of unwanted RNA from any organism, when used in the NEBNext RNA depletion workflow.
<https://depletion-design.neb.com/>


Description: In RNA-seq, highly expressed transcripts with minimal biological interest, such as ribosomal RNA (rRNA) can dominate readouts and mask detection of more informative low-abundance transcripts. This challenge is amplified when working with sample types for which pre-designed RNA depletion kits are not available. The NEBNext RNA Depletion Core Reagent set provides a customized approach to deplete unwanted RNA from any organism, using probe sequences designed with the user-friendly NEBNext Custom RNA Depletion Design Tool.

The efficient RNase-H-based workflow, and close tiling of probes designed using the online tool, enables effective depletion from both low- and high-quality RNA, with a broad range of input amounts.


STEP 1: Use the online NEBNext Custom RNA Depletion Design Tool to obtain custom probe sequences, by entering the sequence of your target RNA.

STEP 2: Order ssDNA probe oligonucleotides from your trusted oligo provider.

STEP 3: Use the probes with the NEBNext Custom RNA Depletion Core Reagent Set or in combination with other NEBNext RNA Depletion Kits



NEBNext Custom RNA Depletion enriches for RNAs of interest by efficiently removing targeted RNA from total RNA across species and a wide range of inputs. The NEBNext Custom RNA Depletion Design Tool was used to design probes against rRNA of the mosquito *Aedes aegypti*, and the archaeal species *Thermococcus kodakarensis* and *Pyrococcus furiosus*. Total RNA (1 µg, 100 ng or 10 ng) was used as input for rRNA depletion using the Core RNA Depletion Reagent Set with the designed probes. RNA-seq libraries were prepared using the NEBNext Ultra II Directional RNA Library Prep Kit for Illumina followed by paired-end sequencing (2 x 75 bp). 20 million reads were sampled (seqtk) from each library. Read pairs were identified as ribosomal using mirabait (6 or more shared 25-mers), and levels of rRNA remaining were calculated by dividing matched reads by the total number of reads passing instrument quality filtering. The data represents an average of 3 replicates. The method efficiently depletes targeted rRNA across species and a wide range of total RNA input amount (1 µg–10 ng).



Combined probe pools efficiently deplete human rRNA and mitochondrial mRNA using NEBNext Custom RNA Depletion. The NEBNext Custom RNA Depletion Design Tool was used to design probes against human mitochondrial mRNA. The probes were used in combination with the probe pool from the NEBNext rRNA Depletion Kit v2 (Human/Mouse/Rat). 1 µg of total Universal Human Reference RNA (Agilent®) was depleted of mitochondrial RNA and rRNA using the Core RNA Depletion Reagent Set. RNA-seq libraries were prepared using the NEBNext Ultra II Directional RNA Library Prep Kit for Illumina followed by paired-end sequencing (2 x 75 bp). 20 million reads were sampled (seqtk) from each library.

A. Read pairs were identified as ribosomal and mitochondrial using mirabait (6 or more, 25-mers), and levels of rRNA and mtRNA remaining were calculated by dividing matched reads by the total number of reads passing instrument quality filtering. Both rRNA and mitochondrial RNA are efficiently depleted.

B. Integrative Genome Viewer (IGV) visualization of read coverage across the human mitochondrial genes.

NEBNext Single Cell/Low Input RNA Library Prep

NEBNext Single Cell/Low Input RNA Library
Prep Kit for Illumina

#E6420S 24 reactions \$1,195
#E6420L 96 reactions \$3,965

NEBNext Single Cell/Low Input cDNA
Synthesis & Amplification Module

#E6421S 24 reactions \$650
#E6421L 96 reactions \$2,210

NEBNext Single Cell Lysis Module


#E5530S 96 reactions \$92

Description: The unique workflow of the NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina meets the demand for a highly sensitive, yet robust method that consistently generates high-quality, full-length transcript sequencing data from a single cell or ultra-low input RNA.


Optimized cDNA synthesis and amplification steps incorporate template switching, and utilize a unique protocol and suite of reagents.

cDNAs are generated directly from single cells or 2 pg–200 ng RNA, and even low-abundance transcripts are represented in the high yields of cDNA obtained. This is followed by library construction that incorporates the Ultra II FS enzymatic DNA fragmentation/end repair/dA-tailing mix in a simple and efficient workflow.

- Generate the highest yields of high-quality full-length transcript sequencing libraries from single cells, or as little as 2 pg–200 ng total RNA
- Experience unmatched detection of low abundance transcripts
- Rely on consistent transcript detection for a wide range of input amounts and sample types
- Obtain full-length, uniform transcript coverage, regardless of input amount or sample type
- Use with cultured or primary cells, or total RNA
- Save time with a fast, streamlined workflow, minimal handling steps and hands-on time
 - Single-tube protocol from cell lysis to cDNA
 - Enzymatic DNA fragmentation, end repair and dA-tailing reagents in a single enzyme mix, with a single protocol, regardless of GC content
- Available with or without library construction reagents



Increased transcript detection with the NEBNext Single Cell/Low Input RNA Library Prep Kit Sequencing libraries were generated from Jurkat single cells (6 replicates) using the NEBNext Single Cell/Low Input RNA Library Prep Kit, or the SMART-Seq v4 Ultra Low Input RNA Kit for Sequencing (Clontech #634891) plus the Nextera® XT DNA Library Prep Kit (Illumina #FC-131-1096). Libraries were sequenced on an Illumina NextSeq 500 using paired-end mode (2 x 76 bp). TPM = Transcripts per Kilobase Million. Each dot represents the number of transcripts identified at the given TPM range, and each box represents the median, first and third quartiles per replicate and method. Salmon 0.6 was used for read mapping and quantification of all GENCODE v25 transcripts. Panels show the number of transcripts detected within the following TPM ranges: 1–5, 5–10, 10–50 and > 50 TPM. Increased identification of low abundance transcripts is observed with the NEBNext libraries.



Higher library yields with the NEBNext Single Cell/Low Input RNA Library Prep Kit Sequencing libraries were generated from HeLa, Jurkat and M1 single cells or 10 pg of Universal Human Reference (UHR) RNA (Agilent #740000) with recommended amounts of ERCC RNA Spike-In Mix I (Thermo Fisher Scientific® #4456740). The NEBNext Single Cell/Low Input RNA Library Prep Kit, or the SMART-Seq v4 Ultra Low Input RNA Kit for Sequencing (Clontech #634891) plus the Nextera XT DNA Library Prep Kit (Illumina #FC-131-1096) were used. Error bars indicate standard deviation for 6–11 replicates. For the NEBNext workflow ~80% of the cDNA was used as input into sequencing library preparation, and libraries were amplified with 8 PCR cycles. For the SMART-Seq v4/Nextera XT workflow, as recommended, 125 pg of cDNA was used as input in sequencing library preparation and 12 PCR cycles were used for amplification. Error bars indicate standard deviation for 6–11 replicates.

NEBNext Small RNA Library Prep Kits

NEBNext Multiplex Small RNA

Library Prep Set for Illumina (Set 1)

#E7300S	24 reactions\$1,470
#E7300L	96 reactions\$4,998

NEBNext Multiplex Small RNA

Library Prep Set for Illumina (Set 2)

#E7580S	24 reactions\$1,470
#E7580L	96 reactions\$4,998

NEBNext Multiplex Small RNA

Library Prep Kit for Illumina (Index Primers 1-48)

#E7560S	96 reactions\$4,998
---------	--------------	-------------




NEBNext Small RNA Library Prep Set for Illumina (Multiplex Compatible)

#E7330S	24 reactions\$1,200
#E7330L	96 reactions\$4,080

For adenylation of custom ssDNA adaptors, the 5' DNA Adenylation Kit is available (NEB #E2610).

- Minimized adaptor-dimer formation
- High yields
- Input RNA can be total RNA
- Suitable for methylated small RNAs (e.g., piRNAs) as well as unmethylated small RNAs

The novel NEBNext Small RNA workflow has been optimized to minimize adaptor-dimers while producing high-yield, high-diversity libraries. Adaptors and primers are included in the Small RNA kits, and multiplexing options are available. The Multiplex kit contains index primers, and the Multiplex-Compatible kit enables use with your own barcode primers.

Input 100 ng – 1 µg						
3' Adaptor Ligation	Primer Hybridization	5' Adaptor Ligation	First Strand Synthesis	PCR Enrichment	Size Selection	Total Workflow
NEBNext Multiplex Small RNA Library Prep Set for Illumina (Set 1 NEB #E7300, Set 2 NEB #E7580)						
<ul style="list-style-type: none"> 3' Ligation Enzyme Mix 3' Ligation Reaction Buffer (2X) 3' SR Adaptor 	<ul style="list-style-type: none"> SR RT Primer 	<ul style="list-style-type: none"> 5' Ligation Enzyme Mix 5' Ligation Reaction Buffer (10X) 5' SR Adaptor Nuclease-Free Water 	<ul style="list-style-type: none"> RNase Inhibitor, Murine M-MuLV Reverse Transcriptase (RNase H⁻) First Strand Synthesis Reaction Buffer 	<ul style="list-style-type: none"> LongAmp[®] Taq 2X Master Mix SR Primer Index Primers 1–12 (Set 1) Index Primers 13–24 (Set 2) 	<ul style="list-style-type: none"> Gel Loading Dye, Blue (6X) Quick-Load[®] pBR322 DNA-MspI Digest DNA Gel Elution Buffer (1X) Linear Acrylamide TE Buffer 	 Hands-On Time 30 min Total Time 6 hrs
NEBNext Multiplex Small RNA Library Prep Kit for Illumina (Index Primers 1-48) (NEB #E7560)						
<ul style="list-style-type: none"> 3' Ligation Enzyme Mix 3' Ligation Reaction Buffer (2X) 3' SR Adaptor 	<ul style="list-style-type: none"> SR RT Primer 	<ul style="list-style-type: none"> 5' Ligation Enzyme Mix 5' Ligation Reaction Buffer (10X) 5' SR Adaptor Nuclease-Free Water 	<ul style="list-style-type: none"> RNase Inhibitor, Murine M-MuLV Reverse Transcriptase (RNase H⁻) First Strand Synthesis Reaction Buffer 	<ul style="list-style-type: none"> LongAmp[®] Taq 2X Master Mix SR Primer NEBNext Index 1-48 Primers for Illumina 	<ul style="list-style-type: none"> Gel Loading Dye, Blue (6X) Quick-Load[®] pBR322 DNA-MspI Digest DNA Gel Elution Buffer (1X) Linear Acrylamide TE Buffer 	 Hands-On Time 30 min Total Time 6 hrs
NEBNext Small RNA Library Prep Set for Illumina (Multiplex Compatible) (NEB #E7330)						
<ul style="list-style-type: none"> 3' Ligation Enzyme Mix 3' Ligation Reaction Buffer (2X) 3' SR Adaptor 	<ul style="list-style-type: none"> SR RT Primer 	<ul style="list-style-type: none"> 5' Ligation Enzyme Mix 5' Ligation Reaction Buffer (10X) 5' SR Adaptor Nuclease-Free Water 	<ul style="list-style-type: none"> RNase Inhibitor, Murine M-MuLV Reverse Transcriptase (RNase H⁻) First Strand Synthesis Reaction Buffer 	<ul style="list-style-type: none"> LongAmp[®] Taq 2X Master Mix SR Primer Index Primer 1 	<ul style="list-style-type: none"> Gel Loading Dye, Blue (6X) Quick-Load[®] pBR322 DNA-MspI Digest DNA Gel Elution Buffer (1X) Linear Acrylamide TE Buffer 	 Hands-On Time 30 min Total Time 6 hrs

NEBNext Magnetic Separation Rack

#S1515S	24 tubes (0.2 ml)\$475
---------	-------------------	-----------

- Fast separations in purification and size-selection steps in next generation sequencing workflows
- Small-scale separation of magnetic particles
- Anodized aluminum rack with Neodymium Iron Boron (NdFeB) rare earth magnets, the most powerful commercially available
- 24 tube capacity:
8- and 12-strip 0.2 ml PCR tubes
or individual 0.2 ml PCR tubes

Description: Next generation sequencing library preparation workflows include magnetic bead-based purification and size-selection steps and it is important for library yield and quality that bead separation be highly efficient and fast.

The NEBNext Magnetic Separation Rack was designed for this application and contains rare earth Neodymium Iron Boron (NdFeB) magnets, the most powerful commercially available magnets, in an anodized aluminium rack. The rack holds 24 (0.2 ml) tubes, and is compatible with single tubes or strip tubes.



NEBNext Adaptors & Primers for Illumina

NEW

NEBNext Multiplex Oligos for Illumina
(Unique Dual Index UMI Adaptors DNA Set 1)
#E7395S 96 rxns (96 indices)\$672
#E7395L 384 rxns (96 indices)\$2,553

NEW

NEBNext Multiplex Oligos for Illumina
(Unique Dual Index Adaptors RNA Set 1)
#E7416S 96 rxns (96 indices)\$648
#E7416L 384 rxns (96 indices)\$2,460

NEBNext Multiplex Oligos for Illumina
(96 Unique Dual Index Primer Pairs)
#E6440S 96 rxns (96 indices)\$535
#E6440L 384 rxns (96 indices)\$1,925

NEW

NEBNext Multiplex Oligos for Illumina
(96 Unique Dual Index Primer Pairs Set 2)
#E6442S 96 rxns (96 indices)\$535
#E6442L 384 rxns (96 indices)\$1,925

NEW

NEBNext Multiplex Oligos for Illumina
(96 Unique Dual Index Primer Pairs Set 3)
#E6444S 96 rxns (96 indices)\$535
#E6444L 384 rxns (96 indices)\$1,925

NEW

NEBNext Multiplex Oligos for Illumina
(96 Unique Dual Index Primer Pairs Set 4)
#E6446S 96 rxns (96 indices)\$535
#E6446L 384 rxns (96 indices)\$1,925

NEBNext Multiplex Oligos for Enzymatic
Methyl-seq (Unique Dual Index Primer Pairs)
#E7140S 24 rxns\$135
#E7140L 96 rxns\$535

NEBNext Multiplex Oligos for Illumina
(Dual Index Primers Set 1)
#E7600S 96 rxns (8 x 12 indices)\$460

NEBNext Multiplex Oligos for Illumina
(Dual Index Primers Set 2)
#E7780S 96 rxns (8 x 12 indices)\$460

NEBNext Multiplex Oligos for Illumina
(Index Primers Set 1)
#E7335S 24 rxns (12 indices)\$105
#E7335L 96 rxns (12 indices)\$385

NEBNext Multiplex Oligos for Illumina
(Index Primers Set 2)
#E7500S 24 rxns (12 indices)\$105
#E7500L 96 rxns (12 indices)\$385

NEBNext Multiplex Oligos for Illumina
(Index Primers Set 3)
#E7710S 24 rxns (12 indices)\$105
#E7710L 96 rxns (12 indices)\$385

NEBNext Multiplex Oligos for Illumina
(Index Primers Set 4)
#E7730S 24 rxns (12 indices)\$105
#E7730L 96 rxns (12 indices)\$385

NEBNext Multiplex Oligos for Illumina
(96 Index Primers)
#E6609S 96 rxns (96 indices)\$650
#E6609L 384 rxns (96 indices)\$2,600

NEBNext Multiplex Oligos for Illumina
(Methylated Adaptor, Index Primers Set 1)
#E7535S 24 rxns (12 indices)\$115
#E7535L 96 rxns (12 indices)\$420

NEBNext Adaptor Dilution Buffer
#B1430S 1 x 9.6 ml\$34





Description: Designed for use in library prep for DNA, ChIP DNA and RNA (but not Small RNA), the NEBNext Multiplex Oligos for Illumina are an essential component of the NGS sample prep workflow. Offering a range of indexing strategies, from Unique Dual Index UMI Adaptors to the truncated, hairpin-loop NEBNext Adaptor, meant for use with UDI, dual, and single index primers, the NEBNext Multiplex Oligos can support NGS across a wide range of formats. Optimized for performance in recommended applications, there's an NEBNext indexing option tailored to you. NEBNext Oligos can be used with NEBNext products, and with other standard Illumina-compatible library preparation protocols.

Single or dual barcode primer options are available. Unique dual index primer pairs are available to address the "index hopping" seen with certain Illumina sequencing instruments.

Unique Dual Index UMI Adaptors (available for both DNA and RNA library prep) offer a ready-to-ligate adaptor for correction of PCR duplicates and errors, while improving the detection of single-nucleotide variants (SNVs). When read without the UMI sequence, the full-length adaptor enables PCR-free DNA library prep.

- Indexing strategies optimized by application
 - Index Primers: For NGS Library Prep workflows that include an amplification step
 - Index Adaptors: Enablement of PCR-free workflows and incorporation of UMIs for error correction/deduplication
- Extensively QC'd for purity and increased library yields
- Flexibility for use with NEBNext library preparation kits and other standard, Illumina-compatible library preparation methods
- Provided with index-pooling guidelines and sample sheets

Multiplex Oligos Selection Chart

	 SINGLE INDEX	 DUAL INDEX	 UNIQUE DUAL INDEX	 UNIQUE DUAL INDEX UMIs
NEB PRODUCTS	NEB #E7335 NEB #E7500 NEB #E7710 NEB #E7730 NEB #E6609	NEB #E7600 NEB #E7780	NEB #E6440 NEB #E6442 NEB #E6444 NEB #E6446 NEB #E7140	NEB #E7395 NEB #E7416
Contains UMI	No	No	No	Yes
Addresses Index Hopping	No	No	Yes	Yes
Indexing Strategy	Index Primer	Index Primer	Index Primer	Index Adaptor
Applications	DNA-seq, RNA-seq (except small RNA)	DNA-seq, RNA-seq (except small RNA)	DNA-seq, RNA-seq (except small RNA)	PCR-free DNA-seq, RNA-seq (except small RNA)
Number of Indices for Multiplexing	up to 144	up to 384	up to 384	up to 96
Compatible with EM-seq*	Yes*	Yes*	Yes*	No
Compatible with EpiMark® Bisulfite Sequencing	Yes**	Yes**	Yes**	No
Number of Sets Available; Formats and Indices Available	Five; Sets 1-4 (12 indices/set): Individual vials 96 Index: premixed plate	Two; Individual vials containing 8 i5 primers and 12 i7 primers for combinatorial mixing	Four; 96 indices in premixed, foil-sealed 96-well plates, including a version for EM-seq (up to 120 indices, either 96-well plate or 24 vial format)	One; 96 indices in premixed, foil-sealed 96-well plate (DNA-seq OR RNA-seq) and primers

* Requires the use of the EM-seq Adaptor; Single, dual and unique dual index are all compatible; NEB recommends using the Unique Dual Index Primers found in the NEBNext Enzymatic Methyl-seq Kit (NEB #E7120) or the NEBNext Multiplex Oligos for EM-seq (NEB #E7140), both supplied with the NEBNext EM-seq Adaptor; For higher levels of multiplexing, Unique Dual Index Primers Sets 3 and 4 (NEB #E6444 and #E6446) are also validated for EM-seq.

** Requires use of NEBNext EM-seq adaptor from NEBNext Multiplex Oligos for Enzymatic Methyl-seq (Unique Dual Index Primer Pairs, #E7140S/L or NEBNext methylated adaptor from NEBNext Multiplex Oligos for Illumina® (Methylated Adaptor, Index Primers Set 1, #E7535S/L).

NEBNext ARTIC Products for SARS-CoV-2 Sequencing

NEW

NEBNext ARTIC SARS-CoV-2 FS Library Prep Kit (Illumina)		
#E7658S	24 reactions\$408
#E7658L	96 reactions\$1,632

NEW

NEBNext ARTIC SARS-CoV-2 Library Prep Kit (Illumina)		
#E7650S	24 reactions\$384
#E7650L	96 reactions\$1,536

NEW

NEBNext ARTIC SARS-CoV-2 Companion Kit (Oxford Nanopore Technologies)		
#E7660S	24 reactions\$192
#E7660L	96 reactions\$768


- Streamlined, high-efficiency protocol
- Ample amplicon yields from a wide range of viral genome inputs
- Improved SARS-CoV-2 genome coverage depth with a more balanced primer pool
- Available for Illumina and Oxford Nanopore Technologies sequencing platforms
- No requirement for amplicon normalization prior to Illumina library preparation

The NEBNext ARTIC kits were developed in response to the critical need for reliable and accurate methods for sequencing viral pathogens, specifically SARS-CoV-2. These kits, for long and short read sequencing, were based on the original work of the ARTIC Network (1). The ARTIC SARS-CoV-2 sequencing workflow is a multiplexed amplicon-based whole-viral-genome sequencing approach.

NEBNext ARTIC Companion kits include primers and reagents for RT-PCR from SARS-CoV-2 gRNA and downstream library preparation for Illumina and Oxford Nanopore Technologies sequencing.

The optimized primers and reagents for RT-PCR deliver uniform, ample amplicon yields from gRNA across a wide copy number range, and library prep and sequencing can be performed downstream of a single RT-PCR procedure.

For Illumina applications, a novel DNA polymerase formulation for the enrichment of next-generation sequencing libraries eliminates the need to normalize amplicon concentrations prior to library preparation. Two library prep options are available: The NEBNext ARTIC SARS-CoV-2 FS Library Prep Kit (Illumina) incorporates enzymatic cDNA fragmentation, and generates libraries with inserts in the 150 bp range. The NEBNext ARTIC SARS-CoV-2 Library Prep Kit (Illumina) does not include DNA fragmentation and library inserts are in the 400 bp range.



(1) Josh Quick 2020. nCoV-2019 sequencing protocol v2 (GunIt). protocols.io <https://dx.doi.org/10.17504/protocols.io.bdp7i5rn>

Reagents for Oxford Nanopore Technologies® Sequencing

NEW

NEBNext Companion Module for Oxford Nanopore Technologies Ligation Sequencing
#E7180S 24 reactions\$920

NEW

NEBNext ARTIC SARS-CoV-2 Companion Kit (Oxford Nanopore Technologies)
#E7660S 24 reactions\$192
#E7660L 96 reactions\$768

For adenylation of custom ssDNA adaptors, the 5' DNA Adenylation Kit is available (NEB #E2610).

Companion Products:

Monarch HMW DNA Extraction Kit for Cells & Blood
#T3050S 5 preps\$70
#T3050L 50 preps\$400

Monarch HMW DNA Extraction Kit for Tissue
#T3060S 5 preps\$75
#T3060L 50 preps\$450

Monarch Genomic DNA Purification Kit
#T3010S 50 preps\$156
#T3010L 150 preps\$402

- Component volumes tailored for use with many SQK-LSK109 and SQK-LSK110 workflows
- Simplified ordering and inventory management
- Compatible with all devices: MinION®, GridION®, PromethION®, Flongle®
- No unnecessary buffers or excess reagents

Many NEBNext and NEB products are recommended for use in multiple sample prep workflows for Oxford Nanopore Technologies sequencing, for a range of sample types and applications.

The NEBNext Companion Module for Oxford Nanopore Technologies Ligation Sequencing includes the NEBNext DNA repair, end repair and ligation reagents recommended in Oxford Nanopore Ligation library preparation. These are provided at volumes designed for use in several protocols alongside Oxford Nanopore Technologies SQK-LSK109 and SQK-LSK110.

The NEBNext Companion Module for Oxford Nanopore Technologies Ligation Sequencing Includes:


- NEBNext FFPE DNA Repair Mix (0.048 ml)
- NEBNext FFPE DNA Repair Buffer (0.084 ml)
- NEBNext Ultra II End Prep Enzyme Mix (0.072 ml)
- NEBNext Ultra II End Prep Reaction Buffer (0.084 ml)
- Quick T4 DNA Ligase (0.240 ml)

Also Available:

Monarch DNA Extraction for Oxford Nanopore Sequencing

Long read sequencing technologies, including Oxford Nanopore sequencing, require high quality extracted DNA. For the longest reads, the Monarch® HMW DNA Extraction kits enable isolation of DNA in the Mb range. The HMW DNA Extraction Kit for Tissue (NEB #T3060) is effective with a variety of tissues, bacteria and other samples (yeast, insect, amphibian), and the HMW DNA Extraction Kit for Cells & Blood (NEB #T3050) isolates HMW DNA from cultured cells and whole blood. When reads < 80 kb are required, the Monarch Genomic DNA Purification Kit (NEB #T3010) produces genomic DNA with a typical peak size of > 50 kb.

For more information see page 136.



Use of varying agitation speeds during lysis produces tunable fragment length of extracted HMW genomic DNA from cells and blood.

DNA extracted with Monarch HMW DNA Extraction Kit for Cells & Blood. 1 x 10⁶ fresh HEK293 cells and 500 µl fresh human blood were used as inputs and for preps performed according to the kit instructions using the agitation speed indicated above the gel lanes. 500 ng of DNA from the replicates was resolved by PFGE (1% agarose gel, 6 V/cm, 13°C for 20 hours, switch times ramped from 0.5–94 seconds on a BioRad® CHEF-DR® III System). Yield and purity ratios of the individual preps are shown in the accompanying tables. Lambda PFG Ladder (NEB #N0341) was used as molecular weight standard.

Target Enrichment for NGS


Target enrichment is used to describe a variety of strategies to selectively isolate specific genomic regions of interest for sequencing analysis.

NEBNext Direct® offers a novel approach to NGS Target Enrichment. Using a unique, capture-based protocol, NEBNext Direct offers a range of solutions for highly specific, efficient targeting of genomic loci of interest. Optimized, 1-day protocols allow for rapid profiling of targeted regions for a variety of applications, including germline and somatic variant calling for translational research and genotyping of genomic markers for applications in agricultural biology and sample identification.

- Define your own targeted sequencing panel by selecting only the genes you want from our expansive list of genes using NEBNext Direct Custom Ready Panels (workflow shown below).
- For genotyping applications, the NEBNext Genotyping Solution enables pre-capture pooling of up to 96 samples for cost-effective and scalable genotyping of up to 5,000 genomic markers.

The NEBNext Microbiome DNA Enrichment Kit facilitates enrichment of microbial DNA from samples containing methylated host DNA (including human), by selective binding and removal of the CpG-methylated host DNA. Importantly, microbial diversity remains intact after enrichment

With the NEBNext Immune Sequencing Kits, available for human and mouse, sequence the full-length immune gene repertoires of B cells and T cells. Profile somatic mutations across all relevant contexts with improved sequence accuracy.



NEBNext Direct Custom Ready Panels employ a fast hybridization-based workflow that couples capture with library preparation



View the
NEBNext Direct
Workflow.

NEBNext Direct Custom Ready Panels

#E6631S	8 reactions
#E6631L	24 reactions
#E6631X	96 reactions

Visit www.neb.com/E6631 to learn more and request a quote.

- Choose from a single gene to hundreds of genes.
- Experience unmatched specificity and coverage uniformity
- Eliminate synthesis and optimization steps, for a faster turnaround
- Improve sensitivity with our Unique Molecular Identifier (UMI)
- Generate results in one day with our automation-friendly workflow


Panels can be designed and ordered by visiting www.neb.com/CustomReadyPanelForm

Description: Employing the unique NEBNext Direct hybridization-based enrichment method, NEBNext Direct Custom Ready Panels allow rapid customization of targeted gene panels for Illumina sequencing. Select from a list of human genes where baits have been carefully designed and optimized to provide complete coverage of the full coding (exon) regions. High quality panels can be designed by you and rapidly delivered, from any combination of genes.

The genes available through the NEBNext Direct Custom Ready offering are continually updated, and currently include those associated with a variety of translational research areas, including cancer, neurological disorders, cardiometabolic disease, autism, severe combined immunodeficiency, cystic fibrosis and the recommended genes for incidental findings by the American College of Medical Genetics. The full list of genes currently available can be found at www.neb.com/CustomReadyPanelForm.


NEBNext Direct Custom Ready Panels can include anywhere from a single specific gene up to 1.5 megabases of total target territory. There are no limitations on genes that can be combined together in a panel. Each panel is tested prior to shipment, and sequencing results are returned through a custom Performance Report.

Bait sets for each gene included in the panel have undergone a rigorous development and optimization process to maximize specificity and target coverage uniformity.




NEBNext Direct Custom Ready Panels demonstrate optimum performance across a wide range of panel sizes.

Key target enrichment metrics demonstrate consistent performance across a range of panel sizes. 100 ng of DNA was tested against panels of 1, 10, 25, 50 and 100 genes, and sequenced using Illumina paired-end 150 bp sequencing. Larger panels included all genes present in smaller panels.



Sensitivity in detection of variants across panel size and DNA input amount. 24 HapMap samples were blended to create a range of variant allele frequencies (VAF) down to 2%. 25, 50, 100, 200 and 500 ng of this blended DNA was enriched using NEBNext Direct Custom Ready Panels of 1, 10, 25, 50 and 100 genes. Larger panels were inclusive of the genes in smaller panels. Resulting libraries were sequenced using 2 x 150 bp Illumina sequencing and variants were called using Mutect and Vardict variant calling algorithms.



NEBNext Direct Custom Ready Panels demonstrate retention of target behavior across panel sizes.

IGV image of coverage profile for 4 BRAF exons included in panels of 1, 10, 25, 50 and 100 genes, demonstrates consistent target behavior with the addition of gene targets. 100 ng of DNA was used as input for NEBNext enrichment using the 5 panels, including the BRAF gene. Libraries were sequenced on an Illumina 2 x 150 basepair sequencing.

NEW

NEBNext Direct Genotyping Solution

NEBNext Direct Genotyping Solution


#E9530B-S 8 reactions

#E9500B-S 96 reactions


- Ideal solution for genotyping hundreds to thousands of markers
- Reduce costs and streamline workflow through pre-capture pooling of up to 96 samples
- Maximize sequencer efficiency through dual barcode sample indexing plus Unique Molecular Identifier
- Unparalleled target coverage uniformity through unique capture-based enrichment
- Eliminate marker dropouts with finely tuned bait design
- Increase sample throughput using 1-day, automatable workflow

Visit www.neb.com/E9500 to learn more.


Description: The NEBNext Direct Genotyping Solution combines highly multiplexed, capture-based enrichment with maximum efficiency next-generation sequencing to deliver cost-effective, high-throughput, genotyping for a wide variety of applications. Applicable for marker ranges spanning 100-5,000 markers, pre-capture multiplexing of up to 96 samples combined with dual indexed sequencing allowing extremely high levels of sample multiplexing in a single Illumina sequencing run.



NEBNext Direct Genotyping Solution workflow.



NEBNext Direct Genotyping Solution demonstrates similar coverage across 96 pooled samples. Mean SNP coverage of 2,309 SolCAP markers across 96 samples. 25 ng of purified tomato DNA was used for each sample. Samples were index-tagged and pooled prior to hybridization and Libraries were sequenced on an Illumina MiSeq with 20 cycles of Read 1 to sequence the 12 base UMI and 8 base sample index, and 75 cycles of Read 2 to sequence the targets.



Mean Coverage across 2309 markers within a single sample. Histogram of coverage across each of the 2,309 SolCAP markers demonstrates evenness of enrichment across targets and coverage levels sufficient for genotyping calls. These data represent enrichment of a single tomato sample pooled with 95 others prior to hybridization. 25 ng of purified tomato DNA was used for each sample. Samples were index-tagged and pooled prior to hybridization and libraries were sequenced on an Illumina MiSeq with 20 cycles of Read 1 to sequence the 12 base UMI and 8 base sample index, and 75 cycles of Read 2 to sequence the targets.

NEW

NEBNext Immune Sequencing Kits (Human & Mouse)


NEBNext Immune Sequencing Kit (Human)			
#E6320S	24 reactions	\$1,680
#E6320L	96 reactions	\$5,760

NEBNext Immune Sequencing Kit (Mouse)			
#E6330S	24 reactions	\$1,680
#E6330L	96 reactions	\$5,760


- Generation of full-length variable sequences (including isotype information), allowing downstream antibody synthesis and functional characterization not possible with approaches sequencing only the CDR3 region.
- Eliminated use of variable region primers, reducing primer pool complexity and realizing unbiased and simultaneous recovery of B Cell and T Cell receptor transcripts.
- Minimized PCR bias and improved sequencing accuracy by allowing a consensus to be generated from duplicate sequencing reads originating from the same transcript; UMIs enable accurate quantitation of each clone present in the sample.
- Optimized high target-capture efficiency for immune repertoire sequencing and analysis from sub-microgram quantities of total RNA.

Description: The NEBNext Immune Sequencing Kits (Human & Mouse) enable exhaustive profiling of somatic mutations in the full-length immune gene repertoires of B cells and T cells, via the expression of complete antibody chains. A unique, UMI-based mRNA barcoding process allows PCR copies derived from an individual molecule to be converted to a consensus sequence. This improves sequence accuracy and eliminates PCR bias, allowing for quantitative digital molecule counting.

Immune repertoire sequencing is frequently used to analyze immune responses, both current and distant. Areas of particular interest include characterization of autoimmune diseases, oncology, discovery of neutralizing antibodies against infectious disease, tumor-infiltrating lymphocytes and use as a tool to study residual disease. Recent improvements in read lengths and throughputs of next-generation sequencing (NGS) platforms have resulted in a rise in the popularity of immune repertoire sequencing.



Simplified representation of the structure of an antibody or TCR. Simplified representation of the structure of an antibody or TCR showing the outcome of V(D)J recombination in mature lymphocytes.



NEBNext Immune Sequencing Kit Workflow.

NEBNext Microbiome DNA Enrichment Kit

#E2612S	6 reactions	\$218
#E2612L	24 reactions	\$744

- Effective separation of microbial DNA from contaminating host DNA
- Fast, simple protocol
- Compatible with downstream applications including next generation sequencing on all platforms, qPCR and end point PCR
- Suitable for a wide range of sample types
- No requirement for live cells
- Captured host DNA can also be eluted and retained


Description: The NEBNext Microbiome DNA Enrichment Kit facilitates separation of microbial DNA from methylated host DNA (including human) by selective binding and removal of the CpG methylated host DNA (1).

Functional Validation: Each set of reagents are functionally validated by enriching *E. coli* DNA from a mixture of *E. coli* and human DNA. Enrichment is evaluated through library construction and sequencing of the enriched sample on an Illumina sequencer.


The Microbiome Enrichment Kit Includes:

- NEBNext MBD2-Fc Protein
- NEBNext Protein A Magnetic Beads
- NEBNext Bind/Wash Buffer (5X)
- 16S RNA Universal Bacteria Control Primers
- RPL30 Human DNA Control Primers


- (1) Feehery, G.R. et al. (2013) *PLoS One*, 8: e76096.
- (2) Chen, T., et al. (2010) *Database*, Vol. 2010, Article ID baq013, doi: 10.1093/database/baq013
- (3) Langmead, B., et al. (2009) *Genome Biol.* 10:R25 doi:10.1186/gb-2009-10-3-r25



Salivary Microbiome DNA Enrichment. DNA was purified from pooled human saliva DNA (Innovative Research) and enriched using the NEBNext Microbiome DNA Enrichment Kit. Libraries were prepared from unenriched and enriched samples and sequenced on the SOLiD 4 platform. The graph shows percentages of 500 M–537 M SOLiD 4 50 bp reads that mapped to either the Human reference sequence (hg19) or to a microbe listed in Human Oral Microbiome Database (HOMD)[2]. (Because the HOMD collection is not comprehensive, ~80% of reads in the enriched samples do not map to either database.) Reads were mapped using Bowtie 0.12.7[3] with typical settings (2 mismatches in a 28 bp seed region, etc.). SOLiD 4 is a registered trademark of Life Technologies, Inc.



NEBNext Microbiome DNA Enrichment Kit workflow.



Microbiome Diversity is Retained after Enrichment with the NEBNext Microbiome DNA Enrichment Kit. DNA was purified from pooled human saliva DNA (Innovative Research) and enriched using the NEBNext Microbiome DNA Enrichment Kit. Libraries were prepared from unenriched and enriched samples, followed by sequencing on the SOLiD 4 platform. The graph shows a comparison between relative abundance of each bacterial species listed in HOMD[2] before and after enrichment with the NEBNext Microbiome DNA Enrichment Kit. Abundance is inferred from the number of reads mapping to each species as a percentage of all reads mapping to HOMD. High concordance continues even to very low abundance species (inset). We compared 501M 50 bp SOLiD 4 reads in the enriched dataset to 537M 50 bp SOLiD 4 reads in the unenriched dataset. Reads were mapped using Bowtie 0.12.7[3] with typical settings (2 mismatches in a 28 bp seed region, etc.).

* *Niesseria flavescens* – This organism may have unusual methylation density, allowing it to bind the enriching beads at a low level. Other *Niesseria* species (*N. mucosa*, *N. sicca* and *N. elongata*) are represented, but do not exhibit this anomalous enrichment.

NEBNext dsDNA Fragmentase®



#M0348S	50 reactions	\$102
#M0348L	250 reactions	\$408

Companion Product:

NEBNext dsDNA Fragmentase	
Reaction Buffer v2	
#B0349S	6 x 1 ml \$31

- Generation of dsDNA fragments for sequencing on next generation sequencing platforms
- Generation of dsDNA fragments for libraries

Description: NEBNext dsDNA Fragmentase generates dsDNA breaks in a time-dependent manner to yield 50–1,000 bp DNA fragments depending on reaction time. NEBNext dsDNA Fragmentase contains two enzymes, one randomly generates nicks on dsDNA and the other recognizes the nicked site and cuts the opposite DNA strand across from the nick, producing dsDNA breaks. The resulting DNA fragments contain short overhangs, 5'-phosphates, and 3'-hydroxyl groups. The random nicking activity of NEBNext dsDNA Fragmentase has been confirmed by preparing libraries for next-generation sequencing. A comparison of the sequencing results between gDNA prepared with NEBNext dsDNA fragmentase and with mechanical shearing demonstrates that the NEBNext dsDNA Fragmentase does not introduce any detectable bias during the sequencing library preparation and no difference in sequence coverage is observed using the two methods.

Source: NEBNext dsDNA Fragmentase is composed of endonucleases isolated from two different *E. coli* sources: one construct expresses a fusion protein consisting of *E. coli* maltose binding protein and *Vibrio vulnificus* nuclease mutant protein; the other expresses a fusion protein consisting of maltose binding protein and T7 endonuclease mutant protein.

Reaction Conditions: 1X NEBNext dsDNA Fragmentase Reaction Buffer v2, supplemented with 100 µM MgCl₂, when required. Incubate at 37°C.

1X NEBNext dsDNA Fragmentase

Reaction Buffer v2:

20 mM Tris-HCl
10 mM MgCl₂
50 mM NaCl
0.15% Triton X-100
pH 7.5 @ 25°C

Reagents Supplied with Enzyme:

10X NEBNext dsDNA Fragmentase Reaction Buffer v2
200 mM MgCl₂

Heat Inactivation: 65°C for 15 minutes in the presence of 50 mM DTT

MINELUTE® is a registered trademark of the Qiagen Group.




Figure 1: Fragmentation of *E. coli* gDNA. *E. coli* gDNA was fragmented with NEBNext dsDNA Fragmentase for the indicated times and purified on MinElute® columns.



Chrystel works at our subsidiary office in France and has been with the team for over 10 years. She spends her time working in Accounting and Human Resources, and also helps out with Customer Service. Chrystel also enjoys interior design and staging homes.

NEBNext Ultra II Q5[®] Master Mix

#M0544S	50 reactions	\$99
#M0544L	250 reactions	\$395

Additional Products:

NEBNext Q5 Hot Start HiFi PCR Master Mix		
#M0543S	50 reactions	\$90
#M0543L	250 reactions	\$360

NEBNext High-Fidelity 2X PCR Master Mix		
#M0541S	50 reactions (50 µl)	\$90
#M0541L	250 reactions (50 µl)	\$360

- Next generation sequencing library preparation
- High-fidelity amplification
- Uniform GC coverage
- Improves sequencing library coverage of known difficult regions of the human genome


Description: The NEBNext Ultra II Q5 Master Mix is a new formulation of Q5 DNA Polymerase that has been optimized for robust, high-fidelity amplification of next-generation sequencing (NGS) libraries. This formulation further improves the uniformity of amplification of libraries, including superior performance with GC-rich regions.

The polymerase component of the master mix, Q5 High-Fidelity DNA Polymerase, is a novel thermostable DNA polymerase that possesses 3'→5' exonuclease activity, and is fused to a processivity-enhancing Sso7d domain. Q5 also has the highest fidelity available (> 100-fold higher than that of *Taq* DNA Polymerase and ~12-fold higher than that of *Pyrococcus furiosus* (Pfu) DNA Polymerase), resulting in ultra-low error rates.

The NEBNext Ultra II Q5 Master Mix is an aptamer-based hot start formulation that allows convenient room temperature reaction set up. The convenient 2X master mix format contains dNTPs, Mg⁺⁺ and a proprietary buffer, and requires only the addition of primers and DNA template for robust amplification. NEBNext Ultra II Q5 Master Mix is also included in the NEBNext Ultra II DNA Library Prep Kit for Illumina.

Source: An *E. coli* strain that carries the Q5 High-Fidelity DNA Polymerase gene.


Reaction Conditions: NEBNext Ultra II Q5 Master Mix, DNA template and 1 µM primers in a total reaction volume of 50 µl.



NEBNext Ultra II Q5 Master Mix provides improved coverage of known low-coverage regions of the human genome. Libraries were prepared from Human NA19240 genomic DNA. One library was not amplified. The other two libraries were amplified using 5 cycles of PCR with NEBNext Q5 Hot Start HiFi PCR Master Mix (NEB #M0543) or with NEBNext Ultra II Q5 Master Mix (NEB #M0544). Libraries were sequenced on an Illumina NextSeq 500. 420 million 75 bp reads were randomly extracted from each dataset, representing an average coverage of 10X. Reads were mapped to the GRCh37 reference genome using Bowtie 2.2.4. Reads on each region were counted using bedtools v2.19.1.

A: The number of reads overlapping distinct low coverage regions of the human genome (1) are shown for each library. **B:** From the 420 million 75 bp reads randomly extracted from each dataset, 10X coverage was expected. The % of difficult regions covered at > 10X are shown for each library. The NEBNext Ultra II Q5 Master Mix provides improved coverage of these known low coverage regions, without drop-outs, and shows similar coverage to the unamplified sample.

(1) Popatov, V. and Ong, J.L. (2017). Examining Sources of Error in PCR by Single-Molecule Sequencing. PLoS ONE. 12(1):e0169774.



NEBNext Ultra II Q5 Master Mix provides uniform GC coverage for microbial genomic DNA with a broad range of GC composition. Libraries were made using 100 ng of the genomic DNAs shown and the NEBNext Ultra II DNA Library Prep Kit. Libraries were amplified using the NEBNext Ultra II Q5 Master Mix, and sequenced on an Illumina MiSeq. GC coverage information was calculated using Picard's CollectGCBiasMetrics (v1.117). Expected normalized coverage of 1.0 is indicated by the horizontal grey line, the number of 100 bp regions at each GC% is indicated by the vertical grey bars, and the colored lines represent the normalized coverage for each library. NEBNext Ultra II Q5 Master Mix provides uniform GC coverage regardless of the GC content of the DNA.

NEBNext Library Quant Kit for Illumina

NEBNext Library Quant Kit for Illumina

#E7630S 100 reactions\$99
#E7630L 500 reactions\$420

Companion Product:

NEBNext Library Dilution Buffer
#B6118S 7.5 ml\$34

- Provides more accurate and reproducible quant values than alternative methods and kits
- Compatible with libraries with a broad range of insert sizes and GC content, made by a variety of methods
- Requires only 4 standards, allowing more libraries to be quantitated per kit
- Supplied with a convenient Library Dilution Buffer
- The NEBNext Library Quant Master Mix requires only the addition of primers
- Utilizes a single extension time for all libraries, regardless of insert size
- Library quant values can be easily calculated using NEB's online tool, at NEBioCalculator.neb.com
- ROX is included in the kit, for use with qPCR instruments that require a reference dye for normalization

With NEBNext, optimal cluster density is achieved from quantitated libraries with a broad range of library size and GC content. Libraries of 310–963 bp from the indicated sources were quantitated using the NEBNext Library Quant Kit, then diluted to 8 pM and loaded onto a MiSeq (v2 chemistry; MCS v2.4.1.3). Library concentrations ranged from 7–120 nM, and resulting raw cluster density for all libraries was 965–1300 k/mm² (ave. =1199). Optimal cluster density was achieved using concentrations determined by the NEBNext Library Quant Kit for all library sizes.

Description: Accurate quantitation of next-generation sequencing libraries is essential for maximizing data output and quality from each sequencing run. For Illumina sequencing specifically, accurate quantitation of libraries is critical to achieve optimal cluster densities, a requirement for optimal sequence performance. qPCR is considered to be the most accurate and effective method of library quantitation, providing considerably higher consistency and reproducibility than electrophoresis or spectrophotometry, which measure total nucleic acid concentration. Amplification-based methods quantitate only those molecules that contain both adaptor sequences, thereby providing a more accurate estimate of the concentration of the library molecules that can be sequenced.



The NEBNext Library Quant Kit delivers significant improvements to qPCR-based library quantitation for next gen sequencing. The NEBNext Library Quant Kit

for Illumina contains components that are optimized for qPCR-based quantitation of libraries prepared for Illumina next-generation sequencing platforms. The NEBNext Library Quant Kit contains primers which target the P5 and P7 Illumina adaptor sequences and a set of high-quality, pre-diluted DNA standards to enable reliable quantitation of diluted DNA libraries between 150–1000 bp.

The Library Quant Kit Includes:

- NEBNext Library Quant Master Mix
- NEBNext Library Quant Primer Mix
- NEBNext Library Quant DNA Standards 1–4
- ROX (Low) and ROX (High)
- NEBNext Library Dilution Buffer



	Reagent Preparation	Library Dilution	Set Up	qPCR	Data Analysis	Total Workflow
Hands-On	5 min.	10 min.	25 min.	1 min.	10 min.	51 min.
Total	5 min.	10 min.	25 min.	60 min.	10 min.	1 hr. 45 min.




Three 340–400 bp libraries were quantitated by 4 different users 2–4 times using either the NEBNext or Kapa Library Quantification Kit (Universal). A notable improvement in quantitation consistency was observed for concentrations determined by the NEBNext Kit (orange) versus those from the Kapa kit (gray).

NEBNext Reagents for Ion Torrent™: DNA Library Preparation

NEBNext kits are available for DNA library preparation for Ion Torrent, with or without enzymatic DNA fragmentation. In addition to stringent QC on individual components, the NEBNext DNA kits are functionally validated by library preparation of a genomic DNA library, followed by Ion Torrent sequencing. Reagent lots are reserved specifically for inclusion in NEBNext kits. Most of these reagents are provided in master mix format, reducing the number of vials provided in the kits, and reducing pipetting steps. Adaptors and primers for singleplex libraries are supplied in the kits. For multiplexed libraries, the Ion XPress™ Barcode Adaptors from Thermo Fisher Scientific can be used.

Input 10 ng – 1 µg*				
Fragmentation	End Repair	Adaptor Ligation/Fill-In	PCR Enrichment	Total Workflow
NEBNext Fast DNA Fragmentation & Library Prep Set for Ion Torrent (NEB #E6285)				<div> Hands-On Time 12 min. Total Time 110 min. – 133 min.</div> <div> Hands-On Time 12 min. Total Time 110 min. – 133 min.</div>
<ul style="list-style-type: none">DNA Fragmentation Master MixDNA Fragmentation Reaction Buffer	<ul style="list-style-type: none">End Repair Enzyme MixEnd Repair Reaction Buffer (10X)	<ul style="list-style-type: none">Adaptors for Ion TorrentT4 DNA LigaseT4 DNA Ligase Buffer for Ion Torrent (10X)Bst 2.0 WarmStart® DNA Polymerase	<ul style="list-style-type: none">Primers for Ion TorrentNEBNext Q5 Hot Start HiFi PCR Master Mix	
NEBNext Fast DNA Library Prep Set for Ion Torrent (NEB #E6270)				
	<ul style="list-style-type: none">End Repair Enzyme MixEnd Repair Reaction Buffer (10X)	<ul style="list-style-type: none">Adaptors for Ion TorrentT4 DNA LigaseT4 DNA Ligase Buffer for Ion Torrent (10X)Bst 2.0 WarmStart DNA Polymerase	<ul style="list-style-type: none">Primers for Ion TorrentNEBNext Q5 Hot Start HiFi PCR Master Mix	

*Note that a minimum of 100 ng is recommended when used in conjunction with Ion Express Barcode Adaptors.



Varying GC Content Libraries. 0.5 µg of DNA from 3 different genomes with varying GC content were used to construct 200 bp and 400 bp libraries using the NEBNext Fast DNA Fragmentation and Library Prep Set for Ion Torrent, analyzed by the Agilent Bioanalyzer.

NEBNext Reagents for DNA Library Preparation – Ordering Information

Illumina Platform:

KITS FOR ILLUMINA DNA LIBRARY PREPARATION		NEB #	SIZE	PRICE
DNA & CHIP	NEBNext Ultra II DNA Library Prep Kit for Illumina	E7645S/L	24/96 rxns	\$590/\$2,245
	NEBNext Ultra II DNA Library Prep with Sample Purification Beads	E7103S/L	24/96 rxns	\$645/\$2,465
	NEBNext Ultra II FS DNA Library Prep Kit for Illumina	E7805S/L	24/96 rxns	\$635/\$2,400
	NEBNext Ultra II FS DNA Library Prep with Sample Purification Beads	E6177S/L	24/96 rxns	\$695/\$2,625
	NEBNext Ultra II DNA PCR-free Library Prep Kit for Illumina	E7410S/L	24/96 rxns	\$560/\$2,132
	NEBNext Ultra II DNA PCR-free Library Prep with Sample Purification Beads	E7415S/L	24/96 rxns	\$612/\$2,341
	NEBNext Ultra II FS DNA PCR-free Library Prep Kit for Illumina	E7430S/L	24/96 rxns	\$603/\$2,280
	NEBNext Ultra II FS DNA PCR-free Library Prep with Sample Purification Beads	E7435S/L	24/96 rxns	\$660/\$2,493
	NEBNext Enzymatic Methyl-seq Kit	E7120S/L	24/96 rxns	\$895/\$3,360
	NEBNext Ultra DNA Library Prep Kit for Illumina	E7370S/L	24/96 rxns	\$565/\$2,160
MODULES & ENZYMES		NEB #	SIZE	PRICE
DNA & CHIP	NEBNext Enzymatic Methyl-seq Conversion Module	E7125S/L	24/96 rxns	\$185/\$675
	NEBNext FFPE DNA Repair Mix	M6630S/L	24/96 rxns	\$165/\$575
	NEBNext Microbiome DNA Enrichment Kit	E2612S/L	6/24 rxns	\$218/\$744
	NEBNext Ultra II FS DNA Module	E7810S/L	24/96 rxns	\$275/\$990
	NEBNext Ultra II End Repair/dA-Tailing Module	E7546S/L	24/96 rxns	\$250/\$795
	NEBNext Ultra II Ligation Module	E7595S/L	24/96 rxns	\$395/\$1,270
	NEBNext Ultra II Q5 Master Mix	M0544S/L	50/250 rxns	\$99/\$395
	NEBNext Ultra End Repair/dA-Tailing Module	E7442S/L	24/96 rxns	\$238/\$761
	NEBNext Ultra Ligation Module	E7445S/L	24/96 rxns	\$379/\$1,219
	NEBNext dsDNA Fragmentase	M0348S/L	50/250 rxns	\$102/\$408
	NEBNext End Repair Module	E6050S/L	20/100 rxns	\$90/\$361
	NEBNext dA-Tailing Module	E6053S/L	20/100 rxns	\$106/\$424
	NEBNext Quick Ligation Module	E6056S/L	20/100 rxns	\$318/\$1,272
	NEBNext Q5 Hot Start HiFi PCR Master Mix	M0543S/L	50/250 rxns	\$90/\$360
	NEBNext High-Fidelity 2X PCR Master Mix	M0541S/L	50/250 rxns	\$90/\$360
	NEBNext Q5U Master Mix	M0597S/L	50/250 rxns	\$125/\$495
	NEBNext dsDNA Fragmentase Reaction Buffer v2	B0349S	6 ml	\$31
ADAPTORS & PRIMERS		NEB #	SIZE	PRICE
	NEBNext Multiplex Oligos for Illumina (Unique Dual Index UMI Adaptors DNA Set 1)	E7395S/L	96/384 rxns	\$672/\$2,553
	NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs)	E6440S/L	96/384 rxns	\$535/\$1,925
	NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs Set 2)	E6442S/L	96/384 rxns	\$535/\$1,925
	NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs Set 3)	E6444S/L	96/384 rxns	\$535/\$1,925
	NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs Set 4)	E6446S/L	96/384 rxns	\$535/\$1,925
	NEBNext Multiplex Oligos for Enzymatic Methyl-seq (Unique Dual Index Primer Pairs)	E7140S/L	24/96 rxns	\$135/\$535
	NEBNext Multiplex Oligos for Illumina (Dual Index Primers Set 1)	E7600S	96 rxns	\$460
	NEBNext Multiplex Oligos for Illumina (Dual Index Primers Set 2)	E7780S	96 rxns	\$460
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 1)	E7335S/L	24/96 rxns	\$105/\$385
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 2)	E7500S/L	24/96 rxns	\$105/\$385
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 3)	E7710S/L	24/96 rxns	\$105/\$385
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 4)	E7730S/L	24/96 rxns	\$105/\$385
	NEBNext Multiplex Oligos for Illumina (96 Index Primers)	E6609S/L	96/384 rxns	\$650/\$2,600
	NEBNext Multiplex Oligos for Illumina (Methylated Adaptor, Index Primers Set 1)	E7535S/L	24/96 rxns	\$115/\$420
	NEBNext Adaptor Dilution Buffer	B1430S	1 x 9.6 ml	\$34
TARGET ENRICHMENT		NEB #	SIZE	PRICE
	NEBNext Direct Custom Ready Panels	E6631S/L/X	8/24/96 rxns	request quote
	NEBNext Direct Genotyping Solution	E9500B-S	96 rxns	request quote
	NEBNext Direct Genotyping Solution	E9530B-S	8 rxns	request quote
	NEBNext Immune Sequencing Kit (Human)	E6320S/L	24/96 rxns	\$1,680/\$5,760
	NEBNext Immune Sequencing Kit (Mouse)	E6330S/L	24/96 rxns	\$1,680/\$5,760
LIBRARY QUANTITATION		NEB #	SIZE	PRICE
	NEBNext Library Quant Kit for Illumina	E7630S/L	100/500 rxns	\$99/\$420
	NEBNext Library Dilution Buffer	B6118S	7.5 ml	\$34

Oxford Nanopore Technologies Platform:

PRODUCTS FOR DNA LIBRARY PREPARATION		NEB #	SIZE	PRICE
DNA	NEBNext Companion Module for Oxford Nanopore Technologies Ligation Sequencing	E7180S	24 rxns	\$920

Ion Torrent Platform:

PRODUCTS FOR DNA LIBRARY PREPARATION		NEB #	SIZE	PRICE
DNA	NEBNext Fast DNA Library Prep Set for Ion Torrent	E6270L	50 rxns	\$891
	NEBNext Fast DNA Fragmentation & Library Prep Set for Ion Torrent	E6285L	50 rxns	\$1,061

Suitable for Any Sequencing Platform:

DNA ENRICHMENT		NEB #	SIZE	PRICE
DNA	NEBNext Microbiome DNA Enrichment Kit	E2612S/L	6/24 rxns	\$218/\$744
DNA REPAIR		NEB #	SIZE	PRICE
DNA	NEBNext FFPE DNA Repair Mix	M6630S/L	24/96 rxns	\$165/\$575
MODULES & ENZYMES		NEB #	SIZE	PRICE
DNA	NEBNext dsDNA Fragmentase	M0348S/L	50/250 rxns	\$102/\$408
	NEBNext Ultra II Q5 Master Mix	M0544S/L	50/250 rxns	\$99/\$395
	NEBNext Q5 Hot Start HiFi PCR Master Mix	M0543S/L	50/250 rxns	\$90/\$360
	NEBNext High-Fidelity 2X PCR Master Mix	M0541S/L	50/250 rxns	\$90/\$360
	NEBNext dsDNA Fragmentase Reaction Buffer v2	B0349S	6 ml	\$31
MAGNETIC SEPARATION		NEB #	SIZE	PRICE
	NEBNext Magnetic Separation Rack	S1515S	24 tubes	\$475

Luis has been with NEB for 24 years and works in our Shipping Department. In addition to making sure NEB orders are packed quickly and efficiently, Luis enjoys snowboarding, saltwater fishing and camping.



NEBNext Reagents for RNA Library Preparation – Ordering Information

Illumina Platform:

KITS FOR ILLUMINA RNA LIBRARY PREPARATION		NEB #	SIZE	PRICE
Directional RNA	NEBNext Ultra II Directional RNA Library Prep Kit for Illumina	E7760S/L	24/96 rxns	\$935/\$3,175
	NEBNext Ultra II Directional RNA Library Prep with Sample Purification Beads	E7765S/L	24/96 rxns	\$1,035/\$3,525
	NEBNext Ultra Directional RNA Library Prep Kit for Illumina	E7420L	96 rxns	\$3,954
Non-directional RNA	NEBNext Ultra II RNA Library Prep Kit for Illumina	E7770S/L	24/96 rxns	\$890/\$3,025
	NEBNext Ultra II RNA Library Prep with Sample Purification Beads	E7775S/L	24/96 rxns	\$985/\$3,350
	NEBNext Ultra RNA Library Prep Kit for Illumina	E7530L	96 rxns	\$3,564
Small RNA	NEBNext Multiplex Small RNA Library Prep Set for Illumina (Set 1)	E7300S/L	24/96 rxns	\$1,470/\$4,998
	NEBNext Multiplex Small RNA Library Prep Set for Illumina (Set 2)	E7580S/L	24/96 rxns	\$1,470/\$4,998
	NEBNext Multiplex Small RNA Library Prep Kit for Illumina (Index Primers 1-48)	E7560S	96 rxns	\$4,998
	NEBNext Small RNA Library Prep Set for Illumina (Multiplex Compatible)	E7330S/L	24/96 rxns	\$1,200/\$4,080
Single Cell	NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina	E6420S/L	24/96 rxns	\$1,195/\$3,965
SARS-CoV-2	NEBNext ARTIC SARS-CoV-2 Library Prep Kit (Illumina)	E7650S/L	24/96 rxns	\$384/\$1,536
	NEBNext ARTIC SARS-CoV-2 FS Library Prep Kit (Illumina)	E7658S/L	24/96 rxns	\$408/\$1,632
MODULES & ENZYMES		NEB #	SIZE	PRICE
RNA	NEBNext RNA Depletion Core Reagent Set	E7865S/L/X	6/24/96 rxns	\$254/\$924/\$3,326
	NEBNext RNA Depletion Core Reagent Set with RNA Sample Purification Beads	E7870S/L/X	6/24/96 rxns	\$262/\$963/\$3,465
	NEBNext Globin & rRNA Depletion Kit	E7750S/L/X	6/24/96 rxns	\$340/\$1,235/\$4,450
	NEBNext Globin & rRNA Depletion Kit (Human/Mouse/Rat) with RNA Sample Purification Beads	E7755S/L/X	6/24/96 rxns	\$350/\$1,290/\$4,635
	NEBNext rRNA Depletion Kit v2 (Human/Mouse/Rat)	E7400S/L/X	6/24/96 rxns	\$297/\$1,080/\$3,885
	NEBNext rRNA Depletion Kit v2 (Human/Mouse/Rat) with RNA Sample Purification Beads	E7405S/L/X	6/24/96 rxns	\$305/\$1,125/\$4,050
	NEBNext rRNA Depletion Kit (Bacteria)	E7850S/L/X	6/24/96 rxns	\$297/\$1,080/\$3,885
	NEBNext rRNA Depletion Kit (Bacteria) with RNA Sample Purification Beads	E7860S/L/X	6/24/96 rxns	\$305/\$1,125/\$4,050
	NEBNext rRNA Depletion Kit (Human/Mouse/Rat)	E6310S/L/X	6/24/96 rxns	\$330/\$1,200/\$4,320
	NEBNext rRNA Depletion Kit (Human/Mouse/Rat) with RNA Sample Purification Beads	E6350S/L/X	6/24/96 rxns	\$340/\$1,250/\$4,500
	NEBNext Poly(A) mRNA Magnetic Isolation Module	E7490S/L	24/96 rxns	\$69/\$247
	NEBNext Magnesium RNA Fragmentation Module	E6150S	200 rxns	\$47
	NEBNext Ultra II RNA First Strand Synthesis Module	E7771S/L	24/96 rxns	\$150/\$480
	NEBNext Ultra II Directional RNA Second Strand Synthesis Module	E7550S/L	24/96 rxns	\$390/\$1,245
	NEBNext Ultra II Non-directional RNA Second Strand Synthesis Module	E6111S/L	20/100 rxns	\$295/\$1,180
	NEBNext RNA First Strand Synthesis Module	E7525S/L	24/96 rxns	\$158/\$507
	NEBNext Single Cell/Low Input cDNA Synthesis and Amplification Module	E6421S/L	24/96 rxns	\$650/\$2,210
	NEBNext Single Cell Lysis Module	E5530S	96 rxns	\$92
DNA	NEBNext Ultra II End Repair/dA-Tailing Module	E7546S/L	24/96 rxns	\$250/\$795
	NEBNext Ultra II Ligation Module	E7595S/L	24/96 rxns	\$395/\$1,270
	NEBNext Ultra End Repair/dA-Tailing Module	E7442S/L	24/96 rxns	\$238/\$761
	NEBNext Ultra Ligation Module	E7445S/L	24/96 rxns	\$379/\$1,219
	NEBNext End Repair Module	E6050S/L	20/100 rxns	\$90/\$361
	NEBNext dA-Tailing Module	E6053S/L	20/100 rxns	\$106/\$424
	NEBNext Quick Ligation Module	E6056S/L	20/100 rxns	\$318/\$1,272
	NEBNext Ultra II Q5 Master Mix	M0544S/L	50/250 rxns	\$99/\$395
	NEBNext Q5 Hot Start HiFi PCR Master Mix	M0543S/L	50/250 rxns	\$90/\$360
	NEBNext High-Fidelity 2X PCR Master Mix	M0541S/L	50/250 rxns	\$90/\$360
ADAPTORS & PRIMERS		NEB #	SIZE	PRICE
	NEBNext Multiplex Oligos for Illumina (Unique Dual Index UMI Adaptors RNA Set 1)	E7416S/L	96/384 rxns	\$648/\$2,460
	NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs)	E6440S/L	96/384 rxns	\$535/\$1,925
	NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs Set 2)	E6442S/L	96/384 rxns	\$535/\$1,925
	NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs Set 3)	E6444S/L	96/384 rxns	\$535/\$1,925
	NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs Set 4)	E6446S/L	96/384 rxns	\$535/\$1,925
	NEBNext Multiplex Oligos for Illumina (Dual Index Primers Set 1)	E7600S	96 rxns	\$460
	NEBNext Multiplex Oligos for Illumina (Dual Index Primers Set 2)	E7780S	96 rxns	\$460
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 1)	E7335S/L	24/96 rxns	\$105/\$385
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 2)	E7500S/L	24/96 rxns	\$105/\$385
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 3)	E7710S/L	24/96 rxns	\$105/\$385
	NEBNext Multiplex Oligos for Illumina (Index Primers Set 4)	E7730S/L	24/96 rxns	\$105/\$385
	NEBNext Multiplex Oligos for Illumina (96 Index Primers)	E6609S/L	96/384 rxns	\$650/\$2,600
	NEBNext Adaptor Dilution Buffer	B1430S	1 x 9.6 ml	\$34

Illumina Platform (continued):

LIBRARY QUANTITATION		NEB #	SIZE	PRICE
	NEBNext Library Quant Kit for Illumina	E7630S/L	100/500 rxns	\$99/\$420
	NEBNext Library Dilution Buffer	B6118S	7.5 ml	\$34

Oxford Nanopore Technologies Platform:

PRODUCTS FOR DNA LIBRARY PREPARATION		NEB #	SIZE	PRICE
RNA	NEBNext ARTIC SARS-CoV-2 Companion Kit (Oxford Nanopore Technologies)	E7660S/L	24/96 rxns	\$192/\$768

Suitable for Any Sequencing Platform:

MODULES & ENZYMES		NEB #	SIZE	PRICE
RNA	NEBNext rRNA Depletion Kit (Human/Mouse/Rat)	E6310S/L/X	6/24/96 rxns	\$330/\$1,200/\$4,320
	NEBNext rRNA Depletion Kit (Human/Mouse/Rat) with RNA Sample Purification Beads	E6350S/L/X	6/24/96 rxns	\$340/\$1,250/\$4,500
	NEBNext Poly(A) mRNA Magnetic Isolation Module	E7490S/L	24/96 rxns	\$69/\$247
	NEBNext Magnesium RNA Fragmentation Module	E6150S	200 rxns	\$47
	NEBNext Ultra II RNA First Strand Synthesis Module	E7771S/L	24/96 rxns	\$150/\$480
	NEBNext Ultra II Non-directional RNA Second Strand Synthesis Module	E6111S/L	20/100 rxns	\$295/\$1,180
	NEBNext Globin & rRNA Depletion Kit	E7750S/L/X	6/24/96 rxns	\$340/\$1,235/\$4,450
	NEBNext Globin & rRNA Depletion Kit (Human/Mouse/Rat) with RNA Sample Purification Beads	E7755S/L/X	6/24/96 rxns	\$350/\$1,290/\$4,635
MAGNETIC SEPARATION		NEB #	SIZE	PRICE
	NEBNext Magnetic Separation Rack	S1515S	24 tubes	\$475

Featured Online Tools**NEBNext Selector**

Use this tool to guide you through selection of NEBNext reagents for next generation sequencing sample preparation. Try it out at NEBNextSelector.neb.com

NEBNext Custom RNA Depletion Design Tool

This tool designs probes to be used with the NEBNext RNA Depletion Core Reagent Set or to supplement an existing NEBNext depletion kit for the depletion of unwanted RNA species. Try it out at <https://depletion-design.neb.com/>.





Be Adventurous

Nurturing a sense of adventure can take us out of our comfort zone, as we embrace new challenges and take risks. It is often easier to choose the familiar path, but in doing so we miss the opportunity to explore something new. It is our innate curiosity that pushes us to pursue new challenges, which can often lead to the most exciting and memorable experiences, and perhaps the opportunity to learn something new.

Don Comb was an adventurer who loved to be outdoors. Trout fishing was one of his most passionate pursuits, and entomology one of his greatest fascinations. As an adult, he enjoyed sailing with his children. They love to regale listeners with stories of their maritime escapades — sailing from St Barthélemy Island in the Caribbean to New England, for example, guided only by Don's celestial navigation skills and fearless determination. Don's family embraced his boundless curiosity and enthusiasm when visiting remote villages in Papua New Guinea on uncharted trekking expeditions, or traveling to the Bermuda Biological Station to snorkel and collect sea urchins for use in developmental research in his position at Harvard Medical School.

Don's wonderment of the world compelled him to question everything. He was always the scientist, even when on vacation, and he would often return with new questions and a sample or two, asking NEB researchers to investigate a specific novelty. A striped Bass that he caught in New England waters was in the process of digesting a lobster, so Don asked researchers at NEB to find the enzyme capable of breaking down the lobster shell. There are many stories like this, and while some of these projects may have seemed outlandish, they often resulted in techniques and findings that were subsequently applied to other research projects at NEB.

NEB researchers and staff are encouraged to bring their curiosity to work, where asking big questions is a vital part of the workplace culture.

In the early 1990s, Chris Taron, now a Scientific Director at NEB, recalls a time when Don was inspired to search for samples in the geothermal springs in Yellowstone National Park. This was the location where *Thermus aquaticus*, the temperature tolerant bacteria that produces *Taq* DNA Polymerase, was first isolated. Don reasoned that the *T. aquaticus* ecosystem might also harbor a bacteriophage that encodes a polymerase exhibiting equivalent performance to *Taq*. Chris, also an outdoor enthusiast, embraced the offer to hike, cross-country ski and snowmobile in the middle of winter to the Lower Geiser Basin. Using a golf ball retriever jerry-rigged with water sampling equipment, he collected samples and transported them back to NEB for analysis. The result — researchers identified the bacteriophage, which unfortunately did not possess the polymerase they had sought. But finding a marketable product was not Don's primary motivation, he just had a question that he wanted an answer to.

The drive behind Don's desire to be outdoors seeking adventure and his endless scientific questioning was the same — a voracious inquisitiveness. Scientific inquiry is driven by a desire to better understand the world around us, and to embrace the innovation that results from being guided by curiosity.

Exploring your adventurous side doesn't have to involve scaling rock walls or taking dangerous risks. Only you can define what adventure means to you. It can be anything that feels a little out of your routine — as long as it gets you excited!

Don maneuvers his sailboat while his son Michael watches for coral from the bow, Exuma Cays, Bahama Islands, 1987.
Credit: Dave Comb

Learn more about
Don's adventures
with his family.



Markers & Ladders (DNA, RNA & Protein)

A wide range of ladders and markers to meet your macromolecule quantification needs.

New England Biolabs provides a wide range of ladders and markers, featuring exceptional quality, uniform band intensities, convenient band spacing and easy-to-identify reference bands.

Double-stranded DNA markers are available for conventional electrophoresis. Conventional electrophoresis markers (size range: ~ 10 to 2.3×10^4 bp) include: HindIII and BstEII digests of Lambda DNA, BstNI and MspI digests of pBR322 DNA, and a HaeIII digest of ϕ X174 DNA.

We also supply a series of DNA ladders ranging from 10 bp to 48.5 kb. Our 100 bp DNA ladder, 1 kb DNA Ladder and 1 kb Plus DNA Ladder are available in several formats: Conventional, Quick-Load® using either non-fluorescing purple dye or bromophenol blue as a tracking dye, and TriDye™ containing three dyes to track gel migration. Our Quick-Load Purple DNA ladders utilize our purple loading dye that improves band visibility by casting no UV shadow and are supplied with a vial of purple dye.

Our RNA ladders and markers have a size range of 17 to 9,000 bases. Both ssRNA ladders are supplied with 2X sample buffer and feature a higher intensity fragment to serve as a reference band. The dsRNA ladders are suitable for use as a size standard in dsRNA and RNAi analysis on both denaturing polyacrylamide and agarose gels.

For your protein analysis needs, NEB offers a selection of highly pure protein standards. Sizes range from 10–250 kDa, which is ideal for accurate molecular weight determination for a wide range of expressed proteins.

Featured Products

- 182** Quick-Load Purple 1 kb Plus DNA Ladder
- 185** Color Prestained Protein Standard, Broad Range (10–250 kDa)
- 182** 1 kb Plus DNA Ladder for Safe Stains
- 182** TriDye Ultra Low Range DNA Ladder

Featured Tools & Resources



Visit www.neb.com/DNALadders to find selection charts for NEB's DNA markers and ladders.

DNA MARKERS & LADDERS

Size Ranges of DNA Ladders	180
Gel Loading Dye, Purple (6X)	180
Gel Loading Dye, Purple (6X), No SDS	180

DNA Ladders

1 kb DNA Ladder	181
100 bp DNA Ladder	181
1 kb Plus DNA Ladder	181
50 bp DNA Ladder	181
Low Molecular Weight DNA Ladder	181
PCR Marker	181
Supercoiled DNA Ladder	184

Quick-Load Purple Formats

Quick-Load Purple 1 kb Plus DNA Ladder	182
Quick-Load Purple 1 kb DNA Ladder	182
Quick-Load Purple 100 bp DNA Ladder	182
Quick-Load Purple 50 bp DNA Ladder	182
Quick-Load Purple Low Molecular Weight DNA Ladder	182

Other Ready-to-Load Formats

1 kb Plus DNA Ladder for Safe Stains	182
Fast DNA Ladder	182

TriDye Formats

TriDye DNA Ladders	182
TriDye Ultra Low Range DNA Ladder	182

Quick-Load Formats

Quick-Load DNA Ladders	182
Quick-Load 1 kb Extend DNA Ladder	182

PFG Ladders

Lambda PFG Ladder	182
MidRange PFG Ladder	182
Lambda DNA-Mono Cut Mix	182

Conventional DNA Markers

Lambda DNA-HindIII Digest	183
Lambda DNA-BstEII Digest	183
φX174 DNA-HaeIII Digest	183
pBR322 DNA-BstNI Digest	183
pBR322 DNA-MspI Digest	183


RNA MARKERS & LADDERS

dsRNA Ladder	184
microRNA Marker	184
ssRNA Ladder	184
Low Range ssRNA Ladder	184

PROTEIN STANDARDS

Unstained Protein Standard, Broad Range (10–200 kDa)	185
Color Prestained Protein Standard, Broad Range (10–250 kDa)	185
Blue Prestained Protein Standard, Broad Range (11–250 kDa)	185

Size Ranges of DNA Ladders from NEB




Purple Loading Dye

Gel Loading Dye, Purple (6X)
#B7024S 4.0 ml\$45

Gel Loading Dye, Purple (6X), no SDS
#B7025S 4.0 ml\$45

Our Gel Loading Dye, Purple (6X) (with and without SDS) is now supplied with all unstained DNA Ladders, sharpens bands and eliminates the UV shadow seen with other dyes. These pre-mixed loading buffers contain a combination of two dyes, Dye 1 (pink/red) and Dye 2 (blue). The red dye serves as the tracking dye for both agarose and non-denaturing polyacrylamide gel electrophoresis. The two dyes separate upon electrophoresis; the red band is the major indicator and migrates similarly to Bromophenol Blue on agarose gels. Specifically chosen, this dye does not leave a shadow under UV light. EDTA is also included to chelate magnesium (up to 10 mM) in enzymatic reactions, thereby stopping the reaction. The dyes also contain Ficoll, which creates brighter and tighter bands when compared to glycerol loading dyes. Gel Loading Dye, Purple (6X) contains SDS, which often results in sharper bands, as some restriction enzymes are known to remain bound to DNA following cleavage.



DNA Ladders

1 kb DNA Ladder

#N3232S	200 gel lanes\$55
#N3232L	1,000 gel lanes\$216

100 bp DNA Ladder

#N3231S	100 gel lanes\$59
#N3231L	500 gel lanes\$235

1 kb Plus DNA Ladder

#N3200S	100-200 gel lanes\$59
#N3200L	500-1,000 gel lanes\$235

50 bp DNA Ladder

#N3236S	100-200 gel lanes\$73
#N3236L	500-1,000 gel lanes\$290

Low Molecular Weight DNA Ladder


#N3233S	100 gel lanes\$69
#N3233L	500 gel lanes\$275

PCR Marker

#N3234S	100 gel lanes\$69
#N3234L	500 gel lanes\$276


NEB offers a variety of DNA Ladders with sizes ranging from 10 bp to 48.5 kb for use in agarose gel electrophoresis.

- Stable at room temperature
- Sharp, uniform bands
- Easy-to-identify reference bands
- Supplied with 1 vial of Gel Loading Dye, Purple (6X), no SDS
- Can be used for sample quantification (see www.neb.com for mass values)




1 kb DNA Ladder

1 kb DNA Ladder visualized by ethidium bromide staining on a 0.8% TAE agarose gel.




100 bp DNA Ladder

100 bp DNA Ladder visualized by ethidium bromide staining on a 1.3% TAE agarose gel.




1 kb Plus DNA Ladder

1 kb Plus DNA Ladder visualized by ethidium bromide staining on a 1.0% TBE agarose gel.




50 bp DNA Ladder

1.0 µg of 50 bp DNA Ladder visualized by ethidium bromide staining on a 3% TBE agarose gel.




Low Molecular Weight DNA Ladder

0.5 µg of LMW DNA Ladder visualized by ethidium bromide staining on a 3% TBE agarose gel.




TriDye Ultra Low Range DNA Ladder

0.5 µg of Ultra Low Range DNA Ladder visualized by ethidium bromide staining on a 20% polyacrylamide gel.




Fast DNA Ladder*

1.2% TBE agarose gel.




PCR Marker

0.3 µg of PCR Marker visualized by ethidium bromide staining on a 1.8% TBE agarose gel.



Quick-Load 1 kb Extend DNA Ladder

0.5% TBE agarose gel.



1 kb Plus DNA Ladder for Safe Stains**

1.0% TBE agarose gel.

Usage Notes: Dilution of these markers is recommended in TE or other buffer of minimal ionic strength. DNA may denature if diluted in dH₂O.

Individual DNA fragments, as seen in the 1 kb DNA Ladder, are available on request.

Contact info@neb.com for more information.

* Fast DNA Ladder can be used for fast electrophoresis systems as well as standard electrophoresis. It is in a ready-to-load format with a xylene cyanol dye.

** 1 kb Plus DNA Ladder for Safe Stains, optimized for use with GelRed®, GelGreen® and SYBR® precast gels.

DNA Ladders in Convenient Pre-mixed Formats

Quick-Load® Purple Formats:

Quick-Load Purple 1 kb Plus DNA Ladder*
 #N0550S 125-250 gel lanes\$75
 #N0550L 375-750 gel lanes\$187

Quick-Load Purple 1 kb DNA Ladder*
 #N0552S 125 gel lanes\$64
 #N0552L 375 gel lanes\$131

Quick-Load Purple 100 bp DNA Ladder*
 #N0551S 125 gel lanes\$88
 #N0551L 375 gel lanes\$216

Quick-Load Purple 50 bp DNA Ladder*
 #N0556S 125-250 gel lanes\$89

Quick-Load Purple Low Molecular Weight DNA Ladder*
 #N0557S 125 gel lanes\$90

* Supplied with 1 vial of Gel Loading Dye, Purple (6X), no SDS.

Other Ready-to-Load Formats:

NEW
 1 kb Plus DNA Ladder for Safe Stains
 #N0559S 125-250 gel lanes\$76

Fast DNA Ladder
 #N3238S 50-200 gel lanes\$63

TriDye™ Formats:

TriDye 1 kb Plus DNA Ladder
 #N3270S 125-250 gel lanes\$85

TriDye 1 kb DNA Ladder
 #N3272S 125 gel lanes\$70

TriDye 100 bp DNA Ladder
 #N3271S 125 gel lanes\$94

NEW

TriDye Ultra Low Range DNA Ladder*
 #N0558S 125-250 gel lanes\$102

* TriDye Ultra Low Range DNA Ladder is suitable for use on both native polyacrylamide and agarose gels.

Quick-Load Formats:

Quick-Load 1 kb Plus DNA Ladder
 #N0469S 125-250 gel lanes\$78

Quick-Load 1 kb DNA Ladder
 #N0468S 125 gel lanes\$66
 #N0468L 375 gel lanes\$156

Quick-Load 1 kb Extend DNA Ladder
 #N3239S 125 gel lanes\$75

Quick-Load 100 bp DNA Ladder
 #N0467S 125 gel lanes\$89
 #N0467L 375 gel lanes\$218

- Ready-to-load
- Uniform band intensities
- Easy-to-identify reference bands
- Defined mass profile for sample quantification



Our 1 kb Plus, 1 kb and 100 bp DNA Ladders are offered in four formats. Choose from the conventional ladder, the Quick-Load version using either non-fluorescing, purple dye or bromophenol blue as a tracking dye, or TriDye containing three dyes to facilitate monitoring of gel migration.

PFG Ladders

Lambda PFG Ladder
 #N0341S 50 gel lanes\$159

MidRange PFG Marker
 #N0342S 50 gel lanes\$159


Lambda DNA-Mono Cut Mix
 #N3019S 100 gel lanes\$76

The Lambda PFG Ladder consists of one GelSyringe™ dispenser, sufficient for 50 gel lanes. Successively larger concatemers of lambda DNA (cl857 ind1 Sam7) are embedded in 1% LMP agarose. Size range: 48.5–1,018 kb.

MidRange PFG Marker consists of concatemers of λ DNA isolated from the bacteriophage λ (cl857 ind1 Sam7) mixed with XhoI digested λ DNA embedded in 1% LMP agarose and supplied in a GelSyringe dispenser. XhoI

produces fragments of 15.0 and 33.5 kb. These fragments anneal to and form concatemers with intact λ DNA. Size range: 15–291 kb.

The Lambda DNA-Mono Cut Mix is best separated by pulsed field gel electrophoresis, but can be alternatively used with standard electrophoresis systems. It is supplied in a liquid format. Size range: 1.5–48.5 kb.



Conventional DNA Markers

Lambda DNA-HindIII Digest

#N3012S 150 gel lanes\$69

#N3012L 750 gel lanes\$276

Lambda DNA-BstEII Digest

#N3014S 150 gel lanes\$76

φX174 DNA-HaeIII Digest

#N3026S 50 gel lanes\$82

#N3026L 250 gel lanes\$326

pBR322 DNA-BstNI Digest

#N3031L 250 gel lanes\$308

pBR322 DNA-MspI Digest


#N3032S 50 gel lanes\$77

NEB offers a wide range of double-stranded DNA molecular weight markers for conventional agarose gel electrophoresis. These standards have a size range of approximately 10–23,000 base pairs.

The typical pattern generated by each of the conventional markers is shown below. The number of fragments generated by each marker, as well as the specific fragment sizes for each of the conventional markers can be found online.

Usage Recommendation: Dilution of these markers is recommended in TE or other buffer of minimal ionic strength. DNA may denature if diluted in dH₂O.


The cohesive ends of fragments 1 and 4 of the Lambda DNA-HindIII and Lambda DNA-BstEII Digests can be separated by heating to 60°C for 3 minutes.



Supercoiled DNA Ladder

#N0472S 100 gel lanes\$107

Mass (ng)	Kilobases
45	10.0
45	8.0
45	6.0
136	5.0
45	4.0
45	3.5
45	3.0
45	2.5
45	2.0



Supercoiled DNA Ladder

0.5 µg/lane.
0.8% TAE agarose gel.

The Supercoiled DNA ladder contains 9 proprietary supercoiled plasmids, ranging in size from 2 to 10 kb, that are suitable for use as supercoiled molecular weight standards for agarose electrophoresis. The 5 kb plasmid has an increased intensity to serve as a reference band.

Source: The 9 proprietary plasmids are purified, phenol extracted and equilibrated to 10 mM Tris-HCl (pH 8.0) and 1 mM EDTA.

Concentration: 500 µg/ml.

Notes: This ladder may contain some traces of nicked DNA and dimers above the 10 kb plasmid. To minimize nicking of the supercoiled DNA, always use sterile pipette tips and avoid multiple freeze-thaw cycles. The migration of supercoiled plasmids in agarose gels can change depending on agarose concentration, buffer and electrophoresis conditions. Dilute in TE or other buffer of minimal ionic strength. DNA may denature if diluted in dH₂O.

Usage Recommendation: Centrifuge briefly and mix gently before use. We recommend loading 0.5 µg (1 µl) of the Supercoiled DNA Ladder diluted in sample buffer. This ladder was not designed for precise quantification of DNA mass, but can be used for approximating the mass of DNA in comparably intense samples of similar size. The approximate mass of DNA in each of the bands in our Supercoiled DNA ladder is as follows (assuming a 0.5 µg loading):

Band	Base Pairs	DNA Mass
1	10,000	45 ng
2	8,000	45 ng
3	6,000	45 ng
4	5,000	136 ng
5	4,000	45 ng
6	3,500	45 ng
7	3,000	45 ng
8	2,500	45 ng
9	2,017	45 ng

RNA Markers & Ladders

dsRNA Ladder
#N0363S 25 gel lanes\$96

microRNA Marker
#N2102S 100 gel lanes\$72

ssRNA Ladder
#N0362S 25 gel lanes\$72


Low Range ssRNA Ladder
#N0364S 100 gel lanes\$71

NEB offers several RNA Markers and Ladders with a size range from 17 to 9,000 bases. The Low Range ssRNA Ladder and the ssRNA Ladder are suitable for use as RNA size standards on denaturing or native gels. Both are supplied with RNA Loading Dye (2X) (NEB #B0363) and feature a higher intensity fragment to serve as a reference band. The microRNA Marker, provided in a ready-to-load denaturing solution, is ideally used as a size marker on denaturing polyacrylamide gels or northern blots and is best visualized stained with SYBR®-Gold. It is supplied

with a 3'-biotinylated 21-mer oligonucleotide probe that can be labeled with γ³²P-ATP and T4 PNK (NEB #M0201). The dsRNA Ladder is suitable for use as a size standard in dsRNA and RNAi analysis on both polyacrylamide and agarose gels.

Concentration: Low Range ssRNA Ladder and dsRNA Ladder are supplied at 500 µg/ml. ssRNA Ladder is supplied at 2,000 µg/ml. MicroRNA Marker is supplied at 12 ng/µl.

SYBR® is a registered trademark of Molecular Probes, Inc.



Protein Standards


Unstained Protein Standard, Broad Range (10–200 kDa)		
#P7717S	150 gel lanes	\$133
#P7717L	750 gel lanes	\$531
Color Prestained Protein Standard, Broad Range (10–250 kDa)		
#P7719S	150 gel lanes	\$174
#P7719L	750 gel lanes	\$735
Blue Prestained Protein Standard, Broad Range (11–250 kDa)		
#P7718S	150 gel lanes	\$153
#P7718L	750 gel lanes	\$653

Companion Products:	
Blue Protein Loading Dye	
3X Blue Loading Buffer (8 ml)	
30X Reducing Agent (1 ml)	
#B7703S	\$30

NEB offers a selection of highly pure protein standards available as unstained, blue prestained or color prestained (containing two colored reference bands for easy identification). Sizes range from 10 to 250 kDa which is ideal for calculating molecular weight determination for a wide range of expressed proteins. NEB protein standards provide uniform band intensities, convenient band spacing and easy-to-identify reference bands.

Recommended Load Volume: 3 µl

Note: For calculating molecular weight determinations, use NEB's Unstained Protein Standard, Broad Range.







Advancement of Science

Don Comb knew that the cultivation of scientific knowledge was critical to driving scientific breakthroughs. He took great joy in watching others learn and grow, and prioritized education and the sharing of knowledge and resources. While NEB has generously donated to numerous causes over the years, Don always felt that sharing science was the purest way to help communities and have the most significant impact.

Supporting the next generation of scientists, both locally and globally, has always been central to NEB's philosophy. This is done through hosting molecular biology courses, welcoming international students and scholars to the NEB campus, maintaining a rigorous postdoctoral program, and managing an active summer student intern program for college students. From the very beginning, NEB has donated reagents to support life science research in schools, community colleges, public and private universities, lab teaching programs, and synthetic biology competitions such as iGEM™ and BioBuilder®. From providing answers through technical support and the NEB Catalog to the sharing of online tools and educational content on our website, NEB has always made access to its scientific information a priority.

In the early days of NEB, Senior Scientist Rick Morgan's role was to find new enzymes in nature with exciting functionality to add to NEB's growing portfolio. In doing this, it occurred to both Don and Rick that the methods for screening the new samples involved relatively simple biology and could therefore be conducted in low resource countries. This, in turn, would help support research in communities that would otherwise not have that opportunity. Researchers from five developing countries came to NEB headquarters to learn molecular biology techniques that would allow them to analyze new samples and support their work. NEB also helped these researchers establish labs in their own countries — over a period of three decades these thriving labs discovered many new restriction endonucleases.

For over 35 years, NEB has collaborated with Smith College in Northampton, MA, to offer a two-week molecular biology program for scientists, lawyers, artists, graduate students, medical doctors, and the generally curious. In fact, anyone with an interest in molecular biology is welcome to participate. NEB provides sponsorships for teachers and people from developing countries. Often, the only way to learn is to roll up one's sleeves, and so the course is largely a practical laboratory experience, teaching foundational techniques and a big picture understanding of the molecular biology toolbox. This course has been attended by over 4,000 students to date.

More recently, "Kids Science Days" were introduced at NEB headquarters in Ipswich, MA. Local students (and teachers) are invited to participate in hands on experiments led by NEB scientists and staff. Seeing the curiosity in the eyes of the students, and being a part of that "aha" moment reminds us all why we are here.

Don's willingness to share what he had and what he knew has created a ripple effect of knowledge throughout research communities, including communities with fewer resources. He reminds us that we must all play a part in the advancement of science and supporting future generations of scientists.

Students participating in a demonstration at
Ipswich High School, hosted by NEB.
Credit: Bree Hall, NEB

Learn more about
NEB's educational
course support.



BIOBUILDER® is a registered trademark of Biobuilder, Inc.

Genome Editing



Programmable nucleases for your applications.

Easily changing the sequence specificity of a DNA binding protein enables many new possibilities for the detection and manipulation of DNA, including for genome editing – creating targeted changes in the DNA of living cells. CRISPR (Clustered Regularly Interspaced Short Palindromic Repeats) -associated (Cas) nucleases have been adapted from bacterial immune systems into useful tools for biotechnology. Cas nucleases including Cas9 and Cas12a, are attractive for genome editing because they are easily programmable with a single guide RNA or crRNA/tracrRNA to introduce a double-stranded break at a specific target. The resultant breaks are repaired by cellular machinery, in some applications using an exogenous repair template.

The ease of programming Cas nucleases, and their conversion into nicking endonucleases, or DNA binding proteins without nuclease activity has expanded their use to include delivering a specific cargo to a locus for applications including visualization, activation, repression and base editing. Furthermore, Cas nuclease and their variants are useful tools for the detection and manipulation of DNA *in vitro*.

Featured Products

- 192** EnGen Mutation Detection Kit
- 193** EnGen sgRNA Synthesis Kit, *S. pyogenes*
- 190** EnGen Lba Cas12a (Cpf1)
- 190** EnGen Sau Cas9

Featured Tools & Resources

















Visit www.neb.com/GenomeEditing for more information, including our feature article and latest brochure.



Learn more about genome editing.



FEATURED PRODUCTS SUPPORTING CRISPR WORKFLOWS

 EnGen Spy Cas9 NLS	190, 191	 Q5 Site-Directed Mutagenesis Kit (with or without competent cells)	98
EnGen Mutation Detection Kit	190, 192	 Q5 High-Fidelity DNA Polymerases	66
EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i>	190, 193	 NEBuilder HiFi DNA Assembly Master Mix	94
 EnGen Spy Cas9 Nickase	190, 191	 NEBuilder HiFi DNA Assembly Cloning Kit	94
 EnGen Spy dCas9 (SNAP-tag)	190, 191	 HiScribe T7 ARCA mRNA Kit (with or without tailing)	200
 EnGen Lba Cas12a (Cpf1)	190, 191	 HiScribe T7 High Yield RNA Synthesis Kit	199
 EnGen Sau Cas9	190, 191	 HiScribe T7 Quick High Yield RNA Synthesis Kit	199
<i>Tth</i> Argonaute (TtAgo)	191	 T7 Endonuclease I	122
 Cas9 Nuclease, <i>S. pyogenes</i>	190, 191		
Monarch Total RNA Miniprep Kit	137		
Monarch RNA Cleanup Kit (50 µg)	138		

 Recombinant Enzyme

Featured NEB Products Supporting CRISPR Workflows

New England Biolabs provides reagents to support a broad variety of CRISPR/Cas genome editing approaches. From introduction of Cas and single guide RNA (sgRNA) on plasmids, to direct introduction of Cas ribonucleoprotein (RNP) and detection of edits using next generation sequencing or enzymatic mutation detection, NEB provides reagents that simplify and shorten genome editing workflows.

PRODUCT	CRISPR/Cas9 APPLICATION	NEB #	SIZE	PRICE
EnGen® Spy Cas9 NLS	<i>in vitro</i> cleavage of dsDNA. Genome engineering by direct introduction of active ribonucleoproteins. Recognizes 5'-NGG-3' PAM.	M0646T/M	400/2,000 pmol	\$150/\$600
EnGen Mutation Detection Kit	Determination of the targeting efficiency of genome editing protocols	E3321S	25 rxns	\$208
EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i>	Generation of microgram quantities of custom sgRNA. Recognizes 5'-NGG-3' PAM.	E3322V/S	10/20 rxns	\$215/\$416
EnGen Spy Cas9 Nickase	<i>in vitro</i> nicking of dsDNA. Genome engineering by direct introduction of active nickase complexes. Recognizes 5'-NGG-3' PAM.	M0650S/T	70/400 pmol	\$60/\$150
EnGen Spy dCas9 (SNAP-tag®)	Programmable binding of DNA. Compatible with SNAP-tag substrates for visualization and enrichment. Recognizes 5'-NGG-3' PAM.	M0652S/T	70/400 pmol	\$60/\$150
EnGen Lba Cas12a (Cpf1)	<i>in vitro</i> cleavage of dsDNA. Genome engineering by direct introduction of active nuclease complexes. Signal generation for sequence detection assays. Recognizes 5'-TTN PAM.	M0653S/T	70/2,000 pmol	\$70/\$250
EnGen Sau Cas9	<i>in vitro</i> cleavage of dsDNA. Genome engineering by direct introduction of active nuclease complexes. Recognizes 5'-NNGRRT-3' PAM.	M0654S/T	70/400 pmol	\$60/\$150
Cas9 Nuclease, <i>S. pyogenes</i>	<i>in vitro</i> cleavage of dsDNA. Genome engineering by direct introduction of active ribonucleoproteins.	M0386S/T/M	70/400/2,000 pmol	\$50/\$135/\$540
Monarch® Total RNA Miniprep Kit	Purification of total RNA, with a binding capacity of up to 100 µg	T2010S	50 preps	\$250
Monarch RNA Cleanup Kit (50 µg)	Purification of sgRNA, with a capacity of up to 50 µg	T2040S/L	10/100 preps	\$54/\$290
Q5® Site-Directed Mutagenesis Kit (with or without competent cells)	Insertion of target sequence into a Cas9-sgRNA construct and modification of HDR templates	E0554S/E0552S	10 rxns	\$201/\$139
Q5 High-Fidelity DNA Polymerases	High-fidelity construct generation for use with CRISPR workflows and for sequencing	Multiple	Multiple	see page 66
NEBuilder® HiFi DNA Assembly Master Mix	Single-tube, isothermal generation of the Cas9-sgRNA construct and HDR templates	E2621S/L/X	10/50/250 rxns	\$159/\$630/\$2520
NEBuilder HiFi DNA Assembly Cloning Kit	Single-tube, isothermal generation of the Cas9-sgRNA construct and HDR templates	E5520S	10 rxns	\$188
HiScribe T7 ARCA mRNA Kit (with or without tailing)	Generation of Cas9 mRNA with ARCA cap	E2060S/E2065S	20 rxns	\$384/\$328
HiScribe T7 High Yield RNA Synthesis Kit	Generation of sgRNA and Cas9 mRNA	E2040S	50 rxns	\$232
HiScribe T7 Quick High Yield RNA Synthesis Kit	Generation of sgRNA and Cas9 mRNA	E2050S	50 rxns	\$284
T7 Endonuclease I	Determination of the editing efficiency of genome editing experiments	M0302S/L	250/1,250 units	\$71/\$284

Programmable Nucleases

The highest efficiency strategy for genome editing with CRISPR/Cas nucleases is direct introduction of Cas/guide RNA complexes. This method further simplifies CRISPR/Cas workflows and has been reported to increase on-target editing activity and reduce off-target events. NEB provides purified Cas9 nucleases from *S. pyogenes* and *S. aureus*, and Cas12a (Cpf1) nuclease from *Lachnospiraceae* bacterium ND2006. In addition, NEB provides variants of Cas9 from *S. pyogenes*, including nicking endonuclease and endonuclease deficient versions. NEB also provides *Thermus thermophilus* argonaute (TtAgo), a programmable DNA endonuclease which requires a short 5'-phosphorylated single-stranded DNA guide to target its activity to a specific corresponding sequence on a substrate.

PRODUCT	NEB #	FEATURES	SIZE	PRICE
Cas9 Nuclease, <i>S. pyogenes</i>	M0386S/T/M	<ul style="list-style-type: none"> Ideal for <i>in vitro</i> digestion of dsDNA Compatible with EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i> (NEB #E3322) and the EnGen Mutation Detection Kit (NEB #E3321) For help with oligo design, try EnGen sgRNA Template Oligo Designer 	70/400/2,000 pmol	\$50/\$135/\$540
EnGen Spy Cas9 NLS	M0646T/M	<ul style="list-style-type: none"> Ideal for direct introduction of Cas9/sgRNA complexes Dual NLS for improved transport to the nucleus Compatible with EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i> (NEB #E3322) and the EnGen Mutation Detection Kit (NEB #E3321) 	400/2,000 pmol	\$150/\$600
EnGen Spy Cas9 Nickase	M0650S/T	<ul style="list-style-type: none"> Variant of Cas9 nuclease differing by a point mutation (D10A) in the RuvC nuclease domain Capable of generating nicks, but not cleaving DNA DNA double strand breaks can be generated, with reduced off-target cleavage, by targeting two sites with EnGen Cas9 Nickase in close proximity Compatible with the EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i> (NEB #E3322) 	70/400 pmol	\$60/\$150
EnGen Spy dCas9 (SNAP-tag)	M0652S/T	<ul style="list-style-type: none"> An inactive mutant of Cas9 nuclease that retains programmable DNA binding activity The N-terminal SNAP-tag allows for covalent attachment of fluorophores, biotin, and a number of other conjugates useful for visualization and target enrichment Compatible with the EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i> (NEB #E3322) 	70/400 pmol	\$60/\$150
EnGen Lba Cas12a (Cpf1)	M0653S/T	<ul style="list-style-type: none"> 5'-TTTN-3' PAM sequence opens up additional genomic regions for targeting Shorter, 40-44 base guide RNA Two nuclear localization signals for improved transport to the nucleus 5' overhanging termini on cleavage products Active from 16 to 48°C <i>Lachnospiraceae</i> bacterium ND2006 (Lba) Cas12a maintains activity at lower temperatures than the <i>Acidaminococcus</i> orthologs, permitting editing in ectothermic organisms such as zebra fish and <i>xenopus</i> High concentration enzyme can be used for microinjection, electroporation and lipofection 	70/2,000 pmol	\$70/\$250
EnGen Sau Cas9	M0654S/T	<ul style="list-style-type: none"> 5'-NNGRRT-3' PAM Dual NLS for improved transport to nucleus High concentration enzyme can be used for microinjection electroporation and lipofection Cleaves 3 bases upstream of PAM, blunt-ended cleavage 	70/400 pmol	\$60/\$150
NEW Tth Argonaute (TtAgo)	M0665S	<ul style="list-style-type: none"> Short 16-18 oligonucleotide 5' phosphorylated ssDNA guides are cost effective and can be phosphorylated with T4 Polynucleotide Kinase Guide target RNA sequence selection is not limited by the requirement of adjacent sequence motif Highly active on ssDNA and most dsDNA substrates (generates a nick in dsDNA substrates), with mild activity on ssRNA substrates Recommended for <i>in vitro</i> applications 	50 pmol	\$75

EnGen Mutation Detection Kit

#E3321S 25 reactions \$208

Companion Products:

Gel Loading Dye, Purple (6X), no SDS
#B7025S 4 ml \$45

Quick-Load Purple 1 kb Plus DNA Ladder
#N0550S 250 gel lanes \$75
#N0550L 750 gel lanes \$187

Q5 Hot Start High-Fidelity 2X Master Mix
#M0494S 100 reactions \$217
#M0494L 500 reactions \$862
#M0494X 500 reactions \$862

Monarch PCR & DNA Cleanup Kit (5 µg)
#T1030S 50 preps \$104
#T1030L 250 preps \$475

EnGen Spy Cas9 NLS
#M0646T 400 pmol \$150
#M0646M 2,000 pmol \$600

EnGen Spy Cas9 Nickase
#M0650S 70 pmol \$60
for high (100X) concentration
#M0650T 400 pmol \$150

EnGen Spy dCas9 (SNAP-tag)
#M0652S 70 pmol \$60
for high (20X) concentration
#M0652T 400 pmol \$150

Cas9 Nuclease, *S.pyogenes*
#M0386S 70 pmol \$50
for high (20X) concentration
#M0386T 400 pmol \$135
for high (20X) concentration
#M0386M 2,000 pmol \$540

EnGen Lba Cas12a (Cpf1)
#M0653S 70 pmol \$70
for high (100X) concentration
#M0653T 2,000 pmol \$250

EnGen Sau Cas9
#M0654S 70 pmol \$60
for high (20X) concentration
#M0654T 400 pmol \$150

T7 Endonuclease I
#M0302S 250 reactions \$71
#M0302L 1,250 reactions \$284

■ T7 Endonuclease-based detection of genome editing events


Description: The EnGen Mutation Detection Kit provides reagents for detection of on-target genome editing events. In the first step, regions targeted for genome editing (i.e. CRISPR/Cas9, TALENs, Zinc-finger Nucleases) are amplified using Q5 Hot Start High-Fidelity 2X Master Mix. Upon denaturation and re-annealing, heteroduplexes are formed when mutations from insertions and deletions (indels) are present in the amplicon pool. In the second step, annealed PCR products are digested with EnGen T7 Endonuclease I, a structure-specific enzyme that will recognize mismatches larger than 1 base. Both strands of the DNA are cut when a mismatch is present, which results in the formation of smaller fragments. Analysis of the resulting fragments provides an estimate of the efficiency of the genome editing experiments.

The EnGen Mutation Detection Kit includes a Control Template and Primer Mix that can be used as a control for the PCR reaction and T7 Endonuclease I digestion. The Control Template and Primer Mix provided contains two plasmids and primers that when amplified, denatured and re-annealed will form heteroduplexes that contain a 10-base insertion. This structure is a substrate for T7 Endonuclease I. The digestion of the 600 bp heteroduplex containing amplicon yields products of 200 bp and 400 bp. 600 bp parental homoduplexes are uncleaved, and are easily distinguished from cleaved heteroduplexes when separated and visualized by agarose gel electrophoresis or fragment analysis instrument.

The protocol has been optimized so that PCR products generated by the Q5 Hot Start High-Fidelity 2X Master Mix can be introduced directly into the T7 Endonuclease I digestion without the need for purification. Digestion of the heteroduplex is complete in only 15 minutes, and Proteinase K is included to stop the reaction efficiently. Additional Q5 Hot Start High-Fidelity 2X Master Mix is also included to allow for optimization of target site amplification before digestion.

The EnGen Mutation Kit Includes:

- Q5 Hot Start High-Fidelity 2X Master Mix
- NEBuffer 2
- EnGen T7 Endonuclease I
- Control Template and Primer Mix
- Proteinase K, Molecular Biology Grade
- Quick-Load® Purple 1 kb Plus DNA Ladder
- Gel Loading Dye, Purple (6X), no SDS



Genomic DNA is amplified with primers bracketing the modified locus. PCR products are then denatured and re-annealed yielding three classes of possible structures. Duplexes containing a mismatch greater than one base are digested by T7 Endonuclease I. The DNA is then electrophoretically separated and fragment analysis is used to estimate targeting efficiency.

EnGen sgRNA Synthesis Kit, *S. pyogenes*

EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i>		
#E3322V	10 reactions	\$215
#E3322S	20 reactions	\$416

Companion Products:

EnGen Spy Cas9 NLS		
#M0646T	400 pmol	\$150
#M0646M	2,000 pmol	\$600
EnGen Spy Cas9 Nickase		
#M0650S	70 pmol	\$60
for high (20X) concentration		
#M0650T	400 pmol	\$150
EnGen Spy dCas9 (SNAP-tag)		
#M0652S	70 pmol	\$60
for high (20X) concentration		
#M0652T	400 pmol	\$150
Cas9 Nuclease, <i>S. pyogenes</i>		
#M0386S	70 pmol	\$50
for high (20X) concentration		
#M0386T	400 pmol	\$135
for high (20X) concentration		
#M0386M	2,000 pmol	\$540
EnGen Mutation Detection Kit		
#E3321S	25 rxns	\$208
RNA Loading Dye, (2X)		
#B0363S	4 ml	\$51
Monarch Total RNA Miniprep Kit		
#T2010S	50 preps	\$260
Monarch RNA Cleanup Kit (50 µg)		
#T2040S	10 preps	\$54
#T2040L	100 preps	\$290
DNase I (RNase-free)		
#M0303S	1,000 units	\$70
#M0303L	5,000 units	\$282

- Rapid generation of microgram quantities of sgRNAs in less than one hour



Configure target-specific DNA oligos design for use with the EnGen sgRNA Synthesis Kit, *S. pyogenes* with our oligo design tool accessible at sgrna.neb.com


Description: The EnGen sgRNA Synthesis Kit, *S. pyogenes* provides a simple and quick method for transcribing high yields of sgRNA in a single 30 minute reaction, using the supplied reagents and target-specific DNA oligos designed by the user.

In nature, *S. pyogenes* Cas9 is programmed with two separate RNAs, the crRNA and tracrRNA. The crRNA, or CRISPR RNA sequence contains approximately 20 nucleotides of homology complementary to the strand of DNA opposite and upstream of a PAM (Protospacer Adjacent Motif) (NGG) sequence. The tracrRNA, or transactivating crRNA, contains partial complementary sequence to the crRNA as well as the sequence and secondary structure that is recognized by Cas9. These sequences have been adapted for use in the lab by combining the tracrRNA and crRNA into one long single guide RNA (sgRNA) species capable of complexing with Cas9 to recognize and cleave the target DNA.

The EnGen sgRNA Synthesis Kit, *S. pyogenes* combines an *S. pyogenes* Cas9-specific Scaffold Oligo (included in the EnGen 2X sgRNA Reaction Mix) that is partially complementary to the target-specific oligos designed by the user. The two oligos anneal at the overlapping region and are filled in by the DNA polymerase, creating a double-stranded DNA (dsDNA) template for transcription. Synthesis of the dsDNA template and transcription of RNA occur in a single reaction, resulting in the generation of a functional sgRNA.

The EnGen sgRNA Synthesis Kit, *S. pyogenes* Includes:

- EnGen sgRNA Enzyme Mix
- EnGen 2X sgRNA Reaction Mix, *S. pyogenes*
- DNase I (RNase-free)
- EnGen sgRNA Control Oligo, *S. pyogenes*
- DTT (0.1 M)



A. The target-specific oligo contains the T7 promoter sequence, ~20 nucleotides of target-specific sequence and a 14 nucleotide overlap sequence complementary to the *S. pyogenes* Cas9 Scaffold Oligo supplied in the reaction mix. Target-specific oligos (or EnGen sgRNA Control Oligo, *S. pyogenes*) are mixed with the EnGen 2X sgRNA Reaction Mix (NTPs, dNTPs, *S. pyogenes* Cas9 Scaffold Oligo), 0.1 M DTT and the EnGen sgRNA Enzyme Mix (DNA and RNA polymerases) at room temperature. **B.** At 37°C the two oligos anneal at the 14 nucleotide overlap region of complementarity. **C.** The DNA polymerase extends both oligos from their 3' ends creating a double-stranded DNA template. **D.** The RNA polymerase recognizes the double-stranded DNA of the T7 promoter and initiates transcription. The resulting sgRNA contains the target-specific/crRNA sequence as well as the tracrRNA. All steps occur in a single reaction during a 30 minute incubation at 37°C.



Recognizing Those Making a Difference

In 1974, Don Comb established New England Biolabs with the goal of providing research tools for the life science community. His devotion to the advancement of science, stewardship of the environment and altruistic philanthropy has been a priority since the beginning, and is reflected in NEB's three core values: passion, humility and being genuine. For over 45 years, these values have served as our inspiration, and they continue to guide us in both our science and our business practices.

In 2014, NEB was coming upon its 40th anniversary, and had been considering ways to recognize this milestone. It didn't feel right to celebrate ourselves — rather, we felt it was important to celebrate our customers who shared the values we so deeply embrace. This was the inspiration for our Passion in Science Awards®. We introduced these awards as a way to recognize members of the scientific community who are committed to making a difference through their science, service, stewardship or artistic endeavors. Based on NEB's foundational values, the Passion in Science Awards were offered in four categories: Scientific Mentorship and Advocacy, Humanitarian Duty, Environmental Stewardship, and Arts and Creativity.

The Passion in Science Awards were first hosted in 2014, and then again in 2016 and 2019 — with over 40 exceptional scientists being recognized to date. Award winners visited NEB's campus in Ipswich, MA for two days of engaging discussions, inspiring presentations, tours of NEB's facilities and an awards celebration. No one anticipated how memorable these events would be.

In 2014, awardee Paul McDonald shared insight into Army policy and practices that were resulting in high-risk behavior, including suicide. Also that year, Shelly Xie told stories of people affected by neglected tropical diseases through a touching sand art presentation. In 2016, Chris Martine described how he is sharing his passion for botany with the next generation of scientists, and Lisa Anderson taught us how laboratory gloves can be upcycled into new materials. In 2019, we learned about new ways science can be shared from several of our winners — Sarah Fankhauser established a peer-reviewed publication for middle- and high-school students, Sarah McAnulty founded Skype a Scientist, Garfield Kwan introduced us to SquidToons, scientifically-accurate comics designed to engage a broader audience, and Will Murta told the story of the discovery of DNA through a musical. And the list goes on....

These events have been most inspiring and humbling for the NEB community — hearing the stories of these “unsung heroes” have taught us how we can find ways, big and small, to give back to our communities and the world. The contributions of each winner aligned with the values that Don has inspired in us all, and helped to create a true celebration of the principles that guide New England Biolabs.

To learn more about our Passion in Science Winners, visit www.neb.com/PassionInScience.

Science and art come together in this sculpture displayed outside the NEB laboratory facility. Sculpture by Dale Rogers. Credit: Kari Goodwin, NEB.

Learn more about the Passion in Science Awards.



RNA Reagents



A broad portfolio of reagents to support RNA research.

RNA molecules play a multitude of cellular roles in all kingdoms of life, perhaps reflecting the hypothetical, prehistoric RNA world. RNA is an essential carrier of genetic information, can scaffold molecular interactions, catalyze chemical reactions, and influences gene expression.

In the last several years, our understanding of RNAs as regulatory molecules in the cell has dramatically changed. Many new classes of small and large non-coding RNAs, with largely unexplored functions, have been reported, ushering in a renaissance of RNA-focused research in biology.

Small RNAs play a major role in the post-transcriptional regulation of gene expression in eukaryotic and prokaryotic organisms. Small RNAs play central roles in CRISPR pathways of adaptive immunity. Researchers have capitalized on these pathways to enable analyses of gene function not previously possible.



Introducing RNA, and in particular mRNA, into cells is an exciting new platform for inducing phenotypic changes in cells with implications for therapy, vaccines, and research applications.

New England Biolabs continues its strong tradition of providing high quality reagents to support RNA research. Our expanding range of products includes tools for the synthesis, processing, cleanup, isolation, analysis, amplification, copying and cloning of RNA molecules. Further, all NEB products pass stringent quality control assays to ensure the highest level of functionality and purity.

Featured Products

- 200** HiScribe™ T7 ARCA mRNA Kit (with tailing)
- 199** HiScribe T7 High Yield RNA Synthesis Kits
- 208** Template Switching RT Enzyme Mix
- 205** mRNA Decapping Enzyme

Featured Tools & Resources

- 198** Avoiding Ribonuclease Contamination
- 210** RNA Ligase Activity Chart
-  View our video tutorial describing high yield *in vitro* synthesis of both capped and uncapped mRNA.
-  View our video for avoiding ribonuclease contamination.



Learn more
about RNA
modifications.



Avoiding Ribonuclease Contamination	198
-------------------------------------	-----

RNA Synthesis

HiScribe T7 Quick High Yield RNA Synthesis Kit	199
HiScribe T7 High Yield RNA Synthesis Kit	199
HiScribe SP6 RNA Synthesis Kit	199
HiScribe T7 ARCA mRNA Kit	200
HiScribe T7 ARCA mRNA Kit (with tailing)	200
EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i>	193
Recommended HiScribe Kits by Application	201
T3 RNA Polymerase	202
T7 RNA Polymerase (Std. & High Conc.)	202
SP6 RNA Polymerase	202
Hi-T7 RNA Polymerase (Std. & High Conc.)	202
<i>E. coli</i> Poly(A) Polymerase	202
Poly(U) Polymerase	202
<i>E. coli</i> RNA Polymerase, (Core and Holoenzyme)	203
Pyrophosphatase, Inorganic (<i>E. coli</i> , yeast)	203
Thermostable Inorganic Pyrophosphatase	203
NudC Pyrophosphatase	203
Ribonucleotide Solution Set	203
Ribonucleotide Solution Mix	203
Sce Pst1	203
Vaccinia Capping System	204
mRNA Cap 2'-O-Methyltransferase	204
RNA Cap Analog Selection Chart	204
3'-Desthiobiotin-GTP	205
3'-Biotin-GTP	205
yDcpS	205
mRNA Decapping Enzyme	205

cDNA Synthesis, RT-PCR & RT-qPCR

cDNA Synthesis Selection Chart	206
ProtoScript II Reverse Transcriptase	206
LunaScript RT SuperMix Kit	73
LunaScript RT SuperMix Kit (Primer-free)	73
M-MuLV Reverse Transcriptase	207
AMV Reverse Transcriptase	207
WarmStart RTx Reverse Transcriptase	207
Template Switching RT Enzyme Mix	208
Primers for cDNA Synthesis	208
ProtoScript II First Strand cDNA Synthesis Kit	209
ProtoScript First Strand cDNA Synthesis Kit	209
Luna Probe One-Step RT-qPCR Kit (No ROX)	74
Luna Probe One-Step RT-qPCR 4X Mix with UDG	74
Luna Cell Ready One-Step RT-qPCR Kit	75
Luna Cell Ready Probe One-Step RT-qPCR Kit	75
Luna Universal One-Step RT-qPCR Kit	74
Luna Universal Probe One-Step RT-qPCR Kit	74
OneTaq One-Step RT-PCR Kit	81
OneTaq RT-PCR Kit	81
Luna SARS-CoV-2 RT-qPCR Multiplex Assay Kit	76

RNA Ligases & Modifying Enzymes

RNA Ligase Activity Chart	210
RNA Ligase Selection Chart	211
T4 RNA Ligase 1	212
T4 RNA Ligase 2	212
T4 RNA Ligase 2, truncated	212
T4 RNA Ligase 2, truncated K227Q and KQ	213

RtcB Ligase	213
Thermostable 5' AppDNA/RNA Ligase	213
5' DNA Adenylation Kit	214
SplintR Ligase	214
RNA 5' Pyrophosphohydrolase (RppH)	214
5' Deadenylase	214
RNase I ₁	215
RNase H	215
Thermostable RNase H	215
RNase HII	215
Quick CIP	105
Antarctic Phosphatase	106
Shrimp Alkaline Phosphatase (rSAP)	106
T4 Polynucleotide Kinase	104
ShortCut RNase III	216
XRN-1	216

RNase Control

Exonuclease T	217
Nucleoside Digestion Mix	217
DNase I (RNase-free)	217
RNase Inhibitor, Murine	217
RNase Inhibitor, Human Placenta	218
Ribonucleoside Vanadyl Complex	218

NEBNext Reagents for RNA Library Preparation


NEBNext Ultra II Directional RNA Library Prep Kit for Illumina	153
NEBNext Ultra II Directional RNA Library Prep with Sample Purification Beads	153
NEBNext Ultra II RNA Library Prep Kit for Illumina	153
NEBNext Ultra II RNA Library Prep with Sample Purification Beads	153
NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina	158
NEBNext Single Cell/Low Input cDNA Synthesis & Amplification Module	158

RNA Cleanup & Isolation

Magnetic mRNA Isolation Kit	219
Oligo d(T) ₂₅ Magnetic Beads	219
Streptavidin Magnetic Beads	219
EpiMark N6-Methyladenosine Enrichment Kit	220
Monarch Total RNA Miniprep Kit	137
Monarch RNA Cleanup Kit (10 µg)	138
Monarch RNA Cleanup Kit (50 µg)	138
Monarch RNA Cleanup Kit (500 µg)	138

RNA Markers & Ladders

dsRNA Ladder	221
microRNA Marker	221
ssRNA Ladder	221
Low Range ssRNA Ladder	221
RNA Loading Dye (2X)	221
Universal miRNA Cloning Linker	221

 Recombinant Enzyme

Avoiding Ribonuclease Contamination

Maintaining the integrity of RNA is a critical aspect of nearly all applications that use it. RNA is more susceptible to degradation than DNA, due to the ability of the 2' hydroxyl groups to act as nucleophiles. Many ribonucleases (RNases) bypass the need for metal ions by taking advantage of the 2' hydroxyl group as a reactive species. These RNases are resistant to metal chelating agents and, some of them, like RNase A family enzymes, can survive prolonged boiling or autoclaving. RNase A-type enzymes rely on active site histidine residues for catalytic activity (1) and can be inactivated by the histidine-specific alkylating agent diethyl pyrocarbonate (DEPC).

Sources of RNase Contamination:

RNases are found in all cell types and organisms from prokaryotes to eukaryotes. These enzymes generally have very high specific activity, meaning miniscule amounts of contamination in an RNA sample is sufficient to destroy the RNA. The major sources of RNase contamination in a typical laboratory include:

- Aqueous solutions, reagents used in experiments
- Exposure to RNase from environmental sources (lab surfaces, aerosols from pipetting, ungloved hands, etc.)
- Contaminated reagents

Laboratory Precautions (2,3):

New England Biolabs' enzymes certified RNase-free work have been purified to be free of ribonucleases. However, it is possible to reintroduce RNases during the course of experimentation from various sources. RNase contamination can be prevented by following a few common sense laboratory procedures:

- Always wear gloves during an experiment and change them often, especially after contact with skin, hair or other potentially RNase-contaminated surfaces, such as doorknobs, keyboards and animals.
- Use RNase-free solutions. Use RNase-free certified, disposable plasticware and filter tips whenever possible.
- Maintain a separate area for RNA work. Carefully clean the surfaces.
- Decontaminate glassware by baking at 180°C or higher for several hours or by soaking in freshly prepared 0.1% (v/v) DEPC in water or ethanol for 1 hour, followed by draining and autoclaving. Autoclaving will destroy any unreacted DEPC which can otherwise react with other proteins and RNA.
- Decontaminate polycarbonate or polystyrene materials (e.g., electrophoresis tanks) by soaking in 3% hydrogen peroxide for 10 minutes. Remove peroxide by extensively rinsing with RNase-free water prior to use.

Preparation of Solutions (2,3):

Preparation of solutions using the following suggestions can help prevent RNase contamination:

- As an alternative to the historic use of DEPC, which can inhibit enzymatic reactions if not completely removed, we have found that Milli-Q® (Millipore) purified water is sufficiently free of RNases for most RNA work. NEB also offers Nuclease-free Water (NEB #B1500).
- DEPC treatment of solutions is accomplished by adding 1 ml DEPC (Sigma) per liter of solution, stirring for 1 hour, and autoclaving for one hour to remove any remaining DEPC. [Note: Compounds with primary amine groups (e.g., Tris) which will react with DEPC, cannot be DEPC-treated. Other compounds, which are not stable during autoclaving, cannot be DEPC-treated].
- Solutions and buffers (e.g. DTT, nucleotides, manganese salts) should be prepared by dissolving the solid (highest available purity) in autoclaved DEPC-treated or Milli-Q water and passing the solution through a 0.22 µm filter to sterilize.

Inhibitors of Ribonucleases:

RNA can also be protected from RNase activity by using one of the following RNase inhibitors:

- RNase Inhibitor, Murine, (NEB #M0314) is a recombinant protein RNase inhibitor of murine origin. Like human and porcine RNase, it is a specific inhibitor of RNases A, B and C, but is more stable due to improved resistance to oxidation (4). The inhibitor requires low concentration of DTT (< 1 mM) to maintain activity, making it ideal for reactions where low DTT concentration is required (e.g., real-time RT-PCR).
- RNase Inhibitor, Human Placenta, (NEB #M0307), a recombinant protein of human placental origin, is a specific inhibitor for RNases A, B and C. Similar to the Murine RNase Inhibitor, it is compatible with many enzymatic reactions involving RNA (e.g., *in vitro* transcription, RT-PCR, ligation, etc.).
- Ribonucleoside Vanadyl Complex (NEB #S1402) is a transition-state analog inhibitor of RNase A-type enzymes with $K_i = 1 \times 10^{-5}$ M. This complex is compatible with many RNA isolation procedures, but it should not be used in the presence of EDTA. The complex also inhibits many other enzymes used in RNA work (5).

References:

- (1) Fersht, A.R. (1977) *Enzyme Structure and Mechanism* Freeman, Reading, PA, 325–329.
- (2) Blumberg, D.D. (1987) *Methods Enzymol.*, 152, 20–24.
- (3) Sambrook, J., Fritsch, E.F. and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual*, (2nd ed.), (pp. 7.3–7.5). Cold Spring Harbor: Cold Spring Harbor Laboratory Press.
- (4) Kim, B.M. et al. (1999) *Protein Science*, 8, 430–434.
- (5) Berger, S.L. (1987) *Methods Enzymol.*, 152, 227–234.

MILLI-Q® is a registered trademark of Millipore, Inc.



Find tips for
avoiding RNA
contamination.

HiScribe™ T7 High Yield RNA Synthesis Kits

HiScribe T7 High Yield RNA Synthesis Kit
#E2040S 50 reactions\$232

HiScribe T7 Quick High Yield RNA Synthesis Kit
#E2050S 50 reactions\$284

- Synthesis of long and short RNA transcripts
- Incorporation of modified nucleotides
- Incorporation of labeled nucleotides
- Generation of capped RNA using cap analogs
- Synthesis of radioactively labeled probes with high or low specific activity

Description: NEB's HiScribe T7 High Yield RNA Synthesis Kits offer robust *in vitro* RNA transcription of many kinds of RNA, including internally labeled and co-transcriptionally capped transcripts. Utilizing T7 RNA Polymerase, these kits achieve efficient transcription with small amounts of template, and can generate up to 180 µg per reaction, or up to 30–40 µg of capped RNA using cap analog. RNA generated can be used in a variety of applications, including RNA structure/function studies, ribozyme biochemistry, probes for RNase protection assays and hybridization-based blots, anti-sense RNA and RNAi experiments, microarray analysis, microinjection, and *in vitro* translation and RNA vaccines.

The HiScribe T7 High Yield RNA Synthesis Kit is an extremely flexible system, with separate NTPs included for flexible reaction setup. The HiScribe T7 Quick High Yield RNA Synthesis Kit is provided in master mix format for fast reaction setup. Each kit can yield up to 9 mg of RNA.



The HiScribe T7 High Yield RNA Synthesis Kit Includes:

- T7 RNA Polymerase Mix
- 10X T7 Reaction Buffer
- ATP, GTP, UTP, CTP (100 mM)
- FLuc Control Template

The HiScribe T7 Quick High Yield RNA Synthesis Kit Includes:

- T7 RNA Polymerase Mix
- NTP Buffer Mix
- FLuc Control Template
- DNase I (RNase-free)
- LiCl Solution

HiScribe SP6 RNA Synthesis Kit

#E2070S 50 reactions\$284

- Synthesis of long and short RNA transcripts
- Incorporation of modified nucleotides
- Incorporation of labeled nucleotides
- Generation of capped RNA using cap analogs
- Synthesis of radioactively labeled probes with high or low specific activity

Description: The HiScribe SP6 RNA Synthesis Kit is designed for the *in vitro* transcription of RNA using SP6 RNA Polymerase. This kit is suitable for synthesis of high yield RNA transcripts and allows for incorporation of cap analogs (not included) or modified nucleotides (not included) to obtain capped, biotin-labeled or dye-labeled RNA. The kit is also capable of synthesizing high specific activity radiolabeled RNA for use as probes or targets.

RNA synthesized from this kit is suitable for many applications including RNA structure and function studies, ribozyme biochemistry, probes for RNase protection or gel shift assays, hybridization-based blots, anti-sense RNA or RNAi experiments, microarray analysis, microinjection and *in vitro* translation studies.

This kit contains sufficient reagents for 50 reactions of 25 µl each. Each standard reaction yields ≥ 80 µg of RNA from 1 µg SP6 Control Template DNA. Each kit can yield ≥ 4 mg of RNA.



The HiScribe SP6 RNA Synthesis Kit Includes:

- SP6 Reaction Buffer (10X)
- ATP (Tris), GTP (Tris), UTP (Tris), CTP (Tris), (50 mM)
- SP6 Control Template
- SP6 RNA Polymerase Mix
- DNase I (RNase-free)
- LiCl Solution

Use Monarch RNA Cleanup Kits to purify your synthesized RNA, see page 137–138.

Companion Products:

DNase I (RNase-Free)
#M0303S 1,000 units\$70
#M0303L 5,000 units\$282

RNA Loading Dye (2X)
#B0363S 4 ml\$51

Vaccinia Capping System
#M2080S 400 units\$140

Q5® Hot Start High-Fidelity DNA Polymerase
#M0493S 100 units\$140
#M0493L 500 units\$554

Monarch® PCR & DNA Cleanup Kit
#T1030S 50 preps\$104
#T1030L 250 preps\$475

Monarch RNA Cleanup Kit (10 µg)
#T2030S 10 preps\$54
#T2030L 100 preps\$290

Monarch RNA Cleanup Kit (50 µg)
#T2040S 10 preps\$54
#T2040L 100 preps\$290

Monarch RNA Cleanup Kit (500 µg)
#T2050S 10 preps\$60
#T2050L 100 preps\$480

Low Range ssRNA Ladder
#N0364S 25 gel lanes\$71

E. coli Poly(A) Polymerase
#M0276S 100 units\$72
#M0276L 500 units\$292

mRNA Cap 2'-O-Methyltransferase
#M0366S 2,000 units\$58

3'-O-Me-m⁷G(5')ppp(5')G RNA Cap Structure Analog
#S1411S 1 µmol\$155
#S1411L 5 µmol\$617

G(5')ppp(5')A RNA Cap Structure Analog
#S1406S 1 µmol\$144
#S1406L 5 µmol\$574

G(5')ppp(5')G RNA Cap Structure Analog
#S1407S 1 µmol\$124
#S1407L 5 µmol\$497

m⁷G(5')ppp(5')G RNA Cap Structure Analog
#S1404S 1 µmol\$161
#S1404L 5 µmol\$641

m⁷G(5')ppp(5')A RNA Cap Structure Analog
#S1405S 1 µmol\$161
#S1405L 5 µmol\$641

HiScribe T7 ARCA mRNA Kits

HiScribe T7 ARCA mRNA Kit
#E2065S 20 reactions\$328

HiScribe T7 ARCA mRNA Kit (with tailing)
#E2060S 20 reactions\$384

Companion Products:

DNase I (RNase-Free)
#M0303S 1,000 units\$70
#M0303L 5,000 units\$282
RNA Loading Dye (2X)
#B0363S 4 ml\$51

- Synthesis of capped and tailed mRNA
- Incorporation of modified nucleotides
- Template removal and mRNA purification reagents included

Description: Most eukaryotic mRNAs require a 7-methyl guanosine (m7G) cap structure at the 5' end and a Poly(A) tail at the 3' end to be efficiently translated. By using a DNA template encoding a poly(A) tail, the HiScribe T7 ARCA mRNA Kit can be used to synthesize capped and tailed mRNAs. The cap structure is added to the mRNA by co-transcriptional incorporation of Anti-Reverse Cap Analog (ARCA, NEB #S1411) using T7 RNA Polymerase. The transcription reaction can be set up easily by combining the ARCA/NTP Mix, T7 RNA Polymerase Mix and a suitable DNA template. The kit also allows for partial incorporation of 5mCTP, Pseudo-UTP and other modified nucleotides into mRNA. mRNAs synthesized with the kit can be used for cell transfection, microinjection, *in vitro* translation and RNA vaccines. Poly(A) tail is incorporated during the transcription reaction. The kit also includes DNase I and LiCl for DNA template removal and quick mRNA purification.

The HiScribe T7 ARCA mRNA Kit (with tailing) is designed for quick production of ARCA capped and poly(A) tailed mRNA *in vitro* from templates without encoded poly(A) tails.



The HiScribe T7 ARCA mRNA Kit Includes:

- T7 RNA Polymerase Mix
- ARCA/NTP Mix
- DNase I (RNase-free)
- LiCl Solution
- CLuc Control Template

The HiScribe T7 ARCA mRNA Kit (with tailing) Includes:

- T7 RNA Polymerase Mix
- ARCA/NTP Mix
- DNase I (RNase-free)
- *E. coli* Poly(A) Polymerase
- Poly(A) Polymerase Reaction Buffer
- LiCl Solution
- CLuc Control Template

Advantages:

- Quick reaction setup and streamlined protocol
- Enables partial incorporation of 5mCTP, Pseudo-UTP and other modified CTP and UTP
- Ultra-high quality components ensure mRNA integrity

EnGen® sgRNA Synthesis Kit, *S. pyogenes*

See page 179 for more information.

#E3322V 10 reactions\$215
#E3322S 20 reactions\$416

The EnGen sgRNA Synthesis Kit, *S. pyogenes* provides a simple and quick method for transcribing high yields of sgRNA in a single 30 minute reaction, using the supplied reagents and target-specific DNA oligos designed by the user.

Bijoyita has been with NEB for four years. In her role as Senior Scientist, Bijoyita leads a research group focused on synthetic mRNA synthesis and modification. In her spare time, Bijoyita likes to travel and try new cuisines.



Recommended HiScribe RNA Synthesis Kits by Application

	APPLICATION	T7 KITS				SP6 KITS
		HiScribe T7 High Yield RNA Synthesis Kit (#E2040)	HiScribe T7 Quick High Yield RNA Synthesis Kit (#E2050)	HiScribe T7 ARCA mRNA Kit (#E2065)	HiScribe T7 ARCA mRNA (with Tailing) (#E2060)	HiScribe SP6 RNA Synthesis Kit (#E2070)
Probe labeling	Fluorescent labeling: FAM, Cyanine (Cy) dyes, etc. • Fluorescent <i>in situ</i> hybridization (FISH)		✓			✓
	Non-fluorescent labeling: Biotin, Digoxigenin • <i>In situ</i> hybridization • Blot hybridization with secondary detection • Microarray		✓			✓
	High specific activity radiolabeling • Blot hybridization • RNase protection	✓				✓
mRNA & RNA for transfection	Streamlined mRNA synthesis with ARCA co-transcriptional capping and enzymatic poly(A) tailing • Transfection • Microinjection • <i>In vitro</i> translation				✓	
	Streamlined ARCA capped RNA synthesis • Template encoded poly(A) tails • Non polyadenylated transcripts • Transfection • Microinjection • <i>In vitro</i> translation			✓		
	Co-transcriptional capping with alternate cap analogs • Transfection • Microinjection • <i>In vitro</i> translation		✓			✓
	Post-transcriptional capping with Vaccinia Capping System • Transfection • Microinjection • <i>In vitro</i> translation	✓	✓			✓
	Complete substitution of NTPs: 5-mC, pseudouridine, etc. • Post-transcriptional capping with Vaccinia mRNA Capping System	✓				✓
	Partial substitution of NTPs: 5-mC, pseudouridine, etc.		✓	✓	✓	✓
	Unmodified RNA		✓			✓
	Hairpins, short RNA, dsRNA • Gene knockdown		✓			✓
Structure, function, & binding studies	Complete substitution of NTPs • Aptamer selection • Isotopic labeling	✓				✓
	Partial substitution of one or more NTPs • Aptamer selection • Structure determination		✓			✓
	Unmodified RNA • SELEX • Structure determination		✓			✓

RNA Polymerases

T3 RNA Polymerase

#M0378S 5,000 units \$71

T7 RNA Polymerase

#M0251S 5,000 units \$70

#M0251L 25,000 units \$282

SP6 RNA Polymerase

#M0207S 2,000 units \$68

- Radiolabeled RNA probe preparation
- RNA generation for *in vitro* translation
- RNA generation for studies of RNA structure, processing and catalysis

NEW

T7 RNA Polymerase (High Concentration)

for high (20X) concentration

#M0460T 50,000 units \$564

Hi-T7 RNA Polymerase

#M0658S 5,000 units \$85

NEW

Hi-T7 RNA Polymerase (High Concentration)

for high (20X) concentration

#M0470T 50,000 units \$680

Description: Initiation of transcription with T3, T7 and SP6 RNA Polymerases is highly specific for the T3, T7 and SP6 phage promoters, respectively. Cloning vectors have been developed which direct transcription from the T3, T7 or SP6 promoter through polylinker cloning sites. These vectors allow *in vitro* synthesis of defined RNA transcripts from a cloned DNA sequence. T7 RNA Polymerase (High Concentration) is offered at a 20-fold higher concentration than our standard T7 RNA Polymerase and is ideal for experienced users interested in building and optimizing their own *in vitro* transcription reactions.

Hi-T7 RNA Polymerase is an engineered, thermoactive T7 RNA Polymerase. Hi-T7 uses T7 RNA Polymerase



Promoters. It can increase capping efficiency and eliminate dsRNA by-product formation during synthesis. Hi-T7 RNA Polymerase (High Concentration) is offered at a 20-fold higher concentration than our standard Hi-T7 RNA Polymerase and is ideal for experienced users interested in building and optimizing their own *in vitro* transcription reactions.

Reaction Conditions: 1X RNAPol Reaction Buffer. Supplement with 0.5 mM each ATP, UTP, GTP, CTP (not included) and DNA template containing the appropriate promoter. Incubate at 37°C (T3, T7 and SP6) or 50°C (Hi-T7). Protocols involving high concentration T7 and Hi-T7 RNA Polymerases are to be designed and optimized by the user.

Unit Definition: One unit is defined as the amount of enzyme required to incorporate 1 nmol ATP into an acid-insoluble material in 1 hour at 37°C or 50°C for Hi-T7. Unit assay conditions can be found at www.neb.com.

Concentration: T3 RNA Polymerase: 50,000 units/ml. T7 RNA Polymerase: 50,000 units/ml. T7 RNA Polymerase (High Concentration): 1,000,000 units/ml. SP6 RNA Polymerase: 20,000 units/ml. Hi-T7 RNA Polymerase: 50,000 units/ml. Hi-T7 RNA Polymerase (High Concentration): 1,000,000 units/ml.

E. coli Poly(A) Polymerase

#M0276S 100 units \$72

#M0276L 500 units \$292

Companion Products:

Adenosine-5' Triphosphate (ATP)

#P0756S 1 ml \$37

#P0756L 5 ml \$143

RNase Inhibitor, Murine

#M0314S 3,000 units \$73

#M0314L 15,000 units \$292

- Labeling of RNA with ATP or cordycepin 5'-triphosphate
- Poly(A) tailing of RNA for cloning or affinity purification
- Enhances translation of RNA transferred into eukaryotic cells

Description: *E. coli* Poly(A) Polymerase catalyzes the template independent addition of AMP from ATP to the 3' end of RNA.

Reaction Conditions: 1X Poly(A) Polymerase Reaction Buffer. Supplement with 1 mM ATP. Incubate at 37°C.



Reagents Supplied with Enzyme:

10X Poly(A) Polymerase Reaction Buffer
10 mM ATP

Unit Definition: One unit is defined as the amount of enzyme that will incorporate 1 nmol of AMP into RNA in a 20 µl volume in 10 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 5,000 units/ml

Poly(U) Polymerase

#M0337S 60 units \$92

Companion Products:

RNase Inhibitor, Murine

#M0314S 3,000 units \$73

#M0314L 15,000 units \$292

Ribonucleotide Solution Set

#N0450S 10 µmol of each \$72

#N0450L 50 µmol of each \$289

- Labeling of RNA with UTP
- Poly(U) tailing of RNA for cloning

Description: Poly(U) Polymerase catalyzes the template-independent addition of UMP from UTP or AMP from ATP to the 3' end of RNA.

Reaction Conditions: 1X NEBuffer 2. Supplement with 0.5 mM UTP (not supplied). Incubate at 37°C. May heat inactivate at 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme that incorporates 1 nmol of UMP into RNA in a



50 µl volume in 10 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 2,000 units/ml

Notes: Poly(U) Polymerase in NEBuffer 2 will incorporate UMP or AMP from UTP or ATP into RNA. Tailing length of poly(U) varies with UTP and primer concentration. Poly(U) Polymerase is highly processive under low primer concentrations (< 100 pmol).

E. coli RNA Polymerase, Core Enzyme & Holoenzyme

E. coli RNA Polymerase, Core Enzyme
#M0550S 100 units\$184

E. coli RNA Polymerase, Holoenzyme
#M0551S 50 units\$139

- RNA synthesis from *E. coli* promoter
- Transcription initiation studies
- In vitro translation with PURExpress

Description: *E. coli* RNA Polymerase Core Enzyme consists of 5 subunits designated α , α' , β , β' , and ω . The enzyme is free of sigma factor and does not initiate specific transcription from bacterial and phage DNA promoters. The enzyme retains the ability to transcribe RNA from nonspecific initiation sequences. Addition of sigma factors will allow the enzyme to initiate RNA synthesis from specific bacterial and phage promoters. The core enzyme has a molecular weight of ~400 kDa.

E. coli RNA Polymerase Holoenzyme is the core enzyme saturated with sigma factor 70. The Holoenzyme initiates RNA synthesis from sigma 70 specific bacterial and phage promoters.

37°

Reaction Conditions: 1X *E. coli* RNA Polymerase Reaction Buffer, 0.5 mM of each rNTP and DNA template. Incubate at 37°C.

Unit Definition: One unit is the amount of enzyme required to incorporate 1 nmol of NTP into RNA in 10 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 1,000 units/ml

Pyrophosphatases

Pyrophosphatase, Inorganic (*E. coli*)
#M0361S 10 units\$64
#M0361L 50 units\$257

Pyrophosphatase, Inorganic (yeast)
#M2403S 10 units\$69
#M2403L 50 units\$278

Thermostable Inorganic Pyrophosphatase
#M0296S 250 units\$74
#M0296L 1,250 units\$299


NEW

NudC Pyrophosphatase
#M0607S 250 pmol\$72

- Increasing RNA yield in transcription reactions
- Enhancing DNA replication

Thermostable Inorganic Pyrophosphatase retains 100% activity after incubation at 100°C for 4 hours.

Description: Inorganic pyrophosphatase (PPase) catalyzes the hydrolysis of inorganic pyrophosphate to form orthophosphate.



Source: Pyrophosphatase, Inorganic (*E. coli*) is prepared from a clone of the *E. coli* inorganic pyrophosphatase gene.

Pyrophosphatase, Inorganic (yeast) is an *E. coli* strain containing a genetic fusion of the *Saccharomyces cerevisiae* ppa gene and the gene coding for *Mycobacterium xenopi* GyrA intein. Developed by BioHelix Corporation, now a wholly owned subsidiary of Quidel Corporation, and produced at New England Biolabs.

Thermostable Inorganic Pyrophosphatase is an *E. coli* strain carrying a plasmid encoding a pyrophosphatase from the extreme thermophile *Thermococcus litoralis*.

NudC is a NUDIX pyrophosphatase that efficiently hydrolyzes NAD⁺- and NADH-capped RNA, generating a ligatable 5' monophosphate on the RNA (NAD⁺ decapping or deNADding). Deletion of the nudC gene has been shown to increase the fraction of NAD⁺-capped RNA in *E. coli*.

RR

Unit Definition: One unit is defined as the amount of enzyme that will generate 1 μ mol of phosphate per minute from inorganic pyrophosphate under standard reaction conditions. Unit assay conditions can be found at www.neb.com.

1 μ M of NudC hydrolyzes 200 μ M or more NAD⁺ into NMN⁺ and AMP in 1X NEBuffer r3.1 and 5 mM DTT at 37°C for 30 min.

Concentration: Pyrophosphatase, Inorganic (*E. coli*) and Pyrophosphatase, Inorganic (yeast) 100 units/ml. Thermostable Inorganic Pyrophosphatase 2,000 units/ml. NudC Pyrophosphatase: 10 μ M

Note: This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

NudC Pyrophosphatase is an **Enzyme for Innovation** (EFI). EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit www.neb.com/EnzymesforInnovation to view the full list.

Ribonucleotides

Ribonucleotide Solution Set
#N0450S 10 μ mol of each\$72
#N0450L 50 μ mol of each\$289

Ribonucleotide Solution Mix
#N0466S 10 μ mol of each\$75
#N0466L 50 μ mol of each\$302

Description:

Ribonucleotide Solution Set:
The set consists of four separate solutions of ATP, GTP, CTP and UTP, pH 7.5, as sodium salts. Each nucleotide is supplied at 100 mM.

Ribonucleotide Solution Mix:
A buffered equimolar solution of ribonucleotide triphosphates: rATP, rCTP, rGTP and rUTP, pH 7.5, as sodium salts. Each nucleotide is supplied at a concentration of 25 mM (total rNTP concentration equals 100 mM).

Note: To ensure maximum activity upon long-term storage, aliquot and store at -80°C.

NEW

Sce Pus1

#M0526S 5,000 pmol\$75

- Sequence-specific pseudouridine modification is an alternative to randomly incorporated modified nucleosides by RNA polymerases

Description: Sce Pseudouridine Synthase I (Sce PUS1) converts Uridine to Pseudouridine in single-stranded RNA, with a preference for Uridines in single-stranded RNA regions over Uridines in double-stranded RNA. The optimal substrate is an unstructured RNA that is 15 nt long or longer.

Reaction Conditions: 1X NEBuffer r1.1. Incubate at 30°C.

RR

Concentration: 100 pmol/ μ l

Note: This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

Sce PUS1 is an **Enzyme for Innovation** (EFI). EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit www.neb.com/EnzymesforInnovation to view the full list.

Vaccinia Capping System

#M2080S 400 units\$140

Companion Product:

RNase Inhibitor, Murine

#M0314S 3,000 units\$73

#M0314L 15,000 units\$292

- Capping mRNA prior to *in vivo* or *in vitro* translation
- Labeling 5' end of mRNA

Description: Based on the Vaccinia Virus Capping Enzyme, the Vaccinia Capping System provides the necessary components to add 7-methylguanylate cap structures (Cap 0) to the 5' end of RNA. In eukaryotes, these terminal cap structures are involved in stabilization, transport and translation of mRNAs. Enzymatic production of capped RNA is an easy way to improve the stability and translational competence of RNA used for *in vitro* translation, transfection and microinjection. Alternatively, use of labeled GTP in a reaction provides a convenient way to label any RNA containing a 5' terminal triphosphate.

This single enzyme is composed of two subunits (D1 and D12) and has three enzymatic activities (RNA triphosphatase and guanylyltransferase by the D1 subunit and guanine methyltransferase by the D12 subunit), all

RI 37°

necessary for addition of a complete Cap 0 structure, m⁷Gppp(5')N. All capped structures are added in the proper orientation, unlike co-transcriptional addition of some cap analogs.

Reaction Conditions: 1X Capping Buffer. Supplement with 0.5 mM GTP and 0.1 mM SAM. Incubate at 37°C.

Reagents Supplied:

Vaccinia Capping Enzyme
Capping Buffer (10X)
GTP Solution (10 mM)
SAM Solution (32 mM)

Unit Definition: One unit is defined as the amount of enzyme required to incorporate 10 pmol of (α-³²P) GTP into an 80 nucleotide transcript in 1 hour at 37°C.

Concentration: 10,000 units/ml

mRNA Cap 2'-O-Methyltransferase

#M0366S 2,000 units\$58

- Enhances translation of RNA
- Improving mRNA expression during microinjection and transfection

Description: mRNA Cap 2'-O-Methyltransferase adds a methyl group at the 2'-O position of the first nucleotide adjacent to the cap structure at the 5' end of the RNA. The enzyme utilizes S-adenosylmethionine (SAM) as a methyl donor to methylate capped RNA (Cap 0) resulting in a Cap 1 structure.

Reaction Conditions: 1X Capping Buffer. Supplement with 0.2 mM SAM (supplied). Incubate at 37°C.

RI 37°

Reagents Supplied:

Capping Buffer (10X)
SAM (32 mM)

Unit Definition: One unit is defined as the amount of enzyme required to methylate 10 pmol of 80 nt long capped RNA transcript in 1 hour at 37°C.

Concentration: 50,000 units/ml

RNA Cap Analog Selection Chart

The 5' terminal m⁷G cap present on most eukaryotic mRNAs is required for translation, *in vitro*, at the initiation level. For most RNAs, the cap structure increases stability, decreases susceptibility to exonuclease degradation, and promotes the formation of mRNA initiation complexes. Certain prokaryotic mRNAs with 5' terminal cap structures are translated as efficiently as eukaryotic mRNA in a eukaryotic cell-free protein synthesizing system. Splicing of certain eukaryotic substrate RNAs has also been observed to require a cap structure.

PRODUCT	APPLICATION	SIZE	PRICE
(ARCA) Anti-Reverse Cap Analog 3'-O-Me-m⁷G(5')ppp(5')G (#S1411S/L)	<ul style="list-style-type: none"> • Produces 100% translatable capped transcripts • Co-transcriptional capping with T7 (NEB #M0251), SP6 (NEB #M0207) and T3 RNA polymerases • Synthesis of m⁷G capped RNA for <i>in vitro</i> splicing assays • Synthesis of m⁷G capped RNA for transfection or microinjection • Ensures incorporation in correct orientation 	1/5 μmol	\$155/ \$617
Standard Cap Analog m⁷G(5')ppp(5')G (#S1404S/L)	<ul style="list-style-type: none"> • Co-transcriptional capping with T7, SP6 and T3 RNA polymerases • Synthesis of m⁷G capped RNA for <i>in vitro</i> splicing assays • Synthesis of m⁷G capped RNA for transfection or microinjection 	1/5 μmol	\$161/ \$641
Unmethylated Cap Analog G(5')ppp(5')G (#S1407S/L)	<ul style="list-style-type: none"> • Co-transcriptional capping with T7, SP6 and T3 RNA polymerases • Synthesis of unmethylated G capped RNA 	1/5 μmol	\$124/ \$497
Methylated Cap Analog for A +1 sites m⁷G(5')ppp(5')A (#S1405S/L)	<ul style="list-style-type: none"> • Co-transcriptional capping with T7 RNA polymerase from the phi2.5 promoter that contains an A at the transcription initiation site • Synthesis of m⁷G capped RNA for <i>in vitro</i> splicing assays • Synthesis of m⁷G capped RNA for transfection or microinjection 	1/5 μmol	\$161/ \$641
Unmethylated Cap Analog for A +1 sites G(5')ppp(5')A (#S1406S/L)	<ul style="list-style-type: none"> • Co-transcriptional capping with T7 RNA polymerase from the phi2.5 promoter that contains an A at the transcription initiation site • Synthesis of unmethylated G capped RNA • Synthesis of A capped RNA 	1/5 μmol	\$144/ \$574

3'-Desthiobiotin-GTP & 3'-Biotin-GTP

3'-Desthiobiotin-GTP

#N0761S 0.5 µmol\$323

3'-Biotin-GTP

#N0760S 0.5 µmol\$323

Description: 3'-Desthiobiotin-GTP or 3'-Biotin-GTP are guanosine triphosphate (GTP) analogs which are modified at their 3' position with desthiobiotin or biotin, respectively. When used with the Vaccinia Capping System, (NEB #M2080) these reagents enable affinity tagging of RNA triphosphate ends. Tagged RNAs are enriched by binding to Hydrophilic Streptavidin Magnetic Beads

(NEB #S1421). Desthiobiotin-tagged RNAs can be eluted with free biotin. This approach is used in Cappable-seq, a method developed at NEB for directly enriching the 5'-ends of primary transcripts (1).

Reference:

(1) Ettwiller, L. et al. (2016) *BMC Genomics*, 17,199.

yDcpS

#M0463S 4,000 units\$261

- mRNA decapping, enabling recapping with tagged-GTP analogs
- Biotinylation of 5' ends of primary transcripts
- Recappable-seq

Description: yDcpS decapping enzyme from *S. cerevisiae* hydrolyzes the phosphodiester bond between the gamma and beta phosphates of m⁷G capped mRNA, leaving behind a diphosphorylated 5' end and m⁷GMP. yDcpS is capable of decapping full length mRNAs and the diphosphorylated 5' end it leaves behind is suitable for recapping using Vaccinia Capping Enzyme.

Reaction Conditions: 1X yDcpS Reaction Buffer. Incubate at 37°C.

Concentration: 200,000 units/ml



NEW

mRNA Decapping Enzyme

#M0608S 2,000 units\$95

- Efficient replacement for Tobacco Acid Pyrophosphatase
- Cap0 and Cap1 are removed with equal efficiency
- Suitable for 5' RLM-RACE and RNA-seq

Description: mRNA Decapping Enzyme catalyzes the removal of 7-methylguanosine cap (m⁷G) from the 5' end of mRNA, producing 5' monophosphate and releasing m⁷GDP. mRNA Decapping Enzyme is capable of decapping mRNAs of various lengths and removes both Cap0 and Cap1 structures with similar efficiency. mRNA Decapping Enzyme also converts 5' triphosphate ends to 5' monophosphate, albeit with reduced efficiency.

Reaction Conditions: 1X mRNA Decapping Enzyme Reaction Buffer. Incubate at 37°C.

Concentration: 100,000 units/ml



Sharon has been the Director of Information Technology for over 7 years, making sure that all of NEB's IT processes run smoothly, both internally and externally. When not at work, Sharon enjoys mountain hiking and creating fabric art.



cDNA Synthesis Selection Chart

cDNA SYNTHESIS	FEATURES	SIZE	PRICE
KITS AND MIXES			
LunaScript® RT SuperMix Kit (NEB #E3010)	<ul style="list-style-type: none"> Ideal for cDNA synthesis in a two-step RT-qPCR workflow Single tube supermix contains random hexamer and oligo-dT primers, dNTPs, Murine RNase Inhibitor, and Luna Reverse Transcriptase Visible blue tracking dye for easy reaction setup Fast 13-minute protocol Primer-free master mix available (NEB #E3025) 	25/100 rxns	\$141/\$437
ProtoScript® II First Strand cDNA Synthesis Kit (NEB #E6560)	<ul style="list-style-type: none"> Generates cDNA at least 10 kb in length Contains ProtoScript II Reverse Transcriptase, an enzyme with increased thermostability and reduced RNase H activity Convenient 2-tube kit includes dNTPs, Oligo-dT primer and Random Primer Mix 	30/150 rxns	\$171/\$685
ProtoScript First Strand cDNA Synthesis Kit (NEB #E6300)	<ul style="list-style-type: none"> Generates cDNA at least 5 kb in length Contains M-MuLV Reverse Transcriptase Convenient 2-tube kit includes dNTPs, Oligo-dT primer and Random Primer Mix 	30/150 rxns	\$162/\$645
Template Switching RT Enzyme Mix (NEB #M0466)	<ul style="list-style-type: none"> Incorporates a universal adaptor sequence at the 3' end of cDNA during the RT reaction Enzyme mix and buffer are optimized for efficient template switching RT enzyme mix includes RNase Inhibitor High sensitivity for cDNA amplification – enables transcriptome analysis by RNA-seq from single cells or as low as 2 pg of human total RNA Robust and simple workflow for 5' Rapid Amplification of cDNA Ends (RACE) Retains the complete 5' end of transcripts for 2nd Strand cDNA Synthesis 	20/100 rxns	\$90/\$360
STANDALONE REAGENTS			
ProtoScript II Reverse Transcriptase (NEB #M0368) An alternative to SuperScript® II	<ul style="list-style-type: none"> RNase H⁻ mutant of M-MuLV Reverse Transcriptase with increased thermostability and reduced RNase H activity Increased reaction temperatures (37–50°C) 	4,000/10,000/40,000 units	\$92/\$184/\$656
M-MuLV Reverse Transcriptase (NEB #M0253)	<ul style="list-style-type: none"> Robust reverse transcriptase for a variety of templates Standard reaction temperatures (37–45°C) 	10,000/50,000 units	\$73/\$294
AMV Reverse Transcriptase (NEB #M0277)	<ul style="list-style-type: none"> Robust reverse transcriptase for a broad temperature range (37–52°C) Can be used for templates requiring higher reaction temperatures 	200/1,000 units	\$74/\$298
WarmStart® RTx Reverse Transcriptase (NEB #M0380)	<ul style="list-style-type: none"> Permits room temperature reaction setup Increased reaction temperatures (50–65°C) Optimized for RT-LAMP isothermal detection 	50/250 rxns	\$69/\$277

For RT-PCR kits see page 81.

For RT-qPCR kits see page 74.

SUPERSRIPT® is a registered trademark of Thermo Fisher, Inc.

ProtoScript® II Reverse Transcriptase



#M0368S	4,000 units	\$92
#M0368L	10,000 units	\$184
#M0368X	40,000 units	\$656

Companion Products:

RNase H			
#M0297S	250 units	\$72
#M0297L	1,250 units	\$289
Monarch Total RNA Miniprep Kit			
#T2010S	50 preps	\$260

- Efficient reverse transcription from different starting RNA amounts
- Increased thermostability
- Generates cDNA up to 10 kb or more

Description: ProtoScript II Reverse Transcriptase is a recombinant M-MuLV reverse transcriptase with reduced RNase H activity and increased thermostability. It can be used to synthesize first strand cDNA at higher temperatures than the wild type M-MuLV. The enzyme is active up to 50°C, providing higher specificity, higher yield of cDNA and more full length cDNA product, up to 12 kb in length. This product was formerly known as M-MuLV Reverse Transcriptase RNase H⁻.

Reaction Conditions: 1X ProtoScript II Reverse Transcriptase Reaction Buffer, 10 mM DTT, 200 units M-MuLV (RNase H⁻), supplemented with 0.5 mM dNTPs (not included) and 5 μM dT₂₃VN (not included).

Incubate at 42°C for 50 minutes. If random primers are used, a 10 minute incubation at room temperature is recommended before transferring to 42°C. May heat inactivate at 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme that will incorporate 1 nmol of dTTP into acid-insoluble material in a total reaction volume of 50 μl in 10 minutes at 37°C using poly(rA)•oligo(dT)₁₈ as template. Unit assay conditions can be found at www.neb.com.

Concentration: 200,000 units/ml

LunaScript® RT SuperMix Kit & LunaScript® RT Master Mix Kit (Primer-free)



See page 73 for more information.

LunaScript RT SuperMix Kit	
#E3010S	25 reactions\$141
#E3010L	100 reactions\$437

NEW

LunaScript RT Master Mix Kit (Primer-free)	
#E3025S	25 reactions\$117
#E3025L	100 reactions\$364

LunaScript RT SuperMix Kit is an optimized master mix containing all the necessary components for first strand cDNA synthesis in the context of a two-step RT-qPCR workflow. It features the thermostable Luna Reverse Transcriptase, which supports cDNA synthesis at elevated temperatures. Murine RNase Inhibitor is also included to protect template RNA from degradation. The LunaScript RT SuperMix Kit contains random hexamer and poly-dT primers, allowing for even coverage across the length of the RNA targets.

LunaScript RT Master Mix Kit (Primer-free) is the same formulation as the LunaScript RT SuperMix Kit but does not contain primers. It includes all components needed for carrying out first strand cDNA synthesis except for RNA template and user-supplied primers. It provides flexible options for using different primers (dT primers, random primers, gene-specific primers, modified primers, etc.) for first strand cDNA synthesis.

M-MuLV Reverse Transcriptase



#M0253S	10,000 units\$73
#M0253L	50,000 units\$294

Companion Product:

Monarch Total RNA Miniprep Kit	
#T2010S	50 preps\$260

Description: Moloney Murine Leukemia Virus (M-MuLV) Reverse Transcriptase is an RNA-directed DNA polymerase. This enzyme can synthesize a complementary DNA strand initiating from a primer using either RNA (cDNA synthesis) or single-stranded DNA as a template. M-MuLV Reverse Transcriptase lacks 3' → 5' exonuclease activity.

Reaction Conditions: 1X M-MuLV Reverse Transcriptase Reaction Buffer. Supplement with dNTPs (not included). Incubate at 37–42°C. May heat inactivate at 65°C for 20 minutes.

Unit Definition: One unit is the amount of enzyme required to incorporate 1 nmol of dTTP into an acid-insoluble material in 10 minutes at 37°C using poly(rA)•oligo(dT) as template primer. Unit assay conditions can be found at www.neb.com.

Concentration: 200,000 units/ml

- cDNA synthesis
- RNA Sequencing
- RT-PCR

AMV Reverse Transcriptase



#M0277S	200 units\$74
#M0277L	1,000 units\$298

Companion Product:

Monarch Total RNA Miniprep Kit	
#T2010S	50 preps\$260

Description: Avian Myeloblastosis Virus (AMV) Reverse Transcriptase is an RNA-directed DNA polymerase. This enzyme can synthesize a complementary DNA strand initiating from a primer using RNA (cDNA synthesis) or single-stranded DNA as a template.

Source: Avian Myeloblastosis Virus (AMV)

Reaction Conditions: 1X AMV Reverse Transcriptase Reaction Buffer. Supplement with dNTPs (not included). Incubate at 37–42°C. May heat inactivate at 85°C for 5 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to incorporate 1 nmol of dTTP into an acid-insoluble form in 10 minutes at 37°C using poly(rA)•oligo(dT) as template primer. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 and 25,000 units/ml

Storage Note: Once thawed, store at –20°C. Repeated freeze thaw cycles will inactivate the enzyme. Aliquots can be stored for longer periods at –70°C.

- cDNA synthesis
- RNA Sequencing
- RT-PCR

WarmStart® RTx Reverse Transcriptase



#M0380S	50 reactions\$69
#M0380L	250 reactions\$277

Companion Product:

Monarch Total RNA Miniprep Kit	
#T2010S	50 preps\$260

Description: WarmStart RTx Reverse Transcriptase is a unique *in silico*-designed, RNA-directed DNA polymerase coupled with a reversibly-bound aptamer that inhibits RTx activity below 40°C. This enzyme can synthesize a complementary DNA strand initiating from a primer using RNA (cDNA synthesis) or single-stranded DNA as a template. RTx is a robust enzyme for RNA detection in amplification reactions and is particularly well suited for use in Loop-mediated Isothermal Amplification (LAMP). The WarmStart property enables high-throughput applications, room-temperature setup, and increases the consistency and specificity of amplification reactions. RTx contains intact RNase H activity.

Reaction Conditions: 1X Isothermal Amplification Buffer, template, primer, dNTPs and 0.25–0.5 µl of WarmStart RTx Reverse Transcriptase in a reaction volume of 25 µl. Incubate at 50–55°C for cDNA synthesis or directly at 65°C for One-step RT-LAMP. May heat inactivate at 80°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme that will incorporate 1 nmol of dTTP into acid-insoluble material in a total reaction volume of 50 µl in 20 minutes at 50°C using poly(rA)•oligo(dT)18 as template.

Concentration: 15,000 units/ml

- RT-LAMP
- cDNA Synthesis
- RT reactions requiring room temperature setup

For more information on products for LAMP, see pages 76–79.

Template Switching RT Enzyme Mix



#M0466S	20 reactions\$90
#M0466L	100 reactions\$360

Companion Products:

NEBNext High-Fidelity 2X PCR Master Mix		
#M0541S	50 reactions\$90
#M0541L	250 reactions\$360

Q5 Hot Start High-Fidelity 2X Master Mix		
#M0494S	100 rxns (2 x 1.25 ml)\$217
#M0494L	500 rxns (10 x 1.25 ml)\$862
#M0494X	500 rxns (1 x 12.5 ml)\$862

NEB PCR Cloning Kit		
#E1202S	20 reactions\$289

Monarch Total RNA Miniprep Kit		
#T2010S	50 preps\$260

LongAmp Hot Start <i>Taq</i> 2X Master Mix		
#M0533S	100 reactions\$179
#M0533L	500 reactions\$719

- Incorporates a universal adaptor sequence at the 3' end of cDNA during the RT reaction
- High sensitivity for cDNA amplification – enables transcriptome analysis by RNA-seq from single cells or as low as 2 pg of human total RNA
- Robust and simple workflow for 5'-Rapid Amplification of cDNA Ends (RACE)
- Retains the complete 5' end of transcripts for 2nd Strand cDNA Synthesis

Description: Template switching reverse transcription (RT) incorporates a universal adaptor sequence to the 3'-end of cDNA. This convenient feature can be utilized in several downstream applications:


- cDNA synthesis and amplification in a one-tube reaction
- 5' Rapid Amplification of cDNA Ends (RACE)
- 2nd strand cDNA synthesis that keeps the 5' end of the transcripts intact

The Template Switching RT Enzyme Mix is optimized for efficient template switching during the RT reaction. The enzyme mix contains RNase Inhibitor in a specially formulated buffer, making reactions easy to setup with no additives needed. It is highly sensitive and specific and can generate RNA-seq libraries from as little as 2 pg of human total RNA or 5'-RACE from 10 ng of total RNA, both with minimal background.

Reagents Supplied:

Template Switching RT Enzyme Mix (10X)

Template Switching RT Buffer (4X)



Primers for cDNA Synthesis

Oligo d(N)_n primers are used for the priming and sequencing of mRNA adjacent to the 3'-poly A tail or tailed cDNA. Note: #S1316 does not contain a 5'-phosphate.

PRODUCT	NEB #	SIZE	PRICE
Random Primer 6 (5' d(N) ₉ 3') ~14.6 nmol	S1230S	1.0 A ₂₆₀ unit	\$101
Random Primer 9 (5' d(N) ₉ 3') ~11.6 nmol	S1254S	1.0 A ₂₆₀ unit	\$101
Oligo d(T) ₂₃ VN	S1327S	1.0 A ₂₆₀ unit	\$101
Random Primer Mix	S1330S	100 µl (60 µM)	\$101
Oligo d(T) ₁₈ mRNA Primer	S1316S	5.0 A ₂₆₀ units	\$95

ProtoScript II First Strand cDNA Synthesis Kit

#E6560S	30 reactions\$171
#E6560L	150 reactions\$685

Companion Products:

Magnetic mRNA Isolation Kit		
#S1550S	25 isolations\$334
RNase Inhibitor, Murine		
#M0314S	3,000 units\$73
#M0314L	15,000 units\$292
Monarch Total RNA Miniprep Kit		
#T2010S	50 preps\$260
Random Primer Mix		
#S1330S	100 µl (60 µM)\$101
Oligo d(T) ₂₃ VN*		
#S1327S	1.0 A ₂₆₀ unit\$101

*Note: V = A, G or C, and N = A, G, C or T

- Enzyme and reaction mixes add flexibility to reaction setup
- Suitable for any PCR format
- Efficient reverse transcription from different starting RNA amounts
- Generates cDNA at least 10 kb

Description: ProtoScript II First Strand cDNA Synthesis Kit features two optimized mixes, ProtoScript II Enzyme Mix and ProtoScript II Reaction Mix. The enzyme mix combines ProtoScript II Reverse Transcriptase and Murine RNase Inhibitor, while the reaction mix contains dNTPs and an optimized buffer. ProtoScript II Reverse Transcriptase is a recombinant M-MuLV reverse transcriptase with reduced RNase H activity and increased thermostability. It can be used to synthesize first strand cDNA at higher temperatures than the wild type M-MuLV. The enzyme is active up to 50°C, providing higher specificity and higher yield of cDNA.

The kit also provides two optimized primers for reverse transcription and nuclease-free water. An anchored Oligo-dT primer [d(T)₂₃ VN] forces the primer to anneal to the beginning of the poly(A) tail. The optimized Random Primer Mix provides random and consistent priming sites covering the entire RNA templates including both mRNAs and non-polyadenylated RNAs. The first strand cDNA product generated is up to 10 kb.



The ProtoScript II First Strand cDNA Synthesis Kit Includes:

- 10X ProtoScript II Enzyme Mix
- 2X ProtoScript II Reaction Mix
- Random Primer Mix (60 µM), Oligo d(T)₂₃ VN Primer (50 µM)** and Nuclease-free Water

**Oligo d(T)₂₃ VN and Random Primer Mix contain 1 mM dNTP

For robust amplification of a wide range of DNA templates, we recommend OneTaq® or Q5® High-Fidelity DNA Polymerases.

ProtoScript First Strand cDNA Synthesis Kit

#E6300S	30 reactions\$162
#E6300L	150 reactions\$645

Companion Products:

Magnetic mRNA Isolation Kit		
#S1550S	25 isolations\$334
Monarch Total RNA Miniprep Kit		
#T2010S	50 preps\$260


- Enzyme and reaction mixes add flexibility to reaction setup
- Suitable for any PCR format
- Efficient reverse transcription from different starting RNA amounts
- Generates cDNA at least 5 kb

Description: ProtoScript First Strand cDNA Synthesis Kit features two optimized mixes, ProtoScript Enzyme Mix and ProtoScript Reaction Mix. ProtoScript Enzyme Mix combines M-MuLV Reverse Transcriptase and RNase Inhibitor, Murine, while ProtoScript Reaction Mix contains dNTPs and an optimized buffer. The kit also contains two optimized primers for reverse transcription and nuclease-free water. An anchored oligo-dT primer [d(T)₂₃ VN] forces the primer to anneal to the beginning of the poly(A) tail. The optimized Random Primer Mix provides random and consistent priming sites covering the entire RNA template including both mRNAs and non-polyadenylated RNAs. The first strand cDNA product generated is over 13.0 kb. This product was formerly known as M-MuLV First Strand cDNA Synthesis Kit.

The ProtoScript First Strand cDNA Synthesis Kit Includes:

- 10X ProtoScript Enzyme Mix
- 2X ProtoScript Reaction Mix
- Random Primer Mix (60 µM), Oligo d(T)₂₃ VN Primer (50 µM)** and Nuclease-free Water

**Oligo d(T)₂₃ VN and Random Primer Mix contain 1 mM dNTP



First Strand cDNA Synthesis with the ProtoScript Kit.

Reactions were performed at 42°C using 2 µg of human spleen total RNA. Negative control reactions (-RT) were carried out with 1X ProtoScript Reaction Mix. A fraction of the first strand cDNA product was used to amplify sequences specific for three different messenger RNAs using 1X LongAmp® Taq 2X Master Mix (NEB #M0287). Lane 1: 1.1 kb of beta-actin gene. Lane 2: -RT control of 1.1 kb of beta-actin gene. Lane 3: 4.7 kb of Xrn-1 gene. Lane 4: -RT control of 4.7 kb of Xrn-1 gene. Lane 5: 9.8 kb of guanine nucleotide exchange factor p532. Lane 6: -RT control of 9.8 kb of guanine nucleotide exchange factor p532. Marker M is 1 kb Plus DNA Ladder (NEB #N3200).

RT-PCR & RT-qPCR Kits



See pages 74–77 for more information.

NEW

Luna® Probe One-Step RT-qPCR Kit (No ROX)
#E3007E 2,500 reactions \$1,908

NEW

Luna Probe One-Step RT-qPCR 4X Mix with UDG
#M3019S 200 reactions \$310
#M3019L 500 reactions \$698
#M3019X 1,000 reactions \$1,241
#M3019E 2,000 reactions \$2,190

NEW

Luna Cell Ready One-Step RT-qPCR Kit
#E3030S 100 reactions \$990

NEW

Luna Cell Ready Probe One-Step RT-qPCR Kit
#E3031S 100 reactions \$939

Luna Universal One-Step RT-qPCR Kit
#E3005S 200 reactions \$244
#E3005L 500 reactions \$549
#E3005X 1,000 reactions \$952
#E3005E 2,500 reactions \$2,116

Luna Universal Probe One-Step RT-qPCR Kit
#E3006S 200 reactions \$220
#E3006L 500 reactions \$495
#E3006X 1,000 reactions \$867
#E3006E 2,500 reactions \$1,908

OneTaq® One-Step RT-PCR Kit
#E5315S 30 reactions \$168

OneTaq RT-PCR Kit
#E5310S 30 reactions \$149


NEW

Luna SARS-CoV-2 RT-qPCR Multiplex Assay Kit
#E3019S 96 reactions \$295
#E3019L 480 reactions \$1,288

RNA Ligase Activity Chart

NEB offers a variety of ligases for DNA and RNA research. Many of these enzymes are recombinant and all offer the quality and value you have come to expect from our products. The chart below highlights reported activities of our T4 ligases ranked by application. A substrate-based selection chart for DNA ligases can be found on page 100.

Reported Activities and Applications for T4 Ligases



<p>T4 RNA Ligase 1</p>	<p>T4 RNA Ligase 2</p>	<p>T4 RNA Ligase 2 Truncated T4 RNA Ligase 2 Truncated K227Q T4 RNA Ligase 2 Truncated KQ</p>
<p>T4 DNA Ligase</p>	<p>Thermostable 5' App DNA/RNA Ligase</p>	<p>T3 DNA Ligase</p>
<p>SplintR Ligase</p>	<p>RtcB Ligase</p>	<p>5' Adenylation Kit</p>

The ligation activities depicted have been reported, but may require optimized reaction conditions.

RNA Ligase Selection Chart

RNA Ligase Selection Chart

	T4 RNA Ligase 1	T4 RNA Ligase 2	T4 RNA Ligase 2 Truncated	T4 RNA Ligase 2, Truncated K227Q	T4 RNA Ligase 2, Truncated KQ	Thermostable 5' App DNA/RNA Ligase	5' Adenylation Kit	SplintR® Ligase	RtcB Ligase
RNA APPLICATIONS									
Ligation of nicks in dsRNA		★★★							
Labeling of 3' termini of RNA	★★★		★	★	★	★			
Ligation of ssRNA to ssRNA	★★★								
Ligation of preadenylated adaptors to RNA	★★		★★	★★	★★★	★★			
5' Adenylation							★★★		
Ligation of 3' P and 5' OH of ssRNA									★★★
DNA APPLICATIONS									
Ligation of preadenylated adaptors to ssDNA						★★★			
DNA/RNA APPLICATIONS									
Joining of RNA & DNA in a ds-structure		★★							
ssDNA Ligation with RNA Splint		★★						★★★	
Ligation of RNA and DNA with 3' P and 5' OH									★★
NGS APPLICATIONS									
NGS Library Prep ssRNA-ssDNA (ligation)	▲		▲	▲	▲				
NGS Library Prep ssRNA-ds-Adaptor splinted ligation		▲							
FEATURES									
Thermostable						✓	✓		
Recombinant	✓	✓	✓	✓	✓	✓	✓	✓	✓

- ★★★ Optimal, recommended ligase for selected application
 ★★ Works well for selected application
 ★ Will perform selected application, but is not recommended
 ▲ Please consult the specific NGS protocol to determine the optimal enzyme for your needs

Michael is a Development Scientist in our Organic Synthesis Division, developing and producing nucleotides, oligonucleotides and beaded resins. When he is not in the lab, Michael enjoys the outdoors, fitness, reading, and music.



T4 RNA Ligase 1 (ssRNA Ligase)



T4 RNA Ligase 1 (ssRNA Ligase)		
#M0204S	1,000 units\$66
#M0204L	5,000 units\$264
T4 RNA Ligase 1 (ssRNA Ligase) High Concentration		
#M0437M	5,000 units\$264

Companion Products:

Adenosine 5'-Triphosphate (ATP)		
#P0756S	1 ml\$37
#P0756L	5 ml\$143
Universal miRNA Cloning Linker		
#S1315S	5 µg\$155

- *Ligation of ss-RNA and DNA*
- *Labeling of 3'-termini of RNA with 5'-[³²P] pCp*
- *Inter- and intramolecular joining of RNA and DNA molecules*
- *Synthesis of ss-oligodeoxynucleotides*
- *Incorporation of unnatural amino acids into proteins*

Description: Catalyzes ligation of a 5' phosphoryl-terminated nucleic acid donor to a 3' hydroxyl-terminated nucleic acid acceptor through the formation of a 3'→5' phosphodiester bond with hydrolysis of ATP to AMP and PPi. Substrates include single-stranded RNA and DNA as well as dinucleoside pyrophosphates.

Reaction Conditions: 1X T4 RNA Ligase Reaction Buffer. Supplement with 1 mM ATP. Incubate at 25°C. May heat inactivate at 65°C for 15 minutes or boiling for 2 minutes.

Notes on Use: Addition of DMSO to 10% (v/v) is required for pCp ligation.

Reagents Supplied with Enzyme:

10X T4 RNA Ligase Reaction Buffer
10 mM ATP (with NEB #M0204) or 100 mM ATP (with NEB #M0437) 50% PEG 8000

Unit Definition: One unit is defined as the amount of enzyme required to convert 1 nmol of 5'-[³²P] rA16 into a phosphate resistant form in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 or 30,000 units/ml

T4 RNA Ligase 2 (dsRNA Ligase)



#M0239S	150 units\$80
#M0239L	750 units\$319

- *Cohesive-end adaptor ligation*
- *Best choice for ligating nicks in dsRNA*
- *Suitable for ligating 3' OH of RNA to 5' phosphate of DNA in a DNA/RNA hybrid*

Description: T4 RNA Ligase 2, also known as T4 Rnl2 (gp24.1), has both intermolecular and intramolecular RNA strand-joining activity. Unlike T4 RNA Ligase 1 (NEB #M0204), T4 RNA Ligase 2 is much more active joining nicks on double stranded RNA than on joining the ends of single stranded RNA. The enzyme requires an adjacent 5' phosphate and 3' OH for ligation. The enzyme can also ligate the 3' OH of RNA to the 5' phosphate of DNA in a double stranded structure.

Reaction Conditions: 1X T4 RNA Ligase 2 Reaction Buffer. Incubate at 37°C. May heat inactivate at 80°C for 5 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to ligate 0.4 µg of an equimolar mix of a 23-mer and 17-mer RNAs in a total reaction volume of 20 µl in 30 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 10,000 units/ml

T4 RNA Ligase 2, truncated



#M0242S	2,000 units\$72
#M0242L	10,000 units\$289

Companion Product:

Universal miRNA Cloning Linker		
#S1315S	5 µg\$155

- *Ligate a pre-adenylated DNA or RNA sequence tag to any RNA 3' end*
- *Join a single-stranded adenylated primer to small RNAs for cDNA library creation*
- *Join a single-stranded adenylated primer to RNA for strand-specific cDNA library construction*

Description: T4 RNA Ligase 2, truncated (T4 Rnl2tr) specifically ligates the pre-adenylated 5' end of DNA or RNA to the 3' end of RNA. The enzyme does not require ATP, but does need the pre-adenylated substrate. T4 Rnl2tr is expressed from a plasmid in *E. coli* which encodes the first 249 amino acids of the full-length T4 RNA Ligase 2. Unlike the full-length ligase, T4 Rnl2tr cannot ligate the phosphorylated 5' end of RNA or DNA to the 3' end of RNA. This enzyme, also known as Rnl2 (1–249), has been used for optimized linker ligation for the cloning of microRNAs. This enzyme reduces background ligation, because it can only use pre-adenylated linkers.

Reaction Conditions: 1X T4 RNA Ligase Reaction Buffer. Incubate at 25°C. May heat inactivate at 65°C for 20 minutes.

Reagents Supplied with Enzyme:

10X T4 RNA Ligase Reaction Buffer
50% PEG 8000

Unit Definition: 200 units is defined as the amount of enzyme required to give 80% ligation of a 31-mer RNA to the pre-adenylated end of a 17-mer [DNA Universal miRNA Cloning Linker (NEB #S1315)] in a total reaction volume of 20 µl in 1 hour at 25°C. Unit assay conditions can be found at www.neb.com.

Concentration: 200,000 units/ml

T4 RNA Ligase 2, truncated K227Q and truncated KQ

T4 RNA Ligase 2, truncated K227Q

#M0351S	2,000 units\$71
#M0351L	10,000 units\$284

T4 RNA Ligase 2, truncated KQ

#M0373S	2,000 units\$71
#M0373L	10,000 units\$284

Companion Product:

Universal miRNA Cloning Linker	
#S1315S	5 µg\$155

- *Ligate a pre-adenylated DNA or RNA sequence tag to any RNA 3' end*
- *Join a single stranded adenylated primer to small RNAs for cDNA library creation*
- *Join a single stranded adenylated primer to RNA for strand-specific cDNA library construction*

Description: T4 RNA Ligase 2, K227Q and truncated KQ (T4 Rnl2tr KQ) specifically ligate the pre-adenylated 5' end of DNA or RNA to the 3' OH end of RNA. The enzymes do not use ATP for ligation, but requires pre-adenylated linkers.

Mutation of K227 in T4 RNA Ligase 2 reduces enzyme lysyl adenylation. K227Q reduces the formation of undesired ligation products (concatemers and circles) by T4 Rnl2tr, by reducing the trace activity of T4 Rnl2tr in transfer of adenylyl groups from linkers to the 5'-phosphates of input RNAs. T4 Rnl2tr KQ is a double-point mutant of T4 RNA Ligase 2, truncated (NEB #M0242). Mutation of R55K in T4 Rnl2tr K227Q increases the ligation activity of the enzyme to levels similar to T4 Rnl2tr.

The exclusion of ATP, use of pre-adenylated linkers, and the reduced enzyme lysyl adenylation activity provide the lowest possible background in ligation reactions. These enzymes have been used for optimized linker ligation for high-throughput sequencing library construction of small RNAs.

Source: Expressed as an MBP fusion from a plasmid in *E. coli* which encodes the first 249 amino acids of the full length T4 RNA Ligase 2. T4 Rnl2tr K227Q has a lysine to glutamine mutation at position 227. T4 Rnl2tr KQ has an arginine to lysine and lysine to glutamine mutation at positions 55 and 227, respectively.

Reaction Conditions: 1X T4 RNA Ligase Reaction Buffer. Incubate at 25°C. May heat inactivate at 65°C for 20 minutes.

Reagents Supplied with Enzyme:

10X T4 RNA Ligase Reaction Buffer
50% PEG 8000

Unit Definition: 200 units is defined as the amount of enzyme required to give 80% ligation of a 31-mer RNA to the pre-adenylated end of a 17-mer DNA [Universal miRNA Cloning Linker (NEB #S1315)] in a total reaction volume of 20 µl in 1 hour at 25°C. Unit assay conditions can be found at www.neb.com.

Concentration: 200,000 units/ml

RtcB Ligase

#M0458S	25 reactions\$74
---------	--------------	-----------

- *Ligate ssRNA or ssDNA with a 3'-phosphate or a 2',3'-cyclic phosphate to the 5'-OH of ssRNA*
- *Circularization of ssRNA with compatible ends*

Description: RtcB Ligase from *E. coli* joins single stranded RNA with a 3'-phosphate or 2',3'-cyclic phosphate to another RNA with a 5'-hydroxyl. Ligation requires both GTP and MnCl₂ and proceeds through a 3'-guanylate intermediate. With substrates having a 2',3'-cyclic phosphate end, hydrolysis to a 3'-phosphate precedes 3' end activation with GMP and ligation.

Source: RtcB Ligase is expressed as His-tagged fusion in *E. coli*.

Reaction Conditions: 1X RtcB Reaction Buffer. Supplement with 0.1 mM GTP and 1 mM MnCl₂ (supplied). Incubate at 37°C.

Reagents Supplied with Enzyme:

10X RtcB Reaction Buffer
MnCl₂ (10 mM)
GTP (10 mM)

Concentration: 15 µM

This is an **Enzyme for Innovation (EFI)**. EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit www.neb.com/EnzymesforInnovation to view the full list.

Thermostable 5' AppDNA/RNA Ligase

#M0319S	10 reactions\$71
#M0319L	50 reactions\$284

Companion Product:

Universal miRNA Cloning Linker	
#S1315S	5 µg\$155

- *Ligation of ssDNA to an adenylated DNA linker for NGS library construction*
- *Ligation of an adenylated linker to RNA at elevated temperatures for small RNA NGS library construction*

Description: Thermostable 5' App DNA/RNA Ligase is a point mutant of catalytic lysine of RNA ligase from *Methanobacterium thermoautotrophicum*. This enzyme is ATP independent. It requires a 5' pre-adenylated linker for ligation to the 3'-OH end of either RNA or single stranded DNA (ssDNA). The enzyme is also active in ligation of RNA with 2'-O-methylated 3' end to 5'-adenylated linkers. The optimal temperature for ligation reaction is 60–65°C. The mutant ligase is unable to adenylate the 5'-phosphate of RNA or ssDNA, which reduces the formation of undesired ligation products (concatemers and circles).

The ability of the ligase to function at 65°C might reduce the constraints of RNA secondary structure in RNA ligation experiments.

Reaction Conditions: 1X NEBuffer 1. Incubate at 65°C.

Reagents Supplied with Enzyme:

10X NEBuffer 1
10X MnCl₂

Concentration: 20 µM

Usage Note: For optimal ligation of ssDNA to preadenylated linkers, we recommend using NEBuffer 1 supplemented with manganese (supplied).

5' DNA Adenylation Kit

#E2610S	10 reactions	\$119
#E2610L	50 reactions	\$476

- Enzymatic 5' adenylation of ss-DNA linkers for next gen sequencing
- One-step reaction gives quantitative adenylation
- Simpler than existing chemical and enzymatic methods
- Reduces need for purification of reaction product

Description: The 5' DNA Adenylation Kit is a simple and efficient enzymatic method for generating 5'-adenylated DNA oligonucleotides using *Mth* RNA ligase, ATP and single stranded 5'-phosphorylated DNA. The kit is optimized to produce the adenylated DNA intermediate with or without 3' terminator. The 5' DNA adenylation kit routinely generates greater than 95% conversion of pDNA to AppDNA. This highly efficient process eliminates the need for gel isolation of the product and increases overall yield.

 **RE** **65°**

The 5' DNA Adenylation Kit Includes:

- *Mth* RNA Ligase (Recombinant)
- 5' DNA Adenylation Reaction Buffer (10X)
- 1 mM ATP

Note: The low turnover of the enzyme requires an approximately equimolar concentration of the enzyme and the oligonucleotide substrate. Adenylated DNA linkers can be used for 3'-end ligation of RNA in cDNA library preparation for next generation sequencing protocols.

SplintR® Ligase

#M0375S	1,250 units	\$98
#M0375L	6,250 units	\$439

- Ligation of adjacent, single-stranded DNA splinted by a complementary RNA
- Characterization of miRNAs and mRNAs, including SNPs

Description: SplintR Ligase, also known as PBCV-1 DNA Ligase or *Chlorella* virus DNA Ligase, efficiently catalyzes the ligation of adjacent, single-stranded DNA oligonucleotides splinted by a complementary RNA strand. This previously unreported activity may enable novel approaches for characterization of miRNAs and mRNAs, including SNPs. SplintR is ideally suited for many target enrichment workflows with applications in next-generation sequencing and molecular diagnostics. The robust activity of the enzyme and its affinity for RNA-splinted DNA substrates (apparent $K_m = 1$ nM) enable sub-nanomolar detection of unique RNA species within a complex mixture, making SplintR ligase a superior choice for demanding RNA detection technologies.

 **RE** **25°** **65°**

Reaction Conditions: 1X SplintR Ligase Reaction Buffer. Incubate at 25°C. May heat inactivate at 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme needed to ligate (to 50% completion) 2 picomoles of a tripartite FAM-labeled DNA:RNA hybrid substrate in a 20 µl reaction at 25°C in 15 minutes in 1X SplintR Ligase Reaction Buffer. Unit assay conditions can be found at www.neb.com.

Concentration: 25,000 units/ml

RNA 5' Pyrophosphohydrolase (RppH)

#M0356S	200 units	\$95
---------	-----------------	------

- Conversion of 5'-triphosphate RNA to monophosphate RNA
- Preparation of 5'-phosphate RNA for ligation
- Characterization of RNA 5' ends

Description: The bacterial RNA 5' Pyrophosphohydrolase (RppH) removes pyrophosphate from the 5' end of triphosphorylated RNA to leave a 5' monophosphate RNA. The RppH protein was also known as NudH/YgdP which can split diadenosine penta-phosphate to ADP and ATP.

Reaction Conditions: 1X NEBuffer 2. Incubate at 37°C.

 **RE** **37°**

Reagents Supplied with Enzyme:

10X NEBuffer 2

Unit Definition: One unit is defined as the amount of enzyme that converts 1 µg 300-mer RNA transcript into an XRN-1 digestible RNA in 30 minutes at 37°C.

Concentration: 5,000 units/ml

5' Deadenylase

#M0331S	2,500 units	\$73
---------	-------------------	------

- Deadenylation of 5' end of DNA and RNA
- Aprataxin-dependent DNA repair assay
- Analysis of dinucleoside tetraphosphate

Description: Yeast 5' Deadenylase is a member of the HIT (histidine triad) family of proteins and specifically a member of the Hint branch. It is the yeast orthologue of aprataxin. Mutations in human aprataxin have been known to be involved in the neurological disorder known as ataxia oculomotor apraxia-1. The human protein has been shown to resolve abortive ligation intermediates by removing the AMP at the 5' end of DNA (AMP-DNA hydrolase activity). It also repairs DNA damage at 3' ends by removing 3'-phosphate and 3'-phosphoglycolate. Human aprataxin acts on small molecules, such as nucleotide polyphosphates diadenosine tetraphosphate (AppppA) and lysyl-AMP.

 **RE** **30°** **70°**

The 5' Deadenylase is encoded by the *HNT3* gene of *S. cerevisiae*. NEB has shown this protein is capable of deadenylation from 5' end of DNA and RNA, leaving the phosphate at 5' end. It also cleaves AppppA into ATP and AMP. Its activity on lysyl-AMP is not detectable.

Reaction Conditions: 1X NEBuffer 1 and 5–50 pmol adenylated DNA (AMP-DNA) in 20 µl. Incubate at 30°C. May heat inactivate at 70°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to remove 10 pmoles of AMP from a 5'-adenylated DNA oligo in 10 minutes at 30°C.

Concentration: 50,000 units/ml

RNase I_f

#M0243S	5,000 units\$69
#M0243L	25,000 units\$276

- Eliminates RNA from DNA and protein preparations
- Degradation of single-stranded RNA to mono-, di- and trinucleotides
- Used in ribonuclease protection assays

Description: Ribonuclease I_f (RNase I_f) is an RNA endonuclease which will cleave at all RNA dinucleotide bonds leaving a 5' hydroxyl and 2', 3' cyclic monophosphate. It has a preference for single-stranded RNA over double-stranded RNA. RNase I_f is a recombinant protein fusion of RNase I (from *E. coli*) and maltose-binding protein. It has identical activity to RNase I.

Reaction Conditions: 1X NEBuffer 3. Incubate at 37°C. May heat inactivate at 70°C for 20 minutes.



Notes on Use: RNase I_f will not degrade DNA. It has a strong preference for single-stranded RNA over double-stranded RNA.

Unit Definition: One unit is defined as the amount of enzyme required to fully digest 1 picomole of synthetic ssRNA 33-mer in a total reaction volume of 10 µl in 15 minutes in 1X NEBuffer 3 as visualized on a 20% acrylamide gel (40:1 Bis) stained with SYBR® Gold. Unit assay conditions can be found at www.neb.com.

Concentration: 50,000 units/ml

SYBR® is a registered trademark of Molecular Probes, Inc.

RNase H

#M0297S	250 units\$72
#M0297L	1,250 units\$289

- Removal of poly(A) tails of mRNA hybridized to poly(dT)
- Removal of mRNA during second strand cDNA synthesis

Description: Ribonuclease H (RNase H) is an endoribonuclease which specifically hydrolyzes the phosphodiester bonds of RNA which is hybridized to DNA. This enzyme does not digest single- or double-stranded DNA.

Reaction Conditions: 1X RNase H Reaction Buffer. Incubate at 37°C. May heat inactivate at 65°C for 20 minutes.



Unit Definition: One unit is defined as the amount of enzyme required to produce 1 nmol of ribonucleotides from 20 picomoles of a fluorescently labelled 50 base pair RNA-DNA hybrid in a total reaction volume of 50 µl in 20 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 5,000 units/ml

Thermostable RNase H

#M0523S	250 units\$144
---------	-----------	------------

- Higher stringency RNA structure mapping and site-specific RNA cleavage
- Removal of poly(A) tails from mRNA hybridized to oligo(dT)
- Removal of mRNA during second strand cDNA synthesis
- Component of isothermal amplification methods

Description: Thermostable RNase H specifically recognizes and cleaves the phosphodiester bonds of an RNA strand in an RNA-DNA hybrid while leaving the DNA strand intact. This thermostable nuclease exhibits the same enzymatic properties as *E. coli* RNase H, but is active at much higher temperatures.

Reaction Conditions: 1X RNase H Reaction Buffer. Incubate at 50°C.



Unit Definition: One unit is defined as the amount of enzyme required to produce 1 nmol of ribonucleotides from 40 picomoles of a fluorescently labelled 25 base pair RNA-DNA hybrid in a total reaction volume of 50 µl in 20 minutes at 50°C. Unit assay conditions can be found at www.neb.com.

Concentration: 5,000 units/ml

RNase HII

#M0288S	250 units\$75
#M0288L	1,250 units\$302

- Nicking of products generated with a polymerase that will incorporate ribonucleotides
- Generation of a double-stranded break at the site of an incorporated ribonucleotide when used with T7 Endo I
- Degradation of the RNA portion of Okazaki fragments or other RNA-DNA hybrids

Description: Ribonuclease HII (RNase HII) is an endoribonuclease that preferentially nicks 5' to a ribonucleotide within the context of a DNA duplex. The enzyme leaves 5' phosphate and 3' hydroxyl ends. RNase HII will also nick at multiple sites along the RNA portion of an Okazaki fragment.

Source: An *E. coli* strain containing a genetic fusion of the RNase HII gene (rnhB) from *E. coli* and the gene coding for maltose binding protein (MBP). Following affinity chromatography, RNase HII is cleaved from the fusion construct by Factor Xa and then purified away from both MBP and Factor Xa. RNase HII cleaved from MBP has four additional amino acids at its N-terminus (Ile-Ser-Glu-Phe).



Reaction Conditions: 1X ThermoPol Reaction Buffer. Incubate at 37°C.

Unit Definition: One unit is defined as the amount of enzyme required to yield a fluorescence signal consistent with the nicking of 100 pmol of synthetic double-stranded DNA substrate containing a single ribonucleotide near the quencher of a fluorophore/quencher pair in 30 minutes at 37°C in 1X ThermoPol Reaction Buffer. Unit assay conditions can be found at www.neb.com.

Concentration: 5,000 units/ml

Note: Incubation with 0.1% SDS is sufficient to inactivate RNase HII.

Phosphorylation and Dephosphorylation

See pages 104–106 for more information.

NEW

Quick CIP

#M0525S	1,000 units\$90
#M0525L	5,000 units\$430

Antarctic Phosphatase

#M0289S	1,000 units\$73
#M0289L	5,000 units\$295

Shrimp Alkaline Phosphatase (rSAP)

#M0371S	500 units\$63
#M0371L	2,500 units\$253

T4 Polynucleotide Kinase

#M0201S	500 units\$57
#M0201L	2,500 units\$228

ShortCut® RNase III

#M0245S	200 units\$128
#M0245L	1,000 units\$515

- Generates siRNAs for RNA interference studies
- Gene silencing
- Target validation
- Removal of long dsRNAs

Description: ShortCut RNase III converts long double-stranded RNA into a heterogeneous mix of short (18–25 bp) interfering RNAs (siRNA) suitable for RNA interference in mammalian cells. 1.5 units (1 µl) of ShortCut RNase III is sufficient to convert 1 µg of dsRNA into siRNA.

Source: An *E. coli* strain containing a genetic fusion of the *E. coli* RNase III gene (*rnc*) and the gene coding for maltose binding protein (MBP).

Reaction Conditions: 1X ShortCut RNase III Reaction Buffer. Supplement with 20 mM MnCl₂ (supplied). Incubate at 37°C.


Unit Definition: One unit is the amount of enzyme required to digest 1 µg of dsRNA to siRNA in a total reaction volume of 50 µl in 20 minutes at 37°C. Unit assay conditions can be found at www.neb.com.

Concentration: 2,000 units/ml

RR 37°

Advantages of the ShortCut RNase III:

- Make effective siRNAs against any gene target
- Heterogeneous population of siRNA ensures silencing of target gene
- From DNA template to transfection in just 1 day
- Eliminates trial and error approach of synthetic siRNA



siRNA production by ShortCut RNase III: A) Varying amounts of ShortCut RNase III were incubated with 2 µg of a 500 bp dsRNA for 20 minutes. B) dsRNA fragments (1 kb and 175 bp) were digested with ShortCut RNase III. Digests were analyzed by 20% TBE polyacrylamide gel electrophoresis.

XRN-1

#M0338S	20 units\$92
#M0338L	100 units\$366

- Removal of RNA containing 5' monophosphate from an RNA mixture

Description: XRN-1 is a highly processive 5' to 3' exoribonuclease, requiring a 5' monophosphate. It also acts on 5' monophosphate ssDNA with reduced efficiency.

Reaction Conditions: 1X NEBuffer 3. Incubate at 37°C. May heat inactivate at 70°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme that digests 1 µg of phosphorylated yeast RNA in 60 minutes at 37°C.

Concentration: 1,000 units/ml

RR 37°

Exonuclease T



#M0265S	250 units\$74
#M0265L	1,250 units\$299

- *Generate blunt ends in DNA or RNA with 3' overhangs*

Description: Exonuclease T (Exo T), also known as RNase T, is a single-stranded RNA or DNA specific nuclease that requires a free 3' terminus and removes nucleotides in the 3'→5' direction. Exo T can be used to generate blunt ends from RNA or DNA having 3' extensions.

Source: Exonuclease T is overexpressed and purified as a C-terminal fusion to maltose-binding protein (MBP). MBP is removed from Exo T by Factor Xa cleavage and then purified. When cleaved from MBP has an additional amino acid on the N-terminus and a Phe instead of a Met (i.e. Glu-Phe-Exo T instead of Met-Exo T).

Reaction Conditions: 1X NEBuffer 4. Incubate at 25°C. May heat inactivate at 65°C for 20 minutes.

Unit Definition: One unit is defined as the amount of enzyme required to produce 0.1 nmol of TCA soluble nucleotides from 1 nmol of [³H]-labeled polythymidine in a total reaction volume of 100 µl in 30 minutes at 25°C in 1X NEBuffer 4 with 1 nmol [³H]-labeled polythymidine DNA.

Concentration: 5,000 units/ml

Usage Note: Exo T has different activity on RNA vs. DNA. For RNA, 1 unit of Exo T is required to completely digest 1.0 pmol of rA20 under standard reaction condition as measured by gel electrophoresis.

Nucleoside Digestion Mix



#M0649S	50 reactions\$121
---------	--------------	------------

- *Convenient one-step protocol*
- *Digests both DNA and RNA to single nucleosides*
- *Low-glycerol formulation significantly reduces glycerol-induced ion suppression during MS analysis*

Description: The Nucleoside Digestion Mix is a mixture of enzymes that provides a convenient one-step method to generate single nucleosides from DNA or RNA. Optimized for quantitative analysis by liquid chromatography-mass spectrometry (LC-MS), this reagent eliminates the need for sequential multi-step, time-consuming digestion protocols. The Nucleoside Digestion Mix digests ssDNA, dsDNA, DNA/RNA hybrids and RNA (except mRNA cap structures) containing epigenetically modified (m5C,

hm5C, f5C, ca5C, m4C, m6A, etc.), unnatural, or damaged bases. Moreover, the low-glycerol formulation (< 1%) significantly reduces glycerol-induced ion suppression during mass spectrometry analysis.

Reaction Conditions: 1X Nucleoside Digestion Mix Reaction Buffer. Incubate at 37°C. May heat inactivate at 65°C for 20 minutes.

DNase I (RNase-free)



#M0303S	1,000 units\$70
#M0303L	5,000 units\$282

- *Degradation of DNA template in transcription reactions*
- *Removal of contaminating genomic DNA from RNA samples*
- *DNase I footprinting*
- *Nick Translation*

Description: DNase I (RNase-free) is an endonuclease that nonspecifically cleaves DNA to release di-, tri- and oligonucleotide products with 5'-phosphorylated and 3'-hydroxylated ends. DNase I acts on single- and double-stranded DNA, chromatin and RNA:DNA hybrids.

Reaction Conditions: 1X DNase I Reaction Buffer. Incubate at 37°C. May heat inactivate at 75°C for 10 minutes.

Unit Definition: One unit is defined as the amount of enzyme which will completely degrade 1 µg of pBR322 DNA in a total reaction volume of 50 µl in 10 minutes at 37°C.

Complete degradation is defined as the reduction of the majority of DNA fragments to tetranucleotides or smaller.

Concentration: 2,000 units/ml

Note: EDTA should be added to a final concentration of 5 mM to protect RNA from being degraded during enzyme inactivation.

RNase Inhibitor, Murine



#M0314S	3,000 units\$73
#M0314L	15,000 units\$292

- *Inhibits common eukaryotic RNases*
- *Compatible with Taq Polymerase, AMV or M-MuLV Reverse Transcriptases*
- *cDNA synthesis & RT-PCR*
- *In vitro transcription/translation*
- *Enzymatic RNA labeling reaction*

RNase Inhibitor, Murine has significantly improved resistance to oxidation compared to human & porcine RNase inhibitors.

Description: RNase Inhibitor, Murine is a 50 kDa recombinant protein of murine origin. It specifically inhibits RNases A, B and C by binding noncovalently in a 1:1 ratio with high affinity. It is not effective against RNase 1, RNase T1, S1 Nuclease, RNase H or RNase from *Aspergillus*. No inhibition of polymerase activity is observed when used with Taq DNA Polymerase, AMV or M-MuLV Reverse Transcriptases, or Phage RNA Polymerases (SP6, T7, or T3).

Recombinant RNase Inhibitor, Murine does not contain the pair of cysteines identified in the human version that are very sensitive to oxidation and lead to inactivation of the inhibitor. As a result, RNase Inhibitor, Murine has significantly improved resistance to oxidation compared

to the human/porcine RNase inhibitors, and is stable at low DTT concentrations (< 1 mM). This makes it ideal for reactions where high concentration DTT is adverse to the reaction (e.g., RT-qPCR).

Unit Definition: One unit is defined as the amount of RNase Inhibitor, Murine required to inhibit the activity of 5 ng of RNase A by 50%. Activity is measured by the inhibition of hydrolysis of cytidine 2', 3'-cyclic monophosphate by RNase A.

Concentration: 40,000 units/ml

RNase Inhibitor, Human Placenta



#M0307S	2,000 units\$83
#M0307L	10,000 units\$330

- Inhibits common eukaryotic RNases
- Compatible with Taq Polymerase, AMV or M-MuLV Reverse Transcriptases
- Active over a broad pH range (pH 5–8)
- cDNA synthesis reactions
- In vitro transcription/translation

Description: RNase Inhibitor, Human Placenta is a recombinant human placental protein which specifically inhibits ribonucleases (RNases) A, B and C. It is not effective against RNase 1, RNase T1, S1 Nuclease, RNase H or RNase from *Aspergillus*. In addition, no inhibition of polymerase activity is observed when RNase Inhibitor, Human Placenta is used with Taq DNA Polymerase, AMV or M-MuLV Reverse Transcriptases, or Phage RNA Polymerases (SP6, T7, or T3).

The 50 kDa protein inhibits RNases by binding noncovalently in a 1:1 ratio with an association constant greater than 10^{14} .

Unit Definition: One unit is defined as the amount of RNase Inhibitor, Human Placenta required to inhibit the activity of 5 ng of RNase A by 50%. Activity is measured by the inhibition of hydrolysis of cytidine 2', 3'-cyclic monophosphate by RNase A.

Concentration: 40,000 units/ml

Ribonucleoside Vanadyl Complex

#S1402S	10 ml (200 mM)\$88
---------	----------------	-----------

Ribonucleoside Vanadyl Complex is an equimolar mixture of all four ribonucleosides, complexed with oxovanadium IV by a modification of the procedures by Berger (1). Each lot of the complex is assayed for oxovanadium V content and inhibition of ribonuclease activity.

Vanadium complexes are used in mRNA purifications as exogenous ribonuclease inhibitors. They are compatible with cell lysis techniques and with sucrose gradient fractionation of cytoplasmic components.

Reference:

- (1) Berger, S.L. and Birkenmeier, C.S. (1979) *Biochemistry* 18, 5143–5149.

NEBNext® Reagents for RNA Library Preparation

NEBNext Ultra II Directional RNA Library Prep Kit for Illumina		
#E7760S	24 rxns\$935
#E7760L	96 rxns\$3,175

NEBNext Ultra II Directional RNA Library Prep with Sample Purification Beads		
#E7765S	24 rxns\$1,035
#E7765L	96 rxns\$3,525

NEBNext Ultra II RNA Library Prep Kit for Illumina		
#E7770S	24 rxns\$890
#E7770L	96 rxns\$3,025

NEBNext Ultra II RNA Library Prep with Sample Purification Beads		
#E7775S	24 rxns\$985
#E7775L	96 rxns\$3,350


NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina		
#E6420S	24 rxns\$1,195
#E6420L	96 rxns\$3,965

NEBNext Single Cell/Low Input cDNA Synthesis & Amplification Module		
#E6421S	24 rxns\$650
#E6421L	96 rxns\$2,210

NEBNext Kits for RNA sample preparation for next generation sequencing keep pace with the use of ever-decreasing input amounts and sub-optimal sample quality, along with the need for superior performance, reliability, and automation compatibility. The fast and streamlined Ultra II Workflow is at the heart of our RNA library prep kits, including our NEBNext Single Cell/Low Input Library Prep Kit for Illumina, and these are all

available in flexible, user-friendly formats. Our NEBNext RNA products also include kits for Small RNA Library Prep, and multiple options for removal of abundant RNAs as well as a growing range of adaptors and primers.

KAPA™ is a trademark of Kapa Biosystems. ILLUMINA® and TRUSEQ® are registered trademarks of Illumina, Inc. AGILENT® and BIOANALYZER® are registered trademarks of Agilent Technologies, Inc.



NEBNext Ultra II Directional RNA produces the highest yields, from a range of input amounts.

Poly(A)-containing mRNA was isolated from 10 ng, 100 ng and 1 µg of Universal Human Reference RNA (Agilent #740000) and libraries were made using the NEBNext Ultra II Directional RNA kit, Kapa™ Stranded mRNA-Seq kit, Kapa mRNA HyperPrep kit and Illumina® TruSeq® Stranded mRNA Kit. The input RNA amount and number of PCR cycles are indicated. Library yields from an average of three replicates are shown. Error bars indicate standard deviation. Library yields were assessed using the Agilent® Bioanalyzer®.



Kits for rRNA depletion from a variety of samples, including human, mouse, rat and bacteria are available, as well as a customizable option for any sample type. Kits for Small RNA library construction are also available. More information on NEBNext reagents for library preparation can be found on pages 144–175.

Magnetic mRNA Isolation Kit

#S1550S 25 isolations\$334

Component Sold Separately:

Oligo d(T)₂₅ Magnetic Beads
#S1419S 25 mg\$307

Companion Products:


6-tube Magnetic Separation Rack
#S1506S 6 tubes (1.5 ml)\$215
12-tube Magnetic Separation Rack
#S1509S 12 tubes (1.5 ml)\$340

- Suitable for automated high-throughput applications
- Eliminates need for organic solvents
- No need to precipitate poly(A)⁺ transcripts in eluent
- Obtain intact poly(A)⁺ RNA in < 1 hour
- Negligible gDNA contamination

Description: The New England Biolabs Magnetic mRNA Isolation Kit is designed to isolate intact poly(A)⁺ RNA from cells and tissue without requiring phenol or other organic solvents. The technology is based on the coupling of Oligo d(T)₂₅ to 1 µm paramagnetic beads which are then used as the solid support for the direct binding of poly(A)⁺ RNA. Thus, the procedure permits the manual processing of multiple samples and can be adapted for automated high-throughput applications. Additionally, magnetic separation technology permits elution of intact mRNA in small volumes eliminating the need for precipitating the poly(A)⁺ transcripts in the eluent. Intact poly(A)⁺ RNA which is fully representative of the mRNA population of the original sample can be obtained in less than one hour. Oligo d(T)₂₅ Magnetic Beads can be reused up to three times and the researcher has the option of eluting the isolated mRNA or using the bound dT DNA as a primer in a first-strand cDNA reaction.

The Magnetic mRNA Isolation Kit Includes:

- Oligo d(T)₂₅ Magnetic Beads
- Lysis/Binding Buffer
- Wash Buffer I, II, III, Elution Buffer



Consistency and wide isolation range are demonstrated by poly(A)⁺ RNA isolation from duplicate samples of decreasing numbers of HEPG2 cells (5×10^5 to 1×10^3) by direct lysis/binding in microtiter plates followed by mRNA isolation with the magnetic method. 1/10th of isolated mRNA is converted to oligo (dT) primed cDNA using ProtoScript M-MuLV First Strand cDNA Synthesis Kit (NEB #E6300) and qPCR done with validated primers for the peptidylpropyl isomerase, a low-abundance housekeeping gene.

Oligo d(T)₂₅ Magnetic Beads

#S1419S 25 mg\$307

Companion Product:

96-well Microtiter Plate Magnetic Separation Rack
#S1511S 96 well\$509

Description: An affinity matrix for the small-scale isolation of mRNA from crude cell lysates and tissue. The isolation occurs through the hybridization of covalently coupled oligo d(T)₂₅ to the poly(A) region present in most eukaryotic mRNAs. Applications include direct mRNA isolation following lysis and second-round purification of previously isolated total RNA. The magnetic separation technology is scalable and permits elution of intact mRNA in small volumes eliminating the need for precipitation of the isolated mRNA. Beads can be reused up to three times, and the researcher has the option of eluting the isolated mRNA or using the bound d(T)₂₅ as a primer in a first-strand cDNA reaction.

Beads are supplied as a 5 mg/ml suspension in phosphate buffer (PBS) (pH 7.4), containing 0.05% Tween-20 and 0.05% Na₃N.

Support Matrix: 1 µm nonporous superparamagnetic microparticles.

Binding Capacity: 1 mg of Oligo d(T)₂₅ Beads will bind 10 µg of poly(A)⁺ RNA.

Streptavidin Magnetic Beads

#S1420S 5 ml (20 mg)\$334

Companion Products:

6-Tube Magnetic Separation Rack
#S1506S 6 tubes (1.5 ml)\$215
96-well Microtiter Plate Magnetic Separation Rack
#S1511S 96 well\$509

Description: Streptavidin Magnetic Beads are 1 µm superparamagnetic particles covalently coupled to a highly pure form of streptavidin. The beads can be used to capture of biotin labeled substrates including antigens, antibodies and nucleic acids. The strength of the biotin-streptavidin interaction coupled with low non-specific binding permits captured substrates to be useful as ligands in subsequent experiments including mRNA isolation and the capture of primary or secondary antibodies.

Beads are supplied as a 4 mg/ml suspension in phosphate buffer (PBS) (pH 7.4) containing 0.1% BSA, 0.05% Tween-20 and 0.05% Na₃N.

Support Matrix: 1 µm non-porous superparamagnetic microparticle.

Binding Capacity: The beads will bind greater than 1000 pmol of free biotin per mg and greater than 500 pmol of single-stranded 25 bp biotinylated oligonucleotide per mg.

EpiMark® N6-Methyladenosine Enrichment Kit

#E1610S 20 reactions\$431

- Enrichment for m6A modified RNA in immunoprecipitation protocols
- Enriched RNA can be used directly for next gen sequencing or RT-qPCR

Description: The EpiMark N6-Methyladenosine Enrichment Kit contains a rabbit monoclonal antibody specific for N6-Methyladenosine (m6A). The kit also contains two control RNAs, one with m6A modification (*Gaussia* luciferase) and one without (*Cypridina* luciferase) to monitor enrichment and depletion. The GLuc RNA control was transcribed in the presence of 20% m6ATP and 80% ATP.

This kit can be used to enrich m6A modified RNA in immunoprecipitation protocols for downstream analysis by next-generation RNA sequencing or RT-qPCR. Modified RNA is isolated from a fragmented RNA sample by binding to the N6-Methyladenosine

Epi

antibody attached to Protein G Magnetic Beads. After multiple wash and clean-up steps, the enriched RNA is eluted in nuclease-free water and is ready for further analysis.

The EpiMark N6-Methyladenosine Enrichment Kit Includes:

- N6-Methyladenosine Antibody
- m6A Control RNA (100 nM)
- Unmodified Control RNA (100 nM)

N6-Methyladenosine Antibody is produced by Cell Signaling Technology, Inc. and sold by New England Biolabs, Inc.

Monarch Kits for Cleanup & Isolation

See pages 137–138 for more information.

Monarch Total RNA Miniprep Kit
#T2010S 50 preps\$260

Monarch RNA Cleanup Kit (10 µg)
#T2030S 10 preps\$54
#T2030L 100 preps\$290

Monarch RNA Cleanup Kit (50 µg)
#T2040S 10 preps\$54
#T2040L 100 preps\$290

Monarch RNA Cleanup Kit (500 µg)
#T2050S 10 preps\$60
#T2050L 100 preps\$480

The Monarch Total RNA Miniprep Kit is a comprehensive solution for sample preservation, cell lysis, gDNA removal, and purification of total RNA from a wide variety of biological samples, including cultured cells, blood, and mammalian tissues. Additionally, tough-to-lyse samples, such as bacteria, yeast, and plant, can be processed with additional steps that enhance lysis. Total RNA, including viral RNA, can also be extracted from clinically-relevant samples like saliva, buccal swabs and nasopharyngeal swabs. Cleanup of enzymatic reactions or purification of RNA from TRIzol®-extracted samples is also possible using this kit. Purified RNA has high quality metrics, including $A_{260/280}$ and $A_{260/230}$ ratios ≥ 1.8 , high RIN scores, and minimal residual gDNA. Captured RNA ranges in size from full-length rRNAs down to intact miRNAs. Additionally, differential binding conditions allow selective capture or exclusion of the sub-200 nucleotide RNA pool that includes miRNA, 5S rRNA, and tRNA. Purified RNA is suitable for downstream applications, such as RT-qPCR, cDNA synthesis, RNA-seq, Northern blot analysis, etc.

The Monarch RNA Cleanup Kits provide a fast and simple silica spin column-based solution for RNA cleanup and concentration after any enzymatic reaction (including *in vitro* transcription, DNase I treatment, capping and labeling) and after other purification methods such as phenol/chloroform extraction. These kits can also be used to extract total RNA from cells, saliva and buccal/nasopharyngeal swabs. The Monarch RNA Cleanup Kits are available in 3 different binding capacities: 10 µg, 50 µg and 500 µg. Each kit contains unique columns, all designed to prevent buffer retention and ensure no


carryover of contaminants, enabling low-volume elution of highly-pure RNA. Following the standard protocol, RNA ≥ 25 nt is purified with this kit; however, a modified protocol is available to enable the binding of RNA as small as 15 nt (including miRNAs).

The Monarch Total RNA Miniprep Kit Includes:

- gDNA Removal Columns
- RNA Purification Columns
- Collection Tubes II
- DNA/RNA Protection Reagent (2X)
- RNA Lysis Buffer
- Proteinase K and associated buffers
- DNase I & associated reaction buffer
- RNA Priming Buffer
- RNA Wash Buffer (5X)
- Nuclease-free Water

The Monarch RNA Cleanup Kits Include:

- RNA Cleanup Columns (10, 50 or 500 µg)
- RNA Cleanup Binding Buffer
- RNA Cleanup Wash Buffer (5X)
- Collection Tubes II
- Nuclease-free Water



The Monarch Total RNA Miniprep Kit successfully purifies synthetic SARS-CoV-2

viral RNA from saliva samples: The Monarch Total RNA Miniprep Kit Proteinase K and Quick Protocols were used to isolate total RNA from saliva samples containing 10-fold serial dilutions of synthetic SARS-CoV-2 N-gene RNA. Purified RNA was eluted in 100 µl nuclease-free water to yield 50 to 500,000 copies of viral RNA/µl. Using the Luna Universal Probe One-Step RT-qPCR Kit (NEB #E3006), titers as low as 50 copies (the lowest input tested) were detected and linear, quantitative recovery of the SARS-CoV-2 N-gene was observed over a 5-Log range.


RNA Markers & Ladders

dsRNA Ladder	
#N0363S	25 gel lanes\$96
microRNA Marker	
#N2102S	100 gel lanes\$72
ssRNA Ladder	
#N0362S	25 gel lanes\$72
Low Range ssRNA Ladder	
#N0364S	100 gel lanes\$71

NEB offers several RNA Markers and Ladders with a size range from 17 to 9,000 bases. The ssRNA ladders are suitable for use as RNA size standards on denaturing or native gels. Both are supplied with 2X Loading Buffer and feature a higher intensity fragment to serve as a reference band. The microRNA Marker, supplied in a ready-to-load denaturing solution, is ideally used as a size marker on polyacrylamide gels or Northern blots and is best visualized stained with ssRNA fluorescent dyes. It is

supplied with a 3'-biotinylated 21-mer oligonucleotide probe that can also be labeled with [γ - 32 P] ATP and T4 PNK (NEB #M0201). The dsRNA Ladder is suitable for use as a size standard in dsRNA analysis on both polyacrylamide and agarose gels.

Concentration: Low Range ssRNA Ladder and dsRNA Ladder are supplied at 500 μ g/ml. ssRNA Ladder is supplied at 2,000 μ g/ml. MicroRNA Marker is supplied at 12 ng/ μ l.



RNA Loading Dye (2X)

#B0363S	4 ml\$51
---------	----------------

Description: The RNA Loading Dye, (2X) is a premixed loading dye for use with denaturing and non-denaturing PAGE/agarose gels.

RNA Loading Dye Composition:
1X RNA Loading Dye:
47.5% Formamide, 0.01% SDS, 0.01% Bromophenol Blue, 0.005% xylene cyanol and 0.5 mM EDTA.

Universal miRNA Cloning Linker

#S1315S	0.83 nmol\$155
Companion Product:	
T4 RNA Ligase 2, truncated KQ	
#M0373S	2,000 units\$71
#M0373L	10,000 units\$284

This 5'-adenylated, 3'-blocked oligoribonucleotide can be used for cloning short RNAs according to the procedure of Bartel (1). RNA ligase recognizes the "activated" adenylated oligo and covalently links (ligates) its 5' end to the 3' OH of a second single stranded sequence in the absence of ATP. In a mixture of nucleic acids use of the 5' adenylated, 3' blocked oligo with T4 RNA Ligase 2,

truncated, T4 RNA Ligase 2, truncated K227Q or T4 RNA Ligase 1 (w/o ATP) results in ligation of the target oligo only.

The sequence of the adenylated DNA oligo is 5'-rAppCTGTAGGCACCATCAAT-NH₂ 3'.

Reference:
(1) Lau et al. (2001) *Science*, 294, 858-856.

NEB

ENGLAND

BIO
LABS



Aim AR camera at this page for a slideshow
of artwork on the NEB campus.

The Intersection of Art and Science

Science and art are both driven by curiosity, discovery and a thirst for knowledge. They are both attempts to communicate truth and understanding of ourselves and the world around us. And, even though the preparation and process leading to scientific discovery and artistic expression is actually surprisingly similar, they are often viewed as separate disciplines.

When we are young, we are both scientists and artists, experimenting and creating. We are less likely to think of ourselves as just one or the other. As we get older, our formal education directs us down a narrower path — one can become either an artist or a scientist, but rarely both. The separation in communication between the artistic and scientific cultures has become a great chasm since the Industrial Revolution due to intellectual specialization. However, there is great potential for both cultures to inform and inspire each other.

Don Comb always felt a strong connection between science and art; he once said “Science and art are both trying to say something — science is trying to tell you the truth, and artists are trying to present what they feel is the truth.” Don celebrated his love of art by filling the NEB campus with sculptures, paintings, ceramics and artifacts from all over the world — an eclectic collection that was curated over many years. When approaching the main building at NEB headquarters, one immediately notices the enormous mosaic incorporated into the patina copper exterior, which was selected by an international art design competition. The main lobby includes access to an art gallery, which in addition to some of the work mentioned above, also features revolving art exhibitions from local and regional artists.

NEB also supports many local art-based organizations and schools. One example is Raw Art Works (RAW), an art therapy-based organization for youth — helping kids share their stories. Teen muralists from Raw Art spent time with Don and were treated to his indomitable enthusiasm, warmth and joy. After learning about NEB, they brought their artistic visions to a painted mural that translated their understanding of the company’s science and culture. Additionally, NEB has long supported the youth of RAW with scholarships to attend Monserrat College of Art, another organization supported by NEB.

Art brings joy and can make you happy — when we are happy, we are more creative and productive. Neurobiologist Semir Zeki from the University College in London found that looking at a work of art increased dopamine levels and heightened activity in the brain’s frontal cortex. He concluded that viewing art can give as much pleasure as being in love. Did Don realize that by sharing his love of art, he created an environment that brought happiness and inspiration to the workplace and instigated lively conversations, especially around his more provocatively chosen pieces?

When we allow science to inspire art and art to inform science, we can approach a challenge with a diverse set of intellectual and creative problem-solving skills. It is in this space, where the lines between art and science are blurred, that innovation can lead to unique, groundbreaking discoveries.

Mural created by students at
Raw Art Works (RAW).
Credit: Mileidy Rodriguez, NEB

View a slideshow
of artwork on the
NEB campus.



Protein Expression & Purification



NEB offers an array of solutions for robust expression of your target protein.

At first glance, recombinant protein expression looks quite simple. Essentially, DNA encoding a target protein is cloned downstream of a promoter in an expression vector. This vector is then introduced into a host cell, and the cell's protein synthesis machinery produces the desired protein. In practice, however, protein expression can be very challenging, because so many factors may influence the process. For example, each protein folds in its own unique manner, a process that may be influenced by the choice of expression host. Similarly, some proteins require post-translational modifications or proper insertion into a biological membrane. Finally, some proteins may have an activity that is detrimental to the host. Thus, no single solution exists for successful production of all recombinant proteins. Instead, it is beneficial to have access to a wide range of expression tools, and a willingness to explore multiple approaches to better one's chances for success.

NEB offers an array of expression systems offering different advantages, enabling you to choose the strategy that best suits your protein expression and purification needs. Many share a compatible polylinker, enabling the gene of interest to be easily shuffled between systems. Additionally, a selection of competent cells is available for *in vitro* expression of difficult-to-express proteins.

Featured Products

- 228** PURExpress® *In Vitro* Protein Synthesis Kit
- 227** NEBExpress® Cell-free *E. coli* Protein Synthesis System
- 234** NEBExpress Ni Spin Columns
- 229** TEV Protease

Featured Tools & Resources

- 233** Purification Beads, Columns & Resins
- 350** Enhancing Efficiency
- 351** Protein Expression with T7 Express Strains



To learn more about NEB's portfolio of products for protein expression and purification, Visit www.neb.com/ProteinExpression to learn more.

Protein Expression Overview	226
-----------------------------	-----

EXPRESSION SYSTEMS

Cell-free Expression

NEBExpress Cell-free <i>E. coli</i> Protein Synthesis System	227
NEBExpress GamS Nuclease Inhibitor	227
PURExpress <i>In Vitro</i> Protein Synthesis Kit	228
PURExpress Δ Ribosome Kit	228
PURExpress Δ (aa, tRNA) Kit	228
PURExpress Δ RF123 Kit	228

PURExpress Companion Products:

PURExpress Disulfide Bond Enhancer	228
<i>E. coli</i> Ribosome	228

E. coli

NEBExpress MBP Fusion and Purification System	229
---	-----

NEBExpress Companion Products:

Amylose Resin	229
Anti-MBP Monoclonal Antibody	229
pMAL-c6T Vector	229
TEV Protease	229

IMPACT Kit	230
------------	-----

IMPACT Vector Comparison	230
--------------------------	-----

IMPACT Companion Products:

Chitin Resin	230
Anti-CBD Monoclonal Antibody	230
IMPACT Vectors	230

Yeast

<i>K. lactis</i> Protein Expression Kit	231
---	-----

K. lactis Companion Products:

SacII	231
Yeast Carbon Base Medium Powder	231
<i>K. lactis</i> GG799 Competent Cells	231

Competent Cells

Competent Cells for Protein Expression	232, 245–249
--	--------------

PURIFICATION

Purification Beads, Columns & Resins	
Selection Chart	233

Polyhistidine-tagged Protein Purification

NEBExpress Ni-NTA Magnetic Beads	234
NEBExpress Ni Resin	234
NEBExpress Ni Spin Column	234
TEV Protease	234

Maltose Binding Protein (MBP) Purification

Amylose Resin	234
Amylose Resin High Flow	234
Amylose Magnetic Beads	234
Anti-MBP Magnetic Beads	234
Anti-MBP Monoclonal Antibody	234

Chitin Binding Domain (CBD) Purification

Chitin Resin	234
Chitin Magnetic Beads	234
Anti-CBD Monoclonal Antibody	234

MAGNETIC BEAD PURIFICATION PRODUCTS

Oligo d(T) ₂₅ Magnetic Beads	235
Magnetic mRNA Isolation Kit	235
Streptavidin Magnetic Beads	235
Hydrophilic Streptavidin Magnetic Beads	235
Protein A Magnetic Beads	235
Protein G Magnetic Beads	235
Goat Anti-Mouse IgG Magnetic Beads	235
Goat Anti-Rabbit IgG Magnetic Beads	235
Goat Anti-Rat IgG Magnetic Beads	235
6-Tube Magnetic Separation Rack	235
50 ml Magnetic Separation Rack	235
12-Tube Magnetic Separation Rack	235
96-Well Microtiter Plate Magnetic Separation Rack	235
NEBNext Magnetic Separation Rack	235

Experts in protein expression and purification.

Protein expression can be a very complex, multi-factorial process. Each protein requires a specific intracellular environment to correctly achieve its secondary and tertiary structures. Proteins may also require post-translational modifications or insertion into a cellular membrane for proper function. Other proteins, once expressed, may be toxic to the host. Thus, no single solution exists for the successful production of all recombinant proteins, and a broad range of expression tools is needed to ensure the successful expression of your target protein.

Our NEBExpress portfolio of products includes solutions for expression and purification, and is supported by access to scientists with over 40 years of experience in developing and using recombinant protein technologies in *E. coli*. We use these solutions in our own research and manufacturing processes, and know that quality and performance are critical – all of our products are stringently tested so that you can be sure they will work optimally for your solution, just as we rely on them to work in ours.

Generate analytical amounts of protein with our cell-free expression systems

- *The PURExpress® In Vitro Protein Synthesis Kit* our novel cell-free transcription/translation system, enables protein expression in approximately two hours and is ideal for high throughput technologies
- *NEBExpress Cell-free E. coli Protein Synthesis System*, our novel *E. coli* cell-extract based transcription/translation system, is designed to synthesize proteins under the control of T7 RNA Polymerase at high yields for a wide variety of proteins

Generate and purify high yields of recombinant proteins


- *The NEBExpress MBP Fusion and Purification System* utilizes a pMAL vector and the *malE* gene for the expression of MBP-fusion proteins, which can be isolated by a two-step affinity purification
- *The IMPACT™ Kit* utilizes engineered protein splicing elements (inteins) to purify recombinant proteins in a single step
- *For analysis, try our protein standards, which are available unstained, prestained, or with two colors for easy identification*

Express a variety of proteins with our competent cells

- *Our popular BL21 and BL21(DE3) Competent E. coli strains* are available for routine expression
- *Lemo21(DE3) Competent E. coli* offers tunable T7 expression for difficult targets
- *For expression of His-tagged proteins, we offer NiCo21(DE3) Competent E. coli*
- *SHuffle® strains* are available for the expression of proteins with multiple disulfide bonds

Purify tagged proteins with our magnetic beads, columns and resins

- *Nickel spin columns, magnetic beads and resin* enable rapid purification of His-tagged proteins
- *Amylose resins* for purification of MBP-tagged proteins are available in a variety of formats (standard, high flow and magnetic)
- *Chitin resin* allows for rapid purification of CBD-tagged proteins
- *Remove affinity tags* following your purification with TEV Protease



NEW

NEBExpress® Cell-free *E. coli* Protein Synthesis System

#E5360S	10 reactions	\$210
#E5360L	100 reactions	\$1,700

Companion Product:

NEW

NEBExpress GamS Nuclease Inhibitor	
#P0774S	75 µg\$90


- Synthesize high yields of protein (typically 0.5 mg/ml)
- Protein can be synthesized and visualized in approximately 2–4 hours
- Synthesize target proteins ranging from 17 to 230 kDa
- Templates can be plasmid DNA, linear DNA, or mRNA
- RNase contamination can be inhibited by the supplied RNase inhibitor, eliminating clean-up steps
- Flexible reaction conditions achieve maximum yield; protein synthesis can be sustained for 10 hours at 37°C or up to 24 hours at lower temperatures
- Reactions can be miniaturized or scaled up to yield milligram quantities of protein

Description: The NEBExpress Cell-free *E. coli* Protein Synthesis System is a coupled transcription/translation system designed to synthesize proteins encoded by a DNA or mRNA template under the control of a T7 RNA Polymerase promoter. The system offers high expression levels, the ability to produce high molecular weight proteins, scalability, and is cost-effective for high throughput expression applications. The speed and robustness of the system facilitates protein synthesis in applications such as protein engineering, mutagenesis studies and enzyme screening.


The NEBExpress Cell-free *E. coli* Protein Synthesis System contains all the components required for protein synthesis, except for the target template DNA. It is a combination of a highly active cell extract from a genetically engineered strain, a reaction buffer, and an optimized T7 RNA Polymerase, which together yield robust expression of a wide variety of protein targets ranging from 17 to 230 kDa.

Applications

- Quickly generate analytical amounts of protein for further characterization
- High throughput screening and liquid handling
- Epitope mapping and protein folding
- Expression of toxic proteins



The NEBExpress Cell-free *E. coli* Protein Synthesis System can be used to express a wide range of proteins. 50 µl reactions containing 250 ng template DNA were incubated at 37°C for 3 hours. The red dot indicates the protein of interest. M = Unstained Protein Standard, Broad Range (NEB #P7717), "neg" = negative control, no DNA.



NEBExpress GamS Nuclease Inhibitor enhances synthesis of linear DNA. GamS inhibits Exonuclease V (RecBCD) activity and stabilizes linear DNA templates in *E. coli* based in vitro protein synthesis reactions. 50 µl reactions containing 100 ng linear template DNA, the components of the NEBExpress Cell-free *E. coli* Protein Synthesis System and 1.5 µg NEBExpress GamS Nuclease Inhibitor incubated for 5 hours at 37°C were monitored for activity as determined by fluorescence signal.



PURExpress® *In Vitro* Protein Synthesis Kits

PURExpress <i>In Vitro</i> Protein Synthesis Kit	
#E6800S 10 reactions \$261
#E6800L 100 reactions \$2,292
PURExpress Δ Ribosome Kit	
#E3313S 10 reactions \$363
PURExpress Δ (aa, tRNA) Kit	
#E6840S 10 reactions \$361
PURExpress Δ RF123 Kit	
#E6850S 10 reactions \$376

Companion Product:

PURExpress Disulfide Bond Enhancer	
#E6820S 50 reactions \$247
<i>E. coli</i> Ribosome	
#P0763S 1 mg \$186

- Generation of analytical amounts of proteins for further characterization
- Confirmation of open reading frames
- Generation of truncated proteins to identify active domains and functional residues
- Introduction of modified, unnatural or labeled amino acids (NEB #E6840, #E6850)
- tRNA structure and function studies (NEB #E6840)
- Ribosome structure and function studies (NEB #E3313, #P0763)
- Release factor function studies/ ribosome display (NEB #E6850)
- Epitope mapping


PURExpress is based on the PURE system technology originally developed by Dr. Takuya Ueda at the University of Tokyo and commercialized as the **PURESYSYSTEM™** by Biocomber (Tokyo, Japan).
PURESYSYSTEM™ is a trademark of Post Genome Institute.

Description: A rapid method for gene expression analysis, PURExpress is a novel cell-free transcription/translation system reconstituted from the purified components necessary for *E. coli* translation. The nuclease-free and protease-free nature of the PURExpress platform preserves the integrity of DNA and RNA templates/complexes and results in proteins that are free of modification and degradation. Transcription and translation are carried out in a one-step reaction, and require the mixing of only two tubes. With results available in a few hours, PURExpress saves valuable laboratory time and is ideal for high throughput technologies.

Advantages:

- Suitable for circular or linear DNA template
- Visualize synthesized protein directly on a Coomassie stained gel
- Protein expression in approximately 2 hours
- Transcription/translation components can be removed by affinity chromatography

PURExpress Disulfide Bond Enhancer: This proprietary blend of proteins and buffer components is designed to correctly fold target proteins with multiple disulfide bonds produced in PURExpress reactions or *E. coli* S30 extracts. Added at the beginning of a reaction, the components assist with the oxidation of cysteine thiols and correcting mis-oxidized substrates, increasing yield of soluble and functionally active protein.




Protein expression using the PURExpress *In Vitro* Protein Synthesis Kit. 25 µl reactions containing 250 ng template DNA and 20 units RNase Inhibitor were incubated at 37°C for 2 hours. 2.5 µl of each reaction was analyzed by SDS-PAGE using a 10–20% Tris-glycine gel. The red dot indicates the protein of interest. Marker M is the Protein Ladder.

PURExpress Δ Ribosome Kit: The ribosomes are removed from this kit, allowing users to add their own modified ribosomes for studies on protein translation. Control ribosomes (sufficient for 2 reactions) are provided in a separate tube.

PURExpress Δ RF123 Kit: Release factors are involved in termination of protein translation by recognizing the stop codons in an mRNA sequence. In a ribosome display experiment using PURExpress, lack of release factors could stabilize the ternary complex of mRNA-ribosome-nascent protein. As a result, the cDNA recovery could be higher. In this kit, the three release factors are supplied separately, allowing users to perform a protein synthesis reaction/ribosome display experiment with/without release factors of their choice.

PURExpress Δ (aa, tRNA) Kit: The tRNA and amino acids are supplied separately in this kit, allowing users to run a protein synthesis reaction by adding modified amino acids and tRNA mixtures to the reaction.

***E. coli* Ribosome:** The 70S *E. coli* Ribosome consists of a small subunit (30S) and a large subunit (50S). This preparation of ribosomes is highly active in NEB's PURExpress Protein Synthesis Kit (NEB #E6800), and can be used in ribosome structure and function studies, as a target for drug screening and as starting material for isolation of native ribosomal RNAs (5S, 16S, 23S). *E. coli* Ribosome is supplied as a 33.3 mg/ml solution.



PURExpress Disulfide Bond Enhancer. (PDBE) promotes proper folding of active vTPA. Reactions were set up according to PURExpress specifications with the vTPA template DNA. After a 2 hour incubation at 37°C, 5 µl of each reaction was used in an activity assay and cleavage of the chromogenic substrate was monitored for one hour.

PURExpress Kit Components

IN VITRO PROTEIN SYNTHESIS (E6800S)	Δ RIBOSOME (E3313S)	Δ RF 123 (E6850S)	Δ (aa, tRNA) (E6840S)
Solution A	Solution A	Solution A	Solution A (minus aa and tRNA)
Solution B	Factor Mix	Solution B (minus RF1, RF2 and RF3)	Solution B
Control (DHFR) template	Control (DHFR) template	Control (DHFR) template	Control (DHFR) template
	Control Ribosomes	RF1, RF2 and RF3	Amino Acid Mixture
			<i>E. coli</i> tRNA

NEW NEBExpress MBP Fusion and Purification System

#E8201S\$691

Components Sold Separately:

Amylose Resin
#E8021S 15 ml\$233
#E8021L 100 ml\$1,216

Anti-MBP Monoclonal Antibody
#E8032S 0.05 ml\$181
#E8032L 0.25 ml\$727

NEW pMAL-c6T Vector
#N0378S 10 µg\$141

TEV Protease
#P8112S 1,000 units\$105

Description: In the NEBExpress MBP Fusion and Purification System, the pMAL-c6T vector provides a method for expressing and purifying a protein produced from a cloned gene or open reading frame. The cloned gene is inserted downstream from and in frame with the *malE* gene of *E. coli*, which encodes maltose-binding protein (MBP); this construct results in the expression of an MBP fusion protein. The pMAL-c6T vector expresses the *N*-terminal hexahistidine tagged *malE* gene (lacking its secretory signal sequence and engineered for tighter binding to amylose) followed by a multiple cloning site containing a TEV protease recognition sequence and stop codons in all three frames. The pMAL-c6T vector expresses the MBP fusion in the cytoplasm. The method uses the strong “*lac*” promoter and the *malE* translation initiation signals to yield high-level expression of the cloned sequences. The fusion protein is then purified by a one-step purification method using amylose resin and MBP’s affinity for maltose.

Following amylose purification, the target protein can be cleaved from the MBP-tag using TEV Protease, without adding any vector-derived residues to the protein. Both

the MBP-tag and TEV Protease are polyhistidine-tagged for easy removal from the reaction. Loading the digest onto NEBExpress Ni Resin (NEB #S1428) sequesters both the MBP-tag and TEV Protease, thereby isolating the target protein in the column flow through. The target protein yield can be up to 100 mg/L, with typical yields in the range of 10–40 mg/L.

References: References for properties and applications of this product can be found at www.neb.com.

The NEBExpress MBP Fusion and Purification System includes:

- pMAL-c6T Vector
- Amylose Resin
- TEV Protease, TEV Protease Reaction Buffer
- Anti-MBP Monoclonal Antibody
- MBP6 Protein
- MBP6-TEV-paramyosin-ΔSal
- *E. coli* ER2523 (NEB Express) (Glycerol Stock)

- *Reliable E. coli* expression: substantial yields (up to 100 mg/L)
- Fusion to MBP significantly enhances proper folding of target proteins
- Two-step purification: amylose elution followed by TEV Protease cleavage and Ni resin isolation results in a highly pure tag-free target protein
- Gentle elution with maltose: no detergents or harsh denaturants required




Figure 1: Schematic illustration of the NEBExpress MBP Fusion and Purification System

IMPACT™ Kit

#E6901S\$427

Components Sold Separately:

pTXB1 Vector
#N6707S 10 µg\$136

pTYB21 Vector
#N6709S 10 µg\$136

Chitin Resin
#S6651S 20 ml\$114
#S6651L 100 ml\$459

Anti-CBD Monoclonal Antibody
#E8034S 0.05 ml\$73

pTWIN1 Vector
#N6951S 10 µg\$129
(not included with the kit)

- Single-column purification without the use of proteases
- Produce target protein without vector-derived amino acids
- Fusion to either N- or C-terminus of target protein
- Ligation and labeling of recombinant proteins
- Isolation of proteins with or without N-terminal methionine

For a more detailed description and a restriction map of the pTXB1 and pTYB21 vectors, including the MCS, see page 363–364; visit www.neb.com for sequence files.

Description: The IMPACT (Intein Mediated Purification with an Affinity Chitin-binding Tag) Kit utilizes engineered protein splicing elements (inteins) to purify recombinant proteins by a single column (Figure 1). This kit distinguishes itself from other protein fusion systems by its ability to separate a recombinant protein from the affinity tag without the use of a protease.

The IMPACT Kit allows fusion of a tag consisting of an intein and the chitin binding domain (CBD) from *Bacillus circulans*, to either the C-terminus (pTXB1) or the N-terminus (pTYB21) of a target protein (Figure 2). In the presence of thiols, such as DTT, the intein undergoes specific self-cleavage which releases the target protein. The pTXB1 vector can also be used to express a protein with a C-terminal thioester for use in Intein-mediated Protein Ligation (IPL). The IPL reaction, also referred to as expressed protein ligation, allows for ligation of a peptide or a protein with a N-terminal cysteine to a bacterially expressed protein with a C-terminal thioester through a native peptide bond for use in protein labeling and semisynthesis. For more information on the IMPACT System and IPL, visit www.neb.com.

pTXB1 is a *E. coli* expression vector that utilizes a mini-intein from the *Mycobacterium xenopi gyrA* gene [*Mxe* GyrA intein; 22 kDa]. This intein has been modified and combined with the CBD to create an affinity tag which can be bound to chitin beads (NEB #S6651). Release of the target protein is induced by thiol-reagents such as DTT or 2-mercaptoethanesulfonic acid (for ligation).

The pTYB21 vector allows for the fusion of the intein tag containing the *Saccharomyces cerevisiae* (*Sce*) VMA1 intein and CBD to the N-terminus of the target protein.

pTWIN1 Vector is available separately and enables isolation of proteins with an N-terminal cysteine and/or a C-terminal thioester. The polylinker is designed for the in-frame fusion of a target gene between the modified *Ssp* DnaB and *Mxe* GyrA inteins. The presence of the CBD facilitates purification.

References: References for properties and applications of this product can be found at www.neb.com.

The IMPACT Kit Includes:

- *E. coli* expression vectors: pTXB1, pTYB21 and pMXB10 Control Vector
- Chitin Beads: Affinity matrix (binding capacity=2 mg/ml)
- Anti-CBD Monoclonal Antibody
- 1, 4-Dithiothreitol (DTT) & SDS Loading Buffer




Figure 1: Purification of Maltoe Binding Protein (MBP) in a single affinity purification step: Lane 1: uninduced cell extract. Lane 2: induced cell extract showing expressed fusion protein. Lane 3: MBP fraction eluted after inducing cleavage overnight at 4°C. Lane 4: MBP ligated to a peptide containing an N-terminal cysteine. Marker M is the protein ladder.




Figure 2: Schematic of the IMPACT System.

Table 1: Guide to IMPACT vectors and applications.

VECTORS	SITE OF TARGET PROTEIN FUSION	INTEIN TAG (KDA)	RECOMMENDED CLONING SITES ^a	PREFERRED RESIDUES AT CLEAVAGE SITE ^b	METHOD OF CLEAVAGE ^{c, d}	APPLICATIONS
pTXB1	C-terminus	Mxe GyrA intein (28)	NdeI-SapI/Spel	Y, F, O, N, T, K, A, H, M (Unfavorable residues: S, P, D, G)	DTT (or MESNA) pH 8.0–8.5, 4°C	Purification; C-terminal thioester for ligation and modification
pTYB21	N-terminus	Sce VMA1 intein (56)	SapI/BsmI/NdeI-PstI	A, Q, M, G, L, N, W, F, Y (Unfavorable residues: P, S, C, T, R)	DTT ^d pH 8.0–8.5, 25°C	Purification
pTWIN1	C-terminus (Intein 2)	Mxe GyrA intein (28)	NdeI-SapI/Spel	M, Y, F, LEM (Unfavorable residues: S, P, E, D)	DTT (or MESNA) pH 8.0–8.5, 4°C	Purification; C-terminal thioester for ligation and modification

^aNEBuilder HiFi DNA Assembly Cloning Kit (NEB #E5520) can be used to generate construct without the use of restriction enzymes. ^bActual cleavage efficiency is dependent on the adjacent residues as well as the folding of the fusion protein. ^cDithiothreitol (DTT) is used only for protein purification. 2-mercaptoethanesulfonic acid (MESNA) is used for isolation of proteins possessing a C-terminal thioester for ligation, labeling and cyclization.

^dCysteine can be used in the place of DTT.

K. lactis Protein Expression Kit

#E1000S \$935

Components Sold Separately:

SacII
 #R0157S 2,000 units \$65
 #R0157L 10,000 units \$266
 Yeast Carbon Base Medium Powder
 #B9017S 12 g \$46
 K. lactis GG799 Competent Cells
 #C1001S 5 transformation reactions \$380

- Clone and express genes toxic to *E. coli*
- Simultaneous expression of multiple genes
- No expensive antibiotics or methanol required
- Easy-to-use protocols for those inexperienced with yeast systems
- Attractive commercial sublicensing

Restriction map for pKLAC2 can be found on page 358; sequence files are available at www.neb.com.

Description: The *K. lactis* Protein Expression Kit provides an easy method for expressing a gene of interest in the yeast *Kluyveromyces lactis*. The gene is cloned into the integrative pKLAC-series of vectors, and may be expressed intracellularly or secreted. The *K. lactis* system offers several advantages over other yeast and bacterial expression systems. Abundant overexpression of protein is achieved through high culture densities as well as the ability to integrate multiple copies of the vector. The pKLAC-series of vectors use a strong *LAC4* promoter, which has been modified to lack expression in *E. coli*, making this system useful for expressing toxic genes. For the selection of transformants, no expensive antibiotics are required. In addition, no methanol is required in growth media. Finally, the *K. lactis* system can express post-translationally modified proteins, making it a useful alternative to bacterial expression systems.

pKLAC2 is a general purpose expression vector. Using this vector, proteins may be produced intracellularly or may be fused to the *K. lactis* α -mating factor sequence for secreted expression. Vector pKLAC2 contains an MCS that is compatible with other expression systems available from NEB.


GG799 competent cells are provided in the *K. lactis* Protein Expression Kit. GG799 cells are characterized by very high cell density growth and efficient expression of foreign proteins. GG799 cells have no genetic modifications.

The *K. lactis* Protein Expression Kit Includes:


- pKLAC2 Vector and pKLAC1-*malE* Control Plasmid
- SacII
- Integration Primer Set
- rCutSmart™ Buffer
- *K. lactis* GG799 Competent Cells and Transformation Reagent
- Yeast Medium Powder & Acetamide Solution

Note: This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

References: References for properties and applications of these products can be found at www.neb.com.



Secreted protein processing. In the nucleus, an integrated expression vector encoding a fusion between the α -MF domain (blue) and a desired protein (black) is expressed. A signal peptide in the α -MF domain directs entry of the fusion protein into the endoplasmic reticulum (ER) and is removed by signal peptidase (SP). The fusion protein is transported to the Golgi apparatus where the Kex protease removes the α -MF domain. The protein of interest is then secreted from the cell.



Protein Expression in *K. lactis*. SDS-polyacrylamide gel electrophoresis separation of secreted recombinant maltose-binding protein (MBP) detected directly in peptone rich growth medium by Coomassie staining. Lane 1: Protein Molecular Weight Markers. Lane 2: spent culture medium (15 µl) from wild-type *K. lactis* cells. Lane 3: spent culture medium (15 µl) from *K. lactis* cells harboring an integrated expression cassette containing the *E. coli malE* gene.

Companion Products:

Enterokinase, light chain
 #P8070S 480 units \$122
 #P8070L 2,560 units \$487
 BstXI
 #R0113S 1,000 units \$70
 #R0113L 5,000 units \$287

OneTaq DNA Polymerase
 #M0480S 200 units \$44
 #M0480L 1,000 units \$174
 #M0480X 5,000 units \$684

OneTaq Hot Start DNA Polymerase
 #M0481S 200 units \$91
 #M0481L 1,000 units \$362
 #M0481X 5,000 units \$1,444

Deoxynucleotide (dNTP) Solution Mix
 #N0447S 8 µmol of each* \$64
 #N0447L 40 µmol of each \$255

*Available in 4 x 0.2 ml aliquots

Competent Cells for Protein Expression

- Free of animal products
- T1 phage resistance (*thiA2*)
- B strains are deficient in proteases *Lon* and *OmpT*
- A variety of convenient formats
- Bulk sales capabilities with custom packaging

NEB offers a wide selection of competent cell strains designed for the expression of a variety of proteins. Proteins with multiple disulfide bonds are correctly oxidized to significantly higher yields with SHuffle® strains. Tunable T7 expression is achieved with Lemo21(DE3), an ideal strain for difficult targets including membrane proteins. NiCo21(DE3) can be used for expression and purification of His-tagged proteins. NEB Express and T7 Express are

offered with varying levels of expression control. Several strains are available with added control by *lacI^q*. Only NEB offers exceptional control of T7 expression from the *lysY* gene, which is ideal for proteins that are difficult to express or toxic to the cell. Each strain is provided with a protocol for optimal expression.

For more information see pages 245–249.

Expression Strains

CHARACTERISTICS	STRAIN	NEB #	SIZE
<ul style="list-style-type: none"> • Versatile non-T7 expression strain • Protease deficient 	NEBExpress Competent <i>E. coli</i> *	C2523H/I	20 x 0.05 ml/6 x 0.2 ml
<ul style="list-style-type: none"> • Control of IPTG induced expression from <i>P_{lac}</i>, <i>P_{lac}</i> and <i>P_{trc}</i> • Protease deficient 	NEBExpress <i>I^q</i> Competent <i>E. coli</i>	C3037I	6 x 0.2 ml
<ul style="list-style-type: none"> • Most popular T7 expression strain • Protease deficient 	T7 Express Competent <i>E. coli</i>	C2566H/I	20 x 0.05 ml/6 x 0.2 ml
<ul style="list-style-type: none"> • T7 expression • Protease deficient • Better reduction of basal expression 	T7 Express <i>lysY</i> Competent <i>E. coli</i>	C3010I	6 x 0.2 ml
<ul style="list-style-type: none"> • T7 expression • Protease deficient • Highest level of expression control 	T7 Express <i>lysY/I^q</i> Competent <i>E. coli</i>	C3013I	6 x 0.2 ml
<ul style="list-style-type: none"> • Enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm • Protease deficient/B strain 	SHuffle® Express Competent <i>E. coli</i>	C3028J	12 x 0.05 ml
<ul style="list-style-type: none"> • Enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm • T7 expression • Protease deficient/B strain 	SHuffle T7 Express Competent <i>E. coli</i>	C3029J	12 x 0.05 ml
<ul style="list-style-type: none"> • T7 expression • Protease deficient/B strain • Tightly controlled expression of toxic proteins • Enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm 	SHuffle T7 Express <i>lysY</i> Competent <i>E. coli</i>	C3030J	12 x 0.05 ml
<ul style="list-style-type: none"> • Enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm • T7 expression/K12 strain 	SHuffle T7 Competent <i>E. coli</i>	C3026J	12 x 0.05 ml
<ul style="list-style-type: none"> • Routine expression for non-T7 Vectors • Protease deficient 	BL21 Competent <i>E. coli</i>	C2530H	20 x 0.05 ml
<ul style="list-style-type: none"> • Routine T7 Expression • Protease deficient 	BL21(DE3) Competent <i>E. coli</i>	C2527H/I	20 x 0.05 ml/6 x 0.2 ml
<ul style="list-style-type: none"> • Tunable T7 Expression for difficult targets • Protease deficient 	Lemo21(DE3) Competent <i>E. coli</i>	C2528J	12 x 0.05 ml
<ul style="list-style-type: none"> • Expression and purification of His-tagged proteins • Protease deficient 	NiCo21(DE3) Competent <i>E. coli</i>	C2529H	20 x 0.05 ml

Note: Store Competent Cells at –80°C. Once thawed, do not refreeze. Storage at –20°C will result in a significant decrease in transformation efficiency. Cells lose efficiency whenever they are warmed above –80°C, even if they do not thaw.

* NEB Express is the recommended strain for the NEBExpress MBP Fusion and Purification System.

Purification Beads, Columns and Resins

Isolation of pure substrates or proteins for downstream experiments is a common, yet time consuming, task. New England Biolabs offers a variety of resins and magnetic beads that are easy-to-use, highly specific, and available in several different formats for rapid isolation and purification of proteins, nucleic acids and immunoglobulins. NEB's magnetic beads are ideally suited for applications involving high-throughput proteomic screening, small-scale protein isolation, immunomagnetic isolations or cell separation experiments. With magnetic beads, affinity

purification of tagged proteins, antigens, antibodies and nucleic acids can be done conveniently and quickly. Immobilized substrates remain biologically active and can be eluted in small volumes or serve as ligands in subsequent pull-down or target interaction experiments involving DNA or proteins. NEB's resins enable simple, one-step purification strategies for tagged proteins, and result in a high yield of highly pure substrate. For the full list of products available for protein expression and purification, visit www.neb.com/ProteinExpression.

Product Selection Chart

	PROTEIN PURIFICATION	LARGE-SCALE PURIFICATIONS	USE IN AUTOMATED CHROMATOGRAPHY	HIGH-THROUGHPUT	BIOTINYLATED SUBSTRATE BINDING	PROTEIN PULL-DOWN	NUCLEIC ACID PULL-DOWN	mRNA PURIFICATION/PULL-DOWN	IMMUNOPRECIPITATION	CELL SEPARATION/CELL SORTING
NEW NEBExpress Ni-NTA Magnetic Beads (NEB #S1423)	● (His-tag)			●		●				
NEW NEBExpress Ni Spin Columns (NEB #S1427)	● (His-tag)			●		●				
NEW NEBExpress Ni Resin (NEB #S1428)	● (His-tag)	●	●			●				
Amylose Resin (NEB #E8021)	● (MBP)	●				●				
Amylose Resin High Flow (NEB #E8022)	● (MBP)	●	●			●				
Amylose Magnetic Beads (NEB #E8035)	● (MBP)			●		●				
Anti-MBP Magnetic Beads (NEB #E8037)	● (MBP)			●		●				
Chitin Resin (NEB #S6651)	● (intein-CBD tag)	●				●				
Chitin Magnetic Beads (NEB #E8036)	● (intein-CBD tag)			●		●				
Oligo d(T) ₂₅ Magnetic Beads (NEB #S1419)				●			●	●		
Streptavidin Magnetic Beads (NEB #S1420)				●	●	● (biotinylated bait)	● (biotinylated bait)			
Hydrophilic Streptavidin Magnetic Beads (NEB #S1421)				●	●	● (biotinylated bait)	● (biotinylated bait)			
Protein A Magnetic Beads (NEB #S1425)				●					●	
Protein G Magnetic Beads (NEB #S1430)				●					●	
Goat Anti-Mouse IgG Magnetic Beads (NEB #S1431)				●					● (Mouse IgGs)	●
Goat Anti-Rabbit IgG Magnetic Beads (NEB #S1432)				●					● (Rabbit IgGs)	●
Goat Anti-Rat IgG Magnetic Beads (NEB #S1433)				●					● (Rat IgGs)	●
Magnetic mRNA Isolation Kit (NEB #S1550)				●				●		

Polyhistidine-tagged Protein Purification

NEW

NEBExpress Ni-NTA Magnetic Beads		
#S1423S	1 ml\$190
#S1423L	5 ml\$899

NEW

NEBExpress Ni Resin		
#S1428S	25 ml\$317

NEW

NEBExpress Ni Spin Columns		
#S1427S	10 columns\$102
#S1427L	25 columns\$204

NEW

TEV Protease		
#P8112S	1,000 units\$105

NEBExpress Ni-NTA Magnetic Beads: An affinity matrix for the small-scale isolation and purification of polyhistidine-tagged (His-tagged) fusion proteins in manual or automated formats. High specific binding yields purities of > 95% in a single-purification step. Matrix tolerates a wide range of conditions, including the presence of denaturants and detergents. Compatible with commercially available detergent-based cell lysis reagents. Elution can be achieved by protonation, ligand exchange (with imidazole) or extraction of the metal ion by a strong chelator (e.g., EDTA).

Support Matrix: Spherical, agarose based superparamagnetic microparticles ranging in size from 20-100 μ m.

Binding Capacity: Varies with target, typically ≥ 7.5 mg His-tagged fusion protein/ml bed volume.

NEBExpress Ni Resin: NEBExpress Ni Resin is an affinity matrix for the isolation and purification of polyhistidine-tagged (His-tagged) fusion proteins. It is intended for use in gravity or pressure flow columns and batch purifications, and high specific binding yields purities of > 95% in a single-purification step. NEBExpress Ni Resin is comprised of a highly uniform and chemical-tolerant resin that is pre-charged with nickel ions on the matrix surface. It is resistant to a wide range of chemicals, including NaOH, EDTA, and commonly used reducing agents such as TCEP, DTT, and β -mercaptoethanol. Can be used under native or denaturing conditions.

Support Matrix: Spherical, agarose based microparticles ranging in size from 10-100 μ m.

Binding Capacity: 1 ml of NEBExpress Ni Resin will bind ≥ 10 mg of His-tagged fusion protein.

NEBExpress Ni Spin Columns: NEBExpress Ni Spin columns are pre-packed with agarose-based microparticles ranging in size from 10-100 μ m for the small-scale isolation and purification of polyhistidine-tagged (His-tagged) fusion proteins. Purification can be performed under native or denaturing conditions, including conditions in which EDTA or reducing reagents are required, yielding highly pure target protein in a single purification step. This enables screening of expression conditions and streamlines the functional and structural characterization of the target protein.

Support Matrix: Spherical, agarose based microparticles ranging in size from 10-100 μ m.

Binding Capacity: Varies with target, ≥ 1 mg His-tagged fusion protein per column.

TEV Protease: TEV Protease, also known as Tobacco Etch Virus (TEV) Protease, is a highly specific cysteine protease that recognizes the amino-acid sequence Glu-Asn-Leu-Tyr-Phe-Gln-(Gly/Ser) and cleaves between the Gln and Gly/Ser residues. It is often used for the removal of affinity purification tags such as maltose-binding protein (MBP) or poly-histidine from fusion proteins. TEV Protease has a 7XHis-tag for easy removal from a reaction using nickel affinity resins and has been engineered to improve thermal stability and decrease autolysis.

Maltose Binding Protein (MBP) Purification

Amylose Resin		
#E8021S	15 ml\$233
#E8021L	100 ml\$1,216

Amylose Resin High Flow		
#E8022S	15 ml\$313
#E8022L	100 ml\$1,817

Amylose Magnetic Beads		
#E8035S	25 mg\$272

Anti-MBP Magnetic Beads		
#E8037S	1 ml\$275

Anti-MBP Monoclonal Antibody		
#E8032S	0.05 ml\$181
#E8032L	0.25 ml\$727

Amylose Resin and Amylose Resin High Flow:

Affinity matrix used for isolation of proteins fused to maltose-binding protein. It is a composite amylose/agarose bead. Amylose Resin High Flow is a more rigid bead, suitable for use in automated chromatography systems.

Binding Capacity: Amylose Resin and Amylose Resin High Flow: > 4 mg MBP5-paramyosin Δ Sal fusion protein/ml of bed volume.

Amylose Magnetic Beads and Anti-MBP Magnetic Beads:

Affinity matrices for the small-scale isolation and purification of MBP-fusion proteins. Amylose or monoclonal Anti-MBP are covalently coupled to a paramagnetic particle through a linkage that is stable and leak resistant over a wide pH range.

Support Matrix: Amylose Magnetic Beads – 10 μ M superparamagnetic particles; Anti-MBP Magnetic Beads – 1 μ M nonporous super paramagnetic particle.

Binding Capacity: 1 mg of Amylose Magnetic Beads will bind ≥ 10 μ g of MBP-fusion protein. 1 mg of Anti-MBP Magnetic Beads will bind 5 μ g of MBP-paramyosin Δ Sal fusion protein.

Anti-MBP Monoclonal Antibody: Anti-MBP Monoclonal Antibody is a murine anti-maltose-binding protein antibody, isotype IgG2a. It is purified from tissue culture supernatant by protein A affinity chromatography

Chitin Binding Domain (CBD) Purification

Chitin Resin		
#S6651S	20 ml\$114
#S6651L	100 ml\$459

Chitin Magnetic Beads		
#E8036S	5 ml\$129

Anti-CBD Monoclonal Antibody		
#E8034S	0.05 ml\$73

Chitin Resin: An affinity matrix for the isolation of target proteins fused on an intein-chitin binding domain (CBD). Strong specific binding enables purification of highly pure protein from crude lysates in one step. Removal of CBD-tag during elution typically yields highly pure, native protein without the use of a protease.

Support Matrix: Approximately 50-70 μ m microparticles

Binding Capacity: 2.0 mg maltose-binding protein/ ml bed volume released from the resin after cleavage of the fusion protein expressed from pMYB5.

Chitin Magnetic Beads: An affinity matrix for the small-scale isolation of target proteins fused to a chitin binding domain (CBD). Chitin beads have been prepared

with encapsulated magnetite, thereby permitting the magnetic isolation of CBD-fusion proteins from cell culture supernatants. Removal of CBD-tag during elution typically yields highly pure, native protein.

Support Matrix: Approximately 50-70 μ m paramagnetic microparticles

Binding Capacity: 2 mg chitin binding domain protein / ml bed volume released

Anti-CBD Monoclonal Antibody: Anti-CBD Monoclonal Antibody is a murine anti-chitin binding domain (CBD) antibody, isotype IgG1. It has high purity and specificity for chitin binding domain tag, and is verified for use in both Western blotting and ELISA.

Magnetic Bead Purification Products

Oligo d(T)₂₅ Magnetic Beads		
#S1419S	5 ml\$307
Magnetic mRNA Isolation Kit		
#S1550S	25 isolations\$334
Streptavidin Magnetic Beads		
#S1420S	5 ml\$334
Hydrophilic Streptavidin Magnetic Beads		
#S1421S	5 ml\$334
Protein A Magnetic Beads		
#S1425S	1 ml\$184
Protein G Magnetic Beads		
#S1430S	1 ml\$207
Goat Anti-Mouse IgG Magnetic Beads		
#S1431S	1 ml\$192
Goat Anti-Rabbit IgG Magnetic Beads		
#S1432S	1 ml\$192
Goat Anti-Rat IgG Magnetic Beads		
#S1433S	1 ml\$190

- *Small-scale purification or immunoprecipitation of IgG species*
- *No centrifugation required*
- *Regenerate matrix without binding capacity loss*

Oligo d(T)₂₅ Magnetic Beads: These beads enable small-scale isolations of mRNA from a variety of samples, including *in vitro* transcribed mRNA, total RNA, crude cell lysates and tissue. The selectivity for mRNA results from the annealing of bead-linked oligo d(T)₂₅ to the poly(A) region present in most eukaryotic mRNAs.

Support Matrix: 1 µm nonporous superparamagnetic microparticles

Binding Capacity: ≥ 5 µg rA₃₀ per mg of beads

Magnetic mRNA Isolation Kit: The Magnetic mRNA Isolation Kit is designed to isolate intact poly(A)⁺ RNA from cells and tissue without requiring phenol or other organic solvents. The technology is based on the coupling of Oligo d(T)₂₅ to 1 µm paramagnetic beads, which is then used as the solid support for the direct binding of poly(A)⁺ RNA.

Streptavidin Magnetic Beads: The beads provide fast magnetic response times and reaction kinetics, and they have high binding capacity and sensitivity while retaining their physical integrity. They can be used to capture biotin-labeled substrates including DNA, RNA, peptides, antigens, antibodies and other proteins of interest in manual or automated workflows. These beads typically exhibit lower non-specific binding of proteins.

Support Matrix: 1 µm nonporous superparamagnetic microparticles

Binding Capacity: ≥ 30 µg biotinylated antibody per mg of beads or > 500 pmol of single-stranded 25 bp biotinylated oligonucleotide per mg of beads

Hydrophilic Streptavidin Magnetic Beads: The beads provide rapid magnetic response times and reaction kinetics, and they have high binding capacity and sensitivity while retaining their physical integrity. They can be used to capture biotin-labeled substrates including DNA, RNA, peptides, antigens, antibodies and other proteins of interest in manual or automated workflows. These beads typically exhibit lower non-specific binding of nucleic acids.

Support Matrix: 2 µm non-porous superparamagnetic microparticles

Binding Capacity: > 400 pmol of single-stranded 25 bp biotinylated oligonucleotide per mg of beads

Protein A and Protein G Magnetic Beads: The beads allow for isolation of most mammalian immunoglobulins (IgGs) and are amenable to immunoprecipitation. Predominant Fc-binding allows optimal IgG orientation upon binding to the outer surface of the Protein A and Protein G Magnetic Beads allowing Fab regions to efficiently bind antigen. These beads can be used to immunoprecipitate target proteins from crude cell lysates using a selected primary antibody. In addition, specific antibodies can be chemically cross-linked to the Protein A- or Protein G- coated surface to create a reusable immunoprecipitation bead, thereby avoiding the co-elution of antibody with the target antigen.

Support Matrix: 2 µm nonporous superparamagnetic microparticles

Binding Capacity: > 280 µg of Human IgG per ml of beads

Goat Anti-Mouse IgG Magnetic Beads: An affinity matrix for the small-scale immunomagnetic separation and purification of mouse IgG. Specifically, the beads consist of Anti-Mouse IgG that is covalently coupled to a 1 µm nonporous superparamagnetic particle.

Support Matrix: 1 µm nonporous superparamagnetic microparticle

Binding Capacity: 5 µg mouse IgG/mg Goat Anti-Mouse IgG Magnetic Beads

Goat Anti-Rabbit IgG Magnetic Beads: An affinity matrix for the small-scale immunomagnetic separation and purification of rabbit IgG. Specifically, the beads consist of Goat Anti-Rabbit IgG that is covalently coupled to a 1 µm nonporous superparamagnetic particle.

Support Matrix: 1 µm nonporous superparamagnetic microparticle

Binding Capacity: 5 µg rabbit IgG/mg Goat Anti-Rabbit IgG Magnetic Beads

Goat Anti-Rat IgG Magnetic Beads: An affinity matrix for the small-scale immunomagnetic separation and purification of rat IgG. Specifically, the beads consist of anti-Rat IgG that is covalently coupled to a 1 µm nonporous superparamagnetic particle.

Support Matrix: 1 µm nonporous superparamagnetic microparticle

Binding Capacity: 5 µg rat IgG/mg Goat Anti-Rat IgG Magnetic Beads

Magnetic Separation Racks

	APPLICATION	MAGNETS	CAPACITY	CONVENIENCE	PRICE
6-Tube Magnetic Separation Rack (NEB #S1506)	Designed for small-scale separations using magnetic particles	Neodymium rare earth permanent magnets	6 tubes (1.5 ml)	Use with magnetic particle-based affinity purification for rapid, small-scale purifications	\$215
50 ml Magnetic Separation Rack (NEB #S1507)	Designed for small-scale separations using magnetic particles	Neodymium rare earth permanent magnets	4 tubes (50 ml)	Use with magnetic particle-based affinity purification for rapid, streamlined purifications	\$304
12-Tube Magnetic Separation Rack (NEB #S1509)	Designed for small-scale separations using magnetic particles	Neodymium rare earth permanent magnets	12 tubes (1.5 ml)	Use with magnetic particle-based affinity purification for rapid, small-scale purifications	\$340
96-Well Microtiter Plate Magnetic Separation Rack (NEB #S1511)	Designed for use with commercially available high-flanged 100 µl to 300 µl flat-bottom 96-well microplates	24 side-pull magnetic pins attract magnetic beads from solution to the side walls of four adjacent wells	96-well	The orientation of the magnetic field ensures complete removal of the magnetic beads from solution during pipetting steps, thereby minimizing sample loss	\$509
NEBNext® Magnetic Separation Rack (NEB #S1515)	Designed for rapid and effective small-scale separations of magnetic particles	Anodized aluminum rack with Neodymium Iron Boron (NdFeB) rare earth magnets	24 tubes (0.2 ml)	Next generation sequencing library preparation workflows include magnetic bead-based purification and size-selection steps. It is important for library yield and quality that bead separation be highly efficient and fast, and this is enabled by the powerful fixed magnet cores in this rack.	\$475



Sustainable Growth

What does “sustainable growth” mean? Originally conceptualized by Robert C. Higgins, sustainable growth rate refers to the growth rate that can be supported by the financial means of the entity itself, without additional equity or debt. Don Comb was a firm believer in the principles of sustainable growth. Throughout its history, New England Biolabs has grown steadily and sustainably, at times purposefully tempering growth to maintain the company for future generations. And over 45 years later, employees remain confident that the company will be there to support them for years to come — as well as its customers and the scientific community.

Sustainable growth applies to more than just companies and institutions though, it also applies to the populations. In 1972, scientists and economists at MIT published *Limits to Growth*, which outlined that the fate of humanity depends on the preservation of biodiversity and the slowing of economic and population growth. Economic growth means greater industrial output: more consumer goods, factories, oil wells, mines and agricultural clearing — and thus less protection of natural resources and ecosystems. While striving for greater economic growth may lead to wealth, it does not necessarily lead to a better world.

So how do we bring the world's populations onto a more sustainable path? The first step is to curb wasteful resource consumption and reduce the disparity in living standards between nations. NEB has supported these efforts through contributions to non-profit organizations dedicated to improving lives by helping communities protect, conserve and manage the natural resources upon which their long-term well-being depends. One such organization is Trees, Water & People, which works with communities across Central America to practice watershed protection, sustainable agriculture and reforestation, as well as to distribute clean cookstoves — all in an effort to protect natural resources.

In addition to supporting efforts like Trees, Water & People, NEB houses the Creative Action Institute (CAI), whose mission is to “catalyze community-driven solutions that advance gender equality and build a sustainable planet.” CAI is working around the world to build the capacity of leaders and organizations for innovation, collaboration and resilience to advance conservation, health and human rights.

NEB has always embraced a “small is beautiful” approach, continuing to make scientific and technological advances — while growing sustainably. Additionally, NEB operates its manufacturing facilities and headquarters efficiently, minimizes waste, subsidizes energy usage with solar panels, and recycles aggressively to protect natural resources. But no company can save the world alone, it will take many companies and individuals adopting similar practices. With the combined efforts of many to protect and preserve natural resources, and to slow the growth rate, we can all enjoy this planet for centuries to come.

Sprout spruce seedling emerging
from the ground.
Credit: Олег Копьев, Adobe Stock

Learn more about
sustainability at
NEB.



Competent Cells

NEB® has a competent cell strain for your needs.

Choose the right cells for your cloning and protein expression applications from NEB's portfolio of high efficiency competent cell strains.

For cloning experiments, choose from a variety of formats, including chemical and electrocompetent. These *E. coli* strains are T1 phage resistant and are Endonuclease-deficient for high-quality plasmid preparations. Additionally, all competent cells from NEB are free of animal products.

NEB also offers a wide variety of competent cell strains ideal for many protein expression applications. These strains address the needs of difficult protein expression control, toxic protein expression and cytoplasmic disulfide bond formation. NEBExpress®, T7 Express and SHuffle® strains are available with varying levels of control. *lacI^q* strains feature added control from increased supply of Lac repressor (*lacI^q*). Only NEB offers the exceptional control of expression from the *lysY* gene that reduces basal expression from T7 strains without inhibiting IPTG-induced expression. Our Lemo21(DE3) strain features tunable T7 expression for difficult targets. Each strain is provided with a detailed protocol for optimal expression.



NO DRY ICE CHARGES
with Competent Cells from NEB



Find tips
for successful
transformation.

Featured Products


- 242** NEB Cloning Competent *E. coli* Sampler
- 244** NEB Stable Competent *E. coli*
- 245** BL21 Competent *E. coli*
- 248** SHuffle Express Competent *E. coli*

Featured Tools & Resources

- 350** Enhancing Transformation Efficiency
- 351** Protein Expression with T7 Express Strains
- 336** Troubleshooting Guide for Cloning



Visit www.neb.com to find additional online tools, including our Competitor Cross-reference Tool for comparing NEB strains to other commercially available strains.



Strain Properties Chart	240	Protein Expression Strains	
Cloning Strains		BL21 Competent <i>E. coli</i>	245
NEB Cloning Competent <i>E. coli</i>		BL21(DE3) Competent <i>E. coli</i>	245
Sampler	242	Lemo21(DE3) Competent <i>E. coli</i>	245
NEB Turbo Competent <i>E. coli</i>		NiCo21(DE3) Competent <i>E. coli</i>	246
(High Efficiency)	242	NEBExpress Competent <i>E. coli</i>	
NEB 10-beta Competent <i>E. coli</i>		(High Efficiency)	246
(High Efficiency)	243	NEBExpress <i>I^q</i> Competent <i>E. coli</i>	
Electrocompetent	243	(High Efficiency)	246
NEB 5-alpha Competent <i>E. coli</i>		T7 Express Competent <i>E. coli</i>	
(High Efficiency)	243	(High Efficiency)	247
(Subcloning Efficiency)	243	T7 Express <i>lysY</i> Competent <i>E. coli</i>	
Electrocompetent	243	(High Efficiency)	247
NEB 5-alpha F' <i>I^q</i> Competent <i>E. coli</i>		T7 Express <i>lysY/I^q</i> Competent <i>E. coli</i>	
(High Efficiency)	243	(High Efficiency)	247
NEB Stable Competent <i>E. coli</i>		SHuffle Express Competent <i>E. coli</i>	248
(High Efficiency)	244	SHuffle T7 Express Competent <i>E. coli</i>	248
<i>dam⁻/dcm⁻</i> Competent <i>E. coli</i>	244	SHuffle T7 Competent <i>E. coli</i>	249
NEB Tube Opener	244	SHuffle T7 Express <i>lysY</i> Competent <i>E. coli</i>	249
<i>Companion Products:</i>		How do I express my protein in SHuffle cells?	249
SOC Outgrowth Medium	242		
NEB 10-beta/Stable Outgrowth Medium	243		

Strain Properties

There are many properties to consider when choosing a strain for your experiments.

Requirements such as plasmid preparation, blue/white screening, proper disulfide bond formation and fast colony growth necessitate specific strain choices. The following selection chart highlights the characteristics of NEB's strains to help select the optimal strain for a particular experiment.

CAUTION: Chemically Competent *E. coli* contain DMSO, a hazardous material. Review the MSDS before handling.

Cloning Strain Properties

STRAIN PROPERTIES	FEATURES	CHEMICAL TRANSFORMATION EFFICIENCY (cfu/μg)	ELECTRO-COMPETENT TRANSFORMATION EFFICIENCY (cfu/μg)	AVAILABLE FORMATS ⁽¹⁾	OUTGROWTH MEDIUM & CONTROL PLASMID INCLUDED?	STRAIN BACKGROUND	LIBRARY CONSTRUCTION
<i>dam^r/dcm^r</i>	• Dam/Dcm methyltransferase free plasmid growth	1-3 x 10 ⁶	n/a	50, 200	●	K12	
NEB Turbo (High Efficiency)	• Fastest growth – colonies visible after 6.5 hours • Plasmid preparation after 4 hours	1-3 x 10 ⁹	n/a	50, 200	●	K12	●
NEB 5-alpha (High Efficiency)	• Versatile cloning strain • DH5α [™] derivative	1-3 x 10 ⁹ ⁽⁸⁾	> 1 x 10 ¹⁰	50, 200, 96, 384, Strips	●	K12	●
NEB 5-alpha F' I ^q (High Efficiency)	• Toxic gene cloning • F' strain with extremely high transformation efficiency	1-3 x 10 ⁹	n/a	50, 200	●	K12	●
NEB 10-beta (High Efficiency)	• Large plasmid and BAC cloning • DH10B [™] derivative	1-3 x 10 ⁹ ⁽⁹⁾	> 2 x 10 ¹⁰	50, 200, 96	●	K12	●
NEB Stable (High Efficiency)	• Cloning unstable inserts • Isolating and propagating retroviral/lentiviral clones	1-3 x 10 ⁹	n/a	50, 200	●	K12	●
NEB 5-alpha (Subcloning Efficiency)	• Ideal for subcloning efficiency transformations, such as plasmid transformation or routine subcloning	> 1 x 10 ⁶	n/a	400	—	K12	

Protein Expression Strain Properties

STRAIN PROPERTIES	FEATURES	CHEMICAL TRANSFORMATION EFFICIENCY (cfu/μg)	AVAILABLE FORMATS ⁽¹⁾	OUTGROWTH MEDIUM & CONTROL PLASMID INCLUDED?	STRAIN BACKGROUND	LIBRARY CONSTRUCTION
NEBExpress	• Versatile non-T7 expression strain • Protease deficient	0.6-1 x 10 ⁹	50, 200	●	B	●
BL21(DE3)	• Routine T7 expression	1-5 x 10 ⁷	50, 200	●	B	
Lemo21(DE3)	• Tunable T7 expression for difficult targets	1-3 x 10 ⁷	50	● ⁽¹⁰⁾	B	
NiCo21(DE3)	• Improved purity of target proteins isolated by IMAC	1-5 x 10 ⁷	50	●	B	
BL21	• Routine non-T7 expression	1-5 x 10 ⁷	50	●	B	
T7 Express	• Most popular T7 expression strain • Protease deficient	0.6-1 x 10 ⁹	50, 200	●	B	●
T7 Express <i>lysY</i>	• T7 expression • Protease deficient • Better reduction of basal expression	0.6-1 x 10 ⁹	200	—	B	●
T7 Express <i>lysY/I^q</i>	• T7 expression • Protease deficient • Highest level of expression control	0.6-1 x 10 ⁹	200	—	B	●
SHuffle T7	• T7 expression/K12 strain • Enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm	1 x 10 ⁶	50	—	K12	
SHuffle Express	• Protease deficient/B strain • Enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm	1 x 10 ⁷	50	—	B	
SHuffle T7 Express	• T7 expression • Protease deficient/B strain • Enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm	1 x 10 ⁷	50	—	B	
SHuffle T7 Express <i>lysY</i>	• T7 expression • Protease deficient/B strain • Tightly controlled expression of toxic proteins • Enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm	1 x 10 ⁷	50	—	B	
NEBExpress I ^q	• Control of IPTG induced expression from <i>P_{lac}</i> , <i>P_{tac}</i> , <i>P_{trc}</i> and <i>T5_{lac}</i> • Protease deficient	0.6-1 x 10 ⁹	200	—	B	●

Competitor Cross Reference

Using another competent cell strain? Try our competitor cross reference tool to find out which NEB strain is compatible.

Cloning Strain Properties

BLUE/WHITE SCREENING	<i>lacZ</i> ⁽¹⁾	<i>endA</i> ⁽²⁾	F'	<i>recA</i> ⁽³⁾	DRUG RESISTANCE ⁽⁴⁾	METHYLATION PHENOTYPE
—	—	●	—	—	cam, str, nit	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
●	●	●	●	—	nit	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
●	—	●	—	●	none	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
●	●	●	●	●	tet	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
●	—	●	—	●	str	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
●	●	●	●	●	tet, str	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
●	—	●	—	●	none	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺

- (1) Rhamnose solution is provided instead of SOC; control plasmid is included.
 (2) Important for high-quality plasmid preparation.
 (3) Lacks Lon and OmpT protease activity.
 (4) Constitutively expresses a chromosomal copy of the disulfide bond isomerase DsbC.
 (5) nit = nitrofurantoin, tet = tetracycline, cam = chloramphenicol, str = streptomycin, spec = spectinomycin
 (6) Resistance to low levels of streptomycin may be observed.
 (7) Legend 50 = 50 µl tubes; 200 = 200 µl tubes; 96 = 96 well plate; 384 = 384 well plate; strips = 96 tube strips (50 µl/tube); 400 = 400 µl tubes
 (8) 1-5 x 10⁸ for R-format.
 (9) 1-3 x 10⁸ for P-format.



For help with choosing the right competent cell strain, try NEBcloner at NEBcloner.neb.com

Protein Expression Strain Properties

<i>lacZ</i> ⁽¹⁾	<i>lysY</i>	<i>endA</i> ⁽²⁾	PROTEASE DEFICIENT ⁽³⁾	F'	T7 RNA POLYMERASE	CYTOPLASMIC DISULFIDE BOND FORMATION ⁽⁴⁾	DRUG RESISTANCE ⁽⁵⁾	METHYLATION PHENOTYPE
—	—	●	●	—	—	—	nit	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
—	—	—	●	—	●	—	none	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
—	●	—	●	—	●	—	cam	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
—	—	—	●	—	●	—	none	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
—	—	—	●	—	—	—	none	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
—	—	●	●	—	●	—	nit	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
—	●	●	●	—	●	—	cam, nit	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
●	●	●	●	—	●	—	cam, nit	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
●	—	—	—	●	●	●	str, spec, nit	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
●	—	●	●	—	—	●	spec ⁽⁶⁾ , nit	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
●	—	●	●	—	●	●	spec ⁽⁶⁾ , nit	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
●	●	●	●	—	●	●	cam, spec ⁽⁶⁾ , nit	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺
●	—	●	●	—	—	—	cam, nit	<i>Dam</i> ⁺ , <i>Dcm</i> ⁺ , M. EcoKI ⁺

Learn how to perform a transformation.



NEB Cloning Competent *E. coli* Sampler

#C1010S 8 vials x 0.05 ml\$122

Companion Product:

SOC Outgrowth Medium
#B9020S 4 x 25 ml\$85

- Outgrowth medium and control plasmid included
- Value pricing
- Free of animal products

Description: A sample pack of four cloning strains of *E. coli* suitable for high efficiency transformation.

Please refer to the individual datacards for each reagent's recommended use and storage conditions.


Transformation Efficiency:

1–3 × 10⁹ cfu/μg pUC19 DNA

Resistance: T1 phage (*thiA2*)

The Sampler Includes:

- NEB Turbo Competent *E. coli* (High Efficiency)
- NEB 5-alpha Competent *E. coli* (High Efficiency)
- NEB 10-beta Competent *E. coli* (High Efficiency)
- NEB Stable Competent *E. coli* (High Efficiency)
- SOC Outgrowth Medium
- NEB 10-beta/Stable Outgrowth Medium
- pUC19 Vector



Benefit from high transformation efficiencies: Transformation efficiencies were compared using manufacturers' recommended protocols. Values shown are the average of triplicate experiments.

NEB Turbo Competent *E. coli*

NEB Turbo Competent *E. coli*
(High Efficiency)

#C2984H 20 x 0.05 ml\$270
#C2984I 6 x 0.2 ml\$205

Companion Product:

SOC Outgrowth Medium
#B9020S 4 x 25 ml\$85

- Tight expression control (*lac^H*)
- Colonies visible after 6.5 hours
- Isolate DNA after 4 hrs growth
- 5 minute transformation protocol with Amp^R plasmids
- Free of animal products

Description: *E. coli* cells featuring fast colony growth (6.5 hours) and tight expression control.

Genotype: F' *proA⁺B⁺ lac^H ΔlacZM15 / thiA2 Δ(lac-proAB) glnV galK16 galE15 R(zgb-210::Tn10)^{Tet^S} endA1 thi-1 Δ(hsdS-mcrB)5*

Features:


- Suitable for blue/white screening
- Activity of nonspecific endonuclease I (*endA1*) eliminated for highest quality plasmid preparations

Transformation Efficiency:

High Efficiency: 1–3 × 10⁹ cfu/μg pUC19 DNA

Resistance: T1 phage (*thiA2*), Nit

Sensitivity: Amp, Cam, Kan, Spec, Str, Tet



Miniprep DNA can be prepared from a single overnight colony after inoculation and only 3 hours growth. DNA yield doubles after an additional hour of growth.

Reagents Supplied:

SOC Outgrowth Medium
pUC19 Control DNA



NEB Turbo Transformation: With NEB Turbo, colonies are visible after only 8 hours. Ligation products were transformed into 50 μl of NEB Turbo Competent *E. coli* and plated on LB/Amp. Plates were incubated for 8 hours, 10 hours and 12 hours at 37°C. NEB Turbo features fast colony growth and blue/white selection to simplify cloning experiments.

ONESHOT® and MAX EFFICIENCY® are registered trademarks of Thermo Fisher Scientific.
DH5A®, DH10B® and MACH1® are trademarks of Thermo Fisher Scientific.



What is the difference between chemical transformation and electroporation?

NEB 10-beta Competent *E. coli*

NEB 10-beta Competent *E. coli*
(High Efficiency)

#C3019H	20 x 0.05 ml\$233
#C3019I	6 x 0.2 ml\$179
#C3019P	1 x 96 well plate\$530

NEB 10-beta Electrocompetent <i>E. coli</i>		
#C3020K	6 x 0.1 ml\$221

Companion Product:

NEB 10-beta/Stable Outgrowth Medium		
#B9035S	4 x 25 ml\$85

- Clone large plasmids and BACs
- DH10B derivative
- Free of animal products

Description: A DH10B derivative suitable for a wide range of applications, including large plasmid and BAC cloning.

Genotype: $\Delta(ara-leu)$ 7697 *araD139* *thiA* $\Delta lacX74$ *galK16* *galE15* *e14-* $\phi 80dlacZ\Delta M15$ *recA1* *relA1* *endA1* *nupG* *rpsL* (Str^R) *rph* *spoT1* $\Delta(mrr-hsdRMS-mcrBC)$

Features:

- Efficient transformation of methylated DNA derived from eukaryotic sources or unmethylated DNA derived from PCR, cDNA and many other sources
- Suitable for blue/white screening without IPTG
- Activity of nonspecific endonuclease I (*endA1*) eliminated for highest quality plasmid preparations
- Reduced recombination of cloned DNA (*recA1*)

Transformation Efficiency:

High Efficiency: $1-3 \times 10^9$ cfu/ μ g pUC19 DNA (NEB #C3019H, #C3019I); $1-3 \times 10^8$ cfu/ μ g pUC19 DNA (NEB #C3019P)


Electrocompetent: $> 2 \times 10^{10}$ cfu/ μ g pUC19 DNA

Resistance: T1 phage (*thiA*), Str

Sensitivity: Amp, Cam, Kan, Nit, Spec, Tet

Reagents Supplied:

NEB 10-beta/Stable Outgrowth Medium
pUC19 Control DNA



Effect of Plasmid Size on Transformation Efficiency:

NEB 10-beta chemically competent cells are more efficiently transformed with large plasmids than NEB 5-alpha cells. The difference in TE between the two cell lines increases with the size of the plasmid being transformed.

NEB 5-alpha Competent *E. coli*

NEB 5-alpha Competent *E. coli*
(High Efficiency)

#C2987H	20 x 0.05 ml\$199
#C2987I	6 x 0.2 ml\$153
#C2987P	1 x 96 well plate\$503
#C2987R	1 x 384 well plate\$1,100
#C2987U	96 x 50 μ l/tube\$748

NEB 5-alpha Competent *E. coli*
(Subcloning Efficiency)*

#C2988J	6 x 0.4 ml\$73
---------	------------	-----------

NEB 5-alpha Electrocompetent <i>E. coli</i>		
#C2989K	6 x 0.1 ml\$188

Companion Product:

SOC Outgrowth Medium		
#B9020S	4 x 25 ml\$85

- DH5 α derivative
- Free of animal products

Description: A DH5 α derivative and versatile *E. coli* cloning strain.

Genotype: *thiA2* Δ (*argF-lacZ*)*U169* *phoA* *glnV44* $\phi 80\Delta$ (*lacZ*)*M15* *gyrA96* *recA1* *relA1* *endA1* *thi-1* *hsdR17*

Features:

- Efficient transformation of unmethylated DNA derived from PCR, cDNA and many other sources (*hsdR*)
- Suitable for blue/white screening
- Activity of nonspecific endonuclease I (*endA1*) eliminated for highest quality plasmid preparations
- Reduced recombination of cloned DNA (*recA1*)

Transformation Efficiency:

High Efficiency: $1-3 \times 10^9$ cfu/ μ g pUC19 DNA (NEB #C2987H, #C2987I, #C2987P, #C2987U); $1-5 \times 10^8$ cfu/ μ g pUC19 DNA (NEB #C2987R)

Subcloning Efficiency: $> 1 \times 10^6$ cfu/ μ g pUC19 DNA

Electrocompetent: $> 1 \times 10^{10}$ cfu/ μ g pUC19 DNA

Resistance: T1 phage (*thiA2*)

Sensitivity: Amp, Cam, Kan, Nit, Spec, Str, Tet

Reagents Supplied:

SOC Outgrowth Medium
pUC19 Control DNA

* NEB 5-alpha Competent *E. coli* (Subcloning Efficiency) is not supplied with SOC Outgrowth Medium or pUC19 Control DNA.

NEB 5-alpha F' I^q Competent *E. coli*

NEB 5-alpha F' I^q Competent *E. coli*
(High Efficiency)

#C2992H	20 x 0.05 ml\$202
#C2992I	6 x 0.2 ml\$158

Companion Product:

SOC Outgrowth Medium		
#B9020S	4 x 25 ml\$85

- Tight expression control (*lacI*^q)
- F' Strain with extremely high TE
- DH5 α derivative
- Free of animal products

Description: An F' *E. coli* strain with extremely high transformation efficiency suitable for toxic gene cloning.

Genotype: F' *proA*⁺*B*⁺ *lacI*^q Δ (*lacZ*)*M15* *zff::Tn10* (Tet^R) / *thiA2* Δ (*argF-lacZ*)*U169* *phoA* *glnV44* $\phi 80\Delta$ (*lacZ*)*M15* *gyrA96* *recA1* *relA1* *endA1* *thi-1* *hsdR17*

Features:

- Efficient transformation of unmethylated DNA derived from PCR, cDNA and many other sources (*hsdR*)
- Suitable for blue/white screening
- Activity of nonspecific endonuclease I (*endA1*) eliminated for highest quality plasmid preparations
- Reduced recombination of cloned DNA (*recA1*)
- Suitable for propagation of M13 clones

Transformation Efficiency:

High Efficiency: $1-3 \times 10^9$ cfu/ μ g pUC19 DNA

Resistance: T1 phage (*thiA2*), Tet

Sensitivity: Amp, Cam, Kan, Nit, Spec, Str

Reagents Supplied:

SOC Outgrowth Medium
pUC19 Control DNA

NEB Stable Competent *E. coli*

NEB Stable Competent *E. coli*
(High Efficiency)

#C3040H 20 x 0.05 ml\$330
#C3040I 6 x 0.2 ml\$254

Companion Product:

NEB 10-beta/Stable Outgrowth Medium
#B9035S 4 x 25 ml\$85

- T1 phage resistance (*fhuA*)
- Free of animal products
- Carries *endA* mutation (isolated plasmids are free of *EndoI*)

Description: Chemically competent *E. coli* cells suitable for high efficiency transformation and isolation of plasmid clones containing repeat elements.

Genotype: F⁻ *proA*⁻ *B*⁻ *lacZ*^Δ(*lacZ*)M15 *zrf::Tn10* (Tet^R)/Δ(*ara-leu*) 7697 *araD139 fhuA ΔlacX74 galK16 galE15 e14- φ80dlacZΔM15 recA1 relA1 endA1 nupG rpsL* (Str^R) *rph spoT1 Δ(mrr-hsdRMS-mcrBC)*

Features:

- Activity of nonspecific endonuclease I (*endA1*) abolished for highest quality plasmid preparations
- Reduced recombination of cloned DNA (*recA1*)

Applications:

- Cloning unstable inserts

- Isolating and propagating retroviral/lentiviral clones
- Compatible with DNA Assembly reactions, as well as ligation reactions

Transformation Efficiency:

High Efficiency: 1–3 × 10⁹ cfu/μg pUC19 DNA

Resistance: T1 phage (*fhuA*), Str, Tet

Sensitivity: Amp, Cam, Kan, Nit, Spec

Reagents Supplied:

NEB 10-beta/Stable Outgrowth Medium
pUC19 Control DNA

dam⁻ / *dcm*⁻ Competent *E. coli*

#C2925H 20 x 0.05 ml\$251
#C2925I 6 x 0.2 ml\$193

Companion Product:

SOC Outgrowth Medium
#B9020S 4 x 25 ml\$85

- Isolate plasmids free of *Dam* and *Dcm* methylation
- Free of animal products

Description: Methyltransferase deficient *E. coli* cells suitable for growth of plasmids free of *Dam* and *Dcm* methylation.

Genotype: *ara-14 leuB6 fhuA31 lacY1 tsx78 glnV44 galK2 galT22 mcrA dcm-6 hisG4 rtdD1 R(zgb210::Tn10) Tet^S endA1 rslP136* (Str^R) *dam13::Tn9* (Cam^R) *xylA-5 mtl-1 thi-1 mcrB1 hsdR2*

Features:

- Allows for propagation of plasmids free of *Dam* and *Dcm* methylation
- Activity of nonspecific endonuclease I (*endA1*) abolished for highest quality plasmid preparations

Transformation Efficiency:

1–3 × 10⁶ cfu/μg pUC19 DNA

Resistance: T1 phage (*fhuA31*), Cam, Nit, Str

Sensitivity: Amp, Kan, Spec, Tet


Reagents Supplied:

SOC Outgrowth Medium
pUC19 Control DNA

NEB Tube Opener

#C1008S 2 units\$21

Description: Use to open a variety of microcentrifuge tubes. Can be used for snap-on caps or screw-on caps.



BL21 Competent *E. coli*

#C2530H 20 x 0.05 ml\$209

Companion Product:

SOC Outgrowth Medium
#B9020S 4 x 25 ml\$85

- Ideal for P_{lac} , P_{tac} , P_{trc} , *ParaBAD* expression vectors
- Resistance to phage T1 (*fhuA2*)
- Protease deficient B strain
- Free of animal products

Description: Widely used non-T7 expression *E. coli* strain. Suitable for transformation and protein expression. This strain does not express the T7 RNA Polymerase.

Genotype: *fhuA2 [lon] ompT gal [dcm] ΔhsdS*

Features:

- Deficient in proteases Lon and OmpT

Transformation Efficiency:

1–5 x 10⁷ cfu/μg pUC19 DNA

Resistance: T1 phage (*fhuA2*)

Sensitivity: Amp, Cam, Kan, Nit, Spec, Str, Tet

Reagents Supplied:

SOC Outgrowth Medium
pUC19 Control DNA

BL21(DE3) Competent *E. coli*

#C2527H 20 x 0.05 ml\$209

#C2527I 6 x 0.2 ml\$158

Companion Product:

SOC Outgrowth Medium
#B9020S 4 x 25 ml\$85

- Routine T7 expression
- Resistance to phage T1 (*fhuA2*)
- Protease deficient B strain
- Free of animal products

Description: Widely used T7 expression *E. coli* strain.

Genotype: *fhuA2 [lon] ompT gal (λ DE3) [dcm] ΔhsdS*
λ DE3 = λ. sBamHlo ΔEcoRI-B int::(lacI::PlacUV5::T7 gene1) i21 Δnin5

Features:

- Deficient in proteases Lon and OmpT
- Resistant to phage T1 (*fhuA2*)

Transformation Efficiency:

1–5 x 10⁷ cfu/μg pUC19 DNA

Resistance: T1 phage (*fhuA2*)

Sensitivity: Amp, Cam, Kan, Nit, Spec, Str, Tet

Reagents Supplied:

SOC Outgrowth Medium
pUC19 Control DNA

Lemo21(DE3) Competent *E. coli*

#C2528J 12 x 0.05 ml\$233

- Expression of difficult targets
- Membrane protein expression
- Ideal for periplasmic expression
- Expression of toxic proteins
- Proteins with solubility issues

Description: Lemo21(DE3) Competent *E. coli* is a tunable T7 expression strain designed for the expression of challenging proteins. A derivative of BL21(DE3), Lemo21(DE3) offers the host features of this popular expression strain, with the added benefit of being able to control expression levels by varying the level of T7 lysozyme (*lysY*), the natural inhibitor of T7 RNA Polymerase. The fine control of expression makes Lemo21(DE3) ideal for membrane proteins, toxic proteins, secreted proteins and proteins prone to insoluble expression.

Genotype: *fhuA2 [lon] ompT gal (λ DE3) [dcm] ΔhsdS/pLemo(Cam^R) λ DE3 = λ. sBamHlo ΔEcoRI-B int::(lacI::PlacUV5::T7 gene1) i21 Δnin5*
pLemo = pACYC184-*PrhaBAD-lysY*

Features:

- Enhanced BL21(DE3) derivative
- Fine control of expression
- Greatest range of expression of any T7 strain (0–2,000 μM rhamnose)
- Potential elimination of inclusion body formation

Transformation Efficiency:

High Efficiency: 1–3 x 10⁷ cfu/μg pUC19 DNA

Resistance: T1 phage (*fhuA2*), Cam

Sensitivity: Amp, Kan, Nit, Spec, Str, Tet

Reagents Supplied:

L-rhamnose solution
pUC19 control DNA

NiCo21(DE3) Competent *E. coli*

#C2529H 20 x 0.05 ml\$388

Companion Product:

SOC Outgrowth Medium
#B9020S 4 x 25 ml\$85

- Superior alternative to BL21(DE3) for routine protein expression
- Improved purity of target proteins isolated by IMAC
- Free of animal products

Description: Poly-histidine tagged recombinant proteins that are isolated by immobilized metal affinity chromatography (IMAC) are often contaminated with significant amounts of endogenous *E. coli* metal binding proteins. The protein expression strain NiCo21(DE3) has been engineered to minimize *E. coli* protein contamination of IMAC fractions: GlnS is mutated to eliminate binding to IMAC resins and three other proteins (SlyD, ArnA and Can) are tagged to enable rapid removal by chitin affinity chromatography.

Genotype: *can::CBD thuA2 [lon] ompT gal (λ DE3) [dcm] araA::CBD slyD::CBD glmS6Ala ΔhsdS λ DE3 = λ sBamHI ΔEcoRI-B int::(lac::PlacUV5::T7 gene1) i21 Δnin5*

Features:

- Identical growth characteristics as BL21(DE3)
- Deficient in proteases Lon and OmpT

Transformation Efficiency:

1–5 x 10⁷ cfu/μg pUC19 DNA

Resistance: T1 phage (*thuA2*)

Sensitivity: Amp, Cam, Kan, Nit, Spec, Str, Tet

Reagents Supplied:

SOC Outgrowth Medium
pUC19 Control DNA

NEBExpress Competent *E. coli*

NEBExpress Competent *E. coli*
(High Efficiency)

#C2523H 20 x 0.05 ml\$209
#C2523I 6 x 0.2 ml\$158

Companion Product:

SOC Outgrowth Medium
#B9020S 4 x 25 ml\$85

- Enhanced BL21 derivative ideal for P_{lac} , $P_{lac'}$, P_{trc} expression vectors
- Fast growth from colonies
- Free of animal products
- Protease deficient

Description: A versatile non-T7 expression *E. coli* strain. NEBExpress is the recommended host strain for the NEBExpress MBP Protein Fusion and Purification System (NEB #E8201).

Genotype: *thuA2 [lon] ompT gal sulA11 R(mcr-73::miniTn10--Tet^s)2 [dcm] R(zgb-210::Tn10--Tet^s) endA1 Δ(mcrC-mrr)114::IS10*

Features:

- Deficient in proteases Lon and OmpT
- Does not restrict methylated DNA

Transformation Efficiency:

High Efficiency: 0.6–1 x 10⁹ cfu/μg pUC19 DNA

Resistance: T1 phage (*thuA2*), Nit

Sensitivity: Amp, Cam, Kan, Tet, Spec, Str

Reagents Supplied:

SOC Outgrowth Medium
pUC19 Control DNA

NEBExpress I^q Competent *E. coli*

NEBExpress I^q Competent *E. coli*
(High Efficiency)

#C3037I 6 x 0.2 ml\$177

Companion Product:

SOC Outgrowth Medium
#B9020S 4 x 25 ml\$85

- Enhanced BL21 derivative ideal for P_{lac} , $P_{lac'}$, P_{trc} , P_{T5} expression vectors
- Better control of IPTG induced expression with non-T7 plasmids
- Fast growth from colonies
- *lacI^q* reduces basal expression
- Protease deficient
- Free of animal products

Description: *E. coli* cells featuring control of IPTG induced expression with non-T7 plasmids.

Genotype: *MiniF lacI^q (Cam^R) / thuA2 [lon] ompT gal sulA11 R(mcr-73::miniTn10--Tet^s)2 [dcm] R(zgb-210::Tn10--Tet^s) endA1 Δ(mcrC-mrr)114::IS10*

Features:

- Deficient in proteases Lon and OmpT
- Does not restrict methylated DNA
- Ideal for controlled protein expression from pUC19 and pUC19 derivatives

Transformation Efficiency:

High Efficiency: 0.6–1 x 10⁹ cfu/μg pUC19 DNA

Resistance: T1 phage (*thuA2*), Cam, Nit

Sensitivity: Amp, Kan, Spec, Str, Tet

T7 Express Competent *E. coli*

T7 Express Competent *E. coli*
(High Efficiency)

#C2566H 20 x 0.05 ml\$209
#C2566I 6 x 0.2 ml\$158

Companion Product:

SOC Outgrowth Medium
#B9020S 4 x 25 ml\$85

- Enhanced BL21 derivative
- Popular T7 expression strain
- Fast growth from colonies
- Free of animal products

Description: Enhanced BL21 *E. coli* derivative for T7 expression.

Genotype: *fhuA2 lacZ::T7 gene1 [lon] ompT gal sulA11 R(mcr-73::miniTn10--Tet^S)2 [dcm] R(zgb-210::Tn10--Tet^S) endA1 Δ(mcrC-mrr)114::IS10*

Features:

- T7 RNA Polymerase in the *lac* operon - no λ prophage
- Deficient in proteases Lon and OmpT
- Does not restrict methylated DNA

Transformation Efficiency:

High Efficiency: 0.6–1 x 10⁹ cfu/μg pUC19 DNA

Resistance: T1 phage (*fhuA2*), Nit

Sensitivity: Amp, Cam, Kan, Spec, Str, Tet

Reagents Supplied:

SOC Outgrowth Medium
pUC19 Control DNA

T7 Express *lysY* Competent *E. coli*

T7 Express *lysY* Competent *E. coli*
(High Efficiency)

#C3010I 6 x 0.2 ml\$177

Companion Product:

SOC Outgrowth Medium
#B9020S 4 x 25 ml\$85

- Enhanced BL21 derivative
- T7 Lysozyme for expression control
- Clone toxic genes
- Fast growth from colonies
- Free of animal products

Description: Enhanced BL21 *E. coli* derivative for T7 expression with enhanced reduction of basal expression.

Genotype: MiniF *lysY* (Cam^R) / *fhuA2 lacZ::T7 gene1 [lon] ompT gal sulA11 R(mcr-73::miniTn10--Tet^S)2 [dcm] R(zgb-210::Tn10--Tet^S) endA1 Δ(mcrC-mrr)114::IS10*

Features:

- T7 RNA Polymerase in the *lac* operon - no λ prophage
- Control of T7 RNA Polymerase by T7 lysozyme allows potentially toxic genes to be expressed
- LysY is a variant of T7 lysozyme lacking amidase activity, thus cells are not susceptible to lysis during induction
- Deficient in proteases Lon and OmpT
- Does not restrict methylated DNA
- No Cam requirement

Transformation Efficiency:

High Efficiency: 0.6–1 x 10⁹ cfu/μg pUC19 DNA

Resistance: T1 phage (*fhuA2*), Cam, Nit

Sensitivity: Amp, Kan, Spec, Str, Tet

T7 Express *lysY/I^q* Competent *E. coli*

T7 Express *lysY/I^q* Competent *E. coli*
(High Efficiency)

#C3013I 6 x 0.2 ml\$177

Companion Product:

SOC Outgrowth Medium
#B9020S 4 x 25 ml\$85

- Enhanced BL21 derivative
- Tight control of expression (*lacI^q*)
- Highest level of expression control
- Clone toxic genes
- Fast growth from colonies
- Free of animal products

Description: Enhanced BL21 *E. coli* derivative with highest level of T7 expression control.

Genotype: MiniF *lysY lacI^q* (Cam^R) / *fhuA2 lacZ::T7 gene1 [lon] ompT gal sulA11 R(mcr-73::miniTn10--Tet^S)2 [dcm] R(zgb-210::Tn10--Tet^S) endA1 Δ(mcrC-mrr)114::IS10*

Features:

- T7 RNA Polymerase in the *lac* operon - no λ prophage
- Tight control of expression by *lacI^q* allows potentially toxic genes to be cloned
- Control of T7 RNA Polymerase by T7 lysozyme allows toxic genes to be expressed
- LysY is a variant of T7 lysozyme lacking amidase activity, thus cells are less susceptible to lysis during induction
- Deficient in proteases Lon and OmpT
- Does not restrict methylated DNA
- No Cam requirement

Transformation Efficiency:

High Efficiency: 0.6–1 x 10⁹ cfu/μg pUC19 DNA

Resistance: T1 phage (*fhuA2*), Cam, Nit

Sensitivity: Amp, Kan, Spec, Str, Tet

Disulfide Bonds

Features of Shuffle® Strains:


- Engineered *E. coli* K12 or B strains promote disulfide bond formation in the cytoplasm
- Constitutively expresses a chromosomal copy of the disulfide bond isomerase DsbC
- DsbC promotes the correction of mis-oxidized proteins into their correct form
- The cytoplasmic DsbC is a chaperone that can also assist in the folding of proteins that do not require disulfide bonds
- Alternative expression strain for proteins that do not fold in wild-type *E. coli*, independent of redox state

SHuffle strains from NEB are engineered *E. coli* strains capable of expressing proteins with increasing disulfide bond complexity in the cytoplasm. SHuffle strains express the disulfide bond isomerase DsbC within the cytoplasm. DsbC isomerizes mis-oxidized substrates into their correctly folded state greatly enhancing the fidelity of disulfide bond formation. Cytoplasmic expression also results in significantly higher protein yields of

disulfide bonded proteins when compared to periplasmic expression. SHuffle strains are sensitive to kan, amp, tet and in most cases, cam, which makes them able to express proteins from a wide variety of expression vectors offering greater versatility in experimental design.

References:

References for properties and applications for these products can be found at www.neb.com.



SHuffle Express Competent *E. coli*

#C3028J 12 x 0.05 ml\$233

Description: *E. coli* cells with enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm.

Genotype: *fhuA2 [lon] ompT ahpC gal λatt::pNEB3-r1-cDsbC (Spec^R, lac^R) ΔtrxB sulA11 R(mcr-73::miniTn10-Tet^S)2 [dcm] R(zgb-210::Tn10-Tet^S) endA1 Δgor Δ(mcrC-mrr)114::IS10*

Transformation Efficiency:

1 x 10⁷ cfu/μg pUC19 DNA

Resistance: T1 phage (*fhuA2*), Nit, Str^R, Spec

Sensitivity: Amp, Cam, Kan, Tet

*Resistance to low levels of streptomycin may be observed.

- Folds disulfide bonded proteins in the cytoplasm
- Protease deficient
- Enhanced BL21 derivative
- Free of animal products

SHuffle T7 Express Competent *E. coli*

#C3029J 12 x 0.05 ml\$233

Description: T7 Expression *E. coli* strain with enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm.

Genotype: *fhuA2 lacZ::T7 gene1 [lon] ompT ahpC gal λatt::pNEB3-r1-cDsbC (Spec^R, lac^R) ΔtrxB sulA11 R(mcr-73::miniTn10-Tet^S)2 [dcm] R(zgb-210::Tn10-Tet^S) endA1 Δgor Δ(mcrC-mrr)114::IS10*

Transformation Efficiency: 1 x 10⁷ cfu/μg pUC19 DNA

Resistance: T1 phage (*fhuA2*), Nit, Str^R, Spec

Sensitivity: Amp, Cam, Kan, Tet

*Resistance to low levels of streptomycin may be observed.

- Folds disulfide bonded proteins in the cytoplasm
- T7 expression
- Protease deficient B strain
- Enhanced BL21 derivative
- Free of animal products



What is a disulfide bond?

SHuffle T7 Competent *E. coli*

#C3026J 12 x 0.05 ml\$233

- Folds disulfide bonded proteins in the cytoplasm
- T7 expression
- K12 strain
- Free of animal products

Description: T7 Expression *E. coli* strain with enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm.

Genotype: F' *lac, pro, lac^R / Δ(ara-leu)7697 araD139 fhuA2 lacZ::T7 gene1 Δ(phoA)PvuII phoR ahpC* galE (or U) galK λatt::pNEB3-r1-cDsbC (Spec^R, lac^R) ΔtrxB rpsL150(Str^R) Δgor Δ(malF)3*

Transformation Efficiency: 1 x 10⁶ cfu/μg pUC19 DNA

Resistance: T1 phage (*fhuA2*), Nit, Spec, Str

Sensitivity: Amp, Cam, Kan, Tet

SHuffle T7 Express *lysY* Competent *E. coli*

#C3030J 12 x 0.05 ml\$233

- Express toxic proteins (*lysY*)
- Folds disulfide bonded proteins in the cytoplasm
- T7 expression
- Protease deficient B strain
- Free of animal products
- Enhanced BL21 derivative

Description: *E. coli* strain with tight T7 Expression control and enhanced capacity to correctly fold proteins with multiple disulfide bonds in the cytoplasm.

Genotype: MiniF *lysY* (Cam^R) / *fhuA2 lacZ::T7 gene1 [lon] ompT ahpC gal λatt::pNEB3-r1-cDsbC (Spec^R, lac^R) ΔtrxB sulA11 R(mcr-73::miniTn10--Tet^S)2 [dcm] R(zgb-210::Tn10--Tet^S) endA1 Δgor Δ(mcrC-mrr)114::IS10*

Transformation Efficiency:

1 x 10⁷ cfu/μg pUC19 DNA

Resistance: T1 phage (*fhuA2*), Cam, Nit, Str*, Spec

Sensitivity: Amp, Kan, Tet

*Resistance to low levels of streptomycin may be observed.

How do I express my protein in SHuffle cells?

Currently there are two SHuffle cell lines available from NEB; SHuffle (NEB #C3026) based on *E. coli* K12, and SHuffle Express (NEB #C3028, #C3029, #C3030) based on *E. coli* B.

We recommend testing both B and K12 expression strains, as we do see variability in expression depending on the protein of interest (Table 1). If T7 expression is not necessary, then we recommend comparing NEB #C3026 and #C3028. If T7 expression is necessary, test NEB #C3026 and #C3029. If T7-driven expression of a protein is toxic, switch to a non-leaky *lysY* version (NEB #C3030). Once the strain is chosen expression conditions should be optimized. This can include temperature as well as auto expression (1).

View our online tutorial for tips on setting up reactions with SHuffle.

Table 1. Percentage of relative solubility of various proteins using SHuffle (K12 and B strains):

PROTEIN	RELATIVE % SOLUBLE		# CYSTEINES
	K12	B	
Gluc	65	100	10
Urokinase	60	100	24
vtPA	5	100	12
BSA	100	0	35
Polymerase	100	0	0
Nuclease	100	10	4

Results are determined based on protein levels detected by SDS-PAGE (not shown)

Reference:

- (1) Ke, N. and Berkmen, M. (2014) *Current Protocols Molecular Biology* 16.1B.21.



Aim AR camera at this page
for a musical 3D experience.

Mostly Mozart

It is often said that music is good for the soul. It has a powerful effect on the mind and can influence our moods and emotions. Music can be social and shared collectively with our peers, or it can be deeply personal. In a diverse world, music provides us with a universal language; while musical interpretations can be culture-specific, around the globe people of all ages embrace and are unified with music. And it is well known that musical education has a long-lasting positive impact on children's social and cognitive development. Music activates the region in the frontal cortex that is also activated when performing cooperative altruistic tasks.

"Listening to Mozart makes you smarter," or so the saying goes. As it turns out, any music genre can positively impact our brain, emotional well-being and creativity. But research indicates that listening to Mozart can improve spatial reasoning skills. Studies have shown that exposure to melodic music, notably Mozart's Sonata for Two Pianos, increases dopamine levels and serotonin release — hormones associated with happiness and positive feelings.

Don Comb was a big fan of Mozart, and he was never happier than when he was working at the lab bench with Mozart playing in the background. This became a familiar soundtrack at NEB for those passing by Don's lab or office.

Don's love of Mozart was the inspiration behind an annual "Mostly Mozart" concert — arranged in his honor as a birthday gift from the NEB community. It was so named for two reasons: the program was composed solely of Mozart arrangements, and the concert was held in Mostly Hall — the historic Jacobean summer mansion located on NEB's campus in Ipswich, MA, which was lovingly restored to its original glory and houses NEB's IT, Finance and Marketing Departments. Legend has it, the manor was named Mostly Hall in 1895: When Mattina Riker, the bride of the original estate owner, James Howe Proctor, entered the large hall through the main entry she remarked, "why, it is mostly hall!" — the name has affectionately endured ever since.

Organized through Rockport Music, a Gloucester, MA-based music organization, these concerts have been held annually for more than ten years. As Mozart's music filled the hall, friends and employees in attendance could witness the joy that music brought to Don's life — he would sit in the front row, with his eyes closed, and enjoy the music with a big smile on his face.


We store musical memories in our implicit memory; they do not require conscious retrieval. This is why we sometimes conjure memories from decades gone by upon hearing just one refrain from a particular song that we attribute to a certain period of our life. Not only do memories come flooding back, but the emotions we associate with those memories can also resurface. For many at NEB, Mozart's music will forever bring a memory of Don Comb and a smile as we go about our day, honoring his legacy and working to continue the prodigious work that he began.

Piano keyboard in the
golden light.
Credit: furtseff, Adobe Stock

Experience
Mozart's Sonata
for Two Pianos.



Glycobiology & Protein Analysis Tools



Trust NEB's expertise in enzymology when you need reagents for glycobiology and protein analysis.

Glycobiology

Proteomics, the systematic study of proteins in biological systems, has expanded the knowledge of protein expression, modification, interaction and function. However, in eukaryotic cells, the majority of proteins are post-translationally modified (1). A common post-translational modification, essential for cell viability, is the attachment of glycans. Glycosylation defines the adhesive properties of glycoconjugates, and it is largely through glycan–protein interactions that cell–cell and cell–pathogen (including viruses) contacts occur, a fact that accentuates the importance of glycobiology.

Glycomics, the study of glycan expression in biological systems, relies on effective enzymatic and analytical techniques for the correlation of glycan structure with function. Glycobiology is a small but rapidly growing field in biology, with relevance to biomedicine, biotechnology, biofuels and basic research. Glycan molecules modulate many other processes important for cell and tissue differentiation, metabolic and gene regulation, protein activity, protein clearance, transport and more (2-9).

Protein Analysis Tools

Not only are proteins a major structural component of living systems, they can also be effector molecules whose states determine downstream activities. Therefore, studying the protein complement within a cell can reveal the mechanisms behind many of the cell's responses to its environment. Given the vast number of applications for protein analysis, several tools and methods for its study exist; determining the correct method for your application is paramount to success.

Phage display technology is an *in vitro* screening technique for identifying ligands for proteins and other macromolecules. At the crux of phage display technology is the ability to express peptide or protein sequences as fusions to the coat proteins of a bacteriophage. Libraries of phage-displayed peptides or proteins are thereby physically linked to their encoding nucleic acid, allowing selection of binding partners for myriad target types by iterative rounds of *in vitro* panning and amplification, followed by DNA sequencing.

All NEB products pass stringent quality control assays to ensure the highest level of functionality and purity.

Featured Products

257 Rapid™ PNGase F

268 Thermolabile Proteinase K

Featured Tools & Resources



Protein Analysis Tools & Glycomics Overview




PNGase F Overview & Selection Chart



Visit www.NEBglycosidase.com to view our online tutorial on *N*- and *O*-linked glycosylation.

(1) Khoury G.A. et al. (2011) *Scientific Reports* 1: 90. PMID: 22034591.
(2) Varki A. (1993) *Glycobiology*, 3, 97–130. PMID: 8490246.
(3) Zhao Y.Y. et al. (2008) *Cancer Sci.* 99, 1304–1310. PMID: 18492092.
(4) Zhao Y. et al. (2008) *FEBS J.* 275, 1939–1948. PMID: 18384383.
(5) Skropeta D. (2009) *Bioorg. Med. Chem.* 17, 2645–2653. PMID: 19285412.
(6) Neu U. et al. (2011) *Curr. Opin. Struct. Biol.* 21, 610–618. PMID: 21917445.
(7) Cerliani J.P. et al. (2011) *J. Clin. Immunol.* 31, 10–21. PMID: 21184154.
(8) Aarnoudse C.A. et al. (2008) *Curr. Opin. Immunol.* 18, 105–111. PMID: 16303292.
(9) Arnold J.N. (2006) *Immunol. Lett.* 106, 103–110. PMID: 16814399.



Endoglycosidases

Endo F2	254
Endo F3	254
Endo D	255
Endo H, Endo H ₁	255
Endo S	255
PNGase A	256
PNGase F and PNGase F, recombinant	256
RNase B	256
Rapid PNGase F and	
Rapid PNGase F (non-reducing format)	257
Rapid PNGase F Antibody Standard	257
Remove-IT PNGase F	257
O-Glycosidase	258
Protein Deglycosylation Mix II	258
Fetuin	258

Exoglycosidases

α -N-Acetylglactosaminidase	259
β -N-Acetylglucosaminidase S	259
β -N-Acetylhexosaminidase,	259
α 1-2 Fucosidase	259
α 1-3,4 Fucosidase	260
α 1-2,3,4,6 Fucosidase	260
α 1-2,4,6 Fucosidase O	260
α 1-3,6 Galactosidase	260
α 1-3,4,6 Galactosidase	261
β 1-3 Galactosidase	261
β 1-3,4 Galactosidase	261
β 1-4 Galactosidase S	261
α 1-2,3 Mannosidase	262
α 1-2,3,6 Mannosidase	262
α 1-6 Mannosidase	262
α 2-3,6,8 Neuraminidase	262
α 2-3,6,8,9 Neuraminidase A	263
α 2-3 Neuraminidase S	263

Heparin Lyases

<i>Bacteroides</i> Heparinase I, II, III	264
--	-----

Proteases

Endoglycoceramidase I (EGCase I)	265
IdeZ Protease	265
O-Glycoprotease	265
Trypsin-ultra, Mass Spectrometry Grade	266
α -Lytic Protease	266
Endoproteinase LysC	267
Endoproteinase GluC	267
Endoproteinase AspN	267
Trypsin-digested BSA MS Standard (CAM Modified)	267
Proteinase K, Molecular Biology Grade	268
Thermolabile Proteinase K	268
TEV Protease	268
Factor Xa Protease	269
Enterokinase, light chain	269
Furin	269

Protein Phosphatases & Kinases

Lambda Protein Phosphatase (Lambda PP)	270
<i>p</i> -Nitrophenylphosphate (PNPP)	270
Sodium Orthovanadate	270
cAMP-dependent Protein Kinase (PKA)	271
Casein Kinase II (CK2)	271
Adenosine 5' Triphosphate (ATP)	271

Phage Display

Ph.D.-7 Phage Display Peptide Library Kit	272
Ph.D.-12 Phage Display Peptide Library Kit	272
Ph.D.-C7C Phage Display Peptide Library Kit	272
Ph.D.-12 Phage Display Peptide Library	272
Ph.D. Peptide Display Cloning System	273
Anti-M13 pIII Monoclonal Antibody	273

 Recombinant Enzyme

Glycosidases

- Enabling Novel Technologies
- Unique Specifications
- Exceptional Value
- High Purity

NEB offers a selection of endoglycosidases and exoglycosidases for glycobiology research. Many of these reagents are recombinant, and all undergo several quality control assays, enabling us to provide products with lower unit cost, high purity, and reduced lot-to-lot variation.

All of our glycosidases are tested for contaminants. Since *p*-nitrophenyl-glycosides are not hydrolyzed by some exoglycosidases, we use only fluorescently-labeled oligosaccharides to screen for contaminating glycosidases.


NEB's glycosidases are provided with 10X buffer to ensure optimal activity. Using more than one glycosidase simultaneously is a common timesaving procedure. Selecting the best buffer to provide reaction conditions that optimize enzyme activity is an important consideration.

Reaction Buffer Compositions:

Visit www.NEBGlycosidase.com for details.

Endo F2

#P0772S 480 units\$215



- Removal of complex biantennary N-linked glycans from glycoproteins and glycopeptides
- Useful for determining N-glycosylation sites

Description: Endo F2 is a highly specific recombinant endoglycosidase which cleaves within the chitobiose core of asparagine-linked complex biantennary and high mannose oligosaccharides from glycoproteins and glycopeptides. Endo F2 cleaves biantennary glycans at a rate approximately 20 times greater than high mannose glycans. The activity of Endo F2 is identical on biantennary structures with and without core fucosylation. However, Endo F2 is not active on hybrid or tri- and tetra-antennary oligosaccharides. Endo F2 is tagged with a chitin binding domain (CBD) for easy removal from a reaction and is supplied glycerol free for optimal performance in HPLC and MS intensive methods.



Source: Cloned from *Elizabethkingia miricola* (formerly *Flavobacterium meningosepticum*) and expressed in *E. coli*.

Reagents Supplied:

10X GlycoBuffer 4

Molecular Weight: 39,800 daltons


Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the carbohydrate from 10 µg Porcine Fibrinogen in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 8,000 units/ml

Heat Inactivation: 65°C for 10 minutes

Endo F3

#P0771S 240 units\$195



- Removal of complex biantennary and triantennary N-linked glycans from glycoproteins and glycopeptides
- Useful for determining N-glycosylation sites

Description: Endo F3 is a highly specific recombinant endoglycosidase which cleaves within the chitobiose core of asparagine-linked fucosylated-biantennary and triantennary complex oligosaccharides from glycoproteins. Endo F3 is tagged with a chitin binding domain (CBD) for easy removal from a reaction and is supplied glycerol free for optimal performance in HPLC and MS intensive methods.

Source: Cloned from *Elizabethkingia miricola* (formerly *Flavobacterium meningosepticum*) and expressed in *E. coli*.



Reagents Supplied:

10X GlycoBuffer 4

Molecular Weight: 38,800 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the carbohydrate from 10 µg Porcine Fibrinogen in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 8,000 units/ml

Heat Inactivation: 65°C for 10 minutes




Find an overview of glycobiology.

Endo D



#P0742S	1,500 units\$195
#P0742L	7,500 units\$775



X= (H or oligosaccharide)

- Removal of paucimannose N-linked glycans from glycoproteins and glycopeptides
- Useful for determining N-glycosylation sites

Description: Endo D, also known as Endoglycosidase D, is a recombinant glycosidase, which cleaves within the chitobiose core of paucimannose N-linked glycans, with or without extensions in the antennae.

Endo D is tagged with a chitin binding domain (CBD) for easy removal from a reaction, and is supplied glycerol-free for optimal performance in HPLC and MS intensive methods.

Source: A truncated Endo D gene cloned from *Streptococcus pneumoniae* and expressed in *E. coli* as a fusion to chitin binding domain.

Reagents Supplied:

10X DTT
10X GlycoBuffer 2

Molecular Weight: 140,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to remove > 95% of the carbohydrate from 10 µg of glycosidase-trimmed (trimannosyl core) Fetuin in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 50,000 units/ml


Heat Inactivation: 65°C for 10 minutes

Endoglycosidase H



Endo H		
#P0702S	10,000 units\$75
#P0702L	50,000 units\$302

Endo H ₁		
#P0703S	100,000 units\$75
#P0703L	500,000 units\$302



Endo H and Endo H₁ cleave only high mannose structures (n = 2–150, x = (Man)_{2–20}, y = H) and hybrid structures (n = 2, x and/or y = AcNeu-Gal-GlcNAc).

- Removal of high mannose N-glycans from glycoproteins

Description: Endoglycosidase H is a recombinant glycosidase which cleaves within the chitobiose core of high mannose and some hybrid oligosaccharides from N-linked glycoproteins.

Endo H₁ is a recombinant protein fusion of Endoglycosidase H and maltose binding protein. It has identical activity to Endo H.

Source: Endo H and Endo H₁ have been cloned from *Streptomyces plicatus* and overexpressed in *E. coli*.

Reaction Conditions:

Denature glycoprotein in 1X Glycoprotein Denaturing Buffer at 100°C for 10 minutes. Incubate in 1X GlycoBuffer 3 at 37°C. Heat inactivation: 65°C for 10 minutes.

Reagents Supplied:

10X Glycoprotein Denaturing Buffer
10X GlycoBuffer 3

Molecular Weight:

Endo H: 29,000 daltons
Endo H₁: 70,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to remove > 95% of the carbohydrate from 10 µg of denatured RNase B in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration:

Endo H concentration: 500,000 units/ml
Endo H₁ concentration: 1,000,000 units/ml


Usage Notes: Enzymatic activity is not affected by SDS.

To deglycosylate a native glycoprotein, longer incubation time as well as more enzyme may be required.

Endo S



#P0741S	6,000 units\$195
#P0741L	30,000 units\$777



- Removal of N-glycans from native IgG
- Useful for determining N-glycosylation sites

Description: Endo S is an endoglycosidase with a uniquely high specificity for removing N-linked glycans from the chitobiose core of the heavy chain of native IgG. Endo S is tagged with a chitin binding domain (CBD) for easy removal from a reaction and is supplied glycerol-free for optimal performance in HPLC- and MS-intensive methods.

Source: Endo S is cloned from *Streptococcus pyogenes* and overexpressed as a fusion to the chitin binding domain in *E. coli*.

Reagents Supplied:

10X GlycoBuffer 1

Molecular Weight: 136,000 daltons


Unit Definition: 5 µg of IgG in 1X GlycoBuffer 1 are incubated with two-fold dilutions of Endo S for 1 hour at 37°C. Separation of reaction products is visualized by SDS-PAGE.

Concentration: 200,000 units/ml

Heat Inactivation: 55°C for 10 minutes

PNGase A

#P0707S	150 units\$207
#P0707L	750 units\$1,080



PNGase A hydrolyzes *N*-glycan chains from glycoproteins/peptides regardless of the presence of xylose or fucose. [x = H or Man or GlcNAc].

- Removal of *N*-linked glycans from glycoproteins

Description: PNGase A is a recombinant amidase, which cleaves between the innermost GlcNAc and asparagine residues of high mannose, hybrid, and short complex oligosaccharides such as those found in plant and insect cells from *N*-linked glycoproteins and glycopeptides. PNGase A differs from PNGase F in that it cleaves *N*-linked glycans with or without $\alpha(1,3)$ -linked core fucose residues.

Source: Cloned from *Oryza sativa* (rice) and expressed in *Pichia pastoris*.

Reaction Conditions:

Denature 1 μ g of recombinant Avidin in 1X Glycoprotein Denaturing Buffer at 100°C for 10 minutes. After the addition of NP-40 and GlycoBuffer 3, two-fold dilutions of PNGase A are added and the reaction mix is incubated for 1 hour at 37°C. Heat inactivation: 65°C for 10 minutes.



Reagents Supplied:

10X Glycoprotein Denaturing Buffer
10X GlycoBuffer 3
10% NP-40

Molecular Weight: 63,800 daltons

Unit Definition: One unit is defined as the amount of enzyme required to remove > 95% of the carbohydrate from 1 μ g of denatured recombinant Avidin produced in Maize in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 5,000 units/ml

Usage Notes: PNGase A is active on both glycoproteins and glycopeptides.

PNGase A cannot cleave larger *N*-glycans such as those from Fetuin, Fibrinogen, IgG, Lactoferrin and Transferrin.

PNGase A is able to cleave high mannose *N*-glycan structures from Man 3 up to Man 9.


PNGase F & PNGase F, Recombinant

PNGase F		
#P0704S	15,000 units\$199
#P0704L	75,000 units\$792

PNGase F (Glycerol-free)		
#P0705S	15,000 units\$199
#P0705L	75,000 units\$792

PNGase F, Recombinant		
#P0708S	15,000 units\$163
#P0708L	75,000 units\$654

PNGase F (Glycerol-free), Recombinant		
#P0709S	15,000 units\$163
#P0709L	75,000 units\$654



PNGase F hydrolyzes nearly all types of *N*-glycan chains from glycopeptides/proteins [x = H or oligosaccharide].

- Removal of *N*-linked glycans from glycoproteins

Description: Peptide-*N*-Glycosidase F, also known as PNGase F, is an amidase which cleaves between the innermost GlcNAc and asparagine residues of high mannose, hybrid and complex oligosaccharides from *N*-linked glycoproteins. A glycerol-free version of PNGase F is also offered for HPLC methods.

Source: NEB #P0704 and NEB #P0705 are purified from *Flavobacterium meningosepticum*.

NEB #P0708 and NEB #P0709 are cloned from *Elizabethkingia miricola* (formerly *Flavobacterium meningosepticum*) and expressed in *E. coli*.

Reaction Conditions:

Denature glycoprotein in 1X Glycoprotein Denaturing Buffer at 100°C for 10 minutes. Heat inactivation: 75°C for 10 minutes.



Reagents Supplied:

10X Glycoprotein Denaturing Buffer
10X GlycoBuffer 2
10% NP-40

Molecular Weight: 36,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to remove > 95% of the carbohydrate from 10 μ g of denatured RNase B in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 500,000 units/ml

Usage Notes: Since PNGase F activity is inhibited by SDS, it is essential to have NP-40 present in the reaction mixture.

To deglycosylate a native glycoprotein, longer incubation time as well as more enzyme may be required.

After addition of NP-40 and Glycobuffer 2, two-fold dilutions of PNGase F are added and the reaction mix is incubated for 1 hour at 37°C.

Companion Products:

RNase B (control substrate)	
#P7817S	250 μ g \$69

Endoglycosidase Reaction Buffer Pack	
10X GlycoBuffer 2, 10X GlycoBuffer 3,	
10X Glycoprotein Denaturing Buffer,	
10% NP-40 (1 ml of each)	
#B0701S \$24

Rapid™ PNGase F & Rapid PNGase F (non-reducing format)




Rapid PNGase F
#P0710S 50 reactions\$457

Rapid PNGase F (non-reducing format)
#P0711S 50 reactions\$457

Companion Product:

Rapid PNGase F Antibody Standard
#P6043S 250 µg\$328



- Complete deglycosylation of antibodies and fusion proteins in minutes
- Release of all N-glycans rapidly and without bias
- Optimal activity is ensured for 12 months, if stored properly
- Purified to > 95% homogeneity, as determined by SDS-PAGE

Description: Rapid PNGase F is an improved reagent that allows the complete and rapid deglycosylation of antibodies and fusion proteins in minutes. All N-glycans are released rapidly and without bias, and are ready to be prepared for downstream chromatography or mass spectrometry analysis. Rapid PNGase F creates an optimized workflow which reduces processing time without compromising sensitivity or reproducibility.


Developed for proteomic applications, Rapid PNGase F (non-reducing format) enables complete and rapid deglycosylation while preserving disulfide bonds. This facilitates high throughput proteomics applications and methods for antibody characterization by mass spectrometry such as intact mass analysis. Rapid PNGase F (non-reducing format) combines the advantages of Rapid PNGase F (fast processing time), with non-reducing conditions, preserving quaternary structure.

Heat inactivation: 75°C for 10 minutes.


Reagents Supplied:

Rapid PNGase F
Rapid PNGase F Reaction Buffer (5X)
Rapid PNGase F (non-reducing format)
Rapid PNGase F (non-reducing format) Buffer (5X)

Specificity: Rapid PNGase F cleaves all complex, hybrid and high-mannose type glycans from antibodies and related proteins. Core α1-3 fucosylation (found in immunoglobulins expressed in plant or insect cells) is resistant to both PNGase F and Rapid PNGase F.



ESI-TOF analysis of an antibody before and after treatment with Rapid PNGase F.



ESI-TOF analysis of an antibody before and after treatment with Rapid PNGase F (non-reducing format).

Remove-iT® PNGase F




#P0706S 6,750 units\$198
#P0706L 33,750 units\$788

Companion Products:

Chitin Magnetic Beads
#E8036S 5 ml\$129

6-Tube Magnetic Separation Rack
#S1506S 6 tubes\$215

12-Tube Magnetic Separation Rack
#S1509S 12 tubes\$340



Remove-iT PNGase F hydrolyzes nearly all types of N-glycan chains from glycopeptides/proteins [x = H or oligosaccharide].

- Removal of N-linked glycans from glycoproteins

Description: Remove-iT PNGase F is an amidase which cleaves between the innermost GlcNAc and asparagine residues of high mannose, hybrid, and complex oligosaccharides from N-linked glycoproteins. Remove-iT PNGase F is tagged with a chitin binding domain (CBD) for easy removal from a reaction and is supplied glycerol free for optimal performance in HPLC and MS intensive methods.

Source: Purified from *Flavobacterium meningosepticum*.

Reaction Conditions:

Denature glycoprotein in 1X DTT (40 mM) at 55°C for 10 minutes. Incubate in 1X GlycoBuffer 2 for 1 hour at 37°C. Heat inactivation: 75°C for 10 minutes.

Reagents Supplied:

10X DTT
10X GlycoBuffer 2

Molecular Weight: 41,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to remove > 95% of the carbohydrate from 5 µg of DTT denatured RNase B in 1 hour at 37°C in a total reaction volume of 10 µl

Concentration: 225,000 units/ml

Usage Notes: Using typical RNase B denaturing conditions with NEB Glycoprotein Denaturing Buffer containing SDS and DTT, Remove-iT PNGase F yields a higher concentration of 500,000 U/ml.

If using Remove-iT PNGase F under typical PNGase F denaturing conditions, it is essential to have NP-40 in the reaction mixture, as Remove-iT PNGase F is inhibited by SDS. It is not known why this non-ionic detergent counteracts the SDS inhibition.

Removal of Remove-iT PNGase F from the deglycosylation reaction can be scaled up linearly with larger volumes of chitin magnetic beads.


O-Glycosidase

#P0733S	2,000,000 units\$141
#P0733L	10,000,000 units\$561

Companion Products:

O-Glycosidase & Neuraminidase Bundle		
#E0540S	1 bundle\$195

α 2-3,6,8 Neuraminidase		
#P0720S	2,000 units\$71
#P0720L	10,000 units\$284



- Removal of Core 1 and Core 3 O-linked disaccharide glycans from glycoproteins

Description: O-Glycosidase, also known as Endo- α -N-Acetylgalactosaminidase, catalyzes the removal of Core 1 and Core 3 O-linked disaccharides from glycoproteins.

Source: Cloned from *Enterococcus faecalis* and expressed in *E. coli*.

Reagents Supplied:

10X Glycoprotein Denaturing Buffer
10X GlycoBuffer 2
10% NP-40

Molecular Weight: 147,000 daltons



Unit Definition: One unit is defined as the amount of enzyme required to remove 0.68 nmol of O-linked disaccharide from 5 mg of neuraminidase-digested, non-denatured fetuin in 1 hour at 37°C in a total reaction volume of 100 μ l (1 unit of both O-Glycosidase and PNGase F will remove equivalent molar amounts of O-linked disaccharides and N-linked oligosaccharides, respectively).

Concentration: 40,000,000 units/ml

Heat Inactivation: 65°C for 10 minutes

Usage Note: Remove sialic acids if present

Protein Deglycosylation Mix II

#P6044S	20 reactions\$529
---------	--------------	------------

- Fast reaction setup
- Enzyme mixture ensures effective deglycosylation of N- and O-linked glycans
- Can be used under native and reducing conditions
- Enzymatic deglycosylation leaves intact core structures suitable for mass spectrometry analysis


Description: The Protein Deglycosylation Mix II contains all of the enzymes, reagents and controls necessary to remove almost all N-linked and simple O-linked glycans, as well as some complex O-linked glycans. This mix contains enzymes sufficient for 20 reactions or the cleavage of as much as 2 mg of glycoprotein.

Reagents Supplied:

Deglycosylation Mix Buffer 1 (10X)
Deglycosylation Mix Buffer 2 (10X)

Deglycosylation Enzyme Mix II:

PNGase F (Glycerol-free), Recombinant: 10,000 units/vial
O-Glycosidase: 80,000 units/vial
 α 2-3,6,8,9 Neuraminidase A: 400 units/vial
 β 1-4 Galactosidase S: 960 units/vial
 β -N-Acetylhexosaminidase: 300 units/vial



Enzymatic Deglycosylation of Bovine Fetuin under both native (10X Deglycosylation Mix Buffer 1) and reducing (10X Deglycosylation Mix Buffer 2) conditions. 20 μ g reactions were loaded onto a 10-20% Tris-glycine SDS-PAGE gel. Lane 1: Color Prestained Protein Standard, Broad Range (11-245 kDa), Lane 2: 20 μ g untreated Fetuin control, Lane 3: 20 μ g Fetuin deglycosylated under native conditions with Deglycosylation Mix Buffer 1, Lane 4: 20 μ g Fetuin deglycosylated under reducing conditions with Deglycosylation Mix Buffer 2, Lane 5: 5 μ l Protein Deglycosylation Mix II

Fetuin

#P6042S	500 μ g\$76
---------	-------------	-----------

Description: Fetuin is a glycoprotein containing sialylated N-linked and O-linked glycans that can be used as a positive control for endoglycosidase enzymes.

Source: Fetal Calf Serum


Molecular Weight: 64 kDa

Concentration: 10 mg/ml

Note: 500 μ g is enough for approximately 20 reactions. Due to heterogeneous glycosylation, Fetuin runs as a doublet on an SDS-PAGE gel.

α -N-Acetylgalactosaminidase

#P0734S 3,000 units\$141



Description: α -N-Acetylgalactosaminidase is a highly specific exoglycosidase that catalyzes the hydrolysis of α -linked D-N-Acetylgalactosamine residues from oligosaccharides and N-glycans attached to proteins.

Source: Cloned from *Chryseobacterium meningosepticum* and expressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 μ g/ml BSA. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

**Reagents Supplied:**

10X GlycoBuffer 1

100X BSA

Molecular Weight: 47,000 daltons


Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal α -D-N-Acetylgalactosamine from 1 nmol of (GalNAc α 1-3)(Fuc α 1-2)Gal β 1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 20,000 units/ml

β -N-Acetylglucosaminidase S

#P0744S 100 units\$136

#P0744L 500 units\$539



■ Removal of bisecting β -GlcNAc residues

Description: β -N-Acetylglucosaminidase S is a highly specific exoglycosidase that catalyzes the hydrolysis of terminal, non-reducing β -N-Acetylglucosamine residues from oligosaccharides.

Source: Cloned from *Streptococcus pneumoniae* and expressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Incubate at 37°C.

**Reagents Supplied:**

10X GlycoBuffer 1


Molecular Weight: 125,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal, non-reducing β -N-Acetylglucosamine from 1 nmol GlcNAc β 1-4GlcNAc β 1-4GlcNAc-AMC, in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 4,000 units/ml

β -N-Acetylhexosaminidase_f

#P0721S 500 units\$69



■ Active only on linear substrates

Description: β -N-Acetylhexosaminidase_f is a recombinant protein fusion of β -N-Acetylhexosaminidase and maltose binding protein with identical activity to β -N-Acetylhexosaminidase. It catalyzes the hydrolysis of terminal β -N-Acetylglactosamine and glucosamine residues from oligosaccharides.

Source: Cloned from *Streptomyces plicatus* and overexpressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Incubate at 37°C. Heat inactivation: 75°C for 10 minutes.

**Reagents Supplied:**

10X GlycoBuffer 1


Molecular Weight: 100,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal β -N-Acetylglactosamine from 1 nmol of GalNAc β 1-4Glc-AMC, in 1 hour at 37°C in 10 μ l volume.

Concentration: 5,000 units/ml

α 1-2 Fucosidase

#P0724S 1,000 units\$136



■ Active only on linear substrates

Description: α 1-2 Fucosidase is a highly specific exoglycosidase that catalyzes the hydrolysis of linear α 1-2 linked fucose residues from oligosaccharides. In this case, a linear substrate is defined as having no branching on the adjacent residue.

Source: Cloned from *Xanthomonas manihotis* and expressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 μ g/ml BSA. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

**Reagents Supplied:**

10X GlycoBuffer 1

100X BSA


Molecular Weight: 70,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the α -L-fucose from 1 nmol of Fuc α 1-2Gal β 1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 20,000 units/ml

α 1-3,4 Fucosidase

#P0769S 200 units\$134



Description: α 1-3,4 Fucosidase, (also known as AMF) is a broad specificity exoglycosidase that catalyzes the hydrolysis of α 1-3 and α 1-4 linked fucose residues from oligosaccharides and glycoproteins.

Source: Cloned from the sweet almond tree (*Prunus dulcis*) and expressed in *Pichia pastoris*.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 μ g/ml BSA. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

**Reagents Supplied:**

10X GlycoBuffer 1

100X BSA

Molecular Weight: 56,000 daltons


Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the α -fucose from 1 nmol of Gal β 1-4GlcNAc β 1-3(Fuc α 1-3)Gal β 1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 4,000 units/ml

α 1-2,3,4,6 Fucosidase

#P0748S 400 units\$147

#P0748L 2,000 units\$588



Description: α 1-2,3,4,6 Fucosidase is a broad specificity exoglycosidase that catalyzes the hydrolysis of α 1-2, α 1-3, α 1-4 and α 1-6 linked fucose residues from oligosaccharides.

Source: Cloned from bovine kidney and expressed in *E. coli*.

Reaction Conditions: 1X Glycobuffer 1. Supplement with 100 μ g/ml BSA. Incubate at 37°C. Heat inactivation: 100°C for 10 minutes.

**Reagents Supplied:**

10X Glycobuffer 1

100X BSA

Molecular Weight: 51,800 daltons


Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the α -L-fucose from 1 nmol of Fuc α 1-2Gal β 1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 8,000 units/ml

α 1-2,4,6 Fucosidase O

#P0749S 80 units\$134

#P0749L 400 units\$531



Description: α 1-2,4,6 Fucosidase O is a broad specificity exoglycosidase that catalyzes the hydrolysis of terminal α 1-2, α 1-4 and α 1-6 linked fucose residues from oligosaccharides. α 1-2,4,6 Fucosidase O cleaves α 1-6 fucose residues more efficiently than other linkages.

Source: Cloned from *Omnitrophica* bacterium and expressed in *E. coli*.

Reaction Conditions: 1X Glycobuffer 1. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

Reagents Supplied:

10X Glycobuffer 1


**Molecular Weight:** 49,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the fucose from 1 nmol of GOF from human IgG [GlcNAc β 1-2Man α 1-6(GlcNAc β 1-2Man α 1-3)Man β 1-4GlcNAc β 1-4GlcNAc(Fuc α 1-6)-AMAC], in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 2,000 units/ml

α 1-3,6 Galactosidase

#P0731S 100 units\$134



Description: α 1-3, 6 Galactosidase is a highly specific exoglycosidase that catalyzes the hydrolysis of α 1-3, 6 linked D-galactopyranosyl residues from oligosaccharides.

Source: Cloned from *Xanthomonas manihotis* and expressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 μ g/ml BSA. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes. Optimal incubation times and enzyme concentrations must be determined empirically for a particular substrate.

**Reagents Supplied:**

10X GlycoBuffer 1

100X BSA


Molecular Weight: 70,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal, α -D-galactose from 1 nmol Gal α 1-3Gal β 1-4Gal-AMC, in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 4,000 units/ml

α1-3,4,6 Galactosidase

#P0747S 200 units\$136
#P0747L 1,000 units\$539



Description: α1-3,4,6 Galactosidase is a highly specific exoglycosidase that catalyzes the hydrolysis of α1-3, α1-4 and α1-6 linked D-galactopyranosyl residues from oligosaccharides.

Source: Cloned from green coffee bean and expressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 µg/ml BSA. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

Optimal incubation times and enzyme concentrations must be determined empirically for a particular substrate.



Reagents Supplied:

10X GlycoBuffer 1
100X BSA


Molecular Weight: 39,700 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal, α-D-galactose from 1 nmol Galα1-3Galβ1-4Gal-AMC, in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 8,000 units/ml

β1-3 Galactosidase

#P0726S 500 units\$136



Description: β1-3 Galactosidase is a highly specific exoglycosidase that catalyzes the hydrolysis of β1-3 and, at a much lower rate, β1-6 linked D-galactopyranosyl residues from oligosaccharides. The approximate kinetic data show > 100-fold preference for β1-3 over β1-6 linkages and > 500-fold preference from β1-3 over β1-4 linkages.

Source: Cloned from *Xanthomonas manihotis* and expressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 µg/ml BSA. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.



Reagents Supplied:

10X GlycoBuffer 1
100X BSA


Molecular Weight: 66,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal β-D-galactose from 1 nmol of Galβ1-3GlcNAcβ1-3Galβ1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 10,000 units/ml

β1-3,4 Galactosidase

#P0746S 400 units\$136



Description: β1-3,4 Galactosidase, cloned from bovine testis and also known as BTG, is a highly specific exoglycosidase that catalyzes the hydrolysis of terminal β1-3 and β1-4 linked galactose residues from oligosaccharides.

Source: Cloned from bovine testis and expressed in *Pichia pastoris*.

Reaction Conditions: 1X GlycoBuffer 4. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.



Reagents Supplied:

10X GlycoBuffer 4


Molecular Weight: 71,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal, β-D-galactose from 1 nmol Galβ1-4GlcNAcβ1-3Galβ1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 8,000 units/ml

β1-4 Galactosidase S

#P0745S 400 units\$136
#P0745L 2,000 units\$539



Description: β1-4 Galactosidase S is a highly specific exoglycosidase that catalyzes the hydrolysis of β1-4 linked galactose residues from oligosaccharides.

Source: Cloned from *Streptococcus pneumoniae* and expressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Incubate at 37°C. Heat inactivation: 75°C for 10 minutes.



Reagents Supplied:

10X GlycoBuffer 1


Molecular Weight: 231,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal, β-D-galactose from 1 nmol Galβ1-4GlcNAcβ1-3Galβ1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 8,000 units/ml

α1-2,3 Mannosidase

#P0729S 640 units\$136



Description: α1-2,3 Mannosidase is a highly specific exoglycosidase that catalyzes the hydrolysis of α1-2 and α1-3 linked D-mannopyranosyl residues from oligosaccharides.

Source: Cloned from *Xanthomonas manihotis* and expressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 µg/ml BSA. Incubate at 37°C. Heat inactivation: 55°C for 10 minutes.



Reagents Supplied:

10X GlycoBuffer 1
100X BSA


Molecular Weight: 90,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal α-D-Mannose from 1 nmol of Manα1-3Manβ1-4GlcNAc-AMC, in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 32,000 units/ml

α1-2,3,6 Mannosidase

#P0768S 80 units\$136



Description: α1-2,3,6 Mannosidase, cloned from Jack Bean, and also known as JBM, is a broad specificity exoglycosidase that catalyzes the hydrolysis of terminal α1-2, α1-3 and α1-6 linked mannose residues from oligosaccharides. α1-2,3,6 Mannosidase has a slight preference for α1-2 mannose residues over α1-3 and α1-6 mannose residues.

Source: Cloned from *Canavalia ensiformis* (Jack Bean) and expressed in *Pichia pastoris*.

Reaction Conditions: 1X GlycoBuffer 4 and 1X Zinc. Incubate at 37°C. Heat inactivation: 95°C for 10 minutes.



Reagents Supplied:

10X GlycoBuffer 4
10X Zinc


Molecular Weight: 110,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal mannose from 1 nmol of Man(α1,3)-Man(β1,4)-GlcNAc-AMC, in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 2,000 units/ml

α1-6 Mannosidase

#P0727S 800 units\$136



Description: α1-6 Mannosidase is a highly specific exoglycosidase that removes unbranched α1-6 linked D-mannopyranosyl residues from oligosaccharides. When used in conjunction with α1-2,3 Mannosidase, the α1-6 Mannosidase will cleave α1-6 Mannose residues from branched carbohydrate substrates.

Source: Cloned from *Xanthomonas manihotis* and expressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Supplement with 100 µg/ml BSA. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.



Reagents Supplied:

10X GlycoBuffer 1
100X BSA


Molecular Weight: 58,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal α-D-mannose from 1 nmol of Manα1-6Manα1-6Man-AMC, in 1 hour at 37°C in a total reaction volume of 10 µl.

Concentration: 40,000 units/ml

Note: p-nitrophenyl-α-D-mannopyranoside is NOT a substrate for this enzyme.

α2-3,6,8 Neuraminidase

#P0720S 2,000 units\$71
#P0720L 10,000 units\$284

■ Active from pH 4.5 to 8.5

Description: Neuraminidase is the common name for Acetyl-neuraminyl hydrolase (Sialidase). This Neuraminidase catalyzes the hydrolysis of α2-3, α2-6 and α2-8 linked N-acetylneuraminic acid residues from glycoproteins and oligosaccharides.

Source: Cloned from *Clostridium perfringens* and overexpressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

Reagents Supplied:

10X GlycoBuffer 1



Molecular Weight: 43,000 daltons


Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal α-Neu5Ac from 1 nmol of Neu5Acα2-3Galβ1-3GlcNAcβ1-3Galβ1-4Glc-AMC, in 5 minutes at 37°C in a total reaction volume of 10 µl.

Concentration: 50,000 units/ml

Note: This enzyme shows a preference for α2,3 and α2,6 linkages over α2,8 linkages.

α 2-3,6,8,9 Neuraminidase A

#P0722S 800 units \$71
 #P0722L 4,000 units \$284



- Removes branched sialic acid residues that are linked to an internal residue

Description: Neuraminidase is the common name for Acetyl-neuraminyl hydrolase (Sialidase). α 2-3,6,8,9 Neuraminidase A catalyzes the hydrolysis of all linear and branched non-reducing terminal sialic acid residues from glycoproteins and oligosaccharides. The enzyme releases α 2-3 and α 2-6 linkages at a slightly higher rate than α 2-8 and α 2-9 linkages.

Source: Cloned from *Arthrobacter ureafaciens* and expressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Incubate at 37°C. Heat inactivation: 65°C for 10 minutes.

RR 37° 65°

Reagents Supplied:

10X GlycoBuffer 1


Molecular Weight: 100,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal α -Neu5Ac from 1 nmol Neu5Ac α 2-3Gal β 1-3GlcNAc β 1-3Gal β 1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 μ l.

Concentration: 20,000 units/ml

α 2-3 Neuraminidase S

#P0743S 400 units \$69
 #P0743L 2,000 units \$282



Description: Neuraminidase is the common name for Acetyl-neuraminyl hydrolase (Sialidase). α 2-3 Neuraminidase S is a highly specific exoglycosidase that catalyzes the hydrolysis of α 2-3 linked *N*-acetylneuraminic acid residues from glycoproteins and oligosaccharides.

Source: Cloned from *Streptococcus pneumoniae* and expressed in *E. coli*.

Reaction Conditions: 1X GlycoBuffer 1. Incubate at 37°C. Heat inactivation: 75°C for 10 minutes.

RR 37° 75°

Reagents Supplied:

10X GlycoBuffer 1

Molecular Weight: 74,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the α -Neu5Ac from 1 nmol of Neu5Ac α 2-3Gal β 1-3GlcNAc β 1-3Gal β 1-4Glc-AMC, in 1 hour at 37°C in a total reaction volume of 10 μ l.


Concentration: 8,000 units/ml



Ann-Marie has been with NEB for 5 years as a Technical Account Manager for our subsidiary office in the UK. She enjoys working with customers across the south of London in both industry and academia. In her spare time, Ann-Marie likes to travel and spend time exploring the outdoors with her family.

Bacteroides Heparinase I

#P0735S 240 units\$160
 #P0735L 600 units\$319



- Degradation of heparin and heparan sulfate glycosaminoglycans

Description: *Bacteroides* Heparinase I cloned from *Bacteroides eggerthii*, also called Heparin Lyase I, is active on heparin and the highly sulfated domains of heparan sulfate. The reaction yields oligosaccharide products containing unsaturated uronic acids which can be detected by UV spectroscopy at 232 nm.

Source: Cloned from *Bacteroides eggerthii* and expressed in *E. coli*.

Reaction Conditions:

10 µg heparin substrate, 10 µl *Bacteroides* Heparinase Reaction Buffer and H₂O in a total reaction volume of 100 µl. Incubate reaction at 30°C.
 Heat inactivation: 100°C for 1 minute.



Reagents Supplied:

10X *Bacteroides* Heparinase Reaction Buffer


Molecular Weight: 42,000 daltons

Unit Definition: One unit is defined as the amount of enzyme that will liberate 1.0 µmol unsaturated oligosaccharides from porcine mucosal heparin per minute at 30°C and pH 7.0 in a total reaction volume of 100 µl. For unit assay conditions, visit www.neb.com.

Concentration: 12,000 units/ml

Bacteroides Heparinase II

#P0736S 80 units\$241
 #P0736L 200 units\$480



- Degradation of heparin and heparan sulfate glycosaminoglycans

Description: *Bacteroides* Heparinase II cloned from *Bacteroides eggerthii*, also called Heparin Lyase II, is active on heparin and heparan sulfate. The reaction yields oligosaccharide products containing unsaturated uronic acids which can be detected by UV spectroscopy at 232 nm.

Source: Cloned from *Bacteroides eggerthii* and expressed in *E. coli*.

Reaction Conditions:

10 µg heparin substrate, 10 µl *Bacteroides* Heparinase Reaction Buffer and H₂O in a total reaction volume of 100 µl. Incubate reaction at 30°C.
 Heat inactivation: 100°C for 1 minute.



Reagents Supplied:

10X *Bacteroides* Heparinase Reaction Buffer


Molecular Weight: 86,000 daltons

Unit Definition: One unit is defined as the amount of enzyme that will liberate 1.0 µmol unsaturated oligosaccharides from porcine mucosal heparin per minute at 30°C and pH 7.0 in a total reaction volume of 100 µl. For unit assay conditions, visit www.neb.com.

Concentration: 4,000 units/ml

Bacteroides Heparinase III

#P0737S 14 units\$323
 #P0737L 35 units\$645



- Degradation of heparan sulfate glycosaminoglycans

Description: *Bacteroides* Heparinase III cloned from *Bacteroides eggerthii*, also called Heparin Lyase III, is active on heparan sulfate. The reaction yields oligosaccharide products containing unsaturated uronic acids which can be detected by UV spectroscopy at 232 nm.

Source: Cloned from *Bacteroides eggerthii* and expressed in *E. coli*.

Reaction Conditions:

10 µg heparan sulfate substrate, 10 µl *Bacteroides* Heparinase Reaction Buffer and H₂O in a total reaction volume of 100 µl. Incubate reaction at 30°C.
 Heat inactivation: 100°C for 1 minute.



Reagents Supplied:

10X *Bacteroides* Heparinase Reaction Buffer

Molecular Weight: 75,000 daltons


Unit Definition: One unit is defined as the amount of enzyme that will liberate 1.0 µmol unsaturated oligosaccharides from heparan sulfate per minute at 30°C and pH 7.0 in a total reaction volume of 100 µl. For unit assay conditions, visit www.neb.com.

Concentration: 700 units/ml

Endoglycoceramidase I (EGCase I)



#P0773S 150 mU \$73



This is an **Enzyme for Innovation (EFI)**. EFI is a project initiated by NEB to provide unique enzymes to the scientific community in the hopes of enabling the discovery of new and innovative applications. Visit www.neb.com/EnzymesforInnovation to view the full list.

Description: Endoglycoceramidase I (EGCase I) catalyzes the hydrolysis of the β -glycosidic linkage between oligosaccharides and ceramides in various glycosphingolipids. One unit of *R. triatomea* EGCase I is defined as the amount of enzyme required to hydrolyze 1 μ mol of ganglioside GM1a per minute at 37°C.

Source: EGCase I is isolated from a strain of *E. coli*, which contains the cloned EGCase I gene from *Rhodococcus triatomea*.

Reaction Conditions: 1X EGCase I Reaction Buffer. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Reagents Supplied:
10X EGCase I Reaction Buffer

Molecular Weight: 50,000 daltons

Unit Definition: One unit is defined as the amount of enzyme required to hydrolyze 1 μ mol of ganglioside GM1a per minute at 37°C.

Concentration: 6 units/ml

IdeZ Protease (IgG-specific)



#P0770S 4,000 units \$431

human IgG1, IgG3, IgG4: CPAPPELLG ∇ GPSVF
human IgG2: CPAPPVA ∇ GPSVF
murine IgG2a: CPAPNLLG ∇ GPSVF
murine IgG3: CPPGNILG ∇ GPSVF

- Complete fragmentation of antibodies and immunoglobulin fusion proteins in 30 minutes under native conditions

Description: IdeZ Protease (IgG-specific) is a recombinant antibody specific protease cloned from *Streptococcus equi* subspecies *zooepidemicus* that recognizes all human, sheep, monkey, and rabbit IgG subclasses, specifically cleaving at a single recognition site below the hinge region, yielding a homogenous pool of F(ab')₂ and Fc fragments. IdeZ Protease more effectively cleaves murine IgG2a than IdeS.

Source: Cloned from *Streptococcus equi* subspecies *zooepidemicus* and expressed in *E. coli*.


Reaction Conditions:
1X GlycoBuffer 2. Incubate at 37°C.
Heat inactivation: 65°C for 10 minutes.

Reagents Supplied with Enzyme:
10X GlycoBuffer 2

Molecular Weight: 35,578 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of 1 μ g of human IgG, in 15 minutes at 37°C in a total reaction volume of 10 μ l.

Concentration: 80,000 units/ml




Digestion of IgG with IdeZ Protease (IgG-specific), followed by denaturation.

O-Glycoprotease



#P0761S 200 reactions \$695



- Efficiently cleaves glycoproteins with or without sialic acid; no neuraminidase treatment necessary
- 200 reactions is sufficient for the cleavage of up to 2 mg of glycoprotein

Description: O-Glycoprotease is a highly specific protease that cleaves the peptide bonds of a glycoprotein or glycopeptide immediately N-terminal to a serine or threonine residue containing a mucin-type O-linked glycan with or without sialylation.

Source: Cloned from *Pseudomonas aeruginosa* and expressed in *E. coli*.

Reaction Conditions:
20 mM TrisHCl, pH 8.0. Incubate at 37°C.
Heat inactivation: 95°C for 10 minutes.

Molecular Weight: 96,727 daltons

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 90% of 2 μ M FAM-labeled O-glycopeptide in 2 hours at 37°C in a total reaction volume of 20 μ l.

Concentration: 1,000 units/ml

Proteases

Proteins found in nature vary greatly in size from 5 kDa to greater than 400 kDa. While it is possible to study intact proteins and the modifications present on these proteins by mass spectrometry, the most common proteomic approaches utilize digestion with site-specific proteases to generate smaller fragments, called peptides, as a first step in the analyses. Peptides are easier to characterize and can be separated using reverse phase supports, such as C18, combined with high performance liquid chromatography (HPLC). HPLC coupled with tandem Mass Spectrometry is used to obtain fragmentation data of individual peptides.

Trypsin-ultra™, Mass Spectrometry Grade

#P8101S 100 µg \$84

Lys/Arg▼XXX

- Digestion of proteins for proteomic analysis by mass spectrometry
- Protein and peptide identification
- TPCK treatment eliminates chymotryptic activity
- Free of contaminating proteases

Description: Trypsin-ultra, Mass Spectrometry Grade is a serine endopeptidase. It selectively cleaves peptide bonds C-terminal to lysine and arginine residues. Trypsin-ultra, Mass Spectrometry Grade is treated with L-(tosylamido-2-phenyl) ethyl chloromethyl ketone (TPCK) to inactivate any remaining chymotryptic activity. It is modified by acetylation of the ε-amino groups of lysine residues to prevent autolysis. Trypsin-ultra, Mass Spectrometry Grade cleaves at Lys-Pro and Arg-Pro bonds at a much slower rate than other amino acid residues.

Source: Isolated from bovine (*Bos taurus*) pancreas.

Reaction Conditions:
1X Trypsin Ultra Reaction Buffer. Incubate at 37°C.

Reagents Supplied with Enzyme:
2X Trypsin-ultra Reaction Buffer

Molecular Weight: 23,675 daltons

Reconstitution: Trypsin-ultra, Mass Spectrometry Grade should be reconstituted by the addition of 20–200 µl of high purity water. Rapid autolysis is a function of enzyme concentration.

Notes: Can be stored frozen in solution at –20°C for up to 1 week. A decrease in activity will occur if stored in solution. Use only freshly reconstituted protease for best results.

NEW α-Lytic Protease

#P8113S 20 µg \$89
#P8113L 100 µg \$339

XX-T/A/S/V▼XX

- Analyze complex proteomes
- Suitable for both in-gel and solution digests
- Optimal activity and stability for up to 24 months
- Ideal for digestion of proteins for proteomic analysis by mass spectrometry

Description: α-Lytic Protease (aLP) cleaves after Threonine (T), Alanine (A), Serine (S) and Valine (V) residues. Its specificity makes it an orthogonal and alternative protease to others commonly used in proteomics applications, including trypsin and chymotrypsin. Peptides generated by aLP are of similar average length to those of Trypsin.

Source: Purified from *Lysobacter enzymogenes*.

Molecular Weight: 19,800 daltons

Concentration: 0.4 mg/ml

Notes: α-Lytic Protease is stable for at least 2 years at –20°C. No loss of activity is observed after 10 freeze-thaw cycles. To avoid autolysis, α-Lytic Protease should be stored long term in 10 mM Sodium Acetate pH 5.0.

John has been with NEB for 10 years, and serves as Packaging Supervisor for NEB's OEM and Customized Solutions Team. His strict attention to detail ensures that your custom product will be packaged quickly and according to your requirements. When not at work, John enjoys boating and coaching youth sports.



Endoproteinase LysC

#P8109S 20 µg\$164

XX-Lys[▼]XXX

- Ideal for proteomic analysis by mass spectrometry
- Free of contaminating proteases
- Best suited for peptide identification

Description: LysC is a serine endoproteinase, isolated from *Lysobacter enzymogenes*, that cleaves peptide bonds C-terminal to lysine residues. LysC is a sequencing grade enzyme and is suitable for proteomics and glycobiology applications.

Source: Isolated from *Lysobacter enzymogenes*.

Molecular Weight: 30,000 daltons

37°

Reconstitution: Endoproteinase LysC should be reconstituted in 200 µl double-distilled water to make a 100 ng/µl solution in 10 mM Tris-HCl, pH 8.0. Rapid autolysis is a function of enzyme concentration. The solution can be stored at 4°C for several days or in single-use aliquots at -20°C for several months.

Note: LysC is compatible for simultaneous co-digestion with Trypsin-ultra, Mass Spectrometry Grade (NEB #P8101).

Endoproteinase GluC

#P8100S 50 µg\$85

XX-Glu[▼]XX

- Ideal for proteomic analysis by mass spectrometry
- Protein & peptide identification
- Free of contaminating proteases. Produced from a protease-deficient *Bacillus subtilis* strain

Description: Endoproteinase GluC (*Staphylococcus aureus* Protease V8) is a serine proteinase which selectively cleaves peptide bonds C-terminal to glutamic acid residues. Endoproteinase GluC also cleaves at aspartic acid residues at a rate 100–300 times slower than at glutamic acid residues.

Source: *Staphylococcus aureus* Protease V8 gene cloned and expressed in *Bacillus subtilis*.

Reaction Conditions:
1X GluC Reaction Buffer. Incubate at 37°C.

RM 37°

Reagents Supplied with Enzyme:

2X GluC Reaction Buffer

Molecular Weight: 29,849 daltons

Reconstitution: Endoproteinase GluC should be reconstituted by the addition of 50–500 µl of high purity water. Rapid autolysis is a function of enzyme concentration.

Note: Can be stored frozen in solution at -20°C for up to 2 weeks. A decrease in activity will occur if stored in solution. Use only freshly reconstituted protease for best results.

Endoproteinase AspN

#P8104S 50 µg\$84

XX[▼]Asp-XXX

- Ideal for proteomic analysis by mass spectrometry
- Free of contaminating proteases
- Best suited for peptide identification

Description: Endoproteinase AspN (flavastacin) is a zinc metalloendopeptidase which selectively cleaves peptide bonds N-terminal to aspartic acid residues.

Source: Purified from *Flavobacterium meningosepticum*.

Reaction Conditions:
1X AspN Reaction Buffer. Incubate at 37°C.

Reagents Supplied with Enzyme:
2X AspN Reaction Buffer

Molecular Weight: 40,089.9 daltons

RM 37°

Reconstitution: Endoproteinase AspN should be reconstituted by the addition of 50–500 µl of high purity water. Rapid autolysis is a function of enzyme concentration.

Notes: Can be stored frozen in solution at -20°C for up to 2 weeks. A decrease in activity will occur if stored in solution. Use only freshly reconstituted protease for best results.


Trypsin-digested BSA MS Standard (CAM Modified)

#P8108S 500 pmol\$84

- Standardization range: 500 to 2400 Da

Description: A complex mixture of peptides produced by Trypsin digestion of Bovine Serum Albumin (BSA) that was reduced and alkylated with Iodoacetimide (CAM modified). This peptide mixture can be used to test a Matrix-Assisted Laser Desorption/Ionization Time-Of-Flight (MALDI-TOF) or Electrospray Ionization (ESI) mass spectrometer (TOF, Q-TOF or Ion Trap).

Source: BSA (GENBANK P02769) was digested using Trypsin-ultra, Mass Spectrometry Grade, (NEB #P8101).



One hundred fmol of resuspended peptide mix may be analyzed by reverse phase liquid chromatography with on-line MS/MS analysis, for example with a Proxeon EASY-nLC and by Orbitrap Mass Spectrometer. Both analytical methods reveal a range of peptides in the standard. At least sixty percent sequence coverage is seen after database search, with greater than 15 unique peptides being identified.

Proteinase K, Molecular Biology Grade

#P8107S 2 ml\$85

- Isolation of plasmid and genomic DNA
- Isolation of RNA
- Inactivation of RNases, DNases and enzymes in reactions

Description: Proteinase K is a subtilisin-related serine protease that will hydrolyze a variety of peptide bonds.

Source: *Engyodontium album* (*Tritirachium album*)

Reaction Conditions: Proteinase K is active in a wide range of buffers including all NEB specific restriction endonuclease buffers. It is highly active between pH 7.5 and 12.0 and temperatures between 20 and 60°C. Proteinase K is also active in chelating agents such as EDTA and activity is stimulated in up to 2% SDS or 4 M urea.

Calcium is important for thermostability of Proteinase K but it is not required for catalysis, therefore Proteinase K is also active in buffers containing chelating agents such as EDTA.

Molecular Weight: 28,900 daltons

Unit Definition: One unit will digest urea-denatured hemoglobin at 37°C (pH 7.5) per minute to produce equal absorbance as 1.0 μmol of L-tyrosine using Folin & Ciocalteu's phenol reagent.

Concentration: 800 units/ml

Thermolabile Proteinase K

#P8111S 30 units\$157

- Heat inactivated following incubation at 55°C for 10 minutes
- Isolation of plasmid and genomic DNA
- Inactivation of RNases, DNases and enzymes in reactions
- Removal of enzymes from DNA to improve cloning efficiency
- PCR purification

Description: Thermolabile Proteinase K is an engineered, subtilisin-related serine protease that will hydrolyze a variety of peptide bonds.

Source: Cloned from *Engyodontium album* (formerly *Tritirachium album*), mutagenized to increase thermolability of the enzyme and expressed in *K. lactis*.

Molecular Weight: 29,000 daltons



Reaction Conditions: Thermolabile Proteinase K is active in a wide range of buffers. It is highly active between pH 7.0-9.5 and temperatures 20-40°C. It is active in chelating agents such as EDTA up to 10 mM.

Unit Definition: One unit is defined as the amount of enzyme required to release 1.0 μmol of 4-nitroaniline per minute from N-Succinyl-Ala-Ala-Pro-Phe-p-nitroanilide at 25°C, in a total reaction volume of 105 μl.

Concentration: 120 units/ml

TEV Protease

#P8112S 1,000 units\$105

ENLYFQ▼(G/S/M)

- Removal of affinity purification tags such as MBP or poly-histidine from fusion proteins
- Contains a His-tag for easy removal from a reaction using NEBExpress Ni Resin (NEB #S1428), NEBExpress NiSpin Columns (NEB #S1427) or NEBExpress Ni-NTA Magnetic Beads (NEB #S1423)

Description: TEV Protease, also known as Tobacco Etch Virus (TEV) Protease, is a highly specific cysteine protease that recognizes the amino-acid sequence Glu-Asn-Leu-Tyr-Phe-Gln-(Gly/Ser/Met) and cleaves between the Gln and Gly/Ser/Met residues. It is often used for the removal of affinity purification tags such as maltose-binding protein (MBP) or poly-histidine from fusion proteins. TEV Protease has a 7xHis-tag for easy removal from a reaction using nickel affinity resins and has been engineered to improve thermal stability and decrease autolysis.

Source: Cloned from Tobacco Etch Virus and expressed in *E. coli*.

Molecular Weight: 28,000 daltons

Unit Definition: 1 unit of TEV Protease will cleave 2 μg of MBP-fusion protein, MBP5-TEV-paramyosin ΔSal, to 95% completion in a total reaction volume of 10 μl in 1 hour at 30°C in 50 mM Tris-HCl (pH 7.5 @ 25°C) with 0.5 mM EDTA and 1 mM DTT.

Concentration: 10,000 units/ml



Factor Xa Protease

#P8010S	50 µg	\$78
#P8010L	250 µg	\$312

Ile-Glu/Asp-Gly-Arg▼

Description: Factor Xa cleaves after the arginine residue in its preferred cleavage site Ile-(Glu or Asp)-Gly-Arg. It will sometimes cleave at other basic residues, depending on the substrate conformation. The most common secondary site, among those that have been sequenced, is Gly-Arg. There seems to be a correlation between proteins that are unstable in *E. coli* and those that are cleaved by Factor Xa at secondary sites; this may indicate that these proteins are in a partially unfolded state. Factor Xa will not cleave a site followed by proline or arginine.

Source: Purified from bovine plasma and activated by treatment with the activating enzyme from Russell's viper venom.

Molecular Weight: The predominant form of Factor Xa has a molecular weight of approximately 43,000 daltons, consisting of two disulfide-linked chains of approximately 27,000 daltons and 16,000 daltons. On SDS-PAGE, the reduced chains have apparent molecular weights of 30,000 daltons and 20,000 daltons.

Unit Definition: 1 µg of Factor Xa will cleave 50 µg of test substrate to 95% completion in 6 hours or less.

Concentration: 1 mg/ml

Removal: Factor Xa will bind specifically to benzamide-agarose.

Enterokinase, light chain

#P8070S	480 units	\$122
#P8070L	2,560 units	\$487

Asp-Asp-Asp-Lys▼

Description: Enterokinase is a specific protease that cleaves after the lysine at its cleavage site, Asp-Asp-Asp-Lys. It will sometimes cleave at other basic residues, depending on the conformation of the protein substrate. Enterokinase will not cleave a site followed by proline.

Source: Purified from *Pichia pastoris* containing a clone of the light chain of the bovine enterokinase gene.

Molecular Weight: The molecular weight of the light chain of enterokinase is 26,300 daltons. Its apparent molecular weight on SDS-PAGE is 31 kDa.



Unit Definition: One unit is defined as the amount of enzyme required to cleave 25 µg of a MBP-EK-paramyosin-ΔSal substrate to 95% completion in 16 hours at 25°C in a total reaction volume of 25 µl.

Concentration: 16,000 U/ml

Removal: Enterokinase will bind specifically to trypsin inhibitor agarose (e.g., Sigma T-0637).

Furin

#P8077S	50 units	\$150
#P8077L	250 units	\$599

Arg-X-X-Arg▼

Description: Furin is an ubiquitous subtilisin-like proprotein convertase. It is the major processing enzyme of the secretory pathway and is localized in the trans-golgi network. Substrates of Furin include blood clotting factors, serum proteins and growth factor receptors such as the insulin-like growth factor receptor. The minimal cleavage site is Arg-X-X-Arg▼. However, the enzyme prefers the site Arg-X-(Lys/Arg)-Arg▼. An additional arginine at the P6 position appears to enhance cleavage. Furin is inhibited by EGTA, α1-Antitrypsin Portland and polyarginine compounds.

Note: The ability to cleave a particular substrate appears to depend on its tertiary structure as well as on the amino acids immediately surrounding the cleavage site.



Source: Isolated from *Spodoptera frugiperda* (Sf9) cells infected with recombinant baculovirus carrying truncated human furin (kindly provided by R. Fuller).

Molecular Weight: The calculated molecular weight of truncated human furin is 52,700 daltons. Its apparent molecular weight in SDS-PAGE gels is 57,000 daltons.

Unit Definition: One unit is defined as the amount of furin required to cleave 25 µg of a MBP-FN-paramyosin-ΔSal substrate to 95% completion in 6 hours at 25°C in a total reaction volume of 25 µl.

Concentration: 2,000 units/ml

Lambda Protein Phosphatase (Lambda PP)



#P0753S	20,000 units\$158
#P0753L	100,000 units\$632

Companion Products:

p-Nitrophenylphosphate (PNPP)

Non-specific substrate for protein, alkaline and acid phosphatases

#P0757S	1 ml\$35
#P0757L	5 ml\$138

Sodium Orthovanadate (Vanadate)

General inhibitor for protein phosphotyrosyl specific phosphatases

#P0758S	1 ml\$33
#P0758L	5 ml\$126

Description: Lambda Protein Phosphatase (Lambda PP) is a Mn^{2+} -dependent protein phosphatase with activity towards phosphorylated serine, threonine and tyrosine residues. Lambda-PPase is active on phosphorylated histidine residues.

Source: Isolated from a strain of *E. coli* that carries the bacteriophage lambda ORF221 open reading frame (kindly provided by Dr. D. Barford).

Reaction Conditions:

1X NEBuffer for Protein MetalloPhosphatases (PMP). Supplement with 1 mM $MnCl_2$. Incubate at 30°C. Heat inactivation: 65°C for 1 hour in the presence of 50 mM EDTA.

Reagents Supplied with Enzyme:

10X NEBuffer for Protein MetalloPhosphatases (PMP)
10X $MnCl_2$ (10 mM)

Molecular Weight: 25,000 daltons

Unit Definition: One unit is defined as the amount of enzyme that hydrolyzes 1 nmol of *p*-Nitrophenyl Phosphate (50 mM, NEB #P0757) in 1 minute at 30°C in a total reaction volume of 50 μ l.

Concentration: 400,000 units/ml

Lea has been with NEB for 6 years as a Production Scientist in our DNA Polymerase Group. She is involved in all facets of the manufacturing process for some of our most popular products. From planning to production and functional testing, Lea helps to make sure these products are always available for our customers.



Protein Kinases

The reversible addition of phosphate groups to proteins is important for the transmission of signals within eukaryotic cells and, as a result, protein phosphorylation and dephosphorylation regulate many diverse cellular processes. As the number of known protein kinases has increased at an ever-accelerating pace, it has become more challenging to determine which protein kinases interact with which substrates in the cell. The determination of consensus phosphorylation site motifs by amino acid sequence alignment of known substrates has proven useful in this pursuit. These motifs can be helpful for predicting phosphorylation sites for specific protein kinases within a potential protein substrate.

Since the determinants of protein kinase specificity involve complex 3-dimensional interactions, these motifs, short amino-acid sequences describing the primary structure around the phosphoacceptor residue, are a significant oversimplification of the issue. They do not take into account possible secondary and tertiary structural elements, or determinants from other polypeptide chains or from distant locations within the same chain. Furthermore, not all of the residues described in a particular specificity motif may carry the same weight in determining recognition and phosphorylation by the kinase. As a consequence, they should be used with some caution.

On the other hand, many of the residues within these consensus sequences have in fact proven to be crucial recognition elements, and the very simplicity of these motifs has made them useful in the study of protein kinases and their substrates. In addition to the prediction of phosphorylation sites, short synthetic oligopeptides based on consensus motifs are often excellent substrates for protein kinase activity assays.

The table below summarizes the specificity motifs for protein kinases that are available from NEB. Phosphoacceptor residue is indicated in red, amino acids which can function interchangeably at a particular residue are separated by a slash (/), and residues which do not appear to contribute strongly to recognition are indicated by an "X".

PROTEIN KINASE	NEB #	RECOGNITION DETERMINANT	SIZE	PRICE
cAMP-dependent Protein Kinase (PKA), Catalytic Subunit	P6000S/L	R - R - X - S / T Y Y = hydrophobic residue	100,000/500,000 units	\$158/\$632
Casein Kinase II (CK2)	P6010S/L	S - X - X - E / D	10,000/50,000 units	\$158/\$632

D = aspartic acid, E = glutamic acid, R = arginine, S = serine, T = threonine, Y = tyrosine, X = any amino acid

Note: More specific information on recognition determinants for each kinase can be found on the corresponding product page at www.neb.com.

Companion Product:

Adenosine-5' Triphosphate (ATP)

#P0756S 1.0 ml \$37

#P0756L 5.0 ml \$143

Ph.D.[™] Phage Display Peptide Library Kits

Ph.D.-7 Kit
10 Panning Experiments
#E8100S \$566

Ph.D.-12 Kit
10 Panning Experiments
#E8110S \$576

Ph.D.-C7C Kit
10 Panning Experiments
#E8120S \$680

Components Sold Separately:

Ph.D.-12 Phage Display Peptide Library
50 Panning Experiments
#E8111L \$2,306

The Ph.D. Kits Include:

- Sufficient Phage Display Library for 10 separate panning experiments, complexity of 10^9 clones
- 28 gIII Sequencing Primer (100 pmol)
- 96 gIII Sequencing Primer (100 pmol)
- Host *E. coli* strain ER2738
- Control Target (Streptavidin) and Elutant (Biotin)
- Detailed Protocols

Description: Phage display describes a selection technique in which a library of peptide or protein variants is expressed on the outside of a phage virion, while the genetic material encoding each variant resides on the inside. This creates a physical linkage between each variant protein sequence and the DNA encoding it, which allows rapid partitioning based on binding affinity to a given target molecule (antibodies, enzymes, cell-surface receptors, etc.) by an *in vitro* selection process called panning. In its simplest form (Figure 1), panning is carried out by incubating a library of phage-displayed peptides with a plate (or bead) coated with the target, washing away the unbound phage, and eluting the specifically-bound phage. The eluted phage is then amplified and taken through additional binding/amplification cycles to enrich the pool in favor of binding sequences. After 3–4 rounds, individual clones are characterized by DNA sequencing and ELISA.

New England Biolabs offers 3 pre-made random peptide libraries, as well as the cloning vector M13KE for construction of custom libraries. The pre-made libraries consist of linear heptapeptide (Ph.D.-7) and dodecapeptide (Ph.D.-12) libraries, as well as a disulfide-constrained heptapeptide (Ph.D.-C7C) library. All of the libraries have complexities in excess of 2 billion independent clones. The randomized peptide sequences in all three libraries are expressed at the N-terminus of the minor coat protein pIII, resulting in a valency of 5 copies of the displayed peptide per virion.

The Ph.D. libraries have been used for myriad applications, including epitope mapping (Figure 2), identification of protein-protein contacts and enzyme inhibitors and discovery of peptide ligands for GroEL, HIV, semiconductor surfaces and small-molecule fluorophores and drugs.




Figure 1: Routine Phage Display Workflow. Round 1 : Incubate 10^{11} pfu Ph.D. library + target incubation, wash away non-binders, elute bound phage, enrich selected phage with amplification in *E. coli*. Carry out 3-4 rounds of selection and then proceed with sequencing and/or phage-ELISA.

β -endorphin:

Y	G	G	F	M	T	S	E	K	Q	T	P	...
---	---	---	---	---	---	---	---	---	---	---	---	-----

1st round sequences:

Y	G	W	I	S	P	L	H	L	P	T	
Y	Q	P	D	N	P	S	R	Q	I	A	N
Y	W	P	A	H	I	R	A	V	P	M	I
R	L	D	D	I	K	N	T	L	A	F	S
S	S	D	V	Y	S	L	Y	P	F	I	M
E	F	F	P	H	P	M	L	H	N	S	R
D	N	W	P	Y	R	P	S	F	S	L	S
S	H	N	T	Y	S	A	P	R	P	S	A
S	L	L	H	Y	A	S	S	L	S	L	M
F	N	Q	N	A	E	P	F	S	S	R	P
H	P	R	Q	L	L	H	H	P	L	S	P

2nd round sequences:

Y	G	G	F	L	I	G	L	Q	D	A	S
Y	G	G	F	H	Y	K	E	T	G	A	L
Y	Q	P	D	N	P	S	R	Q	I	A	N
V	Y	C	Y	I	N	Q	S	M	I	G	N
H	H	D	T	E	Y	R	T	T	Q	L	S
N	L	K	F	P	T	N	P	K	A	M	W
L	P	N	L	T	W	A	L	M	P	R	A
D	N	W	P	Y	R	P	S	F	S	L	S
S	H	N	T	Y	S	A	P	R	P	S	A
S	L	L	H	Y	A	S	S	L	S	L	M
V	T	M	N	T	K	T	P	G	P	M	P

3rd round sequences:

Y	G	G	F	M	T	T	P	S	H	V	P
Y	G	G	F	M	T	T	P	S	H	V	P
Y	G	G	F	I	S	Q	T	Q	H	Y	S
Y	G	G	F	I	S	Q	T	Q	H	Y	S
Y	G	G	F	G	N	S	L	V	M	P	V
Y	G	G	F	S	M	P	F	L	P	A	L
Y	G	A	F	D	V	T	T	G	V	T	S
Y	G	V	F	N	P	H	Y	L	P	S	L
A	P	S	T	D	K	Q	A	T	M	P	L
A	S	V	A	V	S	S	R	Q	D	A	A

Figure 2: Epitope mapping of an anti- β -endorphin monoclonal antibody with the Ph.D.-12 library.

The Ph.D.-12 library was panned against anti- β -endorphin antibody 3-E7 in solution (10 nM antibody), followed by affinity capture of the antibody-phage complexes onto Protein A-agarose (rounds 1 and 3) or Protein G-agarose (round 2). Bound phage were eluted with 0.2 M glycine-HCl, pH 2.2. Selected sequences from each round are shown aligned with the first 12 residues of β -endorphin; consensus elements are boxed. The results clearly show that the epitope for this antibody spans the first 7 residues of β -endorphin, and that the bulk of the antibody-antigen binding energy is contributed by the first 4 residues (YGGF), with some flexibility allowed in the third position. Additionally, the conserved position of the selected sequences within the 12 residue window indicates that the free α -amino group of the N-terminal tyrosine is part of the epitope.

Companion Product:**Anti-M13 pIII Monoclonal Antibody**

Anti-M13 pIII Monoclonal Antibody (mouse isotype IgG2a) is derived from the A23 hybridoma produced by the fusion of mouse myeloma cells and splenocytes from BALB/c mice immunized with the C-terminal half of M13 coat protein III (residues 259–406).

Note: ELISA against whole phage using this antibody is not recommended since the epitope lies in a region that is not accessible on intact M13 phage virions.

#E8033S 0.1 ml\$290

Ph.D. Peptide Display Cloning System

Ph.D. Peptide Display Cloning System


20 μ g M13KE gIII Cloning Vector

16 μ g Extension Primer

#E8101S\$189

Description: The Ph.D. Peptide Display Cloning System facilitates the display of custom peptide libraries on the surface of bacteriophage M13 as coat protein fusions, creating a physical linkage between each displayed peptide and its encoding DNA sequence. Peptide ligands for a variety of targets can then be selected by the straightforward method of panning. The supplied display vector M13KE is an M13 derivative with cloning sites engineered for N-terminal pIII fusion, resulting in a valency of 5 displayed peptides per virion. The use of a phage vector, rather than a phagemid, simplifies the intermediate amplification steps, since neither antibiotic selection nor helper phage superinfection are required. Since displayed proteins longer than 20–30 amino acids have a deleterious effect on the infectivity function of pIII

in phage vectors, **this vector is suitable only for the display of short peptides.** Included with the cloning vector is an insert extension primer as well as a detailed protocol for cloning a peptide library into M13KE.







Promoting Environmental Sustainability

Don Comb was certainly a lover of nature. Growing up, one of his hobbies was entomology, and a formative insect-related observation contributed to his life-long commitment to protecting nature. As a child, Don noticed his parent's car windshield splattered with insects at the end of a journey. However, as an adult, he saw a marked reduction in the insect population, and knew this was a reflection of the health of entire ecosystems. Don was not one to passively make note of a less-than-ideal situation — instead he was a person who would ask big questions, and then search for solutions that he could implement. He did this in both his personal and professional life.

When Don selected a new location for NEB's headquarters in Ipswich, MA, his priority was to minimize and, where possible, mitigate the impact NEB had on the environment. As such, the company implemented sustainable solutions in every aspect of its operations. The campus lies on a historic property comprised of an idyllic 160-acre expanse of mixed forest, wetlands and agricultural land. The beautiful grounds include public trails for hiking, skiing, bird-watching and horseback riding, and is home to a vast array of wildlife species.

In 2003, NEB hosted a competition, inviting architects to submit innovative designs for a 140,000 sq. ft state-of-the-art laboratory facility. The only rule was that they had to design the building around three Copper Beech trees. In fact, NEB installed underground retaining walls to protect the Copper Beech tree's root system to prevent compaction of surrounding soil by heavy construction machinery as the facility was being built.

The laboratory facility is LEED® (Leadership in Energy and Environmental Design) and ISO 14001 certified, which dictates a framework of strict, environmentally sound regulations, including a specific percentage of recycled materials used in construction and limits as to how far those materials should be transferred. No exotic woods were used — only common, local wood such as Maple, White Oak and domestic Ash. The main hallway is mostly constructed of glass, optimizing use of natural light and enabling a view of the majestic surroundings from almost anywhere.

The campus is situated within the Ipswich River Watershed, a complex and sensitive ecosystem that provides drinking water to 350,000 people and thousands of businesses in 14 communities. All of NEB's water needs are pulled from the Ipswich River, and we felt it was important to protect this resource. For this reason, another pioneering sustainable feature was installed on campus — a Solar Aquatics Wastewater System® that processes wastewater through a series of physical and biological filters, essentially utilizing and accelerating processes found naturally in streams and wetlands. Water is discharged from the facility back into the groundwater cleaner than when it came in.


Don believed that all companies should be responsible for protecting the environment. This commitment helped to make NEB a pioneer in eco-friendly building and operating practices, and set a precedent for other biotechnology companies to follow. As NEB continues to expand, we choose to make decisions with the mindset of protecting the planet for future generations.

View of the Copper Beech Tree
adjacent to NEB's laboratory facility.
Credit: Mileidy Rodriguez, NEB

Take a tour of the
NEB campus.



Epigenetics



Simplify your epigenetics research.

Epigenetics is the study of heritable changes in the phenotype of a cell or organism that are not encoded in the DNA of the genome. The molecular basis of an epigenetic profile arises from covalent modifications of the protein and DNA components of chromatin. The epigenetic profile of a cell often dictates cell memory and cell fate and, thus influences mammalian development.

The epigenetic code is hypothesized to be the combined effects of histone modifications and DNA methylation on gene expression. While the genetic code for an individual is the same in every cell, the epigenetic code is tissue- and cell-specific, and may change over time as a result of aging, disease or environmental stimuli (e.g., nutrition, life style, toxin exposure). Cross-talk between histone modifications, DNA methylation or RNAi pathways are being studied in such areas as cancer, X chromosome inactivation, and imprinting.

For over 40 years, New England Biolabs has been committed to understanding the mechanisms of restriction and methylation of DNA. This expertise in enzymology has led to the development of a suite of validated products for epigenetics research. These unique solutions to study DNA and histone modifications are designed to address some of the challenges of the current methods. All NEB products pass stringent quality control assays to ensure the highest level of functionality and purity.







Find an
interactive tutorial
on epigenetics.







Featured Products

- 278** EpiMark 5-hmC and 5-mC Analysis Kit
- 279** EpiMark Methylated DNA Enrichment Kit
- 283** EpiMark Nucleosome Assembly Kit
- 289** NEBNext Enzymatic Methyl-seq






Featured Tools & Resources

-  Videos of NEB Scientists Discussing Epigenetics
-  Epigenetics-related FAQs
-  Feature Articles
-  Visit www.EpiMark.com to view an interactive tutorial explaining the phenomenon of epigenetics at the molecular level.





Epigenetics

 EpiMark 5-hmC and 5-mC Analysis Kit	278
 T4 Phage β -glucosyltransferase	278
 EpiMark Methylated DNA Enrichment Kit	279
 EpiMark N6-Methyladenosine Enrichment Kit	220, 279
 EpiMark Bisulfite Conversion Kit	279
 EpiMark Hot Start <i>Taq</i> DNA Polymerase	71, 279












Methylation-Dependent Restriction Enzymes

 <i>Aba</i> SI	280
 <i>Fsp</i> EI	280
 <i>Lpn</i> PI	280
 <i>Msp</i> JI	280
 <i>Mcr</i> BC	281














Additional Restriction Enzymes for Epigenetic Analysis

 <i>Dpn</i> I	36
 <i>Dpn</i> II	36
 <i>Hpa</i> II	40
 <i>Msp</i> I	43

Histones

 EpiMark Nucleosome Assembly Kit	283
 Nucleosome Control DNA	283
 Histone H3.1/H4 Tetramer Human, Recombinant	283
 Histone H2A/H2B Dimer Human, Recombinant	283
 Histone H1 ^o Human, Recombinant	283
 Histone H2A Human, Recombinant	284
 Histone H2B Human, Recombinant	284
 Histone H3.1 Human, Recombinant	284
 Histone H3.2 Human, Recombinant	285
 Histone H3.3 Human, Recombinant	285
 Histone H4 Human, Recombinant	285

DNA Methyltransferases


 CpG Methyltransferase (M.SssI)	124
 GpC Methyltransferase (M.CviPI)	124
 <i>Alu</i> I Methyltransferase	124
 <i>Bam</i> HI Methyltransferase	124
 <i>dam</i> Methyltransferase	125
 <i>Eco</i> GII Methyltransferase	125
 <i>Eco</i> RI Methyltransferase	125
 <i>Hae</i> III Methyltransferase	125
 <i>Hha</i> I Methyltransferase	125
 <i>Hpa</i> II Methyltransferase	126
 <i>Msp</i> I Methyltransferase	126
 <i>Taq</i> I Methyltransferase	126
 Human DNA (cytosine-5) Methyltransferase (Dnmt1)	286

Control DNA

HeLa Genomic DNA	287
5-methyl-dCTP	287

NEBNext Reagents for ChIP-Seq Library Prep

NEBNext Ultra II DNA Library Prep Kit for Illumina	149
NEBNext Ultra II DNA Library Prep with Sample Purification Beads	149
NEBNext Multiplex Oligos	160
NEBNext Ultra II End Repair/da-Tailing Module	150
NEBNext Ultra II Ligation Module	150
NEBNext Ultra II Q5 Master Mix	150
NEBNext Enzymatic Methyl-seq Kit	151, 289
NEBNext Enzymatic Methyl-seq Conversion Module	151, 289
NEBNext Q5U Master Mix	289
NEBNext Multiplex Oligos for Enzymatic Methyl-seq (Unique Dual Index Primer Pairs)	289

 Recombinant Enzyme

EpiMark® 5-hmC and 5-mC Analysis Kit



#E3317S 20 reactions\$270

- Reproducible quantitation of 5-hmC and 5-mC within a specific loci
- Easy-to-follow protocols
- Compatible with existing techniques (PCR)
- Amenable to high throughput


Visit EpiMark.com to view a video tutorial for this kit.

Description: The EpiMark 5-hmC and 5-mC Analysis Kit is a simple and robust method for the identification and quantitation of 5-hydroxymethylcytosine (5-hmC) and 5-methylcytosine (5-mC) within a specific DNA locus. This enzymatic approach utilizes the differential methylation sensitivity of the isoschizomers MspI and HpaII in a simple 3-step protocol.

Briefly, genomic DNA is treated with T4 phage β -glucosyltransferase and UDP-glucose, which glucosylates all 5-hmC present. DNA is digested with MspI and HpaII, two isoschizomers with different methylation sensitivity. Endpoint or real time PCR can then be used to identify and quantitate the different methylation states. Designed to simplify methylation analysis, the EpiMark Kit expands the potential for new biomarker discovery.

The EpiMark 5-hmC and 5-mC Analysis Kit Includes:

- T4 Phage β -glucosyltransferase
- UDP-Glucose
- MspI
- HpaII
- Proteinase K
- Control DNA (unmodified, 5-mC and 5-hmC)
- Forward and reverse control primer mix
- NEBuffer 4



Analysis of the different methylation states in Balb/C mouse tissue samples (locus 12) using the EpiMark 5-hmC and 5-mC Analysis Kit. A) Endpoint PCR of the 6 different reactions needed for methylation analysis. The boxed lanes indicate the presence of 5-hmC. B) Real time PCR data was used to determine amounts of 5-hmC and 5-mC present. The results demonstrate a variation in 5-hmC levels in the tissue sources indicated.

T4 Phage β -glucosyltransferase



#M0357S 500 units\$92
#M0357L 2,500 units\$366

- Glucosylation of 5-hydroxymethylcytosine in DNA
- Immunodetection of 5-hydroxymethylcytosine in DNA
- Labeling of 5-hydroxymethylcytosine residues by incorporation of [3 H]- or [14 C]-glucose into 5-hmC-containing DNA acceptor after incubation with [3 H]- or [14 C]-UDP-Glc
- Detection of 5-hydroxymethylcytosine in DNA by protection from endonuclease cleavage

Description: T4 Phage β -glucosyltransferase specifically transfers the glucose moiety of uridine diphosphoglucose (UDP-Glc) to the 5-hydroxymethylcytosine (5-hmC) residues in double-stranded DNA, making beta-glucosyl-5-hydroxymethylcytosine.


Reaction Conditions: 1X NEBuffer 4 and 40 μ M UDP-Glucose. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Reagents Supplied with Enzyme:

- 10X NEBuffer 4
- 50X UDP-Glucose (2 mM)

Unit Definition: One unit is defined as the amount of enzyme required to protect 0.5 μ g T4gt-DNA against cleavage by MfeI restriction endonuclease.

Concentration: 10,000 units/ml



Learn about 5-hmC detection in Balb/C brain tissue.




EpiMark Methylated DNA Enrichment Kit

#E2600S 25 reactions\$488

- *High-affinity binding provides greater sensitivity*
- *Elution in a small volume simplifies downstream applications*
- *Easy-to-follow protocol yields enriched fractions in less than 2 hours*
- *Enriched methylated DNA fragments can be easily ligated to double-stranded adaptors for next generation sequencing*
- *Highly pure product from a wide range of input DNA concentrations*

Description: The EpiMark Methylated DNA Enrichment Kit enables the enrichment of double-stranded CpG methylated DNA based on CpG methylation density. It utilizes the methyl-CpG binding domain of human MBD2a protein as a capture agent. The protein is fused to the Fc tail of human IgG1 (MBD2a-Fc), which is coupled to Protein A Magnetic Beads (MBD2a-Fc/Protein A Bead). This stable complex will selectively bind double-stranded methylated CpG containing DNA. The high binding affinity of the MBD2a-Fc coupled beads and optimized reagents increases sensitivity and accuracy. This kit contains all the individual components necessary to achieve enrichment in less than two hours using a four step process.

Enrichment Workflow



Epi

The EpiMark Methylated DNA Enrichment Kit Includes:

- MBD2-Fc protein
- Protein A Magnetic Beads
- Bind/wash buffer
- High-salt elution buffer
- Fragmented HeLa DNA
- Line element primers for methylated-locus controls
- RPL30 primers for unmethylated locus controls
- MirA locus control primers

EpiMark N6-Methyladenosine Enrichment Kit

#E1610S 20 reactions\$431

Description: The EpiMark N6-Methyladenosine Enrichment Kit can be used to enrich m6A modified RNA in immunoprecipitation protocols for downstream analysis by next-generation RNA sequencing or RT-qPCR. Modified RNA is isolated from a fragmented RNA sample by

Epi

See page 220 for more information.

binding to the N6-Methyladenosine antibody attached to Protein G Magnetic Beads. After multiple wash and clean-up steps, the enriched RNA is eluted in nuclease-free water and is ready for further analysis.

EpiMark Bisulfite Conversion Kit

#E3318S 48 reactions\$175

- *Complete conversion of unmodified cytosine to uracil*
- *Easy-to-follow protocol*
- *Reliable and consistent results*
- *Purification columns included*

Description: This technique can reveal the methylation status of every cytosine residue, and it is amenable to massively parallel sequencing methods. Bisulfite conversion involves the conversion of unmodified cytosines to uracil, leaving the modified bases 5-mC and 5-hmC. The EpiMark Bisulfite Conversion Kit is designed for the detection of methylated cytosine, using a series of alternating cycles of thermal denaturation, followed by incubation with sodium bisulfite. This kit includes all the reagents necessary for complete bisulfite conversion, including spin columns. Amplification of bisulfite-treated samples can then be performed using EpiMark Hot Start *Taq* DNA Polymerase.

Epi

The EpiMark Bisulfite Conversion Kit Includes:

- Sodium metabisulfite
- Solubilization buffer
- Desulphonation reaction buffer
- EpiMark spin columns with 2 ml collection tubes
- Binding buffer
- Wash buffer
- Elution buffer

EpiMark Hot Start *Taq* DNA Polymerase

#M0490S 100 reactions\$66
#M0490L 500 reactions\$265

Description: EpiMark Hot Start *Taq* DNA Polymerase is a mixture of *Taq* DNA Polymerase and a temperature sensitive, aptamer-based inhibitor. The inhibitor binds reversibly to the enzyme, inhibiting polymerase activity at temperatures below 45°C, but releases the enzyme during normal PCR cycling conditions. This permits assembly of reactions at room temperature. An advantage of the

Ri Epi

See page 71 for more information.

aptamer-based hot start mechanism is that it does not require a separate high temperature incubation step to activate the enzyme. The advanced aptamer-based hot-start activity coupled with the supplied optimized reaction buffer makes the EpiMark Hot Start *Taq* DNA Polymerase an excellent choice for use on bisulfite-converted DNA.

Methylation-Dependent Restriction Enzymes

- *Specificity to epigenetically-relevant DNA modifications (5-mC and 5-hmC)*
- *Easy-to-follow protocols*
- *Less harsh than bisulfite conversion*
- *Simplified data analysis*

The EpiMark Suite of products has been validated for use in epigenetics applications. Visit EpiMark.com for more information.

Many restriction enzymes are sensitive to DNA methylation states. Cleavage can be blocked or impaired when a particular base in the recognition site is modified. The MspJI family of restriction enzymes are dependent on methylation and hydroxymethylation for cleavage to occur (1). These enzymes excise 32-base pair DNA fragments containing a centrally located 5-hmC or 5-mC modified residue that can be extracted and sequenced. Due to the known position of this epigenetic modification, bisulfite conversion is not required prior to downstream analysis.


These EpiMark validated, methylation-dependent restriction enzymes expand the potential for mapping epigenetic modifications and simplify the study of DNA methylation. Additionally, they provide an opportunity to better understand the role of 5-hydroxymethylcytosine in the genome.

Reference:

(1) Cohen-Karni, D. et al. (2011) *PNAS*, 108, 11040–11045.

AbaSI

#R0665S 1,000 units \$113



*C = ghmC, hmC, mC or C

Description: AbaSI is a DNA modification-dependent endonuclease that recognizes 5-glucosylhydroxymethylcytosine (ghmC) in double-stranded DNA and cleaves 11–13 bases 3' from the modified C leaving a 2–3 base 3' overhang. The enzyme only cleaves if there is a G residue 20–23 nucleotides 3' from the modified C. AbaSI also recognizes 5-hydroxymethylcytosine (hmC) at a much lower efficiency. It does not recognize DNA with 5-methylcytosine (mC) or unmodified cytosine.

rCutSmart   Epi  25° 


Reaction Conditions: 1X rCutSmart Buffer + 1 mM DTT (supplied), 25°C. Heat inactivation: 65°C for 20 minutes.

Concentration: 10,000 units/ml

Note: This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

FspEI

#R0662S 200 units \$112



Description: FspEI is a modification-dependent endonuclease which recognizes C^mC sites and generates a double-stranded DNA break on the 3' side of the modified cytosine at N₁₂/N₁₆. Recognized cytosine modifications include C5-methylation (5-mC) and C5-hydroxymethylation (5-hmC).

Reaction Conditions: 1X rCutSmart Buffer + 1X Enzyme Activator Solution (supplied), 37°C. Heat inactivation: 80°C for 20 minutes.

rCutSmart   Epi  37° 


Concentration: 5,000 units/ml

Note: Star activity may result from extended digestion.

This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

LpnPI

#R0663S 200 units \$112



Description: LpnPI is a modification-dependent endonuclease which recognizes C^mCDG sites and generates a double-stranded DNA break on the 3' side of the modified cytosine at N₁₀/N₁₄. Recognized cytosine modifications include C5-methylation (5-mC) and C5-hydroxymethylation (5-hmC).

Reaction Conditions: 1X rCutSmart Buffer + 1X Enzyme Activator Solution (supplied), 37°C. Heat inactivation: 65°C for 20 minutes.

rCutSmart   Epi  37° 

Concentration: 5,000 units/ml

Note: Star activity may result from extended digestion.

This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

MspJI

#R0661S 200 units \$112
#R0661L 1,000 units \$459



Description: MspJI is a modification-dependent endonuclease that recognizes C^mCNNR sites and generates a double-stranded DNA break on the 3' side of the modified cytosine at N₉/N₁₃. The recognized cytosine modifications include C5-methylation (5-mC) and C5-hydroxymethylation (5-hmC).

Reaction Conditions: 1X rCutSmart Buffer + 1X Enzyme Activator Solution (supplied), 37°C. Heat inactivation: 65°C for 20 minutes.

rCutSmart   Epi  37° 

Concentration: 5,000 units/ml

Note: Star activity may result from extended digestion.

This product will be supplied with Recombinant Albumin-containing buffer. For details, visit page 20 or www.neb.com/BSA-free.

McrBC



#M0272S	500 units	\$73
#M0272L	2,500 units	\$296

5'...Pu^mC (N₄₀₋₃₀₀₀) Pu^mC... 3'

- *Determination of the methylation state of CpG dinucleotides*
- *Detection of cytosine methylated DNA*

Description: McrBC is an endonuclease that cleaves DNA containing methylcytosine* on one or both strands. McrBC will not act upon unmethylated DNA. Sites on the DNA recognized by McrBC consist of two half-sites of the form (G/A)^mC. These half-sites can be separated by up to 3 kb, but the optimal separation is 55–103 base pairs. McrBC requires GTP for cleavage, but in the presence of a non-hydrolyzable analog of GTP, the enzyme will bind to methylated DNA specifically, without cleavage. McrBC will act upon a pair of Pu^mCG sequence elements, thereby detecting a high proportion of methylated CpGs, but will not recognize HpaII/MspI sites (CCGG) in which the internal cytosine is methylated.

* 5-methylcytosine, 5-hydroxymethylcytosine or N4-methylcytosine

Reaction Conditions: 1X NEBuffer 2 + 200 µg/ml BSA + 1 mM GTP, 37°C. Heat inactivation: 65°C for 20 minutes.

Reagents Supplied with Enzyme:
10X NEBuffer 2
100X BSA
100X GTP (100 mM)
Control Plasmid DNA

Unit Definition: One unit is defined as the amount of enzyme required to cleave 0.5 µg of a plasmid containing multiple McrBC sites in 1 hour at 37°C in a total reaction volume of 50 µl.

Concentration: 10,000 units/ml

Note: McrBC makes one cut between each pair of half-sites, cutting close to one half-site or the other, but cleavage positions are distributed over several base pairs approximately 30 base pairs from the methylated base. Therefore, when multiple McrBC half-sites are present in DNA (as is the case with cytosine-methylated genomic DNA), the flexible nature of the recognition sequence results in an overlap of sites, producing a smeared, rather than a sharp, banding pattern.

Additional Restriction Enzymes for Epigenetic Analysis

See pages 36, 40, 43 for more information.

Visit page 321 for more information about Dam, Dcm and CpG methylation.

Methylation sensitive restriction enzymes can be used to generate fragments for further analysis. When used in conjunction with an isoschizomer that has the same recognition site, but is methylation insensitive, information about methylation status can be obtained.

DpnI		
#R0176S	1,000 units	\$67
#R0176L	5,000 units	\$274
DpnII		
#R0543S	1,000 units	\$72
#R0543L	5,000 units	\$296
for high (5X) concentration		
#R0543T	1,000 units	\$72
#R0543M	5,000 units	\$296

HpaII		
#R0171S	2,000 units	\$67
#R0171L	10,000 units	\$274
for high (5X) concentration		
#R0171M	10,000 units	\$274
MspI		
#R0106S	5,000 units	\$67
#R0106L	25,000 units	\$274
for high (5X) concentration		
#R0106T	5,000 units	\$67
#R0106M	25,000 units	\$274


Histones

- Purification and characterization of enzymes that modify histone proteins
- Octamer and nucleosome modification studies
- Carrier chromatin immunoprecipitation (CChIP)
- High throughput studies

In eukaryotes, nuclear DNA is assembled into chromatin by nucleoprotein complexes. The primary unit of chromatin, a nucleosome core particle (NCP), is an octamer complex made up of two molecules each of Histone H2A, H2B, H3 and H4 and approximately 147 base pairs of nuclear DNA. Histone H1 further condenses the DNA by binding the linker segments between NCP complexes (1,2). Histones undergo diverse post-translational modification including acetylation, phosphorylation, mono-, di- or tri-methylation, ubiquitination, isomerization and ADP-ribosylation. Through their potential combinatorial sequences on a given histone and their reversibility, these modifications dynamically restrict or recruit numerous other proteins or protein complexes onto chromatin (3). The study of their roles in gene regulation (4), cellular stress events (4), aging and DNA repair (5) is revealing the multiple functions of histone modifications in the fate of a cell. Additional variability is incorporated into the system by histone variants. Acting individually or combinatorially in conjunction with DNA modification, histone modifications and histone variants are thought to establish an epigenetic code or epigenetic mechanism of gene regulation (3).

In total, seven human histones, including three histone H3 variants (see alignment below), have been individually cloned and expressed in *E. coli* and then highly purified from cell extracts at NEB. Mass spectrometry analysis demonstrates that these recombinant histones are free of post-translational modifications. These histones are ideal substrates for the purification and characterization of histone modifying enzymes.

To aid in studying intact nucleosomes, we now offer the EpiMark Nucleosome Assembly Kit. The precise mixing of a 2:1 ratio of Histone H2A/H2B Dimer to Histone H3.1/H4 Tetramer generates a recombinant human histone octamer, and in the presence of DNA forms nucleosomes (7,8). Enzymes that are unable to modify individual histones or DNA may be active on these nucleosome core particles, the histone dimer, or the histone tetramer (9,10). The NCPs also may be used as carrier chromatin in CChIP (carrier chromatin immunoprecipitation) assays (11). The recombinant human histone dimer and recombinant human histone tetramer are also available as separate products.



Experience the purity of Histones from NEB. SDS-PAGE analysis of the histones available from NEB.

- 1: Histone H1° (NEB #M2501) 1 µg
- 2: Histone H2A (NEB #M2502) 1 µg
- 3: Histone H2B (NEB #M2505) 1 µg
- 4: Histone H3.1 (NEB #M2503) 1 µg
- 5: Histone H3.2 (NEB #M2506) 1 µg
- 6: Histone H3.3 (NEB #M2507) 1 µg
- 7: Histone H4 (NEB #M2504) 1 µg
- 8: Histone H2A/H2B Dimer (NEB #M2508) 2 µg
- 9: Histone H3.1/H4 Tetramer (NEB #M2509) 2 µg
- 10: NEB Protein Ladder

References:

- (1) Kornberg, R.D. (1977) *Annu. Rev. Biochem.*, 46, 931–954.
- (2) van Holde, K.E. (1989) *Chromatin*, 1–497.
- (3) Kim, J.K., Samaranyake, M. and Pradhan, S. (2009) *Cell. Mol. Life Sci.*, 66, 596–612.
- (4) Huang, J. et al. (2006) *Nature*, 444, 629–632.
- (5) Pahlich, S., Zakaryan, R.P. and Gehring, H. (2006) *Biochim. Biophys. Acta.*, 1764, 1890–1903.
- (6) Hake, S.B. et al. (2006) *J. Biol. Chem.*, 281, 559–568.
- (7) Luger et al. (1999) *Methods in Mol. Biol.*, 119, 1–16.
- (8) Comb, D. and Mersha, F. unpublished observations.
- (9) Li, Yan et al. (2009) *J. Biol. Chem.*, 284, 34283–34295.
- (10) Fung, J. et al. (2002) *Current Biology*, 12, 1086–1099.
- (11) O'Neill et al. (2006) *Nature Genetics*, 38, 835–841.

	10	20	30	40	50	60	70	
ARTKQTARKS	TGGKAPRKQL	ATKAARKSAP	ATGGVKKPHR	YRPGTVALRE	IRRYQKSTEL	LIRKLPPQRL		histone H3.1
ARTKQTARKS	TGGKAPRKQL	ATKAARKSAP	ATGGVKKPHR	YRPGTVALRE	IRRYQKSTEL	LIRKLPPQRL		histone H3.2
ARTKQTARKS	TGGKAPRKQL	ATKAARKSAP	STGGVKKPHR	YRPGTVALRE	IRRYQKSTEL	LIRKLPPQRL		histone H3.3
	80	90	100	110	120	130		
VREIAQDFKT	DLRFQSSAVM	ALQEACEAYL	VGLFEDTNLC	AIHAKRVTIM	PKDIQLARRI	RGERA		histone H3.1
VREIAQDFKT	DLRFQSSAVM	ALQEASEAYL	VGLFEDTNLC	AIHAKRVTIM	PKDIQLARRI	RGERA		histone H3.2
VREIAQDFKT	DLRFQSAATIG	ALQEASEAYL	VGLFEDTNLC	AIHAKRVTIM	PKDIQLARRI	RGERA		histone H3.3

Sequence alignment of Human Histone variants H3.1, H3.2 and H3.3. Human Histone H3.1, 3.2 and 3.3 vary by only a few amino acids (changes are highlighted in red), but are associated with different biological functions (6).

EpiMark Nucleosome Assembly Kit

#E5350S 20 reactions\$640

Components Sold Separately:

Histone H3.1/H4 Tetramer Human, Recombinant
#M2509S 1 nmol\$250

Histone H2A/H2B Dimer Human, Recombinant
#M2508S 2 nmol\$250

Nucleosome Control DNA
#N1202S 0.2 nmol\$124

- Highly pure, recombinant system
- Pre-formed histone dimer and tetramer complexes simplify octamer formation
- Components stable for one year
- Ideal for ChIP Assay, HAT Assay and enzyme modification assays (e.g., methylation studies)


Description: This kit contains the components necessary to form an unmodified recombinant human nucleosome using your own target DNA or the supplied control DNA. The protocol requires the mixing of already formed and purified recombinant human Histone H2A/H2B Dimer and Histone H3.1/H4 Tetramer in the presence of DNA in high salt, followed by dialysis down to low salt, to make nucleosomes. One tetramer associates with two dimers to form the histone octamer on the DNA, generating a nucleosome. A method for assaying nucleosome formation by gel shift assay is also provided. These nucleosomes may serve as a better substrate for enzymes that are inactive on the DNA or one of the core histones alone. Each described reaction creates nucleosomes from ~50 pmol of a 208 bp DNA and may be scaled depending on the experiment.

Histone H2A/H2B Dimer Human, Recombinant is generated by refolding the denatured, purified subunits H2A and H2B, followed by gel filtration. Histone H3.1/H4 Tetramer Human, Recombinant is generated by refolding the denatured, purified subunits H3.1 and H4, followed by gel filtration. Both the dimer and tetramer are highly pure and are available separately for histone modification studies.

The Nucleosome Control DNA is a 208 base pair fragment from *Lytechinus variegatus* 5SrDNA, and can be used for mononucleosome formation.

The EpiMark Nucleosome Assembly Kit Includes:

- Histone H2A/H2B Dimer
- Histone H3.1/H4 Tetramer
- Control DNA



Gel shift assay to visualize nucleosome formation.
Samples from nucleosome assembly reactions were run on 6% polyacrylamide gel in 0.5X TBE. Lane 1: Low Molecular Weight DNA Ladder (NEB #N3233). Lane 2: Nucleosome Control DNA. Lane 3: 0.5:1 ratio of Octamer* to DNA. Lane 4: 1:1 ratio of Octamer* to DNA.

*Octamer = 2:1 mix of Histone H2A/H2B Dimer and Histone H3.1/H4 Tetramer.


Histone H1⁰ Human, Recombinant

#M2501S 100 µg\$85

Description: Histone H1 acts on the linker region of polynucleosome DNA to condense the chromatin into structures of ~30 nm. It is not necessary for octamer or nucleosome core particle formation. Eight different Histone H1 proteins have been identified in the human genome. Histone H1⁰ is a non replication-dependent histone that is highly expressed in terminally differentiated cells.

Synonyms: Histone H1.0, Histone H1(0), Histone H1'

Concentration: 1 mg/ml



Mass Spectrometry Analysis of Histone H1⁰ Human, Recombinant.
The average mass calculated from primary sequence is 20,732 Da.

Histone H2A Human, Recombinant




#M2502S

100 µg \$85

Description: Histone H2A combines with Histone H2B to form the H2A-H2B heterodimer. Two H2A/H2B heterodimers interact with an H3/H4 tetramer to form the histone octamer. Histone H2A is also modified by various enzymes and can act as a substrate for them. These modifications have been shown to be important in gene regulation.

Concentration: 1 mg/ml



Mass Spectrometry Analysis of Histone H2A Human, Recombinant.
The average mass calculated from primary sequence is 13,990 Da.

Histone H2B Human, Recombinant




#M2505S

100 µg \$85

Description: Histone H2B combines with Histone H2A to form the H2A-H2B heterodimer. Two H2A/H2B heterodimers interact with an H3/H4 tetramer to form the histone octamer. Histone H2B is also modified by various enzymes and can act as a substrate for them. These modifications have been shown to be important in gene regulation.

Synonyms: Histone H2B/8, Histone H2B.1, Histone H2B-GL105

Concentration: 1 mg/ml



Mass Spectrometry Analysis of Histone H2B Human, Recombinant.
The average mass calculated from primary sequence is 13,789 Da.

Histone H3.1 Human, Recombinant



#M2503S


100 µg \$85

Description: Histone H3 combines with Histone H4 to form the H3/H4 tetramer. Two H2A/H2B heterodimers interact with an H3/H4 tetramer to form the histone octamer. It is also modified by various enzymes and can act as a substrate for them. These modifications have been shown to be important in gene regulation.

Histone H3.1, an H3 variant that has thus far only been found in mammals, is replication dependent and is associated with gene activation and gene silencing.

Synonyms: Histone H3/a

Concentration: 1 mg/ml



Mass Spectrometry Analysis of Histone H3.1 Human, Recombinant.
The average mass calculated from primary sequence is 15,273 Da.

Histone H3.2 Human, Recombinant




#M2506S 100 µg \$85

Description: Histone H3 combines with Histone H4 to form the H3/H4 tetramer. Two H2A/H2B heterodimers interact with an H3/H4 tetramer to form the histone octamer. It is also modified by various enzymes and can act as a substrate for them. These modifications have been shown to be important in gene regulation.

Histone H3.2, an H3 variant that is found in all eukaryotes except budding yeast, is replication dependent and is associated with gene silencing.

Synonyms: Histone H3/m, H3/o

Concentration: 1 mg/ml



Mass Spectrometry Analysis of Histone H3.2 Human, Recombinant.
The average mass calculated from primary sequence is 15,257 Da.

Histone H3.3 Human, Recombinant




#M2507S 100 µg \$85

Description: Histone H3 combines with Histone H4 to form the H3/H4 tetramer. Two H2A/H2B heterodimers interact with an H3/H4 tetramer to form the histone octamer. It is also modified by various enzymes and can act as a substrate for them. These modifications have been shown to be important in gene regulation.

Histone H3.3, an H3 variant that is found in all eukaryotes from yeast to human, is replication and cell cycle phase independent and is the most common H3 in non-dividing cells. It has been shown to be enriched in covalent modifications associated with gene activation.

Synonyms: Histone H3.3A, H3F3, H3.3B

Concentration: 1 mg/ml



Mass Spectrometry Analysis of Histone H3.3 Human, Recombinant.
The average mass calculated from primary sequence is 15,197 Da.

Histone H4 Human, Recombinant




#M2504S 100 µg \$85

Description: Histone H3 combines with Histone H4 to form the H3/H4 tetramer. Two H2A/H2B heterodimers interact with an H3/H4 tetramer to form the histone octamer. Histone H4 is also modified by various enzymes and can act as a substrate for them. These modifications have been shown to be important in gene regulation.

Synonyms: For HIST2H4 gene: H4/N, H4F2, H4FN

Concentration: 1 mg/ml




Mass Spectrometry Analysis of Histone H4 Human, Recombinant.
The average mass calculated from primary sequence is 11,236 Da.

DNA Methyltransferases



See pages 124-126 for more information.



CpG Methyltransferase (M.SssI)
 #M0226S 100 units\$75
 #M0226L 500 units\$302
 for high (5X) concentration
 #M0226M 500 units\$302

GpC Methyltransferase (M.CviPI)
 #M0227S 200 units\$78
 #M0227L 1,000 units\$314

AluI Methyltransferase
 #M0220S 100 units\$77

BamHI Methyltransferase
 #M0223S 100 units\$77

dam Methyltransferase
 #M0222S 500 units\$79
 #M0222L 2,500 units\$318

EcoGII Methyltransferase
 #M0603S 200 units\$76

EcoRI Methyltransferase
 #M0211S 10,000 units\$69

HaeIII Methyltransferase
 #M0224S 500 units\$77

HhaI Methyltransferase
 #M0217S 1,000 units\$77

HpaII Methyltransferase
 #M0214S 100 units\$72


MspI Methyltransferase
 #M0215S 100 units\$77

TaqI Methyltransferase
 #M0219S 1,000 units\$77

Human DNA (cytosine-5) Methyltransferase (Dnmt1)



#M0230S 50 units\$130
 #M0230L 250 units\$520



Description: Dnmt1 methylates cytosine residues in hemimethylated DNA at 5'...CG...3'. Mammalian Dnmt1 is believed to be involved in carcinogenesis, embryonic development and several other biological functions. The bulk of the methylation takes place during DNA replication in the S-phase of the cell cycle.

Reaction Conditions: 1X Dnmt1 Reaction Buffer. Supplement with 100 µg/ml BSA and 160 µM S-adenosyl-methionine. Incubate at 37°C. Heat inactivation: 65°C for 20 minutes.

Reagents Supplied with Enzyme:

10X Dnmt1 Reaction Buffer
 100X BSA
 32 mM S-adenosylmethionine (SAM)

Unit Definition: One unit is the amount of enzyme required to catalyze the transfer of 1 pmol of methyl group to poly dI.dC substrate in a total reaction volume of 25 µl in 30 minutes at 37°C.

Concentration: 2,000 units/ml

Note: For DNA modification and protection applications, M.SssI (NEB #M0226) is preferred because it efficiently methylates both unmethylated and hemimethylated DNA substrates.

HeLa Genomic DNA

#N4006S 15 µg \$120

- PCR, SNP analysis and southern blotting
- Genomic DNA library construction
- Control DNA for Methylation-specific PCR (MSP), Bisulfite sequencing, Methylation-sensitive Single-Nucleotide Primer Extension (MS-SNuPE,) Combined Bisulfite Restriction Analysis (COBRA), Bisulfite treatment and PCR-single strand Confirmation Polymorphism Analysis (Bisulfite-PCR-SSCP/BI/PS)

Description: HeLa (cervix adenocarcinoma) cells were grown to confluency in DMEM plus 10% fetal bovine serum. Genomic DNA was isolated by a standard genomic purification protocol, phenol extracted and equilibrated to 10 mM Tris-HCl (pH 7.5) and 1 mM EDTA.

Note: NEB is one of the only suppliers of this product.

A_{260/280} Ratio: 1.87

5-methyl-dCTP

#N0356S 1 µmol in 0.1 ml \$73

- Generation of fully methylated cytosine-substituted DNA

Description: Cytosine modification at carbon 5 (C5) represents an important epigenetic modification. It is also believed to be the starting substrate for the Ten-Eleven Translocation (TET) family of enzymes and their associated oxidation pathways. 5-methyl-dCTP offers the ability to enzymatically make defined fully methylated cytosine-substituted DNA, which can be used for a variety of biochemical and cellular applications.

5-methyl-dCTP (2'-deoxy-5-methylcytidine 5'-triphosphate) is supplied as a 10 mM solution at pH 7. Nucleotide concentration is determined by measurements of absorbance at 260 nm.

Formula: C₁₀H₁₅N₃O₁₃P₃ (free acid)

Concentration: 10 mM solution

Molecular weight: 481.1 (acid form)

Diluent Compatibility: Can be diluted using sterile distilled water, preferably Milli-Q® water, or can be diluted using sterile TE [10 mM Tris-HCl, 1 mM EDTA (pH 7.5)]

MILLI-Q® is a registered trademark of Millipore, Inc.

As Supply Chain Manager, Dan is responsible for making sure that NEB has uninterrupted access to the critical components needed to manufacture all its products. His attention to detail ensures that our products are ready to ship when your order is placed. In his spare time, Dan spends as much time as possible mountain biking.




NEBNext® Reagents for ChIP-Seq Library Preparation

Epi

NEBNext reagents are a series of highly pure reagents that facilitate library preparation of DNA or RNA for downstream applications, such as next generation sequencing and expression library construction. These reagents undergo stringent quality controls and functional validation, ensuring maximum yield, convenience and value.

For sample preparation of a ChIP-Seq DNA library, NEB offers kits, oligos and modules that support fast workflows. To decide which products to choose, use the selection chart below. For more information on NEBNext reagents for library preparation, see pages 146–175.



Katherine is a Development Associate for our Organic Synthesis Team. In this role, she works to develop, process and purify oligos for our products that help support ongoing research in genetic sequencing and screening. When not at work, Katherine spends time enjoying local trails and beaches in the summer and skiing in the winter.



NEBNext Enzymatic Methyl-seq

NEBNext Enzymatic Methyl-seq Kit

#E7120S	24 reactions\$895
#E7120L	96 reactions\$3,360

NEBNext Enzymatic Methyl-seq Conversion Module

#E7125S	24 reactions\$185
#E7125L	96 reactions\$675

NEBNext Q5U Master Mix

#M0597S	50 reactions\$125
#M0597L	250 reactions\$495

NEBNext Multiplex Oligos for Enzymatic Methyl-seq (Unique Dual Index Primer Pairs)


#E7140S	24 reactions\$135
#E7140L	96 reactions\$535

- Superior sensitivity of detection of 5-mC and 5-hmC
- Larger library insert sizes
- More uniform GC coverage
- Greater mapping efficiency
- High-efficiency library preparation


NEBNext Enzymatic Methyl-seq is an enzymatic alternative to bisulfite conversion with superior performance. For more information, including extensive performance data, visit NEBNext.com.

Description: While bisulfite sequencing has been the gold standard for the study of DNA methylation, this conversion treatment is damaging to DNA, resulting in DNA fragmentation, loss and GC bias. The NEBNext Enzymatic Methyl-seq Kit (EM-seq™) provides an enzymatic alternative to whole genome bisulfite sequencing (WGBS), combined with high-efficiency streamlined library preparation suitable for Illumina sequencing.

The highly effective EM-seq enzymatic conversion minimizes damage to DNA and, in combination with the supplied NEBNext Ultra II library preparation workflow reagents, results in high quality libraries that enable superior detection of 5-mC and 5-hmC from fewer sequencing reads.



NEBNext Enzymatic Methyl-seq libraries have larger insert sizes 50 ng Human NA12878 genomic DNA was sheared to 300 bp using the Covaris® S2 instrument and used as input into EM-seq and WGBS protocols. For WGBS, NEBNext Ultra II DNA was used for library construction, followed by the Zymo Research EZ DNA Methylation-Gold™ kit for bisulfite conversion. Libraries were sequenced on an Illumina MiSeq (2 x 76 bases) and insert sizes were determined using Picard 2.18.14. The normalized frequency of each insert size was plotted, illustrating that library insert sizes are larger for EM-seq than for WGBS, and indicating that EM-seq does not damage DNA as bisulfite treatment does in WGBS.



EM-seq identifies more CpGs than WGBS, at lower sequencing coverage depth with superior uniformity of GC coverage. 10, 50 and 200 ng Human NA12878 genomic DNA was sheared to 300 bp using the Covaris S2 instrument and used as input into EM-seq and WGBS protocols. For WGBS, NEBNext Ultra II DNA was used for library construction, followed by the Zymo Research EZ DNA Methylation-Gold Kit for bisulfite conversion. Libraries were sequenced on an Illumina NovaSeq® 6000 (2 x 100 bases). Reads were aligned to hg38 using bwa-meth 0.2.2.

A: Coverage of CpGs with EM-seq and WGBS libraries was analyzed using 324 million paired end reads, and each top and bottom strand CpGs were counted independently, yielding a maximum of 56 million possible CpG sites. EM-seq identifies more CpGs at lower depth of sequencing.

B: GC coverage was analyzed using Picard 2.17.2 and the distribution of normalized coverage across different GC contents of the genome (0-100%) was plotted. EM-seq libraries have significantly more uniform GC coverage, and lack the AT over-representation and GC under-representation typical of WGBS libraries.





Giving Back

In these times of uncertainty, it is vital to think of others and to help when we can. Philanthropy, through foundation and corporate giving, offers a way to marshal resources to help address the broad array of challenges facing our planet. Don Comb's love of the natural world and his philanthropic impulse compelled him to establish or support many organizations over the years. NEB continues these efforts, supporting organizations throughout the world that help advance science, protect the environment, promote education, and help the less fortunate. NEB creates many opportunities for employees to give back, which, in turn, has created a generous culture throughout.

In 1982, inspired by his deep love of nature and fascination with the arts and diverse cultures, Don created the New England Biolabs Foundation (NEBF) as a private, charitable foundation to address problems including loss of biodiversity and public health challenges in many regions of the world. Today, NEBF's mission is to foster community-based conservation of landscapes and seascapes and the biocultural diversity found in these places, working in countries of Central America, Andean South America and West Africa, and also locally on the North Shore of Massachusetts. NEBF provides small grants to grassroots organizations in these regions for projects such as coral reef conservation, watershed restoration, supporting food sovereignty through agroecology, and helping communities sustain their traditional ecological knowledge and linguistic diversity.

Don believed in a simple but powerful premise that was well ahead of its time: "if you put small amounts of money in the hands of the people closest to the biodiversity, you can make a big impact." His commitment to invest in small and grassroots organizations working to protect their landscapes, livelihoods and biodiversity shaped the foundation's approach early on and remains its touchstone to this day.

In the early days of the foundation, Don felt it was essential to get out in the field and visit the communities and projects that NEBF supported in many different parts of the world. He was intrigued by Papua New Guinea (PNG) and listened to the local people when they expressed concern for conserving their environment and maintaining their traditions. In response, NEBF funded various innovative projects in PNG, including a forest-friendly means for communities to build housing and canoes, and the establishment of the Wau Biological Station, which aimed to protect local plants and trees and establish a butterfly range. Don's commitment to visiting projects in the field laid the groundwork for NEBF's approach to the present day in which the foundation's staff and board make an effort to visit projects and spend time in the field with NEBF's partners.


Building on these early roots, today NEBF's grant-making extends to places as diverse as Peru, Ghana and Belize. NEBF continues to fund in areas that are central to the crucible of threats facing our planet and humanity today: the intersection of climate change, biodiversity loss, extinction of cultures and languages, social injustice and food insecurity. The foundation complements its grant-making with capacity-building opportunities, such as workshops, enabling grantees to develop new skills and strategies and build relationships with each other – as a result, a growing network of NEBF grantees and partners work closely together to tackle these problems.

Don Comb crosses a shaky vine bridge
in the Southern Highlands of Papua
New Guinea, 1985.
Credit: Dave Comb

Learn more
about the New
England Biolabs
Foundation.



Cellular Analysis



Novel tools to study expression & function of proteins.

Cell imaging analysis can use fluorescent dyes, fluorophore-labeled molecules or recombinant protein plasmid systems. Recombinant protein labeling systems and bioluminescent reporter systems are among the most sensitive fluorescence methods for imaging expression, transport, co-localization and degradation in either fixed or living cells. Protein labeling systems offer many advantages. For example, color changes can be easily implemented by using different substrates. Protein labeling systems can involve the use of tag-specific antibodies or antibodies to separate epitopes engineered into a plasmid tag system for detection. Protein labeling systems can be used with non-cell permeable substrates to enable the specific imaging of cell surface targets. This strategy is not possible with bioluminescent recombinant systems. In living cells, protein labeling substrates can be introduced and followed in cells over time. Two separate cellular targets can also be imaged simultaneously, using protein labeling systems with mutually exclusive, tag-specific fluorescent substrates.

Studies of protein expression, interactions and structure, often use reporter systems to introduce and select for gene targets in cells. Reporter genes confer drug resistance, bioluminescence or fluorescence properties in the cells into which they are introduced. Typical reporter studies link reporter genes directly to a promoter region of interest, the function of which can be monitored by the reporter activity. Protein fusion tagging is used to detect subcellular localization, degradation, protein-protein interactions, etc. Typical fusion tags are fluorescent proteins (e.g., eGFP) or small protein epitopes (e.g., FLAG, Myc HA) which can be detected by fluorescence FACS or western blots. New generations of reporter gene systems expand the range of applications and enhance experimental possibilities.

Featured Products

296 SNAP-Cell 647-SiR

296 SNAP-Surface® 649

Featured Tools & Resources

354 Troubleshooting Guide for SNAP-tag Technology



Application Notes



Videos and Tutorials:
SNAP-tag Technology




Cellular Imaging & Analysis

Overview of SNAP- and CLIP-tag Technologies	294
Comparison of SNAP-tag/CLIP-tag Technology to GFP	295
Fluorescent Substrates	296
Blocking Agents	297
Anti-SNAP-tag Antibody, Polyclonal	297
SNAP-tag Purified Protein	297
Cloning Vectors	298
Biotin Labels	298
SNAP-Capture	299
Building Blocks	299



Cellular Imaging & Analysis



Live HeLa cell transfected with pSNAP₇-tubulin and pCLIP₇-Cox8A (mitochondrial cytochrome oxidase 8A). Cells were labeled with 3 μ M SNAP-Cell TMR-Star (red) and 5 μ M CLIP-Cell 505 (green) for 25 minutes and counterstained with Hoechst 33342 (blue) for nuclei.

Features of SNAP-tag and CLIP-tag:


- Clone and express once, then use with a variety of substrates
- Non-toxic to living cells
- Wide selection of fluorescent substrates
- Highly specific covalent labeling
- Simultaneous dual labeling

Applications of SNAP-tag and CLIP-tag:

- Simultaneous dual protein labeling inside or on the surface of live cells
- Protein localization and translocation
- Pulse-chase experiments
- Receptor internalization studies
- Selective cell surface labeling
- Protein pull down assays
- Protein detection in SDS-PAGE
- Flow cytometry
- High throughput binding assays in microtiter plates
- Biosensor interaction experiments
- FRET-based binding assays
- Single molecule labeling
- Super-resolution microscopy

New England Biolabs offers an innovative technology for studying the function and localization of proteins in live and fixed cells. Covalent protein labeling brings simplicity and versatility to the imaging of mammalian proteins in live cells, as well as the ability to capture proteins *in vitro*. The creation of a single genetic

construct generates a fusion protein which, when covalently attached to a variety of fluorophores, biotin or beads, provides a powerful tool for studying proteins. For added flexibility, NEB offers two systems in which the protein is labeled by a self-labeling fusion protein (SNAP-tag[®] and CLIP-tag[™]).



Protein labeling with SNAP-tag (gold) and CLIP-tag (purple). The SNAP- or CLIP-tag is fused to the protein of interest (blue). Labeling occurs through covalent attachment to the tag, releasing either a guanine or a cytosine moiety.

SNAP-tag and CLIP-tag – Self-Labeling Tag Technology

The SNAP- and CLIP-tag protein labeling systems enable the specific, covalent attachment of virtually any molecule to a protein of interest. There are two steps to using this system: cloning and expression of the protein of interest as a SNAP-tag fusion, and labeling of the fusion with the SNAP-tag substrate of choice. The SNAP-tag is based on the human O⁶-alkylguanine-DNA-alkyltransferase (hAGT), a DNA repair protein. SNAP-tag substrates are fluorophores, biotin or beads conjugated to guanine or chloropyrimidine leaving groups via a benzyl linker. In the labeling reaction, the substituted benzyl group of the substrate is covalently attached to the SNAP-tag. CLIP-tag is a modified version of SNAP-tag, engineered to react with benzylcytosine rather than benzylguanine derivatives. When used in conjunction with SNAP-tag, CLIP-tag enables the orthogonal and complementary labeling of two proteins simultaneously in the same cells.

SNAP-Cell[®]: SNAP-Cell labels are cell-permeant and uniquely suited for the labeling of SNAP-tag fusion proteins inside living or fixed cells, on cell surfaces or *in vitro*. These labels are spread across the visible spectrum, ranging from blue to red. Non-fluorescent cell-permeable blocking agent is also available.

SNAP-Surface[®]: SNAP-Surface labels are non-cell-permeant and routinely used to label SNAP-tag fusion proteins on the surface of living cells, in fixed cells or *in vitro*. These labels are spread across the visible spectrum and include the photostable AlexaFluor[®] dyes and a variety of other commonly used fluorophores. Non-fluorescent, non-cell-permeable blocking agent is also available.

CLIP-Cell[™]: CLIP-Cell labels are cell-permeant and uniquely suited for the labeling of CLIP-tag fusion proteins inside living or fixed cells, on cell surfaces or *in vitro*. The CLIP-tag is a derivative of the SNAP-tag that reacts with orthogonal substrates, allowing simultaneous labeling of two expressed proteins with different fluorophores. Non-fluorescent cell-permeable blocking agent is also available.

CLIP-Surface[™]: CLIP-Surface labels are non-cell-permeant and routinely used to label CLIP-tag fusion proteins on the surface of living cells, in fixed cells or *in vitro*. The CLIP-tag is a derivative of the SNAP-tag that reacts with orthogonal substrates, allowing simultaneous labeling of two expressed proteins with different fluorophores. The labels include fluorophores at commonly used areas of the visible spectrum, such as 488, 547 and 647 nm.



Comparison of SNAP-tag/CLIP-tag Technologies to GFP

While SNAP/CLIP-tag technologies are complementary to GFP, there are several applications for which SNAP- and CLIP-tag self-labeling technologies are advantageous.


APPLICATION	SNAP-tag/CLIP-tag	GFP AND OTHER FLUORESCENT PROTEINS
Time-resolved fluorescence	Fluorescence can be initiated upon addition of label	Color is genetically encoded and always expressed. Photoactivatable fluorescent proteins require high intensity laser light, which may activate undesired cellular pathways (e.g., apoptosis)
Pulse-chase analysis	Labeling of newly synthesized proteins can be turned off using available blocking reagents (e.g., SNAP-Cell® Block)	Fluorescence of newly synthesized proteins cannot be specifically quenched to investigate dynamic processes
Ability to change colors	A single construct can be used with different fluorophore substrates to label with multiple colors	Requires separate cloning and expression for each color
Surface specific labeling	Can specifically label subpopulation of target protein expressed on cell surface using non-cell permeant substrates	Surface subpopulation cannot be specifically visualized
Single molecule detection	Conjugation with high quantum yield and photostable fluorophores	Fluorescent proteins are generally less bright and photobleach quicker than most organic fluorophores
Visualizing fixed cells	Resistant to fixation; strong labeling	Labile to fixation; weak labeling
Pull-down studies	"Bait" proteins can be covalently captured on BG beads	Requires anti-GFP antibody to non-covalently capture "bait" protein, complicating downstream analysis
Live animal imaging	Cell permeable far-red dye available, permitting deep tissue visualization	Signal is easily quenched by fixation (whole-mount specimens or thick sections); limited spectral flexibility and weaker fluorescence



Sam Brown works as a Production Technician at our Rowley Production Facility. In this role, he helps with fermentation of several of the GMP-grade* enzymes. In his spare time, Sam enjoys skiing, sailing and reading about history.

*Learn more about GMP-grade enzymes on page 7.

Fluorescent Substrates for Protein Labeling



Live HeLa cell transfected with pSNAP_{ER} (endoplasmic reticulum) and pCLIP₇-tubulin. Cells were labeled with 3 μ M SNAP-Cell TMR-Star (red) and 5 μ M CLIP-Cell 505 (green) for 25 minutes and counterstained with Hoechst 33342 (blue) for nuclei.

- Fluorescently label SNAP-tag or CLIP-tag fusions for cellular imaging
- Labels span fluorescent imaging spectrum from aqua (430 nm) to far-red (647+ nm) wavelengths
- Cell-permeable and non-cell-permeable labels available

NEB offers a large selection of fluorescent labels (substrates) for SNAP-tag and CLIP-tag fusion proteins. SNAP-tag substrates consist of a fluorophore conjugated to guanine or chloropyrimidine leaving groups via a benzyl linker, while CLIP-tag substrates consist of a fluorophore conjugated to a cytosine leaving group via a benzyl linker. These substrates

will label their respective tags without the need for additional enzymes. Cell-permeant substrates (SNAP- and CLIP-Cell) are suitable for both intracellular and cell-surface labeling, whereas non-cell-permeant substrates (SNAP- and CLIP-Surface) are specific for fusion proteins expressed on the cell surface only.

SELF-LABELING TAG						
	APPLICATIONS	NEB #	EXCITATION*	EMISSION ⁽¹⁾	SIZE	PRICE
SNAP-tag	Cell-Permeable					
	SNAP-Cell 430	S9109S	421	444,484	50 nmol	\$353
	SNAP-Cell 505-Star	S9103S	504	532	50 nmol	\$353
	SNAP-Cell Oregon Green®	S9104S	490	514	50 nmol	\$392
	SNAP-Cell TMR-Star	S9105S	554	580	30 nmol	\$353
	SNAP-Cell 647-SiR	S9102S	645	661	30 nmol	\$353
	Non-cell-permeable					
	SNAP-Surface Alexa Fluor® 488	S9129S	496	520	50 nmol	\$392
	SNAP-Surface 488	S9124S	506	526	50 nmol	\$353
	SNAP-Surface Alexa Fluor 546	S9132S	558	574	50 nmol	\$392
	SNAP-Surface 549	S9112S	560	575	50 nmol	\$353
	SNAP-Surface 594	S9134S	606	626	50 nmol	\$353
	SNAP-Surface Alexa Fluor 647	S9136S	652	670	50 nmol	\$392
	SNAP-Surface 649	S9159S	655	676	50 nmol	\$353
CLIP-tag	Cell-Permeable					
	CLIP-Cell 505	S9217S	504	532	50 nmol	\$353
	CLIP-Cell TMR-Star	S9219S	554	580	30 nmol	\$353
	Non-cell-permeable					
	CLIP-Surface 488	S9232S	506	526	50 nmol	\$353
	CLIP-Surface 547	S9233S	554	568	50 nmol	\$353
	CLIP-Surface 647	S9234S	660	673	50 nmol	\$353

* Excitation and emission values determined experimentally for labeled protein tag.

⁽¹⁾ Colors are based on the electromagnetic spectrum. Actual color visualization may vary.

This table lists all currently available fluorescent substrates for SNAP-tag and CLIP-tag, along with excitation and emission wavelengths (determined from a labeled fusion tag, rather than the unreacted substrate).

Blocking Agents

- *Irreversible blocking*
- *Ideal for pulse-chase applications*

Blocking agents are non-fluorescent substrates that block the reactivity of the SNAP- or CLIP-tag intracellularly (SNAP-Cell Block and CLIP-Cell Block) or on the surface of live cells (SNAP-Surface Block and CLIP-Cell Block). They can be used to generate inactive controls in live cell and *in vitro* labeling experiments performed with SNAP- or CLIP-tag fusion proteins.

SNAP- and CLIP-Cell Block are highly membrane permeant and once inside the cell react with the SNAP- or CLIP-tag, irreversibly inactivating them for subsequent labeling steps.

SNAP-Surface Block also reacts with the SNAP-tag irreversibly, inactivating it for subsequent labeling steps. This blocker is largely membrane impermeant essentially limiting blocking to cell surface-exposed SNAP-tags.

PRODUCT	NEB #	APPLICATION	SIZE	PRICE
SNAP-Cell Block	S9106S	Block SNAP-tag inside live cells and <i>in vitro</i>	100 nmol	\$146
CLIP-Cell Block	S9220S	Block CLIP-tag inside or on the surface of live cells and <i>in vitro</i>	100 nmol	\$146
SNAP-Surface Block	S9143S	Block SNAP-tag on the surface of live cells and <i>in vitro</i>	200 nmol	\$148

Anti-SNAP-tag Antibody, Polyclonal

#P9310S 100 µl\$280

Description: The Anti-SNAP-tag Antibody (Polyclonal) can be used in Western blots with SNAP-tag and CLIP-tag proteins. Polyclonal antibodies are produced from the immunization of rabbit with purified recombinant SNAP-tag protein and affinity purified using SNAP-BG resin.

Sensitivity: 5 ng of SNAP-tag per load in Western blotting.

Recommended Dilution: 1:1000

SNAP-tag Purified Protein

#P9312S 50 µg\$102

Description: SNAP-tag Purified Protein can be used as a positive control for *in vitro* labeling with various SNAP-tag fluorescent substrates. The coding sequence of SNAP-tag was cloned into a pTXB1 derived *E. coli* T7 expression vector. SNAP-tag protein was expressed and purified according to the instructions in the IMPACT™ kit manual (NEB #E6901). The purified SNAP-tag protein was dialyzed into


1X phosphate buffered saline (PBS) solution containing 1 mM DTT at 1 mg/ml (50 µM) and stored at -80°C.

Molecular Weight: 19,694

Lovely works at our subsidiary office in Australia. As a member of Customer Service, she oversees shipments to New Zealand. Lovely also enjoys baking and is a budding organist.



Cloning Vectors



Live HeLa cell transfected with pSNAP_i-tubulin and pCLIP_i-H2B constructs generated using pSNAP_i and pCLIP_i vectors. Cells were labeled with 3 μ M SNAP-Cell TMR-Star (red) and 5 μ M CLIP-Cell 505 (green) for 25 minutes and counterstained with Hoechst 33342 (blue) for nuclei.

- Vectors for mammalian and bacterial expression available

Vectors are available for SNAP-tag and CLIP-tag fusion protein expression and labeling in mammalian and bacterial systems. The mammalian SNAP_i and CLIP_i vectors express faster-reacting variants of the SNAP- and CLIP-tags than previously available vectors. Improved polylinker sequences both upstream and downstream from the tag allow expression of the tag on either end of the protein of interest, under control of the CMV promoter. SNAP-tag and CLIP-tag expression vectors contain a neomycin resistance (NeoR) gene for selection of stable transfectants, together with an IRES element for efficient expression of both the fusion protein and NeoR. Codon usage has been optimized for mammalian expression. Control plasmids encoding fusion proteins that are localized to the nucleus (H2B), mitochondria (Cox8A) and cell surface (ADR β 2, NK1R) are also available through Addgene.

The bacterial expression vector pSNAP-tag(T7)-2 includes cloning sites both upstream and downstream from the SNAP-tag, which is under control of the T7 promoter. Codon usage in the SNAP-tag gene has been optimized for *E. coli* expression.

Source: Isolated from an *E. coli* strain by a standard plasmid purification procedure. Plasmids have been purified free of endotoxins for efficient transfection

Concentration: 500 μ g/ml

Restriction Map: For a more detailed description and restriction map of pSNAP_i Vector, see page 362. See www.neb.com for sequence and map files for all expression and control plasmids.


PRODUCT	NEB #	FEATURES	SIZE	PRICE
pSNAP _i Vector	N9183S	stable and transient mammalian expression	20 μ g	\$172
pSNAP-tag(T7)-2 Vector	N9181S	bacterial expression under T7 control	20 μ g	\$172
pCLIP _i Vector	N9215S	stable and transient mammalian expression	20 μ g	\$172

Biotin Labels

SNAP-Biotin®		
#S9110S	50 nmol\$281
CLIP-Biotin		
#S9221S	50 nmol\$281

- Label SNAP-tag and CLIP-tag fusions with biotin
- Compatible with a variety of streptavidin conjugates
- Attach to streptavidin surfaces on microtiter plates

For optimal flexibility with existing technologies, biotinylated labels are available for studies using streptavidin platforms. Cell-permeant (SNAP-Biotin and CLIP-Biotin) substrates are based on biotin with an amidocaproyl linker. Biotin labels are suitable for applications such as biotinylation of fusion proteins for detection with streptavidin fluorophore conjugates or labeling in solution for analysis by SDS-PAGE/ Western blot. Biotin labels are also used for capture with streptavidin surfaces for binding and interaction studies.




Western blot analysis of biotin labeling reactions using anti-Biotin Antibody (CST #7075). SNAP-tag and CLIP-tag (5 μ M) labeled with a biotin-containing substrate (10 μ M). Marker M is Biotinylated protein ladder (CST #7727).

SNAP-Capture

- *Selectively capture SNAP-tag fusion proteins from solution*
- *Ideal for protein pull-down experiments or proteomic analysis*

SNAP-Capture products are agarose or magnetic agarose beads coupled to a benzylguanine substrate, used to selectively capture and immobilize SNAP-tag fusion proteins from solution. These beads have a

high loading capacity for SNAP-tag fusion proteins and show very low non-specific adsorption of proteins from a complex lysate, making them especially suitable for pull-down applications.



Substrate structure on
SNAP-Capture Pull-Down Resin

PRODUCT	NEB #	SIZE	PRICE
SNAP-Capture Pull-Down Resin	S9144S	2 ml	\$235
SNAP-Capture Magnetic Beads	S9145S	2 ml	\$235

Building Blocks

- *Synthesize new SNAP-tag and CLIP-tag substrates*
- *Make surfaces for protein immobilization*
- *Attach novel molecules or ligands to proteins*
- *Create custom substrates for protein labeling*


For advanced users with novel probes interested in working with SNAP-tag and CLIP-tag labeling technologies, building blocks are available for linkage of the core benzylguanine (BG) and benzylcytosine (BC) moieties to activated esters, primary amines and thiol groups. A variety of functional groups allows the choice of chemical coupling approaches to suit the molecule or surface to be coupled. Couple onto surfaces such as Biacore® chips or microarrays for specific protein immobilization. Couple onto peptides, proteins and DNA oligomers.

Couple onto new fluorophores or affinity reagents for specific protein labeling. Labeling is gentle, precise, and versatile: one label is covalently bound under biological conditions in a defined position.

References:


References for enzyme properties and applications for this product can be found at www.neb.com.

PRODUCT	NEB #	STRUCTURE	APPLICATION	SIZE	PRICE
BG-NH ₂	S9148S		SNAP-tag substrate. Suitable for linkage to NHS esters and other activated carboxylic esters.	2 mg	\$416
BG-PEG-NH ₂	S9150S		SNAP-tag substrate. PEG-linker gives superior flexibility. Particularly suited for immobilization on solid surfaces.	2 mg	\$416
BG-GLA-NHS	S9151S		SNAP-tag substrate. Activated as NHS ester. Reacts with primary amines.	2 mg	\$420
BG-Maleimide	S9153S		SNAP-tag substrate. Activated as maleimide. Reacts with thiols.	2 mg	\$416
BC-NH ₂	S9236S		CLIP-tag substrate. Suitable for linkage to NHS esters and other activated carboxylic esters.	2 mg	\$416



RESTRICTION ENDONUCLEASES
NUCLEIC ACIDS
tRNA's
DNA LIGASE
BACTERIAL CELLS

1975 - 1976
CATALOG





And the NEB Catalog

When one thinks of a catalog, they might think of a list of products that are available for sale. As a user, one might browse a catalog to find a particular product, or just to see what is new from their favorite brand. But to NEB, the catalog is so much more — it is a reflection of who we are, the products we offer, and the resources we can provide to help advance science.

In 1975, the first NEB catalog was released, featuring 7 restriction enzymes, 7 DNA or RNA modifying enzymes and 6 transfer RNAs — meager in content but with promise for more. Over the next few years, NEB contemplated how to add value to this list of products, and the answer was simple — share the science. In the early 1980s, NEB began curating technical information related to restriction enzymes and released the first Technical Reference section. This included information on how enzymes behaved and how to set up reactions properly — information that is now easily accessible through the internet, but back then was found mostly in books. Over time, this content expanded to include other cloning-related applications. This idea of a “Benchtop Reference Tool” was unique at the time, and as a result, the NEB Catalog and Technical Reference became a mainstay on the lab benches of undergraduate and graduate researchers around the world. It was the beginning of a collaborative relationship between scientists and NEB.

The catalog also provides an opportunity to share ideas that Don Comb and NEB felt passionate about. Each edition focuses on a theme and contains a collection of minireviews that addresses various scientific, environmental and/or humanitarian challenges facing the world. This represents a unique opportunity to reflect on ideas that are important to NEB, and to think about how we can do better as scientists, and in life.

Over the years, the catalog has expanded to thousands of products, supporting a wide variety of applications. There are even products that may not have an application yet — Don always loved to make unique enzymes available to scientists, with hopes that they would find new applications for them. Sir Richard Roberts, Nobel Laureate and NEB’s Chief Scientific Officer, took this a step further and spearheaded the “Enzymes for Innovation” initiative — sharing enzymes with unique activities that have the potential to lead to new discoveries.

Is it possible to build a relationship through a catalog? And are catalogs even necessary in this digital age? We think so! A catalog can be much more than a transactional exercise — it can be part of a sustainable relationship based on the sharing of knowledge. Sure, your device can give you instant answers, but there is something comforting about flipping through a book, complete with earmarks and notes in the sidebar about how you modified an experiment, or flagged products to look at later. We have often heard scientists recall their first NEB Catalog, or which cover was their favorite — in fact, some have even saved them!

In our minds, the NEB Catalog is more than a book — it embodies the ethos of NEB. The passion that our scientists put into developing our products, the humility in the way the products and data are presented, and the genuine way the stories reflect who we are as a company and what we feel is important — empowering science, preserving the environment, and collaborating with our peers.

Can you guess
what year these
catalog cover images
came from?



A collection of past
NEB Catalog covers.

Reference Appendix



Technical Support – for scientists, by scientists

As a partner to the scientific community, New England Biolabs is committed to providing top quality tools and scientific expertise. This philosophy still stands, and has led to long-standing relationships with many of our fellow scientists. NEB's commitment to scientists is the same regardless of whether or not they purchase product from NEB: their ongoing research is supported by our catalog, website and technical staff.

NEB's technical support model is unique as it utilizes most of the scientists at NEB. Several of our product lines have designated technical support scientists assigned to servicing customers in those application areas. Any questions regarding a product could be dealt with by one of the technical support scientists, the product manager who manufactures it, the product development scientist who optimizes it, or a researcher who uses the product in their daily research. As such, customers are supported by scientists and often experts in the product or its application.

To access technical support:

- Call 1-800-632-7799 (Monday – Friday: 9:00 am - 6:00 pm EST)
- Submit an online form at www.neb.com/techsupport
- Email info@neb.com
- International customers can contact a local NEB subsidiary or distributor. For more information see inside back cover.



Learn more
about NEB's tech
support program.

Featured Tools & Resources

306 Optimizing Restriction
Enzyme Reactions

309 Performance Chart for
Restriction Enzymes

336 Troubleshooting Guide
for Cloning


322 Guidelines for
PCR Optimization

328 Cleavage Close to the
End of DNA Fragments

304 Online Interactive Tools



Visit the Tools & Resources
tab at www.neb.com to find
additional online tools, video
tech tips and tutorials to help
you in your research.



Online Interactive Tools	304	Nucleic Acid Purification	
Databases, Mobile Apps & Freezer Program Locator	305	Troubleshooting Guide for DNA Cleanup & Plasmid Purification Using Monarch Kits	341
Restriction Endonucleases		Choosing Sample Input Amounts When Using the Monarch Genomic DNA Purification Kit	342
Optimizing Restriction Enzyme Reactions	306	Troubleshooting Guide for Genomic DNA Purification Using Monarch Kit	343
Double Digestion	307	Choosing Sample Input Amounts When Using the Monarch Total RNA Miniprep Kit	344
Types of Restriction Enzymes	307	Troubleshooting Guide for Total RNA Extraction & Purification using Monarch Kits	345
Troubleshooting Guide	308	Guidelines for Choosing Sample Input Amounts When Using the Monarch HMW DNA Kit for Cells & Blood	346
Performance		Guidelines for Choosing Sample Input Amounts When Using the Monarch HMW DNA Extraction Kit for Tissue, Bacteria and Other Samples	347
Performance Chart for Restriction Enzymes	309–314	Transformation	
Activity of Enzymes at 37°C	315	Genetic Markers	348–349
Activity of DNA Modifying Enzymes in CutSmart Buffer	315	Enhancing Transformation Efficiency	350
Tips for Avoiding Star Activity	316	Protein Expression with T7 Express Strains	351
High-Fidelity (HF) Restriction Enzymes	316	NGS Library Preparation	
Properties of Restriction Enzymes		Troubleshooting for Guide for NEBNext DNA Library Prep Kits	352
Isoschizomers	317	Guidelines for NGS Library Prep	353
Survival in a Reaction	317–318	Protein Expression	
Cleavage of Supercoiled DNA	318	Labeling with SNAP-tag Technology	354
Generating New Cleavage Sites		Troubleshooting Guide	355–365
Compatible Cohesive Ends	319–320	Odds and Ends	
Methylation Sensitivity		General Molecular Biology Information	366–369
Dam, Dcm and CpG Methylation	321	Index	370
PCR, qPCR & RT-qPCR			
General Guidelines for PCR Optimization	322		
PCR Troubleshooting Guide	323		
Optimization Tips for Luna qPCR	324		
Optimization Tips for Luna One-Step RT-qPCR	325		
Luna qPCR Troubleshooting Guide	326		
Luna One-Step RT-qPCR Troubleshooting Guide	327		
Cleavage Close to the End of DNA Fragments	328		
Activity of Restriction Enzymes in PCR Buffers	329–330		
Cloning & DNA Assembly			
Getting Started with Molecular Cloning	331		
Optimization Tips for Golden Gate Assembly	332		
Optimization Tips for NEBuilder HiFi DNA Assembly	333		
Traditional Cloning Quick Guide	334–335		
Troubleshooting Guide for Cloning	336–339		
Optimization Tips for Your Cloning Reactions	339–340		

Online Interactive Tools

Use the Tools & Resources tab at www.neb.com to access our growing selection of interactive technical tools. These tools can also be accessed directly by visiting www.neb.com/nebtools.

NEB scientists are often involved in the development of online tools that will aid in their research. We are now making these tools and in some circumstances, the source code, available for you to evaluate. To learn more, visit www.neb.com/NEBetaTools.

Online Tools

Competitor Cross-Reference Tool



Use this tool to select another company's product and find out which NEB product is compatible. Choose either the product name or catalog number from the available selections, and this tool will identify the recommended NEB product.

DNA Sequences and Maps Tool



With the DNA Sequences and Maps Tool, find the nucleotide sequence files for commonly used molecular biology tools, including plasmid, viral and bacteriophage vectors.

Double Digest Finder



Use this tool to guide your reaction buffer selection when setting up double digests, a common timesaving procedure. Choosing the right buffers will help you to avoid star activity and loss of product.

EnGen sgRNA Template Oligo Designer



EnGen sgRNA Template Oligo Designer can be used to design target-specific DNA oligos for use with the EnGen sgRNA Synthesis Kit, *S. pyogenes* (NEB #E3322).

Enzyme Finder



Use this tool to select restriction enzymes by name, sequence, overhang or type. Enter your sequence using single letter code nomenclature, and Enzyme Finder will identify the right enzyme for the job.

Exo Selector



Use this tool to simplify the process of selecting the appropriate exonucleases for use in your nucleic acid digestion workflows. The tool guides you to product recommendations based on your answers to a few simple questions.

Glycan Analyzer



Use this tool to interpret ultra or high pressure liquid chromatography (UPLC/HPLC) N-glycan profiles following exoglycosidase digestions.

Ligase Fidelity Tools



This set of beta tools can be used to design your Golden Gate Assemblies – visualize overhang ligation preferences with the Ligase Fidelity Viewer™, Predict high-fidelity junction sets with GetSet™, and Split DNA sequences for scarless high-fidelity assembly with SplitSet™.

NEB Golden Gate Assembly Tool



Use this tool to assist with in silico DNA construct design for Golden Gate DNA assembly. It enables the accurate design of primers with appropriate Type IIS restriction sites and overlaps, quick import of sequences in many formats and export of the final assembly, primers and settings.

NEB LAMP Primer



NEB LAMP Primer Design Tool can be used to design primers for your Loop-mediated Isothermal Amplification. Fixed primers can be specified for the design of LAMP primers, and subsequent Loop primers are then designed based on LAMP primer selection.

NEBaseChanger®



NEBaseChanger can be used to design primers specific to the mutagenesis experiment you are performing using the Q5® Site-Directed Mutagenesis Kit. This tool will also calculate a recommended custom annealing temperature based on the sequence of the primers by taking into account any mismatches.

NEBcloner®



Use this tool to find the right products and protocols for each step (digestion, end modification, ligation, transformation and mutagenesis) of your next traditional cloning experiment. Also, find other relevant tools and resources to enable protocol optimization.

NEBcutter® V2.1



Identify the restriction sites within your DNA sequence using NEBcutter. Choose between Type II and commercially available Type III restriction enzymes to digest your DNA. NEBcutter will indicate cut frequency and methylation-state sensitivity.

NEBioCalculator®



Use this tool for your scientific calculations and conversions for DNA and RNA. Options include conversion of mass to moles, ligation amounts, conversion of OD to concentration, dilution and molarity. Additional features include sgRNA template oligo design and qPCR library quantification.

NEBNext Custom RNA Depletion Design Tool



This tool designs probes to be used with the NEBNext RNA Depletion Core Reagent Set or to supplement an existing NEBNext Depletion Kit for the depletion of unwanted RNA species.

NEBNext Selector



Use this tool to guide you through the selection of NEBNext reagents for next generation sequencing sample preparation.

Online Tools (continued)

NEBuilder® Assembly Tool

Use this tool to design primers for your DNA assembly reaction, based on the entered fragment sequences and the polymerase being used for amplification.

Read Coverage Calculator

This tool allows for easy calculation of values associated with read coverage in NGS protocols.

PCR Fidelity Estimator

Estimate the percentage of correct DNA copies (those without base substitution errors) per cycle of PCR for selected DNA polymerases.

Thermostable Ligase Reaction Temperature Calculator

This tool will help you estimate an optimal reaction temperature to minimize mismatch for thermostable ligation of two adjacent ssDNA probes annealed to a template.

PCR Selector

Use this tool to help select the right DNA polymerase for your PCR setup. Whether your amplicon is long, complex, GC-rich or present in a single copy, the PCR selection tool will identify the perfect DNA polymerase for your reaction.

Tm Calculator

Determine the optimal annealing temperature for your amplicon with our Tm Calculator. Simply input your DNA polymerase, primer concentration and your primer sequence, and the Tm Calculator will guide you to successful reaction conditions.

Additional Databases


Polbase®

Polbase is a repository of biochemical, genetic and structural information about DNA Polymerases.

REBASE®

Use this tool as a guide to the ever-changing landscape of restriction enzymes. REBASE, the Restriction Enzyme DataBASE, is a dynamic, curated database of restriction enzymes and related proteins.

Mobile Apps

NEB Tools for iPhone®, iPad® or Android®

NEB Tools brings New England Biolabs' most popular web tools to your iPhone, iPad or Android devices.

- Use Enzyme Finder to select a restriction enzyme by category or recognition sequence, or search by name to find information on any NEB enzyme. Sort your results so they make sense to you, then email them to your inbox or connect directly to www.neb.com.
- Use Double Digest Finder to determine buffer and reaction conditions for experiments requiring two restriction enzymes.
- Use Tm Calculator to calculate annealing temperatures for your PCR reaction.
- Also included are several popular calculators from the NEBioCalculator web app.

IPHONE® and IPAD® are registered trademarks of Apple Computers, Inc.
 ANDROID® is a registered trademark of Google, Inc.
 GOOGLE PLAY™ is a trademark of Google, Inc.

NEB Tool for Augmented Reality

Download the NEB Augmented Reality (AR) app for iPhone or iPad at the Apple® App Store or for Android on Google Play™.

NEBnow Locator

NEBnow Freezer Programs are ideally suited for researchers in academics and industry looking for on-site access to the world's finest restriction enzymes and related products. NEB freezers offer you convenience, flexibility and value.

Optimizing Restriction Enzyme Reactions

While standard recommended reaction conditions are a good place to start, in some cases, optimization may be necessary to achieve the best results. Depending on the enzyme(s) being used, variables such as incubation time and number of enzyme units used can be tested to find the optimal reaction conditions for your substrate DNA and enzyme(s) of choice.

Protocol: Restriction Enzyme Reactions

	STANDARD PROTOCOL	TIME-SAVER PROTOCOL
DNA	up to 1 µg	up to 1 µg
10X NEBuffer	5 µl (1X)	5 µl (1X)
Restriction Enzymes	10 units*	1 µl
Total Volume	50 µl	50 µl
Incubation Temperature	Enzyme-dependent	Enzyme-dependent
Incubation Time	60 minutes	5–15 minutes**


*Sufficient to digest all types of DNAs.

**Time-Saver qualified enzymes can also be incubated overnight with no star activity.

TOOLS & RESOURCES

Visit NEBRestrictionEnzymes.com to find:


- Online tutorials for setting up restriction enzyme digests
- Tips to avoid star activity
- Restriction Enzyme Performance Chart
- Troubleshooting guide
- Access to NEB's online tools, including: **Enzyme Finder**, **Double Digest Finder** and **NEBcloner**



TIPS FOR SETTING UP RE DIGESTS

TIPS FOR OPTIMIZATION

Enzyme

- Keep on ice when not in the freezer
- Should be the last component added to reaction
- Mix components by pipetting the reaction mixture up and down, or by “flicking” the reaction tube. Follow with a quick (“touch”) spin-down in a microcentrifuge. Do not vortex the reaction.
- In general, we recommend 5–10 units of enzyme per µg DNA, and 10–20 units per µg of genomic DNA in a 1 hour digest
- Some restriction enzymes require more than one recognition site to cleave efficiently. These are designated with the “multi-site” icon . Please review recommendations on working with these enzymes at www.neb.com.

Star Activity

- Unwanted cleavage that can occur when enzyme is used under sub-optimal conditions, such as:
 - Too much enzyme present
 - Too long of an incubation time
 - Using a non-recommended buffer
 - Glycerol concentrations above 5%
- Star activity can be reduced by using a High-Fidelity (HF®) enzyme, reducing incubation time, using a Time-Saver™ enzyme or increasing reaction volume

DNA

- Should be free of contaminants such as phenol, chloroform, alcohol, EDTA, detergents and salts. Extra wash steps during purification are recommended.
- Methylation of DNA can affect digestion with certain enzymes. For more information about methylation visit www.neb.com/methylation

Buffer

- Use at a 1X concentration
- Recombinant Albumin is included in NEBuffer r1.1, r2.1, r3.1 and rCutSmart™ Buffer. BSA is included in NEBuffer 1.1, 2.1, 3.1 and CutSmart Buffer. In either case, no additional BSA is needed.
- Restriction enzymes that do not require BSA for optimal activity are not adversely affected if BSA is present in the reaction

Reaction Volume

- A 50 µl reaction volume is recommended for digestion of up to 1 µg of substrate. This helps maintain salt levels introduced by miniprep DNA low enough that they don't affect enzyme activity.
- Enzyme volume should not exceed 10% of the total reaction volume to prevent star activity due to excess glycerol
- Additives in the restriction enzyme storage buffer (e.g., glycerol, salt), as well as contaminants found in the substrate solution (e.g., salt, EDTA or alcohol), can be problematic in smaller reaction volumes

	RESTRICTION ENZYME*	DNA	10X NEBUFFER
10 µl rxn**	1 unit	0.1 µg	1 µl
25 µl rxn	5 units	0.5 µg	2.5 µl
50 µl rxn	10 units	1 µg	5 µl

* Restriction enzymes can be diluted using the recommended diluent buffer when smaller amounts are needed

** 10 µl rxns should not be incubated for longer than 1 hour to avoid evaporation.

Incubation Time

- Incubation time for the Standard Protocol is 1 hour. Incubation for the Time-Saver Protocol is 5–15 minutes.
- Visit www.neb.com/timesaver for list of Time-Saver qualified enzymes
- It is possible, with many enzymes, to use fewer units and digest for up to 16 hours. For more information, visit www.neb.com

Storage

- Storage at –20°C is recommended for most restriction enzymes. For a few enzymes, storage at –80°C is recommended. Visit www.neb.com for storage information.
- 10X NEBuffers should be stored at –20°C

Stability

- The expiration date is found on the label
- Long term exposure to temperatures above –20°C should be minimized whenever possible

Double Digestion

Digesting a DNA substrate with two restriction enzymes simultaneously (double digestion) is a common timesaving procedure. Over 210 restriction enzymes are 100% active in rCutSmart or CutSmart Buffer, making double digestion simple. If you are using an enzyme that is not supplied with rCutSmart or CutSmart Buffer, the Performance Chart for Restriction Enzymes (pages 309–314) rates the percentage activity of each restriction endonuclease in the four standard NEBuffers.

Setting up a Double Digestion

- Double digests with CutSmart restriction enzymes can be set up in rCutSmart or CutSmart Buffer. Otherwise, choose an NEBuffer that results in the most activity for both enzymes. If star activity is a concern, consider using one of our High-Fidelity (HF) enzymes.
- Set up reaction according to recommended protocol (see page 306). The final concentration of glycerol in any reaction should be less than 5% to minimize the possibility of star activity (see page 316). For example, in a 50 µl reaction, the total amount of enzyme added should not exceed 5 µl. NEBcloner can also be used to determine recommended double-digest conditions.
- If two different incubation temperatures are necessary, choose the optimal reaction buffer and set up reaction accordingly. Add the first enzyme and incubate at the desired temperature. Then, heat inactivate the first enzyme, if it can be heat inactivated, add the second enzyme and incubate at the recommended temperature.

- Depending on an enzyme's activity rating in a non-optimal NEBuffer, the number of units or incubation time may be adjusted to compensate for the slower rate of cleavage. The Performance Chart for Restriction Enzymes indicates if star activity is an issue with sub-optimal buffers.

Setting up a Double Digestion with a Unique Buffer (designated "U")

- NEB currently supplies three enzymes with unique buffers: EcoRI, SspI and DpnII. In most cases, DpnII requires a sequential digest. Note that EcoRI and SspI have HF versions (NEB #R3101 and NEB #R3132, respectively) which is supplied with rCutSmart or CutSmart Buffer.

Setting up a Sequential Digestion

- If there is no buffer in which the two enzymes exhibit > 50% activity, a sequential digest can be performed.
- Set up a reaction using the restriction endonuclease that has the lowest salt concentration in its recommended buffer and incubate to completion.
- Adjust the salt concentration of the reaction (using a small volume of a concentrated salt solution) to approximate the reaction conditions of the second restriction endonuclease.
- Add the second enzyme and incubate to complete the second reaction.
- Alternatively, a spin column can be used to isolate the DNA prior to the second reaction. NEB recommends using our Monarch Nucleic Acid Purification Kits, see pages 132–141.

TOOLS & RESOURCES

Visit www.neb.com/nebtools for:

- Help choosing double digest conditions using NEB's, **Double Digest Finder** and **NEBcloner®**



TIPS FOR SETTING UP DOUBLE DIGESTS

Types of Restriction Enzymes

Restriction enzymes are traditionally classified into four types on the basis of subunit composition, cleavage position, sequence specificity and cofactor requirements. However, amino acid sequencing has uncovered extraordinary variety among restriction enzymes and revealed that at the molecular level there are many more than four different types.

Type I Enzymes are complex, multisubunit, combination restriction-and-modification enzymes that cut DNA at random far from their recognition sequences. Type I enzymes are of considerable biochemical interest, but they have little practical value since they do not produce discrete restriction fragments or distinct gel-banding patterns.

Type II Enzymes cut DNA at defined positions close to or within their recognition sequences. They produce discrete restriction fragments and distinct gel banding patterns, and they are the predominant class used in the laboratory for DNA analysis and gene cloning. Type IIS restriction enzymes recognize asymmetric DNA sequences and cleave outside of their recognition sequences. They are useful for many applications, including Golden Gate Assembly.

Rather than forming a single family of related proteins, Type II enzymes are a collection of unrelated proteins of many different sorts. Type II enzymes frequently differ greatly in amino acid sequence from one another, and indeed from every other known protein, that they exemplify the class of rapidly evolving proteins that are often indicative of involvement in host-parasite interactions.

Type III Enzymes are also large combination restriction-and-modification enzymes. They cleave outside of their recognition sequences and require two such sequences in opposite orientations within the same DNA molecule to accomplish cleavage; they rarely yield complete digests.

Type IV Enzymes recognize modified, typically methylated DNA and are exemplified by the McrBC and Mrr systems of *E. coli*.

TOOLS & RESOURCES

Visit the video library at www.neb.com to find:

- Tutorials on Type I, II and III restriction enzymes



TYPE I, II AND III RESTRICTION ENZYMES

View
double digest
protocol.



Restriction Enzyme Troubleshooting Guide

PROBLEM	CAUSE	SOLUTION
Few or no transformants	Restriction enzyme(s) didn't cleave completely	<ul style="list-style-type: none"> • Check the methylation sensitivity of the enzyme(s) to determine if the enzyme is blocked by methylation of the recognition sequence • Use the recommended buffer supplied with the restriction enzyme • Clean up the DNA to remove any contaminants that may inhibit the enzyme • When digesting a PCR fragment, make sure to have at least 6 nucleotides between the recognition site and the end of the DNA molecule
	The restriction enzyme(s) is bound to the substrate DNA	<ul style="list-style-type: none"> • Lower the number of units • Add SDS (0.1–0.5%) to the loading buffer to dissociate the enzyme from the DNA or use Gel Loading Dye, Purple (6X) (NEB #B7024)
The digested DNA ran as a smear on an agarose gel	Nuclease contamination	<ul style="list-style-type: none"> • Use fresh, clean running buffer and a fresh agarose gel • Clean up the DNA
	Cleavage is blocked by methylation	<ul style="list-style-type: none"> • Check the methylation sensitivity of the enzyme(s) to determine if the enzyme is blocked by methylation of the recognition sequence • DNA isolated from a bacterial source may be blocked by Dam and Dcm methylation • If the enzyme is inhibited by Dam or Dcm methylation, grow the plasmid in a <i>dam-/dcm-</i> strain (NEB #C2925) • DNA isolated from eukaryotic source may be blocked by CpG methylation
Incomplete restriction enzyme digestion	Salt inhibition	<ul style="list-style-type: none"> • Enzymes that have low activity in salt-containing buffers (NEBuffer r3.1 or 3.1) may be salt sensitive, so clean up the DNA prior to digestion • DNA purification procedures that use spin columns can result in high salt levels, which inhibit enzyme activity. To prevent this, DNA solution should be no more than 25% of total reaction volume.
	Inhibition by PCR components	<ul style="list-style-type: none"> • Clean up the PCR fragment prior to restriction digest
	Using the wrong buffer	<ul style="list-style-type: none"> • Use the recommended buffer supplied with the restriction enzyme
	Too few units of enzyme used	<ul style="list-style-type: none"> • Use at least 3–5 units of enzyme per µg of DNA
	Incubation time was too short	<ul style="list-style-type: none"> • Increase the incubation time
	Digesting supercoiled DNA	<ul style="list-style-type: none"> • Some enzymes have a lower activity on supercoiled DNA. Increase the number of enzyme units in the reaction.
	Presence of slow sites	<ul style="list-style-type: none"> • Some enzymes can exhibit slower cleavage towards specific sites. Increase the incubation time, 1–2 hours is typically sufficient.
	Two sites required	<ul style="list-style-type: none"> • Some enzymes require the presence of two recognition sites to cut efficiently
	DNA is contaminated with an inhibitor	<ul style="list-style-type: none"> • Assay substrate DNA in the presence of a control DNA. Control DNA will not cleave if there is an inhibitor present. Miniprep DNA is particularly susceptible to contaminants. • Clean DNA with a spin column or increase volume to dilute contaminant. We recommend the Monarch PCR & DNA Cleanup Kit (NEB #T1030),
Extra bands in the gel	If larger bands than expected are seen in the gel, this may indicate binding of the enzyme(s) to the substrate	<ul style="list-style-type: none"> • Lower the number of units in the reaction • Add SDS (0.1–0.5%) to the loading buffer to dissociate the enzyme from the substrate, or use Gel Loading Dye, Purple (6X) (NEB #B7024)
	Star activity	<ul style="list-style-type: none"> • Use the recommended buffer supplied with the restriction enzyme • Decrease the number of enzyme units in the reaction • Make sure the amount of enzyme added does not exceed 10% of the total reaction volume. This ensures that the total glycerol concentration does not exceed 5% v/v. • Decrease the incubation time. Using the minimum reaction time required for complete digestion will help prevent star activity. • Try using a High-Fidelity (HF) restriction enzyme. HF enzymes have been engineered for reduced star activity.
	Partial restriction enzyme digest	<ul style="list-style-type: none"> • Enzymes that have low activity in salt-containing buffers (e.g., NEBuffer r3.1 or 3.1) may be salt sensitive. Make sure to clean up the DNA prior to digestion. • DNA purification procedures that use spin columns can result in high salt levels, which inhibit enzyme activity. To prevent this, DNA solution should be no more than 25% of total reaction volume. • Clean-up the PCR fragment prior to restriction digest • Use the recommended buffer supplied with the restriction enzyme • Use at least 5–10 units of enzyme per µg of DNA and digest the DNA for 1–2 hours

Performance Chart for Restriction Enzymes

New England Biolabs supplies > 215 restriction enzymes that are 100% active in rCutSmart Buffer. This results in increased efficiency, flexibility and ease-of-use, especially when performing double digests.

This performance chart summarizes the activity information of NEB restriction enzymes. To help select the best conditions for double digests, this chart shows the optimal (supplied) NEBuffer and approximate activity in the four standard NEBuffers for each enzyme. Note that Recombinant Albumin or BSA is included in all NEBuffers, and is not provided as a separate tube. In addition, this performance chart shows recommended reaction temperature, heat-inactivation temperature, recommended diluent buffer, methylation sensitivity and whether the enzyme is Time-Saver qualified (e.g., cleaves substrate in 5–15 minutes under recommended conditions, and can be used overnight without degradation of DNA), and whether the enzyme works better in a substrate with multiple sites.

Chart Legend

U Supplied with a unique reaction buffer that is different from the four standard NEBuffers. The compatibility with the four standard NEBuffers is indicated in the chart.

RH Recombinant

e Engineered enzyme for maximum performance

TS Time-Saver qualified

2⁺ sites Indicates that the restriction enzyme requires two or more sites for cleavage

SAM Supplied with a separate vial of S-adenosylmethionine (SAM). To obtain 100% activity, SAM should be added to the 1X reaction mix as specified on the product data card.

dcm dcm methylation sensitivity

dam dam methylation sensitivity

CpG CpG methylation sensitivity

Note: Starting April 2021, NEB will begin the transition to supply all restriction enzymes with Recombinant Albumin (rAlbumin)-containing buffers (r1.1, r2.1, r3.1 & rCutSmart). There is no difference in enzyme performance when using rAlbumin- or BSA-containing buffers. We have reported the supplied rAlbumin-containing buffer, but show % activity for both buffer sets.

Activity Notes (see last column)

FOR STAR ACTIVITY


































1. Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.
2. Star activity may result from extended digestion.
3. Star activity may result from a glycerol concentration of > 5%.

* May exhibit star activity in this buffer.

+ For added flexibility, NEB offers an isoschizomer or HF enzyme, supplied with rCutSmart Buffer or CutSmart Buffer.

FOR LIGATION AND RECUTTING

- a. Ligation is less than 10%
- b. Ligation is 25% – 75%
- c. Recutting after ligation is < 5%
- d. Recutting after ligation is 50% – 75%
- e. Ligation and recutting after ligation is not applicable since the enzyme is either a nicking enzyme, is affected by methylation, or the recognition sequence contains variable sequences.

	ENZYME	SUPPLIED NEBUFFER	% ACTIVITY IN NEBUFFERS				INCUB. TEMP. (°C)	INACTIV. TEMP. (°C)	DIL.	UNIT SUBSTRATE	METHYLATION SENSITIVITY	NOTE(S)
			r1.1 1.1	r2.1 2.1	r3.1 3.1	rCUTSMART CUTSMART						
	AatII	rCutSmart	< 10	50*	50	100	37°	80°	B	Lambda	CpG	
	AbaSI	rCutSmart + DTT	25	50	50	100	25°	65°	C	T4 wt Phage		e
	AccI	rCutSmart	50	50	10	100	37°	80°	A	Lambda	CpG	
	Acc65I	r3.1	10	75*	100	25	37°	65°	A	pBC4	dcm CpG	
	Acil	rCutSmart	< 10	25	100	100	37°	65°	A	Lambda	CpG	d
	AcII	rCutSmart	< 10	< 10	< 10	100	37°	No	B	Lambda	CpG	
	AcuI	rCutSmart	50	100	50	100	37°	65°	B	Lambda		1, b, d
	AfeI	rCutSmart	25	100	25	100	37°	65°	B	pXba	CpG	
	AflII	rCutSmart	50	100	10	100	37°	65°	A	phiX174		
	AflIII	r3.1	10	50	100	50	37°	80°	B	Lambda		
	AgeI	r1.1	100	75	25	75	37°	65°	C	Lambda	CpG	
	AgeI-HF	rCutSmart	100	50	10	100	37°	65°	A	Lambda	CpG	
	AhdI	rCutSmart	25	25	10	100	37°	65°	A	Lambda	CpG	a
	AleI-v2	rCutSmart	< 10	< 10	< 10	100	37°	65°	B	Lambda	CpG	
	AluI	rCutSmart	25	100	50	100	37°	80°	B	Lambda		b
	AlwI	rCutSmart	50	50	10	100	37°	No	A	Lambda dam-	dam	1, b, d
	AlwNI	rCutSmart	10	100	50	100	37°	80°	A	Lambda	dcm	
	Apal	rCutSmart	25	25	< 10	100	25°	65°	A	pXba	dcm CpG	
	ApaLI	rCutSmart	100	100	10	100	37°	No	A	Lambda HindIII	CpG	
	ApeKI	r3.1	25	50	100	10	75°	No	B	Lambda	CpG	
	ApoI	r3.1	10	75	100	75	50°	80°	A	Lambda		
	ApoI-HF	rCutSmart	10	100	10	100	37°	80°	B	Lambda		
	AscI	rCutSmart	< 10	10	10	100	37°	80°	A	Lambda	CpG	
	Asel	r3.1	< 10	50*	100	10	37°	65°	B	Lambda		3
	AsiSI	rCutSmart	100	100	25	100	37°	80°	B	pXba (Xho digested)	CpG	2, b
	AvaI	rCutSmart	< 10	100	25	100	37°	80°	A	Lambda	CpG	
	AvaII	rCutSmart	50	75	10	100	37°	80°	A	Lambda	dcm CpG	
	AvrII	rCutSmart	100	50	50	100	37°	No	B	Lambda HindIII		
	BaeI	rCutSmart + SAM	50	100	50	100	25°	65°	A	Lambda	CpG	e
	BaeGI	r3.1	75	75	100	25	37°	80°	A	Lambda		
	BamHI	r3.1	75*	100*	100	100*	37°	No	A	Lambda		3
	BamHI-HF	rCutSmart	100	50	10	100	37°	No	A	Lambda		
	BanI	rCutSmart	10	25	< 10	100	37°	65°	A	Lambda	dcm CpG	1

Performance Chart for Restriction Enzymes (continued)

	ENZYME	SUPPLIED NEBUFFER	% ACTIVITY IN NEBUFFERS				INCUB. TEMP. (°C)	INACTIV. TEMP. (°C)	DIL.	UNIT SUBSTRATE	METHYLATION SENSITIVITY	NOTE(S)
			r1.1 1.1	r2.1 2.1	r3.1 3.1	rCUTSMART CUTSMART						
	BanII	rCutSmart	100	100	50	100	37°	80°	A	Lambda		2
	BbsI	r2.1	100	100	25	75	37°	65°	B	Lambda		
	BbsI-HF	rCutSmart	10	10	< 10	100	37°	65°	B	Lambda		
	BbvI	rCutSmart	100	100	25	100	37°	65°	B	pBR322		3
	BbvCI	rCutSmart	10	100	50	100	37°	No	B	Lambda	CpG	1, a
	BccI	rCutSmart	100	50	10	100	37°	65°	A	pXba		3, b
	BceAI	r3.1	100*	100*	100	100*	37°	65°	A	pBR322	CpG	1
	BcgI	r3.1	10	75*	100	50*	37°	65°	A	Lambda	dam CpG	e
	BciVI	rCutSmart	100	25	< 10	100	37°	80°	C	Lambda		b
	BclI	r3.1	50	100	100	75	50°	No	A	Lambda dam-	dam	
	BclI-HF	rCutSmart	100	100	10	100	37°	65°	B	Lambda dam-	dam	
	BcoDI	rCutSmart	50	75	75	100	37°	No	B	Lambda	CpG	
	BfaI	rCutSmart	< 10	10	< 10	100	37°	80°	B	Lambda		2, b
	BfuAI	r3.1	< 10	25	100	10	50°	65°	B	Lambda	CpG	3
	BglI	r3.1	10	25	100	10	37°	65°	B	Lambda	CpG	
	BglII	r3.1	10	10	100	< 10	37°	No	A	Lambda		
	BlpI	rCutSmart	50	100	10	100	37°	No	A	Lambda		d
	BmgBI	r3.1	< 10	10	100	10	37°	65°	B	Lambda	CpG	3, b, d
	BmrI	r2.1	75	100	75	100*	37°	65°	B	Lambda HindIII		b
	BmtI	r3.1	100	100	100	100+	37°	65°	B	pXba		2
	BmtI-HF	rCutSmart	50	100	10	100	37°	65°	B	pXba		
	BpmI	r3.1	75	100	100	100*	37°	65°	B	Lambda		2, d
	Bpu10I	r3.1	10	25	100	25	37°	80°	B	Lambda		3, b, d
	BpuEI	rCutSmart	50*	100	50*	100	37°	65°	B	Lambda		d
	BsaI-HFv2	rCutSmart	100	100	100	100	37°	80°	B	pXba	dam CpG	
	BsaAI	rCutSmart	100	100	100	100	37°	No	C	Lambda	CpG	
	BsaBI	rCutSmart	50	100	75	100	60°	80°	B	Lambda dam-	dam CpG	2
	BsaHI	rCutSmart	50	100	100	100	37°	80°	C	Lambda	dcm CpG	
	BsaJI	rCutSmart	50	100	100	100	60°	80°	A	Lambda		
	BsaWI	rCutSmart	10	100	50	100	60°	80°	A	Lambda		
	BsaXI	rCutSmart	50*	100*	10	100	37°	No	C	Lambda		e
	BseRI	rCutSmart	100	100	75	100	37°	80°	A	Lambda		d
	BseYI	r3.1	10	50	100	50	37°	80°	B	Lambda	CpG	d
	BsgI	rCutSmart	25	50	25	100	37°	65°	B	Lambda		d
	BsiEI	rCutSmart	25	50	< 10	100	60°	No	A	Lambda	CpG	
	BsiHKA1	rCutSmart	25	100	100	100	65°	No	A	Lambda		
	BsiWI	r3.1	25	50*	100	25	55°	65°	B	phiX174	CpG	
	BsiWI-HF	rCutSmart	50	100	10	100	37°	No	B	phiX174	CpG	
	BslI	rCutSmart	50	75	100	100	55°	No	A	Lambda	dcm CpG	b
	BsmI	rCutSmart	25	100	< 10	100	65°	80°	A	Lambda		
	BsmAI	rCutSmart	50	100	100	100	55°	No	B	Lambda	CpG	
	BsmBI-v2	r3.1	< 10	50	100	25	55°	80°	B	Lambda	CpG	
	BsmFI	rCutSmart	25	50	50	100	65°	80°	A	pBR322	dcm CpG	1
	BsoBI	rCutSmart	25	100	100	100	37°	80°	A	Lambda		
	Bsp1286I	rCutSmart	25	25	25	100	37°	65°	A	Lambda		3
	BspCNI	rCutSmart	100	75	10	100	37°	80°	A	Lambda		b
	BspDI	rCutSmart	25	75	50	100	37°	80°	A	Lambda	dam CpG	
	BspEI	r3.1	< 10	10	100	< 10	37°	80°	B	Lambda dam-	dam CpG	
	BspHI	rCutSmart	10	50	25	100	37°	80°	A	Lambda	dam	
	BspMI	r3.1	10	50*	100	10	37°	65°	B	Lambda		
	BspQI	r3.1	100*	100*	100	100*	50°	80°	B	Lambda		3

a. Ligation is less than 10%
b. Ligation is 25% – 75%

c. Recutting after ligation is < 5%
d. Recutting after ligation is 50% – 75%

e. Ligation and recutting after ligation is not applicable since the enzyme is either a nicking enzyme, is affected by methylation, or the recognition sequence contains variable sequences.




















































































	ENZYME	SUPPLIED NEBUFFER	% ACTIVITY IN NEBUFFERS				INCUB. TEMP. (°C)	INACTIV. TEMP. (°C)	DIL.	UNIT SUBSTRATE	METHYLATION SENSITIVITY	NOTE(S)
			r1.1 1.1	r2.1 2.1	r3.1 3.1	rCUTSMART CUTSMART						
	BsrI	r3.1	< 10	50	100	10	65°	80°	B	phiX174		b
	BsrBI	rCutSmart	50	100	100	100	37°	80°	A	Lambda	CpG	d
	BsrDI	r2.1	10	100	75	25	65°	80°	A	Lambda		3, d
	BsrFI-v2	rCutSmart	25	25	0	100	37°	No	C	pBR322	CpG	
	BsrGI	r2.1	25	100	100	25	37°	80°	A	Lambda		
	BsrGI-HF	rCutSmart	10	100	100	100	37°	80°	A	Lambda		
	BssHII	rCutSmart	100	100	100	100	50°	65°	B	Lambda	CpG	
	BssSI-v2	rCutSmart	10	25	< 10	100	37°	No	B	Lambda		
	BstAPI	rCutSmart	50	100	25	100	60°	80°	A	Lambda	CpG	b
	BstBI	rCutSmart	75	100	10	100	65°	No	A	Lambda	CpG	
	BstEII	r3.1	10	75*	100	75*	60°	No	A	Lambda		3
	BstEII-HF	rCutSmart	< 10	10	< 10	100	37°	No	A	Lambda		
	BstNI	r3.1	10	100	100	75	60°	No	A	Lambda		a
	BstUI	rCutSmart	50	100	25	100	60°	No	A	Lambda	CpG	b
	BstXI	r3.1	< 10	50	100	25	37°	80°	B	Lambda	dcm	3
	BstYI	r2.1	25	100	75	100	60°	No	A	Lambda		
	BstZ171-HF	rCutSmart	100	100	10	100	37°	No	A	Lambda	CpG	b
	Bsu36I	rCutSmart	25	100	100	100	37°	80°	C	Lambda HindIII		b
	BtgI	rCutSmart	50	100	100	100	37°	80°	B	pBR322		
	BtgZI	rCutSmart	10	25	< 10	100	60°	80°	A	Lambda	CpG	3, b, d
	BtsI-v2	rCutSmart	100	100	25	100	37°	No	A	Lambda		1
	BtsIMutI	rCutSmart	100	50	10	100	55°	80°	A	pUC19		b
	BtsCI	rCutSmart	10	100	25	100	50°	80°	B	Lambda		
	Cac8I	rCutSmart	50	75	100	100	37°	65°	B	Lambda	CpG	b
	CiaI	rCutSmart	10	50	50	100	37°	65°	A	Lambda dam-	dam	CpG
	CspCI	rCutSmart	10	100	10	100	37°	65°	A	Lambda		e
	CviAI	rCutSmart	50	50	10	100	25°	65°	C	Lambda		
	CviKI-1	rCutSmart	25	100	100	100	37°	No	A	pBR322		1, b
	CviQI	r3.1	75	100*	100	75*	25°	No	C	Lambda		b
	DdeI	rCutSmart	75	100	100	100	37°	65°	B	Lambda		
	DpnI	rCutSmart	100	100	75	100	37°	80°	B	pBR322	CpG	b
	DpnII	U	25	25	100*	25	37°	65°	B	Lambda dam-	dam	
	DraI	rCutSmart	75	75	50	100	37°	65°	A	Lambda		
	DraIII-HF	rCutSmart	< 10	50	10	100	37°	No	B	Lambda	CpG	b
	DrdI	rCutSmart	25	50	10	100	37°	65°	A	pUC19	CpG	3
	EaeI	rCutSmart	10	50	< 10	100	37°	65°	A	Lambda	dcm	CpG
	EagI	r3.1	10	25	100	10	37°	65°	B	pXba	CpG	
	EagI-HF	rCutSmart	25	100	100	100	37°	65°	B	pXba	CpG	
	EarI	rCutSmart	50	10	< 10	100	37°	65°	B	Lambda	CpG	b, d
	EciI	rCutSmart	100	50	50	100	37°	65°	A	Lambda	CpG	2
	Eco53kl	rCutSmart	100	100	< 10	100	37°	65°	A	pXba	CpG	3, b
	EcoNI	rCutSmart	50	100	75	100	37°	65°	A	Lambda		b
	EcoO109I	rCutSmart	50	100	50	100	37°	65°	A	Lambda HindIII	dcm	3
	EcoP15I	r3.1 + ATP	75	100	100	100	37°	65°	A	pUC19		e
	EcoRI	U	25	100*	50	50*	37°	65°	C	Lambda	CpG	
	EcoRI-HF	rCutSmart	10	100	< 10	100	37°	65°	C	Lambda	CpG	
	EcoRV	r3.1	10	50	100	10	37°	80°	A	Lambda		
	EcoRV-HF	rCutSmart	25	100	100	100	37°	65°	B	Lambda	CpG	
	Esp3I	rCutSmart	100	100	< 10	100	37°	65°	B	Lambda	CpG	
	FatI	r2.1	10	100	50	50	55°	80°	A	pUC19		
	FauI	rCutSmart	100	50	10	100	55°	65°	A	Lambda	CpG	3, b, d

1. Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

2. Star activity may result from extended digestion.
3. Star activity may result from a glycerol concentration of > 5%.

* May exhibit star activity in this buffer.
+ NEB isoschizomer or HF enzyme supplied with buffer.







































Performance Chart for Restriction Enzymes (continued)

	ENZYME	SUPPLIED NEBUFFER	% ACTIVITY IN NEBUFFERS				INCUB. TEMP. (°C)	INACTIV. TEMP. (°C)	DIL.	UNIT SUBSTRATE	METHYLATION SENSITIVITY	NOTE(S)
			r1.1 1.1	r2.1 2.1	r3.1 3.1	rCUTSMART CUTSMART						
	Fnu4HI	rCutSmart	< 10	< 10	< 10	100	37°	No	A	Lambda		a
	FokI	rCutSmart	100	100	75	100	37°	65°	A	Lambda		3, b, d
	FseI	rCutSmart	100	75	< 10	100	37°	65°	B	pBC4		
	FspI	rCutSmart	10	100	10	100	37°	No	C	Lambda		b
	FspEI	rCutSmart + Enz. Activ.	< 10	< 10	< 10	100	37°	80°	B	pBR322		1, e
	HaeII	rCutSmart	25	100	10	100	37°	80°	A	Lambda		
	HaeIII	rCutSmart	50	100	25	100	37°	80°	A	Lambda		
	HgaI	r1.1	100	100	25	100*	37°	65°	A	phiX174		1
	HhaI	rCutSmart	25	100	100	100	37°	65°	A	Lambda		
	HincII	r3.1	25	100	100	100	37°	65°	B	Lambda		
	HindIII	r2.1	25	100	50	50	37°	80°	B	Lambda		2
	HindIII-HF	rCutSmart	10	100	10	100	37°	80°	B	Lambda		
	HinfI	rCutSmart	50	100	100	100	37°	80°	A	Lambda		
	HinP1I	rCutSmart	100	100	100	100	37°	65°	A	Lambda		
	HpaI	rCutSmart	< 10	75*	25	100	37°	No	A	Lambda		1
	HpaII	rCutSmart	100	50	< 10	100	37°	80°	A	Lambda		
	HphI	rCutSmart	50	50	< 10	100	37°	65°	B	Lambda		1, b, d
	Hpy99I	rCutSmart	50	10	< 10	100	37°	65°	A	Lambda		
	Hpy166II	rCutSmart	100	100	50	100	37°	65°	C	pBR322		
	Hpy188I	rCutSmart	25	100	50	100	37°	65°	A	pBR322		1, b
	Hpy188III	rCutSmart	100	100	10	100	37°	65°	B	pUC19		3, b
	HpyAV	rCutSmart	100	100	25	100	37°	65°		Lambda		3, b, d
	HpyCH4III	rCutSmart	100	25	< 10	100	37°	65°	A	Lambda		b
	HpyCH4IV	rCutSmart	100	50	25	100	37°	65°	A	pUC19		
	HpyCH4V	rCutSmart	50	50	25	100	37°	65°	A	Lambda		
	I-CeuI	rCutSmart	10	10	10	100	37°	65°	B	pBHS Scal-linearized		
	I-SceI	rCutSmart	10	50	25	100	37°	65°	B	pGPS2 NotI-linearized		
	KasI	rCutSmart	50	100	50	100	37°	65°	B	pBR322		3
	KpnI	r1.1	100	75	< 10	50	37°	No	A	pXba		1
	KpnI-HF	rCutSmart	100	25	< 10	100	37°	No	A	pXba		
	LpnPI	rCutSmart + Enz. Activ.	< 10	< 10	< 10	100	37°	65°	B	pBR322		1, e
	MboI	rCutSmart	75	100	100	100	37°	65°	A	Lambda dam-		
	MboII	rCutSmart	100*	100	50	100	37°	65°	C	Lambda dam-		b
	MfeI	rCutSmart	75	50	10	100	37°	No	A	Lambda		2
	MfeI-HF	rCutSmart	75	25	< 10	100	37°	No	A	Lambda		
	MluI	r3.1	10	50	100	25	37°	80°	A	Lambda		
	MluI-HF	rCutSmart	25	100	100	100	37°	No	A	Lambda		
	MluCI	rCutSmart	100	10	10	100	37°	No	A	Lambda		
	MlyI	rCutSmart	50	50	10	100	37°	65°	A	Lambda		b, d
	MmeI	rCutSmart	50	100	50	100	37°	65°	B	phiX174		b, c
	MnII	rCutSmart	75	100	50	100	37°	65°	B	Lambda		b
	MscI	rCutSmart	25	100	100	100	37°	80°	C	Lambda		
	MseI	rCutSmart	75	100	75	100	37°	65°	A	Lambda		
	MsiI	rCutSmart	50	50	< 10	100	37°	80°	A	Lambda		
	MspI	rCutSmart	75	100	50	100	37°	No	A	Lambda		
	MspA1I	rCutSmart	10	50	10	100	37°	65°	B	Lambda		
	MspJI	rCutSmart + Enz. Activ.	< 10	< 10	< 10	100	37°	65°	B	pBR322		1, e
	MwoI	rCutSmart	< 10	100	100	100	60°	No	B	Lambda		
	NaeI	rCutSmart	25	25	< 10	100	37°	No	A	pXba		b
	NarI	rCutSmart	100	100	10	100	37°	65°	A	pXba		
	Nb.BbvCI	rCutSmart	25	100	100	100	37°	80°	A	pUB		e

a. Ligation is less than 10%
b. Ligation is 25% – 75%

c. Recutting after ligation is < 5%
d. Recutting after ligation is 50% – 75%

e. Ligation and recutting after ligation is not applicable since the enzyme is either a nicking enzyme, is affected by methylation, or the recognition sequence contains variable sequences.

ENZYME	SUPPLIED NEBUFFER	% ACTIVITY IN NEBUFFERS			rCUTSMART CUTSMART	INCUB. TEMP. (°C)	INACTIV. TEMP. (°C)	DIL.	UNIT SUBSTRATE	METHYLATION SENSITIVITY	NOTE(S)
		r1.1 1.1	r2.1 2.1	r3.1 3.1							
 Nbsml	r3.1	< 10	50	100	10	65°	80°	A	pBR322		e
 Nbsrdi	rCutSmart	25	100	100	100	65°	80°	A	pUC19		e
 Nbsssi	r3.1	10	100	100	25	37°	No	B	pUC19		e
 Nbtsl	rCutSmart	75	100	75	100	37°	80°	A	phiX174		e
 Ncil	rCutSmart	100	25	10	100	37°	No	A	Lambda	CpG	b
 Ncoi	r3.1	100	100	100	100+	37°	80°	A	Lambda		
 Ncoi-HF	rCutSmart	50	100	10	100	37°	80°	B	Lambda		
 NdeI	rCutSmart	75	100	100	100	37°	65°	A	Lambda		
 NdeI 2'alt	rCutSmart	100	50	10	100	37°	No	A	pXba	CpG	1
 NheI	r2.1	100	100	10	100+	37°	65°	C	Lambda HindIII	CpG	
 NheI-HF	rCutSmart	100	25	< 10	100	37°	80°	C	Lambda HindIII	CpG	
 NlaIII	rCutSmart	< 10	< 10	< 10	100	37°	65°	B	phiX174		
 NlaIV	rCutSmart	10	10	10	100	37°	65°	B	pBR322	dcm CpG	
 NmeAIII	rCutSmart	10	10	< 10	100	37°	65°	B	phiX174		c
 NotI	r3.1	< 10	50	100	25	37°	65°	C	pBC4	CpG	
 NotI-HF	rCutSmart	25	100	25	100	37°	65°	A	pBC4	CpG	
 NruI	r3.1	< 10	10	100	10	37°	No	A	Lambda	dam CpG	b
 NruI-HF	rCutSmart	0	25	50	100	37°	No	A	Lambda	dam CpG	
 NsiI	r3.1	10	75	100	25	37°	65°	B	Lambda		
 NsiI-HF	rCutSmart	< 10	20	< 10	100	37°	80°	B	Lambda		
 NspI	rCutSmart	100	100	< 10	100	37°	65°	A	Lambda		
 NtAlwI	rCutSmart	10	100	100	100	37°	80°	A	pUC101 dam-dcm-	dam	e
 NtBbvCI	rCutSmart	50	100	10	100	37°	80°	A	pUB	CpG	e
 NtBsmAI	rCutSmart	100	50	10	100	37°	65°	A	pBR322	CpG	e
 NtBspQI	r3.1	< 10	25	100	10	50°	80°	B	pUC19		e
 NtBstNBI	r3.1	0	10	100	10	55°	80°	A	T7		e
 NtCviPII	rCutSmart	10	100	25	100	37°	65°	A	pUC19	CpG	e
 PacI	rCutSmart	100	75	10	100	37°	65°	A	pNEB193		
 PaeR7I	rCutSmart	25	100	10	100	37°	No	A	Lambda HindIII	CpG	
 PaeR7I 2'alt	rCutSmart	< 10	100	10	100	37°	65°	B	Lambda	CpG	1
 PciI	r3.1	50	75	100	50*	37°	80°	B	pXba		
 PfiFI	rCutSmart	25	100	25	100	37°	65°	A	pBC4		b
 PfiIMI	r3.1	0	100	100	50	37°	65°	A	Lambda	dcm	3, b, d
 PI-PspI	U	10	10	10	10	65°	No	B	pAKR XmnI		
 PI-Scel	U	10	10	10	10	37°	65°	B	pBSvdeX XmnI		
 PIel 2'alt	rCutSmart	25	50	25	100	37°	65°	A	Lambda	CpG	b, d
 PluTI 2'alt	rCutSmart	100	25	< 10	100	37°	65°	A	pXba	CpG	
 PmeI	rCutSmart	< 10	50	10	100	37°	65°	A	Lambda	CpG	
 PmiI	rCutSmart	100	50	< 10	100	37°	65°	A	Lambda HindIII	CpG	
 PpuMI	rCutSmart	< 10	< 10	< 10	100	37°	No	B	Lambda HindIII	dcm	
 PshAI	rCutSmart	25	50	10	100	37°	65°	A	Lambda	CpG	
 Psil-v2	rCutSmart	25	50	10	100	37°	65°	B	Lambda		3
 PspGI	rCutSmart	25	100	50	100	75°	No	A	T7	dcm	3
 PspOMI	rCutSmart	10	10	< 10	100	37°	65°	B	pXba	dcm CpG	
 PspXI	rCutSmart	< 10	100	25	100	37°	No	B	Lambda HindIII	CpG	
 PstI	r3.1	75	75	100	50*	37°	80°	C	Lambda		
 PstI-HF	rCutSmart	10	75	50	100	37°	No	C	Lambda		
 PvuI	r3.1	< 10	25	100	< 10	37°	No	B	pXba	CpG	
 PvuI-HF	rCutSmart	25	100	100	100	37°	No	B	pXba	CpG	
PvuII	r3.1	50	100	100	100*	37°	No	B	Lambda		
PvuII-HF	rCutSmart	< 10	< 10	< 10	100	37°	No	B	Lambda		

1. Star activity may result from extended digestion, high enzyme concentration or a glycerol concentration of > 5%.

2. Star activity may result from extended digestion.
3. Star activity may result from a glycerol concentration of > 5%.

* May exhibit star activity in this buffer.

* NEB isoschizomer or HF enzyme supplied with buffer.

Performance Chart for Restriction Enzymes (continued)

	ENZYME	SUPPLIED NEBUFFER	% ACTIVITY IN NEBUFFERS				INCUB. TEMP. (°C)	INACTIV. TEMP. (°C)	DIL.	UNIT SUBSTRATE	METHYLATION SENSITIVITY	NOTE(S)
			r1.1 1.1	r2.1 2.1	r3.1 3.1	rCUTSMART CUTSMART						
	RsaI	rCutSmart	25	50	< 10	100	37°	No	A	Lambda	CpG	
	RsrII	rCutSmart	25	75	10	100	37°	65°	C	Lambda	CpG	
	SacI	r1.1	100	50	10	100+	37°	65°	A	Lambda HindIII	CpG	
	SacI-HF	rCutSmart	10	50	< 10	100	37°	65°	A	Lambda HindIII	CpG	
	SacII	rCutSmart	10	100	10	100	37°	65°	A	pXba	CpG	
	SalI	r3.1	< 10	< 10	100	< 10	37°	65°	A	Lambda HindIII	CpG	
	SalI-HF	rCutSmart	10	100	100	100	37°	65°	A	Lambda HindIII	CpG	
	SapI	rCutSmart	75	50	< 10	100	37°	65°	B	Lambda	CpG	
	Sau3AI	r1.1	100	50	10	100+	37°	65°	A	Lambda	CpG	b
	Sau96I	rCutSmart	50	100	100	100	37°	65°	A	Lambda	dcm CpG	
	SbfI	rCutSmart	50	25	< 10	100	37°	80°	A	Lambda	CpG	3
	SbfI-HF	rCutSmart	50	25	< 10	100	37°	80°	B	Lambda	CpG	
	ScaI-HF	rCutSmart	100	100	10	100	37°	80°	B	Lambda	CpG	
	ScrFI	rCutSmart	100	100	100	100	37°	65°	C	Lambda	dcm CpG	2, a
	SexAI	rCutSmart	100	75	50	100	37°	65°	A	pBC4 dcm-	dcm	3, b, d
	SfaNI	r3.1	< 10	75	100	25	37°	65°	B	phiX174	CpG	3, b
	SfiI	rCutSmart	75	50	25	100	37°	65°	B	Lambda	CpG	3
	SfiI	rCutSmart	25	100	50	100	50°	No	C	pXba	dcm CpG	
	SfiI	rCutSmart	50	100	100	100	37°	No	B	Lambda HindIII	dcm CpG	
	SgrAI	rCutSmart	100	100	10	100	37°	65°	A	Lambda	CpG	1
	SmaI	rCutSmart	< 10	< 10	< 10	100	25°	65°	B	Lambda HindIII	CpG	b
	SmlI	rCutSmart	25	75	25	100	55°	No	A	Lambda	CpG	b
	SnaBI	rCutSmart	50*	50	10	100	37°	80°	A	T7	CpG	1
	SpeI	rCutSmart	75	100	25	100	37°	80°	C	Adenovirus-2	CpG	
	SpeI-HF	rCutSmart	25	50	10	100	37°	80°	C	pXba	CpG	
	SphI	r2.1	100	100	50	100+	37°	65°	B	Lambda	CpG	2
	SphI-HF	rCutSmart	50	25	10	100	37°	65°	B	Lambda	CpG	
	SrfI	rCutSmart	10	50	0	100	37°	65°	B	pNEB193-SrfI	CpG	
	SspI	U	50	100	50	50	37°	65°	C	Lambda	CpG	
	SspI-HF	rCutSmart	25	100	< 10	100	37°	65°	B	Lambda	CpG	
	StuI	rCutSmart	50	100	50	100	37°	No	A	Lambda	dcm	
	StyD4I	rCutSmart	10	100	100	100	37°	65°	B	Lambda	dcm CpG	
	StyI	r3.1	10	25	100	10	37°	65°	A	Lambda	CpG	b
	StyI-HF	rCutSmart	25	100	25	100	37°	65°	A	Lambda	CpG	
	Swal	r3.1	10	10	100	10	25°	65°	B	pXba	CpG	b, d
	TaqI-v2	rCutSmart	50	100	50	100	65°	No	B	Lambda	dam	
	TfiI	rCutSmart	50	100	100	100	65°	No	C	Lambda	CpG	
	TseI	rCutSmart	75	100	100	100	65°	No	B	Lambda	CpG	3
	Tsp45I	rCutSmart	100	50	< 10	100	65°	No	A	Lambda	CpG	
	TspMI	rCutSmart	50*	75*	50*	100	75°	No	B	pUCAdeno	CpG	d
	TspRI	rCutSmart	25	50	25	100	65°	No	B	Lambda	CpG	
	Tth111I	rCutSmart	25	100	25	100	65°	No	B	pBC4	CpG	b
	XbaI	rCutSmart	< 10	100	75	100	37°	65°	A	Lambda HindIII dam-	dam	
	XcmI	r2.1	10	100	25	100*	37°	65°	C	Lambda	CpG	2
	XhoI	rCutSmart	75	100	100	100	37°	65°	A	Lambda HindIII	CpG	b
	XmaI	rCutSmart	25	50	< 10	100	37°	65°	A	pXba	CpG	3
	XmnI	rCutSmart	50	75	< 10	100	37°	65°	A	Lambda	CpG	b
	ZraI	rCutSmart	100	25	10	100	37°	80°	B	Lambda	CpG	

a. Ligation is less than 10%
b. Ligation is 25% – 75%

c. Recutting after ligation is < 5%
d. Recutting after ligation is 50% – 75%

e. Ligation and recutting after ligation is not applicable since the enzyme is either a nicking enzyme, is affected by methylation, or the recognition sequence contains variable sequences.

Activity of Enzymes at 37°C

Listed below is the percentage of activity exhibited at 37°C for enzymes that have an optimal incubation temperature higher (thermophiles) or lower (25°C) than 37°C.

ENZYME	OPTIMAL TEMP. (°C)	% ACTIVITY AT 37°C
AbaSI	25°	0
Apal	25°	100*
ApeKI	75°	10
ApoI	50°	50
BaeI	25°	100
BclI	50°	50
BfuAI	50°	25
BsaBI	60°	25
BsaJI	60°	25
BsaWI	60°	50
BsiEI	60°	10
BsiHKAI	65°	10
BsiWI	55°	25
BsII	55°	50
BsmAI	55°	50
BsmBI-v2	55°	10
BsmFI	65°	100
BsmI	65°	10

ENZYME	OPTIMAL TEMP. (°C)	% ACTIVITY AT 37°C
BspQI	50°	50
BsrI	65°	10
BsrDI	65°	50
BssHII	50°	100
BstAPI	60°	25
BstBI	65°	10
BstEII	60°	10
BstNI	60°	25
BstUI	60°	10
BstYI	60°	10
BtgZI	60°	50
BtsCI	50°	25
BtsIMutI	55°	50
CviAII	25°	10
CviQI	25°	25
FatI	55°	100
FauI	55°	50
MwoI	60°	25

ENZYME	OPTIMAL TEMP. (°C)	% ACTIVITY AT 37°C
Nb.BsmI	65°	100
Nb.BsrDI	65°	50
Nt.BspQI	50°	50
Nt.BstNBI	55°	50
PI-PspI	65°	10
PspGI	75°	25
SfiI	50°	10
SmaI	25°	25
SmlI	55°	10
Swal	25°	25
TaqI-v2	65°	10
TfiI	65°	10
TseI	65°	10
Tsp45I	65°	10
TspMI	75°	10
TspRI	65°	10
Tth111I	65°	10

*Apal has 100% activity at 37°C, however the half-life of this enzyme at 37°C is only 30 minutes.

Activity of DNA Modifying Enzymes in rCutSmart Buffer

A selection of DNA modifying enzymes were assayed in rCutSmart or CutSmart Buffer, in lieu of their supplied buffers. Functional activity was compared to the activity in its supplied buffer, plus required supplements. Reactions were set up according to the recommended reaction conditions, with rCutSmart or CutSmart Buffer replacing the supplied buffer.

ENZYME	ACTIVITY IN rCUTSMART OR CUTSMART	REQUIRED SUPPLEMENTS
Antarctic Phosphatase	+++	Requires Zn ²⁺
Bst DNA Polymerase	+++	
CpG Methyltransferase (M. SssI)	+++	
DNA Polymerase I	+++	
DNA Polymerase I, Large (Klenow) Fragment	+++	
DNA Polymerase Klenow Exo ⁻	+++	
DNase I (RNase-free)	+++	Requires Ca ²⁺
<i>E. coli</i> DNA Ligase	+++	Requires NAD
Endonuclease III (Nth), recombinant	+++	
Endonuclease VIII	+++	
Exonuclease I	+++	
Exonuclease III	+++	
Exonuclease VII	+++	
Exonuclease V (Rec BCD)	+++	Requires ATP
GpC Methyltransferase (M. CviPI)	+	Requires DTT
Hi-T4 DNA Ligase	+++	Requires ATP
Lambda Exonuclease	++	
McrBC	+++	
Micrococcal Nuclease	+++	Requires Ca ²⁺

ENZYME	ACTIVITY IN rCUTSMART OR CUTSMART	REQUIRED SUPPLEMENTS
Nuclease Bal-31	+++	
phi29 DNA Polymerase	+++	Requires DTT
Quick CIP	+++	
RecJ ₁	+++	
Salt-T4 DNA Ligase	+++	Requires ATP
Shrimp Alkaline Phosphatase (rSAP)	+++	
T3 DNA Ligase	+++	Requires ATP + PEG
T4 DNA Ligase	+++	Requires ATP
T4 DNA Polymerase	+++	
T4 Phage β-glucosyltransferase (T4-BGT)	+++	
T4 Polynucleotide Kinase	+++	Requires ATP + DTT
T4 PNK (3' phosphatase minus)	+++	Requires ATP + DTT
T5 Exonuclease	+++	
T7 DNA Ligase	+++	Requires ATP + PEG
T7 DNA Polymerase (unmodified)	+++	
T7 Exonuclease	+++	
Thermolabile ExoI	+++	
USER Enzyme, recombinant	+++	

+++ full functional activity
 ++ 50–100% functional activity
 + 0–50% functional activity

Tips for Avoiding Star Activity

Under non-standard reaction conditions, some restriction enzymes are capable of cleaving sequences which are similar but not identical to their defined recognition sequence. This altered specificity has been termed "star activity". Although the propensity for star activity varies, the vast majority of enzymes from New England Biolabs will not exhibit star activity when used under recommended conditions in their supplied NEBuffers. If an enzyme has been reported to exhibit star activity, it will be indicated in the product entry found in the catalog, on the supplied card and on our website.

CONDITIONS THAT CONTRIBUTE TO STAR ACTIVITY	STEPS THAT CAN BE TAKEN TO INHIBIT STAR ACTIVITY
High glycerol concentration (> 5% v/v)	Restriction enzymes are stored in 50% glycerol, therefore the amount of enzyme added should not exceed 10% of the total reaction volume. Use the standard 50 µl reaction volume to reduce evaporation during incubation.
High concentration of enzyme/µg of DNA ratio (varies with each enzyme, usually 100 units/µg)	Use the fewest units possible to achieve digestion. This avoids overdigestion and reduces the final glycerol concentration in the reaction.
Non-optimal buffer	Whenever possible, set up reactions in the recommended buffer. Buffers with differing ionic strengths and pHs may contribute to star activity.
Prolonged reaction time	Use the minimum reaction time required for complete digestion. Prolonged incubation may result in increased star activity, as well as evaporation.
Presence of organic solvents [DMSO, ethanol (4), ethylene glycol, dimethylacetamide, dimethylformamide, sulfolane (5)]	Make sure the reaction is free of any organic solvents, such as alcohols, that might be present in the DNA preparation.
Substitution of Mg ²⁺ with other divalent cations (Mn ²⁺ , Cu ²⁺ , Co ²⁺ , Zn ²⁺)	Use Mg ²⁺ as the divalent cation. Other divalent cations may not fit correctly into the active site of the restriction enzyme, possibly interfering with proper recognition.

Note: The relative significance of each of these altered conditions will vary from enzyme to enzyme.

New England Biolabs recommends setting up restriction enzyme digests in a 50 µl reaction volume. However, different methods may require smaller reaction volumes. When performing restriction enzyme digests in smaller reaction volumes, extra care must be taken to follow the steps listed above to avoid star activity. Alternatively, using our line of **High Fidelity (HF) restriction enzymes** will allow greater flexibility in reaction setup.

Reference:

(1) Nasri, M. and Thomas, D. (1986) *Nucleic Acids Res.* 14, 811.

TOOLS & RESOURCES

Visit NEBRestrictionEnzymes.com to find:

- Online tutorials on how to avoid star activity, and for setting up restriction enzyme digests
- Access to troubleshooting guides & usage guidelines



RESTRICTION ENZYME STAR ACTIVITY

High-Fidelity (HF) Restriction Enzymes

As part of our ongoing commitment to the study and improvement of restriction enzymes, NEB offers a line of High-Fidelity (HF) restriction enzymes. These engineered enzymes have the same specificity as the native enzymes, are all active in rCutSmart or CutSmart Buffer and have reduced star activity. Star activity, or off-target cleavage, is an intrinsic property of restriction enzymes. Most restriction enzymes will not exhibit star activity when used under recommended reaction conditions. However, for enzymes that have reported star activity, extra caution must be taken to set up reactions according to the recommended conditions to avoid unwanted cleavage.

Many techniques such as cloning, genotyping, mutational analysis, mapping, probe preparation, sequencing and methylation detection employ a wide range of reaction conditions and require the use of enzymes under suboptimal conditions. HF enzymes with reduced star activity offer increased flexibility to reaction setup and help maximize results under a wide range of conditions.

In addition to reduced star activity, all of these engineered enzymes work optimally in rCutSmart or CutSmart Buffer, which has the highest level of enzyme compatibility and will simplify double digest reactions. They are all Time-Saver qualified and digest substrate DNA in 5–15 minutes, and can also be incubated overnight without degradation of DNA. HF enzymes are available at the same price as the native enzymes and are supplied with purple loading dye.

Visit www.neb.com/HF to learn more about HF enzymes.

TOOLS & RESOURCES

Visit www.neb.com/HF

- The full list of HF enzymes available
- Online tutorials on how to avoid star activity and setting up digests using the Time-Saver protocol



WHAT IS A HIGH-FIDELITY ENZYME?



Learn how star activity is reduced with HF enzymes.

Isoschizomers

Restriction enzymes that recognize the same sequence are isoschizomers. The first example discovered is called a prototype, and all subsequent enzymes that recognize the same sequence are isoschizomers of the prototype.

NEB provides a list of isoschizomers for commercially-available restriction endonucleases at www.neb.com/isoschizomers. This table also specifies which isoschizomer is available from New England Biolabs.

All recognition sequences are written 5' to 3' using the single letter code nomenclature with the point of cleavage indicated by a "/".

Numbers in parentheses indicate point of cleavage for non-palindromic enzymes. For example, GGTCTC(1/5) indicates cleavage at: 5' ...GGTCTCN/...3' 3' ...CCAGAGNNNNN/...5'

Isoschizomers with alternative cleavage sites (neoschizomers) are indicated with a "A". Enzymes that are not currently commercially available are indicated with a "X".

For more information on isoschizomers, visit REBASE.neb.com

Survival in a Reaction

Restriction enzymes vary with respect to their ability to maintain activity in a reaction over an extended period of time.

+++	Enzyme is active > 8 hours	N/A	Not Available
++	Enzyme is active 4–8 hours		
+	Enzyme is active 2–4 hours		
–	No benefit from digesting over 1 hour		

While most routine restriction digests are incubated for one hour or less at 37°C, there are certain applications that require the addition of less than 1 unit/μg of DNA and increasing the reaction time beyond one hour. The table below can be used as a guide when low levels of enzyme and extended reaction times are needed.

For example, 1 unit of AatII can be used to digest 8 μg of DNA in a 16 hour digest (+++).

Extended activity was determined by performing the restriction endonuclease unit assay, using a 16 hour incubation in place of the standard 1 hour digestion. After the 16 hour digestion, extended activity enzymes (+++) required only 0.13 units to completely digest 1 μg of DNA. Intermediate activity enzymes required either 0.25 (++) or 0.50 (+) units for complete digestion over this extended incubation time. Finally, enzymes marked (–) required 1.0 unit for complete digestion, the same amount of enzyme required for a 1 hour digestion.

Note: Reaction temperature is 37°C, unless otherwise noted.

ENZYME	SURVIVAL
AatII	+++
AbaSI @25°C	N/A
AccI	+++
Acc65I	+
Acil	–
AcII	+
AcuI	–
AfeI	++
AfiIII	+++
AfiIII	+++
AgeI	+
AgeI-HF	++
AhdI	+++
AleI-v2	N/A
AluI	++
AlwI	+
AlwNI	+++
Apal @25°C	+++
ApaLI	+++
ApeKI @75°C	+++
ApoI @50°C	+++
ApoI-HF	+++
AscI	+++
Asel	+++
AsiSI	+++
AvaI	++
Avall	++
AvrII	+++
BaeI @25°C	+
BaeGI	+
BamHI	+
BamHI-HF	+
BanI	+++
BanII	+
BbsI	++
BbsI-HF	–

ENZYME	SURVIVAL
BbvI	–
BbvCI	+++
BccI	+
BceAI	+++
BcgI	+
BciVI	–
BclI @50°C	+
BclI-HF	N/A
BcoDI	+++
Bfal	+
BfuAI @50°C	++
BglI	+++
BglII	++
BlpI	+
BmgBI	–
Bmri	–
BmtI	+
BmtI-HF	+++
Bpml	–
Bpu10I	+
BpuEI	–
BsaI-HFv2	+++
BsaAI	++
BsaBI @60°C	+
BsaHI	+++
BsaJI @60°C	+++
BsaWI @60°C	+++
BsaXI	++
BseRI	+
BseYI	++
BsgI	+
BsiEI @60°C	++
BsiHKAI @65°C	–
BsiWI @55°C	+++
BsiWI-HF	+++
BsII @55°C	+++

ENZYME	SURVIVAL
BsmI @65°C	+
BsmAI @55°C	++
BsmBI-v2 @55°C	+
BsmFI @65°C	+++
BsoBI	+++
Bsp1286I	+
BspCNI @25°C	–
BspDI	++
BspEI	+++
BspHI	+++
BspMI	++
BspQI @50°C	–
BsrI @65°C	++
BsrBI	+
BsrDI @65°C	+
BsrFI-v2	+++
BsrGI	+++
BsrGI-HF	+++
BssHII @50°C	+
BssSI-v2	+++
BstAPI @60°C	++
BstBI @65°C	+++
BstEI @60°C	–
BstEI-HF	–
BstNI @60°C	–
BstUI @60°C	+++
BstXI @55°C	+++
BstYI @60°C	+++
BstZ17I-HF	++
Bsu36I	+++
BtgI	+
BtgZI @60°C	–
BtsI-v2 @55°C	+++
BtsMutI @55°C	+
BtsCI @50°C	+
Cac8I	++

ENZYME	SURVIVAL
ClaI	+
CspCI	+
CviAI @25°C	–
CviKI-1	–
CviQI @25°C	++
DdeI	+++
DpnI	+++
DpnII	+++
DraI	++
DraIII-HF	+++
DrdI	+++
EaeI	+++
EagI	+++
EagI-HF	+++
EarI	+++
EcII	–
Eco53kI	++
EcoNI	+++
EcoO109I	+++
EcoP15I	–
EcoRI	+++
EcoRI-HF	+++
EcoRV	+
EcoRV-HF	+
Esp3I	+++
FatI @55°C	+
FauI @55°C	–
Fnu4HI	++
FokI	–
FseI	–
FspI	+++
FspEI	+++
HaeII	–
HaeIII	++
HgaI	–
HhaI	++

Survival in a Reaction (continued)

ENZYME	SURVIVAL	ENZYME	SURVIVAL	ENZYME	SURVIVAL	ENZYME	SURVIVAL
HincII	+++	MspI	—	PaeR7I	+++	Scal-HF	++
HindIII	+++	MspA1I	++	PaqCI	N/A	ScrFI	++
HindIII-HF	+++	MspJI	+++	PciI	++	SexAI	++
Hinfl	++	MwoI @60°C	+++	PfIFI	+++	SfaNI	+
HinPI	++	NaeI	+	PfIMI	+	Sfcl	—
HpaI	++	NarI	—	PI-PspI @65°C	+++	Sfil @50°C	++
HpaII	++	Nb.BbvCI	+++	PI-SceI	+++	SfoI	—
HphI	++	Nb.BsmI @65°C	++	PleI	+	SgrAI	—
Hpy99I	—	Nb.BsrDI @65°C	++	PluTI	+	SmaI @25°C	+++
Hpy166II	++	Nb.BssSI	+++	PmeI	—	SmlI @55°C	++
Hpy188I	—	Nb.BtsI	++	PmlI	+	SnaBI	+
Hpy188III	++	NciI	+	PpuMI	+++	SpeI	++
HpyAV	—	NcoI	++	PshAI	—	SpeI-HF	+
HpyCH4III	+++	NcoI-HF	++	Psil-v2	+++	SphI	+++
HpyCH4IV	++	NdeI	+++	PspGI @75°C	++	SphI-HF	—
HpyCH4V	++	NgoMIV	+++	PspOMI	+++	SrfI	+++
I-CeuI	++	NheI	++	PspXI	+++	SspI	++
I-SceI	++	NheI-HF	+++	PstI	+	SspI-HF	+
KasI	—	NlaIII	+	PstI-HF	+	StuI	++
KpnI	++	NlaIV	+	PvuI	+++	StyI-HF	++
KpnI-HF	+	NmeAIII	—	PvuI-HF	+++	StyDI	+++
LpnPI	—	NotI	++	PvuII	+++	Swal @25°C	++
MboI	+	NotI-HF	+++	PvuII-HF	—	TaqI-v2 @65°C	+
MbolI	—	NruI	+++	RsaI	++	TfiI @65°C	++
MfeI	++	NruI-HF	+++	RsrII	++	TseI @65°C	+
MfeI-HF	++	NsiI	++	SacI	+++	Tsp45I @65°C	+
MluI	+++	NsiI-HF	N/A	SacI-HF	+++	TspMI @75°C	+++
MluI-HF	+++	NspI	++	SacII	+++	TspRI @65°C	+++
MluCI	—	Nt.A1wI	+++	Sall	+++	Tth111I @65°C	++
MlyI	—	Nt.BbvCI	+++	Sall-HF	++	XbaI	+++
MmeI	—	Nt.BsmAI	+++	SapI	—	XcmI	+++
MnlI	++	Nt.BspQI @50°C	++	Sau3AI	+	XhoI	+++
MscI	+	Nt.BstNBI @55°C	+	Sau96I	++	XmaI	+
MseI	+++	Nt.CviPII	—	SbfI	+	XmnI	++
MsiI	—	Pacl	+++	SbfI-HF	—	ZraI	+

Cleavage of Supercoiled DNA

Restriction enzymes cleave different DNA substrates with varying efficiency. Restriction enzymes were tested for their ability to cleave various plasmids (pBR322, pUC19 and pLITMUS) under standard reaction conditions. Single sites were tested on each of these plasmids, depending on availability, and average values were taken when there was more than one data point available. Lambda DNA was used as the standard (1 unit to cleave in all cases).

ENZYME	UNITS TO CLEAVE PLASMID	ENZYME	UNITS TO CLEAVE PLASMID	ENZYME	UNITS TO CLEAVE PLASMID	ENZYME	UNITS TO CLEAVE PLASMID
AatII	3	BpmI	1	EcoNI	3	Psil-v2	0.25
AccI	4	BsaAI	20	EcoO109I	8	PstI	1
Acc65I	1	BsgWI	3	EcoRI	3	PvuI	2
AflII	2	BsaWI	3	EcoRV	1	PvuII	2
AflIII	1	BsaXI	2	HincII	4	SacI	5
AgeI	1	BsgI	1	HindIII	5	Sall	10
AhdI	1	BsmI	1	KasI	4	SapI	1
AlwNI	2	BspDI	1	KpnI	2	SmaI	1
Apal	1	BspEI	1	MluI	2	SnaBI	1
ApoI	1	BspMI	**	NarI	10	SpeI	1
Asel	0.3	BspQI	3	NcoI	1	SphI	3
AvaI	10	BsrGI	1	NdeI	3	SspI	4
AvrII	1	BssHI	4	NheI	5	StuI	3
BaeI	3	BstZ17I-HF	1	NmeAIII	**	TspMI	1
BamHI	3	Btgi	5	NruI	1	Tth111I	2
BanII	1	ClaI	5	NsiI	1	XbaI	2
BglII	8	EagI	10	PaqCI	1	XhoI	10
				PciI	3	XmnI	5

** BspMI and NmeAIII require 2 copies of its recognition sequence for cleavage to occur. Thus, the single BspMI site in pBR322 and pUC19 as well as the single NmeAIII site in pUC19 are resistant to cleavage. A 100-fold overdigestion with BspMI cuts less than half the DNA present, while a 10-fold overdigestion with NmeAIII cuts less than half the DNA present. Other plasmids may also be resistant to BspMI and NmeAIII cleavage.

Compatible Cohesive Ends

Restriction enzymes that produce compatible cohesive ends often produce recreatable ligation products. The combinations listed were identified by computer analysis, and have not necessarily been confirmed by experimentation.

Where isoschizomers exist, only one member of each set is listed. A selection of enzymes available from New England Biolabs has been listed. For a more complete listing visit our website, www.neb.com

ENZYME	LIGATED TO	RECLEAVED BY
Acc65I (G/GTACC)	BanI (G/GTACC) BsiWI, BsrGI	Acc65I, BanI, KpnI, NlaIV, RsaI RsaI
AccI (GT/CGAC)	Acil, AcII, BsaHI (GR/CGYC), HinP1I, HpaII, NarI	—
(GT/CGAC)	ClaI, BstBI, TaqI	TaqI
Acil (C/CGC)	AccI (GT/CGAC), AcII, ClaI, BstBI, TaqI BsaHI (GR/CGC), HinP1I, NarI HpaII	Acil HpaII
AcII (AA/CGTT)	AccI (GT/CGAC), Acil, ClaI, BstBI, HinP1I, HpaII, NarI, TaqI	—
AgeI (A/CCGGT)	AvaI (C/CCGGG), XmaI BsaWI, BspEI BsrFI (A/CCGGT), SgrAI (CA/CCGGTG) NgoMIV	HpaII, NciI, ScrFI BsaWI, HpaII AgeI, BsaWI, BsrFI, HpaII BsrFI, HpaII
ApaI (GGGCC/C)	BanII (GGGCC/C), Bsp1286I (GGGCC/C)	ApaI, BanII, Bsp120I, Bsp1286I, HaeIII, NlaIV, Sau96I
ApaLI (G/TGCAC)	SfiI (C/TGCAG)	BsgI
ApoI (A/AATY) (G/AATY) (R/AATY)	EcoRI EcoRI MfeI, MluCI	ApoI, MluCI ApoI, EcoRI, MluCI MluCI
AscI (GG/CGGCC)	MluI BssHII	BstUI, HhaI BssHII, BstUI, Cac8I, HhaI
Asel (AT/TAAT)	BfaI, Csp6I, NdeI MseI	— MseI
AsiSI (GCGAT/CGC)	BsiEI (CGAT/CG) PacI PvuI	DpnII, PvuI MseI DpnII, PvuI
AvaI (C/CCGGG) (C/TCGAG) (C/TCGAG) (C/CCGGG)	AgeI, BsaWI, BspEI, BsrFI (R/CCGGY), NgoMIV, SgrAI (CR/CCGGYG) XhoI Sall XmaI	HpaII, NciI, ScrFI AvaI, TaqI, XhoI TaqI AvaI, BsaJI, HpaII, NciI, ScrFI, SmaI
AvaII (G/GWCC)	PpuMI (RG/GACCY) RsrII PpuMI (RG/GTCCY)	AvaI, NlaIV, Sau96I AvaI, Sau96I AvaI, BsmFI, NlaIV, Sau96I
AvrII (C/CTAGG)	NheI, SpeI, XbaI StyI (C/CTAGG)	BfaI AvrII, BfaI, BsaJI, StyI
BamHI (G/GATCC)	BclI, DpnII BglII, BstYI (R/GATCY) BstYI (G/GATCC)	AlwI, DpnII AlwI, BstYI, DpnII AlwI, BamHI, BstYI, DpnII, NlaIV
BanI (G/GTACC) (G/GCGCC) (G/GTACC)	Acc65I KasI BsiWI, BsrGI	Acc65I, BanI, KpnI, NlaIV, RsaI BanI, BsaHI, HaeII, HhaI, KasI, NarI, NlaIV RsaI
BanII (GGGCC/C) (GAGCT/C)	ApaI, Bsp1286I (GGGCC/C) Bsp1286I (GAGCT/C), SacI	ApaI, BanII, Bsp1286I, HaeIII, NlaIV, Sau96I AluI, BanII, BsiHKAI, Bsp1286I, SacI
BclI (T/GATCA)	BamHI, BstYI (R/GATCY) BglII, MboI	AlwI, DpnII DpnII

Enzymes that have degenerate recognition sequences (e.g., recognize more than one sequence) are followed by a specific sequence in parentheses and are only listed if a non-degenerate equivalent does not exist. Be aware that these degenerate enzymes will cleave sequences in addition to the one listed.

A “—” denotes a ligation product that cannot be recleaved.

ENZYME	LIGATED TO	RECLEAVED BY
BfaI (C/TAG)	Asel, Csp6I, MseI, NdeI	—
BglII (A/GATCT)	BamHI, BstYI (R/GATCY) BclI, DpnII	AlwI, BstYI, DpnII DpnII
BsaHI (GR/CGYC) (GA/CGYC) (GG/CGYC) (GG/CGYC) (GA/CGYC) (GG/CGYC)	AccI (GT/CGAC), ClaI, BstBI, TaqI Acil, HinP1I Acil, HinP1I HpaII NarI NarI	— HgaI HhaI Acil BsaHI, HgaI BanI, BsaHI, HaeII, HhaI, NarI, NlaIV
BsaWI (W/CCGGW)	AgeI, BsrFI (R/CCGGY), SgrAI (CR/ CCGGYG) AvaI (C/CCGGG), XmaI BspEI BsrFI (R/CCGGY), NgoMIV NgoMIV	AgeI, BsaWI, BsrFI, HpaII HpaII, NciI, ScrFI BsaWI, BspEI, HpaII BsrFI, HpaII HpaII
BsiEI (CGAT/CG) (CGAT/CG) (CGGC/CG)	PacI PvuI SacI	MseI BsiEI, DpnII, PvuI Acil
BsiHKAI (GTGCA/C)	Bsp1286I (GTGCA/C) Bsp1286I (GAGCA/C) Bsp1286I (GAGCT/C), SacI NsiI PstI, SbfI	BsiHKAI, Bsp1286I BsiHKAI, Bsp1286I AluI, BanII, BsiHKAI, Bsp1286I, SacI — BsgI
BsiWI (C/GTACG)	Acc65I, BanI (G/GTACC), BsrGI	RsaI
Bsp1286I (GGGCC/C) (GTGCA/C) (GGGCC/C) (GAGCT/C) (GWGCW/C) (GTGCA/C) (GTGCA/C)	ApaI, BanII (GGGCC/C) BsiHKAI BanII (GGGCC/C) BanII (GAGCT/C), BsiHKAI, SacI BsiHKAI NsiI PstI, SbfI	ApaI, BanII, Bsp1286I, HaeIII, NlaIV, Sau96I ApaLI, BsiHKAI, Bsp1286I BanII, Bsp1286I AluI, BanII, BsiHKAI, Bsp1286I, SacI BsiHKAI, Bsp1286I — BsgI
BspEI (T/CCGGA)	AgeI, BsaWI, BsrFI (R/CCGGY), SgrAI (CR/CCGGYG) AvaI (C/CCGGG), XmaI BsaWI BsrFI (R/CCGGY), NgoMIV	BsaWI, HpaII HpaII, NciI, ScrFI BsaWI, BspEI, HpaII HpaII
BspHI (T/CATGA)	FatI, NcoI, PciI	FatI, NlaIII
BsrFI (A/CCGGY) (G/CCGGY) (R/CCGGY) (A/CCGGY) (R/CCGGY) (G/CCGGY) (CR/CCGGYG)	AgeI, BsaWI AgeI, BsaWI, NgoMIV AvaI (C/CCGGG), XmaI BsaWI, BspEI BsaWI, BspEI NgoMIV SgrAI	AgeI, BsaWI, BsrFI, HpaII BsrFI, HpaII HpaII, NciI, ScrFI BsaWI, HpaII HpaII BsrFI, Cac8I, HpaII, NaeI BsrFI, HpaII
BsrGI (T/GTACA)	Acc65I, BanI (G/GTACC), BsiWI	RsaI
BssHII (G/CGCGC)	MluI AscI	BstUI, HhaI BssHII, BstUI, Cac8I, HhaI
BstBI (TT/CGAA)	AccI (GT/CGAC), ClaI, TaqI Acil, AcII, BsaHI (GR/CGYC), HinP1I, HpaII, NarI	TaqI —

ENZYME	LIGATED TO	RECLEAVED BY
BstYI (A/GATCY) (G/GATCY) (R/GATCY) (G/GATCY) (A/GATCY)	BamHI, BglII BamHI BclI, DpnII BclI, DpnII BglII	AlwI, BstYI, DpnII AlwI, BamHI, BstYI, DpnII, NlaIV DpnII AlwI, DpnII BglII, BstYI, DpnII
Clal (AT/CGAT)	AccI (GT/CGAC), BstBI, TaqI Acil, AcII, BsaHI (GR/CGYC), HinP1I, HpaII, NarI	TaqI —
DpnII/MboI/ Sau3AI (GATC)	BamHI, BstYI (R/GATCC) BclI, BglII, BstYI (R/GATCY)	AlwI, DpnII DpnII
EaeI (Y/GGCCR) (C/GGCCR) (T/GGCCR) (C/GGCCR) (T/GGCCR)	PspOMI EagI EagI NotI NotI	HaeIII, Sau96I BsiEI, EaeI, EagI, HaeIII EaeI, HaeIII Acil, BsiEI, EaeI, EagI, Fnu4HI, HaeIII Acil, EaeI, Fnu4HI, HaeIII
EagI (C/GGCCG)	PspOMI EaeI (Y/GGCCR) EaeI (C/GGCCG) NotI	HaeIII, Sau96I EaeI, HaeIII BsiEI, EaeI, EagI, HaeIII Acil, BsiEI, EaeI, EagI, Fnu4HI, HaeIII
EcoRI (G/AATTC)	ApoI (G/AATTC) ApoI (R/AATTY) MfeI, MluCI	ApoI, EcoRI, MluCI ApoI, MluCI MluCI
FatI (/CATG)	BspHI, NcoI, PciI	FatI, NlaIII
HinP1I (G/CGC)	AccI (GT/CGAC), AcII, ClaI, BstBI, TaqI Acil, BsaHI (GR/CGCC), NarI BsaHI (GR/CGTC) HpaII	— HhaI HgaI Acil
HpaII/MspI (C/CGG)	AccI (GT/CGAC), AcII, ClaI, BstBI, TaqI Acil, BsaHI (GR/CGCC), HinP1I, NarI	— Acil
KasI (G/GCGCC)	BanI (G/GCGCC)	BanI, BsaHI, HaeII, HhaI, KasI, NarI, NlaIV
MfeI (C/AATTG)	ApoI (R/AATTY), EcoRI, MluCI	MluCI
MluI (A/CGCGT)	AscI, BssHII	BstUI, HhaI
MluCI (/AATT)	ApoI (R/AATTY), EcoRI, MfeI	MluCI
MseI (T/TAA)	Asel BfaI, Csp6I, NdeI	MseI —
NarI (GG/CGCC)	AccI (GT/CGAC), AcII, ClaI, BstBI, TaqI Acil, HinP1I BsaHI (GR/CGCC) BsaHI (GR/CGTC) HpaII	— HhaI BanI, BsaHI, HaeII, HhaI, NarI, NlaIV BsaHI, HgaI Acil
NcoI (C/CATGG)	BspHI, FatI, PciI	FatI, NlaIII
NdeI (CA/TATG)	Asel, BfaI, Csp6I, MseI	—
NgoMIV (G/CCGGC)	AgeI, BsaWI, BsrFI (R/CCGGY), SgrAI AvaI (C/CCGGG), XmaI BsaWI, BspEI BsrFI (R/CCGGC), SgrAI	BsrFI, HpaII HpaII, NciI, ScrFI HpaII BsrFI, Cac8I, HpaII, NaeI
NheI (G/CTAGC)	AvrII, SpeI, Styl (C/CTAGG), XbaI	BfaI
NlaIII (CATG/)	SphI, NspI	NlaIII
NotI (GC/GGCCCC)	PspOMI EagI EaeI (Y/GGCCR)	Acil, EaeI, Fnu4HI, HaeIII Acil, BsiEI, EaeI, EagI, Fnu4HI, HaeIII Acil, BsiEI, EaeI, Fnu4HI, HaeIII

ENZYME	LIGATED TO	RECLEAVED BY
NsiI (ATGCA/T)	BsiHKAI (GTGCA/C), Bsp1286I (GTGCA/C), PstI, SbfI	—
NspI (RCATG/Y)	NlaIII, SphI	NlaIII, NspI
PacI (TTAAT/TAA)	AsiSI BsiEI (CGAT/CG), PvuI	MseI
PciI (A/CATGT)	BspHI, FatI, NcoI	FatI, NlaIII
PluTI (GGCGC/C)	HaeII	HaeII
PpuMI (RG/GWCCY) (GG/GTCCY) (GG/GACCY)	AvaII, RsrII AvaII, RsrII AvaII, RsrII	AvaII, Sau96I AvaII, BsmFI, NlaIV, Sau96I AvaII, NlaIV, Sau96I
PspOMI (G/GCCCC)	EaeI (Y/GGCCR), EagI NotI	HaeIII, Sau96I Acil, Fnu4HI, HaeIII, Sau96I
PspXI (VC/TCGAGB)	XhoI, TiiI SalI	XhoI, TiiI TaqI
PstI (CTGCA/G)	BsiHKAI, Bsp1286I (GTGCA/C) NsiI SbfI	BsgI — PstI
PvuI (CGAT/CG)	AsiSI PacI BsiEI (CGAT/CG)	DpnII, PvuI MseI BsiEI, DpnII, PvuI
RsrII (CG/GWCCG)	AvaII, PpuMI (RG/GACCY) PpuMI (RG/GACCY) PpuMI (RG/GTCCY)	AvaII, Sau96I AvaII, NlaIV, Sau96I AvaII, BsmFI, NlaIV, Sau96I
SacI (GAGCT/C)	BanII (GAGCT/C), BsiHKAI, Bsp1286I (GAGCT/C)AluI, BanII, BsiHKAI, Bsp1286I, SacI	
SacII (CCGC/GG)	BsiEI (CGGC/CG)	Acil
SalI (G/TCGAC)	PspXI, XhoI	TaqI
SbfI (CCTGCA/GG)	BsiHKAI, Bsp1286I (GTGCA/C) NsiI PstI	BsgI — PstI
SfcI (C/TGCAG)	ApaLI	BsgI
SgrAI (CR/CCGGYG)	See BsrFI	
SpeI (A/CTAGT)	AvrII, NheI, Styl (C/CTAGG), XbaI	BfaI
SphI (GCATG/C)	NlaIII, NspI	NlaIII, NspI
StyI (C/CTAGG) (C/CATGG)	AvrII NheI, SpeI, XbaI BspHI NcoI	AvrII, BfaI, BsaJI, Styl BfaI NlaIII BsaJI, NcoI, NlaIII, Styl
TaqI (T/CGA)	AccI (GT/CGAC), ClaI, BstBI Acil, AcII, BsaHI (GR/CGYC), HinP1I, HpaII, NarI	TaqI —
XbaI (T/CTAGA)	AvrII, NheI, SpeI, Styl (C/CTAGG)	BfaI
XhoI (TiiI) (C/TCGAG)	PspXI SalI	XhoI, TiiI TaqI
XmaI (C/CCGGG)	AgeI, BsaWI, BspEI, BsrFI, NgoMIV, SgrAI AvaI (C/CCGGG)	HpaII, NciI, ScrFI AvaI, BsaJI, HpaII, NciI, ScrFI, SmaI, XmaI

Dam (G^mATC), Dcm (C^mCWGG) and CpG (C^mCG) Methylation

DNA methyltransferases (MTases) that transfer a methyl group from S-adenosylmethionine to either adenine or cytosine residues are found in a wide variety of prokaryotes and eukaryotes. Methylation should be considered when digesting DNA with restriction endonucleases because cleavage can be blocked or impaired when a particular base in the recognition site is methylated.

Prokaryotic Methylation

In prokaryotes, MTases have most often been identified as elements of restriction/modification systems that act to protect host DNA from cleavage by the corresponding restriction endonuclease. Most laboratory strains of *E. coli* contain three site-specific DNA methyltransferases.

- Dam methyltransferases—methylation at the N⁶ position of the adenine in the sequence GATC (1,2).
- Dcm methyltransferases—methylation at the C⁵ position of cytosine in the sequences CCAGG and CCTGG (1,3).
- EcoKI methylase—methylation of adenine in the sequences AAC(N⁶A)GTGC and GCAC(N⁶A)GTT.

Some or all of the sites for a restriction endonuclease may be resistant to cleavage when isolated from strains expressing the Dam or Dcm MTase if the methylase recognition site overlaps the endonuclease recognition site. For example, plasmid DNA isolated from *dam⁻* *E. coli* is completely resistant to cleavage by MboI, which cleaves at GATC sites.

Not all DNA isolated from *E. coli* is methylated to the same extent. While pBR322 DNA is fully modified (and is therefore completely resistant to MboI digestion), only about 50% of λ DNA Dam sites are methylated, presumably because the methylase does not have the opportunity to methylate the DNA fully before it is packaged into the phage head. As a result, enzymes blocked by Dam or Dcm modification will yield partial digestion patterns with λ DNA.

Restriction sites that are blocked by Dam or Dcm methylation can be un-methylated by cloning your DNA into a *dam⁻*, *dcm⁻* strain of *E. coli*, such as *dam⁻/dcm⁻* Competent *E. coli* (NEB #C2925).

Restriction sites can also be blocked if an overlapping site is present. In this case, part of the Dam or Dcm sequence is generated by the restriction enzyme sequence, followed by the flanking sequence. This situation should also be considered when designing restriction enzyme digests.

Eukaryotic Methylation

CpG MTases, found in higher eukaryotes (e.g., Dnmt1), transfer a methyl group to the C⁵ position of cytosine residues. Patterns of CpG methylation are heritable, tissue specific and correlate with gene expression. Consequently, CpG methylation has been postulated to play a role in differentiation and gene expression (4).

Note: The effects of CpG methylation are mainly a concern when digesting eukaryotic genomic DNA. CpG methylation patterns are not retained once the DNA is cloned into a bacterial host.

Methylation Sensitivity

Information on methylation sensitivity for NEB restriction enzymes can be found in the Restriction Enzymes Performance chart on page 309–314, as well as at REBASE.neb.com.

References

- (1) Marinus, M.G. and Morris, N.R. (1973) *J. Bacteriol.*, 114, 1143–1150.
- (2) Geier, G.E. and Modrich, P. (1979) *J. Biol. Chem.*, 254, 1408–1413.
- (3) May, M.S. and Hattman, S. (1975) *J. Bacteriol.*, 123, 768–770.
- (4) Siegfried, Z. and Cedar, H. (1997) *Curr. Biol.*, 7, r305–307.



Kelly works as a Production Technician at our Rowley Production Facility, helping with the formulation of enzymes for catalog and custom products. In her spare time, Kelly enjoys surfing, playing soccer and hiking.

General Guidelines for PCR Optimization

New England Biolabs offers a diverse group of DNA Polymerases for PCR-based applications. Specific recommendations for PCR optimization can be found in the product literature or on the individual product webpages. However, these general guidelines will help to ensure success using New England Biolabs' PCR enzymes.

SETUP GUIDELINES

DNA Template

- Use high quality, purified DNA templates whenever possible. Please refer to specific product information for amplification from unpurified DNA (e.g., colony PCR or direct PCR).
- For low complexity templates (e.g., plasmid, lambda, BAC DNA), use 1 pg–10 ng of DNA per 50 µl reaction
- For higher complexity templates (e.g., genomic DNA), use 1 ng–1 µg of DNA per 50 µl reaction
- Higher DNA concentrations tend to decrease amplicon specificity, particularly for high numbers of cycles

Primers

- Primers should typically be 20–30 nucleotides in length, with 40–60% GC Content
- Primer T_m values should be determined with NEB's T_m Calculator (TmCalculator.neb.com)
- Primer pairs should have T_m values that are within 5°C
- Avoid secondary structure (e.g., hairpins) within each primer and potential dimerization between the primers
- Higher than recommended primer concentrations may decrease specificity
- When engineering restriction sites onto the end of primers, 6 nucleotides should be added 5' to the site
- Annealing temperatures should be determined according to specific enzyme recommendations. *Please note that Q5® and Phusion®** annealing temperature recommendations are unique.*
- Final concentration of each primer should be 0.05–1 µM in the reaction. Please refer to the more detailed recommendations for each specific enzyme.
- When amplifying products > 20 kb in size, primers should be ≥ 24 nucleotides in length with a GC content above 50% and matched T_m values above 60°C
- To help eliminate primer degradation and subsequent non-specific product formation, use a hot-start enzyme (e.g., OneTaq® Hot Start DNA Polymerase or Q5 Hot Start High-Fidelity DNA Polymerase)

Magnesium Concentration

- Optimal Mg²⁺ concentration is usually 1.5–2.0 mM for most PCR polymerases
- Most PCR buffers provided by NEB already contain sufficient levels of Mg²⁺ at 1X concentrations.
- NEB offers a variety of Mg-free reaction buffers to which supplemental Mg²⁺ can be added for applications that require complete control over Mg²⁺ concentration
- Further optimization of Mg²⁺ concentration can be done in 0.2–1 mM increments, if necessary. For some specific applications, the enzyme may require as much as 6 mM Mg²⁺ in the reaction.
- Excess Mg²⁺ may lead to spurious amplification; Insufficient Mg²⁺ concentrations may cause reaction failure

Deoxynucleotides

- Ideal dNTP concentration is typically 200 µM of each, however, some enzymes may require as much as 400 µM each. Please refer to specific product literature for more detailed recommendations.
- Excess dNTPs can chelate Mg²⁺ and inhibit the polymerase
- Lower dNTP concentration can increase fidelity, however, yield is often reduced
- The presence of uracil in the primer, template, or deoxynucleotide mix will cause reaction failure when using archaeal PCR polymerases. Use OneTaq or Taq DNA Polymerases for these applications.

Enzyme Concentration

- Optimal enzyme concentration in the reaction is specific to each polymerase. Please see the product literature for specific recommendations.
- In general, excess enzyme can lead to amplification failure, particularly when amplifying longer fragments

Starting Reactions

- Unless using a hot start enzyme (e.g., OneTaq Hot Start DNA Polymerase or Q5 Hot Start High-Fidelity DNA Polymerase), assemble all reaction components on ice
- Add the polymerase last, whenever possible
- Transfer reactions to a thermocycler that has been pre-heated to the denaturation temperature. Please note that pre-heating the thermocycler is not necessary when using a hot start enzyme (e.g., OneTaq Hot Start DNA Polymerase or Q5 Hot Start High-Fidelity DNA Polymerase).

CYCLING GUIDELINES

Denaturation

- Optimal denaturation temperature ranges from 94°–98°C and is specific to the polymerase in the reaction. Please refer to product information for recommended conditions.
- Avoid longer or higher temperature incubations unless required due to high GC content of the template
- For most PCR polymerases, denaturation of 5–30 seconds is recommended during cycling
- NEB's aptamer-based hot start enzymes do not require additional denaturation steps to activate the enzymes

Annealing

- Primer T_m values should be determined using the NEB T_m Calculator (TmCalculator.neb.com)
- For PCR polymerases other than Q5 High-Fidelity DNA Polymerase or Phusion High-Fidelity DNA Polymerase*, annealing temperatures are usually set at 2°–5°C below the lowest T_m of the primer pair
- When using Q5 High-Fidelity DNA Polymerase or Phusion High-Fidelity DNA Polymerase*, annealing temperatures should be set at 0°–3°C above the lowest T_m of the primer pair. Please refer to the product literature for detailed recommendations.
- Non-specific product formation can often be avoided by optimizing the annealing temperature or by switching to a hot start enzyme (e.g., OneTaq Hot Start DNA Polymerase or Q5 Hot Start High-Fidelity DNA Polymerase)
- Annealing temperatures can be optimized by doing a temperature gradient PCR, starting at 5°C below the lowest T_m of the primer pair
- Ideally, primer T_m values should be less than the extension temperature. However, if T_m values are calculated to be greater than the extension temperature, a two-step PCR program (combining annealing and extension into one step) can be employed.

Extension

- Extension temperature recommendations range from 65°–72°C and are specific to each PCR polymerase. Please refer to the product literature for specific recommendations.
- Extension rates are specific to each PCR polymerase. In general, extension rates range from 15–60 seconds per kb. Please refer to the recommendations for each specific product.
- Longer than recommended extension times can result in higher error rates, spurious banding patterns and/or reduction of amplicon yields

* Phusion DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific. Phusion® is a registered trademark and property of Thermo Fisher Scientific.

PCR Troubleshooting Guide

The following guide can be used to troubleshoot PCR reactions. Additional tips for optimizing reactions can be found in the technical reference section of our website, www.neb.com.

PROBLEM	POSSIBLE CAUSE	SOLUTION
Sequence errors	Low fidelity polymerase	• Choose a higher fidelity polymerase such as Q5 High-Fidelity (NEB #M0491) or Phusion (NEB #M0530)* DNA Polymerases
	Suboptimal reaction conditions	• Reduce number of cycles • Decrease extension time
	Unbalanced nucleotide concentrations	• Prepare fresh deoxynucleotide mixes
	Template DNA has been damaged	• Start with a fresh template • Try repairing DNA template with the PreCR® Repair Mix (NEB #M0309) • Limit UV exposure time when analyzing or excising PCR product from the gel
	Desired sequence may be toxic to host	• Clone into a non-expression vector • Use a low-copy number cloning vector
Incorrect product size	Incorrect annealing temperature	• Recalculate primer T _m values using the NEB T _m calculator (TmCalculator.neb.com)
	Mispriming	• Verify that primers have no additional complementary regions within the template DNA
	Improper Mg ²⁺ concentration	• Adjust Mg ²⁺ concentration in 0.2–1 mM increments
	Nuclease contamination	• Repeat reactions using fresh solutions
No product	Incorrect annealing temperature	• Recalculate primer T _m values using the NEB T _m calculator (TmCalculator.neb.com) • Test an annealing temperature gradient, starting at 5°C below the lower T _m of the primer pair
	Poor primer design	• Check specific product literature for recommended primer design • Verify that primers are non-complementary, both internally and to each other • Increase length of primer
	Poor primer specificity	• Verify that oligos are complementary to proper target sequence
	Insufficient primer concentration	• Primer concentration can range from 0.05–1 µM in the reaction. Please see specific product literature for ideal conditions
	Missing reaction component	• Repeat reaction setup
	Suboptimal reaction conditions	• Optimize Mg ²⁺ concentration by testing 0.2–1 mM increments • Thoroughly mix Mg ²⁺ solution and buffer prior to adding to the reaction • Optimize annealing temperature by testing an annealing temperature gradient, starting at 5°C below the lower T _m of the primer pair
	Poor template quality	• Analyze DNA via gel electrophoresis before and after incubation with Mg ²⁺ • Check 260/280 ratio of DNA template
	Presence of inhibitor in reaction	• Further purify starting template by alcohol precipitation, drop dialysis or commercial clean up kit (NEB #T1030) • Decrease sample volume
	Insufficient number of cycles	• Rerun the reaction with more cycles
	Incorrect thermocycler programming	• Check program, verify times and temperatures
	Inconsistent thermocycler block temperature	• Test calibration of heating block
	Contamination of reaction tubes or solutions	• Autoclave empty reaction tubes prior to use to eliminate biological inhibitors • Prepare fresh solutions or use new reagents
	Complex template	• Use Q5 High-Fidelity (NEB #M0491) or OneTaq DNA Polymerase (NEB #M0482) • For GC-rich templates, use OneTaq DNA Polymerase (NEB #M0480) with OneTaq GC Reaction Buffer (plus OneTaq High GC Enhancer, if necessary) or Q5 High-Fidelity DNA Polymerase (NEB #M0491) with the High GC Enhancer • For longer templates, we recommend LongAmp® Taq DNA Polymerase (NEB #M0323), Q5 or Q5 Hot Start High Fidelity DNA Polymerase (NEB #M0493)
Multiple or non-specific products	Premature replication	• Use a hot start polymerase, such as Q5 Hot Start High-Fidelity (NEB #M0493) or OneTaq Hot Start (NEB #M0481) DNA Polymerases • Set up reactions on ice using chilled components and add samples to thermocycler preheated to the denaturation temperature
	Primer annealing temperature too low	• Recalculate primer T _m values using the NEB T _m Calculator (TmCalculator.neb.com) • Increase annealing temperature
	Incorrect Mg ²⁺ concentration	• Adjust Mg ²⁺ in 0.2–1 mM increments
	Poor primer design	• Check specific product literature for recommended primer design • Verify that primers are non-complementary, both internally and to each other • Increase length of primer • Avoid GC-rich 3' ends
	Excess primer	• Primer concentration can range from 0.05–1 µM in the reaction. Please see specific product literature for ideal conditions.
	Contamination with exogenous DNA	• Use positive displacement pipettes or non-aerosol tips • Set-up dedicated work area and pipettor for reaction setup • Wear gloves during reaction setup
	Incorrect template concentration	• For low complexity templates (e.g., plasmid, lambda, BAC DNA), use 1 pg–10 ng of DNA per 50 µl reaction • For higher complexity templates (e.g., genomic DNA), use 1 ng–1 µg of DNA per 50 µl reaction

* Phusion DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific. Phusion® is a registered trademark and property of Thermo Fisher Scientific.

TIPS FOR OPTIMIZATION

New England Biolabs provides Luna products for your qPCR and RT-qPCR experiments. For more information on these products, visit LUNAQPCR.com. The following tips can be used to help optimize qPCR. For RT-qPCR guidelines, please see page 325.

Target Selection

- Short PCR amplicons, ranging from 70 to 200 bp, are recommended for maximum PCR efficiency
- Target sequences should ideally have a GC content of 40–60%
- Avoid highly repetitive sequences when possible

DNA Template

- Use high quality, purified DNA templates whenever possible. Luna qPCR is compatible with DNA samples prepared through typical nucleic acid purification methods.
- Template dilutions should be freshly prepared in either TE or water for each qPCR experiment
- Generally, useful concentrations of standard and unknown material will be in the range of 10^6 copies to 1 copy. For gDNA samples from large genomes, (e.g., human, mouse) a range of 50–1 pg of gDNA is typical. For small genomes, adjust as necessary using 10^6 –1 copy input as an approximate range. Note that for dilutions in the single-copy range, some samples will contain multiple copies and some will have none, as defined by the Poisson distribution.
- To generate cDNA, use of the LunaScript® RT SuperMix Kit (NEB #E3010) is recommended. Up to 1 µg total RNA, 1 µg mRNA or 100 ng specific RNA can be used in a 20 µl reaction.
- cDNA does not need to be purified before addition to the Luna reaction but should be diluted at least 1:20 before addition to qPCR

Primers

- Primers should typically be 15–30 nucleotides in length
- Ideal primer content is 40–60% GC
- Primer T_m should be approximately 60°C
- Primer T_m calculation should be determined with NEB's TmCalculator (TmCalculator.neb.com) using the Hot Start Taq setting
- For best results in qPCR, primer pairs should have T_m values that are within 3°C
- Avoid secondary structure (e.g., hairpins) within each primer and potential dimerization between primers
- G homopolymer repeats ≥ 4 should be avoided

- Optimal primer concentration for dye-based experiments (250 nM) is lower than for probe-based experiments (400 nM). If necessary, the primer concentration can be optimized between 100–500 nM for dye-based qPCR or 200–900 nM for probe-based experiments.
- Higher primer concentrations may increase secondary priming and create spurious amplification products
- When using primer design software, enter sufficient sequence around the area of interest to permit robust primer design and use search criteria that permit cross-reference against relevant sequence databases to avoid potential off-target amplification.
- For cDNA targets, it is advisable to design primers across known exon-exon junctions in order to prevent amplification from genomic DNA
- Primers designed to target intronic regions can ensure amplification exclusively from genomic DNA

Hydrolysis Probes

- Probes should typically be 15–30 nucleotides in length to ensure sufficient quenching of the fluorophore
- The optimal probe concentration is 200 nM but may be optimized between 100 to 500 nM
- Both single or double-quenched probes may be used
- In general, non-fluorescence quenchers result in better signal-to-noise ratio than fluorescence quenchers
- Ideal probe content is 40–60% GC
- The probe T_m should be 5–10°C higher than the T_m of the primers to ensure all targeted sequences are saturated with probe prior to amplification by the primers
- Probes may be designed to anneal to either the sense or antisense strand
- Generally, probes should be designed to anneal in close proximity to either the forward or reverse primer without overlapping
- Avoid a 5'-G base which is known to quench 5'-fluorophores

Multiplexing

- Avoid primer/probe combinations that contain complementary sequences, and ensure target sequences do not overlap
- Probes should be designed such that each amplicon has a unique fluorophore for detection
- Select fluorophores based on the detection capabilities of the available real-time PCR instrument

- The emission spectra of the reporter fluorophores should not overlap
- Test each primer/probe combination in a singleplex reaction to establish a performance baseline. Ensure C_q values are similar when conducting the multiplex qPCR.
- Pair dim fluorescence dyes with high abundance targets and bright dyes with low abundance targets
- Optimization may require lower primer/probe concentrations to be used for high copy targets along with higher concentrations for low copy targets

Cycling Conditions

- Generally, best performance is achieved using the cycling conditions provided in the manual
- Longer amplicons (> 400 bp) can be used but may require optimization of extension times
- Due to the hot start nature of the polymerase, it is not necessary to preheat the thermocycler prior to use
- Select the "Fast" ramp speed where applicable (e.g., Applied Biosystems QuantStudio®)
- Amplification for 40 cycles is sufficient for most applications, but for very low input samples 45 cycles may be used

Reaction Setup

- For best results, keep reactions on ice prior to thermocycling
- A reaction volume of 20 µl is recommended for 96-well plates while a reaction volume of 10 µl is recommended for 384-well plates
- Reactions should be carried out in triplicate for each sample
- For each amplicon, ensure to include no template controls (NTC)
- To prevent carry-over contamination, treat reactions with 0.2 units/µl Antarctic Thermolabile UDG (NEB #M0372) for 10 minutes at room temperature prior to thermocycling
- The Luna reference dye supports broad instrument compatibility (High-ROX, Low-ROX, ROX-independent) so no additional ROX is required for normalization

Assay Performance

- Ensure 90–110% PCR efficiency for the assay over at least three log₁₀ dilutions of template
- Linearity over the dynamic range (R²) should ideally be ≥ 0.99
- Target specificity should be confirmed by product size, sequencing or melt-curve analysis

Optimization Tips for Luna One-Step RT-qPCR

TIPS FOR OPTIMIZATION

New England Biolabs provides Luna products for your qPCR and RT-qPCR experiments. For more information on these products, visit [LUNAqPCR.com](https://www.neb.com/LUNAqPCR.com). The following tips can be used to help optimize your one-step RT-qPCR. For qPCR guidelines (DNA/cDNA starting material), please see page 324.

Target Selection

- Short PCR amplicons, ranging from 70 to 200 bp, are recommended for maximum PCR efficiency
- Target sequences should ideally have a GC content of 40–60%
- Avoid highly repetitive sequences when possible
- Target sequences containing significant secondary structure should be avoided

RNA Template

- Use high quality, purified RNA templates whenever possible. Luna qPCR is compatible with RNA samples prepared through typical nucleic acid purification methods.
- Prepared RNA should be stored in an EDTA-containing buffer (e.g., 1X TE) for long-term stability
- Template dilutions should be freshly prepared in either TE or water for each qPCR experiment
- Treatment of RNA samples with DNase I (NEB #M0303) may minimize amplification from genomic DNA contamination
- Generally, useful concentrations of standard and unknown material will be in the range of 10^6 copies to 10 copies. Note that for dilutions in the single-copy range, some samples will contain multiple copies and some will have none, as defined by the Poisson distribution. For total RNA, Luna One-Step Kits can provide linear quantitation over an 8-order input range of 1 µg–0.1 pg. For most targets, a standard input range of 100 ng–10 pg total RNA is recommended. For purified mRNA, input of ≤ 100 ng is recommended. For *in vitro*-transcribed RNA, input of ≤ 10^9 copies is recommended.

Primers

- Primers should typically be 15–30 nucleotides in length
- Ideal primer content is 40–60% GC
- Primer T_m should be approximately 60°C
- Primer T_m calculation should be determined with NEB's TmCalculator. ([TmCalculator.neb.com](https://www.neb.com/TmCalculator)) using the Hot Start Tag setting.
- For best results in qPCR, primer pairs should have T_m values that are within 3°C
- Avoid secondary structure (e.g., hairpins) within each primer and potential dimerization between primers
- G homopolymer repeats ≥ 4 should be avoided

- The optimal primer concentration for dye-based experiments and probe-based experiments is 400 nM. If necessary, the primer concentration can be optimized between 100–900 nM.
- Higher primer concentrations may increase secondary priming and create spurious amplification products
- When using primer design software, enter sufficient sequence around the area of interest to permit robust primer design and use search criteria that permit cross-reference against relevant sequence databases to avoid potential off-target amplification
- It is advisable to design primers across known exon-exon junctions in order to prevent amplification from genomic DNA

Hydrolysis Probes

- Probes should typically be 15–30 nucleotides in length to ensure sufficient quenching of the fluorophore
- The optimal probe concentration is 200 nM but may be optimized between 100 to 500 nM
- Both single or double-quenched probes may be used
- In general, non-fluorescence quenchers result in better signal-to-noise ratio than fluorescence quenchers
- Ideal probe content is 40–60% GC
- The probe T_m should be 5–10°C higher than the T_m of the primers to ensure all targeted sequences are saturated with probe prior to amplification by the primers
- Probes may be designed to anneal to either the sense or antisense strand
- Generally, probes should be designed to anneal in close proximity to either the forward or reverse primer without overlapping
- Avoid a 5'-G base which is known to quench 5'-fluorophores

Multiplexing

- Avoid primer/probe combinations that contain complementary sequences, and ensure target sequences do not overlap
- Probes should be designed such that each amplicon has a unique fluorophore for detection
- Select fluorophores based on the detection capabilities of the available real-time PCR instrument
- The emission spectra of the reporter fluorophores should not overlap
- Test each primer/probe combination in a singleplex reaction to establish a performance baseline. Ensure C_q values are similar when conducting the multiplex qPCR.
- Pair dim fluorescence dyes with high abundance targets and bright dyes with low abundance targets

- Optimization may require lower primer/probe concentrations to be used for high copy targets along with higher concentrations for low copy targets

Reverse Transcription

- The default reverse transcription temperature is 55°C
- For difficult targets, the temperature of reverse transcription may be increased to 60°C for 10 minutes
- Due to the WarmStart feature of the Luna RT, reverse transcription temperatures lower than 50°C are not recommended

Cycling Conditions

- Generally, best performance is achieved using the cycling conditions provided in the manual
- Longer amplicons (> 400 bp) can be used but may require optimization of extension times
- Due to the dual WarmStart/Hot Start feature of the Luna kits, it is not necessary to preheat the thermocycler prior to use
- Select the "Fast" ramp speed where applicable (e.g., Applied Biosystems QuantStudio).
- Amplification for 40 cycles is sufficient for most applications, but for very low input samples 45 cycles may be used

Reaction Setup

- For best results, keep reactions on ice prior to thermocycling
- A reaction volume of 20 µl is recommended for 96-well plates while a reaction volume of 10 µl is recommended for 384-well plates
- Reactions should be carried out in triplicate for each sample
- For each amplicon, ensure to include no template controls (NTC)
- A no Luna RT control should be conducted to guarantee amplification is specific for RNA input and not due to genomic DNA contamination
- To prevent carry-over contamination, treat reactions with 0.2 units/µl Antarctic Thermolabile UDG (NEB #M0372) for 10 minutes at room temperature prior to thermocycling
- The Luna reference dye supports broad instrument compatibility (High-ROX, Low-ROX, ROX-independent) so no additional ROX is required for normalization

Assay Performance

- Ensure 90–110% PCR efficiency for the assay over at least three log₁₀ dilutions of template.
- Linearity over the dynamic range (R²) should ideally be ≥ 0.99
- Target specificity should be confirmed by product size, sequencing or melt-curve analysis

Luna qPCR Troubleshooting Guide

PROBLEM	PROBABLE CAUSE(S)	SOLUTION(S)
qPCR traces show low or no amplification	Reagent omitted from qPCR assay	• Verify all steps of the protocol were followed correctly
	Reagent added improperly to qPCR assay	
	Incorrect cycling protocol	• Refer to the proper qPCR cycling protocol in product manual
	Incorrect channel selected for the qPCR thermal cycler	• Verify correct optical settings on the qPCR instrument
	DNA template or reagents are contaminated or degraded	• Confirm the expiration dates of the kit reagents • Verify proper storage conditions provided in this user manual • Rerun the qPCR assay with fresh reagents • Confirm template input amount
Inconsistent qPCR traces for triplicate data	Improper pipetting during qPCR assay set-up	• Ensure proper pipetting techniques
	qPCR plate film has lost its seal, causing evaporation in the well. The resulting qPCR trace may show significantly different fluorescence values relative to its replicates	• Ensure the qPCR plate is properly sealed before inserting into the qPCR thermal cycler. • Exclude problematic trace(s) from data analysis.
	Poor mixing of reagents during qPCR set-up	• Make sure all reagents are properly mixed after thawing them
	Bubbles cause an abnormal qPCR trace	• Avoid bubbles in the qPCR plate • Centrifuge the qPCR plate prior to running it in the thermal cycler • Exclude problematic trace(s) from data analysis
DNA standard curve has a poor correlation coefficient/efficiency of the DNA standard curve falls outside the 90–110% range	Presence of outlying qPCR traces	• Omit data produced by qPCR traces that are clearly outliers caused by bubbles, plate sealing issues, or other experimental problems
	Improper pipetting during qPCR assay set-up	• Ensure that proper pipetting techniques are used
	Reaction conditions are incorrect	• Verify that all steps of the protocol were followed correctly
	Bubbles cause an abnormal qPCR trace	• Avoid bubbles in the qPCR plate • Centrifuge the qPCR plate prior to running it in the thermal cycler
	Poor mixing of reagents	• After thawing, make sure all reagents are properly mixed
	Threshold is improperly set for the qPCR traces	• Ensure the threshold is set in the exponential region of qPCR traces • Refer to the real-time instrument user manual to manually set an appropriate threshold
Melt curve shows different peaks for low input samples	Non-template amplification is occurring Infrequently, denaturation of a single species can occur in a biphasic manner, resulting in two peaks	• Compare melt curve of NTC to samples • Redesign primers with a T _m of 60°C or use our T _m calculator to determine the optimal annealing temperature of the primers • Perform a primer matrix analysis to determine optimal primer concentrations
No template control qPCR trace shows amplification, NTC C_q is close to or overlapping lower copy standards	Reagents are contaminated with carried-over products of previous qPCR (melt curve of NTC matches melt curve of higher input standards)	• Replace all stocks and reagents • Clean equipment and setup area with a 10% chlorine bleach • Consider use of 0.2 U/μl Antarctic Thermolabile UDG to eliminate carryover products
	Primers produce non-specific amplification (melt curve of NTC does not match melt curve of higher input standards)	• Redesign primers with a T _m of 60°C or use qPCR primer design software

Brian is a member of our Digital Marketing Team. He is responsible for uploading content to the many pages found on our website. When not working in his home office, Brian can be found tending to his farm and his pet pig.



Luna One-Step RT-qPCR Troubleshooting Guide

PROBLEM	PROBABLE CAUSE(S)	SOLUTION(S)
qPCR traces show low or no amplification	Incorrect RT step temperature or RT step omitted	• For typical use, a 55°C RT step temperature is optimal for the Luna WarmStart Reverse Transcriptase.
	Incorrect cycling protocol	• Refer to the proper RT-qPCR cycling protocol in product manual
	Reagent omitted from RT-qPCR assay	• Verify all steps of the protocol were followed correctly
	Reagent added improperly to RT-qPCR assay	• Verify correct optical settings on the qPCR instrument
	Incorrect channel selected for the qPCR thermal cycler	• Prepare high quality RNA without RNase/DNase contamination • Confirm template input amount • Confirm the expiration dates of the kit reagents • Verify proper storage conditions provided in product manual • Rerun the RT-qPCR assay with fresh reagents
Inconsistent qPCR traces for triplicate data	Improper pipetting during RT-qPCR assay set-up	• Ensure proper pipetting techniques
	qPCR plate film has lost its seal, causing evaporation in the well. The resulting qPCR trace may show significantly different fluorescence values relative to its replicates.	• Ensure the qPCR plate is properly sealed before inserting into the qPCR thermal cycler • Exclude problematic trace(s) from data analysis
	Poor mixing of reagents during RT-qPCR set-up	• Make sure all reagents are properly mixed after thawing them
	Bubbles cause an abnormal qPCR trace	• Avoid bubbles in the qPCR plate • Centrifuge the qPCR plate prior to running it in the thermal cycler • Exclude problematic trace(s) from data analysis
Standard curve has a poor correlation coefficient/efficiency of the standard curve falls outside the 90–110% range	Cycling protocol is incorrect	• Refer to the proper RT-qPCR cycling protocol in product manual • Use a 55°C RT step temperature • For ABI instruments, use a 1 minute 60°C annealing/extension step
	Presence of outlying qPCR traces	• Omit data produced by qPCR traces that are clearly outliers caused by bubbles, plate sealing issues, or other experimental problems
	Improper pipetting during RT-qPCR assay set-up	• Ensure that proper pipetting techniques are used
	Reaction conditions are incorrect	• Verify that all steps of the protocol were followed correctly
	Bubbles cause an abnormal qPCR trace	• Avoid bubbles in the qPCR plate • Centrifuge the qPCR plate prior to running it in the thermal cycler
	Poor mixing of reagents	• After thawing, make sure all reagents are properly mixed
	Threshold is improperly set for the qPCR traces	• Ensure the threshold is set in the exponential region of qPCR traces • Refer to the real-time instrument user manual to manually set an appropriate threshold
Melt curve shows different peaks for low input samples	Non-template amplification is occurring Infrequently, denaturation of a single species can occur in a biphasic manner, resulting in two peaks	• Compare melt curve of NTC to samples • Redesign primers with a T _m of 60°C or use our T _m calculator to determine the optimal annealing temperature of the primers • Perform a primer matrix analysis to determine optimal primer concentrations
No template control qPCR trace shows amplification/NTC C _q is close to or overlapping lower copy standards	Reagents are contaminated with carried-over products of previous qPCR (Melt curve of NTC matches melt curve of higher input standards)	• Replace all stocks and reagents • Clean equipment and setup area with a 10% chlorine bleach • Consider use of 0.2 U/μl Antarctic Thermolabile UDG to eliminate carryover products
	Primers produce non-specific amplification (Melt curve of NTC does not match melt curve of higher input standards)	• Redesign primers with a T _m of 60°C or use qPCR primer design software
Amplification in No-RT control	RNA is contaminated with genomic DNA	• Treat sample with DNase I • Redesign primer to span exon-exon junction

Cleavage Close to the End of DNA Fragments

Annealed 5' FAM-labeled oligos were incubated with the indicated enzyme (10 units/ 1pmol oligo) for 60 minutes at the recommended incubation temperature and NEBuffer. The digest was run on a TBE acrylamide gel and analyzed by fluorescent imaging. The double stranded oligos were designed to have the indicated number of base pairs from the end followed by the recognition sequence and an additional 12 bases. In some cases asymmetric cleavage was observed and interpreted as a negative result. Asymmetric cleavage decreased with increasing base pairs from the end.

Note: As a general rule and for enzymes not listed below, 6 base pairs should be added on either side of the recognition site to cleave efficiently.

The extra bases should be chosen so that palindromes and primer dimers are not formed. In most cases there is no requirement for specific bases.

Chart Legend

— 0% + 0–20%
++ 20–50% +++ 50–100%
NT not tested

ENZYME	BASE PAIRS FROM END				
	1 bp	2 bp	3 bp	4 bp	5 bp
AccI*	—	—	—	—	—
Acil	—	+	+	++	+++
AgeI	+++	+++	+++	+++	+++
AgeI-HF	++	+++	+++	+++	+++
AlaI-v2	+++	+++	+++	+++	+++
AluI	—	+++	+++	+++	+++
ApaI	+++	+++	+++	+++	+++
AscI	+++	+++	+++	+++	+++
AvrII	++	++	+++	+++	+++
BamHI	+	++	+++	+++	+++
BamHI-HF	+	+	+++	+++	+++
BbsI-HF	+++	+++	+++	+++	+++
BclI-HF	—	—	+++	+++	+++
BglII	++	+++	+++	+++	+++
BmtI	+++	+++	+++	+++	+++
BmtI-HF	+++	+++	+++	+++	+++
BsaI-HFv2	+++	+++	+++	+++	+++
BsiWI	++	+++	+++	+++	+++
BsiWI-HF	++	+++	+++	+++	+++
BsmBI-v2	+++	+++	+++	+++	+++
BsrGI	+++	+++	+++	+++	+++
BssHII	+	+++	+++	+++	+++
BstZ17I-HF	+	+++	+++	+++	+++
ClaI	—	—	+	+++	+++
DdeI	+++	+++	+++	+++	+++
DpnI	—	++	++	NT	NT
DraIII-HF	+++	+++	+++	+++	+++
EagI	++	+++	+++	+++	+++
EagI-HF	+	+++	+++	+++	+++
EcoRI	+	+	++	++	+++
EcoRI-HF	+	+	++	+++	+++
EcoRV	++	++	++	++	+++
EcoRV-HF	+	++	++	++	+++
Esp3I	+++	+++	+++	+++	+++
FseI	+	++	+++	+++	+++
HindIII	—	+	+++	+++	+++
HindIII-HF	—	+	+++	+++	+++
HpaI	+++	+++	+++	+++	+++
KpnI	+	+++	+++	+++	+++
KpnI-HF	+	+++	+++	+++	+++
MfeI	+	++	+++	+++	+++
MfeI-HF	+	++	+++	+++	+++
MluI	+	++	+++	+++	+++

* AccI requires at least 13 base pairs beyond the end of its recognition sequence to cleave efficiently.

ENZYME	BASE PAIRS FROM END				
	1 bp	2 bp	3 bp	4 bp	5 bp
MseI	+++	+++	+++	+++	+++
NcoI	—	++	+++	+++	+++
NcoI-HF	+	++	+++	+++	+++
NdeI	+	+	+++	+++	+++
NheI	+	++	+++	+++	+++
NheI-HF	++	++	+++	+++	+++
NlaIII	++	+++	+++	+++	+++
NotI	++	++	++	++	++
NotI-HF	++	++	++	++	++
NsiI	+	+	+++	+++	+++
NspI	—	—	+	+	+++
PacI	+++	+++	+++	+++	+++
PaqCI	++	+++	—	—	—
PciI	+++	+++	+++	+++	+++
PmeI	+++	+++	+++	+++	+++
PstI	+	+++	+++	+++	+++
PstI-HF	++	+++	+++	+++	+++
PvuI	+++	+++	+++	+++	+++
PvuI-HF	++	+++	+++	+++	+++
PvuII	++	++	++	+++	+++
PvuII-HF	—	++	++	+++	+++
RsaI	+	+++	+++	+++	+++
SacI	—	++	+++	+++	+++
SacI-HF	—	+	+++	+++	+++
SacII	+++	+++	+++	+++	+++
Sall	—	++	+++	+++	+++
Sall-HF	—	++	+++	+++	+++
SapI	+++	+++	+++	+++	+++
Sau3AI	+++	+++	+++	+++	+++
SbfI	++	+++	+++	+++	+++
SbfI-HF	++	+++	+++	+++	+++
Scal-HF	+	+++	+++	+++	+++
SfiI	+++	+++	+++	+++	+++
SmaI	+++	+++	+++	+++	+++
SpeI	+	++	++	++	++
SpeI-HF	+	++	++	++	++
SphI	+++	+++	+++	+++	+++
SphI-HF	++	++	+++	+++	+++
SspI	+	+++	+++	+++	+++
SspI-HF	+	+++	+++	+++	+++
StuI	+++	+++	+++	+++	+++
StyI-HF	+	+++	+++	+++	+++
XbaI	++	++	++	++	++
XhoI	++	++	++	+++	+++
XmaI	+++	+++	+++	+++	+++

Activity of Restriction Enzymes in PCR Buffers

Frequently, a PCR product needs to be digested directly restriction enzymes. For convenience, digestion can be performed directly in the PCR mix without any purification of the DNA. This table summarizes the activity of restriction enzymes on the DNA in *Taq*, Q5, Phusion*, *OneTaq* and LongAmp *Taq* PCR mixes. 50 µl reactions containing 5 units of restriction enzyme were incubated at the appropriate temperature for 1 hour in a PCR mix containing the following: 1 µg DNA, 1 unit of DNA Polymerase and 1X ThermoPol Reaction Buffer, Standard *Taq* Reaction Buffer, Phusion HF Buffer, *OneTaq* Standard Reaction Buffer or LongAmp *Taq* Reaction Buffer. Reactions were supplemented with 200 µM dNTPs. Enzyme activity was analyzed by gel electrophoresis.

Notes: The polymerase is still active and can alter the ends of DNA fragments after cleavage, affecting subsequent ligation. Primers containing the restriction

enzyme recognition site can act as competitive inhibitors in the cleavage reaction. The use of restriction enzymes under non-optimal conditions may increase the likelihood of star activity. If any problems are encountered, the DNA should be purified by spin column or phenol/chloroform extraction followed by alcohol precipitation.

* Phusion DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific. Phusion® is a registered trademark and property of Thermo Fisher Scientific.

Chart Legend

Cleavage in extension mix with 5 units of enzyme:

+++ complete cleavage ++ ~50% cleavage
+ ~25% cleavage – no cleavage

** It has been shown that the addition of 1X Restriction Enzyme Buffer may help to improve the ability of some enzymes to cleave.

ENZYME	<i>Taq</i> IN THERMOPOL RXN BUFFER	Q5 IN Q5 BUFFER**	PHUSION IN PHUSION HF BUFFER	ONE <i>Taq</i> IN ONE <i>Taq</i> RXN BUFFER	LONGAMP <i>Taq</i> IN LONGAMP <i>Taq</i> RXN BUFFER
AatII	<+	<+	+	++	+
AccI	<++	<+	<+	+++	+++
Acc65I	+++	<+	<+	<+	+
Acil	++	++	+++	+++	+++
AcII	+++	<+	<+	+++	+++
AcuI	+++	<+	++	+++	+++
AfeI	+++	<+	++	+++	+++
AflII	+	<+	<+	+	<+
AflIII	<+	+++	+	<+	<+
AgeI	+++	+	+++	+++	<+
AgeI-HF	+++	<+	++	+++	+++
AhdI	<+	–	–	<+	<+
AleI-v2	–	–	–	+	+
AluI	+++	+	+++	+++	+++
AlwI	–	<+	<+	<+	<+
AlwNI	<++	+	+++	<+	+
ApaI	+++	<+	<+	<+	–
ApaLI	+++	<+	<+	+++	+++
ApeKI	<++	++	+++	<+	+
ApoI	+++	++	+++	++	+++
ApoI-HF	+++	+	++	+++	+++
AscI	+++	<+	<+	<+	–
Asel	+++	<+	+	++	++
AsiSI	+++	<+	++	+++	+++
AvaI	+++	<+	+++	+++	+
AvaII	+++	<+	++	+++	+++
AvrII	+++	<+	<+	+++	+++
BaeGI	+++	<+	+++	+++	+++
BaeI	–	<+	++	<+	<+
BamHI	+++	<+	+++	+++	+++
BamHI-HF	+++	<+	–	<+	++
BanI	+++	<+	+++	+++	+++
BanII	+++	<+	+++	+++	+++
BbsI	+++	<+	<+	+++	+++
BbsI-HF	+	–	–	–	+
BbvCI	+++	–	–	<+	<+
BbvI	+++	<+	++	+++	+++
BccI	<+	<+	<+	<+	<+
BceAI	<+	<+	++	+	<+
BcgI	<+	<+	+	++	++
BciVI	–	–	–	<+	–
BclI	+++	++	+++	+++	+++
BclI-HF	+++	–	–	+	+
BcoDI	<+	<+	+	+	<+
BfaI	–	<+	–	–	–
BfuAI	<++	–	+	<+	–
BglI	<+	++	+	<+	<+
BglII	<+	+	++	<+	<+
BipI	<++	<+	<+	<+	–
BmgBI	–	++	+	<+	<+
BmrI	<++	<+	+++	+++	+++
BmtI	+++	<+	++	+++	+++
BmtI-HF	++	<+	+	++	+++
BpmI	<+	<+	+++	<++	<++
BpuEI	+++	–	++	<++	<++
Bpu10I	<+	<+	+++	++	+++
BsaAI	+++	++	+++	+++	+++
BsaBI	+	<+	++	++	+++

ENZYME	<i>Taq</i> IN THERMOPOL RXN BUFFER	Q5 IN Q5 BUFFER**	PHUSION IN PHUSION HF BUFFER	ONE <i>Taq</i> IN ONE <i>Taq</i> RXN BUFFER	LONGAMP <i>Taq</i> IN LONGAMP <i>Taq</i> RXN BUFFER
BsaHI	+++	+	+++	+++	+++
BsaI-HFv2	+	<+	+	+	++
BsaJI	+++	<+	++	+++	+++
BsaWI	<++	<+	++	+	+
BsaXI	<++	<+	<+	<+	<+
BseRI	+++	<+	++	++	+
BseYI	+++	++	++	+++	+++
BsgI	<+	<+	+	<+	<+
BsiEI	+++	<+	++	++	++
BsiHKA1	–	++	+	–	–
BsiWI	+++	<+	+++	+++	+++
BsiWI-HF	–	–	–	–	–
BsII	+++	++	+++	+++	+++
BsmAI	+++	++	+++	<+	<+
BsmBI	<++	+	++	<+	<+
BsmFI	<+	+++	++	+	+
BsmI	+++	+	<+	+++	+
BsoBI	+++	+++	+++	++	+++
BspCNI	<+	<+	+	–	–
BspDI	<++	<+	++	+++	+++
BspEI	–	<+	<+	–	–
BspHI	+++	<+	+++	+++	+++
Bsp1286I	<+	<+	<+	<+	<+
BspMI	+++	<+	++	<+	<+
BspQI	+	++	+++	+++	+++
BsrBI	+++	<+	+	+++	+++
BsrDI	<+	<+	+	<+	<+
BsrFI-v2	<+	–	–	–	–
BsrGI	<+	+	+++	<+	+++
BsrI	+++	<+	+++	++	+++
BssHII	+++	<+	+	+++	+++
BssSI-v2	+++	–	+	+++	+++
BstAPI	+++	<+	++	+++	+++
BstBI	+++	++	+++	+++	+++
BstEI	+++	<+	<+	+++	+++
BstEII-HF	+++	<+	<+	++	++
BstNI	+++	<+	<+	<+	<+
BstUI	+++	<+	<+	+++	+
BstXI	<++	+	+	+	<+
BstYI	+++	<+	<+	++	+
BstZ17I-HF	+++	–	+	+++	+++
Bsu36I	<+	<+	<+	<+	+
BtgI	+++	<+	+	<+	<+
BtgZI	+++	+	++	++	++
BtsI-v2	+++	–	+	+++	+++
BtsCI	+++	<+	<+	+++	+++
CacBI	+++	<+	<+	+++	++
Clal	++	<+	<+	<+	++
CspCI	<+	–	+	<+	<+
CviAI	+++	<+	+	+++	+++
CviKI-1	+++	<+	++	+++	+++
CviQI	+++	+	+++	++	+++
DdeI	+++	++	+	+++	+++
DpnI	+++	++	+++	++	++
DpnII	+++	++	+++	+++	++
DraI	+++	<+	+++	+++	+++
DraIII-HF	++	++	+++	++	++
DrdI	+++	<+	+++	+++	+++

ENZYME	Taq IN THERMOPOL RXN BUFFER	Q5 IN Q5 BUFFER**	PHUSION IN PHUSION HF BUFFER	ONE Taq IN ONE Taq RXN BUFFER	LONGAMP Taq IN LONGAMP Taq RXN BUFFER
EaeI	+++	<+	—	<+	<+
EagI	<+	+++	+++	+++	+++
EagI-HF	+	<+	+	++	++
EarI	+++	<+	+++	+	<+
Ecil	<+	++	+++	<++	<++
Eco53kI	+++	<+	<+	+++	+++
EcoNI	+++	<+	+	+++	+++
EcoO109I	+++	<+	—	<+	+
EcoP15I	<+	<+	+	<+	+
EcoRI	+	<+	+++	—	—
EcoRI-HF	+++	<+	+	+++	+++
EcoRV	<+	<+	+	—	<+
EcoRV-HF	+	<+	<+	+	++
Esp3I	+++	—	+++	+	+++
FatI	++	<+	+++	<+	+++
FauI	+	<+	++	+++	++
Fnu4HI	+++	<+	<+	++	+
FokI	+++	+	+	+++	+++
FseI	+	<+	++	+++	—
FspI	<++	<+	+	+	+
HaeII	+++	<+	+++	+++	+++
HaeIII	+++	<+	+++	+++	+++
HgaI	<+	<+	+	<++	<++
HhaI	+++	<+	+++	+++	+++
HincII	+++	<+	<+	+++	+++
HindIII	+++	<+	+	++	+++
HindIII-HF	+++	<+	<+	+++	+++
HinfI	+++	+++	+++	+	+++
HinPII	+++	+	+++	+++	+++
HpaI	+++	<+	+++	+++	+++
HpaII	+++	<+	<+	<+	<+
HphI	<++	<+	<+	<+	<+
HpyAV	+++	—	++	+	++
HpyCH4III	<++	<+	+	<++	<++
HpyCH4IV	+++	<+	<+	+++	+++
HpyCH4V	+++	<+	<+	+++	+++
Hpy99I	+++	—	+	<+	<+
Hpy188I	+++	<+	+	++	++
Hpy166II	+++	+	++	+++	+++
Hpy188III	+	<+	<+	+	<+
KasI	+++	<+	+++	+++	—
KpnI	+++	++	+	++	<+
KpnI-HF	++	—	++	<+	<+
MboI	+++	<+	+++	+++	+++
MboII	+++	+	++	+	+
MfeI	+++	<+	<+	+++	+
MfeI-HF	+	—	—	+++	<+
MluCI	+	<+	<+	++	+
MluI	+++	++	++	++	++
MluI-HF	++	—	++	++	++
MlyI	+++	+	++	<+	+
MmeI	<+	—	++	<+	<+
MnII	+++	+	+	+	+
MscI	<+	<+	+	<+	<+
MseI	<+	<+	<+	<+	<+
MsiI	+++	<+	+	+++	++
MspA1I	+++	<+	+++	++	+++
MspI	+++	<+	+++	++	+++
MwoI	+++	+++	+++	++	+++
NaeI	<+	<+	+	<+	<+
NarI	—	<+	++	+++	+++
NciI	+++	<+	<+	+	<+
NcoI	+++	<+	+	++	++
NcoI-HF	+++	<+	—	++	+
NdeI	<++	++	+++	++	<+
NgoMIV	—	<+	+	<+	<+
NheI	+++	<+	<+	+++	+++
NheI-HF	+++	<+	—	++	++
NlaIII	<+	<+	+	++	<+
NlaIV	+++	<+	+++	+++	+++
NmeAIII	<+	—	+++	<+	<+
NotI	++	<+	+	<+	<+
NotI-HF	+++	<+	<+	<+	+

ENZYME	Taq IN THERMOPOL RXN BUFFER	Q5 IN Q5 BUFFER**	PHUSION IN PHUSION HF BUFFER	ONE Taq IN ONE Taq RXN BUFFER	LONGAMP Taq IN LONGAMP Taq RXN BUFFER
NruI	++	+	+	++	++
NruI-HF	++	—	—	+	—
NsiI	+++	+	+++	++	+
NsiI-HF	+++	++	+++	+++	+++
NspI	+++	<+	<+	+++	++
PacI	+++	<+	<+	++	+++
PaeR7I	+++	<+	<+	+++	+++
PaqCI	—	—	—	+++	++
PciI	<+	<+	—	—	—
PfiFI	+++	<+	<+	<+	+
PfiIMI	+	<+	+++	++	+++
PleI	+++	<+	<+	<+	<+
PluTI	+++	<+	+	+++	+++
PmeI	+++	<+	<+	+++	+++
PmlI	—	—	—	+	<+
PpuMI	+++	<+	+++	+++	+++
PshAI	+++	<+	<+	<+	<+
PsiI-v2	+++	—	—	+++	+++
PspGI	+++	+++	+++	+++	+++
PspOMI	+++	<+	+	+++	+++
PspXI	+++	<+	++	+++	+++
PstI	++	+	+	<+	<+
PstI-HF	+++	<+	++	++	+
PvuI	<+	<+	+++	—	<+
PvuI-HF	+++	<+	+++	++	+++
PvuII	+++	<+	+	+++	+++
PvuII-HF	+	—	—	<+	<+
RsaI	+++	<+	++	+++	+++
RsrII	<++	—	—	<+	<+
SacI	+++	<+	+	++	++
SacI-HF	+++	<+	<+	<+	++
SacII	+++	<+	+++	++	+
SaiI	<+	+	++	—	—
SaiI-HF	+	<+	+++	+	+++
SapI	<++	<+	++	++	++
Sau3AI	+++	<+	<+	<+	<+
Sau96I	<++	+	+	+++	+++
SbfI	<++	<+	+	<+	+++
SbfI-HF	+	—	—	<+	<+
ScaI-HF	+	<+	<+	—	—
ScrFI	+++	+++	+++	+++	+++
SexAI	+++	<+	+++	+++	+++
SfaNI	—	<+	++	<++	<++
Sfcl	+++	<+	<+	+	+
SfiI	+++	—	—	+++	+++
SfoI	+++	<+	+++	+	+++
SgrAI	<++	<+	++	+	+++
SmaI	+++	<+	++	+++	+++
SmlI	<+	<+	+	+	+
SnaBI	<+	<+	<+	+++	+++
SpeI	+++	+	<+	+++	+++
SpeI-HF	+++	—	<+	+++	+++
SphI	+++	+	++	<+	<+
SphI-HF	+++	<+	+	+++	+++
SrfI	<+	<+	+++	+	++
SspI-HF	++	<+	+	+++	+++
StuI	+++	<+	<+	+++	+++
StyD4I	<++	<+	+	<+	<+
StyI	<+	+	<+	<+	<+
StyI-HF	+	<+	<+	++	+++
Swal	<+	<+	<+	<+	+++
TaqI-v2	+++	<+	+	+++	+++
TfiI	<++	<+	<+	++	++
TseI	+++	+++	+++	+++	+++
Tsp45I	+++	—	—	+	<+
TspMI	+++	<+	+	+++	+++
TspRI	+	<+	<+	+++	+++
Tth111I	+++	<+	++	<+	+
XbaI	+++	—	<+	++	++
XcmI	+++	<+	+	+++	+++
XhoI	<+	<+	+++	++	+++
XmaI	+++	<+	+	—	—
XmnI	+++	<+	<+	++	+++
ZraI	+++	<+	<+	++	+

Getting Started with Molecular Cloning

Molecular cloning has traditionally used restriction enzymes to excise a fragment from source DNA, and to linearize a plasmid vector, while creating compatible ends. After purification, insert and vector are ligated to form a recombinant vector, which is transformed into an *E. coli* host. Alternatively, PCR can be used to generate both the vector and insert, which can be joined using a variety of techniques, such as standard DNA ligation, enzymatic joining using a recombinase or topoisomerase, homologous recombination, or synthetic biology (see NEBuilder HiFi DNA Assembly and Gibson Assembly).

Regardless of the method chosen, the process can be made more efficient and successful by following good practices in the lab. The following tips will help improve the success of your cloning experiments.

1. Take the time to plan your experiments

Pay attention to the junction sequences and the effect on reading frames of any translated sequences. Check both the vector and insert for internal restriction sites (we recommend NEBcutter at NEBcutter.neb.com) prior to designing PCR primers that contain similar sites to those used for cloning. Verify that the antibiotic selective marker in the vector is compatible with the chosen host strain.

2. Start with clean DNA at the right concentration

Ensure that your source DNA is free of contaminants, including nucleases and unwanted enzymatic activities. Use commercially-available spin columns to purify starting DNA, (e.g., Monarch Plasmid Miniprep Kit, NEB #T1010 for DNA plasmids, Monarch PCR & DNA Cleanup Kit, NEB #T1030 for DNA Fragments). Completely remove solvents, such as phenol, chloroform and ethanol, prior to manipulation of the DNA. Elute DNA from the spin columns with salt-free buffer to prevent inhibition of the downstream steps, either restriction digestion or PCR amplification.

3. Perform your restriction digests carefully

The reaction volume should be compatible with the downstream step (e.g., smaller than the volume of the well of an agarose gel used to resolve the fragments). For a typical cloning reaction, this is often between 20–50 μ l. The volume of restriction enzyme(s) added should be no more than 10% of the total reaction volume, to ensure that the glycerol concentration stays below 5%; this is an important consideration to minimize star activity (unwanted cleavage).

4. Mind your ends

DNA ends prepared for cloning by restriction digest are ready for ligation without further modification, assuming the ends to be joined are compatible. If the ends are non-compatible, they can be modified using blunting reagents, phosphatases, etc.

DNA ends prepared by PCR for cloning may have a 3' addition of a single adenine (A) residue following amplification using *Taq* DNA Polymerase (NEB #M0273). High-fidelity DNA polymerases, such as Q5 (NEB #M0491),

leave blunt ends. PCR using standard commercial primers produces non-phosphorylated fragments, unless the primers were 5' phosphorylated. The PCR product may need to be kinase treated to add a 5' phosphate prior to ligation with a dephosphorylated vector.

5. Clean up your DNA prior to vector:insert joining

This can be done with gel electrophoresis or column purification (e.g., Monarch PCR & DNA Cleanup Kit, NEB #T1030). Isolating the desired DNA from unwanted parent vectors and/or other DNA fragments can dramatically improve your cloning results.

Confirm digested DNA on an agarose gel prior to ligation. For a single product, run a small amount of the digest, and then column purify to isolate the remainder (e.g., Monarch PCR & DNA Cleanup Kit, NEB #T1030). When multiple fragments are produced and only one is to be used, resolve the fragments on a gel and excise the desired fragment under UV light. Using longwave (365 nm) UV light will minimize any radiation-induced DNA damage to the fragment. Recover the DNA fragment from the agarose slice using a gel extraction kit (e.g., Monarch DNA Gel Extraction Kit, NEB #T1020) or β -Agarase I (NEB #M0392).

6. Quantitate your isolated material

Simple quantitation methods, such as gel electrophoresis with mass standards or spectroscopic quantitation on low-input spectrophotometers (such as a NanoSpec®), ensure that the proper amount of material is used for the downstream reactions.

7. Follow the manufacturer's guidelines for the joining reaction

For traditional cloning, follow the guidelines specified by the ligase supplier. If a 3:1 molar ratio of insert to vector is recommended, try this first for best results. Using a 3:1 mass ratio is not the same thing (unless the insert and vector have the same mass). Ligation usually proceeds quickly and, unless your cloning project requires the generation of a high-complexity library that benefits from the absolute capture of every possible ligation product, long incubation times are not necessary.

Follow the manufacturer's guidelines for the joining reactions in PCR cloning and seamless cloning. If you are performing a cloning protocol for the first time, adhere to the recommended protocol for optimal results.

8. Use competent cells that are suited to your needs

While some labs prepare their own competent cells "from scratch" for transformations, the levels of competence achieved rarely matches the high levels attained with commercially-available competent cells. Commercially-available competent cells save time and resources, and make cloning more reproducible.

NANOSPEC® is a registered trademark of Nanometrics, Inc.

Find an overview of traditional cloning.



TIPS FOR OPTIMIZATION

Looking to assemble multiple DNA fragments in a single reaction? Here are some tips to keep in mind when planning your Golden Gate Assembly experiment.

Check Your Sequences

- Always check your assembly sequences for internal sites before choosing which Type IIS restriction endonuclease to use for your assembly. While single insert Golden Gate assembly has such high efficiencies of assembly that the desired product is obtainable regardless of the presence of an internal site, this is not true for assemblies with multiple inserts. Options include choosing a different Type IIS restriction enzyme to direct your assembly, or eliminating internal sites through mutagenesis. The Q5® Site-Directed Mutagenesis Kit (NEB #E0554) and the NEB web tool NEBaseChanger work well for this purpose. Alternately, a junction point can be created at the internal site's recognition sequence.

Orient Your Primers

- When designing PCR primers to introduce Type IIS restriction enzyme sites, either for amplicon insert assembly or as an intermediate for pre-cloning the insert, remember that the recognition sites should always face inwards towards your DNA to be assembled. Consult the Golden Gate Assembly Kit manuals for further information regarding the placement and orientation of the sites.

Choose the Right Plasmid

- Consider switching to the pGGAselect Destination Plasmid for your Golden Gate assembly. This versatile new destination construct is included in all Golden Gate Assembly kits and can be used for BsaI-HFv2, BsmBI-v2 or BbsI directed assemblies. It also features T7 and SP6 promoter sequences flanking the assembly site, and has no internal BsaI, BsmBI or BbsI sites. The pGGAselect plasmid can also be transformed into any *E. coli* strain compatible with pUC19 for producing your own plasmid preparation if so desired.

Choose the Right Buffer

- T4 DNA Ligase Buffer works best for Golden Gate Assembly with BsaI-HFv2, BsmBI-v2 and PaqCI. However, alternate buffers would be NEBuffer r1.1 or 1.1 for Bsa-HFv2, NEBuffer r2.1 or 2.1 for BsmBI-v2 & rCutSmart or CutSmart Buffer for PaqCI, if these buffers are supplemented with 1 mM ATP and 5-10 mM DTT.

Increase Your Complex Assembly Efficiency by Increasing the Golden Gate Cycling Levels

- T4 DNA Ligase, BsaI-HFv2, BsmBI-v2 and PaqCI are very stable and continue to be active during extended cycling protocols; an easy way to increase assembly efficiencies without sacrificing fidelity is to increase the total cycles from 30 to 45-65, even when using long (5-minute) segments for the temperature steps.

Make Sure Your Plasmid Prep is RNA-free

- For pre-cloned inserts/modules, make sure your plasmid prep is free of RNA to avoid an overestimation of your plasmid concentrations.

Avoid Primer Dimers

- For amplicon inserts/modules, make sure your PCR amplicon is a specific product and contains no primer dimers. Primer dimers, with Type IIS restriction endonuclease sites (introduced in the primers used for the PCR reactions), would be active in the assembly reaction and result in mis-assemblies.

Avoid PCR-induced Errors

- Do not over-cycle and use a proofreading high fidelity DNA polymerase, such as Q5 DNA High-Fidelity Polymerase.

Decrease Insert Amount for Complex Assemblies

- For complex assemblies involving >10 fragments, pre-cloned insert/modules levels can be decreased from 75 to 50 ng each without significantly decreasing the efficiencies of assembly.

Carefully Design EVERY Insert's Overhang

- An assembly is only as good as its weakest junction. Research at NEB has led to an increased understanding of ligase fidelity, including the development of a comprehensive method for profiling end-joining ligation fidelity in order to predict which overhangs will result in improved accuracy. This ligase fidelity information can be used in conjunction with the NEB Golden Gate Assembly Kits (BsaI-HFv2 or BsmBI-v2) to achieve high efficiencies and accurate complex assemblies. Please use the free NEB Golden Gate Assembly Tool to design primers for your Golden Gate Assembly reactions, predict overhang fidelity or find optimal Golden Gate junctions for long sequences. When working with complex assemblies (> 20 fragments), refer to the ligase fidelity tools on the NEBeta Tools site.

Check for a Sequence Error if Your Assembly Becomes Non-functional

- Be aware that occasionally a pre-cloned insert/ module can become corrupted by an error during propagation in *E. coli*, usually a frameshift due to slippage in a run of a single base (e.g., AAAA) by the *E. coli* DNA Polymerase. This should be suspected if previously functional assembly components suddenly become nonfunctional.

For more information on Golden Gate, visit www.neb.com/GoldenGate



Optimization Tips for NEBuilder HiFi DNA Assembly and NEB Gibson Assembly

TIPS FOR OPTIMIZATION

Looking to assemble multiple DNA fragments in a single reaction? Here are some tips to keep in mind when planning your NEBuilder HiFi DNA Assembly or Gibson Assembly experiments.

Decide How You Want to Generate the Linearized Vector. You can Choose from Two Methods:

- 1) Restriction enzyme digestion: good for large plasmids you don't want to amplify; background may be higher if undigested vector is present.
- 2) PCR: achieves lower background versus restriction enzyme digestion, but is limited by the size of the vector. Typically, vectors up to 10 kb can be amplified; for amplicons greater than 10 kb, divide into 2 fragments.

Design the Primers

- Use the NEBuilder Assembly Tool (NEBuilder.neb.com) to design the primers and check the sequence of the final assembly. Primers will contain the overlap sequence. We recommend watching the tutorials before using the tool for the first time. There is one for restriction enzyme digestion and another for PCR. The videos can be found at NEBuilderHiFi.com
- Make sure the overlap is the correct length for the number of fragments in the assembly: Refer to the section titled "Use the Correct Amount of DNA" for more details.

Column Purify the PCR Products

- If you do not purify the PCR products, limit the unpurified PCR products to 20% of the reaction volume (4 µl for a standard 20 µl reaction).
- If PCR produces a single band of the correct size and the yield is good, DNA purification is not necessary.

- If PCR produces multiple products or a smear, it is best to optimize the PCR. If it is not possible to optimize, purify the products using gel extraction. Be careful, however, as gel extraction can introduce guanidine thiocyanate (from the gel dissolving buffer) and can reduce the efficiency of the assembly reaction. To minimize this contamination, trim the gel slice so that a smaller amount of gel dissolving buffer is required. Due to the potential for residual guanidine salt being present in fragments isolated by gel-extraction, PCR or DNA column purification (NEB #T1030) is preferable to gel extraction (NEB #T1020).

Use the Correct Amount of DNA

- Make sure you calculate the optimum ratio of insert(s):vector. If the ratio is not ideal, we recommend using NEBioCalculator (NEBioCalculator.neb.com) to determine molar amounts.

For NEBuilder HiFi DNA Assembly:

2-3 fragments: 15-20 nt overlaps, total DNA = 0.03-0.2 pmol, 2 fold molar excess of each insert:vector

4-6 fragments: 20-30 nt overlaps, total DNA = 0.2-0.5 pmol, 1:1 molar ratio of each insert:vector

For NEB Gibson Assembly:

2-3 fragments: 15-25 nt overlaps, total DNA = 0.02-0.5 pmol, 2-3 fold molar excess of each insert:vector

4-6 fragments: 20-80 nt overlaps, total DNA = 0.2-1.0 pmol, 1:1 molar ratio of each insert:vector

Perform a PCR Assay to Determine if the Assembly is Successful

- Determine if the assembly works *in vitro* by amplifying the assembled product directly from the assembly reaction. Dilute 1 µl of the assembly reaction with 3 µl water then use 1 µl as a template in a 50 µl PCR. Use primers that anneal to the vector and amplify across the insert. Do not use primers that anneal across the assembly junction because this can lead to false positive results. If you can amplify the assembled product but cannot recover clones by transformation, then the problem is either with the transformation step, or the inability of the cells to maintain the transformed construct due to toxicity.

Check the reaction conditions, DNA amounts, overlap sequences and perform the assembly control.

Always use High Competency Cells with a Transformation Efficiency of 10^8 - 10^9 cfu/µg

- We recommend NEB 5-alpha High Efficiency Competent *E. coli* (NEB #C2987) or NEB 10-beta High Efficiency Competent *E. coli* (NEB #C3019).

More information can be found on NEBuilderHiFi.com



Traditional Cloning Quick Guide

PREPARATION OF INSERT AND VECTORS

Insert From a Plasmid Source

- Digest plasmid with the appropriate restriction enzymes to produce a DNA fragment that can be cloned directly into a vector. Unidirectional cloning is achieved with restriction enzymes that produce non-compatible ends.

Insert From a PCR Product

- Design primers with appropriate restriction sites to clone unidirectionally into a vector
- Addition of 6 bases upstream of the restriction site is sufficient for digestion with most enzymes
- If fidelity is a concern, choose a proofreading polymerase such as Q5 High-Fidelity DNA Polymerase (NEB #M0491)
- Visit www.NEBPCRPolymers.com for additional guidelines for PCR optimization
- Purify PCR product by running the DNA on an agarose gel and excising the band or by using a spin column (e.g., Monarch® DNA Gel Extraction Kit, NEB #T1020, Monarch PCR & DNA Cleanup Kit, NEB #T1030)
- Digest with the appropriate restriction enzyme

Standard Restriction Enzyme Protocol

DNA	1 µg
10X NEBuffer	5 µl (1X)
Restriction Enzyme	10 units is sufficient, generally 1 µl is used
Nuclease-free Water	To 50 µl
Incubation Time	1 hour*
Incubation Temperature	Enzyme dependent

* Can be decreased by using a Time-Saver qualified enzyme

Time-Saver Restriction Enzyme Protocol

DNA	1 µg
10X NEBuffer	5 µl (1X)
Restriction Enzyme	1 µl
Nuclease-free Water	To 50 µl
Incubation Time	5–15 minutes*
Incubation Temperature	Enzyme dependent

* Time-Saver qualified enzymes can also be incubated overnight with no star activity

Insert from Annealed Oligos

- Annealed oligos can be used to introduce a fragment (e.g., promoter, polylinker, etc.)
- Anneal two complementary oligos that leave protruding 5' or 3' overhangs for ligation into a vector cut with appropriate enzymes
- Non-phosphorylated oligos can be phosphorylated using T4 Polynucleotide Kinase (NEB #M0201)

Typical Annealing Reaction

Oligo 1	20 µM Final concentration
Oligo 2	20 µM Final concentration
10X NEBuffer r2.1 or 2.1	5 µl
Nuclease-free Water	To 50 µl
Incubation	95°C for 5 minutes, cool slowly to room temp.

Vector

- Digest vector with appropriate restriction enzymes. Enzymes that leave non-compatible ends are ideal as they prevent vector self-ligation

DEPHOSPHORYLATION

- Dephosphorylation is sometimes necessary to prevent self-ligation. NEB offers four products for dephosphorylation of DNA:
- Quick CIP (NEB #M0525), Shrimp Alkaline Phosphatase (rSAP) (NEB #M0371) and Antarctic Phosphatase (AP) (NEB #M0289) are heat-inactivatable phosphatases. They work in all NEBuffers, but AP requires supplementation with Zn²⁺.

Dephosphorylation of 5' ends of DNA Using the Quick CIP

DNA	1 pmol of DNA ends
10X rCutSmart or CutSmart Buffer	2 µl
Quick CIP	1 µl
Nuclease-free Water	To 20 µl
Incubation	37°C for 10 minutes
Heat Inactivation	80°C for 2 minutes

Note: Scale larger reaction volumes proportionally.

BLUNTING

- In some instances, the ends of the insert or vector require blunting
- PCR with a proofreading polymerase will leave a predominantly blunt end
- T4 DNA Polymerase (NEB #M0203) or Klenow (NEB #M0210) will fill in a 5' overhang and chew back a 3' overhang
- The Quick Blunting Kit (NEB #E1201) is optimized to blunt and phosphorylate DNA ends for cloning in less than 30 minutes
- Analyze agarose gels with longwave UV (365 nm) to minimize UV exposure that may cause DNA damage

Blunting with the Quick Blunting Kit

DNA	Up to 5 µg
10X Blunting Buffer	2.5 µl
dNTP Mix (1 mM)	2.5 µl
Blunt Enzyme Mix	1 µl
Nuclease-free Water	To 25 µl
Incubation	room temperature; 15 min for RE-digested DNA; 30 min for sheared/nebulized DNA or PCR products*
Heat Inactivation	70°C for 10 minutes

* PCR-generated DNA must be purified before blunting using a purification kit (NEB #T1030), phenol extraction/ethanol precipitation, or gel extraction (NEB #T1020).

Traditional Cloning Quick Guide (continued)

PHOSPHORYLATION

- For ligation to occur, at least one of the DNA ends (insert or vector) should contain a 5' phosphate
- Primers are usually supplied non-phosphorylated; therefore, the PCR product will not contain a 5' phosphate
- Digestion of DNA with a restriction enzyme will always produce a 5' phosphate
- A DNA fragment can be phosphorylated by incubation with T4 Polynucleotide Kinase (T4 PNK, NEB #M0201). T4 PNK can be inactivated at 65°C for 20 minutes.

Phosphorylation with T4 PNK

DNA (20 mer)	up to 300 pmol of 5' termini
10X T4 PNK Buffer	5 µl
10 mM ATP	5 µl (1 mM final conc.)
T4 PNK	1 µl (10 units)
Nuclease-free Water	To 50 µl
Incubation	37°C for 30 minutes

PURIFICATION OF VECTOR AND INSERT

- Purify the vector and insert by either running the DNA on an agarose gel and excising the appropriate bands or by using a spin column, such as Monarch DNA Gel Extraction Kit or PCR & DNA Cleanup Kit (NEB #T1020 or T1030)
- DNA can also be purified using β -Agarase I (NEB #M0392) with low melt agarose, or an appropriate spin column or resin
- Analyze agarose gels with longwave UV (365 nm) to minimize UV exposure that may cause DNA damage

LIGATION OF VECTOR AND INSERT

- Use a molar ratio of 1:3 vector to insert. Use NEBioCalculator to calculate molar ratios.
- If using T4 DNA Ligase (NEB #M0202) or the Quick Ligation Kit (NEB #M2200), thaw and resuspend the Ligase Buffer at room temp. If using Ligase Master Mixes, no thawing is necessary.
- The Quick Ligation Kit (NEB #M2200) is optimized for ligation of both sticky and blunt ends
- Instant Sticky-end Ligase Master Mix (NEB #M0370) is optimized for instant ligation of sticky/cohesive ends
- Blunt/TA Ligase Master Mix (NEB #M0367) is optimized for ligation of blunt or single base overhangs, which are the more challenging type of ends for T4 DNA Ligase
- Following ligation, chill on ice and transform
- DO NOT heat inactivate when using the Quick Ligation Buffer or Ligase Master Mixes, as this will inhibit transformation
- Electroligase (NEB #M0369) is optimized for ligation of both sticky and blunt ends and is compatible with electroporation (i.e., no cleanup step required)
- Improved Golden Gate Assembly can be achieved by selecting high fidelity overhangs [Potapov, V. et al. (2018) *ACS Synth. Biol.* 7(11), 2665–2674.

Ligation with the Quick Ligation Kit

Vector DNA (3 kb)	50 ng
Insert DNA (1 kb)	To 50 ng
2X Quick Ligation Buffer	10 µl
Quick T4 DNA Ligase	1 µl
Nuclease-free Water	20 µl (mix well)
Incubation	Room temperature for 5 minutes

Ligation with Instant Sticky-end Ligase Master Mix

Vector DNA (3 kb)	50 ng
Insert DNA (1 kb)	50 ng
Master Mix	5 µl
Nuclease-free Water	To 10 µl
Incubation	None

Ligation with Blunt/TA Ligase Master Mix

Vector DNA (3 kb)	50 ng
Insert DNA (1 kb)	50 ng
Master Mix	5 µl
Nuclease-free Water	To 10 µl
Incubation	Room temperature for 15 minutes

TRANSFORMATION

- To obtain transformants in 8 hrs., use NEB Turbo Competent *E. coli* (NEB #C2984)
- If recombination is a concern, then use the *recA*⁻ strains NEB 5-alpha Competent *E. coli* (NEB #C2987), NEB-10 beta Competent *E. coli* (NEB #C3019) or NEB Stable Competent *E. coli* (NEB #C3040)
- NEB 10-beta Competent *E. coli* works well for constructs larger than 5 kb
- NEB Stable Competent *E. coli* (NEB #C3040) can be used for constructs with repetitive sequences such as lentiviral constructs
- If electroporation is required, use NEB 5-alpha (NEB #C2989) or NEB 10-beta (NEB #C3020) Electrocompetent *E. coli*
- Use pre-warmed selection plates
- Perform several 10-fold serial dilutions in SOC or NEB 10-beta/Stable Outgrowth Medium for plating

Transformation with NEB 5-alpha Competent *E. coli*

DNA	1–5 µl containing 1 pg – 100 ng of plasmid DNA
Competent <i>E. coli</i>	50 µl
Incubation	On ice for 30 minutes
Heat Shock	Exactly 42°C for exactly 30 seconds
Incubation	On ice for 5 minutes Add 950 µl room temperature SOC 37°C for 60 minutes, with shaking

Troubleshooting Guide for Cloning

We strongly recommend running the following controls during transformations. These controls may assist in identifying which step(s) in the cloning workflow has failed.

- Transform 100 pg – 1 ng of uncut vector to check cell viability, calculate transformation efficiency and verify the antibiotic resistance of the plasmid.
- Transform the cut vector to determine the amount of background due to undigested plasmid. The number of colonies in this control should be < 1% of the number of colonies in the uncut plasmid control transformation (from control #1).
- Transform a vector only ligation reaction. The ends of the vector should not be able to re-ligate because either they are incompatible (e.g., digested with two restriction enzymes that do not generate compatible ends) or the 5' phosphate group has been removed in a dephosphorylation reaction (e.g., blunt ends treated with rSAP). This control transformation should yield the same number of colonies as control #2.
- Digest vector DNA with a single restriction enzyme, re-ligate and transform. The ends of the vector DNA should be compatible and easily joined during the ligation reaction, resulting in approximately the same number of colonies as control #1.

The cloning workflow often benefits from an accurate quantitation of the amount of DNAs that are being worked with. We recommend quantification of DNAs whenever possible.

PROBLEM	CAUSE	SOLUTION
Few or no transformants	Cells are not viable	<ul style="list-style-type: none"> • Transform an uncut plasmid (e.g., pUC19) and calculate the transformation efficiency of the competent cells. If the transformation efficiency is low (< 10⁴) re-make the competent cells or consider using commercially available high efficiency competent cells.
	Incorrect antibiotic or antibiotic concentration	<ul style="list-style-type: none"> • Confirm antibiotic and antibiotic concentration
	DNA fragment of interest is toxic to the cells	<ul style="list-style-type: none"> • Incubate plates at lower temperature (25–30°C). • Transformation may need to be carried out using a strain that exerts tighter transcriptional control over the DNA fragment of interest (e.g., NEB 5-alpha F' I^a Competent <i>E. coli</i> (NEB #C2992))
	If using chemically competent cells, the wrong heat-shock protocol was used	<ul style="list-style-type: none"> • Follow the manufacturer's specific transformation protocol (Note: going above the recommended temperature during the heat shock can result in competent cell death)
	If using electrocompetent cells, PEG is present in the ligation mix	<ul style="list-style-type: none"> • Clean up DNA prior to transformation with the Monarch PCR & DNA Cleanup Kit (NEB #T1030) • Try NEB's ElectroLigase (NEB #M0369)
	If using electrocompetent cells, arcing was observed or no voltage was registered	<ul style="list-style-type: none"> • Clean up the DNA prior to ligation with the Monarch PCR & DNA Cleanup Kit (NEB #T1030) • Tap the cuvette to get rid of any trapped air bubbles • Be sure to follow the manufacturer's specified electroporation parameters
	Construct is too large	<ul style="list-style-type: none"> • Select a competent cell strain that can be transformed efficiently with large DNA constructs (≥ 10 kb, we recommend trying NEB 10-beta Competent <i>E. coli</i> (NEB #C3019)) or NEB Stable Competent <i>E. coli</i> (NEB #C3040) • For very large constructs (> 10 kb), consider using electroporation
	Construct may be susceptible to recombination	<ul style="list-style-type: none"> • Select a <i>recA</i>- strain such as NEB 5-alpha (NEB #C2987), NEB 10-beta (NEB #C3019) or NEB Stable (NEB #C3040) Competent <i>E. coli</i>
	The insert comes directly from mammalian or plant DNA and contains methylated cytosines, which are degraded by many <i>E. coli</i> strains	<ul style="list-style-type: none"> • Use a strain that is deficient in McrA, McrBC and Mrr, such as NEB 10-beta Competent <i>E. coli</i>
	Too much ligation mixture was used	<ul style="list-style-type: none"> • Use < 5 µl of the ligation reaction for the transformation
	Inefficient ligation	<ul style="list-style-type: none"> • Make sure that at least one fragment being ligated contains a 5' phosphate moiety • Vary the molar ratio of vector to insert from 1:1 to 1:10. Use NEBioCalculator to calculate molar ratios • Purify the DNA to remove contaminants such as salt and EDTA with Monarch PCR & DNA Cleanup Kit (5 µg) (NEB #T1030) • ATP will degrade after multiple freeze-thaws; repeat the ligation with fresh buffer • Heat inactivate or remove the phosphatase prior to ligation • Ligation of single base-pair overhangs (most difficult) may benefit from being carried out with Blunt/TA Master Mix (NEB #M0367), Quick Ligation Kit (NEB #M2200) or concentrated T4 DNA Ligase (NEB #M0202)
	Inefficient phosphorylation	<ul style="list-style-type: none"> • Purify the DNA prior to phosphorylation with Monarch PCR & DNA Cleanup Kit (5 µg) (NEB #T1030). Excess salt, phosphate or ammonium ions may inhibit the kinase. • If the ends are blunt or 5' recessed, heat the substrate/buffer mixture for 10 minutes at 70°C. Rapidly chill on ice before adding the ATP and enzyme, then incubate at 37°C. • ATP was not added. Supplement the reaction with 1 mM ATP, as it is required by T4 Polynucleotide Kinase (NEB #M0201) • Alternatively, use 1X T4 DNA Ligase Buffer (contains 1 mM ATP) instead of the 1X T4 PNK Buffer

Troubleshooting Guide for Cloning (continued)

PROBLEM	CAUSE	SOLUTION
Few or no transformants	Inefficient blunting	<ul style="list-style-type: none"> Heat inactivate or remove the restriction enzymes prior to blunting Clean up the PCR fragment prior to blunting with Monarch PCR & DNA Cleanup Kit (NEB #T1030) Sonicated gDNA should be blunted for at least 30 minutes Do not use > 1 unit of enzyme/μg of DNA Do not incubate for > 15 minutes Do not incubate at temperatures > 12°C (for T4 DNA Polymerase, NEB #M0203) or > 24°C (for Klenow, NEB #M0210) Make sure to add a sufficient amount of dNTPs to the reaction (33 μM each dNTP for DNA Polymerase I, Large (Klenow) Fragment, NEB #M0210 and 100 μM each dNTP for T4 DNA Polymerase, NEB #M0203). When using Mung Bean Nuclease (NEB #M0250), incubate the reaction at room temperature. Do not use > 1 unit of enzyme/μg DNA or incubate the reaction > 30 minutes.
	Inefficient A-Tailing	<ul style="list-style-type: none"> Clean up the PCR prior to A-tailing. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030). High-fidelity polymerases will remove any non-templated nucleotides.
	Restriction enzyme(s) didn't cleave completely	<ul style="list-style-type: none"> Check the methylation sensitivity of the enzyme(s) to determine if the enzyme is blocked by methylation of the recognition sequence Use the recommended buffer supplied with the restriction enzyme Clean up the DNA to remove any contaminants that may inhibit the enzyme. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030). When digesting a PCR fragment, make sure to have at least 6 nucleotides between the recognition site and the end of the DNA molecule
Colonies don't contain a plasmid	Antibiotic level used was too low	<ul style="list-style-type: none"> Increase the antibiotic level on plates to the recommended amount Use fresh plates with fresh antibiotics
	Satellite colonies were selected	<ul style="list-style-type: none"> Choose large, well-established colonies for analysis
Colonies contain the wrong construct	Recombination of the plasmid has occurred	<ul style="list-style-type: none"> Use a <i>recA</i>⁻ strain such as NEB 5-alpha, NEB 10-beta or NEB Stable Competent <i>E. coli</i>
	Incorrect PCR amplicon was used during cloning	<ul style="list-style-type: none"> Optimize the PCR conditions Gel purify the correct PCR fragment. NEB recommends the Monarch DNA Gel Extraction Kit (NEB #T1020).
	Internal recognition site was present	<ul style="list-style-type: none"> Use NEBcutter to analyze insert sequence for presence of an internal recognition site
	DNA fragment of interest is toxic to the cells	<ul style="list-style-type: none"> Incubate plates at lower temperature (25–30°C) Transformation may need to be carried out using a strain that exerts tighter transcriptional control of the DNA fragment of interest (e.g., NEB 5-alpha F' I^q Competent <i>E. coli</i>)
	Mutations are present in the sequence	<ul style="list-style-type: none"> Use a high-fidelity polymerase (e.g., Q5 High-Fidelity DNA Polymerase, NEB #M0491) Re-run sequencing reactions
Too much background	Inefficient dephosphorylation	<ul style="list-style-type: none"> Heat inactivate or remove the restriction enzymes prior to dephosphorylation
	Kinase is present/active	<ul style="list-style-type: none"> Heat inactivate the kinase after the phosphorylation step. Active kinase will re-phosphorylate the dephosphorylated vector.
	Restriction enzyme(s) didn't cleave completely	<ul style="list-style-type: none"> Check the methylation sensitivity of the restriction enzyme(s) to be sure it is not inhibited by methylation of the recognition sequence Use the recommended buffer supplied with the restriction enzyme Clean up the DNA to remove contaminants. (e.g., too much salt). NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030).
	Antibiotic level is too low	<ul style="list-style-type: none"> Confirm the correct antibiotic concentration
Ran the ligation on a gel and saw no ligated product	Inefficient ligation	<ul style="list-style-type: none"> Make sure at least one DNA fragment being ligated contains a 5' phosphate Vary the molar ratios of vector to insert from 1:1 to 1:10. Use NEBioCalculator to calculate molar ratios. Purify the DNA to remove contaminants such as salt and EDTA. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030). ATP will degrade after multiple freeze-thaws; repeat the ligation with fresh buffer Heat inactivate or remove the phosphatase prior to ligation Ligation of single base-pair overhangs (most difficult) may benefit from being carried out with Blunt/TA Master Mix, Quick Ligation Kit or concentrated T4 DNA Ligase
The ligated DNA ran as a smear on an agarose gel	The ligase is bound to the substrate DNA	<ul style="list-style-type: none"> Treat the ligation reaction with Proteinase K (NEB #P8107) prior to running on a gel
The digested DNA ran as a smear on an agarose gel	The restriction enzyme(s) is bound to the substrate DNA	<ul style="list-style-type: none"> Lower the number of units Add SDS (0.1–0.5%) to the loading buffer to dissociate the enzyme from the DNA or use Gel Loading Dye, Purple (6X) (NEB #B7024)
	Nuclease contamination	<ul style="list-style-type: none"> Use fresh, clean running buffer Use a fresh agarose gel Clean up the DNA. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030).

PROBLEM	CAUSE	SOLUTION
Incomplete restriction enzyme digestion	Cleavage is blocked by methylation	<ul style="list-style-type: none"> Check the methylation sensitivity of the enzyme(s) to determine if the enzyme is blocked by methylation of the recognition sequence DNA isolated from a bacterial source may be blocked by Dam and Dcm methylation, which some enzymes are sensitive to If the enzyme is inhibited by Dam or Dcm methylation, grow the plasmid in a dam-/dcm- strain (NEB #C2925) or use PCR DNA DNA isolated from eukaryotic source may be blocked by CpG methylation
	Salt inhibition	<ul style="list-style-type: none"> Enzymes that have low activity in salt-containing buffers (NEBuffer r3.1) may be salt sensitive, so clean up the DNA prior to digestion. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030). DNA purification procedures that use spin columns can result in high salt levels, which inhibit enzyme activity. Monarch kits (NEB #T1010, #T1020, #T1030) use columns that have been designed to minimize salt carry over into the eluted DNA, so using them can minimize this issue. To prevent this, DNA solution should be no more than 25% of total reaction volume.
	Inhibition by PCR components	Clean up the PCR fragment prior to restriction digest. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030).
	Using the wrong buffer	Use the recommended buffer supplied with the restriction enzyme
	Too few units of enzyme used	Use at least 5–10 units of enzyme per µg of DNA
	Incubation time was too short	Increase the incubation time
	Digesting supercoiled DNA	Some enzymes have a lower activity on supercoiled DNA. Increase the number of enzyme units in the reaction.
	Presence of slow sites	Some enzymes can exhibit slower cleavage towards specific sites. Increase the incubation time, 1–2 hours is typically sufficient.
	Two sites required	Some enzymes require the presence of two recognition sites to cut efficiently. For more information, visit the table "Restriction Enzymes Requiring Multi-sites" on neb.com.
	DNA is contaminated with an inhibitor	<ul style="list-style-type: none"> Assay substrate DNA in the presence of a control DNA. Control DNA will not cleave if there is an inhibitor present. Miniprep DNA is particularly susceptible to contaminants. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030). Clean DNA with a spin column, we recommend the Monarch PCR & DNA Cleanup Kit (NEB #T1030), or increase volume to dilute contaminant
Extra bands in the gel	If larger bands than expected are seen in the gel, this may indicate binding of the enzyme(s) to the substrate	<ul style="list-style-type: none"> Lower the number of units in the reaction Add SDS (0.1–0.5%) to the loading buffer to dissociate the enzyme from the substrate or add Gel Loading Dye, Purple (6X) (NEB #B7024)
	Star activity	<ul style="list-style-type: none"> Use the recommended buffer supplied with the restriction enzyme Decrease the number of enzyme units in the reaction Make sure the amount of enzyme added does not exceed 10% of the total reaction volume. This ensures that the total glycerol concentration does not exceed 5% v/v Decrease the incubation time. Using the minimum reaction time required for complete digestion will help prevent star activity. Try using a High-Fidelity (HF) restriction enzyme. HF enzymes have been engineered for reduced star activity.
	Partial restriction enzyme digest	<ul style="list-style-type: none"> Enzymes that have low activity in salt-containing buffers (e.g., NEBuffer r3.1 or 3.1) may be salt sensitive. Make sure to clean up the DNA prior to digestion. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030). DNA purification procedures that use spin columns can result in high salt levels, which inhibit enzyme activity. Monarch kits (NEB #T1010, #T1020, #T1030) use columns that have been designed to minimize salt carry over into the eluted DNA, so using them can minimize this issue. To prevent this, DNA solution should be no more than 25% of total reaction volume Clean-up the PCR fragment prior to restriction digest. NEB recommends the Monarch PCR & DNA Cleanup Kit (NEB #T1030). Use the recommended buffer supplied with the restriction enzyme Use at least 3–5 units of enzyme per µg of DNA Digest the DNA for 1–2 hours
No PCR fragment amplified	Used the wrong primer sequence	Double check the primer sequence
	Incorrect annealing temperature	Use the NEB Tm calculator to determine the correct annealing temperature (www.neb.com/TmCalculator)
	Incorrect extension temperature	Each polymerase type has a different extension temperature requirement. Follow the manufacturer's recommendations.
	Too few units of polymerase	Use the recommended number of polymerase units based on the reaction volume
	Incorrect primer concentration	Each polymerase has a different primer concentration requirement. Make sure to follow the manufacturer's recommendations.
	Mg ²⁺ levels in the reaction are not optimal	Titrate the Mg ²⁺ levels to optimize the amplification reaction. Follow the manufacturer's recommendations.
	Difficult template	With difficult templates, try different polymerases and/or buffer combinations
The PCR reaction is a smear on a gel	If bands are larger than expected it may indicate binding of the enzyme(s) to the DNA	Add SDS (0.1–0.5%) to the loading buffer to dissociate the enzyme from the DNA
Extra bands in PCR reaction	Annealing temperature is too low	Use the NEB Tm calculator to determine the annealing temperature of the primers
	Mg ²⁺ levels in the reaction are not optimal	Titrate the Mg ²⁺ levels to optimize the amplification reaction. Make sure to follow the manufacturer's recommendations.
	Additional priming sites are present	Double check the primer sequence and confirm it does not bind elsewhere in the DNA template
	Formation of primer dimers	Primer sequence may not be optimal. Additional primers may need to be tested in the reaction.
	Incorrect polymerase choice	Try different polymerases and/or buffer combinations

Optimization Tips for Your Cloning Reactions

New England Biolabs offers a wide selection of reagents for your cloning experiments. For more information, visit ClonewithNEB.com. The following tips can be used to help optimize each step in your cloning workflow. Tips for restriction enzyme digestion and amplification can be found on pages 306 and 322, respectively.

CDNA SYNTHESIS

Starting Material

- Intact RNA of high purity is essential for generating cDNA for cloning applications
- Total RNA or mRNA can be used in the reverse transcription reaction. Total RNA is generally sufficient for cDNA synthesis reactions. However, if desired, mRNA can be easily obtained using a PolyA Spin mRNA Isolation Kit (NEB #S1560) or Magnetic mRNA Isolation Kit (NEB #S1550).
- The amount of RNA required for cDNA cloning depends on the abundance of the transcript-of-interest. In general, 1 ng to 1 µg total RNA or 0.1–100 ng mRNA are recommended.

Product Selection

- Streamline your reaction setup by using the ProtoScript II First Strand cDNA Synthesis Kit (NEB #E6560). This kit combines ProtoScript II Reverse Transcriptase (NEB #M0360), a thermostable M-MuLV (RNase H⁻) Reverse

Transcriptase, and recombinant RNase Inhibitor in an enzyme Master Mix, along with a separate Reaction Mix containing dNTPs. Additionally, the kit contains two optimized reverse transcription primer mixes.

Yield

- ProtoScript II Reverse Transcriptase is capable of generating cDNA of more than 10 kb up to 48°C. We recommend 42°C for routine reverse transcription.
- You can increase the yield of a long cDNA product by doubling the amount of enzyme and dNTPs

Additives

- For most RT-PCR reactions, RNase H treatment is not required. But for some difficult amplicons or sensitive assays, add 2 units of *E. coli* RNase H to the reaction and incubate at 37°C for 20 minutes

PHOSPHORYLATION

Enzyme

- T4 Polynucleotide Kinase (NEB #M0201) and T4 DNA Ligase (NEB #M0202) can be used together in the T4 DNA Ligase Buffer
- T4 Polynucleotide Kinase is inhibited by high levels of salt (50% inhibition by 150 mM NaCl), phosphate (50% inhibition by 7 mM phosphate) and ammonium ions (75% inhibited by 7 mM (NH₄)₂SO₄)

- If using T4 Polynucleotide Kinase and working with 5'-recessed ends, heat the reaction mixture for 10 min at 70°C, chill rapidly on ice before adding the ATP (or Ligase Buffer containing ATP) and enzyme, then incubate at 37°C

Additives

- The addition of PEG 8000 (up to 5%) can improve results

DEPHOSPHORYLATION

Enzyme

- When dephosphorylating a fragment following a restriction enzyme digest, a DNA clean up step is required if the restriction enzyme(s) used is NOT heat inactivatable. We recommend the Monarch PCR & DNA Cleanup Kit (NEB #T1030).
- When working with Quick CIP (NEB #M0525), rSAP (NEB #M0371) or AP (NEB #M0289), which are heat-inactivatable enzymes, a DNA clean-up step after dephosphorylation is not necessary prior to the ligation step.

Additives

- AP requires the presence of Zn²⁺ in the reaction, so don't forget to supplement the reaction with 1X Antarctic Phosphatase Reaction Buffer when using other NEBuffers

BLUNTING/END REPAIR

Enzyme

- Make sure that you choose the correct enzyme to blunt your fragment. The Quick Blunting Kit (NEB #E1201), T4 DNA Polymerase (NEB #M0203) and DNA Polymerase I, Large (Klenow) Fragment (NEB #M0210) will fill 5' overhangs and degrade 3' overhangs. Mung Bean Nuclease (NEB #M0250) degrades 5' overhangs.
- T4 DNA Polymerase and DNA Polymerase I, Large (Klenow) Fragment are active in all NEBuffers. Please remember to add dNTPs.

Clean-up

- When trying to blunt a fragment after a restriction enzyme digestion, if the restriction enzyme(s) used are heat inactivatable, then a clean-up step prior to blunting is not needed. Alternatively, if the restriction enzyme(s) used are not heat inactivatable, a DNA clean-up step is recommended prior to blunting.
- When trying to blunt a fragment amplified by PCR, a DNA clean-up step is necessary prior to the blunting step to remove the nucleotides and polymerase

- When trying to dephosphorylate a fragment after the blunting step, you will need to add a DNA clean-up step (e.g., Monarch PCR & DNA Cleanup Kit, NEB #T1030) after the blunting and before the addition of the phosphatase

Temperature

- When trying to blunt a fragment with Mung Bean Nuclease, the recommended temperature of incubation is room temperature, since higher temperatures may cause sufficient breathing of the dsDNA ends that the enzyme may degrade some of the dsDNA sequence. The number of units to be used and time of incubation may be determined empirically to obtain best results.

Heat Inactivation

- Mung Bean Nuclease reactions should not be heat inactivated. Although Mung Bean Nuclease can be inactivated by heat, this is not recommended because the DNA begins to "breathe" before the Mung Bean Nuclease is inactivated and undesirable degradation occurs at breathing sections. Purify DNA by phenol/chloroform extraction and ethanol precipitation or spin column purification [e.g., Monarch PCR & DNA Cleanup Kit (NEB #T1030)].

A-TAILING

- If the fragment to be tailed has been amplified with a high-fidelity polymerase, the DNA needs to be purified prior to the tailing reaction. For this we recommend the Monarch PCR & DNA Cleanup Kit (NEB T1030). Otherwise,

any high-fidelity polymerase present in the reaction will be able to remove any non-templated nucleotides added to the end of the fragments.

DNA LIGATION

Reaction Buffers

- T4 DNA Ligase Buffer (NEB #B0202) should be thawed on the bench or in the palm of your hand, and not at 37°C (to prevent breakdown of ATP)
- Once thawed, T4 DNA Ligase Buffer should be placed on ice
- Ligations can be performed in any of the four standard restriction endonuclease NEBuffers or in T4 Polynucleotide Kinase Buffer (NEB #B0201) supplemented with 1 mM ATP
- When supplementing with ATP, use ribo-ATP (NEB #P0756). Deoxyribo-ATP will inhibit ligation.
- Before ligation, completely inactivate restriction enzyme by heat inactivation, spin column (e.g., Monarch PCR & DNA Cleanup Kit, NEB #T1030) or Phenol/EtOH purification

DNA

- Heat inactivate (AP, SAP, Quick Dephosphorylation Kit) before ligation
- Keep total DNA concentration between 5–10 µg/ml
- Vector:Insert molar ratios between 1:1 and 1:10 are optimal for single insertions. Use NEBioCalculator at NEBioCalculator.neb.com to calculate molar ratios.
- For cloning more than one insert, we recommend the NEBuilder® HiFi DNA Assembly Master Mix (NEB #E2621) or Cloning Kit (NEB #E5520)

- If you are unsure of your DNA concentration, perform multiple ligations with varying ratios

Ligase

- For cohesive-end ligations, standard T4 DNA Ligase, Instant Sticky-end Ligase Master Mix or the Quick Ligation Kit are recommended.
- For blunt and single-base overhangs the Blunt/TA Ligase Master Mix is recommended
- For ligations that are compatible with electroporation, Electroligase is recommended
- Standard T4 DNA Ligase can be heat inactivated at 65°C for 20 minutes
- Do not heat inactivate the Quick Ligation Kit or the ligase master mixes

Transformation

- Add between 1–5 µl of ligation mixture to competent cells for transformation
- Extended ligation with PEG causes a drop off in transformation efficiency
- Electroporation is recommended for larger constructs (> 10,000 bp). Dialyze samples or use a spin column first if you have used the Quick Ligation Kit or ligase master mixes.
- Electroligase is recommended for ligations that will be electroporated

TRANSFORMATION

Thawing

- Cells are best thawed on ice
- DNA should be added as soon as the last trace of ice in the tube disappears
- Cells can be thawed by hand, but warming decreases efficiency

DNA

- Up to 10 µl of DNA from a ligation mix can be used with only a 2-fold loss of efficiency

Incubation & Heat Shock

- Incubate on ice for 30 minutes. Expect a 2-fold loss in transformation efficiency for every 10 minutes this step is shortened.
- Both temperature and time are specific to the transformation volume and vessel. Typically, 30 seconds at 42°C is recommended, except when using BL21 (NEB #C2530) which requires exactly 10 seconds.

Outgrowth

- Outgrowth at 37°C for 1 hour is best for cell recovery and for expression of antibiotic resistance. Expect a 2-fold loss in transformation efficiency for every 15 minutes this step is shortened.
- SOC and NEB 10-beta/Stable Outgrowth Medium give 2-fold higher transformation efficiency than LB medium
- Incubation with shaking or rotation results in 2-fold higher transformation efficiency

Plating

- Selection plates can be used warm or cold, wet or dry with no significant effects on transformation efficiency
- Warm, dry plates are easier to spread and allow for the most rapid colony formation

DNA Contaminants to Avoid

CONTAMINANT	REMOVAL METHOD
Detergents	Ethanol precipitate
Phenol	Extract with chloroform and ethanol precipitate
Ethanol or Isopropanol	Dry pellet before resuspending
PEG	Column purify (e.g., Monarch PCR & DNA Cleanup Kit) or phenol/chloroform extract and ethanol precipitate

Troubleshooting Guide for DNA Cleanup & Plasmid Purification using Monarch® Kits

PROBLEM	PRODUCT	POSSIBLE CAUSE	SOLUTION
No DNA purified	Monarch Plasmid Miniprep Kit (NEB #T1010)	Buffers added incorrectly	<ul style="list-style-type: none"> Add buffers in the correct order so that the sample is bound, washed and eluted in the correct sequence Ensure ethanol was added to Plasmid Wash Buffer 2
		Plasmid loss during culture growth	<ul style="list-style-type: none"> Ensure proper antibiotic and concentration was used to maintain selection during culture growth
	Monarch DNA Gel Extraction Kit (NEB #T1020)	Ethanol not added to wash buffer	<ul style="list-style-type: none"> Ensure the proper amount of ethanol was added to Monarch DNA Wash Buffer
	Monarch PCR & DNA Cleanup Kit (5 µg) (NEB #T1030)		
Low DNA yield	Monarch Plasmid Miniprep Kit (NEB #T1010)	Incomplete lysis	<ul style="list-style-type: none"> Pellet must be completely resuspended before addition of Plasmid Lysis Buffer (B2) – color should change from light to dark pink Avoid using too many cells; this can overload the column. If culture volume is larger than recommended, scale up buffers B1-B3 or consider using two columns.
		Too many cells used	<ul style="list-style-type: none"> The amount of cells used should not exceed 15 O.D. units
		Plasmid loss during culture growth	<ul style="list-style-type: none"> Ensure proper antibiotic and concentration was used to maintain selection during culture growth
		Low-copy plasmid selected	<ul style="list-style-type: none"> Increase amount of cells processed and scale buffers accordingly
		Lysis of cells during growth	<ul style="list-style-type: none"> Harvest culture during transition from logarithmic growth to stationary phase (~12-16 hours)
		Incomplete neutralization	<ul style="list-style-type: none"> Invert tube several times until color changes to a uniform yellow color
		Incomplete elution	<ul style="list-style-type: none"> Deliver Elution Buffer directly to center of column Larger elution volumes and longer incubation times can increase yield For elution of plasmids > 10 kb, heat the DNA Elution Buffer to 50°C and extend incubation time to 5 minutes
	Monarch DNA Gel Extraction Kit (NEB #T1020)	Buffers added incorrectly	<ul style="list-style-type: none"> Be sure that buffers have been reconstituted correctly and that reagents have been added in the correct order
		Gel slice not fully dissolved	<ul style="list-style-type: none"> Undissolved agarose may clog the column and interfere with binding. Incubate in Monarch Gel Dissolving Buffer for proper time and temperature.
		Gel dissolved above 60°C	<ul style="list-style-type: none"> Dissolve gel slice in specified range (37-55°C). Higher temperatures can denature DNA
		Incomplete elution during preparation	<ul style="list-style-type: none"> Deliver Elution Buffer directly to center of column Larger elution volumes and longer incubation times can increase yield For elution of DNA > 10 kb, heat the Elution Buffer to 50°C and extend incubation time to 5 minutes Multiple rounds of elution can also be performed
	Monarch PCR & DNA Cleanup Kit (5 µg) (NEB #T1030)	Buffers added incorrectly	<ul style="list-style-type: none"> Be sure that buffers have been reconstituted correctly and that reagents have been added in the correct order
		Incomplete elution during preparation	<ul style="list-style-type: none"> Deliver Elution Buffer directly to center of column Larger elution volumes and longer incubation times can increase yield For elution of DNA > 10 kb, heat the Elution Buffer to 50°C and extend incubation time to 5 minutes Multiple rounds of elution can also be performed
Low DNA quality	Monarch Plasmid Miniprep Kit (NEB #T1010)	Plasmid degradation	<ul style="list-style-type: none"> Be cautious of strains with high levels of endogenous endonuclease (e.g., HB101 and JM 100 series)
		Plasmid is denatured	<ul style="list-style-type: none"> Limit incubation with Plasmid Lysis Buffer (B2) to two minutes, as NaOH in the buffer can denature the plasmid
		gDNA contamination	<ul style="list-style-type: none"> Use careful inversion mixing after cell lysis to avoid shearing of host cell chromosomal DNA. Do not vortex.
		RNA contamination	<ul style="list-style-type: none"> Incubate sample in neutralization buffer for the full 2 minutes. For cell culture volumes > 3 ml, increase the spin after neutralization to 5 minutes.
		Improper storage	<ul style="list-style-type: none"> Elute DNA in DNA Elution Buffer or nuclease-free water, and store at -20°C. Do not store in solutions containing magnesium.
Low DNA purity	Monarch Plasmid Miniprep Kit (NEB #T1010)	Ethanol has been carried over	<ul style="list-style-type: none"> Centrifuge final wash for 1 minute to ensure complete removal Ensure column tip does not come in contact with flow through
		Excessive salt in sample	<ul style="list-style-type: none"> Use both plasmid wash buffers and do not skip wash steps
		Excessive carbohydrate has been carried over	<ul style="list-style-type: none"> Avoid strains with high amounts of endogenous carbohydrate (e.g., HB101 and JM 100 series). Be sure to follow protocol and include Plasmid Wash Buffer 1 step.
	Monarch DNA Gel Extraction Kit (NEB #T1020)	Guanidine carryover from Gel Dissolving Buffer	<ul style="list-style-type: none"> Use smaller agarose plugs, which will require less Gel Dissolving Buffer Use 3 or 3.5 volumes of Gel Dissolving Buffer instead of 4
		Ethanol has been carried over	<ul style="list-style-type: none"> Centrifuge final wash for 1 minute to ensure complete removal Ensure column tip does not come in contact with flow through
		Trace amounts of salts have been carried over	<ul style="list-style-type: none"> Ensure column tip does not come in contact with new tube for elution
	Monarch PCR & DNA Cleanup Kit (5 µg) (NEB #T1030)	Ethanol has been carried over	<ul style="list-style-type: none"> Centrifuge final wash for 1 minute to ensure complete removal Ensure column tip does not come in contact with flow through
		Trace amounts of salts have been carried over	<ul style="list-style-type: none"> Ensure column tip does not come in contact with new tube

Guidelines for Choosing Sample Input Amounts When Using the Monarch Genomic Purification Kit

Genomic DNA yield, purity and integrity vary immensely based on sample type, input amount and sample condition. Below, we have provided some empirical yield and DIN data from a wide variety of sample types, as well as guidance on the maximal input amounts for each of those samples when using the Monarch Genomic DNA Purification Kit. It is very important not to overload the column and the buffer system when extracting and purifying gDNA, as DNA yields, purity, integrity, and length may suffer.

SAMPLE TYPE	RECOMMENDED INPUT AMOUNT	TYPICAL YIELD (µg)	DIN	MAXIMUM INPUT AMOUNT
TISSUE*				
Tail (mouse)	10 mg	12–20	8.5–9.5	25 mg
Ear (mouse)	10 mg	18–21	8.5–9.5	10 mg
Liver (mouse and rat)	10 mg	15–30	8.5–9.5	15 mg
Kidney (mouse)	10 mg	10–25	8.5–9.5	10 mg
Spleen (mouse)	10 mg	30–70	8.5–9.5	10 mg
Heart (mouse)	10 mg	9–10	8.5–9.5	25 mg
Lung (mouse)	10 mg	14–20	8.5–9.5	15 mg
Brain (mouse and rat)	10 mg	4–10	8.5–9.5	12 mg
Muscle (mouse and rat)	10 mg	4–7	8.5–9.5	25 mg
Muscle (deer)	10 mg	5	8.5–9.5	25 mg
BLOOD**				
Human (whole)	100 µl	2.5–4	8.5–9.5	100 µl
Mouse	100 µl	1–3	8.5–9.5	100 µl
Rabbit	100 µl	3–4	8.5–9.5	100 µl
Pig	100 µl	3.5–5	8.5–9.5	100 µl
Guinea pig	100 µl	3–8	8.5–9.5	100 µl
Cow	100 µl	2–3	8.5–9.5	100 µl
Horse	100 µl	4–7	8.5–9.5	100 µl
Dog	100 µl	2–4	8.5–9.5	100 µl
Chicken (nucleated)	10 µl	30–45	8.5–9.5	10 µl
CELLS				
HeLa	1 x 10 ⁶ cells	7–9	9.0–9.5	5 x 10 ⁶ cells
HEK293	1 x 10 ⁶ cells	7–9	9.0–9.5	5 x 10 ⁶ cells
NIH3T3	1 x 10 ⁶ cells	6–7.5	9.0–9.5	5 x 10 ⁶ cells
BACTERIA				
<i>E. coli</i> (Gram-negative)	2 x 10 ⁹ cells	6–10	8.5–9.0	2 x 10 ⁹ cells
<i>Rhodobacter</i> sp. (Gram-negative)	2 x 10 ⁹ cells	6–10	8.5–9.0	2 x 10 ⁹ cells
<i>B. cereus</i> (Gram-positive)	2 x 10 ⁹ cells	6–9	8.5–9.0	2 x 10 ⁹ cells
ARCHAEA				
<i>T. kodakarensis</i>	2 x 10 ⁹ cells	3–5	8.5–9.0	2 x 10 ⁹ cells
YEAST				
<i>S. cerevisiae</i>	5 x 10 ⁷ cells	0.5–0.6	8.5–9.0	5 x 10 ⁷ cells
SALIVA/BUCCAL CELLS***				
Saliva (human)	200 µl	2–3	7.0–8.0	500 µl
Buccal swab (human)	1 swab	5–7	6.0–7.0	1 swab

* Tissue gDNA yields are shown for frozen tissue powder, frozen tissue pieces and RNAlater-stabilized tissue pieces. Though frozen tissue powder results in highly-intact gDNA, lower yields can be expected than when using frozen or RNAlater-stabilized tissue pieces. Residual nuclease activity in tissue pieces will cut the gDNA, resulting in a slightly smaller overall size; however, this gDNA is optimal for silica-based purification.

** Human whole blood samples stabilized with various anticoagulants (e.g., EDTA, citrate and heparin) and various counter-ions were evaluated and results were comparable in all cases. Additionally, all indicated blood samples were tested both as fresh and frozen samples, yielding comparable results. Human samples were donated by healthy individuals; yields from unhealthy donors may differ.

*** Buccal swabs and saliva samples partially consist of dead cell material with degraded gDNA. Therefore, the purified gDNA from those samples will naturally have lower DIN values.

Troubleshooting Guide for Genomic DNA Purification using the Monarch Kit

PROBLEM	CAUSE	SOLUTION
LOW YIELD		
Cells	Frozen cell pellet was thawed and/or resuspended too abruptly	• Thaw cell pellets slowly on ice and flick tube several times to release the pellet from bottom of tube. Use cold PBS, and resuspend gently by pipetting up and down 5–10 times until pellet is dissolved
	Cell Lysis Buffer was added concurrently with enzymes	• Add Proteinase K and RNase A to sample and mix well before adding the Cell Lysis Buffer
Blood	Blood was thawed, allowing for DNase activity	• Keep blood samples frozen and add Proteinase K, RNase A and Blood Lysis Buffer directly to the frozen samples
	Blood sample is too old	• Fresh (unfrozen) whole blood should not be older than 1 week. Older samples will show progressive DNA degradation and loss of yield.
	Formation of hemoglobin precipitates	• Species with high hemoglobin content (e.g., guinea pig) may accumulate insoluble hemoglobin complexes that clog the membrane. Reduce Proteinase K lysis time from 5 to 3 minutes.
Tissue	Tissue pieces are too large	• Cut starting material to the smallest possible pieces or grind with liquid nitrogen. In large tissue pieces, nucleases will destroy the DNA before the Proteinase K can lyse the tissue.
	Membrane is clogged with tissue fibers	• Proteinase K digestion of fibrous tissues (e.g., muscle, heart, skin, ear clips), brain tissue and all RNAlater-stabilized tissues leads to the release of small indigestible protein fibers, which block the binding sites of the silica membrane. To remove fibers, centrifuge lysate at maximum speed for 3 minutes, as indicated in the protocol. For ear clips and brain tissue, use no more than 12–15 mg input material.
	Sample was not stored properly	• Samples stored for long periods of time at room temperature, 4°C or -20°C, will show degradation and loss of gDNA. Flash freeze tissue samples with liquid nitrogen or dry ice and store at -80°C. Alternatively, use stabilizing reagents to protect the gDNA.
	Genomic DNA was degraded (common in DNase-rich tissues)	• Organ tissues (e.g., pancreas, intestine, kidney, liver) contain significant amounts of nucleases. Store properly to prevent DNA degradation. Keep on ice during sample preparation. Refer to the protocol for the recommended amount of starting material and Proteinase K to use.
	Column is overloaded with DNA	• Some organ tissues (e.g., spleen, kidney, liver) are extremely rich in genomic DNA. Using inputs larger than recommended will result in the formation of tangled, long-fragment gDNA that cannot be eluted from the silica membrane. Reduce the amount of input material.
	Incorrect amount of Proteinase K added	• Most samples are digested with 10 µl Proteinase K, but for brain, kidney and ear clips, use 3 µl.
DNA DEGRADATION		
Tissue	Tissue samples were not stored properly	• Samples stored for long periods of time at room temperature, 4°C or -20°C, will show degradation and loss of gDNA. Flash freeze tissue samples with liquid nitrogen or dry ice and store at -80°C. Alternatively, use stabilizing reagents to protect the gDNA.
	Tissue pieces are too large	• Cut starting material to the smallest possible pieces or grind with liquid nitrogen. In large tissue pieces, nucleases will degrade the DNA before Proteinase K can lyse the tissue.
	High DNase content of soft organ tissue	• Organ tissues (e.g., pancreas, intestine, kidney, liver) contain significant amounts of nucleases. Store properly to prevent DNA degradation. Keep on ice during sample preparation. Refer to the protocol for the recommended amount of starting material and Proteinase K to use.
Blood	Blood sample is too old	• Fresh (unfrozen) whole blood should not be older than 1 week. Older samples will show progressive DNA degradation and loss of yield.
	Blood was thawed, allowing for DNase activity	• Keep frozen blood samples frozen and add enzymes and lysis buffer directly to the frozen samples
SALT CONTAMINATION		
	Guanidine Thiocyanate salt from the binding buffer was carried over into the eluate	<ul style="list-style-type: none"> • When transferring the lysate/binding buffer mix, avoid touching the upper column area with the pipet tip and always pipet carefully onto the silica membrane. • Avoid transferring any foam that may have been present in the lysate; foam can enter into the cap area of the spin column. • Close the caps gently to avoid splashing the mixture into the upper column area and move the samples with care in and out of the centrifuge. • If salt contamination is a concern, invert the columns a few times (or vortex briefly) with gDNA Wash Buffer as indicated in the protocol.
PROTEIN CONTAMINATION		
Tissue	Incomplete digestion	• Cut samples to the smallest possible pieces. Incubate sample in the lysis buffer for an extra 30 minutes to 3 hours to degrade any remaining protein complexes.
	Membrane is clogged with tissue fibers	• Proteinase K digestion of fibrous tissues (e.g., muscle, heart, skin, ear clips), brain tissue and all RNAlater-stabilized tissues leads to the release of small, indigestible protein fibers, which block the binding sites of the silica membrane. To remove fibers, centrifuge the lysate at maximum speed for 3 minutes as indicated in the protocol. For ear clips and brain tissue, use no more than 12–15 mg input material.
Blood	High hemoglobin content	• Some blood samples (e.g., horse) are rich in hemoglobin, evidenced by their dark red color. Extend lysis time by 3–5 minutes for best purity results.
	Formation of hemoglobin precipitates	• Species with high hemoglobin content (e.g., guinea pig) may accumulate insoluble hemoglobin complexes that clog the membrane. Reduce Proteinase K lysis from 5 to 3 minutes.
RNA CONTAMINATION		
Tissue	Too much input material	• DNA-rich tissues (e.g., spleen, liver and kidney) will become very viscous during lysis and may inhibit RNase A activity. Do not use more than the recommended input amount.
	Lysis time was insufficient	• Extend lysis time by 30 minutes to 3 hours after the tissue piece has completely dissolved
TISSUE DIGESTION TAKES TOO LONG		
	Tissue pieces too large	• Cut tissue pieces to the smallest possible size or grind with liquid nitrogen before starting lysis
	Tissue pieces are stuck to bottom of tube	• Vortex to release pieces from the tube bottom, and immediately after adding Proteinase K and Tissue Lysis Buffer
	Too much starting material	• Use recommended input amount
TISSUE LYSATE APPEARS TURBID		
	Formation of indigestible fibers	• Proteinase K digestion of fibrous tissues (e.g., muscle, heart, skin, ear clips), brain tissue and all RNAlater-stabilized tissues leads to the release of small indigestible protein fibers, which block the binding sites of the silica membrane. To remove fibers, centrifuge lysate at maximum speed for 3 minutes, as indicated in the protocol. For ear clips and brain tissue, use no more than 12–15 mg input material.
RATIO $A_{260}/A_{230} > 2.5$		
	Slight variations in EDTA concentration in eluates	• EDTA in elution buffer may complex with cations like Mg^{2+} and Ca^{2+} samples present in genomic DNA, which may lead to higher than usual A_{260}/A_{230} ratio. In some cases, this ratio exceeds a value of 3.0 and is consistent with highly pure samples. In these cases, the elevated value does not have any negative effect on downstream applications.

Guidelines for Choosing Sample Input Amounts

When Using the Monarch Total RNA Miniprep Kit

RNA yield, purity, and integrity vary immensely based on sample type, input amount and sample condition. Below, we have provided some empirical yield and RIN data from a wide variety of sample types, as well as guidance on the maximal input amounts for each of those samples when using the Monarch Total RNA Miniprep Kit. It is very important not to overload the column when extracting and purifying RNA, as yields, purity and integrity may suffer.

SAMPLE TYPE ⁽¹⁾		INPUT	AVERAGE YIELD (µg)	OBSERVED RIN	MAXIMUM STARTING MATERIAL
CULTURED CELLS					
HeLa		1 x 10 ⁶ cells	12–15	9–10	1 x 10 ⁷ cells
HEK 293		1 x 10 ⁶ cells	12–14	9–10	1 x 10 ⁷ cells
NIH3T3		1 x 10 ⁶ cells	8–12	9–10	1 x 10 ⁷ cells
MAMMALIAN BLOOD⁽²⁾					
Human	Fresh	200 µl	0.5–1.0	7–8	3 ml
	Frozen	200 µl	0.5–1.0	7–8	3 ml
	Stabilized	200 µl	0.5–1.0	7–8	3 ml
Rat	Frozen	100 µl	5.6	9	1 ml*
BLOOD CELLS					
PBMC (isolated from 5 ml whole blood)		5 ml	3	7	1 x 10 ⁷ cells
TISSUE					
Rat liver	Frozen pulverized	10 mg	25	8–9	20 mg
	Stabilized solid	10 mg	50–60	8–9	20 mg
Rat spleen (stabilized solid with bead homogenizer)		10 mg	40–50	9	20 mg
Rat kidney (frozen pulverized)		10 mg	7–10	9	50 mg
Rat brain	Frozen pulverized	10 mg	2–3	8–9	50 mg
	Stabilized solid	10 mg	0.5–1.5	8–9	50 mg
	Stabilized solid with bead homogenizer	10 mg	5–8	8–9	50 mg
Rat muscle (frozen pulverized)		10 mg	2–3	8–9	50 mg
Mouse muscle	Frozen pulverized	10 mg	3	8–9	50 mg
	Powder with bead homogenizer	10 mg	5	7–8	50 mg
	Stabilized solid with bead homogenizer	10 mg	8–10	9	50 mg
Mouse heart (stabilized solid w/bead homogenizer)		10 mg	5–6	8–9	50 mg
YEAST					
<i>S. cerevisiae</i>	Frozen with bead homogenizer	1 x 10 ⁷ cells	50	9–10**	5 x 10 ⁷ cells
	Fresh with Zymolyase®	1 x 10 ⁷ cells	60	9**	5 x 10 ⁷ cells
BACTERIA					
<i>E. coli</i>	Frozen	1 x 10 ⁸ cells	5	10	1 x 10 ⁸ cells
	Frozen with bead homogenizer	1 x 10 ⁸ cells	10	10	1 x 10 ⁸ cells
	Frozen with lysozyme	1 x 10 ⁸ cells	70	10	1 x 10 ⁸ cells
<i>B. cereus</i>	Frozen with lysozyme	1 x 10 ⁸ cells	20–30	9	1 x 10 ⁸ cells
	Frozen with bead homogenizer	1 x 10 ⁸ cells	8	9–10	1 x 10 ⁸ cells
PLANT					
Corn leaf (frozen pulverized with bead homogenizer)		100 mg	45	8	100 mg
Tomato leaf (frozen pulverized with bead homogenizer)		100 mg	30	8	100 mg

⁽¹⁾ RNA for other blood samples, including drosophila, zebrafish embryos/larvae, plasma, serum, saliva, buccal swabs and nucleated blood have been successfully purified with this kit; protocols are available in the product manual.

⁽²⁾ A protocol for nucleated blood (e.g., birds, reptiles) is also available.

* Mouse blood also has a maximum input of 1 ml.

** *S. cerevisiae* total RNA was run on an Agilent® Nano 6000 Chip using plant assay.



Troubleshooting Guide for Total RNA Extraction & Purification Using Monarch Kits

PROBLEM	CAUSE	SOLUTION
Clogged column	Insufficient sample disruption or homogenization	<ul style="list-style-type: none"> • Increase time of sample digestion or homogenization • Centrifuge sample after Proteinase K digestion or homogenization to pellet debris and use only supernatant for next steps • Use larger volume of DNA/RNA Protection Reagent (NEB #T2011) and/or RNA Lysis Buffer (NEB #T2012) for sample disruption and homogenization. See sample-specific protocols in the product manual.
	Too much sample	<ul style="list-style-type: none"> • Reduce amount of starting material to match kit specifications to ensure buffer amounts are sufficient and column is not overloaded. See Guidelines for Choosing Sample Input Amounts on page 344.
Low RNA yield	Incomplete elution	<ul style="list-style-type: none"> • After addition of Nuclease-free Water (NEB #B1500) to column matrix, incubate 5-10 min at room temperature and then centrifuge to elute • Perform a second elution (note: this will dilute sample)
	Sample is degraded	<ul style="list-style-type: none"> • Store input sample at -80°C prior to use • Use Monarch DNA/RNA Protection Reagent (NEB #T2011) to maintain RNA integrity during storage
	Insufficient disruption or homogenization	<ul style="list-style-type: none"> • Increase time of sample digestion or homogenization • Centrifuge sample after Proteinase K digestion or homogenization to pellet debris and use only supernatant for next steps • Use larger volume of DNA/RNA Reagent (NEB #T2011) and/or RNA Lysis Buffer (NEB #T2012) for sample disruption and homogenization. See sample specific protocol in the product manual. • For Proteinase K treated samples, doubling Proteinase K (from 5% to 10%) may lead to an increase in RNA yield
	Too much sample	<ul style="list-style-type: none"> • Reduce amount of starting material to match kit specifications to ensure buffer amounts are sufficient and column is not overloaded. See Guidelines for Choosing Sample Input Amounts on page 344.
RNA degradation	Starting material not handled/stored properly	<ul style="list-style-type: none"> • Store input sample at -80°C prior to use. Degradation of RNA may occur if sample is not flash frozen or protected by a preservation reagent. Use Monarch DNA/RNA Protection Reagent (NEB #T2011) to maintain RNA integrity during storage.
	Deviation from the stated protocol may expose RNA to unwanted RNase activities	<ul style="list-style-type: none"> • Refer to the General Guidelines for working with RNA in the product manual
	RNase contamination of eluted materials or kit buffers may have occurred	<ul style="list-style-type: none"> • See General Guidelines for working with RNA in the product manual for advice on reducing risks of contamination
Low OD ratios	Low $A_{260/280}$ values indicate residual protein in the purified sample	<ul style="list-style-type: none"> • Ensure the Proteinase K step was utilized for the recommended time. Ensure samples have no debris prior to addition of ethanol and loading onto RNA Purification Column.
	Low $A_{260/230}$ values indicate residual guanidine salts have been carried over during elution	<ul style="list-style-type: none"> • Ensure wash steps are carried out prior to eluting sample. Use care to ensure the tip of the column does not contact the flow-through after the final wash. If unsure, please repeat centrifugation. When reusing collection tubes, blot rim of tube on a Kimwipe prior to reattachment to the column to remove any residual wash buffer.
DNA contamination	Genomic DNA not removed by column	<ul style="list-style-type: none"> • Perform optional on-column DNase I treatment to remove unwanted gDNA from lysed sample • Perform in-tube/off-column DNase I treatment to remove gDNA
	Too much sample	<ul style="list-style-type: none"> • Reduce amount of starting material to match kit specifications to ensure buffer amounts are sufficient and column is not overloaded. See Guidelines for Choosing Sample Input Amounts on page 344.
Low performance of RNA in downstream steps	Salt and/or ethanol carryover has occurred	<ul style="list-style-type: none"> • Use care to ensure the tip of the RNA Purification Column does not contact the flow-through after the final wash. If unsure, please repeat centrifugation. • Be sure to spin the RNA Purification Column for 2 minutes following the final wash with RNA Wash Buffer • When reusing collection tubes, blot rim of tube on a Kimwipe prior to reattachment to the column to remove any residual wash buffer • Add additional wash step and/or extend spin time for final wash
Unusual spectrophotometric readings	RNA concentration is too low for spectrophotometric analysis	<ul style="list-style-type: none"> • For more concentrated RNA, elute with 30 µl of nuclease-free water • Increase amount of starting material (within kit specifications). See Guidelines for Choosing Sample Input Amounts on page 344.
	Silica fines in eluate	<ul style="list-style-type: none"> • Re-spin eluted samples and pipet aliquot from the top of the liquid to ensure the $A_{260/230}$ is unaffected by possible elution of silica particles

Guidelines for Choosing Sample Input Amounts

When Using the Monarch HMW DNA Extraction Kit for Cells & Blood

The table below provides data on minimum, maximum, and recommended input amounts for various cell lines and blood samples using the Monarch HMW DNA Extraction Kit for Cells & Blood. Data on yield, purity and RNA content is also provided. Samples that were successfully tested in standard ligation-based Oxford Nanopore Technologies® sequencing runs are indicated. RNA content was determined by HPLC analysis of nucleoside content after digestion of 1 µg of eluted nucleic acid with the Nucleoside Digestion Mix (NEB #M0649). Yields from blood samples vary by donor due to different leukocyte content; yield can vary up to 3-fold by donor. Similar yield and purity results were obtained with different anticoagulants (e.g., EDTA, citrate, heparin and PAXgene Blood DNA tubes were tested).

Using input amounts below the recommended minimum will reduce yields drastically. Exceeding maximum input amounts will result in DNA eluates that are highly viscous and difficult to dissolve and will reduce purity of the isolated DNA. Results are shown for samples that were lysed with agitation at 2,000 rpm.

Cells

	MINIMUM INPUT (CELLS)	MAXIMUM INPUT (CELLS)*	RECOMMENDED INPUT AMOUNT (CELLS)	YIELD (µg) FROM 1 x 10 ⁶ CELLS	PURITY RATIOS		RNA CONTENT	VALIDATED FOR ONT SEQUENCING?
					A260/280	A260/230		
HEK293	1 x 10 ⁵	1 x 10 ⁷	1 x 10 ⁶	11.5–13	1.86	2.4	≤ 1%	Yes
HeLa	1 x 10 ⁵	1 x 10 ⁷	1 x 10 ⁶	12.9	1.86	2.4	≤ 1%	Yes
NIH3T3	1 x 10 ⁵	1 x 10 ⁷	1 x 10 ⁶	9.4	1.86	2.4	≤ 1%	Yes
Jurkat	1 x 10 ⁵	1 x 10 ⁷	1 x 10 ⁶	13.7	1.86	2.5	≤ 1%	Yes
K562 (suspension cells)	1 x 10 ⁵	1 x 10 ⁷	1 x 10 ⁶	13.7	1.86	2.4	≤ 1%	Yes
HCT116	1 x 10 ⁵	1 x 10 ⁷	1 x 10 ⁶	16.9	1.86	2.5	≤ 1%	Yes
A549	1 x 10 ⁵	1 x 10 ⁷	1 x 10 ⁶	12.7	1.86	2.3	≤ 1%	Yes
U50s	1 x 10 ⁵	1 x 10 ⁷	1 x 10 ⁶	10.6	1.86	2.4	≤ 1%	Yes
HepG2	1 x 10 ⁵	1 x 10 ⁷	1 x 10 ⁶	13.4	1.81	2.2	≤ 1%	Yes
NCI-460	1 x 10 ⁵	1 x 10 ⁷	1 x 10 ⁶	9.5	1.86	2.4	≤ 1%	Yes
SK-N-SH	1 x 10 ⁵	1 x 10 ⁷	1 x 10 ⁶	9.5	1.86	2.4	≤ 1%	Yes
Aa23	1 x 10 ⁵	1 x 10 ⁷	1 x 10 ⁶	8.7	1.81	2.3	≤ 1%	Yes

Mammalian Blood

		MINIMUM INPUT (µl)	MAXIMUM INPUT (µl)	RECOMMENDED INPUT (µl)	YIELD (µg) for 500 µl**	PURITY RATIOS		RNA CONTENT	VALIDATED FOR ONT SEQUENCING?
						A260/280	A260/230		
Human***	Fresh	100	2,000	500	12–32	1.86	2.4	≤ 1%	Yes
	Frozen	100	2,000	500	9–30	1.86	2.4	≤ 1%	Yes
Mouse	Fresh	100	2,000	500	7–11	1.88	2.4	≤ 1%	Yes
	Frozen	100	2,000	500	16–17	1.88	2.4	≤ 1%	ND
Rat (fresh only)	Fresh	100	2,000	500	29–38	1.87	2.4	≤ 1%	Yes
Rabbit	Fresh	100	500	200	12–15	1.72	1.9	≤ 1%	Yes
	Fresh	100	500	200	200 µl: 4–5	1.89	2.4	≤ 1%	Yes
	Frozen	100	500	200	200 µl: 4–5	1.89	2.4	≤ 1%	Yes
Pig	Fresh	100	2,000	500	up to 42	1.86	2.4	≤ 1%	Yes
	Frozen	100	2,000	500	up to 40	1.86	2.4	≤ 1%	Yes
Horse	Fresh	100	2,000	500	16	1.86	2.3	≤ 1%	Yes
	Frozen	100	2,000	500	22.3	1.86	2.4	ND	ND
Cow	Fresh	200	2,000	500	7	1.86	2.4	≤ 1%	Yes
	Frozen	200	2,000	500	9.1	1.86	2.4	ND	ND
Rhesus monkey	Fresh	100	2,000	500	52	1.86	2.4	≤ 1%	Yes
	Frozen	100	2,000	500	52.6	1.86	2.5	ND	ND
Goat (fresh only)	Fresh	100	2,000	500	24	1.87	2.4	≤ 1%	Yes
Sheep (fresh only)	Fresh	100	2,000	500	15.3	1.87	2.4	ND	ND

Nucleated Blood

		MINIMUM INPUT (µl)	MAXIMUM INPUT (µl)	RECOMMENDED INPUT (µl)	YIELD (µg) per 5 µl	PURITY RATIOS		RNA CONTENT	VALIDATED FOR ONT SEQUENCING?
						A260/280	A260/230		
Chicken	Fresh	2	20	5	33	1.86	2.5	ND	Yes
	Frozen	2	20	5	30	1.86	2.5	ND	ND
Turkey	Fresh	2	20	5	37	1.87	2.4	ND	Yes
	Frozen	2	20	5	28	1.87	2.5	ND	ND

ND = Not determined

* For low agitation speeds, do not exceed 5 x 10⁶ cells

** Unless otherwise stated

*** Compatible with K2-EDTA, Na-citrate, Na-heparin, PAXgene® Blood DNA

Guidelines for Choosing Sample Input Amounts When Using the Monarch HMW DNA Extraction Kit for Tissue, Bacteria and Other Samples

The table below provides guidance on the minimum, maximum, and recommended input amounts for various sample types when using the Monarch HMW DNA Extraction Kit for Tissue. Data on yield, purity and RNA content is also provided. Samples that were successfully tested in standard ligation-based Oxford Nanopore Technologies sequencing runs are indicated. Using input amounts that exceed the maximum will lead to challenges in solubility and viscosity, and purity may be affected. If more starting material is required, splitting the sample and performing multiple preps is recommended. RNA content was determined by HPLC analysis of nucleoside content after digestion of 1 µg of eluted nucleic acid with the Nucleoside Digestion Mix (NEB #M0649). Using input amounts below the recommended minimums will reduce yields drastically.

		MINIMUM INPUT (mg)	MAXIMUM INPUT (mg)	RECOMMENDED INPUT (mg)	YIELD (µg) FOR RECOMMENDED INPUT (YIELD PER mg)	PURITY RATIOS		RNA CONTENT	VALIDATED FOR ONT SEQUENCING?
						A260/280	A260/230		
MAMMALIAN TISSUE									
Mouse brain	Fresh	2**	20	15	12–21	1.87	2.39	ND	YES
	Frozen	2**	20	15	15–21 (1–1.5)	1.86	2.48	ND	YES
Mouse liver	Fresh (w/NaCl)	2	15	10	7	1.84	2.10	1.2%	YES
	Frozen (w/NaCl)	2	15	10	17–19 (1.7–1.9)	1.89	2.50	ND	YES
	Fresh*	2	15	10	20	1.84	1.52*	8.7%	YES
	Frozen*	2	15	10	27–31 (2.7–3.1)	1.89	1.93**	ND	YES
Mouse muscle	Fresh	2**	25	20	8–9	1.87	2.25	2.1%	YES
	Frozen	2**	25	20	12–16 (0.6–0.8)	1.87	2.30	ND	YES
Mouse kidney	Fresh	2	15	10	23–34	1.86	2.44	ND	YES
	Frozen	2	15	10	32–41 (3.2–4.1)	1.86	2.53	0.8%	YES
Mouse tail	Frozen	2**	25	20	20 (1.8–2.1)	1.86	2.43	ND	YES***
Mouse ear punch	Fresh	2**	15	10	15–16 (1.5–1.6)	1.86	2.29	ND	YES
Rat kidney	Frozen	2	15	10	20–25	1.87	2.40	ND	YES
BACTERIA									
<i>E. coli</i> (Gram-negative)	Frozen	5 x 10 ⁸ cells	5 x 10 ⁸ cells	1 x 10 ⁹ cells	8–9	1.89	2.31	1.7%	YES
<i>B. cereus</i> (Gram-positive)	Frozen	2 x 10 ⁸ cells	4 x 10 ⁸ cells	2 x 10 ⁸ cells	4–5	1.86	2.20	3.9%	YES
<i>M. luteus</i> (Gram-positive)	Frozen	ND	ND	1 x 10 ⁸ cells	2.0	1.89	2.09	ND	ND
AMPHIBIAN									
<i>X. laevis</i>	Fresh	ND	ND	3–4	5	1.86	2.51	2.3%	ND
YEAST									
<i>S. cerevisiae</i>	Fresh	ND	ND	20 x 10 ⁷ cells	3–6***	1.90	2.01	ND	ND
INSECT									
<i>A. aegypti</i>	Frozen	ND	ND	15	6	1.84	2.53**	2.7%	ND
NEMATODE									
<i>C. elegans</i> ****	Frozen	ND	ND	2 plates	8.2	1.91	2.5	ND	ND

ND = Not determined

* Standard protocol without recommended NaCl treatment.

** If working with input amounts < 5 mg, refer to the product manual for guidance on reducing buffer volumes.

*** Total nucleic acid yields are 4–10 µg and 6–12 µg for haploid and diploid strains, respectively. Though an RNase A step is included, RNA is co-purified. Yields may vary depending on the strain.

**** Rotor-stator homogenization is recommended.

+ Measured with Nanodrop One; systems that differentiate turbidity in the content profiling will give higher values.

++ Measured with Unchained Labs Lunatic (formerly Trinean DropSense16); devices without content profiling that differentiates turbidity may give lower values.

+++ Size selection is recommended.

Genetic Markers

A *genotype* indicates the genetic state of the DNA in an organism. It is a theoretical construct describing a genetic situation that explains the observed properties (phenotype, see below) of a strain. *E. coli* genotypes list only genes that are defective (1). If a gene is not mentioned, then it is not known to be mutated*. **. Prophages and plasmids that were present in the original K-12 strain (F, λ , e14, *rac*) are normally listed only if absent. However, for simplicity, we have not listed λ except when it is present, and we have listed F and its variants in all cases. Parentheses or brackets surround a prophage or plasmid when listed. Genes are given three-letter, lower-case, italicized names (e.g., *dam*) that are intended to be mnemonics suggesting the function of the gene (here, DNA adenine methylase). If the same function is affected by several genes, the different genes are distinguished with uppercase italic letters (e.g., *recA*, *recB*, *recC*, *recD* all affect recombination). Proper notation omits superscript + or – in a genotype, but these are sometimes used redundantly for clarity, as with F' *lac*⁺ *proA*⁺ *B*⁺. Deletion mutations are noted as Δ , followed by the names of deleted genes in parentheses, [e.g., Δ (*lac-pro*)]. All genes between the named genes are also deleted. Specific mutations are given allele numbers that are usually italic arabic numerals (e.g., *hsdR17*) and may be characterized as *am*=amber (UAG) mutation or *ts*=inactive at high temperature, as appropriate. Some common alleles [e.g., Δ (*lac-pro*)*X111*] break the rules. If two strains' genotypes list a gene with the same allele number, they should carry exactly the same mutation.

The *phenotype* of a strain is an observable behavior, e.g., *Lac*[–] fails to grow on lactose as a sole carbon source. Phenotypes are capitalized and in Roman type, and the letters are always followed by superscript + or – (or sometimes r, resistant, or s, sensitive). Although phenotypes do not, strictly speaking, belong in a genotype, they are sometimes included following the genotype designation when the former is not obvious from the latter [e.g., *rpsL104* (Str^r)—gene name from ribosomal protein, small subunit, S12, confers resistance to streptomycin].

Some common genes of interest are described below and on the next page; a catalogue of genetically defined genes can be found in reference 2 and on the very useful internet site maintained by the *E. coli* Genetic Stock Center (CGSC) at Yale University <<http://cgsc.biology.yale.edu/>>. Additional information from CGSC can be obtained from curator Mary Berlyn by e-mail <cgsc@yale.edu>.

* Most *E. coli* laboratory strains have been heavily mutagenized over forty years of study, and different lines may carry different, so far undiscovered, mutations that may or may not affect your situation. For this reason, it is sometimes useful to try more than one line or strain background in your experiments.

** *E. coli* B and its derivatives are naturally Lon[–] and Dcm[–]. We have listed this in brackets even though it is the wild type state for these strains.

References

- (1) Demerec et al. (1966) *Genetics*, 54, 61–76.
- (2) Berlyn, M.K.B. (1996). In F. C. Niedhardt et al. (Ed.), *Escherichia coli and Salmonella: cellular and molecular biology*, (2nd ed.), Vol. 2, (pp. 1715–1902). ASM Press.
- (3) Raleigh, E.A. et al. (1991) *J. Bacteriol.*, 173, 2707–2709.

dam Endogenous adenine methylation at GATC sequences is abolished. *dam* strains have a high recombination frequency, express DNA repair functions constitutively, and are poorly transformed by Dam-modified plasmids. Used for making DNA susceptible to cleavage by some restriction enzymes (e.g., BclI).

dcm Endogenous cytosine methylation at CCWGG sequences is abolished. Used for making DNA susceptible to cleavage by some restriction enzymes (e.g., AclI).

dnaJ One of several "chaperonins" is inactive. This defect has been shown to stabilize certain mutant proteins expressed in *E. coli*.

dut dUTPase activity is abolished. This mutation, in combination with *ung*, allows incorporation of uracil into DNA. Used for oligonucleotide mutagenesis.

endA Activity of nonspecific Endonuclease I is abolished. DNA preparations are thought to be of higher quality when prepared from *endA* strains.

e14 An excisable prophage-like element, present in K-12 but missing from many derivatives. e14 carries the *mcrA* gene among others, therefore e14[–] strains are McrA[–].

F A low-copy number self-transmissible plasmid. F' factors carry portions of the *E. coli* chromosome, most notably the *lac* operon and *proAB* on F' *lac-proA*⁺ *B*⁺.

fhuA An iron uptake receptor is mutated. This mutation confers resistance to phage T1 (ferric hydroxamate uptake). Former name is *tonA*.

gal The ability to metabolize galactose is abolished.

glnV See *supE*.

gyrA A point mutation in DNA gyrase, subunit A. This mutation confers resistance to the antibiotic nalidixic acid.

hflA This mutation results in high frequency lysogenization by λ .

hsdR, **hsdS** DNA that does not contain methylation of certain sequences is recognized as foreign by EcoKI or EcoBI and restricted (degraded). These enzymes recognize different sequences and are encoded by different alleles of *hsdRMS*. *hsdR* mutations abolish restriction but not protective methylation (*r*[–]*m*⁺), while *hsdS* mutations abolish both (*r*[–]*m*[–]). DNA made in the latter will be restricted when introduced into a wild-type strain.

lacIq The *lac* repressor is overproduced, turning off expression from *P_{lac}* more completely.

lacZ β -galactosidase activity is abolished.

lacZ::T7gene 1 The phage T7 RNA polymerase (= gene 1) is inserted into the *lacZ* gene.

lacY Lactose permease activity is abolished.

Δ (*lac*) = deletion; there are four common deletions involving *lac*:

Δ (*lacZ*)M15 expresses a fragment that complements the *lac* α -fragment encoded by many vectors. These vectors will yield blue color on X-Gal only if the host carries Δ M15.

Δ U169, Δ X111, and Δ X74 all delete the entire *lac* operon from the chromosome, in addition to varying amounts of flanking DNA. Δ X111 deletes *proAB* as well, so that the cell requires proline for growth on minimal medium, unless it also carries F' *lac proA*⁺ *B*⁺.

lon Activity of a protease responsible for degrading aberrant proteins is abolished. Some eukaryotic proteins are stabilized in *lon* strains. *E. coli* B naturally lacks Lon.

Genetic Markers (continued)

lysY	The lysozyme gene from the T7 bacteriophage is mutated. The mutation K128Y eliminates lysozyme activity, but the mutant protein still binds to and inhibits T7 RNA polymerase.	recD	Exonuclease activity of ExoV is abolished, but recombination activity is elevated. Inverted repeat sequences in λ can be propagated in <i>recD</i> strains. Plasmid replication is aberrant.	supE	A glutamine-inserting amber (UAG) suppressor tRNA; required for growth of some phage vectors. Now called <i>glnV</i> .
malB	The <i>malB</i> region encompasses the genes <i>malEFG</i> and <i>malK lamb malM</i> . $\Delta(malB)$ deletes most or all of this region and eliminates expression of Maltose Binding Protein (MalE).	recF	Plasmid-by-plasmid homologous recombination is abolished.	supF	A tyrosine-inserting amber (UAG) suppressor tRNA; required for lytic growth of S7 or S100 λ phage, such as λ gt11. Now called <i>tyrT</i> .
mcrA, mcrBC	A restriction system that requires methyl cytosine is abolished. DNA containing methylcytosine in some sequences is restricted by Mcr ⁺ . <i>dcm</i> -modified DNA is not restricted by Mcr ⁺ . $\Delta(mcrC-mrr)$ deletes six genes: <i>mcrC-mcrB-hsdS-hsdM-hsdR-mrr</i> ; <i>mcrA</i> is lost with <i>e14</i> .	recJ	Plasmid-by-plasmid homologous recombination is abolished.	thi-1	The ability to synthesize thiamine is abolished (vitamin B1).
mrr	A restriction system that requires cytosine or adenine methylation is abolished; however, <i>dam</i> ⁺ , <i>dcm</i> ⁺ or EcoKI-modified DNA is not restricted by Mrr ⁺ . The methylcytosine-dependent activity is also known as McrF (3).	relA1	Lacks ppGpp synthesis during the stringent response to amino acid starvation; activity of ATP:GTP 3'-pyrophosphotransferase (EC2.7.6.5) is abolished.	traD	The self-transmissibility of the F factor is severely reduced.
mtl	The ability to metabolize the sugar alcohol mannitol is abolished.	rtbD	Lacks functional TDP-rhamnose synthetase, and thus does not synthesize the cell surface O-antigen.	tsp	A periplasmic protease that may degrade secreted or cytoplasmically overexpressed proteins after lysis is abolished. Now called <i>prc</i> .
ompT	Activity of outer membrane protease (protease VII) is abolished.	rpoH	(also known as <i>htrP</i>) Lack of this heat-shock transcription factor abolishes expression of some stress-induced protease activities in addition to <i>lon</i> . Some cloned proteins are more stable in <i>rpoHam supCts</i> strains at high temperature.	tsx	Confers resistance to bacteriophage T6.
phoA	Activity of alkaline phosphatase is abolished.	sbcB	Exo I activity is abolished. Strains carrying <i>recB recC</i> and <i>sbcB</i> are usually also <i>sbcC</i> . These quadruple mutant strains are recombination-proficient and propagate inverted repeats in λ , but plasmid replication is aberrant.	tyrT	See <i>supC</i> , <i>supF</i> .
prc	See <i>tsp</i> .	sbcC	Usually found with <i>recB recC sbcB</i> . However, strains carrying <i>sbcC</i> alone are recombination-proficient and stably propagate inverted repeats both in λ and in plasmids.	ung	Uracil N-glycosylase activity is abolished. Uracil incorporated into DNA is removed by Ung ⁺ , leaving baseless site. See <i>dut</i> .
recA	Homologous recombination is abolished; particularly desirable when working with sequences containing direct repeats > 50 bp.	sulA	Mutations in this gene allows cells to divide and recover from DNA damage in a <i>lon</i> mutant background (<i>suppressor of Lon</i>).	xyl	The ability to metabolize the sugar xylose is abolished.
recB, recC	Exonuclease and recombination activity of Exonuclease V is abolished. Homologous recombination is much reduced in <i>recB recC</i> strains that are not also <i>sbcB</i> or <i>sbcA</i> . Stability of inverted repeat sequences is enhanced in <i>recB recC</i> strains, especially if they are also <i>sbcB sbcC</i> . Plasmid replication may be aberrant.	supC(ts)	A thermosensitive tyrosine-inserting ochre (UAA) and amber (UAG) suppressor tRNA. Nonsense mutations in the same strain are suppressed only at low temperatures. Now called <i>tyrT</i> .	(P1)	The cell carries a P1 prophage. Such strains express the P1 restriction system.
				(P2)	The cell carries a P2 prophage. This allows selection against Red ⁺ Gam ⁺ λ (Spi ⁻ selection).
				(ϕ80)	The cell carries the lambdoid prophage ϕ 80. A defective ϕ 80 prophage carrying the <i>lac</i> M15 deletion is present in some strains.
				(Mu)	Mu prophage; Mud means the phage is defective.



Carole has been with NEB for over 33 years and serves as Director of International Business. Carole works closely with our subsidiary offices and distributors to make sure that our global network runs smoothly. In past years, Carole spent much of her time traveling to various locations around the world – now she is spending more time in her home office, which her dog Bonnie much appreciates!

Enhancing Transformation Efficiency

Transformation efficiency is defined as the number of colony forming units (cfu) that would be produced by transforming 1 µg of plasmid into a given volume of competent cells. However, 1 µg of plasmid is rarely transformed. Instead, efficiency is routinely calculated by transforming 100 pg–1 ng of highly purified supercoiled plasmid under ideal conditions. Transformation Efficiency (TE) is calculated as: $TE = \text{Colonies}/\mu\text{g}/\text{Dilution}$. Efficiency calculations can be used to compare cells or ligations. Our recommended protocols and tips are presented here to help you achieve maximum results.

Recommended Protocols

High Efficiency Transformation Protocol

1. Thaw cells on ice for 10 minutes
2. Add 1 pg–100 ng of plasmid DNA (1–5 µl) to cells and mix without vortexing
3. Place on ice for 30 minutes
4. Heat shock at 42°C for 10–30 seconds or according to recommendations. For BL21, use exactly 10 seconds.
5. Place on ice for 5 minutes
6. Add 950 µl of room temperature SOC or NEB 10-beta/Stable Outgrowth Medium.
7. Place at 37°C for 60 minutes. Shake vigorously (250 rpm) or rotate.
8. Mix cells without vortexing and perform several 10-fold serial dilutions in SOC or NEB 10-beta/Stable Outgrowth Medium.
9. Spread 50–100 µl of each dilution onto pre-warmed selection plates and incubate overnight at 37°C (30°C for SHuffle® strains) or according to recommendations

5 Minute Transformation Protocol

(10% efficiency compared to above protocol)

1. Thaw cells in your hand
2. Add 1 pg–100 ng of plasmid DNA (1–5 µl) to cells and mix without vortexing
3. Place on ice for 2 minutes
4. Heat shock at 42°C for 30 seconds or according to recommendations.
5. Place on ice for 2 minutes
6. Add 950 µl of room temperature SOC or NEB 10-beta/Stable Outgrowth Medium. Immediately spread 50–100 µl onto a selection plate and incubate overnight at 37–42°C. (30°C for SHuffle strains) NOTE: Selection using antibiotics other than ampicillin may require some outgrowth prior to plating.

Transformation Tips

Thawing

- Cells are best thawed on ice
- DNA should be added as soon as the last trace of ice in the tube disappears
- Cells can be thawed by hand, but warming above 0°C decreases efficiency

Incubation of DNA with Cells on Ice

- Incubate on ice for 30 minutes. Expect a 2-fold loss in TE for every 10 minutes this step is shortened.

Heat Shock

- Both temperature and time are specific to the transformation volume and vessel. Typically, 30 seconds at 42°C is recommended.

Outgrowth

- Outgrowth at 37°C for 1 hour is best for cell recovery and for expression of antibiotic resistance. Expect a 2-fold loss in TE for every 15 minutes this step is shortened.
- SOC gives 2-fold higher TE than LB medium
- Incubation with shaking or rotation results in 2-fold higher TE

Plating

- Selection plates can be used warm or cold, wet or dry with no significant effects on TE
- Warm, dry plates are easier to spread and allow for the most rapid colony formation

DNA

- DNA should be purified and resuspended in water or TE Buffer
- Up to 10 µl of DNA from a ligation mix can be used with only a 2-fold loss of efficiency
- Purification by either a spin column or phenol/chloroform extraction and ethanol precipitation is ideal
- The optimal amount of DNA is lower than commonly recognized. Using clean, supercoiled pUC19, the efficiency of transformation is highest in the 100 pg–1 ng range. However, the total colonies which can be obtained from a single transformation reaction increase up to about 100 ng.

DNA Contaminants to Avoid

CONTAMINANT	REMOVAL METHOD
Detergents	Ethanol precipitate
Phenol	Extract with chloroform and ethanol precipitate
Ethanol or Isopropanol	Dry pellet before resuspending
PEG	Column purify or phenol/chloroform extract and ethanol precipitate
DNA binding proteins (e.g., ligase)	Column purify or phenol/ chloroform extract and ethanol precipitate

Electroporation Tips

NEB 5-alpha (NEB #C2989) and NEB 10-beta (NEB #C3020) Competent *E. coli* Strains are available as electrocompetent cells. The following tips will help maximize transformation efficiencies.

- Pre-chill electroporation cuvettes and microcentrifuge tubes on ice
- Thaw cells on ice and suspended well by carefully flicking the tubes
- Once DNA is added, electroporation can be carried out immediately. It is not necessary to incubate DNA with cells. The maximum recommended volume of a DNA solution to be added is 2.5 µl. Addition of a large volume of DNA decreases transformation efficiency.
- DNA should be purified and suspended in water or TE. Transformation efficiency is > 10-fold lower for ligation mixtures than the control pUC19 plasmid due to the presence of ligase and salts. If used directly, ligation reactions should be heat-inactivated at 65°C for 20 min and then diluted 10-fold. For optimal results, spin columns are recommended for clean up of ligation reactions.
- Electroporation conditions vary with different cuvettes and electroporators. If you are using electroporators not specified in the protocol, you may need to optimize the electroporation conditions. Cuvettes with 1 mm gap are recommended (e.g., BTX Model 610/613 and Bio-Rad #165-2089). Higher voltage is required for cuvettes with 2 mm gap.
- Arcing may occur due to high concentration of salts or air bubbles
- It is essential to add recovery medium to the cells immediately after electroporation. One minute delay can cause a 3-fold reduction in efficiency.
- Cold and dry selection plates lead to lower transformation efficiency. Pre-warm plates at 37°C for 1 hour. Using 37°C pre-warmed recovery medium increases the efficiency by about 20%.
- Refreeze unused cells in a dry ice/ethanol bath for 5 min and then store at -80°C. Do not use liquid nitrogen. Additional freeze-thaw cycles result in lower transformation efficiency.



Find tips
for successful
transformation.

Protein Expression with T7 Express Strains

T7 Protein Expression

1. Transform expression plasmid into a T7 expression strain. Plate out on antibiotic selection plates and incubate overnight at 37°C (24 hours at 30°C for SHuffle strains).
2. Resuspend a single colony in 10 ml liquid culture with antibiotic
3. Incubate at 37°C until OD₆₀₀ reaches 0.4–0.6
4. Induce with 40 µl of a 100 mM stock of IPTG (final conc. = 0.4 mM) and induce for 2 hours at 37°C (4 hours at 30°C or 16°C overnight for SHuffle strains)
5. Check expression by Coomassie stained protein gel, Western Blot or activity assay. Check expression in the total cell extract (soluble + insoluble) and the soluble fraction alone.
6. For large scale, inoculate 1 L of liquid medium (with antibiotic) with a freshly grown colony or 10 ml of freshly grown culture. Incubate at 37°C (30°C for SHuffle strains) until OD₆₀₀ reaches 0.4–0.6. Add IPTG to 0.4 mM. Induce 2 hours at 37°C or 15°C overnight (4 hours at 30°C or 16°C overnight for SHuffle strains).

Troubleshooting Tips

No Colonies or No Growth in Liquid Culture

- Even though T7 expression is tightly regulated, there may be a low level of basal expression in the T7 Express host. If toxicity of the expressed protein is likely, transformation of the expression plasmid should be carried out in a more tightly controlled expression strain:
 - In *l⁻* strains over-expression of the *Lac^l* repressor reduces basal expression of the T7 RNA polymerase
 - In *lysY* strains, mutant T7 lysozyme is produced which binds to T7 RNA polymerase, reducing basal expression of the target protein. Upon induction, newly made T7 RNA polymerase titrates out the lysozyme and results in expression of the target protein.
- Incubation at 30°C or room temperature may also alleviate toxicity issues
- Check antibiotic concentration (test with control plasmid)

No Protein Visible on Gel or No Activity

- Check for toxicity – the cells may have eliminated or deleted elements in the expression plasmid. If this is the case, test *l⁻* and/or *lysY* strains to reduce basal expression.
- Culture cells for protein induction. Just before induction, plate a sample on duplicate plates with and without antibiotic selection. If toxicity is an issue, significantly fewer colonies will be seen on plates containing antibiotic (indicating that the plasmid has been lost) compared to plates without antibiotic.

Induced Protein is Insoluble

T7 expression often leads to very high production of protein that can result in the target protein becoming insoluble. In this case:

- Induce at lower temperatures (12–15°C overnight)
- Reduce IPTG concentration to 0.01 mM – 0.1 mM
- Induce for less time (as little as 15 minutes)
- Induce earlier in growth (OD₆₀₀ = 0.3 or 0.4)



Ron has worked at NEB for 13 years as an Automation Engineer. He is responsible for the implementation of many of the processes needed to package our products. In his spare time, Ron enjoys the outdoors and riding his motorcycle.

Troubleshooting Guide for NEBNext DNA Library Prep Kits

Troubleshooting guides are available at www.neb.com for NEBNext products including NEBNext RNA Library Prep, NEBNext Ultra II FS DNA Library Prep, NEBNext DNA Library Prep, NEBNext rRNA Depletion Kit (Bacteria) and NEBNext Custom RNA Depletion Design Tool with NEBNext RNA Depletion Core Reagent Set

PROBLEM	CAUSE	SOLUTION
Failed Library Prep. For example: <ul style="list-style-type: none"> You may see nothing on the Bioanalyzer, or similar instrument After amplification, you may see library fragments that are still the same size as the starting input DNA rather than ~120 bp longer than the input DNA 	Input DNA contains an inhibitor	<ul style="list-style-type: none"> Ensure DNA does not contain inhibitor Consider additional cleanup step
	Failed step - Any of the enzymatic steps can fail if a critical reagent is omitted, or if the reagent has become inactive	<ul style="list-style-type: none"> Confirm reagents were added for each step in the protocol
Low Library Yield	Input DNA is damaged	<ul style="list-style-type: none"> Shear input DNA on a Covaris® instrument in 1X TE Buffer, and/ or use the NEBNext FFPE DNA Repair Mix (NEB #M6630) after shearing and prior to library prep
	Adaptor is denatured	<ul style="list-style-type: none"> When diluting NEBNext adaptors, use 10 mM Tris HCl (pH 7.5-8.0) with 10 mM NaCl Keep the adaptor on ice until use
	Insufficient mixing	<ul style="list-style-type: none"> Mix samples well with 80-90% of the total volume in the well or tube by pipetting up and down. Keep the tip in the liquid to avoid the formation of bubbles. For enzymatic steps, follow the manual recommendations (usually 10 mix cycles) Try to avoid losing sample in the pipette tip or on the source tube during transfer
	SPRI beads have dried out before elution	<ul style="list-style-type: none"> Add Elution Buffer and mix before the beads turn lighter brown and start cracking For additional tips about SPRI beads view our video
	Incomplete ethanol removal during SPRI bead wash	<ul style="list-style-type: none"> Quickly spin the tube after the last ethanol wash at each SPRI bead step, keep the tube on the magnet and remove residual ethanol with a p10 tip For additional tips about SPRI beads view our video
	SPRI bead sample loss	<ul style="list-style-type: none"> Mix slowly to avoid droplets clinging to the inside of the tip, which may not combine with the sample before the tip is ejected. Dispense the last mix slowly into the sample tube so that the liquid stays together. Wait 1 second before pushing the pipette to the second stop. When removing the supernatant, take care not to remove any beads. Check your tip over a white piece of paper. If beads are visible, dispense everything back into the tube and allow beads to resettle. For additional tips about SPRI beads view our video
	Sample storage after A-tailing	<ul style="list-style-type: none"> Avoid prolonged storage of sample before moving to ligation. If sample inputs are low, avoid overnight storage and move immediately from end prep to adaptor ligation.
	Adaptor self-ligation (Adaptor dimer formation)	<ul style="list-style-type: none"> Do not add adaptor to the ligation master mix. This can cause increased adaptor dimer formation. For best results, add the adaptor to the sample, mix and then add ligase master mix and ligation enhancer
Adaptor Dimer Formation (sharp 127 bp peak on Bioanalyzer)	Ligation incubation temperature is too warm	<ul style="list-style-type: none"> If ligation incubation occurs above 20°C the DNA ends may breathe, which could reduce ligation
		<ul style="list-style-type: none"> To recover the samples, repeat the bead cleanup using a 0.9 x bead ratio.
	Adaptor concentration too high	<ul style="list-style-type: none"> Optimize adaptor dilution based on your sample input, quality and type using an adaptor titration experiment Adaptor titration may need to be repeated if the source of the sample input changes (e.g., extraction method, tissue type, etc.)
Adaptor or primers remaining after PCR (e.g., visible on Bioanalyzer or similar instrument after PCR)	Adaptor self-ligation (Adaptor Dimer formation)	<ul style="list-style-type: none"> Do not add adaptor to ligation master mix. This can cause increased adaptor dimer formation. For best results, add adaptor to sample, mix and then add ligase master mix and ligation enhancer. Mix again.
	Excess adaptor or primer used or inefficient cleanup	<ul style="list-style-type: none"> Perform another 0.9 x SPRI cleanup
Overamplification (Once PCR primers are depleted, library fragments will become single stranded and/ or form heteroduplexes. These appear as high molecular weight fragments on a Bioanalyzer or similar instrument)	Too many PCR Cycles	<ul style="list-style-type: none"> We recommend starting with the number of PCR cycles recommended in the product manual. The ideal number of PCR cycles for your samples may vary. Reduce the number of PCR cycles if you are seeing overamplification Data quality may be compromised if overamplified libraries are sequenced
	Not enough PCR primer	<ul style="list-style-type: none"> Check primer concentration and ensure that you are adding the primer volume recommended in the manual Store primers at the correct temperature to prevent degradation Data quality may be compromised overamplified libraries are sequenced
	Too much input DNA	<ul style="list-style-type: none"> The higher the input of template for the PCR, the sooner the primers will be depleted NEBNext adaptor and PCR primers require a minimum of 3 PCR cycles. If you cannot further reduce the number of PCR cycles, consider a size selection step, or using only a fraction of the ligated library as input for PCR.

Guidelines for NGS Library Prep

DNA SAMPLE INPUT GUIDELINES

Integrity of DNA

- The quality of the input material directly affects the quality of the library. Absorbance measurements can be used as an indication of DNA purity. Ideally, the ratio of the absorbance at 260 nm to 280 nm should be between 1.8 – 2.0. However, measurements can be affected by the presence of RNA or small nucleic acid fragments. A DNA Integrity Number can be determined using the Agilent TapeStation® or similar instrumentation, and qPCR-based methods can also provide a measurement of DNA integrity.

Quantitation of DNA

- It is important to quantify accurately the DNA sample prior to library construction. Fluorescence-based detection which utilizes dsDNA-specific dyes, such as the Qubit® from Thermo Fisher Scientific®, Qubit is more accurate than UV spectrometer-based measurements, as the presence of RNA or other contaminants can result in overestimation of the amount of the DNA sample by the latter. This can result in use of non-optimal adaptor dilutions and numbers of PCR cycles, compromising library prep efficiency.

RNA SAMPLE INPUT GUIDELINES

Integrity of RNA

- We recommend determining the RNA sample input using the RNA Integrity Number (RIN) estimated by the Agilent TapeStation or similar instrumentation. Ideally the RNA sample will have a RIN value of 7 or higher but NEBNext RNA products are compatible for use even with samples with low RIN values.
- RNA should be completely free of DNA, and DNA digestion of the purified RNA using RNase-free DNase I (such as that provided with the Monarch Total RNA Miniprep Kit) is recommended

Quantitation of RNA

- It is important to quantify accurately the RNA sample prior to library construction. The concentration can be estimated with the Agilent Bioanalyzer or similar instrumentation, using a pico or nano chip. Alternatively, RNA concentration can be determined by measuring the absorbance at 260 nm (A_{260}) in a spectrophotometer such as a NanoDrop®. Note that free nucleotides or organic compounds used in some RNA extraction methods also absorb UV light near 260 nm and will cause an over-estimation of RNA concentration.

LIBRARY QUANTITATION

- Library quantity and quality can be determined using the Agilent TapeStation or similar instrumentation, or qPCR-based methods.

The NEBNext Library Quant Kit for Illumina is an ideal method for highly accurate library quantitation, that is especially important for PCR-free workflows.

BEAD-BASED CLEAN-UPS AND SIZE SELECTION

- Be sure to vortex the beads well just before use. They should form a uniform suspension. If beads have settled for a long time period without being agitated a tight bead sediment can form. When beads have not been used for several weeks, plan for extra time for bead vortexing and agitation.
- Do not over-dry the beads. Beads should still be dark brown and glossy looking when eluting. Over-drying can make resuspension difficult and reduce yield.
- Take care not to remove beads after separation. If beads are accidentally aspirated, dispense everything back, allow the beads to settle and then try again.
- Remove all of the supernatant after the bind step. After removing most of the liquid with a p200 pipette, aspirate any remaining drops with a p10 pipette if necessary. Incomplete supernatant removal can cause leftover adaptor dimer or PCR primers to remain in the libraries.
- Bead ratios for cleanup and size selection after NEBNext ligation steps are appropriate just for the ligation step. Different ratios would apply if size selection is done after PCR or at any other step in the workflow.
- When adding beads to the sample, aspirate slowly to make sure the correct volume of beads is drawn into the tip. Remove any droplets of beads from the outside of the tip and make sure you dispense the full volume into the sample.

INDICES/BARCODES

- When using a subset of the indices supplied in a kit, or using indices from more than one kit, it is important to optimize the combination of indices used, in order to ensure balanced sequencing reads. We provide recommendations for NEBNext index combinations in the manuals for NEBNext Oligos products.
- For index primers provided in vials, open only one vial at a time, to minimize the risk of contamination
- Be sure to change pipette tips for each index primer
- For 96-well plate formats, NEBNext index primers are provided in single-use plates with pierceable foil lids. To avoid risk of contamination, do not pipette from a well more than one time.

Labeling with SNAP-tag® Technology-Troubleshooting Guide

APPLICATION	PROBLEM	CAUSE	SOLUTION
Cellular Labeling	No labeling	Fusion protein not expressed	<ul style="list-style-type: none"> • Verify transfection • Check expression of fusion protein via Western blot or SDS-PAGE with an appropriate fluorescent substrate
	Weak labeling	Poor expression and/or insufficient exposure of fusion protein to substrate	<ul style="list-style-type: none"> • Increase substrate concentration • Increase incubation time
		Rapid turnover of fusion protein	<ul style="list-style-type: none"> • Analyze samples immediately or fix cells directly after labeling • Label at lower temperature (4°C or 16°C)
	High background	Non-specific binding of substrates	<ul style="list-style-type: none"> • Reduce substrate concentration and/or incubation time • Allow final wash step to proceed for up to 2 hours • Include fetal calf serum or BSA during labeling
	Signal strongly reduced after short time	Instability of fusion protein	<ul style="list-style-type: none"> • Fix cells • Switch tag from N-terminus to C-terminus or vice versa
		Photobleaching	<ul style="list-style-type: none"> • Add commercially available anti-fade reagent • Reduce illumination time and/or intensity
Labeling in Solution	Precipitation	Insoluble fusion	<ul style="list-style-type: none"> • Test from pH 5.0 to 10.0 • Optimize salt concentration [50 to 250 mM] • Add 0.05 to 0.1% Tween 20
	Weak or no labeling	Exhaustive labeling has not been achieved	<ul style="list-style-type: none"> • Increase incubation time to 2 hrs at 25°C or 24 hrs at 4°C • Reduce the volume of protein solution labeled • Check expression of fusion protein via SDS-PAGE with an appropriate fluorescent substrate
	Loss of activity	Instability of fusion protein	<ul style="list-style-type: none"> • Reduce labeling time • Decrease labeling temperature (4°C or 16°C)

Cellular Imaging & Analysis FAQs

Q. How does SNAP-tag® labeling differ from using GFP fusion proteins?

A. GFP and SNAP-tag are both valuable technologies used to visualize proteins in live cells. GFP is an intrinsically fluorescent protein derived from *Aequorea victoria* while SNAP-tag is derived from hAGT, a human DNA repair protein. In contrast to GFP, the fluorescence of SNAP-tag fusions can be readily turned on with the addition of a variety of fluorescent probes added directly to the culture media. Substituting different fluorophores or other functionalities (biotin, magnetic beads, blocking agents) requires no new cloning or expression, merely incubation of the appropriate substrate with cells, cell lysates or recombinant proteins.

Q. What is the difference between SNAP- and CLIP-tag™?

A. SNAP-tag and CLIP-tag are both derived from O⁶-alkylguanine-DNA-alkyltransferase (hAGT). SNAP-tag recognizes O⁶-labeled benzylguanine substrates while CLIP-tag recognizes O²-labeled benzylcytosine substrates. Each tag transfers the label from the substrate to itself, resulting in specific covalent labeling. In creating the tags, hAGT has been engineered to no longer interact with DNA, but rather with derivatives of the free benzylguanine or benzylcytosine substrates. The tags exhibit no cross-reactivity with one another, enabling researchers to simultaneously label fusion proteins containing SNAP- and CLIP-tags with different fluorophores in live cells.

Q. Can I clone my protein as a fusion to the N- or C-terminus of the tags?

A. Yes. SNAP- and CLIP-tags can be fused to either the N- or C-terminus of a protein of interest. However, to label surface proteins on the outside of cells, the SNAP-tag or CLIP-tag must be cloned so that it is oriented to the extracellular surface of the plasma membrane. In this orientation, the tag is accessible to its fluorophore conjugated substrate.

Q. Are the substrates toxic to cells?

A. No toxicity has been noted by proliferation or viability assays when using up to 20 µM substrate for 2 hours. Most of the substrates can be incubated with cells for 24 hours up to a concentration of 20 µM without significant toxicity.

Q. How stable is the labeled protein in mammalian cells?

A. The stability of the tagged protein in the cell is dependent upon the stability of protein of interest. Labeled SNAP-tag fusion protein has been detected for up to 2 days in mammalian cells.

Q. Are SNAP-tag substrates stable to fixation?

A. Yes. SNAP-tag substrates are derived from organic fluorophores which are stable to fixation. Fluorescently-labeled SNAP-tag fusion proteins do not lose signal intensity in contrast to some GFP spectral variants. After labeling the SNAP-tag fusion proteins, the cells can be fixed with standard fixation methods such as para-formaldehyde, ethanol, methanol, methanol/acetone etc. without loss of signal.

Q. What conditions are recommended for SNAP-tag labeling in vitro?

A. The SNAP-tag labeling reaction is tolerant of a wide range of buffers. The requirements of the fusion partner should dictate the buffer selected. The following buffer guidelines are recommended: pH between 5.0 and 10.0, monovalent salts (e.g. sodium chloride) between 50 mM and 250 mM and at least 1 mM DTT. Non-ionic detergents can be added to 0.5% v/v if required, but SDS and other ionic detergents should be avoided entirely because they inhibit the activity of the SNAP-tag. Metal chelating reagents (e.g., EDTA and EGTA) also inhibit SNAP-tag activity and should be avoided.

Lambda

48,502 base pairs

GenBank Accession #: NC_001416

See page 128 for ordering information.

There are no restriction sites for the following enzymes: *AbsI*(x), *AsiSI*, *FseI*, *MauBI*(x), *MreI*(x), *NottI*, *PacI*, *SfiI*, *SpeI*, *SrfI*, *SwaI*

(x) = enzyme not available from NEB




We recommend **NEBcutter** at **NEBcutter.com** to identify the restriction sites within your DNA sequence. **NEBcutter** indicates cut frequency and methylation-state sensitivity.

Lambda (λ) is a large, temperate *E. coli* bacteriophage with a linear, largely double-stranded DNA genome (1–5). At each end, the 5' strand overhangs the 3' strand by 12 bases. These single-stranded overhangs are complementary and anneal to form a *cos* site following entry into a host cell. Once annealed, the genome is a circular, completely double-stranded molecule which serves as a template for rolling-circle replication.

Many laboratory strains of lambda are derivatives of the strain λ . *cl857 ind1 Sam7*, which contains four point mutations relative to the wild type strain. The *ind1* mutation in the *cl* gene creates a new *HindIII* site at 37584 not present in the wild type. All lambda products sold by NEB are λ . *cl857 ind1 Sam7*.

Numbering of the genome sequence begins at the first (5'-most) base of the left end (bottom of diagram below) and continues rightward from late genes *nu1* and *A* towards the early genes. The map below shows the positions of all known ORFs larger than 200 codons.

Enzymes with unique restriction sites are shown in **bold type** and enzymes with two restriction sites are shown in regular type. **Coordinates indicate position of cutsite on the top strand. In previous catalogs, coordinates referred to the position of the 5' most base on the top strand, please make note of new numbering system.**



References

- (1) Echols, H. and Murialdo, H. (1978) *Microbiol. Rev.*, 42, 577–591.
- (2) Szybalski, E.H. and Szybalski, W. (1979) *Gene*, 7, 217–270
- (3) Daniels, D.L., de Wet, J.R. and Blattner, F.R. (1980) *J. Virol.*, 33, 390–400.
- (4) Sanger, F., Coulson, A.R., Hong, G.F., Hill, D.F. and Petersen, G.B. (1982) *J. Mol. Biol.*, 162, 729–773.
- (5) Daniels, D.L. et al. (1983). In R.W. Hendrix, J.W. Roberts, F.W. Stahl and R.A. Weisberg (Eds.), *Lambda II: Appendix*, New York: Cold Spring Harbor Press.

APPENDIX

(x) = enzyme not available from NEB


We recommend NEBcutter at NEBcutter.com to identify the restriction sites within your DNA sequence. NEBcutter indicates cut frequency and methylation-state sensitivity.

Coordinates indicate position of cutsite on the top strand. In previous catalogs, coordinates referred to the position of the 5' most base on the top strand, please make note of new numbering system.

M13 origin of replication arrows indicate the direction of synthesis of both the (+) and (-) strands.

Feature	Description	Coordinates
gene II	replication	6848-831 (cw)
gene X	replication	496-831
gene V	replication	843-1106
gene VII	minor coat protein	1108-1209
gene IX	minor coat protein	1206-1304
gene VIII	major coat protein	1301-1522
gene III	minor coat protein	1578-2852
gene VI	minor coat protein	2855-3193
gene I	phage assembly	3195-4241
gene XI (I ⁺)	phage assembly	3915-4241
gene IV	phage assembly	4219-5499
ori	M13 origin (+) of replication	5487-5867
<i>lacZα</i>	for α-complementation	6216-6722
MCS	multiple cloning site	6230-6286

(cw) = clockwise



- (1) Stewart, F.J. (2002) unpublished observations.
- (2) Messing, J. et al. (1977) *Proc. Natl. Acad. Sci. USA*, 74, 3652–3646.
- (3) Yanisch-Perron, C., Vieira, J. and Messing, J. (1985) *Gene*, 33, 103–119.

pBR322

GenBank Accession #: J01749

See page 128 for ordering information.

There are no restriction sites for the following enzymes: AbsI(x), Acc65I, AflII, AgeI, AjuI(x), AleI, Alol(x), ApaI, ArsiI(x), AscI, AsiSI, AvrII, BaeI, BaeI(x), BbvCI, BclI, BglII, BlnI, BmgBI, BplI(x), BsaXI, BseRI, BsiWI, BsrGI, BssHII, BstBI, BstEII, BstXI, Bsu36I, CspCI, DraIII, Eco53KI, Fall(x), FseI, HpaI, KflI(x), KpnI, MauBI(x), MfeI, MluI, MreI(x), MteI(x), NcoI, NotI, NsiI, PacI, PaeR7I, PaqCI, PstI(x), PmeI, PmlI, PstI, PspOMI, PspXI, PstI(x), RsrII, SacI, SacII, SbfI, SexAI, SfiI, SgrDI(x), SmaI, SnaBI, SpeI, SrfI, StuI, SwaI, TspMI, XbaI, XcmI, XhoI, XmaI

(x) = enzyme not available from NEB

pBR322 is an *E. coli* plasmid cloning vector containing the origin of replication from pMB1 (a plasmid in the ColE1 compatibility group; 1–3). The *rop* gene product, which regulates plasmid replication by stabilizing the interaction between RNAI and RNAII transcripts, maintains the copy number at about 20 per cell. However, pBR322 can be amplified with chloramphenicol or spectinomycin (4).

Enzymes with unique restriction sites are shown in **bold** type and enzymes with two restriction sites are shown in regular type. **Coordinates indicate position of cutsite on the top strand. In previous catalogs, coordinates referred to the position of the 5' most base on the top strand, please make note of new numbering system.**

Open reading frame (ORF) coordinates are in the form "translational start – translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons.

Origin of replication coordinates include the region from the -35 promoter sequence of the RNAII transcript to the RNA/DNA switch point. *bla* (Ap^R) gene coordinates include the signal sequence.


Feature	Coordinates	Source
<i>tet</i> (Tc ^R)	86-1276	pSC101
<i>bla</i> (Ap ^R)	4153-3293	<i>Tn3</i>
<i>rop</i>	1915-2106	pMB1
origin	3122-2534	pMB1

ori = origin of replication

Ap = ampicillin, Tc = tetracycline



We recommend NEBcutter at NEBcutter.com to identify the restriction sites within your DNA sequence. NEBcutter indicates cut frequency and methylation-state sensitivity.



References

- (1) Bolivar, F. et al. (1977) *Gene*, 2, 95–113.
- (2) Sutcliffe, J.G. (1979) *Cold Spring Harb. Symp. Quant. Biol.*, 43, 77–90.
- (3) Watson, N. (1988) *Gene*, 70, 399–403.
- (4) Sambrook, J., Fritsch, E.F., and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual*, (2nd ed.), Cold Spring Harbor, Cold Spring Harbor Laboratory Press.

pKLAC2

GenBank Accession #: EU196354
See page 231 for ordering information.

There are no restriction sites for the following enzymes: AatII, AbsI(x), Acc65I, AfeI, AflII, ApaI, AscI, AsiSI, AvrII, BbvCI, BliI, Bpu10I, BsiWI, FseI, FspAI(x), KfiI(x), KpnI, MauBI(x), MluI, MreI(x), MscI, MteI(x), PacI, PaeCI, PaeI(x), PmeI, PmlI, PspOMI, PspXI, PstI(x), RsrII, SfiI, SgrAI, SpeI, SrfI, SwaI, ZraI

(x) = enzyme not available from NEB



We recommend NEBcutter at NEBcutter.com to identify the restriction sites within your DNA sequence. NEBcutter indicates cut frequency and methylation-state sensitivity.

pKLAC2 is an expression vector capable both of replication in *E. coli* and stable integration into the genome of the yeast *Kluyveromyces fragilis* (1). It is designed for high-level expression of recombinant protein in *K. fragilis* using the *K. fragilis* Protein Expression Kit (NEB #E1000). pKLAC2 contains a universal multiple cloning site (MCS) that is compatible with all NEB expression systems.

In *E. coli*, it replicates using the pMB1 origin of replication from pBR322 (although the *rop* gene is missing) and carries the *bla* (Ap^R) marker for selection with ampicillin. Upon transformation of *K. fragilis* GG799 competent cells (NEB #C1001), SacII- or BstXI-linearized pKLAC2 integrates into the *K. fragilis* chromosome at the *LAC4* locus. Yeast transformants can be selected using the acetamidase selectable marker (*amdS*), which is expressed from the yeast *ADH1* promoter. Acetamidase expressed from pKLAC2 permits transformed cells to utilize acetamide as a sole nitrogen source on defined medium (2).

The multiple cloning site (MCS) is positioned to allow translational fusion of the *K. fragilis* α -mating factor secretion domain (α -MF) to the N-terminus of the recombinant target protein. This directs the fusion protein to the general secretory pathway, but the α -MF domain is cleaved off in the Golgi apparatus by the Kex protease, resulting in secretion of the recombinant protein alone.

Expression of the recombinant fusion protein is driven by the *K. fragilis* *LAC4* promoter, which has been modified to be transcriptionally silent in *E. coli* (1). This facilitates the cloning of proteins that are toxic to *E. coli*. This promoter is split such that when pKLAC2 is cleaved with SacII or BstXI, the recombinant protein and selectable marker are flanked by the two halves of the promoter. When these ends recombine with the *LAC4* promoter in the *K. fragilis* chromosome, the result is integration of the recombinant fusion protein (driven by the *LAC4* promoter) and *amdS* upstream of the *LAC4* gene (driven by a duplicate copy of the *LAC4* promoter) (2).

Enzymes with unique restriction sites are shown in bold type and selected enzymes with two restriction sites are shown in regular type. **Coordinates indicate position of cutsite on the top strand. In previous catalogs, coordinates referred to the position of the 5' most base on the top strand, please make note of new numbering system.**

Open reading frame (ORF) coordinates are in the form "translational start – translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons. Components of coordinated regions are indented below the region itself.


pMB1 origin of replication coordinates include the region from the –35 promoter sequence of the RNAII transcript to the RNA/DNA switch point. Promoter and transcription terminator coordinates represent cloned regions and not necessarily the precise functional elements.

Feature	Coordinates	Source
expression region:		
α -mating factor		
leader sequence	14-349	<i>K. fragilis</i>
MCS	257-354	—
<i>LAC4</i> TT region	371-953	<i>K. fragilis</i>
<i>AdH1</i> promoter region	1010-1712	<i>S. cerevisiae</i>
<i>amdS</i>	1713-3359	<i>A. nidulans</i>
<i>LAC4</i> promoter		
region (5' end)	4068-4648	<i>K. fragilis</i>
origin	5102-5690	pMB1
<i>bla</i> (Ap ^R)	6721-5861	Tn3
<i>LAC4</i> promoter		
region (3' end)	7475-9107	<i>K. fragilis</i> (modified)

ori = origin of replication

Ap = ampicillin

TT = transcription terminator



References

- (1) Colussi, P.A. and Taron, C.H. (2005) *Appl. Environ. Microbiol.*, 71, 7092–7098.
- (2) van Ooyen, A.J. et al. (2006) *FEMS Yeast Res.*, 6, 381–392.

pMAL-c6T

Sequence file available at www.neb.com.
See page 229 for ordering information.

Feature	Coordinates	Source
<i>lacI^q</i>	80-1162	<i>E. coli</i>
<i>P_{tac}</i>	1405-1432	—
expression ORF	1527-2761	—
<i>malE</i>	1527-2721	<i>E. coli</i>
MCS	2722-2761	—
<i>bla</i> (Ap ^R)	3101-3961	<i>Tn3</i>
origin	4049-4637	pMB1
<i>rop</i>	5007-5198	pMB1

There are no restriction sites for the following enzymes: AatII, AbsI(x), Acc65I, AflII, AgeI, AjuI(x), AleI, AseI(x), AscI, AsiSI, AvrII, BaeI, BarI(x), BbvCI, BmtI, BplI(x), BsaAI, BseRI, BsmFI, BspDI, BsrGI, BstBI, BstZ17I, ClaI, CspCI, DraIII, EcoNI, Fall(x), FseI, FspAI(x), KfiI(x), KpnI, MauBI(x), MreI(x), MscI, MteI(x), NaeI, NcoI, NdeI, NgoMIV, NheI, NruI, NsiI, PacI, PaeR7I, PaeGI, PaeI(x), PmeI, PmlI, PshAI, PspXI, PstI(x), SacII, SexAI, SfiI, SgrAI, SmaI, SnaBI, SpeI, SphI, SrfI, StuI, StylI, SwaI, TspMI, XbaI, XhoI, XmaI, XmnI, ZraI

(x) = enzyme not available from NEB



We recommend NEBcutter at NEBcutter.com to identify the restriction sites within your DNA sequence. NEBcutter indicates cut frequency and methylation-state sensitivity.

pMAL-c6T is an *E. coli* plasmid cloning vector designed for recombinant protein expression and purification using the NEBExpress MBP Fusion and Purification System (NEB #E8201) (1–3). It contains the pMB1 origin of replication from pBR322 and is maintained at a similar copy number to pBR322.

The multiple cloning site (MCS) is positioned to allow translational fusion of the *E. coli* maltose binding protein (MBP, encoded by the *malE* gene lacking its secretory signal sequence) to the *N*-terminus of the cloned target protein, thus resulting in an MBP-fusion protein localized in the cytoplasm. The pMAL-c6T vector contains a multiple cloning site (MCS) that is compatible with other NEB expression systems and is followed by stop codons in all three reading frames. In this vector, MBP has been engineered for tighter binding to amylose. This allows easy purification of the fusion protein, and the MBP domain can be subsequently removed using TEV Protease (NEB #P8112).

Transcription of the gene fusion is controlled by the inducible “*tac*” promoter (*P_{tac}*). Basal expression from *P_{tac}* is minimized by the binding of the Lac repressor, encoded by the *lacI^q* gene, to the *lac* operator immediately downstream of *P_{tac}*. A portion


of the *rnnB* operon containing two terminators, derived from the vector pKK233-2, prevents transcription originating from *P_{tac}* from interfering with plasmid functions.

Enzymes with unique restriction sites are shown in bold type.

Coordinates indicate position of cutsite on the top strand. In previous catalogs, coordinates referred to the position of the 5' most base on the top strand, please make note of new numbering system.

Open reading frame (ORF) coordinates are in the form “translational start – translational stop”; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons.

The pMB1 origin of replication includes the region from the -35 promoter sequence of the RNAlI transcript to the RNA/DNA switch point (labeled “ori”) and the *rop* gene, which controls expression of the RNAlI transcript. *bla* (Ap^R) gene coordinates include the signal sequence.



pMAL-c6T Polylinker:

5' *malE*...ATG CTG ATG GGC GGC CGC GAT ATC GTC GAC GGA TCC GAA TTC CCT GCA GGT AAT
 HindIII
 TAA ATA AGC TAA ATA AGC TTC AA...

References

- (1) Guan, C. et al. (1987) *Gene*, 67, 21–30.
- (2) Maina, C.V. et al. (1988) *Gene*, 74, 365–373.
- (3) Riggs, P.D. (1992). In F.M. Ausubel, et al. (Eds.), *Current Prot. in Molecular Biol.* New York: John Wiley & Sons, Inc.

pMiniT 2.0

Sequence available at www.neb.com
See page 97 for more information.

Feature	Coordinates	Source
Constitutive promoter	1-214	pNK2138
SP6 promoter	479-496	SP6
Toxic minigene	541-549	—
Synthetic T7 promoter	619-602	T7
<i>bla</i> (Ap ^R)	733-1593	<i>Tn3</i>
origin	1764-2352	pUC19

There are no restriction sites for the following enzymes: AbsI(x), Acc65I, AccI, AflII, AgeI, AjuI(x), AleI, Aol(x), ApaI, Arsl(x), AscI, AsiSI, AvrII, BaeI, BanII, Bari(x), BbsI, BbvCI, BclI, BglII, Bpl(x), BmgBI, BmtI, BplI(x), Bpu10I, BsaI, BsaAI, BsaBI, BseRI, BsgI, BsiWI, BsmFI, BsmI, BspDI, BspEI, BsrGI, BssHII, BstAPI, BstBI, BstEII, BstXI, BstZ17I, Bsu36I, BtgI, ClaI, CspCI, DraIII, Eco53kI, EcoNI, EcoO109I, EcoRV, Fall(x), FseI, FspAI(x), HincII, HindIII, HpaI, Kasi, KflI(x), KpnI, MauBI(x), MfeI, MluI, MreI(x), MscI, MteI(x), NaeI, NarI, NcoI, NgoMIV, NheI, NsiI, PaeI(x), PfiFI, PflMI, PfoI(x), PluTI, PmlI, PpuMI, PshAI, PsiI, PspOMI, PstI(x), PvuII, RsrII, SacI, SacII, SalI, SexAI, SfiI, SfoI, SgrAI, SgrDI(x), SmaI, SnaBI, SphI, SrfI, StuI, Styl, SwaI, TspMI, Tth111I, XbaI, XcmI, XmaI

(x) = enzyme not available from NEB




We recommend NEBcutter at NEBcutter.com to identify the restriction sites within your DNA sequence. NEBcutter indicates cut frequency and methylation-state sensitivity.




pMiniT 2.0 is an *E. coli* plasmid cloning vector designed for cloning blunt-ended or single-base overhang PCR products, or amplicons, using the NEB PCR Cloning Kit (NEB #E1202, #E1203). The pMiniT2.0 also enables *in vitro* transcription using SP6 and T7 promoters. It is compatible with Golden Gate Assembly as the BsaI site has been removed from the Ampicillin resistance gene.

In *E. coli*, it replicates using the pMB1 origin of replication from pUC19 and carries the *bla* (Ap^R) marker for selection with ampicillin. pMiniT2.0 contains a toxic minigene that is under the control of a constitutive promoter. If the pMiniT 2.0 vector recircularizes without an insert, the toxic minigene it will cause lethal inhibition of protein synthesis and no colony will result. If the pMiniT 2.0 Vector carries an insert, a colony will grow.

The map shown below displays the construct formed if no insert is present. Unique restriction sites are shown in **bold**. Additional restriction sites that can be used for subcloning are also shown. Expanded box below shows location of sequencing primers, restriction sites for subcloning or linearization for *in vitro* transcription, RNA Polymerase promoter sequences and placement of insertion site within the toxic minigene. **Coordinates indicate position of cutsite on the top strand. In previous catalogs, coordinates referred to the position of the 5' most base on the top strand, please make note of new numbering system.**



Features within Sequence Flanking the Toxic Minigene/Cloning Site:



pNEB206A

Sequence file available at www.neb.com.

There are no restriction sites for the following enzymes: AbsI(x), Acc65I, AccI, AfeI, AflII, AgeI, AjuI(x), AleI, AolI(x), ApaI, ArsiI(x), AsiSI, Aval, AvrII, BaeI, BclI(x), BbsI, BclI, BfuAI, BglII, BplI, BmgBI, BmtI, BplI(x), BsaAI, BsaBI, BsgI, BsiWI, BsmFI, BsmI, BsoBI, BspDI, BspEI, BspMI, BsrGI, BstBI, BstEII, BstXI, BstZ171, Bsu36I, BtgI, BtgZI, ClaI, CspCI, DraIII, EagI, EcoNI, EcoRV, Fall(x), FseI, FspAI(x), HincII, HpaI, KflI(x), KpnI, MauBI(x), MfeI, MluI, MreI(x), MscI, MteI(x), NaeI, NcoI, NgoMIV, NheI, NotI, NruI, NsiI, PaeR7I, PacCI, PasiI(x), PfiFI, PfiMI, PmlI, PpuMI, PshAI, PstI, PspOMI, PspXI, PstI(x), RsrII, SacII, SalI, SexAI, SfiI, SgrAI, SgrDI(x), SmaI, SnaBI, SpeI, SphI, SrfI, StuI, Styl, SwaI, TspMI, Tth111I, XcmI, XhoI, XmaI

(x) = enzyme not available from NEB

pNEB206A is an *E. coli* plasmid vector designed for fast and efficient cloning of PCR products to be used in conjunction with USER Enzyme (NEB #M5505; 1). It is derived from pNEB193 containing the high-copy pUC19 origin of replication and *lacZα* gene for screening of insertions at the cloning site using α-complementation (2).

The plasmid is supplied in a linearized form 2,706 bp in length (with bp 438-453 excised from the circular form), flanked by two noncomplementary 8-base 3' overhangs at the intended cloning site. Amplification with deoxyuridine-containing primers and subsequent treatment (as defined in the protocol "Cloning with USER Enzyme" found on our website), results in PCR products with 5' overhangs complementary to those in pNEB206A. These products can be directionally cloned into pNEB206A at high efficiency without the use of restriction enzymes or DNA ligase, forming recombinant circular molecules.

Enzymes with unique restriction sites are shown in **bold type**, and enzymes with two restriction sites are shown in regular type. **Coordinates indicate position of cutsite on the top strand. In previous catalogs, coordinates referred to the position of the 5' most base on the top strand, please make note of new numbering system.**

Open reading frame (ORF) coordinates are in the form "translational start – translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons.


Origin of replication coordinates include the region from the -35 promoter sequence of the RNAlI transcript to the RNA/DNA switch point. *bla* (Ap^R) gene coordinates include the signal sequence. Cloning site coordinates include those bases in the circular form that are single-stranded in or missing from the supplied linear form.

Feature	Coordinates	Source
<i>lacZα</i>	505-146	—
cloning site	430-461	—
origin	1491-903	pUC19
<i>bla</i> (Ap ^R)	2522-1662	Tn3


ori = origin of replication
Ap = ampicillin



We recommend NEBcutter at NEBcutter.com to identify the restriction sites within your DNA sequence. NEBcutter indicates cut frequency and methylation-state sensitivity.



pNEB206A (linearized form) cloning site:



References

- (1) Bitinaite, J. and Vaiskunaite, R. (2003) unpublished observations.
- (2) Yanisch-Perron, C. et al. (1985) *Gene*, 33, 103-119.

Sequence file available at www.neb.com.
See page 298 for ordering information.

There are no restriction sites for the following enzymes: *AbstI*(x), *AfeI*, *AfIII*, *AjuI*(x), *AlfI*(x), *Alol*(x), *AsiSI*, *BaeI*, *BamHI*(x), *BbvCI*, *BlpI*, *BplI*(x), *BsiWI*, *BsmBI*, *BspDI*, *BspEI*, *BstAPI*, *BstBI*, *BstEII*, *ClaI*, *EcoNI*, *Esp3I*, *FseI*, *FspAI*(x), *KllI*(x), *MauBI*(x), *MreI*(x), *MteI*(x), *PasI*(x), *PfoI*(x), *PshAI*, *PsrI*(x), *SexAI*, *SgrAI*, *SrlI*, *StuI*, *XcmI*

(x) = enzyme not available from NEB



We recommend NEBcutter at NEBcutter.com to identify the restriction sites within your DNA sequence. NEBcutter indicates cut frequency and methylation-state sensitivity.

pSNAP_f Vector is a mammalian expression plasmid intended for the cloning and stable or transient expression of SNAP-tag[®] protein fusions in mammalian cells. This plasmid encodes SNAP_f, a SNAP-tag protein, which is expressed under control of the CMV promoter. SNAP_f is an improved version of the SNAP-tag which exhibits faster labeling kinetics. The SNAP-tag is a novel tool for protein research, allowing the specific, covalent attachment of virtually any molecule to a protein of interest. The SNAP-tag is a small protein based on human O6-alkyl-guanine-DNA-alkyltransferase (hAGT). SNAP-tag substrates are derivatives of benzyl purines and benzyl pyrimidines. In the labeling reaction, the substituted benzyl group of the substrate is covalently attached to the SNAP-tag. Use of this system involves two steps: sub-cloning and expression of the protein of interest as a SNAP-tag fusion, and labeling of the fusion with the SNAP-tag substrate of choice. Further details are provided with the SNAP-Cell Starter Kit (NEB #E9100) and SNAP-Surface Starter Kit (NEB #E9120).

Codon usage of the gene is optimized for expression in mammalian cells. pSNAP_f contains two multiple cloning sites to allow cloning of the fusion partner as a fusion to the N- or C-terminus of the SNAP-tag. The expression vector has an Internal Ribosome Entry Site (IRES) and a neomycin resistance gene downstream of the SNAP-tag for the efficient selection of stable transfectants.

Enzymes with unique restriction sites are shown in bold type.

Coordinates indicate position of cutsite on the top strand. In previous catalogs, coordinates referred to

the position of the 5' most base on the top strand, please make note of new numbering system.

Open reading frame (ORF) coordinates are in the form "translational start – translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons. Component genes or regions of fusion ORFs are indented below the ORF itself.

pUC19 origin of replication coordinates include the region from the -35 promoter sequence of the RNAlI transcript to the RNA/DNA switch point. *bla* (Ap^R) gene coordinates include the signal sequence.


Feature	Coordinates	Source
CMV promoter	251-818	—
expression region	915-1564	—
MCS1	915-965	—
SNAP _f	969-1514	—
MCS2	1515-1564	—
IRES	1910-2500	ECMV
Neo ^R	2536-3339	Tn5
ori	4094-4682	pUC19
<i>bla</i> (Ap ^R)	4853-5713	Tn3

ori = origin of replication

Ap = ampicillin

Neo = neomycin

IRES = internal ribosomal entry site



MCS1

NheI *AscI* *SwaI* *BsrGI* *EcoRI*
...GCTAGC GATATCGGCG CGCCAGCATT TAAATCTGTA CAGACCGGTG AATTC
CGATCG CTATAGCCGC GCGGTCGTAA ATTTAGACAT GTCTGGCCAC TTAAG...

MCS2

SbfI *BamHI* *PmeI* *XhoI* *PacI* *NotI*
...CCTGCA GGC GGATCCG CGTTTAACT CGAGGTTAAT TAATGAGCGG CCGC
GGACGT CCGCCTAGGC GCAAATTGA GCTCCAATTA ATTACTCGCC GGCG...

pTXB1

Sequence file available at www.neb.com.
See page 230 for ordering information.

Feature	Coordinates	Source
<i>bla</i> (Ap ^R)	140-1000	<i>Tn3</i>
M13 origin	1042-1555	M13
origin	1666-2254	pMB1
<i>rop</i>	2814-2623	pMB1
<i>lacI</i>	4453-3371	<i>E. coli</i>
T7 promoter	5637-5654	T7
expression ORF	5725-6558	—
MCS	5722-5775	—
<i>Mxe</i> GyrA intein	5776-6369	<i>M. xenopi</i>
CBD	6400-6558	<i>B. circulans</i>

ori = origin of replication

Ap = ampicillin

There are no restriction sites for the following enzymes: Acc65I, AfIII, AjuI(x), AleI, AroI(x), AscI, AsiSI, AvrII, BaeI, BbvCI, BglII, BlnI(x), BmgBI, Bpu10I, BseRI, BspDI, BstBI, Bsu36I, ClaI, CspCI, Eco53KI, Fall(x), FseI, FspAI(x), HindIII, KflI(x), KpnI, MauBI(x), MscI, MteI(x), NcoI, NsiI, PacI, PaeCI, PaeI(x), PmlI, PpuMI, PstI(x), RsrII, SacI, SanDI(x), SbfI, SexAI, SfiI, SgrDI(x), SmaI, SnaBI, SrfI, TspMI, XmaI

(x) = enzyme not available from NEB



We recommend NEBcutter at NEBcutter.com to identify the restriction sites within your DNA sequence. NEBcutter indicates cut frequency and methylation-state sensitivity.

pTXB1 is an *E. coli* plasmid cloning vector designed for recombinant protein expression, purification, and ligation using the IMPACT™ Kit (NEB #E6901) (1,2). It contains the pMB1 origin of replication from pBR322 and is maintained at a similar copy number to pBR322; in addition, pTXB1 also contains an M13 origin of replication.

The multiple cloning site (MCS) is positioned to allow translational fusion of the *Mxe* GyrA intein tag to the C-terminus of the cloned target protein (2,3). The chitin binding domain (CBD) from *B. circulans*, fused to the C-terminus of the intein, facilitates purification of the intein-target protein precursor.

Transcription of the gene fusion is controlled by the inducible T7 promoter, requiring *E. coli* strains containing integrated copies of the T7 RNA polymerase gene [e.g., C2566 or BL21(DE3)] for expression. Basal expression from the T7 promoter is minimized by the binding of the Lac repressor, encoded by the *lacI* gene, to the *lac* operator immediately downstream of the T7 promoter (4). Translation of the fusion utilizes the translation initiation signal (Shine Dalgarno sequence) from the strongly expressed T7 gene 10 protein (φ10).



pTXB1 and pTXB3 are identical except for the MCS regions: pTXB1 contains an NdeI site, and pTXB3 an NcoI site, overlapping the initiating methionine codon of the intein fusion gene. The N-terminal cysteine residue ("Cys_N") of the intein is shaded.

Enzymes with unique restriction sites are shown in bold type.

Coordinates indicate position of cutsite on the top strand. In previous catalogs, coordinates referred to the position of the 5' most base on the top strand, please make note of new numbering system.

Open reading frame (ORF) coordinates are in the form "translational start – translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons. Component genes or regions of fusion ORFs are indented below the ORF itself.

pMB1 origin of replication coordinates include the region from the -35 promoter sequence of the RNAlI transcript to the RNA/DNA switch point. For the M13 origin, the arrow shows the direction of synthesis of the (+) strand, which gets packaged into phage particles. *bla* (Ap^R) gene coordinates include the signal sequence.



References

- (1) Chong, S. et al. (1997) *Gene*, 192, 271–281.
- (2) Evans, T.C., Benner and Xu, M.-Q. (1998) *Protein Sci.*, 7, 2256–2264.
- (3) Southworth, M.W. et al. (1999) *Biotechniques*, 27, 110–120.
- (4) Dubendorff, J.W. and Studier, F.W. (1991) *J. Mol. Biol.*, 219, 45–59.

pTYB21

Sequence file available at www.neb.com.
See page 230 for ordering information.

Feature	Coordinates	Source
<i>bla</i> (Ap ^R)	140-1000	<i>Tn3</i>
M13 origin	1042-1555	M13
origin	1666-2254	pMB1
<i>rop</i>	2814-2623	pMB1
<i>lacI</i>	4453-3371	<i>E. coli</i>
T7 promoter	5637-5654	T7
expression ORF	5725-7368	—
MCS	7301-7361	—
<i>Sce</i> VMA intein	5770-7299	<i>S. cerevisiae</i>
CBD	6595-6747	<i>B. circulans</i>

ori = origin of replication
Ap = ampicillin

There are no restriction sites for the following enzymes: AatII, AbsI(x), AflII, AgeI, AjuI(x), AscI, AsiSI, AvrII, BbvCI, BmgBI, BplI(x), BseRI, BsiWI, BsmI, BspDI, Bsu36I, ClaI, CspCI, Fall(x), FseI, FspAI(x), KflI(x), MauBI(x), MreI(x), MteI(x), NruI, NsiI, PacI, PaeR7I, PaqCI, Pasi(x), PpuMI, PspXI, PstI(x), RsrII, SexAI, SfiI, SgrAI, SmaI, SnaBI, SrfI, TspMI, XhoI, XmaI, ZraI

(x) = enzyme not available from NEB



We recommend NEBcutter at NEBcutter.com to identify the restriction sites within your DNA sequence. NEBcutter indicates cut frequency and methylation-state sensitivity.

pTYB21 is an *E. coli* plasmid cloning vector designed for recombinant protein expression and purification using the IMPACT™ Kit (NEB #E6901) (1,2). It contains the pMB1 origin of replication from pBR322 and is maintained at a similar copy number to pBR322; in addition, pTYB21 also contains an M13 origin of replication.

The multiple cloning site (MCS) is positioned to allow translational fusion of the *Sce* VMA intein tag to the N-terminus of the cloned target protein (2). The chitin binding domain (CBD) from *B. circulans*, facilitates purification of the intein-target protein precursor.

Transcription of the gene fusion is controlled by the inducible T7 promoter, requiring *E. coli* strains containing integrated copies of the T7 RNA polymerase gene [e.g., C2566 or BL21(DE3)] for expression. Basal expression from the T7 promoter is minimized by the binding of the Lac repressor, encoded by the *lacI* gene, to the *lac* operator immediately downstream of the T7 promoter (3). Translation of the fusion utilizes the translation initiation signal


(Shine Dalgarno sequence) from the strongly expressed T7 gene 10 protein (ϕ10).

pTYB21 contains a *SapI* site which allows for cloning of a target gene without any extra amino acids. pTYB22 is identical to pTYB21 except for the MCS regions (see below). pTYB22 contains an *NdeI* site overlapping the initiating methionine codon of the intein fusion gene. pTYB21 differs from pTYB11 in that it contains a universal MCS that is compatible with all NEB expression systems.

Enzymes with unique restriction sites are shown in bold type. Coordinates indicate position of cutsite on the top strand. In previous catalogs, coordinates referred to the position of the 5' most base on the top strand, please make note of new numbering system.

Open reading frame (ORF) coordinates are in the form "translational start – translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons. Component genes or regions of fusion ORFs are indented below the ORF itself.

pMB1 origin of replication coordinates include the region from the -35 promoter sequence of the RNAII transcript to the RNA/DNA switch point. For the M13 origin, the arrow shows the direction of synthesis of the (+) strand, which gets packaged into phage particles. *bla* (Ap^R) gene coordinates include the signal sequence.



...TAATACGACTCACTATAGGGGAATTGTG...GAAGACGATTATTATGGGATTACTTTATCTGATGATTCTGATCATCAGTTTTTCTTGATCTCAG
5650 7220 7240 7260 7280

pTYB21 MCS

Scd VMA Intein
GTTGTTGTACAGAAC GGAAGAGCTCATATGTCCATGGGCGCGGATATCGTCGACGGATCCGAATTCCTCGAGGTAATTAATAAC...
V V V Q N G R A H M S M G G R D I V D G S E F P A G N *

pTYB22 MCS

Scd VMA Intein
GTTGTTGTACAGAAAT GCTGTGTCATATGTCCATGGGCGCGGATATCGTCGACGGATCCGAATTCCTCGAGGTAATTAATAAC...
V V V Q N A G H M S M G G R D I V D G S E F P A G N *

References

- (1) Chong et al. (1996) *J. Biol. Chem.*, 271, 22159-22168
- (2) Chong et al. (1998) *NAR*, 26, 5109-5115.
- (3) Dubendorff, J.W. and Studier, F.W. (1991) *J. Mol. Biol.*, 219, 45-59.

pUC19

GenBank Accession #: L09137
See page 128 for ordering information.

Feature	Coordinates	Source
<i>lacZα</i>	469-146	—
origin	1455-867	pMB1 (mutant)
<i>bla</i> (Ap ^R)	2486-1626	<i>Tn3</i>

ori = origin of replication
Ap = ampicillin

There are no restriction sites for the following enzymes: AbsI(x), AfeI, AflII, AgeI, AjuI(x), AleI, AfoI(x), ApaI, AseI(x), AscI, AsiSI, AvrII, BaeI, BclI(x), BbsI, BbvCI, BclI, BglII, BliI, BmgBI, BmtI, BplI(x), Bpu10I, BsaAI, BsaBI, BseRI, BsgI, BsiWI, BsmFI, BsmI, BspDI, BspEI, BsrGI, BssHII, BstBI, BstEII, BstXI, BstZ17I, Bsu36I, BtgI, BtgZI, ClaI, CspCI, DraIII, EagI, EcoNI, EcoRV, Fall(x), FseI, FspAI(x), HpaI, KfiI(x), MauBI(x), MfeI, MluI, MreI(x), MscI, MteI(x), NaeI, NcoI, NgoMIV, NheI, NotI, NruI, NsiI, PacI, PaeR7I, PaeCI, PaeI(x), PfiFI, PfiMI, PmeI, PmlI, PpuMI, PshAI, PstI, PspOMI, PspXI, PstI(x), RsrII, SacII, SexAI, SfiI, SgrAI, SgrDI(x), SnaBI, SpeI, SrfI, StuI, Styl, SwaI, Tth111I, XcmI, XhoI

(x) = enzyme not available from NEB



We recommend NEBcutter at NEBcutter.com to identify the restriction sites within your DNA sequence. NEBcutter indicates cut frequency and methylation-state sensitivity.

pUC19 is a small, high-copy number *E. coli* plasmid cloning vector containing portions of pBR322 and M13mp19 (1). It contains the pMB1 origin of replication from pBR322, but it lacks the *rop* gene and carries a point mutation in the RNAII transcript (G 2975 in pBR322 to A 1308 in pUC19; 2). These changes together result in a temperature-dependent copy number of about 75 per cell at 37°C and > 200 per cell at 42°C (2,3). The multiple cloning site (MCS) is in frame with the *lacZα* gene, allowing screening for insertions using α-complementation.

pUC18 is identical to pUC19 except that the MCS region (nt 397-454) is inverted.

pNEB193 is also identical to pUC19 except for the addition of several restriction endonuclease sites to the MCS. Its total length is 2713 bp.


Enzymes with unique restriction sites are shown in **bold type**, and enzymes with two restriction sites are shown in regular type. **Coordinates indicate position of cutsite on the top strand. In previous catalogs, coordinates referred to the position of the 5' most base on the top strand, please make note of new numbering system.**

Open reading frame (ORF) coordinates are in the form "translational start – translational stop"; numbers refer to positions on the top (clockwise) strand, regardless of the direction of transcription and include the start and stop codons.


Origin of replication coordinates include the region from the -35 promoter sequence of the RNAII transcript to the RNA/DNA switch point. *bla* (Ap^R) gene coordinates include the signal sequence.

References

- (1) Yanisch-Perron, C., Vieira, J. and Messing, J. (1985) *Gene*, 33, 103–119.
- (2) Lin-Chao, S., Chen, W.-T. and Wong, T.-T. (1992) *Mol. Microbiol.*, 6, 3385–3393
- (3) Miki, T. et al. (1987) *Protein Eng.*, 1, 327–332.



pUC19 MCS



pNEB193 MCS



The Genetic Code

	A	R	N	D	C	Q	E	G	H	I	L	K	M	F	P	S	T	W	Y	V
	Ala	Arg	Asn	Asp	Cys	Gln	Glu	Gly	His	Ile	Leu	Lys	Met	Phe	Pro	Ser	Thr	Trp	Tyr	Val
5'	GCA	CGA	AAC	GAC	UGC	CAA	GAA	GGA	CAC	AUA	CUA	AAA	AUG	UUC	CCA	UCA	ACA	UGG	UAC	GUA
	C	C	U	U	U	G	G	C	U	C	C	G	U	U	C	C	C	U	U	C
	G	G	U	U	U	G	G	G	U	U	G	U	U	U	G	G	G	U	U	G
	U	U	U	U	U	G	G	U	U	U	U	U	U	U	U	U	U	U	U	U
		or									or					or				
		AGA									UUA					AGC				
		G									G					U				

Second Position

	U	C	A	G	
U	UUU] Phe UUC] UUA] Leu UUG]	UCU] UCC] Ser UCA] UCG]	UAU] Tyr UAC] Stop UAA] Stop UAG]	UGU] Cys UGC] Stop UGA] Stop UGG] Trp	U C A G
C	CUU] CUC] Leu CUA] CUG]	CCU] CCC] Pro CCA] CCG]	CAU] His CAC] Gln CAA] CAG]	CGU] CGC] Arg CGA] CGG]	U C A G
A	AUU] AUC] Ile AUA] Met AUG]	ACU] ACC] Thr ACA] ACG]	AAU] Asn AAC] Lys AAA] AAG]	AGU] Ser AGC] Arg AGA] AGG]	U C A G
G	GUU] GUC] Val GUA] GUG]	GCU] GCC] Ala GCA] GCG]	GAU] Asp GAC] Glu GAA] GAG]	GGU] GGC] Gly GGA] GGG]	U C A G

Third Position (3' end)

Termination Signals

UAA (Ochre)
UAG (Amber)
UGA (Opal)

Single Letter Code



A = adenosine
C = cytidine
G = guanosine
T = thymidine
U = uridine

B = C or G or T
D = A or G or T
H = A or C or T
K = G or T
M = A or C
N = A or C or G or T
R = A or G
S = C or G
V = A or C or G
W = A or T
Y = C or T




Amino Acid Structures

Each amino acid is accompanied by its three- and one-letter code, residue molecular weight (actual molecular weight minus water) and side-chain pK_a where appropriate.






SMALL






NUCLEOPHILIC





HYDROPHOBIC



AROMATIC





ACIDIC






APPENDIX

AMIDE




BASIC



DNA Base Pairs

The structures of the adenosine:thymidine and guanosine:cytidine base pairs are shown in the context of the ribose phosphodiester backbones. The numbering schemes of the ribose and nucleotide moieties are indicated. Arrows indicate the polarity of each strand from 5' to 3'.



Nucleic Acid Data

Average weight of a DNA basepair (sodium salt) = 650 daltons

1.0 A_{260} unit ds DNA = 50 $\mu\text{g}/\text{ml}$ = 0.15 mM (in nucleotides)

1.0 A_{260} unit ss DNA = 33 $\mu\text{g}/\text{ml}$ = 0.10 mM (in nucleotides)

1.0 A_{260} unit ss RNA = 40 $\mu\text{g}/\text{ml}$ = 0.11 mM (in nucleotides)

MW of a double-stranded DNA molecule = (# of base pairs) \times (650 daltons/base pair)

Moles of ends of a double-stranded DNA molecule = $2 \times (\text{grams of DNA}) / (\text{MW in daltons})$

Moles of ends generated by restriction endonuclease cleavage:

a) circular DNA molecule: $2 \times (\text{moles of DNA}) \times (\text{number of sites})$

b) linear DNA molecule: $2 \times (\text{moles of DNA}) \times (\text{number of sites}) + 2 \times (\text{moles of DNA})$

1 μg of 1000 bp DNA = 1.52 pmol = 9.1×10^{11} molecules

1 μg of pUC18/19 DNA (2686 bp) = 0.57 pmol = 3.4×10^{11} molecules

1 μg of pBR322 DNA (4361 bp) = 0.35 pmol = 2.1×10^{11} molecules

1 μg of M13mp18/19 DNA (7249 bp) = 0.21 pmol = 1.3×10^{11} molecules

1 μg of λ DNA (48502 bp) = 0.03 pmol = 1.8×10^{10} molecules

1 pmol of 1000 bp DNA = 0.66 μg

1 pmol of pUC18/19 DNA (2686 bp) = 1.77 μg

1 pmol of pBR322 DNA (4361 bp) = 2.88 μg

1 pmol of M13mp18/19 DNA (7249 bp) = 4.78 μg

1 pmol of λ DNA (48502 bp) = 32.01 μg

1.0 kb DNA = coding capacity for 333 amino acids \approx 37,000 dalton protein

10,000 dalton protein \approx 270 bp DNA


50,000 dalton protein \approx 1.35 kb DNA

Isotope Data


Isotope	Particle Emitted	Half Life	
^{14}C	β	5,730 years	1 Ci = 1,000 mCi
^3H	β	12.3 years	1 mCi = 1,000 μCi
^{125}I	γ	60 days	1 μCi = 2.2×10^6 disintegrations/minute
^{32}P	β	14.3 days	1 Becquerel = 1 disintegration/second
^{33}P	β	25 days	1 μCi = 3.7×10^4 Becquerels
^{35}S	β	87.4 days	1 Becquerel = 2.7×10^{-5} μCi

Messenger RNA Modifications


In nature, ribonucleic acid undergoes extensive chemical modification that can result in altered function or stability. The figure below shows examples of base and ribose modifications commonly found in native mRNAs




1-methyladenosine (m¹A) [1A]




N⁶-methyladenosine (m⁶A) [6A]




N⁶,2'-O-dimethyladenosine (m⁶Am) [6Am]




inosine (I) [9A]




2'-O-methylnucleoside (Nm) [ON]




N⁴-acetylcytidine (ac⁴C) [42C]




3-methylcytidine (m³C) [3C]



5-methylcytidine (m⁵C) [5C]



5-hydroxymethylcytidine (hm⁵C) [51C]



pseudouridine (Ψ) [9U]

Key:

Nucleoside Name (short name)
[new nomenclature]

Acids and Bases

Compound	Formula	Molecular Weight	Specific Gravity	% by Weight	Conc Reagent Molarity
Acetic acid, glacial	CH ₃ COOH	60.0	1.05	99.5	17.4
Formic acid	HCOOH	46.0	1.20	90	23.4
Hydrochloric acid	HCl	36.5	1.18	36	11.6
Nitric acid	HNO ₃	63.0	1.42	71	16.0
Perchloric acid	HClO ₄	100.5	1.67	70	11.6
Phosphoric acid	H ₃ PO ₄	98.0	1.70	85	18.1
Sulfuric acid	H ₂ SO ₄	98.1	1.84	96	18.0
Ammonium hydroxide	NH ₄ OH	35.0	0.90	28	14.8
Potassium hydroxide	KOH	56.1	1.52	50	13.5
Sodium hydroxide	NaOH	40.0	1.53	50	19.1
β-mercaptoethanol	HSCH ₂ CH ₂ OH	78.1	1.11	100	14.3

Protein Data

Bacterial Cells: *E. coli* or *Salmonella typhimurium*

	Cell Data	per cell	per liter at 10 ⁹ cells per ml
Theoretical maximum yield for a 1 liter culture (10 ⁹ cells/ml) if protein of interest is:	Wet Weight	9.5 x 10 ⁻¹³ g	0.95 g
0.1% of total protein: 150 µg/liter	Dry Weight	2.8 x 10 ⁻¹³ g	0.28 g
2.0% of total protein: 3 mg/liter	Total Protein	1.55 x 10 ⁻¹³ g	0.15 g
50.0% of total protein: 75 mg/liter	Volume	1.15 µm ³ = 1 femtoliter	
	Protein Conc. in the cell:	135 mg/ml	

Common Plasmid
Gene Products

Gene	Gene Product # of Residues	Molecular Weight (daltons)
<i>tet</i> (pBR322)	401	43,267
<i>amp</i> (pBR322, bla)	286	31,515
<i>kan</i> (pACYC177, nptI)	264	29,047
<i>cam</i> (pACYC184, cat)	219	25,663
<i>lacZα</i> (pUC19)	107	12,232
<i>lacZ</i>	1,023	116,351

Nucleotide
Physical Properties

Compound	Molecular Weight	λ max (pH 7.0)	Absorbance at λ max 1 M solution (pH 7.0)
ATP	507.2	259	15,400
CTP	483.2	271	9,000
GTP	523.2	253	13,700
UTP	484.2	262	10,000
dATP	491.2	259	15,200
dCTP	467.2	271	9,300
dGTP	507.2	253	13,700
dTTP	482.2	267	9,600

pH vs Temperature
for Tris Buffer

pH of Tris Buffer (0.05 M)		
5°C	25°C	37°C
7.76	7.20	6.91
7.89	7.30	7.02
7.97	7.40	7.12
8.07	7.50	7.22
8.18	7.60	7.30
8.26	7.70	7.40
8.37	7.80	7.52
8.48	7.90	7.62
8.58	8.00	7.71
8.68	8.10	7.80
8.78	8.20	7.91
8.88	8.30	8.01
8.98	8.40	8.10
9.09	8.50	8.22
9.18	8.60	8.31
9.28	8.70	8.42

Agarose Gel Resolution

% Gel	Optimum Resolution for Linear DNA (kb)
0.5	30 to 1.0
0.7	12 to 0.8
1.0	10 to 0.5
1.2	7 to 0.4
1.5	3 to 0.2

INDEX

#

6-tube Magnetic Separation Rack	235
50 ml Magnetic Separation Rack	235
12-tube Magnetic Separation Rack	235
96-well Microtiter Plate Magnetic Separation Rack	235
1 kb DNA Ladder	181
1 kb Plus DNA Ladder	181
1 kb Plus DNA Ladder for Safe Stains	182
3'-Biotin-GTP	205
5' Deadenylase	214
7-deaza-dGTP	85
3' Desthiobiotin-GTP	205
5' DNA Adenylation Kit	211, 214
9° N DNA Ligase	104
50 bp DNA Ladder	181
100 bp DNA Ladder	181
5-hydroxymethyluridine DNA Kinase	105
5-methyl-dCTP	85, 287

A

α -N-Acetylgalactosaminidase	259
α 1-2 Fucosidase	259
α 1-2,3,4,6 Fucosidase	260
α 1-2,4,6 Fucosidase O	260
α 1-3,4 Fucosidase	260
α 1-3,6 Galactosidase	260
α 1-3,4,6 Galactosidase	261
α -Lytic Protease	266
α 1-2,3 Mannosidase	262
α 1-2,3,6 Mannosidase	262
α 1-6 Mannosidase	262
α 2-3 Neuraminidase S	263
α 2-3,6,8 Neuraminidase	262
α 2-3,6,8,9 Neuraminidase A	263
AatII	22
AbaSI	280
AccI	22
Acc65I	22
Acil	22
Acids and Bases	369
AcII	22
Activity Chart for Restriction Enzymes	309-314
Activity of DNA Modifying Enzymes in CutSmart Buffer	315
Activity of enzymes at 37°C	315
AcuI	22
Acyclonucleotide Set	85
Adaptors, Sequencing	146, 147, 160
Adenosine 5'-Triphosphate (ATP)	85, 271
AfeI	22
AfiII	22
AfiIII	23
Afu Uracil-DNA Glycosylase (UDG)	121
Agarase (β)	123
Agel	23
Agel-HF	23
AhdI	23
Albumin	58
AleI-v2	23
AluI	23
AluI Methyltransferase	124, 286
AlwI	23
AlwNI	23
Amino Acid Structures	366
Amplification, DNA	62-87
Amylose Magnetic Beads	229, 233, 234
Amylose Resin	233, 234
Amylose Resin High Flow	233, 234
AMV Reverse Transcriptase	86, 206, 207
Antarctic Phosphatase	106, 216
Antarctic Thermolabile UDG	64, 122

Antibodies

Anti-CBD Monoclonal Antibody	230, 234
Anti-M13 pIII Monoclonal Antibody	273
Anti-MBP Magnetic Beads	233, 234
Anti-MBP Monoclonal Antibody	229, 234
Anti-SNAP-tag Antibody, Polyclonal	297
Goat Anti-Mouse IgG Magnetic Beads	233, 235
Goat Anti-Rabbit IgG Magnetic Beads	233, 235
Goat Anti-Rat IgG Magnetic Beads	233, 235
Anti-MBP Magnetic Beads	233, 234
Apal	24
ApalI	24
APE I	118
ApeKI	24
Apol	24
Apol-HF	24
Apyrase	107
ARCA Cap Analog	204
ARTIC, Library Preparation	161
Ascl	24
Asel	24
AsiSI	24
AspN, Endoproteinase	267
ATP	85, 271
Aval	25
Avall	25
Avoiding Ribonuclease Contamination	198
AvrII	25

B

β -Agarase I	123
β 1-3 Galactosidase	261
β 1-3,4 Galactosidase	261
β 1-4 Galactosidase S	261
β -N-Acetylglucosaminidase S	259
β -N-Acetylhexosaminidase	259
<i>Bacteriodes</i> Heparinase I, II, III	264
Bael	25
BaeGI	25
BamHI	25
BamHI-HF	25
BamHI Methyltransferase	124, 286
BanI	25
BanII	26
Basic Research at NEB	2
BbsI	26
BbsI-HF	26
BbvI	26
BbvCI	26
BC-NH2	299
BccI	26
BceAI	26
Bcgl	26
BciVI	27
BclI	27
BclI-HF	27
BcoDI	27
Bfal	27
BfuAI	27
BG-GLA-NHS	299
BG-Maleimide	299
BG-NH2	299
BG-PEG-NH2	299
BglI	27
BglII	27
BioBrick Assembly Kit	97
Biotin Labels	298
Bisulfite Sequencing	151, 279
BL21 Competent <i>E. coli</i>	232, 240, 241, 245
BL21(DE3) Competent <i>E. coli</i>	232, 240, 241, 245
BlpI	28
Blue Gel Loading Dye	59
Blue Protein Loading Dye	185

Blue Prestained Protein Std., Broad range (11–250 kDa)	185
Blunt/TA Ligase Master Mix	99, 100, 101
Blunting, DNA	
Exonuclease T	112, 217
Klenow DNA Polymerase	83
Quick Blunting Kit	98
T4 DNA Polymerase	83
T4 Polynucleotide Kinase	104
BmgBI	28
Bmrl	28
BmtI	28
BmtI-HF	28
Bovine Serine Albumin	20, 58
Bpml	28
Bpu10I	28
BpuEI	28
BSA	20, 58
BsaI-HFv2	29
BsaAI	29
BsaBI	29
BsaHI	29
BsaJI	29
BsaWI	29
BsaXI	29
BseRI	29
BseYI	30
BsgI	30
BsiEI	30
BsiHKA1	30
BsiWI	30
BsiWI-HF	30
BsII	30
BsmI	30
BsmAI	31
BsmBI-v2	31
BsmFI	31
BsoBI	31
Bsp1286I	31
BspCNI	31
BspDI	31
BspEI	31
BspHI	32
BspMI	32
BspQI	32
BsrI	32
BsrBI	32
BsrDI	32
BsrFI-v2	32
BsrGI	32
BsrGI-HF	33
BssHII	33
BssSI-v2	33
Bst 2.0, Bst 2.0 WarmStart DNA Polymerases	64, 79
Bst 3.0 DNA Polymerase	64, 79
BstAPI	33
BstBI	33
Bst DNA Polymerase, Full Length	79
Bst DNA Polymerase, Large Fragment	79
BstEII	33
BstEII-HF	33
BstNI	33
BstUI	34
BstXI	34
BstYI	34
BstZ17I-HF	34
Bsu DNA Polymerase, Large Fragment	64, 84
Bsu36I	34
BtgI	34
BtgZI	34
BtsI-v2	34
BtsIMutI	35
BtsCI	35
Buffers	
Diluent, Restriction Endonucleases	59, 309
Double Digestion, Restriction Enzymes	307
Reaction, Restriction Endonucleases	58, 309

Thermophilic Polymerases, Reaction Buffers	85
Tris Buffer, pH of as a Function of Temperature	369

C

Cac8I	35
cAMP-dependent Protein Kinase	271
Cap Analogs, RNA	204
Capping System, Vaccinia	204
Cas9 Enzymes	190, 191
Cas9 Nuclease, <i>S. pyogenes</i>	190, 191
Casein Kinase II	271
cDNA Synthesis	
AMV Reverse Transcriptase	86, 206, 207
LunaScript RT Supermix Kit	73, 86, 206
LunaScript RT Master Mix Kit (Primer-free)	73, 86, 206
M-MuLV Reverse Transcriptase	86, 206, 207
Primers	208
ProtoScript II Reverse Transcriptase	86, 206
ProtoScript II First Strand cDNA Synthesis Kit	86, 206, 209
ProtoScript First Strand Synthesis Kit	86, 206, 209
Selection Chart	86, 206
Template Switching RT Enzyme Mix	86, 206, 208
WarmStart RTx Reverse Transcriptase	86, 206, 207
cDNA, Primers	208
Cell-free Expression	227-229
Cellular Analysis	
Cellular Imaging & Analysis	292-299
Labeling with SNAP-tag Technology/Troubleshooting Guide	354
Chitin Resin	230, 233, 234
Chitin Magnetic Beads	233, 234
Clal	35
Cleavage Close to the End of DNA Fragments	328
Cleavage of Supercoiled DNA	318
CLIP-Biotin	298
CLIP-Cell 505	296
CLIP-Cell Block	297
CLIP-Cell TMR-Star	296
CLIP-Surface 488	296
CLIP-Surface 547	296
CLIP-Surface 647	296
CLIP-tag	
Overview	294
Comparison to GFP	295
Vectors	298
Cloning	
BioBrick Assembly Kit	97
Cloning Plasmids	128
Competent Cells for Cloning	129, 242-244
Enhancing Transformation Efficiency	350
Getting Started with Molecular Cloning	331
Gibson Assembly Cloning Kit	95
NEB Golden Gate Assembly Kit (BsaI-HFv2)	96
NEB Golden Gate Assembly Kit (BsmBI-v2)	96
NEB PCR Cloning Kit	97
NEBuilder HiFi DNA Assembly Cloning Kit	94
Optimization Tips	339-340
Q5 Site Directed Mutagenesis Kit	98
Traditional Cloning Quick Guide	334-335
Troubleshooting Guide	336-338
Workflow Comparison	92-93
Colony PCR	
One Taq DNA Polymerase	65, 68
Taq DNA Polymerase	65, 69
Color Prestained Protein Standard, Broad Range (10–250 kDa)	185
Combinatorial Libraries	272-273
Compatible Cohesive Ends	319-320
Competent Cells	129, 232, 238-249
Competent Cell Strain Properties Chart	240-241
Competent Cell Cloning Selection Chart	129
Competent <i>E. coli</i>	
BL21	232, 240, 241, 245
BL21(DE3)	232, 240, 241, 245
dam ⁻ / dcm ⁻	129, 240, 241, 244

Enhancing Transformation Efficiency	350
Lemo21(DE3)	232, 240, 241, 245
NEB 5-alpha (High Efficiency)	129, 240, 241, 243
NEB 5-alpha (Subcloning Efficiency)	129, 240, 241, 243
NEB 5-alpha Electrocompetent	129, 240, 241, 243
NEB 5-alpha F ['] /h (High Efficiency)	129, 240, 241, 243
NEB 10-beta (High Efficiency)	129, 240, 241, 243
NEB 10-beta Electrocompetent	129, 240, 241, 243
NEB Express (High Efficiency)	232, 240, 241, 246
NEB Express /h (High Efficiency)	232, 240, 241, 246
NEB Stable Competent <i>E. coli</i>	129, 240, 241, 244
NEB Turbo (High Efficiency)	129, 240, 241, 242
NiCo21(DE3)	232, 240, 241, 246
Protein Expression with T7 Express Strains	351
SHuffle Express	232, 240, 241, 248
SHuffle T7	232, 240, 241, 249
SHuffle T7 Express	232, 240, 241, 248
SHuffle T7 Express /ysY	232, 240, 241, 249
T7 Express (High Efficiency)	232, 240, 241, 247
T7 Express /ysY (High Efficiency)	232, 240, 241, 247
T7 Express /ysY/h (High Efficiency)	232, 240, 241, 247
Competitor Cross Reference Tool	304
COVID-19 research, supporting	6
CpG Methyltransferase (M.SssI)	124, 286
Cre Recombinase	123
Creative Action Institute	8
CRISPR Products	190
CspCI	35
Customized Solutions	4
CutSmart Buffer	20, 58, 309
CviAll	35
CviKI-1	35
CviQI	35

D

5' Deadenylase	214
3' Desthiobiotin-GTP	205
5' DNA Adenylation Kit	211, 214
7-deaza-dGTP	85
dam ⁻ /dcm ⁻ Competent <i>E. coli</i>	232, 240, 241, 244
Dam, Dcm and CpG Methylation	321
dam Methyltransferase	125, 286
dATP Solution	85
DdeI	36
Deep Vent DNA Polymerase	71
Deep Vent (exo ⁻) DNA Polymerase	71
Deoxynucleotide (dNTP) Solution Set	85
Deoxynucleotide (dNTP) Solution Mix	64, 85
Dephosphorylation	105-106
dGTP Solution	85
Diluent Buffers (Restriction Endonuclease)	59, 309
Distributors	inside back cover
DNA	
Adenylation	214
Assembly	94-96
Base Pair Diagrams	367
Cleanup, Monarch Products for	129, 139
Cloning Plasmids	128
Enrichment	167, 220, 279
Extraction, Monarch Products for	129, 136, 140
Fragmentation	149
Input Guidelines for Library Prep	353
Ladders	180-183
Library Preparation	144-175
Ligases	99-104
Loading Dye	59
Maps	355-365
Markers	183
Methylation Analysis	151, 280-281, 289, 321
Methyltransferases	124-126, 286
Modifying Enzymes	90-129
Phosphatases	105-106

Plasmid Miniprep, Monarch Products for	140
Polymerases	65-71
Reagents for Library Preparation	144-175
Substrates	128
DNA Markers and Ladders	
φX174 DNA-HaeIII Digest	183
50 bp DNA Ladder	181
1 kb Plus DNA Ladder	181
1 kb Plus DNA Ladder for Safe Stains	182
100 bp DNA Ladder	181
1 kb DNA Ladder	181
Fast DNA Ladder	182
Lambda DNA-BstEII Digest	183
Lambda DNA-HindIII Digest	183
Lambda DNA-Mono Cut Mix	182
Lambda PFG Ladder	182
Low Molecular Weight DNA Ladder	181
Midrange PFG Ladder	182
pBR322 DNA-BstNI Digest	183
pBR322 DNA-MspI Digest	183
PCR Marker	181
Quick-Load 1 kb Plus DNA Ladder	182
Quick-Load 1 kb DNA Ladder	182
Quick-Load 1 kb Extend DNA Ladder	182
Quick-Load 100 bp DNA Ladder	182
Quick-Load Purple 1 kb Plus DNA Ladder	182
Quick-Load Purple 1 kb DNA Ladder	182
Quick-Load Purple 100 bp DNA Ladder	182
Quick-Load Purple 50 bp DNA Ladder	182
Quick-Load Purple Low Molecular Weight DNA Ladder	182
Supercoiled DNA Ladder	184
Tridye 1 kb Plus DNA Ladder	182
Tridye 1 kb DNA Ladder	182
Tridye 100 bp DNA Ladder	182
Tridye Ultra Low Range DNA Ladder	182
DNA Modifying Enzymes, Activity in CutSmart Buffer	315
DNA Modifying Enzymes and Cloning Technologies	90-129
DNA Polymerase I	83
DNA Polymerase I, Large (Klenow) Fragment	83
DNA Repair	
Enzymes on Damaged & Non-standard Bases	117
NEBNext FFPE DNA Repair Mix	152
PreCR Repair Mix	82
Properties of DNA Repair Enzymes & Structure Specific Endonucleases	116
Enzymes	116-122
DNA Library Preparation Reagents for Next Generation Sequencing	144-175
DNA Sequences and Maps Tool	304
DNase I (RNase-Free)	110, 217
Dnmt1, Methyltransferase	126, 286
Double Digestion	307
Double Digest Finder	304
DpnI	36, 281
DpnII	36, 281
DraI	36
DraIII-HF	36
DrdI	36
dsRNA Ladder	184, 221
dsRNA Ligase (T4 RNA Ligase 2)	212
dUTP Solution	85

E

EaeI	36
EagI-HF	37
EarI	37
Ecil	37
<i>E. coli</i> DNA Ligase	102
<i>E. coli</i> Ribosome	228
<i>E. coli</i> RNA Polymerase, Core Enzyme	203
<i>E. coli</i> RNA Polymerase, Holoenzyme	203
<i>E. coli</i> Poly(A) Polymerase	202
Eco53KI	37
EcoNI	37
EcoO109I	37

EcoP15I	37
EcoRI	37
EcoRI-HF	38
EcoGII Methyltransferase	125, 286
EcoRI Methyltransferase	125, 286
EcoRV	38
EcoRV-HF	38
EGCasel	265
ElectroLigase	99, 101
EM-Seq	151, 289
Endoglycosidase I	265
Endoglycosidases	
Endo D	255
Endo F2	254
Endo F3	254
Endo H, Endo H _p	255
Endo S	255
O-Glycosidase	258
PNGase A	256
PNGase F (Glycerol Free, Recombinant)	256
Rapid PNGase F (non-reducing format)	257
PNGase F, Remove-iT	257
Rapid PNGase F	257
Endonucleases	
APE I	118
Common Applications	109
DNase I (RNase-Free)	217
Endonuclease III (Nth)	119
Endonuclease IV	118
Endonuclease V	119
Endonuclease VIII	120
Fpg	120
hAAG	120
Homing	57-58
Mung Bean Nuclease	113
Nicking	55-57
Properties of	108
Restriction Endonucleases	18-59
T4 PDG (T4 Endonuclease V)	121
T7 Endonuclease I	122
Thermostable FEN1	112
Tma Endonuclease III	119
Tth Endonuclease IV	118
Endoproteinase AspN	267
Endoproteinase GluC	267
Endoproteinase LysC	267
EnGen Lba Cas12a (Cpf1)	190, 191
EnGen Mutation Detection Kit	190, 192
EnGen Sau Cas9	190, 191
EnGen sgRNA Synthesis Kit, <i>S. pyogenes</i>	190, 193
EnGen sgRNA Template Oligo Designer Tool	304
EnGen Spy Cas9 Nickase	190, 191
EnGen Spy Cas9 NLS	190, 191
EnGen Spy dCas9 (SNAP-tag)	190, 191
Enhancing Transformation Efficiency	350
Enterokinase, light chain	269
Environmental Commitment	2
Enzymatic Methyl-Seq	151, 289
Enzyme Finder	304
Enzymes for Innovation	4
Epigenetics	276-289
EpiMark 5-hmC and 5-mC Analysis Kit	278
EpiMark Bisulfite Conversion Kit	279
EpiMark Hot Start Taq DNA Polymerase	71, 279
EpiMark Methylated DNA Enrichment Kit	279
EpiMark N6-Methyladenosine Enrichment Kit	220, 279
EpiMark Nucleosome Assembly Kit	283
<i>E. coli</i> DNA Polymerase I	83
Esp3I	38
ET SSB (Extreme Thermostable ssBinding Protein)	127
Exo-CIP Rapid PCR Cleanup Kit	87
Exoglycosidases	
α -N-Acetylgalactosaminidase	259
α 1-2 Fucosidase	259
α 1-2,3,4,6 Fucosidase	260

α 1-2,4,6 Fucosidase O	260
α 1-3,4 Fucosidase	260
α 1-3,6 Galactosidase	260
α 1-3,4,6 Galactosidase	261
α 1-2,3 Mannosidase	262
α 1-2,3,6 Mannosidase	262
α 1-6 Mannosidase	262
α 2-3 Neuraminidase S	263
α 2-3,6,8 Neuraminidase	262
α 2-3,6,8,9 Neuraminidase A	263
β -N-Acetylglucosaminidase S	259
β -N-Acetylhexosaminidase	259
β 1-3 Galactosidase	261
β 1-3,4 Galactosidase	261
β 1-4 Galactosidase S	261
Exonucleases	
Common Applications	109
Exonuclease I (<i>E. coli</i>)	110
Exonuclease III (<i>E. coli</i>)	111
Exonuclease T	112, 217
Exonuclease V (RecBCD)	111
Exonuclease VII	111
Exonuclease VIII, truncated	112
Lambda Exonuclease	110
Nuclease BAL-31	113
Properties of	108
RecJ	114
T7 Exonuclease	115
Exo Selector	304

F

Factor Xa Protease	269
Fast DNA Ladder	182
FatI	38
FatII	38
Fetuin	258
Fluorescent LAMP	77, 78
Fluorescent Substrates for Protein Labeling	296
Fnu4HI	38
FokI	38
Fpg	120
Freezer Program Locator, NEBNow	305
Freezer Program, On-site	4
FseI	39
FspI	39
FspEI	280
Fucosidase	259, 260
Furin	269

G

Gel Extraction, Monarch Products for	140
Gel Loading Dye, Blue (6X)	59
Gel Loading Dye, Orange (6X)	59
Gel Loading Dye, Purple (6X)	59, 180
Gel Loading Dye, Purple (6X), No SDS	59, 180
Gene Assembly	94-96
General Guidelines for PCR Optimization	322
General Molecular Biology Data	366-369
Genetic Code	366
Genetic Markers	348-349
Genomic DNA Substrates	
HeLa Genomic DNA	287
Genome Editing	188-193
Genomic DNA Purification, Monarch Kits for	135-136
Genotyping, NEBNext Direct	165
Gibson Assembly Cloning Kit	95
Gibson Assembly Master Mix	95
Gibson Assembly Tips for Optimization	331
Global Development	4
Globin RNA & rRNA Depletion	156
GluC, Endoproteinase	267

INDEX

Glycan Analyzer	304
Glycobiology and Protein Tools	252-273
Glycosidases	
Endoglycosidases	254-258
Exoglycosidases	259-263
Glycosylases	
<i>Afu</i> Uracil-DNA Glycosylase (UDG)	121
Endonuclease III (Nth)	119
Endonuclease IV	118
Endonuclease V	119
Endonuclease VIII	120
Fpg	120
hAAG	120
T4 PDG	121
<i>Tma</i> Endonuclease III	119
<i>Tth</i> Endonuclease IV	118
Uracil-DNA Glycosylase (UDG)	121
Glucosylation (T4-BGT)	278
GMP Capabilities	7
Goat Anti-Mouse IgG Magnetic Beads	233, 235
Goat Anti-Rabbit IgG Magnetic Beads	233, 235
Goat Anti-Rat IgG Magnetic Beads	233, 235
Golden Gate Assembly Kit (Bsal-HFv2)	96
Golden Gate Assembly Kit (BsmBI-v2)	96
Golden Gate Assembly, Tips	332
GpC Methyltransferase (M.CviPI)	124, 286

H

hAAG	120
HaeII	39
HaeIII	39
HaeIII Methyltransferase	125, 286
Half Life of Radioisotopes	367
Heat Inactivation of Restriction Endonucleases	309-314
HeLa Genomic DNA	287
Helper Phage, M13K07	128
Hemo Klen <i>Taq</i>	64, 65, 70
Heparin Lyases	
<i>Bacteriodes</i> Heparinase I	264
<i>Bacteriodes</i> Heparinase II	264
<i>Bacteriodes</i> Heparinase III	264
Hgal	39
Hhal	39
Hhal Methyltransferase	125, 286
HiFi <i>Taq</i> DNA Ligase	103
High-Fidelity PCR	
Deep Vent DNA Polymerase	65, 71
Q5 High-Fidelity DNA Polymerase	64, 65, 66
Q5 Hot Start High-Fidelity DNA Polymerase	64, 65, 66
Phusion High-Fidelity DNA Polymerase	65, 67
Phusion Hot Start Flex	65, 67
Vent DNA Polymerase	65, 71
High-Fidelity (HF) Restriction Enzymes	21, 316
HincII	39
HindIII	39
HindIII-HF	40
Hinfl	40
HinPI	40
HiScribe SP6 RNA Synthesis Kit	199, 201
HiScribe T7 ARCA mRNA Kit (with or without tailing)	190, 200, 201
HiScribe T7 High Yield RNA Synthesis Kit	190, 199, 201
HiScribe T7 Quick High Yield RNA Synthesis Kit	190, 199, 201
Histones	282-285
EpiMark Nucleosome Assembly Kit	283
H1 ⁰ Human, Recombinant	283
H2A Human, Recombinant	284
H2A/H2B Dimer	283
H2B Human, Recombinant	284
Histone H3.1 Human, Recombinant	284
H3.1/H4 Tetramer	283
Histone H3.2 Human, Recombinant	285

Histone H3.3 Human, Recombinant	285
Histone H4 Human, Recombinant	285
Hi-T4 DNA Ligase	99, 101
Hi-T7 RNA Polymerase (std. & high conc.)	202
HMW DNA Extraction, Monarch Kits for	136
Homing Endonucleases	
I-CeuI	57
I-SceI	57
PI-PsPI	58
PI-SceI	58
Hot Start PCR	
EpiMark Hot Start DNA Polymerase	71
Hot Start <i>Taq</i> DNA Polymerase	64, 65, 69
Hot Start <i>Taq</i> 2X Master Mix	64, 69
LongAmp Hot Start DNA Polymerase	65, 70
LongAmp Hot Start <i>Taq</i> 2X Master Mix	70
One <i>Taq</i> Hot Start DNA Polymerase	65, 68
One <i>Taq</i> Hot Start 2X Master Mixes	68
Phusion Hot Start Flex DNA Polymerase	67
Phusion Hot Start Flex 2X Master Mix	67
Q5 Hot Start High-Fidelity DNA Polymerase	65, 66
Q5 Hot Start High-Fidelity 2X Master Mix	66
HpaI	40
HpaII	40, 281
HpaII Methyltransferase	126, 286
HphI	40
Hpy99I	40
Hpy166II	40
Hpy188I	41
Hpy188III	41
HpyAV	41
HpyCH4III	41
HpyCH4IV	41
HpyCH4V	41
hSMUG1	122
Human Alkyladenine Glycosylase (hAAG)	120
Human DNA (cytosine-5) MTase (Dnmt1)	126, 286
Human RNase Inhibitor	218
Hydrophilic Streptavidin Magnetic Beads	233, 285
Hydroxymethylcytosine, Identification	278

I

I-CeuI	57
Icon Descriptions	18, 309
IdeZ Protease (IgG-specific)	265
Immobilized T4 DNA Ligase	99, 101
Immune Sequencing, NEBNext	166
IMPACT Systems	
IMPACT Kit	230
IMPACT Vectors	230
IMPACT Vectors Comparison	230
Inhibitors, RNA	217, 218
Inorganic Pyrophosphatase (<i>E. coli</i>)	107, 203
Inorganic Pyrophosphatase, Thermotable	107, 203
Inorganic Pyrophosphatase, (yeast)	107, 203
Instant Sticky-end Ligase Master Mix	99, 100, 101
Inteins	230
Interactive Tools	304-305
Internet Technical Support	302
<i>In vitro</i> Protein Synthesis, PURExpress	228
<i>In vitro</i> Transcription, Protein	228
I-SceI	57
ISO 9001:2000, ISO 13485, ISO 14001	2
IsoAmp II Universal tHDA Kit	64, 80
Isoschizomers	317
Isothermal Amplification	76-80
Isotope Data	367

K

KasI	41
Kinases	
cAMP-dependent Protein Kinase (PKA)	271
Casein Kinase II (CK2)	271
T4 Polynucleotide Kinase	104
T4 Polynucleotide Kinase (3' phosphatase minus)	104
<i>K. lactis</i> Protein Expression	
<i>K. lactis</i> Competent Cells	231
<i>K. lactis</i> Protein Expression Kit	231
<i>K. lactis</i> Vectors	231
Klenow Fragment	83
Klenow Fragment (3' → 5' exo ⁻)	83
KpnI	41
KpnI-HF	42

L

Labeling with SNAP-tag/Troubleshooting	354
Ladders	
DNA	180-182
RNA	184, 221
Lambda	
BstEII Digest	183
DNA	128
DNA (N ⁶ -methyladenine-free)	128
Exonuclease (I, T, III, λ)	109-113
HindIII Digest	183
Map of	355
Mono Cut Mix	182
PFG Ladders	182
Protein Phosphatase (Lambda PP)	270
LAMP	76-78
LAMP Fluorescent Dye	77
Lemo21(DE3) Competent <i>E. coli</i>	232, 240, 241, 245
Lesion Bypass Polymerase, <i>Sulfolobus</i> DNA Polymerase IV	82
Library Preparation, DNA	144-175
Adaptors and Primers	160
DNA Enrichment, Microbiome Kit	167
DNA Input Guidelines	353
DNA Repair	152
Fragmentation	149, 154
Illumina Compatible	146-160
Ion Torrent Compatible	171
Methylated – DNA	151, 289
Modules	158-160
Oxford Nanopore Technologies	162
Target Enrichment	149-152, 158-159
Library Preparation, RNA	144-175
Illumina Compatible	137, 142-148, 160-161
Modules	150, 154, 159, 160-161, 168-169
rRNA Depletion	155-157
Library Quantitation	170
Ligases	
9°N DNA Ligase	104
<i>E. coli</i> DNA Ligase	102
ElectroLigase	99, 101
HiFi <i>Taq</i> DNA Ligase	103
Hi-T4 DNA Ligase	99, 101
Immobilized T4 DNA Ligase	99, 101
NEBNext Ultra II Ligation Module	145, 146, 150, 154
RtcB Ligase	211, 213
Salt-T4 DNA Ligase	99, 101
SplintR Ligase	104, 211, 214
T3 DNA Ligase	102
T4 DNA Ligase	99, 100, 101
T7 DNA Ligase	102
T4 RNA Ligase 1	210, 211, 212
T4 RNA Ligase 2	210, 211, 212
T4 RNA Ligase 2, truncated	210, 211, 212
T4 RNA Ligase 2, truncated K227Q	210, 211, 213
T4 RNA Ligase 2, truncated KQ	210, 211, 213

<i>Taq</i> DNA Ligase	103
Thermostable 5' App DNA/RNA	211, 213
Ligase Fidelity Tools	304
Ligase Master Mixes	
Blunt/TA Ligase Master Mix	99, 100, 101
Instant Sticky-end Ligase Master Mix	99, 100, 101
Ligase Selection Chart, DNA	99
Ligase Selection Chart, RNA	211
Ligation	
DNA	99-104
Protein	230
RNA	210-214
Ligation Kit, Quick	99, 100, 101
Loading Dye, Blue	59
Loading Dye, Orange	59
Loading Dye, Purple	59
LongAmp Hot Start <i>Taq</i> 2X (Polymerase, Master Mix)	65, 70
LongAmp <i>Taq</i> (Polymerase, Master Mix, Kit)	65, 70
Long PCR	
LongAmp <i>Taq</i>	65, 70
LongAmp Hot Start <i>Taq</i>	65, 70
Phusion High-Fidelity DNA Polymerase	65, 67
Phusion Hot Start Flex DNA Polymerase	65, 67
Q5 High-Fidelity DNA Polymerase	65, 66
Q5 Hot Start High-Fidelity DNA Polymerase	65, 66
Low Range ssRNA Ladder	184, 221
Low Molecular Weight DNA Ladder	181
LpnPI	280
Luna Cell Ready Lysis Module	75
Luna Cell Ready One-Step RT-qPCR Kit	75, 210
Luna Cell Ready Probe One-Step RT-qPCR Kit	75, 210
Luna Probe One-Step RT-pPCR 4X Mix with UDG	64, 74, 210
Luna Probe One-Step RT-qPCR Kit (No ROX)	64, 74, 210
Luna SARS-CoV-2 RT-qPCR Multiplex Assay Kit	64, 77, 210
Luna Universal One-Step RT-qPCR Kit	64, 74, 210
Luna Universal Probe One-Step RT-qPCR Kit	64, 74, 210
Luna Universal Probe qPCR Master Mix	64, 73
Luna Universal qPCR Master Mix	64, 73
Luna qPCR & RT-qPCR Optimization Tips	324-325
Luna qPCR & RT-qPCR Troubleshooting Guide	326-327
LunaScript RT Master Mix Kit (Primer-free)	73
LunaScript RT SuperMix Kit	64, 73

M

5-methyl-dCTP	85, 287
M13K07 Helper Phage	128
M13mp18 DNA	128
M13mp18, Map of	356
Magnetic Matrices (beads)	
Amylose	233, 234
Anti-MBP	233, 234
Chitin	233, 234
Goat Anti-Mouse	233, 235
Goat Anti-Rabbit	233, 235
Goat Anti-Rat	233, 235
Hydrophilic Streptavidin	233, 235
Oligo (dT) ₂₅ Magnetic Beads	219, 233, 235
Protein A	233, 235
Protein G	233, 235
SNAP-Capture	299
Streptavidin	219, 233, 235
Magnetic mRNA Isolation Kit	219, 235
Magnetic Protein Purification	233-235
Magnetic Separation Racks	
6-Tube Magnetic Separation Rack	235
12-Tube Magnetic Separation Rack	235
50 ml Magnetic Separation Rack	235
96-Well Microtiter Plate Magnetic Separation Rack	235
NEBNext Magnetic Separation Rack	159, 235
Maps	
Lambda	355
M13mp18	356

pBR322	357	MluI-HF	42
pKLAC2	358	MluCI	42
pMAL-c6T	359	MlyI	43
pMiniT 2.0	360	MmeI	43
pNEB206A	361	M-MuLV Reverse Transcriptase	86, 206, 207
pSNAP	362	MnII	43
pTXB1	363	Modifying Enzymes, Activity in Cutsmart Buffer	315
pTYB21	364	Molecular Diagnostics	7, 64
pUC19	365	Monarch 2 ml Tubes	136
Maltose Binding Protein (MBP) Purification	233, 234	Monarch Bead Retainers	136
Markers and Ladders (DNA, RNA, Protein)	178-185	Monarch Collection Tubes II	135, 137, 138
Markers, Molecular Weight		Monarch DNA Capture Beads	136
DNA	180-182, 184	Monarch DNA Cleanup Binding Buffer	139
RNA	184, 221	Monarch DNA Cleanup Columns (5 µg)	139, 140
Protein	185	Monarch DNA Elution Buffer	139, 140
Mass Spectrometry Standards, Trypsin Digested BSA	267	Monarch DNA Gel Extraction Kit	140
Master Mixes, PCR, qPCR, LAMP & RT-LAMP		Monarch DNA/RNA Protection Reagent	137
Hot Start <i>Taq</i> 2X Master Mix	64, 69	Monarch DNA Wash Buffer	139, 140
LongAmp Hot Start <i>Taq</i> 2X Master Mix	70	Monarch Gel Dissolving Buffer	140
LongAmp <i>Taq</i> 2X Master Mix	70	Monarch gDNA Binding Buffer	135
Luna Probe One-Step RT-qPCR 4X Mix with UDG	64, 74	Monarch gDNA Blood Lysis Buffer	135
LunaScript RT Master Mix Kit (Primer-free)	73	Monarch gDNA Cell Lysis Buffer	135
Luna Universal Probe qPCR Master Mix	64, 74	Monarch gDNA Elution Buffer	135
Luna Universal qPCR Master Mix	64, 73	Monarch gDNA Elution Buffer II	136
Multiplex PCR 5X Master Mix	69	Monarch gDNA Kit Sample Input Amounts	342
NEBNext Ultra II Q5 Master Mix	146, 147, 150, 169, 172-174, 288	Monarch gDNA Nuclei Prep & Lysis Buffer Pack	136
One <i>Taq</i> 2X Master Mix with GC Buffer	68	Monarch gDNA Purification Columns	135
One <i>Taq</i> 2X Master Mix with Standard Buffer	68	Monarch gDNA Removal Column	137
One <i>Taq</i> Hot Start 2X Master Mix with GC Buffer	68	Monarch gDNA Tissue Lysis Buffer	135
One <i>Taq</i> Hot Start 2X Master Mix with Standard Buffer	68	Monarch gDNA Troubleshooting Guide	343
One <i>Taq</i> Hot Start Quick-Load 2X Master Mix with GC Buffer	68	Monarch gDNA Wash Buffer	135
One <i>Taq</i> Hot Start Quick-Load 2X Master Mix with Standard Buffer	68	Monarch Genomic DNA Purification Kit	135, 162
One <i>Taq</i> Quick-Load 2X Master Mix with GC Buffer	68	Monarch HMW DNA Extraction Kit for Cells & Blood	136, 162
One <i>Taq</i> Quick-Load 2X Master Mix with Standard Buffer	68	Monarch HMW DNA Extraction Kit for Tissue	136, 162
Phusion High-Fidelity PCR Master Mix with GC Buffer	67	Monarch HMW DNA Tissue Lysis Buffer	135
Phusion High-Fidelity PCR Master Mix with HF Buffer	67	Monarch HMW DNA, Sample Input Amounts	346-347
Phusion Hot Start Flex 2X Master Mix	67	Monarch Microfuge Tube EcoRack	141
Q5 High-Fidelity 2X Master Mix	64, 66	Monarch Nucleic Acid Purification Kits	129, 132-141
Q5 Hot Start High-Fidelity 2X Master Mix	64, 66	Monarch PCR & DNA Cleanup Kit (5 µg)	87, 139
Quick-Load <i>Taq</i> 2X Master Mix	69	Monarch Pestle Set	136
<i>Taq</i> 2X Master Mix	69	Monarch Plasmid Lysis Buffer (B2)	141
<i>Taq</i> 5X Master Mix	69	Monarch Plasmid Miniprep Columns	141
WarmStart Colorimetric LAMP/RT-LAMP 2X Master Mix	64, 78	Monarch Plasmid Miniprep Kit	141
WarmStart Colorimetric LAMP/RT-LAMP 2X Master Mix (with UDG)	64, 78	Monarch Plasmid Neutralization Buffer (B3)	141
WarmStart Multi-Purpose LAMP/RT-LAMP 2X Master Mix (with UDG)	64, 77	Monarch Plasmid Resuspension Buffer (B1)	141
Mbol	42	Monarch Plasmid Wash Buffer 1	141
MbolI	42	Monarch Plasmid Wash Buffer 2	141
McrBC	281	Monarch Precipitation Enhancer	136
Messenger RNA Modifications	368	Monarch Protein Separation Solution	136
Methyl-seq	151, 289	Monarch RBC Lysis Buffer	136
Methylation-Dependent Restriction Enzymes	280-281	Monarch RNA Cleanup Binding Buffer	138
Methylation Sensitivity (dam, dcm, CpG)	309-314, 321	Monarch RNA Cleanup Columns (50 µg, 500 µg)	138
Methyltransferases		Monarch RNA Cleanup Kit (10 µg, 50 µg, 500 µg)	138, 190, 193, 220
AluI	124	Monarch RNA Cleanup Wash Buffer	138
BamHI	124	Monarch RNA Lysis Buffer	137
CpG (M.SssI)	124	Monarch RNA Priming Buffer	137
dam	125	Monarch RNA Purification Column	137
EcoGII	125	Monarch RNase A	135
EcoRI	125	Monarch Total RNA Miniprep Enzyme Pack	137
GpC (M.CviPI)	124	Monarch Total RNA Miniprep Kit	137, 190, 193, 220
HaeIII	125	Monarch RNA Wash Buffer	137
HhaI	125	Monarch Total RNA Kit Sample Input Amounts	344
HpaII	126	Monarch Troubleshooting Guide for DNA Cleanup & Plasmid Purification	341
Human DNA (cytosine-5) MTase (Dnmt1)	286	Monarch Troubleshooting Guide for gDNA Purification	343
MspI	126	Monarch Troubleshooting Guide for Total RNA Extraction & Purification	345
<i>Taq</i> I	126	mRNA Cap 2'-O-Methyltransferase	204
MfeI	42	mRNA Decapping Enzyme	205
MfeI-HF	42	mRNA Isolation	219, 220
Microbiome Enrichment	167, 173	MscI	43
Micrococcal Nuclease	113	MseI	43
microRNA Marker	184, 221	MslI	43
Midrange PFG Ladder	182	MspI	43, 281
miRNA Cloning Linker	221	MspI Methyltransferase	126, 286
MluI	42	MspA1I	43

MspJI	280
Multiplex PCR	
Hemo Klen <i>Taq</i>	70
Multiplex PCR 5X Master Mix	69
One <i>Taq</i> DNA Polymerase	68
Phusion High-Fidelity DNA Polymerase	67
Q5 High-Fidelity DNA Polymerase	66
Mung Bean Nuclease	113
Murine RNase Inhibitor	217
Mutation Detection	190, 192
Mwol	44

N

9°N DNA Ligase	104
N6-Methyladenosine Enrichment	220
NaeI	44
NarI	44
Nb.BbvCI	55
Nb.BsmI	55
Nb.BsrDI	56
Nb.BssSI	56
Nb.BtsI	56
NciI	44
NcoI	44
NcoI-HF	44
NdeI	44
NEB 5-alpha (High Efficiency) Competent <i>E. coli</i>	129, 240, 241, 243
NEB 5-alpha (Subcloning Efficiency) Competent <i>E. coli</i>	129, 240, 241, 243
NEB 5-alpha Electrocompetent <i>E. coli</i>	129, 240, 241, 243
NEB 5-alpha F ⁺ (High Efficiency) Competent <i>E. coli</i>	129, 240, 241, 243
NEB 10-beta (High Efficiency) Competent <i>E. coli</i>	129, 240, 241, 243
NEB 10-beta Electrocompetent <i>E. coli</i>	129, 240, 241, 243
NEB 10-beta/Stable Outgrowth Medium	243
NEB Cloning Competent <i>E. coli</i> Sampler	242
NEB Custom RNA Depletion Design Tool	304
NEB Golden Gate Assembly Kit (Bsal-HFv2)	96
NEB Golden Gate Assembly Kit (BsmBI-v2)	96
NEB Golden Gate Assembly Tool	304
NEB LAMP Primer	304
NEB PCR Cloning Kit	97
NEB Online Tools	304-305
NEB Stable Competent <i>E. coli</i>	129, 240, 241, 244
NEB Tool for Augmented Reality	305
NEB Tube Opener	58, 244
NEB Turbo (High Efficiency) Competent <i>E. coli</i>	129, 240, 241, 242
NEBaseChanger	304
NEBcloner	304
NEBcutter	304
NEBExpress Cell-free <i>E. coli</i> Protein Synthesis System	227
NEBExpress (High Efficiency) Competent <i>E. coli</i>	232, 240, 241, 246
NEBExpress ⁺ (High Efficiency) Competent <i>E. coli</i>	232, 240, 241, 246
NEBExpress GamS Nuclease Inhibitor	227
NEBExpress MBP Fusion and Purification System	229
NEBExpress Ni-NTA Magnetic Beads	233, 234
NEBExpress Ni Resin	233, 234
NEBExpress Ni Spin Columns	233, 234
NEBioCalculator	304
NEBNext Adaptors & Primers for Illumina	146, 147, 160, 172, 174
NEBNext Adaptor Dilution Buffer	146, 147, 160, 172, 174
NEBNext ARTIC SARS-CoV-2 Companion Kit (Oxford Nanopore Technologies)	161, 162, 175
NEBNext ARTIC SARS-CoV-2 FS Library Prep Kit (Illumina)	161, 174
NEBNext ARTIC SARS-CoV-2 Library Prep Kit (Illumina)	161, 174
NEBNext Companion Module for Oxford Nanopore Technologies Ligation Sequencing	162, 173
NEBNext Custom Ready Panels	164, 172
NEBNext Custom RNA Depletion Design Tool	304
NEBNext dA-Tailing Module	172
NEBNext Direct Genotyping Solution	165, 172
NEBNext DNA Library Prep, Troubleshooting	352
NEBNext dsDNA Fragmentase	146, 168, 173
NEBNext dsDNA Fragmentase Reaction Buffer v2	168, 173
NEBNext End Repair Module	172
NEBNext Enzymatic Methyl-seq Kit	146, 151, 172, 289

NEBNext Enzymatic Methyl-seq Conversion Module	146, 151, 172, 289
NEBNext FFPE DNA Repair Mix	152, 173
NEBNext Globin RNA & rRNA Depletion Kit	147, 156, 174, 175
NEBNext Globin & rRNA Depletion Kit w/RNA Sample Purif. Beads	147, 156, 174, 175
NEBNext High-Fidelity 2X PCR Master Mix	169, 172
NEBNext Immune Sequencing Kit (Human)	166, 172
NEBNext Immune Sequencing Kit (Mouse)	166, 172
NEBNext Library Dilution Buffer	170, 172, 175
NEBNext Library Quant Kit for Illumina	170, 172
NEBNext Magnetic Separation Rack	159, 173, 175
NEBNext Microbiome DNA Enrichment Kit	167, 173
NEBNext Magnesium RNA Fragmentation Module	154, 174, 175
NEBNext Multiplex Oligos (Dual Index Primers Set 1)	160, 172, 174
NEBNext Multiplex Oligos (Dual Index Primers Set 2)	160, 172, 174
NEBNext Multiplex Oligos (Methylated Adaptor, Index Primers Set 1)	160, 172
NEBNext Multiplex Oligos for Enzymatic Methyl-seq (Unique Dual Index Primer Pairs)	151, 172, 289
NEBNext Multiplex Oligos for Illumina (96 Index Primers)	160, 172, 174
NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs)	160, 172, 174
NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs Set 2)	160, 172, 174
NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs Set 3)	160, 172, 174
NEBNext Multiplex Oligos for Illumina (96 Unique Dual Index Primer Pairs Set 4)	160, 172, 174
NEBNext Multiplex Oligos for Illumina (Index Primers Set 1, 2, 3, 4)	160, 172, 174
NEBNext Multiplex Oligos for Illumina (Unique Dual Index UMI Adaptors DNA Set 1)	160, 172
NEBNext Multiplex Oligos for Illumina (Unique Dual Index UMI Adaptors RNA Set 1)	160, 174
NEBNext Multiplex Small RNA Library Prep Kit for Illumina (Index Primers 1-48)	159, 174
NEBNext Multiplex Small RNA Library Prep Set for Illumina (Set 1)	159, 174
NEBNext Multiplex Small RNA Library Prep Set for Illumina (Set 2)	159, 174
NEBNext Poly(A) mRNA Magnetic Isolation Module	154, 174
NEBNext Q5U Master Mix	151, 172, 289
NEBNext Quick Ligation Module	172
NEBNext Quick Ligation Reaction Buffer	172
NEBNext Reagents for Next Generation Sequencing	
Illumina Compatible	146-160, 161, 163-170
Ion Torrent	152
Oxford Nanopore Technologies	162
Selection Charts	146-147
NEBNext RNA Depletion Core Reagent Set	
(with or without RNA Sample Purification Beads)	154, 157, 174
NEBNext RNA First Strand Synthesis Module	154, 174
NEBNext rRNA Depletion Kit (Bacteria)	
(with or without RNA Sample Purification Beads)	155, 174
NEBNext rRNA Depletion Kit (Human/Mouse/Rat)	
(with or without RNA Sample Purification Beads)	174
NEBNext rRNA Depletion Kit v2 (Human/Mouse/Rat)	
(with or without RNA Sample Purification Beads)	154, 155, 174
NEBNext Selection Charts	146-147
NEBNext Selector	175, 304
NEBNext Single Cell/Low Input cDNA Synthesis and Amplification Module	158, 174, 218
NEBNext Single Cell/Low Input RNA Library Prep Kit for Illumina	158, 174, 218
NEBNext Small RNA Library Prep Set for Illumina (Multiplex Compatible)	159, 174
NEBNext Ultra Directional RNA Library Prep Kit for Illumina	174
NEBNext Ultra DNA Library Prep Kit for Illumina	172
NEBNext Ultra End Repair/dA-Tailing Module	172, 174
NEBNext Ultra Ligation Module	172, 174
NEBNext Ultra RNA Library Prep Kit for Illumina	153, 174
NEBNext Ultra II Directional RNA Library Prep Kit for Illumina	
(with or without Sample Purification Beads)	147, 153, 174, 218
NEBNext Ultra II Directional RNA Second Strand Synthesis Module	147, 154, 174
NEBNext Ultra II DNA Library Prep Kit for Illumina	
(with or without Sample Purification Beads)	146, 149, 150, 172, 288
NEBNext Ultra II DNA PCR-free Library Prep Kit for Illumina	
(with or without Sample Purification Beads)	149, 150, 172
NEBNext Ultra II End Repair/dA-Tailing Module	146, 150, 154, 172, 174, 288
NEBNext Ultra II FS DNA Library Prep Kit	149, 150, 172
NEBNext Ultra II FS DNA Library Prep with Sample Purification Beads	149, 150, 172
NEBNext Ultra II FS DNA Module	150, 172
NEBNext Ultra II FS DNA PCR-free Library Prep Kit for Illumina	
(with or without Sample Purification Beads)	146, 149, 150, 172
NEBNext Ultra II Ligation Module	146, 150, 154, 172, 174, 288
NEBNext Ultra II Non-directional RNA Second Strand Synthesis Module	147, 154, 174
NEBNext Ultra II Q5 Master Mix	146, 147, 150, 169, 172, 173, 174, 288
NEBNext Ultra II RNA First Strand Synthesis Module	147, 154, 174
NEBNext Ultra II RNA Library Prep Kit for Illumina	147, 153, 154, 174, 218
NEBNext Ultra II RNA Library Prep with Sample Purification Beads	147, 153, 154, 174, 218

NEBnow Locator	305
NEBuffers for Restriction Enzymes	58, 309
NEBuilder Assembly Tool	304
NEBuilder HiFi DNA Assembly Bundle for Large Fragments	94, 95
NEBuilder HiFi DNA Assembly Cloning Kit	94, 95, 190
NEBuilder HiFi DNA Assembly Master Mix	94, 95, 190
NEBuilder Tips for Optimization	333
Neuraminidase	262, 263
New England Biolabs Foundation	8
NgoMIV	44
NGS Library Prep, Guidelines	353
NheI	45
NheI-HF	45
Nicking Endonucleases	
Nb.BbvCI	55
Nb.BsmI	55
Nb.BsrDI	56
Nb.BssSI	56
Nb.BtsI	56
Nt.AlwI	56
Nt.BbvCI	56
Nt.BsmAI	56
Nt.BspQI	56
Nt.BstNBI	56
Nt.CviPII	57
NiCo21(DE3) Competent <i>E. coli</i>	232, 240, 241, 246
Ni-NTA Magnetic Beads	233, 234
NlaIII	45
NlaIV	45
NmeAIII	45
NotI	45
NotI-HF	45
NruI	45
NruI-HF	46
NsiI	46
NsiI-HF	46
NspI	46
Nt.AlwI	56
Nt.BbvCI	56
Nt.BsmAI	56
Nt.BspQI	56
Nt.BstNBI	56, 64
Nt.CviPII	57
Nucleases	
BAL-31 Nuclease	113
Cas9	190, 191
DNase I (RNase-Free)	110, 217
Endonuclease III (Nth)	119
Endonuclease VIII	120
Exonuclease I (<i>E. coli</i>)	110
Exonuclease III (<i>E. coli</i>)	111
Exonucleases and Endonucleases, Properties of	108
Exonuclease T	112, 217
Exonuclease V	111
Exonuclease VII	111
Exonuclease VIII, truncated	112
Homing	57-58
Lambda Exonuclease	110
Micrococcal	113
Mung Bean Nuclease	113
Nicking Endonucleases	55-57
Nuclease P1	114
RecJ	114
Restriction Endonucleases	18-59
T5 Exonuclease	114
T7 Exonuclease	115
Nuclease-free Water	58
Nucleic Acid Data	367
Nucleic Acid Purification	132-141
Nucleic Acid Purification Troubleshooting Guide for DNA Cleanup & Plasmid Prep	341
Nucleic Acid Purification Troubleshooting Guide for gDNA Purification	343
Nucleic Acid Purification Troubleshooting Guide for Total RNA Extraction & Purification	345
Nucleosome Control DNA	283
Nucleoside Digestion Mix	115, 217

Nucleosome Assembly	283
Nucleotides	
7-deaza-dGTP	85
ATP	85
Acyclonucleotide Set	85
dATP/dUTP/dGTP Solutions	85
Deoxynucleotide (dNTP) Solution Set	85
Deoxynucleotide (dNTP) Solution Mix	85
5-methyl-dCTP	85
Physical Properties	369
Ribonucleotide Solution Mix	85, 203
Ribonucleotide Solution Set	85, 203
NudC Pyrophosphatase	107, 203

O

Ocean Genome Legacy Center	8
OEM	4
O-Glycoprotease	265
O-Glycosidase	258
Oligos, NEBNext (Multiplex)	160, 172, 174
Oligo d(T) ₂₅ Magnetic Beads	219, 233, 235
Oligo d(T) ₂₃ VN	208
Online Tools	304-305
OneTaq DNA Polymerase (Enzyme, Master Mix, Kit)	65, 68
OneTaq Hot Start DNA Polymerase (Enzyme, Master Mix)	65, 68
OneTaq One-Step RT-PCR Kit	81, 210
OneTaq RT-PCR Kit	81, 210
Optimization Tips	
Cloning	339-340
Golden Gate Assembly	332
Luna qPCR	324
Luna RT-qPCR	325
NEBuilder HiFi and NEB Gibson Assembly	333
NGS Library Prep	353
PCR	322
Restriction Endonucleases	306
Orange Gel Loading Dye	59
Ordering Information	8

P

Pacl	46
PaeR7I	46
PaqCI	46
Partnering with NEB	4
pBR322	
DNA	128
Map of	357
BstNI Digest	183
MspI Digest	183
PciI	46
pCLIP ₁ Vector	298
PCR	
Cleanup	87, 127
Cloning Kit	97
Deep Vent DNA Polymerase	71
Deep Vent (exo ⁻) DNA Polymerase	71
Deoxynucleotides	85
EpiMark Hot Start Taq	65, 71
Fidelity Estimator	305
Guidelines for Optimization	322
Hemo Klen Taq	64, 65, 70
Hot Start Taq DNA Polymerase	65, 69
Hot Start Taq 2X Master Mix	69
LongAmp Hot Start Taq 2X Master Mix	70
LongAmp Hot Start Taq DNA Polymerase	65, 70
LongAmp Taq DNA Polymerase	65, 70
LongAmp Taq 2X Master Mix	70
LongAmp Taq PCR Kit	70
Marker	181
Multiplex PCR 5X Master Mix	69

One <i>Taq</i> 2X Master Mix with Standard Buffer	68	<i>p</i> -Nitrophenylphosphate (PNPP)	270
One <i>Taq</i> DNA Polymerase	65, 68	pNEB206A, Map of	361
One <i>Taq</i> Hot Start 2X Master Mix with GC Buffer	68	PNase A	256
One <i>Taq</i> Hot Start 2X Master Mix with Standard Buffer	68	PNase F	256
One <i>Taq</i> Hot Start DNA Polymerase	65, 68	PNase F (Glycerol Free, Recombinant, Rapid, non-reducing)	256-257
One <i>Taq</i> Hot Start Quick-Load 2X Master Mix with GC Buffer	68	PNase F, Rapid Antibody Standard	257
One <i>Taq</i> Hot Start Quick-Load 2X Master Mix with Standard Buffer	68	PNase F, Remove-IT	257
One <i>Taq</i> Quick-Load 2X Master Mix with Standard Buffer	68	Polbase	305
PCR Buffers, Activity of Restriction Enzymes in	329-330	Poly(A) Polymerase, <i>E. coli</i>	202
PCR Optimization, Guidelines	322	Poly(U) Polymerase	202
Phusion High-Fidelity DNA Polymerase	65, 67	Polyhistidine Tagged Protein Purification	233-234
Phusion High-Fidelity PCR Kit	67	Polymerases, DNA	
Phusion High-Fidelity PCR Master Mix with GC Buffer	67	Activity of Restriction Enzymes in PCR Buffers	329-330
Phusion High-Fidelity PCR Master Mix with HF Buffer	67	<i>Bst</i> 2.0 DNA Polymerase	79
Phusion Hot Start Flex 2X Master Mix	67	<i>Bst</i> 3.0 DNA Polymerase	79
Phusion Hot Start Flex DNA Polymerase	67	<i>Bst</i> DNA Polymerase (Large Fragment, Full Length)	79
Polymerase Selection Chart	65	<i>Bst</i> 2.0 WarmStart DNA Polymerase	79
PreCR Repair Mix	82	<i>Bsu</i> DNA Polymerase Large Fragment	84
Q5 High-Fidelity 2X Master Mix	66	Deep Vent DNA Polymerases	71
Q5 High-Fidelity DNA Polymerase	65, 66	DNA Polymerase I (<i>E. coli</i>)	83
Q5 High-Fidelity PCR Kit	66	DNA Polymerase I, Large (Klenow) Fragment	83
Q5 Hot Start High-Fidelity 2X Master Mix	66	Epimark Hot Start <i>Taq</i> DNA Polymerase	71
Q5 Hot Start High-Fidelity DNA Polymerase	66	Guidelines for Optimization	322
qPCR	72-73	Hemo KlenTaq	65, 70
RT-qPCR	72-77	Hot Start <i>Taq</i>	65, 69
Selection Chart	65	Klenow Fragment (3'→5' exo-)	83
<i>Taq</i> Master Mixes (2X, 5X, Quick-Load)	69	LongAmp <i>Taq</i> DNA Polymerases	65, 70
<i>Taq</i> DNA Polymerase with Standard <i>Taq</i> Buffer	65, 69	One <i>Taq</i> DNA Polymerases	65, 68
<i>Taq</i> DNA Polymerase with Standard <i>Taq</i> (Mg-free) Buffer	69	phi29 DNA Polymerase	80
<i>Taq</i> DNA Polymerase with ThermoPol Buffer	69	Phusion High-Fidelity DNA Polymerases	65, 67
<i>Taq</i> PCR Kit	69	Q5 High-Fidelity DNA Polymerases	65, 66
Troubleshooting Guide	323	Reaction Buffers	85
Vent DNA Polymerase	71	Selection Chart	65
Vent (exo-) DNA Polymerase	71	<i>Sulfolobus</i> DNA Polymerase IV	82
PCR Fidelity Estimator	305	T4 DNA Polymerase	83
PCR from Blood		T7 DNA Polymerase	84
Hemo KlenTaq	70	<i>Taq</i> DNA Polymerases	65, 69
PCR Selector	305	Terminal Transferase	84
Peptide Display Cloning System, Ph.D.	273	Therminator DNA Polymerase	82
Peptide Ligation	230	Troubleshooting guide	323
Performance Chart for Restriction Enzymes	309-314	Vent DNA Polymerase	71
PFG Ladders	182	Polymerases, RNA	
PfiFI	47	<i>E. coli</i> Poly(A) Polymerase	202
PfiMI	47	Hi-T7 RNA Polymerase (Std & High Conc.)	202
Phage Display Libraries, Ph.D. Kits	272-273	Poly(U) Polymerase	202
Ph.D. Peptide Display Cloning System	273	SP6 RNA Polymerase	202
Ph.D.-7 Phage Display Peptide Library Kit	272	T3 RNA Polymerase	202
Ph.D.-12 Phage Display Peptide Library	272	T7 RNA Polymerase	202
Ph.D.-12 Phage Display Peptide Library Kit	272	T7 RNA Polymerase (Std & High Conc.)	202
Ph.D.-C7C Phage Display Peptide Library Kit	272	Polymerases and Amplification Technologies	62-87
phi29 DNA Polymerase	80	PpuMI	47
phiX174		PreCR Repair Mix	82
HaeIII Digest	183	Primers for cDNA Synthesis	
RF I, RF II DNA	128	Oligo d(T)18 mRNA Primer	208
Virion DNA	128	Oligo d(T)23 VN	208
Phosphatases		Random Primer 6	208
Phosphatase, Antarctic	106	Random Primer 9	208
Phosphatase, Shrimp Alkaline (rSAP)	106	Random Primer Mix	208
Quick CIP	106	Primers, Sequencing	146, 147, 160
Phosphatases, Protein		Properties of Exonucleases and Non-specific Endonucleases	108
Lambda Protein Phosphatase (λ-PPase)	270	Properties of DNA Repair Enzymes & Structure Specific Endos	116
Phosphorylation	104	Proteases	265-269
Phusion High-Fidelity DNA Polymerase (Enzyme, Master Mix, Kit)	67	Protein A Magnetic Beads	233, 235
Phusion Hot Start Flex DNA Polymerase (Enzyme, Master Mix)	67	Protein Data	369
pKLAC2, Map of	358	Protein Deglycosylation Mix II	258
PI-PspI	58	Protein Expression & Purification Technologies	224-235
PI-SceI	58	Protein Expression with T7 Express Strains	351
Plasmid Miniprep, Monarch Products for	140	Protein G Magnetic Beads	233, 235
PleI	47	Protein Kinases	271
PluTI	47	Proteinase K	64, 268
pMAL-c6T Vector	229	Protein Ligation	230
pMAL-c6T, Map of	359	Protein Markers and Ladders	185
pMiniT 2.0, Map of	360	Blue Prestained Protein Std., Broad range (10–200 kDa)	185
PmeI	47	Color Prestained Protein Std., Broad Range (10–250 kDa)	185
PmlI	47	Unstained Protein Std., Broad Range (10–200 kDa)	185

Protein Phosphatases	270
Protein Standards	185
Proteome Analysis	266-267
ProtoScript II First Strand cDNA Synthesis Kit	86, 206, 209
ProtoScript II Reverse Transcriptase	86, 206
ProtoScript First Strand cDNA Synthesis Kit	86, 206, 209
PshAI	47
Psil-v2	48
pSNAP, Vector	298
pSNAP, Vector, Map of	362
pSNAP-tag(T7)-2 Vector	298
PspGI	48
PspOMI	48
PspXI	48
PstI	48
PstI-HF	48
pTWIN1 Vector	230
pTXB1 Vector, Map of	363
pTXB1 Vector	230
pTYB21 Vector	230
pTYB21 Vector, Map of	364
pUC19 DNA	128
pUC19, Map of	365
PURExpress <i>In Vitro</i> Protein Synthesis Kit	228
PURExpress Δ (aa, tRNA) Kit	228
PURExpress Δ RF123 Kit	228
PURExpress Δ Ribosome Kit	228
PURExpress Disulfide Bond Enhancer	228
Purple Gel Loading Dye	59, 180
PvuI	48
PvuI-HF	48
PvuII	49
PvuII-HF	49
Pyrophosphatase, Inorganic (yeast)	107, 203
Pyrophosphatase, Inorganic (<i>E. coli</i>)	107, 203
Pyrophosphatase, NudC	107, 203
Pyrophosphatase, Thermostable Inorganic	107, 203
Pyrophosphohydrolase, RNA 5' (RppH)	214

Q

Quality	2
Q5 High-Fidelity DNA Polymerases (Enzyme, Master Mix, Kit)	64, 65, 66, 190
Q5 Hot Start High-Fidelity DNA Polymerases (Enzyme, Master Mix)	64, 65, 66
Q5 Site Directed Mutagenesis Kit (with or without comp cells)	98, 190
Q5U Hot Start High-Fidelity DNA Polymerase	66
qPCR	
Luna Universal qPCR Master Mix	73
Luna Universal Probe qPCR Master Mix	73
Optimization Tips	324
Troubleshooting Guide	326
Quick Blunting Kit	98
Quick CIP	105, 216
Quick Ligation Kit	99, 100, 101
Quick-Load DNA Ladders	182
Quick-Load Purple DNA Ladders	182
Quick-Load <i>Taq</i> 2X Master Mix	69

R

Random Primers	208
Random Primer Mix	208
Rapid PNGase F Antibody Standard	257
Rapid PNGase F, (non-reducing format)	257
Reaction Buffers (Polymerases)	85
Reaction Buffers (Restriction Enzymes)	58, 309
Read Coverage Calculator	305
Reagents for DNA Library Preparation (DNA & RNA)	144-175
REBASE	305
RecA, RecA ₁	127
RecJ	114

Recombinant Albumin, Molecular Biology Grade	58
Recombinase, Cre	123
Reference Appendix	302-369
Remove-IT PNGase F	257
Repair Proteins	82, 116-122
Restriction Endonucleases	
Activity Chart	309-314
Activity of DNA Modifying Enzymes in CutSmart Buffer	315
Activity of Enzymes at 37°C	315
Activity of Restriction Enzymes in PCR Buffers	329-330
Cleavage of Supercoiled DNA	318
Cleavage Close to the End of DNA Fragments	328
Compatible Cohesive Ends	319-320
Diluent Buffers	59, 309
Double Digestion	307
Extended Digests	328-329
Heat Inactivation	309-314
High Fidelity (HF)	21, 316
Isoschizomers	317
Methylation-Dependent	280-281
Methylation Sensitivity (Dam, Dcm, CpG)	321
Optimizing Reactions	306
Performance Chart	309-314
Reaction Buffers	58
Star Activity	316
Survival in a Reaction	317-318
Thermophiles at 37°C, Activity of	315
Troubleshooting Guide	308
Types	307
Restriction Maps	355-365
Reverse Transcription	
AMV Reverse Transcriptase	86, 206, 207
M-MuLV Reverse Transcriptase	86, 206, 207
One <i>Taq</i> One-Step RT-PCR Kit	68, 81
One <i>Taq</i> RT-PCR Kit	68, 81
ProtoScript First Strand cDNA Synthesis Kit	86, 206, 209
ProtoScript II First Strand cDNA Synthesis Kit	86, 206, 209
ProtoScript II Reverse Transcriptase	86, 206
Template Switching RT Enzyme Mix	86, 206, 208
WarmStart RTx Reverse Transcriptase	86, 206, 207
Ribonuclease Contamination	198
Ribonucleoside Vanadyl Complex	218
Ribonucleotide Solution (Mix, Set)	85, 203
Ribosomal RNA Depletion	155-157
Ribosome, <i>E. coli</i>	228
RNA 5' Pyrophosphohydrolase (RppH)	214
RNA Cap Structure Analogs	204
RNA Cleanup, Monarch Kits for	138, 220
RNA Depletion	155, 157
RNA Library Preparation for NGS	144-175, 218
RNA Ligases	
Reported Activities and Applications	210
RtcB Ligase	210, 211, 213
Selection Chart	211
SplintR Ligase	211, 214
T4 RNA Ligase 1	210, 211, 212
T4 RNA Ligase 2	210, 211, 212
T4 RNA Ligase 2, truncated	210, 211, 212
T4 RNA Ligase 2, truncated K227Q, KQ	210, 211, 213
Thermostable 5' App DNA/RNA Ligase	210, 211, 213
RNA Loading Dye, 2X	221
RNA Markers and Ladders	
dsRNA Ladder	184, 221
Low Range ssRNA Ladder	184, 221
microRNA Marker	184, 221
ssRNA Ladder	184, 221
RNA Polymerases	
<i>E. coli</i> Poly(A) Polymerases	202
<i>E. coli</i> (Core & Holo.)	203
Hi-T7 RNA Polymerase (Std & High Conc.)	202
Poly(U) Polymerase	202
SP6 RNA Polymerase	202
T3 RNA Polymerase	202
T7 RNA Polymerase (Std & High Conc.)	202

RNA Purification, Monarch Kits for	137-138, 220
RNA Library Preparation	144-175
RNA Polymerases	202-203
RNA Purification, Monarch Kits	137-138, 220
RNA Reagents	196-221
RNA Synthesis	
EnGen sgRNA Synthesis Kit	190, 193
HiScribe SP6 RNA Synthesis Kit	199, 201
HiScribe T7 ARCA mRNA Kit (with or without Tailing)	190, 200, 201
HiScribe T7 Quick High Yield RNA Synthesis Kit	190, 199, 201
HiScribe T7 High Yield RNA Synthesis Kit	190, 199, 201
RNase I	215
RNase III, Shortcut	216
RNase B	256
RNase Contamination, Avoiding	198
RNase H	215
RNase HII	215
RNase Inhibitor, Human Placenta	218
RNase Inhibitor, Murine	217
RNase T	112, 217
Ribonucleotides	85, 203
Routine PCR	
LongAmp Taq DNA Polymerase	65, 70
OneTaq DNA Polymerase	65, 68
Phusion DNA Polymerase	65, 67
Q5 High-Fidelity DNA Polymerase	65, 66
Taq DNA Polymerase	65, 69
RppH	211, 214
RsaI	49
RsrII	49
RtcB Ligase	210, 211, 213
RT-LAMP	78
RT-PCR	
OneTaq One-Step RT-PCR Kit	81, 210
OneTaq RT-PCR Kit	81, 210
RT-qPCR	
Luna Cell Ready One-Step RT-qPCR Kit	75, 210
Luna Cell Ready Probe One-Step RT-qPCR Kit	75, 210
Luna Probe One-Step RT-qPCR 4X Mix with UDG	74, 210
Luna Probe One-Step RT-qPCR Kit (No ROX)	74, 210
Luna SARS-CoV-2 RT-qPCR Multiplex Assay Kit	76, 210
Luna Universal One-Step RT-qPCR Kit	74, 210
Luna Universal Probe One-Step RT-qPCR Kit	74, 210
Optimization Tips for Luna RT-qPCR	325
Troubleshooting Guide for Luna RT-qPCR	327

S

7-deaza-dGTP	85
SacI	49
SacI-HF	49
SacII	49
S-adenosylmethionine (SAM)	58
Sall	49
Sall-HF	50
Salt-T4 DNA Ligase	99, 101
Sample Preparation	
Next Gen Sequencing	144-175
Nucleic Acid Purification	132-141
Target Enrichment	163-165
DNA Enrichment	167, 220, 279
SapI	50
SARS-CoV-2 Detection	6, 76, 77
SARS-CoV-2 Rapid Colorimetric LAMP Assay Kit	64, 76
SARS-CoV-2 Positive Control (N gene)	76
Sau3AI	50
Sau96I	50
SbfI	50
SbfI-HF	50
Scal-HF	50
ScrFI	50
Sce PstI	203
SexAI	51

SfiNI	51
SfiI	51
SfiII	51
SfiI	51
SgrAI	51
sgRNA Synthesis	190, 193
Shipping Information	8
Shortcut RNase III	216
Shrimp Alkaline Phosphatase (rSAP)	106, 216
SHuffle Competent <i>E. coli</i> (Express, T7 Express)	232, 240, 241, 248
SHuffle T7 Competent <i>E. coli</i> (Express, lysY)	232, 240, 241, 249
Single-stranded DNA Binding Protein	127
Site-Directed Mutagenesis Kit, Q5	98
SmaI	51
SmlI	51
SnaBI	52
SNAP Analysis, Substrates for	296
SNAP-Biotin	298
SNAP-Capture Magnetic Beads	299
SNAP-Capture Pull-Down Resin	299
SNAP-Cell 430	296
SNAP-Cell 505-star	296
SNAP-Cell 647-SiR	296
SNAP-Cell Block	297
SNAP-Cell Oregon Green	296
SNAP-Cell TMR-Star	296
SNAP-Surface 488	296
SNAP-Surface 549	296
SNAP-Surface 594	296
SNAP-Surface 649	296
SNAP-Surface Alexa Fluor 488	296
SNAP-Surface Alexa Fluor 546	296
SNAP-Surface Alexa Fluor 647	296
SNAP-Surface Block	297
SNAP-tag	
Comparison to GFP	295
Labeling with SNAP-tag/Troubleshooting Guide	354
Overview	294
SNAP-tag Purified Protein	297
Vectors	298
SOC Outgrowth Medium	242
Social Responsibility	2
Sodium Orthovanadate (Vanadate)	270
SP6 RNA Polymerase	202
SpeI	52
SpeI-HF	52
SphI	52
SphI-HF	52
SplintR Ligase	104, 211, 214
SrfI	52
SspI	52
SspI-HF	52
ssRNA Ladder	184, 221
ssRNA Ligase (T4 RNA Ligase 1)	210, 211, 212
SssI (CpG) Methyltransferase	124, 286
Standards, Protein	185
Star Activity	316
Streptavidin Magnetic Beads	219, 233, 235
StuI	53
Styl-HF	53
StyD4I	53
<i>Sulfolobus</i> DNA Polymerase IV (Lesion Bypass Polymerase)	82
Supercoiled DNA, Cleavage of	318
Supercoiled DNA Ladder	184
Survival in a Reaction	317-318
Swal	53
Synthetic Biology	
BioBrick Assembly Kit	97
Gibson Assembly	95
Golden Gate Assembly	96
NEBuilder HiFi DNA Assembly Kit	94
USER Enzyme	97

T

T3 DNA Ligase	102
T3 RNA Polymerase	202
T4 DNA Ligase	99, 100, 101
T4 DNA Polymerase	83
T4 Gene 32 Protein	64, 127
T4 PDG (T4 Endonuclease V)	121
T4 Phage β -glucosyltransferase (T4-BGT)	278
T4 Polynucleotide Kinases	104
T4 RNA Ligase 1 (ssRNA Ligase)	210, 211, 212
T4 RNA Ligase 2 (dsRNA Ligase)	210, 211, 212
T4 RNA Ligase 2, truncated	210, 211, 212
T4 RNA Ligase 2, truncated K227Q and KQ	210, 211, 213
T4 RNA Ligase Selection Chart	211
T7 DNA Ligase	102
T7 DNA Polymerase	84
T7 Endonuclease I	122, 190, 192
T5 Exonuclease	114
T7 Exonuclease	115
T7 Express Competent <i>E. coli</i> , SHuffle	232, 240, 241, 249
T7 Express (High Efficiency), Competent <i>E. coli</i>	232, 240, 241, 247
T7 Express <i>lysY</i> Competent <i>E. coli</i> , SHuffle	232, 240, 241, 249
T7 Express, <i>lysY</i> (High Efficiency), Competent <i>E. coli</i>	232, 240, 241, 247
T7 Express, <i>lysY/h</i> (High Efficiency), Competent <i>E. coli</i>	232, 240, 241, 247
T7 Express Strains, Protein Expression with	351
T7 RNA Polymerase (Std. & High Conc.)	202
<i>TaqI</i> -v2	53
<i>TaqI</i> Methyltransferase	126, 286
<i>Taq</i> 2X Master Mix	69
<i>Taq</i> 2X Master Mix, Quick-Load	69
<i>Taq</i> 5X Master Mix	69
<i>Taq</i> DNA Ligase	103
<i>Taq</i> DNA Polymerase with Standard <i>Taq</i> Buffer	69
<i>Taq</i> DNA Polymerase with Standard <i>Taq</i> (Mg-free) Buffer	69
<i>Taq</i> DNA Polymerase with ThermoPol Buffer	69
<i>Taq</i> PCR Kit	69
<i>Taq</i> Reaction Buffers	85
Target Enrichment	163-165
Technical Support	302
TelN Proteolomerase	123
Template Switching RT Enzyme Mix	86, 206, 208
Terminal Transferase	84
TEV Protease	229, 234, 268
TfiI	53
Therminator DNA Polymerase	82
Thermolabile Exonuclease I	110
Thermolabile USER II Enzyme	97
Thermolabile Proteinase K	64, 268
Thermophiles at 37°C, Activity of	315
Thermophilic Helicase Dependent Amplification	80
Thermophilic Polymerases, Diluent Buffers	85
Thermophilic Polymerases, Reaction Buffers	85
Thermostable 5' App DNA/RNA Ligase	210, 211, 213
Thermostable DNA Ligase (<i>Taq</i>)	103
Thermostable FEN1	112
Thermostable Inorganic Pyrophosphatase	107, 203
Thermostable Ligase Reaction Temperature Calculator	305
Thermostable RNase H	215
Tm Calculator	305
<i>Tma</i> Endonuclease III	119
Topoisomerase	
Cre Recombinase	123
Topoisomerase I	123
Total RNA Miniprep, Monarch Kits for	137
Transformation Efficiencies, Enhancing	350
TriDye DNA Ladders	182
TriDye Ultra Low Range DNA Ladder	182
Troubleshooting Guides	
Cloning	336-338
Genomic DNA Purification	343
Luna qPCR	326
Luna RT-qPCR	327
DNA Cleanup & Plasmid Purification	341
NEBNext DNA Library Prep Kits	352

PCR	323
Protein Expression	351
Restriction Enzymes	308
SNAP-tag Technology	354
Total RNA Extraction & Purification	345
Trypsin-digested BSA MS Standard (CAM-modified)	267
Trypsin-ultra, Mass Spectrometry Grade	266
TseI	53
Tsp45I	53
TspMI	54
TspRI	54
<i>Tte</i> UvrD Helicase	64, 107
Tth111I	54
<i>Tth</i> Argonaute (TiAgo)	128, 191
<i>Tth</i> Endonuclease IV	118
TubeOpener	58

U

UDG	121, 122
Universal miRNA Cloning Linker	221
Unstained Protein Standard, Broad Range (10-200 kDa)	185
Uracil-DNA Glycosylase (UDG)	121
Uracil-DNA Glycosylase (UDG), Afu	121
USER Enzyme	97

V

Vaccinia Capping System	204
(Vanadate) Sodium Orthovanadate	270
Vectors	
Cloning	128
IMPACT	230
<i>K. lactis</i>	231
NEBExpress	229
SNAP/CLIP	298
Vent DNA Polymerase	71
Vent (exo-) DNA Polymerase	71

W

WarmStart Colorimetric LAMP/RT-LAMP 2X Master Mix	64, 78
WarmStart Colorimetric LAMP/RT-LAMP 2X Master Mix (with UDG)	64, 78
WarmStart Fluorescent LAMP/RT-LAMP Kit	64, 78
WarmStart Fluorescent LAMP/RT-LAMP Kit (with UDG)	64, 77
WarmStart Multi-Purpose LAMP/RT-LAMP 2X Master Mix (with UDG)	64, 77
WarmStart RTx Reverse Transcriptase	64, 86, 206
Water, Nuclease-free	58

X

XbaI	54
XcmI	54
XhoI	54
XmaI	54
XmnI	54
XRN-1	216

Y

Yeast Carbon Base Medium Powder	231
Yeast Expression	231
yDcpS	205

Z

Zral	55
------	----

For over 45 years, the New England Biolabs catalog has been a resource for scientists around the world. Each edition contains a collection of mini-reviews that addresses various scientific, environmental and/or humanitarian topics. This year, we are dedicating the Catalog to our Founder, Donald G. Comb, and sharing some of the values that he was passionate about. As part of this effort, NEB will offer support to an organization working to conserve biodiversity throughout the world.

The Marine Conservation Action Fund

The New England Aquarium's Marine Conservation Action Fund (MCAF) builds hope for the future of the ocean by supporting conservation leaders who tackle critical marine research and conservation needs. From funding for urgent, community-based initiatives to mentorship, networking, and professional development, MCAF has supported more than 170 projects in over 60 countries since it was founded in 1999. The program focuses its support on country-national scientists and grassroots leaders in low-and middle-income countries – areas that host much of the world's biodiversity but have fewer financial resources for conservation efforts.

MCAF, which is housed in the Aquarium's Anderson Cabot Center for Ocean Life and is a grantee of the New England Biolabs Foundation, has supported community leaders in bringing about landmark conservation victories through pioneering scientific research, outreach and collaboration with local stakeholders. These landmark successes include full legal protection of giant manta rays in Peru; the launch of the first marine mammal stranding network in Iran; the rescue of hundreds of entangled sea turtles in collaboration with local fishers in Ghana; and, the designation or strengthening of marine protected areas (MPA's) in countries such as Costa Rica, Bangladesh, the Philippines, Indonesia and Panama.

MCAF builds enduring two-way relationships with conservation leaders, connects them to a supportive global community of their peers and helps to ensure their long-term conservation impact. Together we are striving towards our vision of a healthy and vibrant ocean stewarded by local leaders and their communities for generations to come.

To learn more, please visit: www.neaq.org/mcaf



Planeta Océano



Nature Conservation Foundation



Misión Tiburón



Planeta Océano



Bird's Head Leatherback Program

Find out what's trending in science



Hear from scientists as they discuss a wide range of molecular biology topics and techniques. Examples of past webinar topics include updates in SARS-CoV-2 detection, high molecular weight DNA extraction, rRNA depletion and cell-free protein expression.

www.neb.com/NEBTVwebinars



Our podcast series features conversations with scientific colleagues on science, careers and backstories that impacted their lives. Be sure to check out our COVID-19 Researcher Spotlight Series, featuring scientists from around the world discussing publications, protocols and detection methods related to COVID-19.

www.neb.com/nebpodcast



Our science blog is designed to share inspirational stories about trends in the life sciences, lab tips to help you save time and life lessons to reflect on. Learn about topics such as large scale mRNA synthesis, loop-mediated isothermal amplification, bacterial art and lab safety during COVID-19.

www.neb.com/NEBinspired



Our video series explores cutting edge topics and trends in molecular biology. Find interviews with NEB scientists, highlights from recent publications, and overviews of new products and online tools for your research. Recent episodes address COVID-19 related topics such as vaccine development and methods for SARS-CoV-2 detection.

www.neb.com/NEBTV

Environmental Philosophy & the NEB Catalog

The NEB Catalog & Technical Reference is printed with sustainability in mind.

- The cover is printed on Forest Stewardship Council® (FSC®) Certified paper
- The interior pages are printed on Forest Stewardship Council® (FSC®) Certified, 30% recycled paper
- The printing facility is FSC® Certified
- The printing facility uses environmentally-friendly inks and coatings

Once again, NEB has neutralized the carbon emissions associated with producing and delivering this catalog by working with Carbonfund.org, a world-leading provider of carbon reduction solutions. This means that the negative environmental impact of the unavoidable CO₂ generated by printing and distributing this catalog has been mitigated through verified carbon offset projects. Offsets will be used to support the Blandin Improved Forestry Management Project in Minnesota, which has undertaken an approach that combines reduced impact logging with the establishment of conservation easements. In 2015, the NEB catalog was the first carbon neutral-certified catalog produced in the United States. More information can be found at Carbonfund.org.

Because the catalog pages, excluding the cover, use 30% post-consumer fiber instead of virgin fiber, the following savings were incurred:


- 1,270 trees saved
- 5,000 lbs of solid waste not generated
- 280,000 gallons of water conserved
- 720,000 pounds of greenhouse gas emissions avoided
- 570,000,000 British Thermal Units (BTUs) of energy not consumed

Note: Statistics and calculations are based on the Environmental Paper Network Paper Calculator v3.2. For more information, visit www.papercalculator.org.

FOREST STEWARDSHIP COUNCIL® is a registered trademark of The Forest Stewardship Council.



NEB's Global Distribution Network



8 SUBSIDIARY
OFFICES

>60 DISTRIBUTION
PARTNERS*

To find your local subsidiary or distributor, visit www.neb.com/international-support

USA

New England Biolabs, Inc.
240 County Road
Ipswich, MA 01938-2723
Telephone (978) 927-5054
Toll Free (USA Orders) 1-800-632-5227
Toll Free (USA Tech) 1-800-632-7799
Fax (978) 921-1350
email: info@neb.com
www.neb.com

AUSTRALIA/NEW ZEALAND

New England Biolabs (Australia) PTY
Telephone: 1800 934 218
info.au@neb.com
NEB New Zealand Ltd
Telephone: 0800 437 209
info.nz@neb.com

CANADA

New England Biolabs, Ltd.
Telephone (905) 665-4632
Toll Free 1-800-387-1095
Fax (905) 665-4630
Fax Toll Free 1-800-563-3789
email: info.ca@neb.com
www.neb.ca

CHINA

New England Biolabs (Beijing), Ltd.
Telephone 010-82378265/82378266
Free Call 400-811-2220
Free Call 400-690-3366
Fax 010-82378262
email: info@neb-china.com
www.neb-china.com

FRANCE

New England Biolabs France
Telephone +33 (0) 1 69 90 87 15
Free Call 0800 100 632
Free Fax 0800 100 610
email: info.fr@neb.com
www.neb-online.fr

GERMANY/AUSTRIA

New England Biolabs GmbH
Telephone +49 (0) 69/305 23140
Free Call 0800/246 5227 (Germany)
Free Call 00800/246 52277 (Austria)
Fax +49 (0) 69/305 23149
Free Fax 0800/246 5229 (Germany)
email: info.de@neb.com
www.neb-online.de

JAPAN

New England Biolabs Japan, Inc.
Telephone +81 (0) 3 5669 6191
Fax +81 (0) 3 5669 6192
email: info.jp@neb.com
www.nebj.jp

SINGAPORE

New England Biolabs Pte Ltd.
Telephone +65 638 59623
Fax +65 638 59617
email: sales.sg@neb.com
www.neb.sg

UNITED KINGDOM

New England Biolabs (UK), Ltd.
Telephone +44 (0) 1462 420616
Free Call 0800 318486
Fax +44 (0) 1462 421057
Free Fax 0800 435682
email: info.uk@neb.com
www.neb.uk.com

**NEB also works with additional agents and customers globally. For more information contact info@neb.com.*



LOOKING TO PLACE AN ORDER?

Telephone 1-800-632-5227
9:00 AM to 8:00 PM EST
fax: 1-800-632-7440
email: orders@neb.com
online: www.neb.com



NEED TECHNICAL SUPPORT?

Telephone 1-800-632-7799
email: info@neb.com
online: www.neb.com/techsupport

WANT TO EXPERIENCE THE NEB CATALOG IN AUGMENTED REALITY?



DOWNLOAD

the NEB Augmented Reality (AR) app
for iPhone® or iPad® at the Apple® App
Store or for Android™ on Google Play™.



EXPLORE

the NEB Catalog & Technical
Reference and keep an eye out
for the augmented reality icon.



DISCOVER

videos, tutorials and immersive
experiences by scanning the
page using the AR app.

iPhone®, iPad® and Apple® are registered
trademarks of Apple Computer, Inc.
Android™ and Google Play™ are trademarks
of Google, Inc.

Swimming sea turtle and sandy seabed.
Underwater video from scuba diving.
Credit: blueseac/Pond5, Adobe Stock

