




Now includes EnGen® Spy Cas9 HF1  
& EnGen SpRY Cas9

# Genome Editing

TOOLS FOR CRISPR/CAS APPLICATIONS




# Genome Editing: Tools for CRISPR/Cas Applications

Genome editing is enabled by the development of tools to make precise, targeted changes to the genome of living cells. Recent approaches to targeted genome modification – zinc-finger nucleases (ZFNs) and transcription-activator like effector nucleases (TALENs) – enable researchers to generate mutations by introducing double-stranded breaks to activate repair pathways. These approaches are costly and time consuming to engineer, limiting their widespread use, particularly for large scale, high-throughput studies. Recently, methods based on a bacterial CRISPR-associated protein-9 nuclease (Cas9) from *Streptococcus pyogenes* have generated considerable excitement.

CRISPR (Clustered Regularly Interspaced Short Palindromic Repeats) and CRISPR-associated (Cas) genes are essential for adaptive immunity in select bacteria and archaea, enabling the organisms to respond to and eliminate invading genetic material.

## CRISPR/Cas *in vivo*: Bacterial Adaptive Immunity



*In the acquisition phase, foreign DNA is incorporated into the bacterial genome at the CRISPR locus. The CRISPR locus is then transcribed and processed into crRNA during crRNA biogenesis. During interference, Cas endonuclease complexed with crRNA cleaves foreign DNA containing a crRNA complementary sequence adjacent to the PAM sequence. (Figure not drawn to scale.)*

## CRISPR/Cas Genome Editing

The simplicity of the CRISPR nuclease system (nuclease and guide RNA), makes this system attractive for laboratory use. Breaks activate repair through error prone Non-Homologous End Joining (NHEJ) or Homology Directed Repair (HDR). In the presence of a donor template with homology to the targeted locus, the HDR pathway may operate, allowing for precise mutations to be made. In the absence of a template, NHEJ is activated, resulting in insertions and/or deletions (indels), which disrupt the target locus (1,2).

### TOOLS & RESOURCES


Visit [www.neb.com/GenomeEditing](http://www.neb.com/GenomeEditing) to find our up-to-date listing of products and protocols to support this application.

# Direct Introduction of Cas RNP Complexes

The highest efficiency strategy for genome engineering with CRISPR/Cas is direct introduction of Cas9/guide RNA complexes (3–8) or Cas12a (Cpf1)/guide RNA complexes. This method further simplifies CRISPR/Cas workflows and has been reported to increase mutagenic activity (3–5) and reduce off-target editing events (3,4).

NEB® provides purified Cas9 Nuclease, *S. aureus*, *S. equinus* and *S. pyogenes* variants, and Cas12a nuclease (Cpf1), *Lachnospiraceae* bacterium ND2006 with nuclear localization signals as standalone enzymes to support direct introduction of Cas RNP complexes.

## Increased genome editing efficiency using Cas9 RNP delivery



*Cas9* and *sgRNA* targeting a human gene were delivered to HEK293 cells by transfection. Transfected plasmid DNA contained expression cassettes for 2x NLS (N- and C-terminal) *Cas9* and *sgRNA*. Plasmid DNA was delivered using TransIT-X2 (Mirus). Transfected mRNA was modified with pseudouridine and 5-methylcytosine and encoded 2x NLS (N- and C-terminal) *Cas9*. *sgRNA* was co-transfected with the mRNA using TransIT-mRNA. *Cas9* RNPs were delivered in reverse transfections using Lipofectamine RNAiMAX (Life Technologies) using 10 nanomolar final concentration of ribonucleoprotein (RNP). *Cas9* has no NLS in the protein sequence. EnGen *Cas9* has N- and C-terminal NLSs. The efficiency of editing was determined using T7 Endonuclease I assay and is expressed as % modification.

## Ordering Information

PRODUCT	NEB #
EnGen Spy Cas9 HF1	<a href="#">M0667T/M</a>
EnGen SpRY Cas9	<a href="#">M0669TM</a>
EnGen Spy Cas9, NLS	<a href="#">M0646T/M</a>
Cas9 Nuclease, <i>S. pyogenes</i>	<a href="#">M0386S/T/M</a>
EnGen Spy Cas9 Nickase	<a href="#">M0650S/T</a>
EnGen Spy dCas9 (SNAP-tag®)	<a href="#">M0652S/T</a>
EnGen Sau Cas9	<a href="#">M0654S/T</a>
EnGen Seq1 Cas9	<a href="#">M0668T</a>
EnGen Lba Cas12a (Cpf1)	<a href="#">M0653S/T</a>

# NEB's CRISPR Nuclease toolbox (selected products)

**EnGen Spy Cas9 HF1, *Streptococcus pyogenes***  
High-fidelity, quadruple substitution variant of EnGen Spy Cas9 Nuclease from *Streptococcus pyogenes* with reduced non-specific DNA cleavage

**EnGen Lba Cas12a, *Lachnospiraceae* bacterium ND2006**  
AT-rich PAM, extended temperature range

**EnGen Seq1 Cas9, *Streptococcus equinus***  
5'- NAGA -3' PAM sequence allows targeting of additional genomic regions.

**EnGen Spy Cas9 Nickase**  
*In vitro* nicking of dsDNA. Increased genome editing specificity *in vivo*, requiring adjacent targets.


# Reducing Off-Target DNA Cleavage

EnGen Spy Cas9 HF1 is a high-fidelity, quadruple substitution (N497A/R661A/Q695A/Q926A) variant of EnGen Spy Cas9 NLS from *Streptococcus pyogenes* with reduced non-specific DNA cleavage. Spy Cas9 is an RNA-guided endonuclease that catalyzes site-specific cleavage of double-stranded DNA. The single guide RNA (sgRNA) targets Cas9 to the region immediately upstream of a 5'-NGG-3' protospacer adjacent motif (PAM) producing a double-stranded break 3 bases upstream of the PAM (9). EnGen Spy Cas9 HF1 contains Simian virus 40 (SV40) T antigen nuclear localization sequence (NLS) on the N- and C-termini of the protein.

## Ordering Information

PRODUCT	NEB #
EnGen Spy Cas9 HF1	<a href="#">M0667M/T</a>
EnGen sgRNA Synthesis Kit	<a href="#">E3322V/S</a>

## EnGen Spy Cas9 HF1 demonstrates increased sensitivity to mismatches between guide RNA and DNA targets *in vitro*



Comparison of the tolerance of mismatches between the guide RNA sequence and target DNA sequence of EnGen Spy Cas9 NLS, EnGen Spy Cas9 HF1, and other commercially available high fidelity Cas9 variants. One of several guide RNAs encoding a single, double, or triple mismatch with a fluorescently labeled dsDNA substrate were allowed to form a ribonucleoprotein (RNP) complex with each of five Cas9 variants. A fully matched guide RNA was included as a control. The RNPs were incubated with the substrate at a 2:1 ratio at 37°C for 5 minutes. The percent substrate cleavage for each RNP complex was measured by capillary electrophoresis. Results were graphed as a heat map with white representing no cleavage and increasing intensity of blue indicating increasing percent cleavage. The guide RNA sequence is indicated in each row, with mismatches denoted in green. The DNA protospacer sequence is 5' - AGAAGTGGCAGAGGAGGAGG - 3' and the protospacer adjacent motif (PAM) is 5' - TGG - 3'. EnGen Spy Cas9 HF1 demonstrates increased sensitivity to mismatches by showing the greatest ratio of on-target cleavage to average cleavage of off-targets.


# Rapid Generation of sgRNA for Spy Cas9

The EnGen sgRNA Synthesis Kit simplifies the generation of microgram quantities of custom sgRNAs in an hour or less by combining template synthesis and transcription. The single-tube reaction is easy to set up and requires a single ~55 nt ssDNA target-specific oligonucleotide, which is combined with the Reaction Mix and Enzyme Mix included in the kit. sgRNAs are suitable for use in downstream applications, including CRISPR/Cas9-based genome editing and *in vitro* DNA cleavage. This single-reaction format offers ease-of-use and eliminates separate DNA amplification and template clean up steps. This kit is compatible with EnGen Spy Cas9 NLS, Cas 9 Nuclease, *S. pyogenes*, EnGen Spy Cas9 Nickase, EnGen Spy dCas9 (SNAP-tag) and EnGen Spy Cas9 HF1.

 **Need help configuring target-specific DNA oligos?**

Try our **EnGen sgRNA Template Oligo Designer** (accessible through NEBioCalculator<sup>®</sup> at [NEBiocalculator.neb.com](http://NEBiocalculator.neb.com))

## EnGen sgRNA Synthesis Kit overview



“ This kit is really easy to use and will save us plenty of time in making sgRNAs! Thanks for the streamlined method! ”

– Postdoctoral Researcher,  
Harvard University

**A.** The target-specific oligo contains the T7 promoter sequence, ~20 nucleotides of target-specific sequence and a 14-nucleotide overlap sequence complementary to the *S. pyogenes* Cas9 Scaffold Oligo, supplied in the reaction mix. Target-specific oligos are mixed with the EnGen 2X sgRNA Reaction Mix and the EnGen sgRNA Enzyme Mix at room temperature.

**B.** At 37°C the two oligos anneal at the 14-nucleotide overlap region of complementarity.

**C.** The DNA polymerase contained in the EnGen sgRNA Enzyme Mix extends both oligos from their 3' ends, creating a dsDNA template.

**D.** The RNA polymerase contained in the EnGen sgRNA Enzyme Mix recognizes the dsDNA of the T7 promoter and initiates transcription, resulting in a target-specific sgRNA.



All steps occur in a single reaction during a 30-minute incubation at 37°C.

## Evaluating Targeting Efficiency with Mutation Detection Kits and Assays

A widely used method to identify mutations is the T7 Endonuclease I mutation detection assay (10,11). This assay detects heteroduplex DNA that results from the annealing of a DNA strand, including desired mutations, with a wild-type DNA strand. Multiple options are available for performing robust detection of genome editing events including the T7 Endonuclease-based EnGen Mutation Detection Kit and Authenticase, which provides improved detection of single- and two-basepair mutations.

### Authenticase Mismatch Detection Assay

**MISMATCH DETECTION ASSAY**  
to estimate genome editing efficiency



Authenticase can replace T7 Endonuclease I in the mismatch detection assay used to assess the efficiency of genome editing (S1 is the starting material. P1 and P2 are products of Authenticase digestion).

Try Monarch RNA Cleanup Kits for purification of sgRNA after synthesis.

Learn more at

[www.neb.com/MonarchRNACleanup](http://www.neb.com/MonarchRNACleanup).

### Ordering Information

PRODUCT	NEB #
Authenticase <sup>®</sup>	<a href="#">M0689S/L</a>
EnGen Mutation Detection Kit	<a href="#">E3321S</a>
T7 Endonuclease I	<a href="#">M0302S/L</a>
Q5 <sup>®</sup> Hot Start High-Fidelity 2X Master Mix	<a href="#">M0494S/L/X</a>

Need to determine targeting efficiencies over 50%?

Visit [www.neb.com/Cas9locusmod](http://www.neb.com/Cas9locusmod) to find out how.




Learn more about genome editing vocabulary through the NIST Genome Editing Lexicon  
[www.nist.gov/programs-projects/nist-genome-editing-lexicon](http://www.nist.gov/programs-projects/nist-genome-editing-lexicon)

# EnGen SpRY Cas9 Nuclease

EnGen SpRY Cas9 from *Streptococcus pyogenes* is an engineered RNA-guided DNA endonuclease that catalyzes site-specific cleavage of double-stranded DNA (dsDNA). Targeting requires an ~100 nucleotide single guide RNA (sgRNA) with complementarity to the 20 nucleotide region immediately upstream of a protospacer adjacent motif (PAM) on the dsDNA substrate. Unlike the canonical 5'-NGG-3' PAM of wild-type Spy Cas9, SpRY Cas9 is essentially PAMless *in vitro*, requiring a PAM of 5'-NNN-3' (1,2). DNA cleavage by EnGen SpRY Cas9 produces a double-stranded break occurring 3 nucleotides upstream of the PAM. EnGen SpRY Cas9 contains Simian virus 40 (SV40) T antigen nuclear localization sequence (NLS) on the C-terminus of the protein.

## EnGen SpRY Cas9 has no sequence constraints for dsDNA targeting *in vitro*



EnGen SpRY Cas9 is a variant of Cas9 nuclease from *S. pyogenes* with several point mutations within the PAM interacting domain (1). Unlike wildtype Cas9, EnGen SpRY Cas9 is not constrained by the presence of an NGG PAM and can produce a double stranded break three nucleotides upstream of any trinucleotide sequence in *in vitro* applications.

## BENEFITS

- Eliminate sequence constraints for dsDNA targeting with non-specific PAM (5'-NNN-3' PAM)
- Digest large plasmids in cloning workflows successfully
- Use in conjunction with the EnGen sgRNA Synthesis Kit, *S. pyogenes* (NEB #E3322), EnGen Mutation Detection Kit (NEB #E3321) and NEBuilder® HiFi DNA Assembly Master Mix (NEB #E2621)

## Ordering Information


### PRODUCT

EnGen SpRY Cas9

### NEB #

M0669T/M

## Streamline large construct cloning workflows with EnGen SpRY Cas9 and NEBuilder HiFi DNA Assembly



1 µg of pXba (22,563 bp) was linearized downstream of the start codon of the pVII ORF with 50 nM EnGen SpRY Cas9 and 50 nM sgRNA in 1X NEBuffer™ r3.1 for 1 hour at 37°C. The reactions were either spin column purified or left unpurified before proceeding to DNA assembly. The NEBuilder HiFi DNA Assembly kit was used to insert an oligonucleotide encoding a nuclear localization signal (NLS) tag to pVII according to the recommended protocol. The number of colonies grown after transformation with DNA assembly reactions and no insert controls were measured to qualitatively assess how many transformants arise from undigested plasmid following EnGen SpRY Cas9 digest. Though not a requirement, purification of linearized pXba prior to DNA assembly reduced the percentage of background colonies.

## References:

1. Walton, R.T., et al. (2020) *Science*. 368(6488):290–6.
2. Christie, K.A., et al. (2023) *Nat. Biotechnol.* 41(3):409–16.

# Featured NEB Products Supporting CRISPR Workflows

CRISPR NUCLEASES									
NEB #	PRODUCT NAME	VIAL SIZE	CONC.	VOLUME PER VIAL	REACTION BUFFER	PAM SEQUENCE	NLS	TAGS	MUTATIONS (IF APPLICABLE)
<a href="#">M0386S</a>	Cas9 Nuclease, <i>S. pyogenes</i>	90 pmol	1 μM	0.090 ml	Buffer r3.1	5'- NGG -3'	No NLS	N-terminal 6X His tag	N/A (wild-type)
<a href="#">M0386M</a>	Cas9 Nuclease, <i>S. pyogenes</i>	2500 pmol	20 μM	0.125 ml					
<a href="#">M0386T</a>	Cas9 Nuclease, <i>S. pyogenes</i>	500 pmol	20 μM	0.025 ml					
<a href="#">M0646M</a>	EnGen® Spy Cas9 NLS	500 pmol	20 μM	0.125 ml	Buffer r3.1	5'- NGG -3'	SV40 NLS Sequence on N & C termini	N-terminal 6X His tag	N/A (wild-type)
<a href="#">M0646T</a>	EnGen Spy Cas9 NLS	2500 pmol	20 μM	0.025 ml					
<a href="#">M0650S</a>	EnGen Spy Cas9 Nickase	90 pmol	1 μM	0.090 ml	Buffer r3.1	5'- NGG -3'	SV40 NLS Sequence on N & C termini	N-terminal 6X His tag	D10A in the RuvC nuclease domain
<a href="#">M0650T</a>	EnGen Spy Cas9 Nickase	500 pmol	20 μM	0.025 ml					
<a href="#">M0652S</a>	EnGen Spy dCas9 (SNAP-tag®)	90 pmol	1 μM	0.090 ml	Buffer r3.1	5'- NGG -3'	SV40 NLS Sequence on N & C termini	"N-terminal 6X His tag N-terminal SNAP-tag"	D10A in the RuvC and H840A in the HNH nuclease domains
<a href="#">M0652T</a>	EnGen Spy dCas9 (SNAP-tag)	500 pmol	20 μM	0.025 ml					
<a href="#">M0653S</a>	EnGen Lba Cas12a (Cpf1)	70 pmol	1 μM	0.070 ml	Buffer r2.1	5'- TTTV -3'	SV40 NLS Sequence on N & C termini	N-terminal 6X His tag	N/A (wild-type)
<a href="#">M0653T</a>	EnGen Lba Cas12a (Cpf1)	2000 pmol	100 μM	0.020 ml					
<a href="#">M0654T</a>	EnGen Sau Cas9	500 pmol	20 μM	0.025 ml	Buffer r3.1	5'- NNGRRT -3'	SV40 NLS Sequence on N & C termini	C-terminal 6X His tag	N/A (wild-type)
<a href="#">M0667M</a>	EnGen Spy Cas9 HF1	500 pmol	20 μM	0.125 ml	Buffer r3.1	5'- NGG -3'	SV40 NLS Sequence on N & C termini	N-terminal 6X His tag	N497A/R661A/Q695A/Q926A
<a href="#">M0667T</a>	EnGen Spy Cas9 HF1	2500 pmol	20 μM	0.025 ml					
<a href="#">M0668T</a>	EnGen Seq1 Cas9	500 pmol	20 μM	0.025 ml	Buffer r3.1	5'- NAGA -3'	SV40 NLS Sequence on N & C termini	C-terminal 6X His tag	N/A (wild-type)
<a href="#">M0669M</a>	EnGen SpRY Cas9	500 pmol	20 μM	0.125 ml	Buffer r3.1	5'- NNN -3'	SV40 NLS Sequence on N & C termini	C-terminal 6X His tag	A61R, L1111R, D1135L, S1136W, G1218K, E1219Q, N1317R, A1322R, R1333P, R1335Q, T1337R
<a href="#">M0669T</a>	EnGen SpRY Cas9	2500 pmol	20 μM	0.025 ml					

## ADDITIONAL PRODUCTS SUPPORTING CRISPR WORKFLOWS

PRODUCT NAME	CRISPR/CAS9 APPLICATION	NEB #	SIZE
Q5 Site-directed Mutagenesis Kit (with or without competent cells)	Insertion of target sequence into the Cas9-sgRNA construct and modification of HDR templates	<a href="#">E0554S/E0552S</a>	10 rxns
Q5 High-fidelity DNA Polymerases	High-fidelity construct generation for use with CRISPR workflows	Multiple*	Multiple*
NEBuilder HiFi DNA Assembly Master Mix	Single-tube, isothermal generation of the Cas9-sgRNA construct and HDR templates	<a href="#">E2621S/L/X</a>	10/50/250 rxns
NEBuilder HiFi DNA Assembly Cloning Kit	Single-tube, isothermal generation of the Cas9-sgRNA construct and HDR templates	<a href="#">E5520S</a>	10 rxns
HiScribe T7 ARCA mRNA Kit (with or without tailing)	Generation of Cas9 mRNA with ARCA cap	<a href="#">E2060S/E2065S</a>	20 rxns
HiScribe T7 High Yield RNA Synthesis Kit	Generation of sgRNA and Cas9 mRNA	<a href="#">E2040S</a>	50 rxns
HiScribe T7 Quick High Yield RNA Synthesis Kit	Generation of sgRNA and Cas9 mRNA	<a href="#">E2050S</a>	50 rxns
HiScribe T7 mRNA Kit with CleanCap® Reagent AG	Generation of Cas9 mRNA with CleanCap Reagent AG	<a href="#">E2080S</a>	20 rxns
T7 Endonuclease I	Determination of the targeting efficiency of genome editing protocols	<a href="#">M0302S/L</a>	250/1,250 units
<b>NEW</b> Authenticase	Determination of the targeting efficiency of genome editing protocols	<a href="#">M0689S/L</a>	250/1,250 units
Monarch RNA Cleanup Kit	Cleanup of sgRNA and Cas9 mRNA	<a href="#">T2040S/L</a>	10/100 preps

\* Visit [Q5PCR.com](https://www.neb.com) for ordering information.

## USA

New England Biolabs, Inc.  
Telephone (978) 927-5054  
Toll Free (USA Orders) 1-800-632-5227  
Toll Free (USA Tech) 1-800-632-7799  
Fax (978) 921-1350  
info@neb.com  
www.neb.com

## Australia & New Zealand

New England Biolabs (Australia) PTY  
Telephone: +61 (1800) 934218  
info.au@neb.com

## Canada

New England Biolabs, Ltd.  
Toll Free: 1-800-387-1095  
info.ca@neb.com

## China

New England Biolabs (Beijing), Ltd.  
Telephone: 010-82378265/82378266  
info@neb-china.com

## France

New England Biolabs France SAS  
Telephone: 0800 100 632  
info.fr@neb.com

## Germany & Austria

New England Biolabs GmbH  
Free Call: 0800/246 5227 (Germany)  
Free Call: 00800/246 52277 (Austria)  
info.de@neb.com

## Japan

New England Biolabs Japan, Inc.  
Telephone: +81 (0)3 4545 1422  
cs.jp@neb.com

## Singapore

New England Biolabs, PTE. Ltd.  
Telephone: +65 638 59623  
sales.sg@neb.com

## United Kingdom

New England Biolabs (U.K.) Limited  
Call Free: 0800 318486  
customersupport.uk@neb.com

[www.neb.com](http://www.neb.com)



Products and content are covered by one or more patents, trademarks and/or copyrights owned or controlled by New England Biolabs, Inc. (NEB). The use of trademark symbols does not necessarily indicate that the name is trademarked in the country where it is being read; it indicates where the content was originally developed. See [www.neb.com/trademarks](http://www.neb.com/trademarks). The use of these products may require you to obtain additional third-party intellectual property rights for certain applications. For more information, please email [busdev@neb.com](mailto:busdev@neb.com).

Your purchase, acceptance, and/or payment of and for NEB's products is pursuant to NEB's Terms of Sale at [www.neb.com/support/terms-of-sale](http://www.neb.com/support/terms-of-sale). NEB does not agree to and is not bound by any other terms or conditions, unless those terms and conditions have been expressly agreed to in writing by a duly authorized officer of NEB.

CleanCap® and products incorporating it are sold under license from TriLink Biotechnologies, LLC and may be used for research use only, not for diagnostics, therapeutic procedures or for use in humans. For additional information, please see the products Limited Use License at [www.trilinkbiotech.com/cleancap-research-license](http://www.trilinkbiotech.com/cleancap-research-license).

© Copyright 2024, New England Biolabs, Inc.; all rights reserved

Gene\_Editing – Version 6.0 – 12/24



For help with configuring target-specific DNA oligos, try our **EnGen sgRNA Template Oligo Designer** (accessible through **NEBioCalculator®** at [NEBiocalculator.neb.com](http://NEBiocalculator.neb.com))



For help with designing primers for DNA assembly, try **NEBuilder® DNA Assembly Tool** ([NEBuilder.neb.com](http://NEBuilder.neb.com))



Did you know that many of these products can be purchased in large volumes and custom formats? Learn more at [www.neb.com/customizedsolutions](http://www.neb.com/customizedsolutions)