

# mRNA Fragmentation Protocol (E6100)

## Overview

Starting Material: Purified mRNA (50–250 ng)

## Protocol

### 1. mRNA Fragmentation Protocol

1. Mix the following components in a sterile PCR tube:

Purified mRNA: 1–18  $\mu$ l

10X RNA Fragmentation Reaction Buffer: 2  $\mu$ l

Nuclease-Free Water: variable

-----  
Total volume: 20  $\mu$ l

2. Incubate in a preheated thermal cycler for 5 minutes at 94°C. This is the optimal condition for eukaryotic mRNA ([see Figure 1 on the product page](#)). Other types of mRNA may require optimizing incubation time to obtain desired fragment size distribution.
3. Transfer tube to ice.
4. Add 2  $\mu$ l 10X RNA Fragmentation Stop Solution.

### 2. Clean Up Fragmented RNA Using RNeasy MinElute Spin Columns

1. Add 78  $\mu$ l of the Nuclease-Free Water to the 22  $\mu$ l fragmented RNA from step 4. Purify sample using RNeasy MinElute Cleanup Kit (Qiagen #74204) following manufacture instructions. Elute in 15.5  $\mu$ l Nuclease-Free Water. The recovered volume should be ~14.5  $\mu$ l.

**Note: column purification removes short RNA Fragments and enriches the sample for RNA fragments longer than 200 nucleotides.**

### 3. Alternatively, Clean Up Fragmented RNA Using Ethanol Precipitation

1. Mix the following components in a sterile 1.5 ml microcentrifuge tube:

Fragmented RNA from Step 4: 22  $\mu$ l

3M Sodium Acetate, pH 5.2: 2  $\mu$ l

Linear Acrylamide, 10 mg/ml: 1–2  $\mu$ l

100% Ethanol: 60  $\mu$ l

-----  
Total volume: 85–86  $\mu$ l

2. Incubate at -80°C for 30 minutes.
3. Centrifuge at 14,000 rpm for 25 minutes at 4°C in a microcentrifuge.
4. Carefully remove ethanol.
5. Wash pellet with 300  $\mu$ l of 70% ethanol.
6. Centrifuge and carefully remove 70% ethanol.
7. Air dry pellet for up to 10 minutes at room temperature to remove residual ethanol.
8. Resuspend in 14.5  $\mu$ l Nuclease-Free Water.