

# RNA Circularization using T4 RNA Ligase 1 (ssRNA Ligase) (NEB #M0204)

## Materials Required but not Supplied

### T4 RNA Ligase 1 (ssRNA Ligase)

- Nuclease-free Water (NEB #B1500)
- RNase Inhibitor, Murine (NEB #M0314) (optional)

## Overview

This is a general protocol that serves as a starting point to circularize single-stranded RNA using T4 RNA Ligase 1 (ssRNA Ligase). Increasing the concentration, length of incubation, and adding PEG can increase the number of intramolecular ligations. Blocking one end of the molecule with a dideoxy terminator will prevent the molecule from forming a circle.

## Protocol:

1. Set up the following reaction on ice:

COMPONENTS	20 $\mu$ l REACTION
10X T4 RNA Ligase Reaction Buffer	2 $\mu$ l (1X)
ssRNA with 5'P and 3'OH ends	X $\mu$ l (200 ng - 1 $\mu$ g)
10,000 units/ml T4 RNA Ligase 1	1 $\mu$ l (10 units)
50% PEG 8000	4 $\mu$ l (10%)
10 mM ATP*	0.1 $\mu$ l (50 $\mu$ M)*
<i>Optional</i> - RNase Inhibitor, Murine	0.5 $\mu$ l (20 units)
Nuclease-free Water	to 20 $\mu$ l

\*Note: A working solution of 1mM ATP can be prepared so that 1  $\mu$ l ATP [1mM] can be pipetted into a 20  $\mu$ l reaction volume.

2. Incubate at 25°C for 1-2 hours. For longer oligos, overnight incubation at 16°C may improve yield.
3. Heat inactivate at 65°C for 15 minutes.

## Linked Resources:

- [RNA Ligation Overview](#)

- RNA Cloning Overview