

Immunoprecipitation using Protein A/G Magnetic Beads

Overview

Protocols.io also provides an [interactive version of this protocol](#) where you can discover and share optimizations with the research community.

Introduction

Use 25 μ l of [Protein A](#) or [Protein G Magnetic Beads](#) per 200 μ l of crude cell lysate containing 200-500 μ g of total protein in a standard immunoprecipitation protocol. It is important to increase the volume of beads proportionately for larger cell lysate volumes.

Protocol

1. Cell Lysis

- (1) Rinse a 60 mm culture dish of confluent cells with PBS.
- (2) Lyse the cells with 0.5 ml cold Immunoprecipitation Buffer (150 mM NaCl, 10 mM Tris-HCl (pH 7.4), 1 mM EDTA, 1 mM EGTA (pH 8.0), 0.2 mM sodium ortho-vanadate, 0.2 mM PMSF, 1% Triton X-100, 0.5% NP-40).
- (3) Maintain constant agitation for 30 minutes at 4°C.
- (4) Scrape the cells from the dish. Sonicate on ice for 5 seconds; repeat 4 times. Centrifuge for 5 minutes at 4°C. The supernatant is the crude cell lysate. Assay for total protein then adjust concentration to approximately 1 mg/ml with Immunoprecipitation Buffer.

2. Immunoprecipitation

- (1) This step pre-clears crude cell extract of proteins which can bind non-specifically to the beads. In a 1.5 ml microcentrifuge tube, add 25 μ l Protein A/G Magnetic Beads to 200 μ l of crude cell extract. Gently vortex and incubate at 4°C for 1 hour. Apply magnetic field for 30 seconds to pull beads to the side of the tube. Pipette supernatant to a clean 1.5 ml microcentrifuge tube. Discard beads.
- (2) Add 1-5 μ g of antibody to crude cell lysate vortex and incubate at 4°C for 1 hour. (If monoclonal antibodies are used, add 5 μ g rabbit anti-mouse IgG antibody. Vortex and incubate an additional 30 minutes at 4°C).*
*Alternatively, Protein G Magnetic Beads ([NEB #S1430S](#)) can be used for immunoprecipitations with monoclonal antibodies.
- (3) Add 25 μ l of Protein A/G Magnetic Beads suspension. Gently vortex and incubate with agitation for 1 hour at 4°C.
- (4) Apply magnetic field to pull beads to the side of the tube. Carefully pipette to remove supernatant.
- (5) Wash with 500 μ l of Immunoprecipitation Buffer by gentle vortex. Apply magnetic field then remove supernatant and discard. Repeat wash 2 times.
- (6) Resuspend bead pellet in 30 μ l of 3X SDS Sample Loading Buffer (187.5 mM Tris-HCl (pH 6.8), 6% (w/v) SDS, 30% glycerol, 150 mM DTT, 0.03% (w/v) bromophenol blue, 2% β -mercaptoethanol).
- (7) Incubate sample at 70°C for 5 minutes.
- (8) Apply magnetic field to sample then load supernatant on SDS-PAGE gel and electrophorese.