

NEBNext™ dsDNA
Fragmentase™
Reaction Buffer



1-800-632-7799
info@neb.com
www.neb.com



B0348S 008140116011

B0348S

6 x 1.0 ml (6.0 ml) rxn buffer Lot: 0081401

Store at -20°C Exp: 1/16

Description: New England Biolabs supplies a 10X reaction buffer with all of its enzymes. At a 1X concentration this reaction buffer assures optimal activity of the enzyme.

1X NEBNext dsDNA Fragmentase Reaction Buffer:

20 mM Tris-HCl
10 mM MgCl₂
50 mM NaCl
0.15% Triton X-100
pH 7.5 @ 25°C

Quality Control Assay

16-Hour Incubation: A 50 µl reaction containing this reaction buffer at a 1X concentration and 1 µg of HaeIII digested φX174 RF I DNA incubated for 16 hours resulted in no detectable non-specific nuclease degradation.

Endonuclease Activity: Incubation of this reaction buffer at a 1X concentration with 1 µg φX174 RF I DNA for 4 hours at 37°C in 50 µl reactions resulted in less than 5% conversion to RF II.

Protease Assay: Incubation of at least 1X dsDNA Fragmentase buffer with 0.2 nmol of a standardized mixture of proteins, for 20 hours at 37°C, resulted in no proteolytic activity detected by SDS-PAGE.

Phosphatase Assay: Incubation of 10 µl of at least 1X dsDNA Fragmentase Buffer in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl₂) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

CERTIFICATE OF ANALYSIS

NEBNext™ dsDNA
Fragmentase™
Reaction Buffer



1-800-632-7799
info@neb.com
www.neb.com



B0348S 008140116011

B0348S

6 x 1.0 ml (6.0 ml) rxn buffer Lot: 0081401

Store at -20°C Exp: 1/16

Description: New England Biolabs supplies a 10X reaction buffer with all of its enzymes. At a 1X concentration this reaction buffer assures optimal activity of the enzyme.

1X NEBNext dsDNA Fragmentase Reaction Buffer:

20 mM Tris-HCl
10 mM MgCl₂
50 mM NaCl
0.15% Triton X-100
pH 7.5 @ 25°C

Quality Control Assay

16-Hour Incubation: A 50 µl reaction containing this reaction buffer at a 1X concentration and 1 µg of HaeIII digested φX174 RF I DNA incubated for 16 hours resulted in no detectable non-specific nuclease degradation.

Endonuclease Activity: Incubation of this reaction buffer at a 1X concentration with 1 µg φX174 RF I DNA for 4 hours at 37°C in 50 µl reactions resulted in less than 5% conversion to RF II.

Protease Assay: Incubation of at least 1X dsDNA Fragmentase buffer with 0.2 nmol of a standardized mixture of proteins, for 20 hours at 37°C, resulted in no proteolytic activity detected by SDS-PAGE.

Phosphatase Assay: Incubation of 10 µl of at least 1X dsDNA Fragmentase Buffer in protein phosphatase assay buffer (1 M diethanolamine @ pH 9.8 and 0.5 mM MgCl₂) containing 2.5 mM *p*-nitrophenyl phosphate at 37°C for 4 hours yields no detectable *p*-nitrophenylene anion as determined by spectrophotometric analysis at 405 nm.

CERTIFICATE OF ANALYSIS