E. coli K12 ER2420/pACYC184







E4152S

Lot: 0111703

Exp: 3/22

Store at -20°C

Description: A suspension of *E. coli* strain ER2420 containing the plasmid pACYC184 which has been grown in rich medium and brought to 50% glycerol. ER2420 background: *Eco*K r⁻ m⁻ McrBC⁻ Mrr⁻ Dam⁺ Dcm⁺.

Genotype: F⁻ ara-14 leu fhuA2 Δ (gpt-proA)62 lacY1 glnV44 galK2 rpsL20(Str^R) xyl-5 mtl-1 Δ (mcrC-mrr)_{HR101}

Antibiotic Resistance: The plasmid pACYC184 carries a gene encoding resistance to tetracycline (tet) and a gene encoding resistance to chloramphenicol (chloramphenicol acetyltransferase). To maintain the plasmid, cells should be grown with $25-50 \mu g/ml$ of tet and $35-50 \mu g/ml$ of cam. The plasmid, which is present at a low copy number, cannot be amplified using a standard chloramphenicol procedure but can be amplified using spectinomycin ($50 \mu g/ml$ of growth media).

Notes: Storage at -70°C is recommended for periods longer than 30 days. Avoid repeated freeze/thaw cycles.

References:

- 1. Chang, A. and Cohen, S. (1978) *J. Bacteriol.* 134, 1114–1156.
- 2. Rose, R.E. (1988) Nucleic Acids. Res. 16, 355.

nga.
ISO 9001
Registered
Quality
Management
Medical Device

NEW ENCLAND BIOLARS is a register

NEW ENGLAND BIOLABS $^{\!\circ}$ is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

CERTIFICATE OF ANALYSIS

E. coli K12 ER2420/pACYC184



1-800-632-7799 info@neb.com www.neb.com



E4152S

Lot: 0111703

Exp: 3/22

Store at -20°C

Description: A suspension of *E. coli* strain ER2420 containing the plasmid pACYC184 which has been grown in rich medium and brought to 50% glycerol. ER2420 background: *Eco*K r⁻ m⁻ McrBC⁻ Mrr⁻ Dam⁺ Dcm⁺.

Genotype: F^- ara-14 leu fhuA2 Δ (gpt-proA)62 lacY1 glnV44 galK2 rpsL20(Str^R) xyl-5 mtl-1 Δ (mcrC-mrr)_{HB101}

Antibiotic Resistance: The plasmid pACYC184 carries a gene encoding resistance to tetracycline (tet) and a gene encoding resistance to chloramphenicol (chloramphenicol acetyltransferase). To maintain the plasmid, cells should be grown with 25–50 μ g/ml of tet and 35–50 μ g/ml of cam. The plasmid, which is present at a low copy number, cannot be amplified using a standard chloramphenicol procedure but can be amplified using spectinomycin (50 μ g/ml of growth media).

Notes: Storage at -70°C is recommended for periods longer than 30 days. Avoid repeated freeze/ thaw cycles.

References:

- 1. Chang, A. and Cohen, S. (1978) *J. Bacteriol.* 134, 1114–1156.
- 2. Rose, R.E. (1988) Nucleic Acids. Res.16, 355.







NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.