

## Anti-MBP Monoclonal Antibody



1-800-632-7799  
info@neb.com  
www.neb.com



E8032S 010150817081

# E8032S

**0.05 ml**      **Lot: 0101508**      **Exp: 8/17**  
**1 mg/ml**      **Store at -20°C**

**Description:** Anti-MBP Monoclonal Antibody is a murine anti-maltose binding protein antibody, isotype IgG2a. It is purified from tissue culture supernatant by Protein A affinity chromatography.

**Source:** Tissue culture supernatant from cell line B48.

LABORATORY REAGENT – FOR RESEARCH USE ONLY

Supplied in: 10 mM HEPES pH 7.5, 150 mM NaCl and 50% glycerol.

**Recommended Dilution:** 1/10,000.

**Quality Assurance:** In an ELISA assay, a dilution of 1/10,000 gives a signal of at least 20% of the maximum signal using high concentrations of antibody. The same 1/10,000 dilution gives a strong signal when used to detect maltose-binding protein in Western blots developed with a variety of detection systems. This antibody does not cross-react with other *E. coli* proteins.

### Western Transfer Protocol

Materials:

Transfer apparatus and associated buffers  
Nitrocellulose or PVDF membrane  
TBST (20 mM Tris-Cl, 150 mM NaCl, 0.1% Tween 20)

Blocking Buffer (TBST + 5% Nonfat Dry Milk)  
Anti-MBP Monoclonal Antibody NEB #E8032  
anti-mouse antibody conjugated to peroxidase  
Detection reagent

For a 10 cm x 10 cm gel:

1. Transfer protein from the gel to a nitrocellulose or PVDF membrane following the directions of the transfer apparatus manufacturer. Mark the wells of the gel on the filter with a blunt pencil before removing and discarding the gel.
2. Rinse the membrane with TBST.
3. Incubate the membrane with 25 ml Blocking Buffer for 1 hour at room temperature (or overnight at 4°C) with gentle shaking.
4. Wash the membrane in 25 ml TBST for 5 minutes with gentle shaking, 3 times for 5 minutes each.
5. Add 1 µl of the Anti-MBP Monoclonal Antibody to 10 ml Blocking Buffer (a 1:10,000 dilution). Cover the membrane with the antibody dilution and incubate for 1 hour at room temperature with gentle shaking.
6. Wash the membrane in 25 ml TBST at room temperature with gentle shaking, 3 times for 5 minutes each.

7. Make a dilution of the anti-mouse IgG-peroxidase conjugate in 10 ml Blocking Buffer according to the manufacturer's recommendation, and incubate the membrane in the solution for 1 hour.
8. Wash the membrane in 25 ml TBST at room temperature with gentle shaking, 3 times for 5 minutes each.
9. Follow the manufacturer's directions for detection.

**Note:** Store at -20°C undiluted. May be stored at 4°C diluted in buffer containing 1 mM Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub> or an equivalent antimicrobial agent



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

CERTIFICATE OF ANALYSIS

## Anti-MBP Monoclonal Antibody



1-800-632-7799  
info@neb.com  
www.neb.com



E8032S 010150817081

# E8032S

**0.05 ml**      **Lot: 0101508**      **Exp: 8/17**  
**1 mg/ml**      **Store at -20°C**

**Description:** Anti-MBP Monoclonal Antibody is a murine anti-maltose binding protein antibody, isotype IgG2a. It is purified from tissue culture supernatant by Protein A affinity chromatography.

**Source:** Tissue culture supernatant from cell line B48.

LABORATORY REAGENT – FOR RESEARCH USE ONLY

Supplied in: 10 mM HEPES pH 7.5, 150 mM NaCl and 50% glycerol.

**Recommended Dilution:** 1/10,000.

**Quality Assurance:** In an ELISA assay, a dilution of 1/10,000 gives a signal of at least 20% of the maximum signal using high concentrations of antibody. The same 1/10,000 dilution gives a strong signal when used to detect maltose-binding protein in Western blots developed with a variety of detection systems. This antibody does not cross-react with other *E. coli* proteins.

### Western Transfer Protocol

Materials:

Transfer apparatus and associated buffers  
Nitrocellulose or PVDF membrane  
TBST (20 mM Tris-Cl, 150 mM NaCl, 0.1% Tween 20)

Blocking Buffer (TBST + 5% Nonfat Dry Milk)  
Anti-MBP Monoclonal Antibody NEB #E8032  
anti-mouse antibody conjugated to peroxidase  
Detection reagent

For a 10 cm x 10 cm gel:

1. Transfer protein from the gel to a nitrocellulose or PVDF membrane following the directions of the transfer apparatus manufacturer. Mark the wells of the gel on the filter with a blunt pencil before removing and discarding the gel.
2. Rinse the membrane with TBST.
3. Incubate the membrane with 25 ml Blocking Buffer for 1 hour at room temperature (or overnight at 4°C) with gentle shaking.
4. Wash the membrane in 25 ml TBST for 5 minutes with gentle shaking, 3 times for 5 minutes each.
5. Add 1 µl of the Anti-MBP Monoclonal Antibody to 10 ml Blocking Buffer (a 1:10,000 dilution). Cover the membrane with the antibody dilution and incubate for 1 hour at room temperature with gentle shaking.
6. Wash the membrane in 25 ml TBST at room temperature with gentle shaking, 3 times for 5 minutes each.

7. Make a dilution of the anti-mouse IgG-peroxidase conjugate in 10 ml Blocking Buffer according to the manufacturer's recommendation, and incubate the membrane in the solution for 1 hour.
8. Wash the membrane in 25 ml TBST at room temperature with gentle shaking, 3 times for 5 minutes each.
9. Follow the manufacturer's directions for detection.

**Note:** Store at -20°C undiluted. May be stored at 4°C diluted in buffer containing 1 mM Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub> or an equivalent antimicrobial agent



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

CERTIFICATE OF ANALYSIS