

pNEB206A Linearized Vector



1-800-632-7799
info@neb.com
www.neb.com



N5502S 005161118111

N5502S

1.0 µg **Lot: 0051611** **Exp: 11/18**

20 µg/ml **Store at -20°C**

Description: Linearized pNEB206A Vector is supplied with the USER™ Friendly Cloning Kit (NEB #E5500). pNEB206A has been linearized within the Multiple Cloning Site to produce 8-nucleotide, 3' single-stranded extensions on both vector ends. The single-stranded extensions are not complementary; this prevents the vector termini from re-annealing to form transformable circular DNA and also allows control over the orientation of the inserted PCR product.

Source: pNEB206A is isolated from *E. coli* by a standard purification procedure, digested to completion with Xba I and nicked with N.BbvC IB. The DNA is phenol extracted and resuspended in 10 mM Tris-HCl (pH 8.0), 1 mM EDTA

Supplied in: 10 mM Tris-HCl (pH 8.0), 1 mM EDTA.

Quality Control Assays

A standard USER reaction was performed as described in Appendix V of the USER Friendly Cloning Kit (NEB #E5500) manual [20 ng linearized pNEB206A, 1 µl USER Enzyme and 10 µl (100 ng) of a 950 bp control PCR product amplified using *Taq* DNA Polymerase and primers containing uracil, designed as recommended in the USER Friendly cloning Kit manual]. After transformation into chemically-competent cells (NEB #ER2267 at 5×10^6 c.f.u./µg pNEB206A), 50 µl of the 1 ml outgrowth was spread on Amp + Xgal + IPTG plates. A minimum of 200 colonies were obtained and > 90% of these were white (i.e., contained recombinant molecules).



U.S. Publication No. US-2009-0042258

NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

USER™ is a trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

CERTIFICATE OF ANALYSIS

pNEB206A Linearized Vector



1-800-632-7799
info@neb.com
www.neb.com



N5502S 005161118111

N5502S

1.0 µg **Lot: 0051611** **Exp: 11/18**

20 µg/ml **Store at -20°C**

Description: Linearized pNEB206A Vector is supplied with the USER™ Friendly Cloning Kit (NEB #E5500). pNEB206A has been linearized within the Multiple Cloning Site to produce 8-nucleotide, 3' single-stranded extensions on both vector ends. The single-stranded extensions are not complementary; this prevents the vector termini from re-annealing to form transformable circular DNA and also allows control over the orientation of the inserted PCR product.

Source: pNEB206A is isolated from *E. coli* by a standard purification procedure, digested to completion with Xba I and nicked with N.BbvC IB. The DNA is phenol extracted and resuspended in 10 mM Tris-HCl (pH 8.0), 1 mM EDTA

Supplied in: 10 mM Tris-HCl (pH 8.0), 1 mM EDTA.

Quality Control Assays

A standard USER reaction was performed as described in Appendix V of the USER Friendly Cloning Kit (NEB #E5500) manual [20 ng linearized pNEB206A, 1 µl USER Enzyme and 10 µl (100 ng) of a 950 bp control PCR product amplified using *Taq* DNA Polymerase and primers containing uracil, designed as recommended in the USER Friendly cloning Kit manual]. After transformation into chemically-competent cells (NEB #ER2267 at 5×10^6 c.f.u./µg pNEB206A), 50 µl of the 1 ml outgrowth was spread on Amp + Xgal + IPTG plates. A minimum of 200 colonies were obtained and > 90% of these were white (i.e., contained recombinant molecules).



U.S. Publication No. US-2009-0042258

NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

USER™ is a trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

CERTIFICATE OF ANALYSIS