pTWIN1



1-800-632-7799 info@neb.com www.neb.com

10 µg Lot: 0011512 Exp: 12/17 200 µg/ml Store at -20°C

Description: pTWIN1 is an *E. coli* expression vector which can be used with the IMPACT™ Kit (NEB #E6901). pTWIN vectors are designed for protein purification or for the isolation of proteins with an N-terminal cysteine and/or a C-terminal thioester (1). A polylinker in the vector is designed for the in-frame fusion of a target gene between the modified Ssp DnaB (2) and Mxe GyrA inteins (3). The presence of the chitin binding domain from *Bacillus circulans* (4,5) facilitates purification. The double-stranded vector is 7.375 base pairs in length.

N6951S

derived from the *Synechocystis sp* DnaB intein (154 amino acids) (6) and the other from the Mycobacterium xenopi GyrA intein (198 amino acids) (7).

Supplied in: 10 mM Tris-HCI (pH 8.0), 1 mM EDTA.

Source: pTWIN1 contains two mini-inteins, one

Features of pTWIN1:

- A pBR322 derivative
- The Sapl sites should be used for directional cloning of both the 5' and 3' ends of an insert.

5'...AC TGG GAC TCC ATC GTT TCT ATT ACG GAG ACT GGA GTC GAA GAG GTT TTT

GAT TTG ACT GTG CCA GGA CCA CAT AAC TTT GTC GCG AAT GAC ATC ATT GTA CAC AAC

GGA AGA GCC ATG GGC GGC CGC GAA TTC CTC GAG GGC TCT TCC TGC ATC ACG GGA GAT

GCA CTA GTT GCC CTA CCC GAG GGC GAG TCG GTA CGC ATC GCC GAC ATC GTG CC...3

XhoI

Gly Arg Ala Met Gly Gly Arg Glu Phe Leu Glu Gly Ser Ser Cys Ile Thr Gly

EcoRI

...Ssp DnaB Intein... Val Ala Asn Asp Ile Ile Val His Asn

SapI

NruI

Polylinker Region: pTWIN1

NcoI

...Mxe GyrA Intein...

Ssp DnaB Intein Forward Primer →

- Expression of the fusion gene is under the control of the T7 promotor (8) and is regulated by IPTG due to the presence of a lac gene.
- Expression requires an *E. coli* host that carries the T7 RNA Polymerase gene [e.g., T7 Express Competent E. coli (High Efficiency), (NEB #C2566) or BL21(DE3) Competent E. coli, (NEB #C2527) and derivatives].
- · Origin of DNA replication from the bacteriophage M13 allows for the production of singlestranded DNA by helper phage superinfection of cells bearing the plasmid.

← Intein_

- Thiol-induced cleavage of the Mxe GyrA intein is dependent on the amino acids adjacent to the intein. The amino acid residues M or Y at the C-terminus of the target protein is recommended for use with this intein.
- Controllable cleavage of the Ssp DnaB intein is dependent on the amino acids adjacent to the intein. The amino acid residues CRA or GRA at the N-terminus of the target protein is recommended for use with this intein.
- Ampicillin resistance.

Recommended Buffers

- Cell Lvsis Buffer: 50 mM Tris-HCl (pH 8.5) containing 500 mM NaCl.
- Ssp DnaB Intein Cleavage Buffer: 50 mM Tris-HCl (pH 6.0) containing 500 mM NaCl.
- Mxe GyrA Intein Cleavage Buffer: 50 mM Tris-HCI (pH 8.5) containing 500 mM NaCl and 50 mM 2-mercaptoethanesulfonic acid.

(see other side)

CERTIFICATE OF ANALYSIS

pTWIN1



1-800-632-7799 info@neb.com www.neb.com

Source: pTWIN1 contains two mini-inteins, one derived from the Synechocystis sp DnaB intein (154 amino acids) (6) and the other from the Mycobacterium xenopi GyrA intein (198 amino acids) (7).

NoI

Supplied in: 10 mM Tris-HCI (pH 8.0), 1 mM EDTA.

Features of pTWIN1:

SapI

SpeI

- A pBR322 derivative
- The Sapl sites should be used for directional cloning of both the 5' and 3' ends of an insert.

- Expression of the fusion gene is under the control of the T7 promotor (8) and is regulated by IPTG due to the presence of a lac gene.
- Expression requires an *E. coli* host that carries the T7 RNA Polymerase gene [e.g., T7 Express Competent E. coli (High Efficiency), (NEB #C2566) or BL21(DE3) Competent E. coli, (NEB #C2527) and derivatives].
- Origin of DNA replication from the bacteriophage M13 allows for the production of singlestranded DNA by helper phage superinfection of cells bearing the plasmid.

N6951S

Lot: 0011512 10 µg Exp: 12/17 Store at -20°C 200 μg/ml

Description: pTWIN1 is an E. coli expression vector which can be used with the IMPACT™ Kit (NEB #E6901). pTWIN vectors are designed for protein purification or for the isolation of proteins with an N-terminal cysteine and/or a C-terminal thioester (1). A polylinker in the vector is designed for the in-frame fusion of a target gene between the modified Ssp DnaB (2) and Mxe GyrA inteins (3). The presence of the chitin binding domain from *Bacillus circulans* (4,5) facilitates purification. The double-stranded vector is 7,375 base pairs in length.

Polylinker Region: pTWIN1

SpeI

5'...<u>AC TGG GAC TCC ATC GTT TCT</u> ATT ACG GAG ACT GGA GTC GAA GAG GTT TTT Ssp DnaB Intein Forward Primer →

← Intein _ ...Ssp DnaB Intein... Val Ala Asn Asp Ile Ile Val His Asn GAT TTG ACT GTG CCA GGA CCA CAT AAC TTT GTC GCG AAT GAC ATC ATT GTA CAC AAC NruI

Gly Arg Ala Met Gly Gly Arg Glu Phe Leu Glu Gly Ser Ser Cys Ile Thr Gly GGA AGA GCC ATG GGC GGC CGC GAA TTC CTC GAG GGC TCT TCC TGC ATC ACG GGA GAT NcoI SapI EcoRI XhoI ... Mxe GyrA Intein... GCA CTA GTT GCC CTA CCC GAG GGC GAG TCG GTA CGC ATC GCC GAC ATC GTG CC...3'

at the C-terminus of the target protein is recommended for use with this intein. · Controllable cleavage of the Ssp DnaB intein is dependent on the amino acids adjacent to the intein. The amino acid residues CRA or GRA at

Thiol-induced cleavage of the Mxe GvrA intein

the intein. The amino acid residues M or Y

is dependent on the amino acids adjacent to

- the N-terminus of the target protein is recommended for use with this intein.
- Ampicillin resistance.

Recommended Buffers

- Cell Lysis Buffer: 50 mM Tris-HCl (pH 8.5) containing 500 mM NaCl.
- Ssp DnaB Intein Cleavage Buffer: 50 mM Tris-HCI (pH 6.0) containing 500 mM NaCl.
- · Mxe GyrA Intein Cleavage Buffer: 50 mM Tris-HCI (pH 8.5) containing 500 mM NaCl and 50 mM 2-mercaptoethanesulfonic acid.

(see other side)

CERTIFICATE OF ANALYSIS

References:

- 1. Evans, T.C., Benner, J., and Xu, M.-Q. (1999) The cyclization and polymerization of bacterially expressed proteins using modified self-splicing inteins. J. Biol. Chem. 274, 18359-18363.
- 2. Mathys, S., Evans, T.C., Chute, I.C., Wu, H., Chong, S., Benner, J., Liu, X.-Q. and Xu, M.-Q. (1999). Characterization of a self-splicing miniintein and its conversion into autocatalytic N-and C-terminal cleavage elements: facile production of protein building blocks for protein ligation. Gene 231, 1-13.
- 3. Evans, T.C., Benner, J. and Xu, M.-Q. (1998) Semisynthesis of cytotoxic proteins using a modified protein splicing element. *Protein Sci.* 7. 2256-2264.
- 4. Chong. S., Mersha, F.B., Comb. D.G., Scott. M. E., Landry, D., Vence, L.M., Perler, F.B., Benner, J., Kucera, R.B., Hirvonen, C.A., Pelletier, J.J., Paulus, H. and Xu, M.-Q. (1997) Single-column purification of free recombinant proteins using a self-cleavable affinity tag derived from a protein splicing element. Gene 192, 271-281.

1. Evans, T.C., Benner, J., and Xu, M.-Q. (1999)

inteins. J. Biol. Chem. 274, 18359-18363.

2. Mathys, S., Evans, T.C., Chute, I.C., Wu, H.,

The cyclization and polymerization of bacterially

expressed proteins using modified self-splicing

Chong, S., Benner, J., Liu, X.-Q. and Xu, M.-Q.

(1999). Characterization of a self-splicing mini-

intein and its conversion into autocatalytic N-and

C-terminal cleavage elements: facile production

of protein building blocks for protein ligation.

3. Evans, T.C., Benner, J. and Xu, M.-Q. (1998)

Semisynthesis of cytotoxic proteins using a

4. Chong, S., Mersha, F.B., Comb, D.G., Scott, M.

E., Landry, D., Vence, L.M., Perler, F.B., Benner,

J., Kucera, R.B., Hirvonen, C.A., Pelletier, J.J.,

Paulus, H. and Xu. M.-Q. (1997) Single-column purification of free recombinant proteins using a

self-cleavable affinity tag derived from a protein

splicing element. Gene 192, 271-281.

modified protein splicing element. Protein Sci. 7,

Page 2 (N6951S)

References:

Gene 231, 1-13,

2256-2264.

- 5. Watanabe, T., Ito, Y., Yamada, T., Hashimoto, M., Sekine, S. and Tanaka, H. (1994). The role of the C-terminal domain and type III domains of chitinase A1 from Bacillus circulans WL-12 in chitin degradation. J. Bacteriol. 176, 4465-4472.
 - 6. Wu, H., Xu, M.-Q. and Liu, X.-Q. (1998) Protein trans-splicing and functional mini-inteins of a cyanobacterial DnaB intein. Biochem. Biophys. Acta 1387, 422-432.
 - 7. Telenti, A., Southworth, M., Alcaide, F., Daugelat, S., Jacobs, W.R. Jr. and Perler, F.B. (1997) The Mycobacterium xenopi GyrA protein splicing element: Characterization of a minimal intein. J. Bacteriol. 179, 6378-6382.
 - 8. Dubendorff, J.W. and Studier, F.W. (1991) Controlling basal expression in an inducible T7 expression system by blocking the target T7 promoter with lac repressor. J. Mol. Biol. 219, 45-59.

Additional information such as vector sequences and frequently asked questions, are available at www.neb.com.

5. Watanabe, T., Ito, Y., Yamada, T., Hashimoto, M., Sekine, S. and Tanaka, H. (1994). The role of the C-terminal domain and type III domains of chitinase A1 from Bacillus circulans WL-12 in chitin degradation. J. Bacteriol. 176, 4465-

- 6. Wu, H., Xu, M.-Q. and Liu, X.-Q. (1998) Protein trans-splicing and functional mini-inteins of a cyanobacterial DnaB intein. Biochem. Biophys, Acta 1387, 422-432.
- Daugelat, S., Jacobs, W.R. Jr. and Perler, F.B. (1997) The Mycobacterium xenopi GyrA protein splicing element: Characterization of a minimal intein. J. Bacteriol. 179, 6378-6382.
- 8. Dubendorff, J.W. and Studier, F.W. (1991) Controlling basal expression in an inducible T7 expression system by blocking the target T7 promoter with lac repressor. J. Mol. Biol. 219, 45-59.

and frequently asked questions, are available at www.neb.com.








U.S. Patent Nos. 6,569,669

NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

IMPACT™ is a trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals

Notice to Buyer/User: The buyer/user has a non-exclusive license for Research Purposes Only. See RESEARCH USE ASSURANCE STATE-MENT on neb.com





7. Telenti, A., Southworth, M., Alcaide, F.,

Additional information such as vector sequences







U.S. Patent Nos. 6,569,669

NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

IMPACT™ is a trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals

Notice to Buyer/User: The buyer/user has a non-exclusive license for Research Purposes Only. See RESEARCH USE ASSURANCE STATE-MENT on neh com