

β 1-4 Galactosidase S



P0745S



400 units 8,000 U/ml Lot: 0011511

RECOMBINANT Store at -20°C Exp: 11/17

Description: β 1-4 Galactosidase S is a highly specific exoglycosidase that catalyzes the hydrolysis of β 1-4 linked galactose residues from oligosaccharides.

Specificity:

Source: Cloned from *Streptococcus pneumoniae*



Gal R = any sugar

and expressed in *E. coli* (1).

Supplied in: 50 mM NaCl, 20 mM Tris-HCl (pH 7.5 @ 25°C) and 1 mM EDTA.

Reagents Supplied with Enzyme:

10X GlycoBuffer 1
(0.5 M Sodium Acetate, pH 5.5 @ 25°C and 50 mM CaCl₂)

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal, β -D-galactose from 1 nmol Gal β 1-4GlcNAc β 1-3Gal β 1-4Glc-7-amino-4-methyl-coumarin (AMC), in 1 hour at 37°C in a total reaction volume of 10 μ l.

Unit Definition Assay: Two fold dilutions of β 1-4 Galactosidase S are incubated with 1 nmol AMC-labeled substrate and 1X GlycoBuffer 1 in a 10 μ l reaction. The reaction mix is incubated at 37°C for 1 hour. Separation of reaction products are visualized via thin layer chromatography (2).

Specific Activity: ~123,000 units/mg.

Molecular Weight: 231,000 daltons.

Quality Assurance: No contaminating exoglycosidase or endoglycosidase F1, F2 or F3 activity could be detected. No contaminating proteolytic activity could be detected.

Quality Controls

Glycosidase Assays: 80 units of β 1-4 Galactosidase S were incubated with 0.1 mM of fluorescently-labeled oligosaccharides and glycopeptides, in a 10 μ l reaction for 20 hours at 37°C. The reaction products were analyzed by TLC for digestion of substrate.

No other glycosidase activities were detected (ND) with the following substrates:

β -N-Acetylglucosaminidase:
GlcNAc β 1-4GlcNAc β 1-4GlcNAc-AMC ND

β -N-Acetylgalactosaminidase:
GalNAc β 1-4Gal β 1-4Glc-AMC ND

α -N-Acetylgalactosaminidase:
GalNAc α 1-3(Fuc α 1-2)Gal β 1-4Glc-AMC ND

α -Fucosidase:
Fuc α 1-2Gal β 1-4Glc-AMC ND

β -Galactosidase:
Gal β 1-3GlcNAc β 1-4Gal β 1-4Glc-AMC ND

α -Galactosidase:
Gal α 1-3Gal β 1-4Gal-AMC ND
Gal α 1-6Gal α 1-6Glc α 1-2Fru-AMC ND

α -Neuraminidase:
Neu5Aca2-3Gal β 1-3GlcNAc β 1-3Gal β 1-4Glc-AMC ND

α -Mannosidase:
Man α 1-3Man β 1-4GlcNAc-AMC ND
Man α 1-6Man α 1-6(Man α 1-3)Man-AMC ND

α -Glucosidase:
Glc α 1-6Glc α 1-4Glc-AMC ND

β -Xylosidase:
Xyl β 1-4Xyl β 1-4Xyl β 1-4Xyl-AMC ND

β -Mannosidase:
Man β 1-4Man β 1-4Man-AMC ND

Endo F₁, F₂, H:
Dansylated invertase high mannose. ND

(see other side)

CERTIFICATE OF ANALYSIS

β 1-4 Galactosidase S



P0745S



400 units 8,000 U/ml Lot: 0011511

RECOMBINANT Store at -20°C Exp: 11/17

Description: β 1-4 Galactosidase S is a highly specific exoglycosidase that catalyzes the hydrolysis of β 1-4 linked galactose residues from oligosaccharides.

Specificity:

Source: Cloned from *Streptococcus pneumoniae*



Gal R = any sugar

and expressed in *E. coli* (1).

Supplied in: 50 mM NaCl, 20 mM Tris-HCl (pH 7.5 @ 25°C) and 1 mM EDTA.

Reagents Supplied with Enzyme:

10X GlycoBuffer 1
(0.5 M Sodium Acetate, pH 5.5 @ 25°C and 50 mM CaCl₂)

Unit Definition: One unit is defined as the amount of enzyme required to cleave > 95% of the terminal, β -D-galactose from 1 nmol Gal β 1-4GlcNAc β 1-3Gal β 1-4Glc-7-amino-4-methyl-coumarin (AMC), in 1 hour at 37°C in a total reaction volume of 10 μ l.

Unit Definition Assay: Two fold dilutions of β 1-4 Galactosidase S are incubated with 1 nmol AMC-labeled substrate and 1X GlycoBuffer 1 in a 10 μ l reaction. The reaction mix is incubated at 37°C for 1 hour. Separation of reaction products are visualized via thin layer chromatography (2).

Specific Activity: ~123,000 units/mg.

Molecular Weight: 231,000 daltons.

Quality Assurance: No contaminating exoglycosidase or endoglycosidase F1, F2 or F3 activity could be detected. No contaminating proteolytic activity could be detected.

Quality Controls

Glycosidase Assays: 80 units of β 1-4 Galactosidase S were incubated with 0.1 mM of fluorescently-labeled oligosaccharides and glycopeptides, in a 10 μ l reaction for 20 hours at 37°C. The reaction products were analyzed by TLC for digestion of substrate.

No other glycosidase activities were detected (ND) with the following substrates:

β -N-Acetylglucosaminidase:
GlcNAc β 1-4GlcNAc β 1-4GlcNAc-AMC ND

β -N-Acetylgalactosaminidase:
GalNAc β 1-4Gal β 1-4Glc-AMC ND

α -N-Acetylgalactosaminidase:
GalNAc α 1-3(Fuc α 1-2)Gal β 1-4Glc-AMC ND

α -Fucosidase:
Fuc α 1-2Gal β 1-4Glc-AMC ND

β -Galactosidase:
Gal β 1-3GlcNAc β 1-4Gal β 1-4Glc-AMC ND

α -Galactosidase:
Gal α 1-3Gal β 1-4Gal-AMC ND
Gal α 1-6Gal α 1-6Glc α 1-2Fru-AMC ND

α -Neuraminidase:
Neu5Aca2-3Gal β 1-3GlcNAc β 1-3Gal β 1-4Glc-AMC ND

α -Mannosidase:
Man α 1-3Man β 1-4GlcNAc-AMC ND
Man α 1-6Man α 1-6(Man α 1-3)Man-AMC ND

α -Glucosidase:
Glc α 1-6Glc α 1-4Glc-AMC ND

β -Xylosidase:
Xyl β 1-4Xyl β 1-4Xyl β 1-4Xyl-AMC ND

β -Mannosidase:
Man β 1-4Man β 1-4Man-AMC ND

(see other side)

CERTIFICATE OF ANALYSIS

Endo F₂, F₃:
Dansylated fibrinogen biantennary. ND

Protease Assay: After incubation of 400 units of β 1-4 Galactosidase S with 0.2 nmol of a standard mixture of proteins in a 20 μ l reaction, for 20 hours at 37°C, no proteolytic activity could be detected by SDS-PAGE.

Physical Purity: Purified to > 95% homogeneity as determined by SDS-PAGE analysis using Coomassie Blue detection.

Heat Inactivation: 75°C for 10 minutes.

Reaction Conditions: Optimal incubation times and enzyme concentrations must be determined empirically for a particular substrate. Typical reaction

conditions are as follows:

1. Combine 1 μ g of glycoprotein or 100 nM of oligosaccharide and H₂O (if necessary) to make a 9 μ l total reaction volume.
2. Add 1 μ l of 10X GlycoBuffer 1 to make a 10 μ l total reaction volume.
3. Add 1 μ l of β 1-4 Galactosidase S.
4. Incubate at 37°C for 1 hour.

Notes on Use:

- Reactions may be scaled-up linearly to accommodate larger reaction volumes.
- The amount of exoglycosidase enzyme required varies when different substrates are used. Start with 1–2 μ l for 1 μ g of glycoprotein or 100 nM of oligosaccharide for one hour in a 10–25 μ l reaction. If there is still undigested material, let the reaction go overnight.

References:

1. Chen, M. New England Biolabs, Inc., unpublished results.

Page 2 (P0745)

Endo F₁, F₂, H:
Dansylated invertase high mannose. ND

Endo F₂, F₃:
Dansylated fibrinogen biantennary. ND

Protease Assay: After incubation of 400 units of β 1-4 Galactosidase S with 0.2 nmol of a standard mixture of proteins in a 20 μ l reaction, for 20 hours at 37°C, no proteolytic activity could be detected by SDS-PAGE.

Physical Purity: Purified to > 95% homogeneity as determined by SDS-PAGE analysis using Coomassie Blue detection.

Heat Inactivation: 75°C for 10 minutes.

Reaction Conditions: Optimal incubation times and enzyme concentrations must be determined empirically for a particular substrate. Typical reaction conditions are as follows:

1. Combine 1 μ g of glycoprotein or 100 nM of oligosaccharide and H₂O (if necessary) to make a 9 μ l total reaction volume.
2. Add 1 μ l of 10X GlycoBuffer 1 to make a 10 μ l total reaction volume.
3. Add 1 μ l of β 1-4 Galactosidase S.
4. Incubate at 37°C for 1 hour.

Notes on Use:

- Reactions may be scaled-up linearly to accommodate larger reaction volumes.
- The amount of exoglycosidase enzyme required varies when different substrates are used. Start with 1–2 μ l for 1 μ g of glycoprotein or 100 nM of oligosaccharide for one hour in a 10–25 μ l reaction. If there is still undigested material, let the reaction go overnight.

2. Wong-Madden, S. T. and Landry, D. (1995) *Glycobiology* 5, 19–28.



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

References:

1. Chen, M. New England Biolabs, Inc., unpublished results.
2. Wong-Madden, S. T. and Landry, D. (1995) *Glycobiology* 5, 19–28.



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

Page 2 (P0745)