

240 County Road Ipswich, MA 01938-2723 Tel 978-927-5054 Fax 978-921-1350 www.neb.com info@neb.com

## New England Biolabs Product Specification

Product Name:	ProtoScript® II Reverse Transcriptase
Catalog #:	M0368S/L/X
Concentration:	200,000 units/ml
Unit Definition:	One unit is defined as the amount of enzyme that will incorporate 1 nmol of dTTP into an acid-insoluble form in 10 minutes at 37°C.
Shelf Life:	24 months
Storage Temp:	-20°C
Storage Conditions:	20 mM Tris-HCl, 100 mM NaCl, 1 mM DTT, 0.1 mM EDTA, 0.01 % IGEPAL® CA-630, 50 % Glycerol, (pH 7.5 @ 25°C)
Specification Version:	PS-M0368S/L/X v2.0
Effective Date:	12 Feb 2020

## Assay Name/Specification (minimum release criteria)

Endonuclease Activity (Nicking) - A 50  $\mu$ l reaction in ProtoScript® II Reverse Transcriptase Reaction Buffer containing 1  $\mu$ g of supercoiled PhiX174 DNA and a minimum of 200 units of ProtoScript® II Reverse Transcriptase incubated for 4 hours at 37°C results in <10% conversion to the nicked form as determined by agarose gel electrophoresis.

**Exonuclease Activity (Radioactivity Release)** - A 50  $\mu$ l reaction in ProtoScript® II Reverse Transcriptase Reaction Buffer containing 1  $\mu$ g of a mixture of single and double-stranded [<sup>3</sup>H] *E. coli* DNA and a minimum of 200 units of ProtoScript® II Reverse Transcriptase incubated for 4 hours at 37°C releases <0.1% of the total radioactivity.

Non-Specific DNase Activity (16 Hour) - A 50 µl reaction in NEBuffer 2 containing 1 µg of T3 or T7 DNA in addition to a reaction containing Lambda-HindIII DNA and a minimum of 200 units of ProtoScript® II Reverse Transcriptase incubated for 16 hours at 37°C results in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis.

**Protein Purity Assay (SDS-PAGE)** - ProtoScript<sup>®</sup> II Reverse Transcriptase is  $\geq$  99% pure as determined by SDS-PAGE analysis using Coomassie Blue detection.

**qPCR DNA Contamination (***E. coli* **Genomic)** - A minimum of 200 units of ProtoScript® II Reverse Transcriptase is screened for the presence of *E. coli* genomic DNA using SYBR® Green qPCR with primers specific for the *E. coli* 16S rRNA locus. Results are quantified using a standard curve generated from purified *E. coli* genomic DNA. The measured level of *E. coli* genomic DNA contamination is  $\leq 1 E$ . *coli* genome.



PS-M0368S/L/X v2.0 Page 1 of 2



240 County Road Ipswich, MA 01938-2723 Tel 978-927-5054 Fax 978-921-1350

www.neb.com info@neb.com

## New England Biolabs Product Specification

Assay Name/Specification (minimum release criteria)

**RNAse Activity Assay (4 Hour Digestion)** - A 10  $\mu$ l reaction in ProtoScript® II Reverse Transcriptase Reaction Buffer containing 40 ng of a 300 base single-stranded RNA and a minimum of 1  $\mu$ l of ProtoScript® II Reverse Transcriptase is incubated at 37°C. After incubation for 4 hours, >90% of the substrate RNA remains intact as determined by gel electrophoresis using fluorescent detection.

One or more products referenced in this document may be covered by a 3rd-party trademark. Please visit <u>www.neb.com/trademarks</u> for additional information.

Date 12 Feb 2020

Derek Robinson Director, Quality Control



PS-M0368S/L/X v2.0 Page 2 of 2