240 County Road Ipswich, MA 01938-2723 Tel 978-927-5054 Fax 978-921-1350

Date

31 Mar 2022

www.neb.com info@neb.com

New England Biolabs Product Specification

Product Name: BbvCI
Catalog #: R0601S/L
Concentration: 2,000 units/ml

Unit Definition:

One unit is defined as the amount of enzyme required to digest 1 µg of Lambda DNA in rCutSmart™ Buffer in 1 hour at

37°C in a total reaction volume of 50 μ l.

Shelf Life: 12 months
Storage Temp: -20°C

Storage Conditions: 10 mM Tris-HCl, 300 mM NaCl, 1 mM DTT, 0.1 mM EDTA, 50% Glycerol, 500 µg/ml rAlbumin (pH 7.4 @

25°C,

Specification Version: PS-R0601S/L v3.0

Effective Date: 31 Mar 2022

Assay Name/Specification (minimum release criteria)

Ligation and Recutting (Terminal Integrity) - After a 2-fold over-digestion of Lambda DNA with BbvCI, <10% of the DNA fragments can be ligated with T4 DNA ligase in 16 hours at 16°C. Of these ligated fragments, >95% can be recut with BbvCI.

Protein Purity Assay (SDS-PAGE) - BbvCI is ≥ 95% pure as determined by SDS-PAGE analysis using Coomassie Blue detection.

qPCR DNA Contamination (*E. coli* Genomic) - A minimum of 2 units of BbvCI is screened for the presence of *E. coli* genomic DNA using SYBR® Green qPCR with primers specific for the *E. coli* 16S rRNA locus. Results are quantified using a standard curve generated from purified *E. coli* genomic DNA. The measured level of *E. coli* genomic DNA contamination is ≤ 1 *E. coli* genome.

Exonuclease Activity (Radioactivity Release) - A 50 μl reaction in rCutSmartTM Buffer containing 1 μg of a mixture of single and double-stranded [³H] *E. coli* DNA and a minimum of 10 units of BbvCI incubated for 4 hours at 37°C releases <0.1% of the total radioactivity.

Non-Specific DNase Activity (16 hour) - A 50 μ l reaction in rCutSmartTM Buffer containing 1 μ g of Lambda DNA and a minimum of 2 units of BbvCI incubated for 16 hours at 37°C results in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis. NOTE: although no nuclease degradation is detected under these conditions, extended incubations and/or high concentrations of this enzyme may result in star activity. See the product FAQ for recommended reaction conditions for this enzyme.

One or more products referenced in this document may be covered by a 3rd-party trademark. Please visit www.neb.com/trademarks for additional information.

Derek Robinson

Director, Quality Control





