

New England Biolabs Product Specification

| | |
|-------------------------------|--|
| <i>Product Name:</i> | <i>Magnesium Sulfate (MgSO₄) Solution</i> |
| <i>Catalog #:</i> | <i>B1003S</i> |
| <i>Concentration:</i> | <i>100 mM</i> |
| <i>Shelf Life:</i> | <i>60 months</i> |
| <i>Storage Temp:</i> | <i>-20°C</i> |
| <i>Composition (1X):</i> | <i>100 mM MgSO₄</i> |
| <i>Specification Version:</i> | <i>PS-B1003S v1.0</i> |
| <i>Effective Date:</i> | <i>10 Aug 2016</i> |

Assay Name/Specification (minimum release criteria)

Conductivity (buffers/solutions) - The conductivity of 100 mM Magnesium Sulfate (MgSO₄) Solution is between 8.5 and 10.5 mS/cm at 25°C.

Endonuclease Activity (Nicking) - A 50 µl reaction in NEBuffer 2 containing 1 µg of supercoiled PhiX174 DNA and a minimum of 5 µl of Magnesium Sulfate (MgSO₄) Solution incubated for 4 hours at 37°C results in <10% conversion to the nicked form as determined by agarose gel electrophoresis.

Non-Specific DNase Activity (16 Hour) - A 50 µl reaction in NEBuffer 2 containing 1 µg of T3 DNA in addition to a reaction containing Lambda-HindIII DNA and a minimum of 5 µl of Magnesium Sulfate (MgSO₄) Solution incubated for 16 hours at 37°C results in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis.

PCR Amplification (5.0 kb Lambda DNA, Mg²⁺) - A 50 µl reaction in ThermoPol II® (Mg-free) Reaction Buffer containing 2 mM Magnesium Sulfate (MgSO₄) Solution in the presence of 200 µM dNTPs and 0.2 µM primers containing 5 ng Lambda DNA with 1.25 units of *Taq* DNA Polymerase for 25 cycles of PCR amplification results in the expected 5.0 kb product.

pH (buffers/solutions) - The pH of 100 mM Magnesium Sulfate (MgSO₄) Solution is between pH 5.3 and 5.7 at 25°C.

Phosphatase Activity (pNPP, Buffer) - A 200 µl reaction in 1M Diethanolamine @ pH 9.8 and 0.5 mM MgCl₂ containing 2.5 mM *p*-Nitrophenyl Phosphate (pNPP) and a minimum of 20 µl Magnesium Sulfate (MgSO₄) Solution incubated for 4 hours at 37°C yields <0.0001 unit of alkaline phosphatase activity as determined by spectrophotometric analysis.

qPCR DNA Contamination (*E. coli* Genomic) - A minimum of 1 µl of Magnesium Sulfate (MgSO₄) Solution is screened for the presence of *E. coli* genomic DNA using SYBR® Green qPCR with primers specific for the *E. coli* 16S rRNA locus. Results are quantified using a standard curve generated from purified *E. coli* genomic DNA. The measured level of *E. coli* genomic DNA contamination is ≤ 1 *E. coli* genome.



New England Biolabs Product Specification

| Assay Name/Specification (minimum release criteria) |
|---|
|---|

| |
|--|
| RNase Activity (Extended Digestion) - A 10 µl reaction in NEBuffer 4 containing 40 ng of a 300 base single-stranded RNA and a minimum of 1 µl of Magnesium Sulfate (MgSO ₄) Solution is incubated at 37°C. After incubation for 16 hours, >90% of the substrate RNA remains intact as determined by gel electrophoresis using fluorescent detection. |
|--|



Date 10 Aug 2016

Derek Robinson
Director of Quality Control

