

## New England Biolabs Product Specification

*Product Name:* SHuffle<sup>®</sup> T7 Express lysY Competent *E. coli*  
*Catalog #:* C3030J  
*Shelf Life:* 12 months  
*Storage Temp:* -80°C  
*Specification Version:* PS-C3030J v2.0  
*Effective Date:* 02 May 2023

### Assay Name/Specification (minimum release criteria)

**Antibiotic Resistance (Chloramphenicol)** - 15 µl of untransformed SHuffle<sup>®</sup> T7 Express lysY Competent *E. coli* streaked onto a LB or Rich Broth plate containing Chloramphenicol will form colonies after incubation for 16 hours at 37°C.

**Antibiotic Resistance (Nitrofurantoin)** - 15 µl of untransformed SHuffle<sup>®</sup> T7 Express lysY Competent *E. coli* streaked onto a LB or Rich Broth plate containing Nitrofurantoin will form colonies after incubation for 16 hours at 37°C.

**Antibiotic Resistance (Spectinomycin)** - 15 µl of untransformed SHuffle<sup>®</sup> T7 Express lysY Competent *E. coli* streaked onto a LB or Rich Broth plate containing Spectinomycin will form colonies after incubation for 16 hours at 37°C.

**Antibiotic Resistance (Streptomycin)** - 15 µl of untransformed SHuffle<sup>®</sup> T7 Express lysY Competent *E. coli* streaked onto a LB or Rich Broth plate containing Streptomycin will form colonies after incubation for 16 hours at 37°C.

**Antibiotic Sensitivity (Ampicillin)** - 15 µl of untransformed SHuffle<sup>®</sup> T7 Express lysY Competent *E. coli* streaked onto a LB or Rich Broth plate containing Ampicillin will not form colonies after incubation for 16 hours at 37°C.

**Antibiotic Sensitivity (Kanamycin)** - 15 µl of untransformed SHuffle<sup>®</sup> T7 Express lysY Competent *E. coli* streaked onto a LB or Rich Broth plate containing Kanamycin will not form colonies after incubation for 16 hours at 37°C.

**Antibiotic Sensitivity (Tetracycline)** - 15 µl of untransformed SHuffle<sup>®</sup> T7 Express lysY Competent *E. coli* streaked onto a LB or Rich Broth plate containing Tetracycline will not form colonies after incubation for 16 hours at 37°C.

**Phage Resistance (Φ 80)** - 15 µl of untransformed SHuffle<sup>®</sup> T7 Express lysY Competent *E. coli* streaked onto a Rich Broth plate does not support plaque formation by phage Φ 80 after incubation for 16 hours at 37°C.

**Functional Testing (Disulfide Bond Formation)** - The nuclease NucA requires disulfide bonds for its stability. When expressed at 30°C in *E. coli*, NucA is toxic to cells only in its oxidized disulfide-bonded state. Transformation of SHuffle<sup>®</sup> T7 Express lysY Competent *E. coli* using 100 pg of plasmid that expresses a MBP-NucA fusion results in < 1% of the colonies when compared to a control transformation of its wild type parent strain NEB 10-beta after 0.4 mM IPTG induction at 30°C.



---

## New England Biolabs Product Specification

<b>Assay Name/Specification (minimum release criteria)</b>
--

<b>Transformation Efficiency</b> - 50 µl of SHuffle <sup>®</sup> T7 Express lysY Competent <i>E. coli</i> cells were transformed with 100 pg of pUC19 DNA using the transformation protocol provided. Incubation overnight on LB-Ampicillin plates at 37°C resulted in >1 x 10 <sup>7</sup> cfu/µg of DNA.
--

*One or more products referenced in this document may be covered by a 3rd-party trademark.  
Please visit [www.neb.com/trademarks](http://www.neb.com/trademarks) for additional information.*



Date 02 May 2023

Lauren Brown  
Quality Approver

