

New England Biolabs Product Specification

| | |
|------------------------|--|
| <i>Product Name:</i> | <i>NEBNext[®] Immune Sequencing Kit (Mouse)</i> |
| <i>Catalog #:</i> | <i>E6330L</i> |
| <i>Kit Components:</i> | <i>NEBNext[®] IS RT Buffer (4X) (E6321) — Store at -20°C</i> <i>NEBNext[®] Cell Lysis Buffer (E6322) — Store at -20°C</i> <i>dNTP Solution Mix (E6323) — Store at -20°C</i> <i>NEBNext[®] IS RT Primer (E6324) — Store at -20°C</i> <i>NEBNext[®] IS TS Oligos (E6326) — Store at -20°C</i> <i>NEBNext[®] First Strand Synthesis Enzyme Mix (E7761) — Store at -20°C</i> <i>Q5[®] Hot Start High-Fidelity DNA Polymerase (E6327) — Store at -20°C</i> <i>Q5[®] Reaction Buffer (5X) (E6328) — Store at -20°C</i> <i>NEBNext[®] IS PCR2 Universal Primer (E6329) — Store at -20°C</i> <i>TE (0.1X) (E6331) — Store at -20°C</i> <i>Nuclease-free Water (E6332) — Store at -20°C</i> <i>NEBNext[®] i701 Primer (E6333) — Store at -20°C</i> <i>NEBNext[®] i702 Primer (E6334) — Store at -20°C</i> <i>NEBNext[®] i703 Primer (E6335) — Store at -20°C</i> <i>NEBNext[®] i704 Primer (E6336) — Store at -20°C</i> <i>NEBNext[®] i705 Primer (E6337) — Store at -20°C</i> <i>NEBNext[®] i706 Primer (E6338) — Store at -20°C</i> <i>NEBNext[®] i707 Primer (E6347) — Store at -20°C</i> <i>NEBNext[®] i708 Primer (E6348) — Store at -20°C</i> <i>NEBNext[®] i709 Primer (E6349) — Store at -20°C</i> <i>NEBNext[®] i710 Primer (E6352) — Store at -20°C</i> <i>NEBNext[®] i711 Primer (E6353) — Store at -20°C</i> <i>NEBNext[®] i712 Primer (E6354) — Store at -20°C</i> <i>NEBNext[®] i501 Primer (E6339) — Store at -20°C</i> <i>NEBNext[®] i502 Primer (E6340) — Store at -20°C</i> <i>NEBNext[®] i503 Primer (E6341) — Store at -20°C</i> <i>NEBNext[®] i504 Primer (E6342) — Store at -20°C</i> <i>NEBNext[®] i505 Primer (E6343) — Store at -20°C</i> <i>NEBNext[®] i506 Primer (E6344) — Store at -20°C</i> <i>NEBNext[®] i507 Primer (E6345) — Store at -20°C</i> <i>NEBNext[®] i508 Primer (E6346) — Store at -20°C</i> |



New England Biolabs Product Specification

Streptavidin Magnetic Beads (E6355) — Store at 4°C
NEBNext[®] IS Bind/Wash Buffer (2X) (E6356) — Store at 4°C
0.1% Tween (E6357) — Store at 4°C
NEBNext[®] Sample Purification Beads (E6358) — Store at 25°C
NEBNext[®] IS BCR Primers (Mouse) (E2624) — Store at -20°C
NEBNext[®] IS TCR Primers (Mouse) (E2625) — Store at -20°C

Shelf Life: 18 months
Storage Temp: Multi-temp
Specification Version: PS-E6330L v1.0
Effective Date: 14 Jun 2021

Assay Name/Specification (minimum release criteria)

Functional Testing (Library Construction) - NEBNext[®] Immune Sequencing Kit (Mouse) is functionally validated using commercially available mouse spleen RNA. Libraries are constructed with the new and previous kit lots at the highest and lowest RNA input using both BCR and TCR primers. All libraries are sequenced on the same flowcell and compared with expected thresholds and between lots across a variety of metrics. Library assessment is based on metrics including library yield, read quality, distribution of reads per UMI, primer usage, assembled sequence length, and distribution of clonotypes.

* **Individual Product Component Note** - Standard Quality Control Tests are performed for each component included in NEBNext[®] Immune Sequencing Kit (Mouse) and meet the designated specifications.

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 14 Jun 2021

Derek Robinson
Director, Quality Control

