

---

## New England Biolabs Product Specification

*Product Name:* Monarch<sup>®</sup> DNA Capture Beads  
*Catalog #:* T3005L  
*Shelf Life:* 60 months  
*Storage Temp:* 25°C  
*Specification Version:* PS-T3005L v2.0  
*Effective Date:* 31 Mar 2022

<b>Assay Name/Specification (minimum release criteria)</b>
--

<b>Functional Testing (Bead Binding Capacity)</b> - Eight (8) pairs of beads are tested with 10 µg of Lambda DNA using a modified High Molecular Weight DNA Extraction protocol, resulting in ≥80% recovery in ≥80% of samples.
---

One or more products referenced in this document may be covered by a 3rd-party trademark.  
Please visit [www.neb.com/trademarks](http://www.neb.com/trademarks) for additional information.



Date 31 Mar 2022

---

Derek Robinson  
Director, Quality Control

