

Phusion® GC Buffer Pack



1-800-632-7799
info@neb.com
www.neb.com



B0519S 001140717071

B0519S

6.0 ml **Lot: 0011407** **Exp: 7/17**
5X Concentration **Store at -20°C**

Description: New England Biolabs supplies a unique 5X reaction buffer and supplemental reagents for use with Thermo Scientific® Phusion High-Fidelity DNA Polymerase and associated products, all of which are manufactured and quality controlled at New England Biolabs.

Reagents Supplied:

| | |
|-------------------------------|------------|
| 5X Phusion GC Reaction Buffer | 4 x 1.5 ml |
| 50 mM MgCl ₂ | 2 x 1.5 ml |
| 100% DMSO | 0.5 ml |

Reaction Conditions: At a 1X concentration the reaction buffer contains 1.5 mM MgCl₂ and assures optimal activity for Phusion DNA Polymerase. Additional MgCl₂ (50 mM) is provided for those amplicons that require higher Mg⁺⁺ concentrations.

Quality Assurance

The reaction buffer and supplements are free of detectable nucleases.

Phusion® DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific.

PHUSION® is a registered trademark and property of Thermo Fisher Scientific.

THERMO SCIENTIFIC® is a registered trademark of Thermo Fisher Scientific.

CERTIFICATE OF ANALYSIS

Phusion® GC Buffer Pack



1-800-632-7799
info@neb.com
www.neb.com



B0519S 001140717071

B0519S

6.0 ml **Lot: 0011407** **Exp: 7/17**
5X Concentration **Store at -20°C**

Description: New England Biolabs supplies a unique 5X reaction buffer and supplemental reagents for use with Thermo Scientific® Phusion High-Fidelity DNA Polymerase and associated products, all of which are manufactured and quality controlled at New England Biolabs.

Reagents Supplied:

| | |
|-------------------------------|------------|
| 5X Phusion GC Reaction Buffer | 4 x 1.5 ml |
| 50 mM MgCl ₂ | 2 x 1.5 ml |
| 100% DMSO | 0.5 ml |

Reaction Conditions: At a 1X concentration the reaction buffer contains 1.5 mM MgCl₂ and assures optimal activity for Phusion DNA Polymerase. Additional MgCl₂ (50 mM) is provided for those amplicons that require higher Mg⁺⁺ concentrations.

Quality Assurance

The reaction buffer and supplements are free of detectable nucleases.

Phusion® DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific.

PHUSION® is a registered trademark and property of Thermo Fisher Scientific.

THERMO SCIENTIFIC® is a registered trademark of Thermo Fisher Scientific.

CERTIFICATE OF ANALYSIS