

Detergent-free
Phusion® HF
Buffer Pack



1-800-632-7799
info@neb.com
www.neb.com



B0520S 001130916091

B0520S

6.0 ml **Lot: 0011309** **Exp: 9/16**
5X Concentration **Store at -20°C**

Description: New England Biolabs supplies a unique 5X reaction buffer and supplemental reagents for use with Thermo Scientific® Phusion High-Fidelity DNA Polymerase and associated products, all of which are manufactured and quality controlled at New England Biolabs.

Reagents Supplied:

5X Detergent-free Phusion HF Reaction Buffer	4 x 1.5 ml
50 mM MgCl ₂	2 x 1.5 ml
100% DMSO	0.5 ml

Reaction Conditions: Use at a 1X final concentration in the Phusion reaction. The 1X buffer contains 1.5 mM MgCl₂. Additional MgCl₂ (50 mM) is provided for those amplicons that require higher Mg⁺⁺ concentrations.

Quality Assurance

The reaction buffer and supplements are free of detectable nucleases.

Phusion® DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific.

PHUSION® is a registered trademark and property of Thermo Fisher Scientific.

THERMO SCIENTIFIC® is a registered trademark of Thermo Fisher Scientific.

CERTIFICATE OF ANALYSIS

Detergent-free
Phusion® HF
Buffer Pack



1-800-632-7799
info@neb.com
www.neb.com



B0520S 001130916091

B0520S

6.0 ml **Lot: 0011309** **Exp: 9/16**
5X Concentration **Store at -20°C**

Description: New England Biolabs supplies a unique 5X reaction buffer and supplemental reagents for use with Thermo Scientific® Phusion High-Fidelity DNA Polymerase and associated products, all of which are manufactured and quality controlled at New England Biolabs.

Reagents Supplied:

5X Detergent-free Phusion HF Reaction Buffer	4 x 1.5 ml
50 mM MgCl ₂	2 x 1.5 ml
100% DMSO	0.5 ml

Reaction Conditions: Use at a 1X final concentration in the Phusion reaction. The 1X buffer contains 1.5 mM MgCl₂. Additional MgCl₂ (50 mM) is provided for those amplicons that require higher Mg⁺⁺ concentrations.

Quality Assurance

The reaction buffer and supplements are free of detectable nucleases.

Phusion® DNA Polymerase was developed by Finnzymes Oy, now a part of Thermo Fisher Scientific. This product is manufactured by New England Biolabs, Inc. under agreement with, and under the performance specifications of Thermo Fisher Scientific.

PHUSION® is a registered trademark and property of Thermo Fisher Scientific.

THERMO SCIENTIFIC® is a registered trademark of Thermo Fisher Scientific.

CERTIFICATE OF ANALYSIS