

XRN-1



1-800-632-7799
info@neb.com
www.neb.com



M0338S 003160118011

M0338S



20 units 1,000 U/ml Lot: 0031601

RECOMBINANT Store at -20°C Exp: 1/18

Description: XRN-1 is highly processive 5' → 3' exoribonuclease, requiring 5' monophosphate. It also acts on 5' monophosphate ssDNA with greatly reduced efficiency. XRN-1 is not efficient in removing RNAs with recessed 5'P, like tRNA.

Source: Purified from *E. coli* carrying a plasmid overexpressing the yeast XRN-1 gene (1).

Supplied in: 20 mM Tris-HCl (pH 7.5), 500 mM NaCl, 0.1 mM EDTA, 2 mM dithiothreitol, 0.1% Triton X-100 and 50% glycerol.

Applications:

- Removal of RNA containing 5' monophosphate from an RNA mixture.

Reagents Supplied with Enzyme:
10X NEBuffer 3.

Reaction Conditions: 1X NEBuffer 3.
Incubate at 37°C.

1X NEBuffer 3:

100 mM NaCl
50 mM Tris-HCl
10 mM MgCl₂
1 mM dithiothreitol
pH 7.9 @ 25°C

Unit Definition: One unit is defined as the amount of enzyme that digests 1 µg of phosphorylated yeast RNA in 60 minutes at 37°C.

Quality Control Assays

Protein Purity Assay (SDS-PAGE): XRN-1 is ≥ 95% pure as determined by SDS-PAGE analysis using Coomassie Blue detection.

Endonuclease Activity: dsDNA: Incubation of 2 units of XRN-1 with 300 ng of supercoiled plasmid in 10 µl for 4 hours at 37°C produced less than 10% nicked or linear molecules as determined by agarose gel electrophoresis.

Endonuclease Activity: ssDNA: Incubation of 2 units of XRN-1 with 250 ng of M13 SSDNA in 10 µl for 4 hours at 37°C resulted in no detectable degradation as determined by agarose gel electrophoresis.

Heat Inactivation: 70°C for 10 minutes.

Notes on use: We recommend adding Murine RNase Inhibitor (NEB #M0314) in your reaction at 1 unit/µl to prevent non-specific RNA degradation.

Reference:

1. Stevens, A. (1980) *J. Biol. Chem.* 255, 3080-3085

Companion Products:

Murine RNase Inhibitor
#M0314S 3,000 units
#M0314L 15,000 units

RNase Inhibitor
#M0307S 2,000 units
#M0307L 10,000 units



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

CERTIFICATE OF ANALYSIS

XRN-1



1-800-632-7799
info@neb.com
www.neb.com



M0338S 003160118011

M0338S



20 units 1,000 U/ml Lot: 0031601

RECOMBINANT Store at -20°C Exp: 1/18

Description: XRN-1 is highly processive 5' → 3' exoribonuclease, requiring 5' monophosphate. It also acts on 5' monophosphate ssDNA with greatly reduced efficiency. XRN-1 is not efficient in removing RNAs with recessed 5'P, like tRNA.

Source: Purified from *E. coli* carrying a plasmid overexpressing the yeast XRN-1 gene (1).

Supplied in: 20 mM Tris-HCl (pH 7.5), 500 mM NaCl, 0.1 mM EDTA, 2 mM dithiothreitol, 0.1% Triton X-100 and 50% glycerol.

Applications:

- Removal of RNA containing 5' monophosphate from an RNA mixture.

Reagents Supplied with Enzyme:
10X NEBuffer 3.

Reaction Conditions: 1X NEBuffer 3.
Incubate at 37°C.

1X NEBuffer 3:

100 mM NaCl
50 mM Tris-HCl
10 mM MgCl₂
1 mM dithiothreitol
pH 7.9 @ 25°C

Unit Definition: One unit is defined as the amount of enzyme that digests 1 µg of phosphorylated yeast RNA in 60 minutes at 37°C.

Quality Control Assays

Protein Purity Assay (SDS-PAGE): XRN-1 is ≥ 95% pure as determined by SDS-PAGE analysis using Coomassie Blue detection.

Endonuclease Activity: dsDNA: Incubation of 2 units of XRN-1 with 300 ng of supercoiled plasmid in 10 µl for 4 hours at 37°C produced less than 10% nicked or linear molecules as determined by agarose gel electrophoresis.

Endonuclease Activity: ssDNA: Incubation of 2 units of XRN-1 with 250 ng of M13 SSDNA in 10 µl for 4 hours at 37°C resulted in no detectable degradation as determined by agarose gel electrophoresis.

Heat Inactivation: 70°C for 10 minutes.

Notes on use: We recommend adding Murine RNase Inhibitor (NEB #M0314) in your reaction at 1 unit/µl to prevent non-specific RNA degradation.

Reference:

1. Stevens, A. (1980) *J. Biol. Chem.* 255, 3080-3085

Companion Products:

Murine RNase Inhibitor
#M0314S 3,000 units
#M0314L 15,000 units

RNase Inhibitor
#M0307S 2,000 units
#M0307L 10,000 units



NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

CERTIFICATE OF ANALYSIS