

pNEB193 Vector



1-800-632-7799
info@neb.com
www.neb.com



N3051S 013120514051

N3051S

50 µg Lot: 0131205 Exp: 5/14

500 µg/ml Store at -20°C

Description: pNEB193 is a pUC19 derivative that carries unique *Ascl*, *PmeI*, and *PacI* sites. The only differences between pUC19 and pNEB193 are in the polylinker region. A unique *Ascl* site (GGCGGCC) is located between the *Bam*HI site and the *Sma*I site. A unique *PacI* site (TTAATTAA) is located between the *Bam*HI site and the *Xba*I site, and a unique *PmeI* site (GTTTAAAC), is located between the *Pst*I site and the *Sal*I site.

The polylinker does not interrupt the *lacZ* reading frame. All or part of the polylinker can be easily excised with the appropriate restriction endonucleases and placed in your own DNA vector.

Supplied in: 10 mM Tris-HCl (pH 8.0), 1 mM EDTA.

Preparation: pNEB193 is isolated from *E. coli* TB1 (*hsd M*⁺) by a standard plasmid purification procedure.

Notes: The *Ascl* site is also a unique *Bss*HI site in pNEB193.

E. coli strain TB1 is recommended for transformations.

pNEB193 Polylinker:

EcoRI KpnI Ascl PacI SalI PstI HindIII
GAATTCGAGCTCGGTACCCGGGGGCGCGCCGGATCCTTAATTAAGTCTAGAGTCGACTGTTTAAACCTGCAGGCATGCAAGCTT
SacI SmaI BamHI XbaI PmeI SphI

CERTIFICATE OF ANALYSIS

pNEB193 Vector



1-800-632-7799
info@neb.com
www.neb.com



N3051S 013120514051

N3051S

50 µg Lot: 0131205 Exp: 5/14

500 µg/ml Store at -20°C

Description: pNEB193 is a pUC19 derivative that carries unique *Ascl*, *PmeI*, and *PacI* sites. The only differences between pUC19 and pNEB193 are in the polylinker region. A unique *Ascl* site (GGCGGCC) is located between the *Bam*HI site and the *Sma*I site. A unique *PacI* site (TTAATTAA) is located between the *Bam*HI site and the *Xba*I site, and a unique *PmeI* site (GTTTAAAC), is located between the *Pst*I site and the *Sal*I site.

The polylinker does not interrupt the *lacZ* reading frame. All or part of the polylinker can be easily excised with the appropriate restriction endonucleases and placed in your own DNA vector.

Supplied in: 10 mM Tris-HCl (pH 8.0), 1 mM EDTA.

Preparation: pNEB193 is isolated from *E. coli* TB1 (*hsd M*⁺) by a standard plasmid purification procedure.

Notes: The *Ascl* site is also a unique *Bss*HI site in pNEB193.

E. coli strain TB1 is recommended for transformations.

pNEB193 Polylinker:

EcoRI KpnI Ascl PacI SalI PstI HindIII
GAATTCGAGCTCGGTACCCGGGGGCGCGCCGGATCCTTAATTAAGTCTAGAGTCGACTGTTTAAACCTGCAGGCATGCAAGCTT
SacI SmaI BamHI XbaI PmeI SphI

CERTIFICATE OF ANALYSIS

References:

1. Guan, C. unpublished observations.

References:

1. Guan, C. unpublished observations.