

# $\alpha$ 1-2,3 Mannosidase



1-800-632-7799  
info@neb.com  
www.neb.com



P0729S 017151216121

## P0729S



**640 units**    **32,000 U/ml**    **Lot: 0171512**  
**RECOMBINANT**    **Store at 4°C**    **Exp: 12/16**


**Description:**  $\alpha$ 1-2,3 Mannosidase is a highly specific exoglycosidase that catalyzes the hydrolysis of  $\alpha$ 1-2 and  $\alpha$ 1-3 linked D-mannopyranosyl residues from oligosaccharides (1).

**Specificity:**  
↓  
**Man  $\alpha$  1-2 R**  
 **$\alpha$  1-3 R**


**Note:** p-nitrophenyl- $\alpha$ -D-mannopyranoside is NOT a substrate for this enzyme.

**Detailed Specificity:** Specificity can vary depending on incubation time and concentration of substrate (Figure 1).


### A. 0.1 nm/ $\mu$ l substrate, 1 hour incubation




### B. 0.1 nm/ $\mu$ l substrate, 1 hour incubation




### C. 0.1 nm/ $\mu$ l substrate, 18 hour incubation



### D. 0.1 nm/ $\mu$ l substrate, 18 hour incubation



### E. 0.045 nm/ $\mu$ l substrate, 18 hour incubation



**Figure 1:** Detailed specificity of  $\alpha$ 1,2-3 Mannosidase. All reactions contained 32 units of  $\alpha$ 1,2-3 Mannosidase, 1X GlycoBuffer 1 and 1X BSA in a total reaction volume of 10  $\mu$ l. Reactions were incubated at 37°C. The substrate depicted in (E) will not cut to completion.

**Source:** Cloned from *Xanthomonas manihotis* and expressed in *E. coli* (2).

Supplied in: 50 mM NaCl, 20 mM Tris-HCl (pH 7.5 @ 25°C) and 1 mM Na<sub>2</sub>EDTA.

**Reagents Supplied with Enzyme:**  
10X GlycoBuffer 1  
100X BSA

### Reaction Conditions:

1X GlycoBuffer 1:  
50 mM Sodium Acetate (pH 5.5 @ 25°C)  
and 5 mM CaCl<sub>2</sub>. Supplement with 100  $\mu$ g/ml BSA. Incubate at 37°C.

Optimal incubation times and enzyme concentrations must be determined empirically for a particular substrate.

**Unit Definition:** One unit is defined as the amount of enzyme required to cleave > 95% of the non-reducing terminal  $\alpha$ -D-Mannose from 1 nmol Man $\alpha$ 1-3Man $\beta$ 1-4GlcNAc-7-amino-4-methylcoumarin (AMC), in 1 hour at 37°C in a total reaction volume of 10  $\mu$ l.

(see other side)  
CERTIFICATE OF ANALYSIS

# $\alpha$ 1-2,3 Mannosidase



1-800-632-7799  
info@neb.com  
www.neb.com



P0729S 017151216121

## P0729S



**640 units**    **32,000 U/ml**    **Lot: 0171512**  
**RECOMBINANT**    **Store at 4°C**    **Exp: 12/16**


**Description:**  $\alpha$ 1-2,3 Mannosidase is a highly specific exoglycosidase that catalyzes the hydrolysis of  $\alpha$ 1-2 and  $\alpha$ 1-3 linked D-mannopyranosyl residues from oligosaccharides (1).

**Specificity:**  
↓  
**Man  $\alpha$  1-2 R**  
 **$\alpha$  1-3 R**


**Note:** p-nitrophenyl- $\alpha$ -D-mannopyranoside is NOT a substrate for this enzyme.

**Detailed Specificity:** Specificity can vary depending on incubation time and concentration of substrate (Figure 1).


### A. 0.1 nm/ $\mu$ l substrate, 1 hour incubation




### B. 0.1 nm/ $\mu$ l substrate, 1 hour incubation




### C. 0.1 nm/ $\mu$ l substrate, 18 hour incubation



### D. 0.1 nm/ $\mu$ l substrate, 18 hour incubation



### E. 0.045 nm/ $\mu$ l substrate, 18 hour incubation



**Figure 1:** Detailed specificity of  $\alpha$ 1,2-3 Mannosidase. All reactions contained 32 units of  $\alpha$ 1,2-3 Mannosidase, 1X GlycoBuffer 1 and 1X BSA in a total reaction volume of 10  $\mu$ l. Reactions were incubated at 37°C. The substrate depicted in (E) will not cut to completion.

**Source:** Cloned from *Xanthomonas manihotis* and expressed in *E. coli* (2).

Supplied in: 50 mM NaCl, 20 mM Tris-HCl (pH 7.5 @ 25°C) and 1 mM Na<sub>2</sub>EDTA.

**Reagents Supplied with Enzyme:**  
10X GlycoBuffer 1  
100X BSA

### Reaction Conditions:

1X GlycoBuffer 1:  
50 mM Sodium Acetate (pH 5.5 @ 25°C)  
and 5 mM CaCl<sub>2</sub>. Supplement with 100  $\mu$ g/ml BSA. Incubate at 37°C.

Optimal incubation times and enzyme concentrations must be determined empirically for a particular substrate.

**Unit Definition:** One unit is defined as the amount of enzyme required to cleave > 95% of the non-reducing terminal  $\alpha$ -D-Mannose from 1 nmol Man $\alpha$ 1-3Man $\beta$ 1-4GlcNAc-7-amino-4-methylcoumarin (AMC), in 1 hour at 37°C in a total reaction volume of 10  $\mu$ l.

(see other side)  
CERTIFICATE OF ANALYSIS

**Specific Activity:** ~ 80,000 units/mg

**Molecular Weight:** 90,000 daltons.

**Unit Definition Assay:** Two fold serial dilutions of  $\alpha$ 1-2,3 Mannosidase are incubated with 1 nmol AMC-labeled substrate in 1X GlycoBuffer 1, supplemented with 100  $\mu$ g/ml BSA, in a 10  $\mu$ l reaction. The reaction mix is incubated for 1 hour at 37°C. Separation of reaction products are visualized via thin layer chromatography (1).

**Quality Assurance:** No contaminating exoglycosidase or proteolytic activity could be detected (ND).

### Quality Controls

#### Glycosidase Assays:

32 units of  $\alpha$ 1-2,3 Mannosidase were incubated with 0.1 mM of fluorescently-labeled oligosaccharides and glycopeptides, in a 10  $\mu$ l reaction for 20 hours at 37°C. The reaction products were analyzed by TLC for digestion of substrate.

Page 2 (P0727)

**Specific Activity:** ~ 80,000 units/mg

**Molecular Weight:** 90,000 daltons.

**Unit Definition Assay:** Two fold serial dilutions of  $\alpha$ 1-2,3 Mannosidase are incubated with 1 nmol AMC-labeled substrate in 1X GlycoBuffer 1, supplemented with 100  $\mu$ g/ml BSA, in a 10  $\mu$ l reaction. The reaction mix is incubated for 1 hour at 37°C. Separation of reaction products are visualized via thin layer chromatography (1).

**Quality Assurance:** No contaminating exoglycosidase or proteolytic activity could be detected (ND).

### Quality Controls

#### Glycosidase Assays:

32 units of  $\alpha$ 1-2,3 Mannosidase were incubated with 0.1 mM of fluorescently-labeled oligosaccharides and glycopeptides, in a 10  $\mu$ l reaction for 20 hours at 37°C. The reaction products were analyzed by TLC for digestion of substrate.

Page 2 (P0727)

No other glycosidase activities were detected (ND) with the following substrates:

**$\beta$ -N-Acetylglucosaminidase:**  
GlcNAc $\beta$ 1-4GlcNAc $\beta$ 1-4GlcNAc-AMC ND

**$\alpha$ -N-Acetylgalactosaminidase:**  
GalNAc $\alpha$ 1-3(Fuc $\alpha$ 1-2)Gal $\beta$ 1-4Glc-AMC ND

**$\alpha$ -Fucosidase:**  
Fuc $\alpha$ 1-2Gal $\beta$ 1-4Glc-AMC ND  
Gal $\beta$ 1-4 (Fuc $\alpha$ 1-3)GlcNAc $\beta$ 1-3Gal $\beta$ 1-4Glc-AMC ND

**$\beta$ -Galactosidase:**  
Gal $\beta$ 1-3GlcNAc $\beta$ 1-4Gal $\beta$ 1-4Glc-AMC ND  
Gal $\beta$ 1-4GlcNAc $\beta$ 1-2Man $\alpha$ 1-6Man $\beta$ 1-4GlcNAc-AMC ND

**$\alpha$ -Galactosidase:**  
Gal $\alpha$ 1-3Gal $\beta$ 1-4Gal $\alpha$ 1-3Gal-AMC ND

**$\alpha$ -Mannosidase:**  
Man $\alpha$ 1-6Man $\alpha$ 1-6(Man $\alpha$ 1-3)Man-AMC ND

**$\alpha$ -Neuraminidase:**

No other glycosidase activities were detected (ND) with the following substrates:

**$\beta$ -N-Acetylglucosaminidase:**  
GlcNAc $\beta$ 1-4GlcNAc $\beta$ 1-4GlcNAc-AMC ND

**$\alpha$ -N-Acetylgalactosaminidase:**  
GalNAc $\alpha$ 1-3(Fuc $\alpha$ 1-2)Gal $\beta$ 1-4Glc-AMC ND

**$\alpha$ -Fucosidase:**  
Fuc $\alpha$ 1-2Gal $\beta$ 1-4Glc-AMC ND  
Gal $\beta$ 1-4 (Fuc $\alpha$ 1-3)GlcNAc $\beta$ 1-3Gal $\beta$ 1-4Glc-AMC ND

**$\beta$ -Galactosidase:**  
Gal $\beta$ 1-3GlcNAc $\beta$ 1-4Gal $\beta$ 1-4Glc-AMC ND  
Gal $\beta$ 1-4GlcNAc $\beta$ 1-2Man $\alpha$ 1-6Man $\beta$ 1-4GlcNAc-AMC ND

**$\alpha$ -Galactosidase:**  
Gal $\alpha$ 1-3Gal $\beta$ 1-4Gal $\alpha$ 1-3Gal-AMC ND

**$\alpha$ -Mannosidase:**  
Man $\alpha$ 1-6Man $\alpha$ 1-6(Man $\alpha$ 1-3)Man-AMC ND

**$\alpha$ -Neuraminidase:**

Neu5Ac $\alpha$ 2-3Gal $\beta$ 1-3GlcNAc $\beta$ 1-3Gal $\beta$ 1-4Glc-AMC ND

**$\beta$ -Glucosidase:**  
Glc $\beta$ 1-4Glc $\beta$ 1-4Glc-AMC ND

**$\alpha$ -Glucosidase:**  
Glc $\alpha$ 1-6Glc $\alpha$ 1-4Glc-AMC ND

**$\beta$ -Xylosidase:**  
Xyl $\beta$ 1-4Xyl $\beta$ 1-4Xyl $\beta$ 1-4Xyl-AMC ND

**$\beta$ -Mannosidase:**  
Man $\beta$ 1-4Man $\beta$ 1-4Man-AMC ND

**Endo F<sub>1</sub>, F<sub>2</sub>, H:**  
Dansylated invertase high mannose. ND

**Endo F<sub>2</sub>, F<sub>3</sub>:**  
Dansylated fibrinogen biantennary. ND

**PNase F:**  
Fluoresceinated fetuin triantennary. ND

Neu5Ac $\alpha$ 2-3Gal $\beta$ 1-3GlcNAc $\beta$ 1-3Gal $\beta$ 1-4Glc-AMC ND

**$\beta$ -Glucosidase:**  
Glc $\beta$ 1-4Glc $\beta$ 1-4Glc-AMC ND

**$\alpha$ -Glucosidase:**  
Glc $\alpha$ 1-6Glc $\alpha$ 1-4Glc-AMC ND

**$\beta$ -Xylosidase:**  
Xyl $\beta$ 1-4Xyl $\beta$ 1-4Xyl $\beta$ 1-4Xyl-AMC ND

**$\beta$ -Mannosidase:**  
Man $\beta$ 1-4Man $\beta$ 1-4Man-AMC ND

**Endo F<sub>1</sub>, F<sub>2</sub>, H:**  
Dansylated invertase high mannose. ND

**Endo F<sub>2</sub>, F<sub>3</sub>:**  
Dansylated fibrinogen biantennary. ND

**PNase F:**  
Fluoresceinated fetuin triantennary. ND

**Protease Assay:** After incubation of 220 units of  $\alpha$ 1-2,3 Mannosidase with 0.2 nmol of a standard mixture of proteins for 20 hours at 37°C, no proteolytic activity could be detected by SDS-PAGE.

### References:

1. Wong-Madden, S.T. and Landry, D. (1995) *Glycobiology* 5, 19–28.
2. Guthrie, E.P., Taron, C.H., New England Biolabs, Inc. unpublished results.

U.S. Patent No. 7,094,563

NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.



**Protease Assay:** After incubation of 220 units of  $\alpha$ 1-2,3 Mannosidase with 0.2 nmol of a standard mixture of proteins for 20 hours at 37°C, no proteolytic activity could be detected by SDS-PAGE.

### References:

1. Wong-Madden, S.T. and Landry, D. (1995) *Glycobiology* 5, 19–28.
2. Guthrie, E.P., Taron, C.H., New England Biolabs, Inc. unpublished results.

U.S. Patent No. 7,094,563

NEW ENGLAND BIOLABS® is a registered trademark of New England Biolabs, Inc.

This product is intended for research purposes only. This product is not intended to be used for therapeutic or diagnostic purposes in humans or animals.

