



Discontinuation Notice/New Product Introduction

Effective June 15th, the HiScribe™ T7 *In Vitro* Transcription Kit (NEB #E2030) will be discontinued.

As an improved alternative, NEB now offers the **T7 Quick High Yield RNA Synthesis Kit**. This new kit provides the same level of performance and yield as the HiScribe T7 *In Vitro* Transcription Kit, but with several protocol enhancements.

Easier Workflow

- Streamlined reaction setup reducing hands-on time
- Fewer pipetting steps reducing chances for technical error

Expanded protocols

- Supporting co-transcriptional capping
- Supporting incorporation of modified nucleotides

Complete Solution

- Including template removal and RNA cleanup reagents

Questions or Concerns?

If you have been using the HiScribe T7 *In Vitro* Transcription Kit, you can expect the same level of performance with an easier bench experience. For questions or concerns, please contact: rna@neb.com.