
New England Biolabs Product Specification

Product Name: Anti-MBP Monoclonal Antibody
Catalog #: E8032S/L
Recommended Dilution: 1:10,000
Shelf Life: 24 months
Storage Temp: -20°C
Storage Conditions: 150 mM NaCl, 10 mM HEPES, 50 % Glycerol, (pH 7.5 @ 25°C)
Specification Version: PS-E8032S/L v1.0
Effective Date: 14 Feb 2018

Assay Name/Specification (minimum release criteria)

Functional Testing (Western Blot, Antibody) - A 1:10,000 dilution of the antibody is tested for its ability to detect 10 ng of MBP by Western Blot.



Date 14 Feb 2018

Derek Robinson
Director of Quality Control

