

---

## New England Biolabs Product Specification

*Product Name:* SNAP-Cell<sup>®</sup> Starter Kit  
*Catalog #:* E9100S  
*Kit Components:* pSNAPf Vector (N9183)  
pSNAPf-Cox8A Control Plasmid (N9185)  
SNAP-Cell<sup>®</sup> 505-Star (S9103)  
SNAP-Cell<sup>®</sup> TMR-Star (S9105)  
SNAP-Cell<sup>®</sup> Block (S9106)  
*Shelf Life:* 24 months  
*Storage Temp:* -20°C  
*Specification Version:* PS-E9100S v1.0  
*Effective Date:* 22 Dec 2017

Assay Name/Specification (minimum release criteria)
---

* <b>Individual Product Component Note</b> - Standard Quality Control Tests are performed for each component included in SNAP-Cell <sup>®</sup> Starter Kit and meet the designated specifications.
---



Date 22 Dec 2017

---

Derek Robinson  
Director of Quality Control

