

## New England Biolabs Product Specification

*Product Name:* Low Range PFG Marker  
*Catalog #:* N0350S  
*Concentration:* 25 µg/ml  
*Unit Definition:* N/A  
*Shelf Life:* 24 months  
*Storage Temp:* -20°C  
*Storage Conditions:* 10 mM Tris-HCl (pH 8.0), 1 mM EDTA, 50 % Glycerol, 1 % LMP Agarose  
*Specification Version:* PS-N0350S v1.0  
*Effective Date:* 16 May 2014

Assay Name/Specification (minimum release criteria)
---

<b>Electrophoretic Pattern (PFGE Marker)</b> - The banding pattern of Low Range PFG Marker on a 1% CHEF PFG gel shows discrete, clearly identifiable bands at each band of the marker, as determined by gel electrophoresis using Ethidium Bromide.
---



Derek Robinson  
Director of Quality Control

Date 16 May 2014

