

New England Biolabs Product Specification

Product Name: MidRange I PFG Marker
Catalog #: N3551S
Concentration: 50 µg/ml
Unit Definition: N/A
Shelf Life: 24 months
Storage Temp: -20°C
Storage Conditions: 10 mM Tris-HCl (pH 8.0), 1 mM EDTA, 50 % Glycerol, 1 % LMP Agarose
Specification Version: PS-N3551S v1.0
Effective Date: 15 May 2014

Assay Name/Specification (minimum release criteria)

Electrophoretic Pattern (PFGE Marker) - The banding pattern of MidRange I PFG Marker on a 1% CHEF PFG gel shows discrete, clearly identifiable bands at each band of the marker, as determined by gel electrophoresis using Ethidium Bromide.
--



Derek Robinson
Director of Quality Control

Date 15 May 2014

