
New England Biolabs Product Specification

Product Name: ssRNA Ladder
Catalog #: N0362S
Shelf Life: 12 months
Storage Temp: -20°C
Composition: 20 mM Sodium Citrate , 1 mM EDTA, (pH 6.0 @ 25°C)
Specification Version: PS-N0362S v2.0
Effective Date: 20 Mar 2019

Assay Name/Specification (minimum release criteria)

Electrophoretic Pattern (Marker) - The banding pattern of ssRNA Ladder on a 1% agarose-TBE gel shows discrete, clearly identifiable bands at each fragment of the marker, when stained with Ethidium Bromide at a concentration of 0.5 µg/ml.
--



Date 20 Mar 2019

Derek Robinson
Director of Quality Control

