

New England Biolabs Product Specification

Product Name: Oligo (dT)₂₅ Cellulose Beads
Catalog #: S1408S
Concentration: 80 mg/ml
Shelf Life: 24 months
Storage Temp: 4°C
Storage Conditions: 0.02% NaN₃
Specification Version: PS-S1408S v1.0
Effective Date: 26 Nov 2018

Assay Name/Specification (minimum release criteria)

Binding Capacity (Poly rA) - Oligo (dT)₂₅ Cellulose Beads (0.15 ml) were equilibrated and incubated with 0.5 ml of PolyA substrate solution for 5 minutes at 25°C, then washed and the PolyA eluted. Binding capacity was determined to be >400 OD per gram of cellulose.



Date 26 Nov 2018

Derek Robinson
Director of Quality Control

