

New England Biolabs Product Specification

Product Name: SNAP-Surface[®] Alexa Fluor[®] 546
Catalog #: S9132S
Shelf Life: 60 months
Storage Temp: -20°C
Specification Version: PS-S9132S v3.0
Effective Date: 11 Apr 2019

Assay Name/Specification (minimum release criteria)

Identity (Mass Spectrometry) - The observed molecular mass of SNAP-Surface[®] Alexa Fluor[®] 546 is 1211.2 Da +/- 1 Da as determined by mass spectrometry analysis.

In Vitro Protein Labeling - A 50 µl reaction in 1X PBS and 1 mM DTT containing 5 µM of SNAP-tag[®] Purified Protein and a minimum of 10 µM of SNAP-Surface[®] Alexa Fluor[®] 546 is incubated for 1 hour at 37°C results in the expected labeled product that is visualized on SDS-PAGE by fluorescent detection.

Physical Purity (HPLC) - SNAP-Surface[®] Alexa Fluor[®] 546 is ≥ 90% pure as determined by HPLC analysis.



Date 11 Apr 2019

Derek Robinson
Director of Quality Control

