
New England Biolabs Product Specification

Product Name: Monarch[®] Bead Retainers
Catalog #: T3004L
Storage Temp: 25°C
Specification Version: PS-T3004L v1.0
Effective Date: 16 Dec 2020

Assay Name/Specification (minimum release criteria)

Integrity (Physical Inspection) - A representative selection of each production batch of Monarch [®] Bead Retainers is visually inspected for signs of damage or manufacturing defect prior to batch release.

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 16 Dec 2020

Derek Robinson
Quality Approver

