
New England Biolabs Product Specification

Product Name: Monarch[®] gDNA Elution Buffer
Catalog #: T3016L
Shelf Life: 18 months
Storage Temp: 25°C
Composition: Proprietary
Specification Version: PS-T3016L v1.0
Effective Date: 18 Mar 2019

Assay Name/Specification (minimum release criteria)
--

Functional Testing (DNA Recovery and Purity) - Twenty-four Monarch [®] gDNA Purification Columns are tested with 5 µg of HindIII digested Lambda resulting in ≥75% recovery in ≥80% of the samples. OD 260/280 and 260/230 are ≥1.75 in ≥80% of the samples.
--



Date 18 Mar 2019

Derek Robinson
Director of Quality Control

