

NEB #	DNA Modifying Enzymes and Cloning Technologies	Size
M0202S	T4 DNA Ligase	20,000 units
M0201S	T4 Polynucleotide Kinase	500 units
NEB #	Nucleic Acid Purification	Size
T1120S	Monarch® Spin DNA Gel Extraction Kit	50 preps
T1130S	Monarch® Spin PCR & DNA Cleanup Kit (5 µg)	50 preps
T1110S	Monarch® Spin Plasmid Miniprep Kit	50 preps
T2110S	Monarch® Spin RNA Isolation Kit (Mini)	50 preps
NEB #	PCR, Polymerases & Amplification Technologies	Size
N0447S	Deoxynucleotide (dNTP) Solution Mix	8 µmol
M0530S	Phusion™ High-Fidelity DNA Polymerase	100 units
M0491S	Q5® High-Fidelity DNA Polymerase	100 units
M0267S	Taq DNA Polymerase with ThermoPol® Buffer	400 units
NEB #	Protein Tools	Size
P7719S	Color Prestained Protein Standard, Broad Range (10-250 kDa)	0.45 ml
NEB #	Restriction Endonucleases	Size
R3142S	KpnI-HF®	4,000 units
R3198S	MluI-HF®	1,000 units
R3138S	SalI-HF®	2,000 units

NEB Representative: Alexandre Rochette, Ph.D.

arochette@neb.com

(514) 233-4557

Freezer Support: 1-800-632-5227

Email: info.ca@neb.com