



University of Saskatchewan - Health Science Supply Centre (HSSC)
HSC Supply Centre, Julie Zacharias
hsc.supplycentre@usask.ca
306-966-6166

NEB #	DNA Modifying Enzymes and Cloning Technologies	Size
M0202S	T4 DNA Ligase	20,000 units
NEB #	Nucleic Acid Purification	Size
T1120S	Monarch [®] Spin DNA Gel Extraction Kit	50 preps
T1110S	Monarch [®] Spin Plasmid Miniprep Kit	50 preps
NEB #	PCR, Polymerases & Amplification Technologies	Size
M0530S	Phusion [™] High-Fidelity DNA Polymerase	100 units
M0491S	Q5 [®] High-Fidelity DNA Polymerase	100 units
M0273S	Taq DNA Polymerase with Standard Taq Buffer	400 units
M0267S	Taq DNA Polymerase with ThermoPol [®] Buffer	400 units
NEB #	Restriction Endonucleases	Size
R0176S	DpnI	1,000 units
R3189S	NotI-HF [®]	500 units

NEB Representative: Julie Wong, Ph.D.
jwong@neb.com
(604) 767-5078

Freezer Support: 1-800-632-5227
Email: info.ca@neb.com