
New England Biolabs Product Specification

Product Name: NEBuffer™ Set (1.1, 2.1, 3.1 and CutSmart®)

Catalog #: B7200S

Kit Components: NEBuffer™ 1.1 (B7201)
NEBuffer™ 2.1 (B7202)
NEBuffer™ 3.1 (B7203)
CutSmart® Buffer (B7204)

Shelf Life: 36 months

Storage Temp: -20°C

Specification Version: PS-B7200S v1.0

Effective Date: 16 Jul 2019

Assay Name/Specification (minimum release criteria)
--

* Individual Product Component Note - Standard Quality Control Tests are performed for each component included in NEBuffer™ Set (1.1, 2.1, 3.1 and CutSmart®) and meet the designated specifications.
--



Date 16 Jul 2019

Derek Robinson
Director of Quality Control

