240 County Road Ipswich, MA 01938-2723

Tel 978-927-5054 Fax 978-921-1350 www.neb.com info@neb.com

New England Biolabs Product Specification

Product Name: NEBufferTM 1.1

Catalog #: B7201S

Concentration: 10X Concentrate

Shelf Life: 36 months
Storage Temp: -20°C

Composition (1X): 10 mM Bis-Tris Propane, 10 mM MgCl₂, 100 µg/ml BSA, (pH 7.0 @ 25°C)

Specification Version: PS-B7201S v1.0
Effective Date: 14 Sep 2017

Assay Name/Specification (minimum release criteria)

Conductivity (buffers/solutions) - The conductivity of 10X NEBuffer 1.1 is between 23 and 26 mS at 25°C.

Endonuclease Activity (Nicking, Buffer) - A 50 μ l reaction in 1X NEBuffer 1.1 containing 1 μ g of supercoiled PhiX174 DNA incubated for 4 hours at 37°C results in <10% conversion to the nicked form as determined by agarose gel electrophoresis.

Functional Testing (Restriction Digest, Buffer) - A 50 µl reaction in 1X NEBuffer 1.1 containing 1 µg of Lambda-HindIII DNA and 1 unit of SacI incubated for 1 hour at 37°C results in complete digestion of the substrate DNA as determined by agarose gel electrophoresis.

Functional Testing (Restriction Digest, Buffer) - A 50 μ l reaction in 1X NEBuffer 1.1 containing 1 μ g of pXba DNA and 1 unit of KpnI incubated for 1 hour at 37°C results in complete digestion of the substrate DNA as determined by agarose gel electrophoresis.

Non-Specific DNase Activity (16 hour, Buffer) - A 50 μ l reaction in 1X NEBuffer 1.1 containing 1 μ g of PhiX174-HaeIII DNA incubated for 16 hours at 37°C results in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis.

pH (buffers/solutions) - The pH of 10X NEBuffer 1.1 is between pH 6.9 and 7.1 at 25°C.

RNase Activity (Buffer) - A 10 μ l reaction in 1X NEBuffer 1.1 containing 40 ng of a 300 base single-stranded RNA is incubated at 37°C. After incubation for 16 hours, >90% of the substrate RNA remains intact as determined by fluorescent detection.

Derek Robinson

Director of Quality Control







14 Sep 2017

Date

^{*} The BSA in this product has been granted an EDQM "Certificate of Suitability" from the European Directorate for the Quality of Medicines (# R1-CEP-2003-204-Rev00) and has been granted a USDA Certificate for Export of Bovine Blood Plasma/Serum for Manufacture into Pharmaceutical Products.