

240 County Road Ipswich, MA 01938-2723 Tel 978-927-5054 Fax 978-921-1350 www.neb.com info@neb.com

New England Biolabs Product Specification

| Product Name: | NEBNext® Ultra™ II End Repair/dA-Tailing Module |
|-----------------|-------------------------------------------------------------------------------------------------------|
| Catalog #: | E7546S/L |
| Kit Components: | NEBNext® Ultra™ II End Prep Enzyme Mix (E7646) NEBNext® Ultra™ II End Prep Reaction Buffer (E7647) |
| Shelf Life: | 12 months |

Storage Temp:-20°CSpecification Version:PS-E7546S/L v1.0Effective Date:10 May 2018

Assay Name/Specification (minimum release criteria)

Functional Testing (Library Construction, DNA Modules) - Each set of reagents is functionally validated and compared to a previous lot through construction of DNA libraries. Libraries made from previous and current lots for the recommended input DNA amounts are sequenced together on the same Illumina® flow cell and compared across various metrics including library yield, fraction of reads aligning to the reference, GC bias, and insert size.

* Individual Product Component Note - Standard Quality Control Tests are performed for each component included in NEBNext® Ultra™ II End Repair/dA-Tailing Module and meet the designated specifications.

Date 10 May 2018

Derek Robinson Director of Quality Control



PS-E7546S/L v1.0 Page 1 of 1