
New England Biolabs Product Specification

Product Name: Anti-M13 pIII Monoclonal Antibody
Catalog #: E8033S
Recommended Dilution: 1:1,000
Shelf Life: 24 months
Storage Temp: -20°C
Storage Conditions: 150 mM NaCl, 10 mM HEPES, 50 % Glycerol, (pH 7.5 @ 25°C)
Specification Version: PS-E8033S v1.0
Effective Date: 31 Oct 2018

Assay Name/Specification (minimum release criteria)
Functional Testing (Western Blot, Antibody, Exp) - A 1:1,000 dilution of the antibody is tested for its ability to detect 10 ⁷ pfu M13 by Western Blot.



Date 31 Oct 2018

Derek Robinson
Director of Quality Control

