

240 County Road Ipswich, MA 01938-2723 Tel 978-927-5054 Fax 978-921-1350 www.neb.com info@neb.com

## New England Biolabs Product Specification

Product Name:	Luna® Universal qPCR Master Mix
Catalog #:	M3003S/L/G
Concentration:	2X Concentrate
Shelf Life:	24 months
Storage Temp:	-20°C
Composition (1X):	Proprietary
Specification Version:	PS-M3003S/L/G v1.0
Effective Date:	01 Dec 2016

Assay Name/Specification (minimum release criteria)

Functional Testing (qPCR) - Luna® Universal qPCR Master Mix is functionally tested in qPCR with human cDNA template, resulting in a standard curve with a calculated qPCR efficiency of 90-110%, and a dynamic range of 5 orders of magnitude.

Non-Specific DNase Activity (16 hour, Master Mix) - A 50  $\mu$ l reaction in 1X Luna® Universal qPCR Master Mix containing 1  $\mu$ g of T3 DNA in addition to a reaction containing Lambda-HindIII DNA incubated for 16 hours at 37°C results in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis.

**qPCR DNA Contamination (***E. coli* **Genomic)** - A minimum of 1  $\mu$ l of Luna® Universal qPCR Master Mix is screened for the presence of *E. coli* genomic DNA using SYBR® Green qPCR with primers specific for the *E. coli* 16S rRNA locus. Results are quantified using a standard curve generated from purified *E. coli* genomic DNA. The measured level of *E. coli* genomic DNA contamination is  $\leq$  1 *E. coli* genome.

**RNase Activity Assay (4 Hour Digestion)** - A 10  $\mu$ l reaction in NEBuffer 4 containing 40 ng of a 300 base single-stranded RNA and a minimum of 1  $\mu$ l of Luna® Universal qPCR Master Mix is incubated at 37°C. After incubation for 4 hours, >90% of the substrate RNA remains intact as determined by gel electrophoresis using fluorescent detection.

Date 01 Dec 2016

Derek Robinson Director of Quality Control



PS-M3003S/L/G v1.0 Page 1 of 1