

240 County Road Ipswich, MA 01938-2723 Tel 978-927-5054 Fax 978-921-1350

www.neb.com info@neb.com

New England Biolabs **Product Specification**

Product Name: p-Nitrophenyl Phosphate (PNPP)

Catalog #: P0757S/L Concentration: 500 mM Shelf Life: 24 months Storage Temp: -20°C

Storage Conditions: Supplied in sterile purified water

Specification Version: PS-P0757S/L v1.0 Effective Date: 22 Oct 2018

Assay Name/Specification (minimum release criteria)

Functional Testing (PNPP) - A 50 µl reaction in 1X Phosphatase Buffer containing 50 mM p-Nitrophenyl Phosphate (PNPP) and 8 units of Lambda Protein Phosphatase incubated for 5 minutes at 30°C results in an OD of ≥ 0.40 as determined by UV absorption at 405 nm.

Date 22 Oct 2018

Derek Robinson Director of Quality Control





