

240 County Road Ipswich, MA 01938-2723 Tel 978-927-5054 Fax 978-921-1350

www.neb.com info@neb.com

New England Biolabs Product Specification

Product Name: SNAP-tag® Purified Protein

Catalog #: P9312S
Concentration: 50 µM
Shelf Life: 24 months
Storage Temp: -20°C

Storage Conditions: 1 X PBS; 1 mM DTT

Specification Version: PS-P9312S v1.0

Effective Date: 15 Sep 2017

Assay Name/Specification (minimum release criteria)

In Vitro Protein Labeling - A 50 μ l reaction in 1X PBS and 1 mM DTT containing 10 μ M of a fluorescent SNAP-tag substrate and a minimum of 5 μ M of SNAP-tag® Purified Protein is incubated for 1 hour at 37°C results in the expected labeled product that is visualized on SDS-PAGE by agarose gel electrophoresis.

Ruh Cotum

Date 15 Sep 2017

Derek Robinson

Director of Quality Control





