

240 County Road Ipswich, MA 01938-2723 Tel 978-927-5054 Fax 978-921-1350

Date

02 Nov 2021

www.neb.com info@neb.com

New England Biolabs Product Specification

Product Name: $m^7G(5')ppp(5')G$ RNA Cap Structure Analog

Catalog #: S1404S/LConcentration: $1 \mu mol$ Shelf Life: 36 monthsStorage Temp: $-20^{\circ}C$

Storage Conditions: Supplied as dried sodium salt.

Specification Version: PS-S1404S/L v2.0
Effective Date: 02 Nov 2021

Assay Name/Specification (minimum release criteria)

Functional Testing (Incorporation using RNA Polymerase) - A 20 μ l reaction in RNA Polymerase Reaction Buffer in the presence of 4 mM NTPs +/- m⁷G(5')ppp(5')G RNA Cap Structure Analog containing 2 μ g of template DNA and 50 units of T7 RNA Polymerase incubated for 2 hours at 37°C results in the expected product as determined by polyacrylamide gel electrophoresis.

Molecular Weight Determination (Mass Spectrometry) - The molecular weight of m⁷G(5')ppp(5')G RNA Cap Structure Analog is between 802.4 and 804.4 as determined by mass spectrometry analysis.

Physical Purity (HPLC) - $m^7G(5')ppp(5')G$ RNA Cap Structure Analog is $\geq 95\%$ pure as determined by HPLC analysis.

One or more products referenced in this document may be covered by a 3rd-party trademark. Please visit www.neb.com/trademarks for additional information.

Derek Robinson

Director, Quality Control





