
New England Biolabs Product Specification

Product Name: Monarch[®] Plasmid Miniprep Columns
Catalog #: T1017L
Shelf Life: 18 months
Storage Temp: 25°C
Specification Version: PS-T1017L v1.0
Effective Date: 03 Dec 2017

Assay Name/Specification (minimum release criteria)
--

Functional Testing (DNA Recovery and Purity) - Twenty-four Monarch [®] Plasmid Miniprep Columns are tested with 1.5 µg of 1 Kb Ladder resulting in ≥75% recovery in ≥80% of the samples. OD 260/280 and 260/230 are ≥1.75 in ≥80% of the samples.



Date 03 Dec 2017

Derek Robinson
Director of Quality Control

