
New England Biolabs Product Specification

Product Name: Monarch[®] Plasmid Miniprep Columns
Catalog #: T1017L
Shelf Life: 24 months
Storage Temp: 25°C
Specification Version: PS-T1017L v2.0
Effective Date: 01 Nov 2021

| |
|--|
| Assay Name/Specification (minimum release criteria) |
|--|

| |
|---|
| Functional Testing (DNA Recovery and Purity) - Twenty-four Monarch [®] Plasmid Miniprep Columns are tested with 1.5 µg of 1 Kb Ladder resulting in ≥75% recovery in ≥80% of the samples. OD 260/280 and 260/230 are ≥1.75 in ≥80% of the samples. |
|---|

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 01 Nov 2021

Derek Robinson
Director, Quality Control

