240 County Road Ipswich, MA 01938-2723 Tel 978-927-5054 Fax 978-921-1350 www.neb.com info@neb.com

## New England Biolabs Product Specification

Product Name: Monarch® PCR & DNA Cleanup Kit (5 µg)

*Catalog #: T1030S/L* 

Kit Components: Monarch® DNA Elution Buffer (T1016)

Monarch® DNA Cleanup Binding Buffer (T1031)

Monarch® DNA Wash Buffer (T1032)

Monarch® DNA Cleanup Columns (5 μg) (T1034) Monarch® Collection Tubes (2 ml) (T1018)

Shelf Life: 24 months
Storage Temp: 25°C

Specification Version: PS-T1030S/L v2.0
Effective Date: 21 Mar 2019

## Assay Name/Specification (minimum release criteria)

Functional Testing (DNA Recovery and Purity) - Twenty-four Monarch® DNA Cleanup Columns are tested with 1  $\mu$ g of 1 Kb Ladder resulting in  $\geq$ 75% recovery in  $\geq$ 80% of the samples. OD 260/280 and 260/230 are  $\geq$ 1.75 in  $\geq$ 80% of the samples.

Functional Testing (Fragment Purification and Restriction Digestion) - A linear pLA3 fragment is added to 8 Monarch® DNA Cleanup Columns, eluted and digested with DraIII-HF® resulting in ≥90% conversion to two-fragment pattern observed by gel electrophoresis in ≥75% of the samples.

Derek Robinson

Director of Quality Control





