

New England Biolabs Product Specification

Product Name: Monarch[®] gDNA Blood Lysis Buffer
Catalog #: T3013S/L
Shelf Life: 24 months
Storage Temp: 25°C
Composition: Proprietary
Specification Version: PS-T3013S/L v2.0
Effective Date: 16 May 2022

Assay Name/Specification (minimum release criteria)

Functional Testing (Blood Cell Lysis Buffer, Monarch[®]) - Genomic DNA is purified from 8 samples of whole pig blood using the Protocol for Extraction and Purification of Genomic DNA from Mammalian Whole Blood (non-nucleated), with $\geq 2 \mu\text{g}$ of gDNA being recovered in $\geq 80\%$ of the samples. OD 260/280 and 260/230 are ≥ 1.75 in $\geq 80\%$ of the samples.

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 16 May 2022

Derek Robinson
Director, Quality Control

