

NEB #	DNA Modifying Enzymes and Cloning Technologies	Size	Stockroom #
M0289S	Antarctic Phosphatase	1,000 units	M0289S
M0289L	Antarctic Phosphatase	5,000 units	M0289L
M0226S	CpG Methyltransferase (M.SssI)	100 units	M0226S
M2200S	Quick Ligation™ Kit	30 reactions	M2200S
M0371L	Shrimp Alkaline Phosphatase (rSAP)	2,500 units	M0371L
M0371S	Shrimp Alkaline Phosphatase (rSAP)	500 units	M0371S
M0202S	T4 DNA Ligase	20,000 units	M0202S
M0202T	T4 DNA Ligase	20,000 units	M0202T
M0201S	T4 Polynucleotide Kinase	500 units	M0201S
NEB #	Markers & Ladders	Size	Stockroom #
N3232L	1 kb DNA Ladder	1,000 gel lanes (1 ml)	N3232L
N3231L	100 bp DNA Ladder	500 gel lanes (0.5 ml)	N3231L
NEB #	PCR, Polymerases & Amplification Technologies	Size	Stockroom #
M0210S	DNA Polymerase I, Large (Klenow) Fragment	200 units	M0210S
M0323S	LongAmp® Taq DNA Polymerase	500 units	M0323S
M0480S	OneTaq® DNA Polymerase	200 units	M0480S
M0481S	OneTaq® Hot Start DNA Polymerase	200 units	M0481S
M0530S	Phusion™ High-Fidelity DNA Polymerase	100 units	M0530S
M0491L	Q5® High-Fidelity DNA Polymerase	500 units	M0491L
M0491S	Q5® High-Fidelity DNA Polymerase	100 units	M0491S
M0273L	Taq DNA Polymerase with Standard Taq Buffer	2,000 units	M0273L
M0267L	Taq DNA Polymerase with ThermoPol® Buffer	2,000 units	M0267L
NEB #	Restriction Endonucleases	Size	Stockroom #
R0117S	AatII	500 units	R0117S
R0520S	AflII	2,000 units	R0520S
R0114S	Apal	5,000 units	R0114S
R0558S	Ascl	500 units	R0558S
R0174S	AvrII	100 units	R0174S
R0136S	BamHI	10,000 units	R0136S
R3136S	BamHI-HF®	10,000 units	R3136S
R0160S	BclI	3,000 units	R0160S
R0144S	BglII	2,000 units	R0144S
R0113S	BstXI	1,000 units	R0113S

R0197S	Clal	1,000 units	R0197S
R0175S	Ddel	1,000 units	R0175S
R0176S	Dpnl	1,000 units	R0176S
R0101S	EcoRI	10,000 units	R0101S
R3101S	EcoRI-HF®	10,000 units	R3101S
R0195S	EcoRV	4,000 units	R0195S
R3195S	EcoRV-HF®	4,000 units	R3195S
R0104S	HindIII	10,000 units	R0104S
R0193S	NcoI	1,000 units	R0193S
R3193S	NcoI-HF®	1,000 units	R3193S
R0111S	NdeI	4,000 units	R0111S
R3131S	NheI-HF®	1,000 units	R3131S
R0189L	NotI	2,500 units	R0189L
R0189S	NotI	500 units	R0189S
R3189S	NotI-HF®	500 units	R3189S
R0547S	PacI	250 units	R0547S
R0560S	PmeI	500 units	R0560S
R0532S	PmlI	2,000 units	R0532S
R0140S	PstI	10,000 units	R0140S
R0151S	PvuII	5,000 units	R0151S
R3151S	PvuII-HF®	5,000 units	R3151S
R3156S	SacI-HF®	2,000 units	R3156S
R0157S	SacII	2,000 units	R0157S
R0138S	Sall	2,000 units	R0138S
R3138S	Sall-HF®	2,000 units	R3138S
R0141S	SmaI	2,000 units	R0141S
R0182S	SphI	500 units	R0182S
R3182S	SphI-HF®	500 units	R3182S
R0145S	XbaI	3,000 units	R0145S
R0146S	XhoI	5,000 units	R0146S
R0180S	XmaI	500 units	R0180S
R0194S	XmnI	1,000 units	R0194S

NEB Representative: James DeBoe
jdeboe@neb.com
(978) 778-6011

NEBnow Support: 1-800-632-5227
Email: freezers@neb.com