
New England Biolabs Product Specification

Product Name: *K. lactis GG799 Competent Cells*
Catalog #: *C1001S*
Shelf Life: *12 months*
Storage Temp: *-80°C*
Specification Version: *PS-C1001S v1.0*
Effective Date: *08 Aug 2018*

Assay Name/Specification (minimum release criteria)
Transformation Efficiency (<i>K. lactis</i>) - One microgram of linearized pKLAC1-malE was used to transform one tube of <i>K. lactis</i> GG799 competent cells. A minimum of 1×10^2 colonies were formed after incubation for 3 days at 30°C.



Date 08 Aug 2018

Derek Robinson
Director of Quality Control

