

New England Biolabs Product Specification

Product Name: EpiMark[®] Nucleosome Assembly Kit
Catalog #: E5350S
Kit Components: Histone H2A/H2B Dimer Human, Recombinant (M2508)
Histone H3.1/H4 Tetramer Human, Recombinant (M2509)
Nucleosome Control DNA (N1202)
Shelf Life: 12 months
Storage Temp: -20°C
Specification Version: PS-E5350S v1.0
Effective Date: 11 Dec 2017

Assay Name/Specification (minimum release criteria)

Functional Testing (EpiMark[®] Nucleosome Assembly Kit) - The EpiMark[®] Nucleosome Assembly Kit has been functionally tested using the kit components provided as described in the Dilution Assembly Protocol section of the EpiMark[®] Nucleosome Assembly Kit manual.

* **Individual Product Component Note** - Standard Quality Control Tests are performed for each component included in EpiMark[®] Nucleosome Assembly Kit and meet the designated specifications.



Date 11 Dec 2017

Derek Robinson
Director of Quality Control

